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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of JIANGSU KONSUNG BIO-MEDICAL SCIENCE AND TECHNOLOGY CO.,LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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证明网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

APPROVED BY"/ УТВЕРЖДАЮ
General Director/Генеральный директор
Jiangsu Konsung Bio-Medical Science and Technology
Co., Ltd
Wang Qiang

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User's manual for medical device:

Medical pulseoxymetry AORTIS® device in the options

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- Medical pulse oximeter AORTIS® F010 -1 pc;
- Operation manual -1pc.

Accessories:

1. Battery - from 1 to 100 pcs.

2. Lace - from 1 to 100 pcs.

3. Cover - from 1 to 100 pcs.

4. Bag - from 1 to 100 pcs.

5. Mini USB cable - from 1 to 100 pcs.

6. Cuff - from 1 to 100 pcs.

6. Medical pulse oximeter AORTIS® F011 consisting of:

- Medical pulse oximeter AORTIS® F011 -1 pc;
- Operation manual -1pc.

Accessories:

1. Battery - from 1 to 100 pcs.

2. Lace - from 1 to 100 pcs.

3. Cover - from 1 to 100 pcs.

4. Bag - from 1 to 100 pcs.

5. Mini USB cable - from 1 to 100 pcs.

6. Cuff - from 1 to 100 pcs.

7. Medical pulse oximeter AORTIS® F012 consisting of:

- Medical pulse oximeter AORTIS® F012 -1 pc;
- Operation manual -1pc.

Accessories:

1. Battery - from 1 to 100 pcs.

2. Lace - from 1 to 100 pcs.

3. Cover - from 1 to 100 pcs.

4. Bag - from 1 to 100 pcs.

5. Mini USB cable - from 1 to 100 pcs.

6. Cuff - from 1 to 100 pcs.

8. Medical pulse oximeter AORTIS® F013 consisting of:

- Medical pulse oximeter AORTIS® F013 -1 pc;
- Operation manual -1pc.

Accessories:

1. Battery - from 1 to 100 pcs.

2. Lace - from 1 to 100 pcs.

3. Cover - from 1 to 100 pcs.

4. Bag - from 1 to 100 pcs.

5. Mini USB cable - from 1 to 100 pcs.

6. Cuff - from 1 to 100 pcs.

9. Medical pulse oximeter AORTIS® F014 consisting of:

- Medical pulse oximeter AORTIS® F014 -1 pc;
- Operation manual -1pc.

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4. Bag - from 1 to 100 pcs.
5. Mini USB cable - from 1 to 100 pcs.
6. Cuff - from 1 to 100 pcs.

10. Medical pulse oximeter AORTIS® F015 consisting of:

- Medical pulse oximeter AORTIS® F015 -1 pc;
- Operation manual -1pc.

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11. Medical pulse oximeter AORTIS® F016 consisting of:

- Medical pulse oximeter AORTIS® F016 -1 pc;
- Operation manual -1pc.

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12. Medical pulse oximeter AORTIS® F017 consisting of:

- Medical pulse oximeter AORTIS® F017 -1 pc;
- Operation manual -1pc.

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13. Medical pulse oximeter AORTIS® F020 consisting of:

- Medical pulse oximeter AORTIS® F020 -1 pc;
- Operation manual -1pc.

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14. Medical pulse oximeter AORTIS® F021 consisting of:

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6. Cuff - from 1 to 100 pcs.

19. Medical pulse oximeter AORTIS® F026 consisting of:
- Medical pulse oximeter AORTIS® F026 -1 pc;
- Operation manual -1pc.

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20. Medical pulse oximeter AORTIS® F027 consisting of:
- Medical pulse oximeter AORTIS® F027 -1 pc;
- Operation manual -1pc.

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21. Medical pulse oximeter AORTIS® F030 consisting of:
- Medical pulse oximeter AORTIS® F030 -1 pc;
- Operation manual -1pc.

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22. Medical pulse oximeter AORTIS® F031 consisting of:
- Medical pulse oximeter AORTIS® F031 -1 pc;
- Operation manual -1pc.

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23. Medical pulse oximeter AORTIS® F032 consisting of:
- Medical pulse oximeter AORTIS® F032 -1 pc;
- Operation manual -1pc.

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6. Cuff - from 1 to 100 pcs.

28. Medical pulse oximeter AORTIS® F037 consisting of:

- Medical pulse oximeter AORTIS® F037 -1 pc;
- Operation manual -1pc.

Accessories:

1. Battery - from 1 to 100 pcs.
2. Lace - from 1 to 100 pcs.
3. Cover - from 1 to 100 pcs.
4. Bag - from 1 to 100 pcs.
5. Mini USB cable - from 1 to 100 pcs.
6. Cuff - from 1 to 100 pcs.

29. Medical pulse oximeter AORTIS® K01T consisting of:

- Medical pulse oximeter AORTIS® K01T -1 pc;
- Operation manual -1pc.

Accessories:

1. Battery - from 1 to 100 pcs.
2. Lace - from 1 to 100 pcs.
3. Cover - from 1 to 100 pcs.
4. Bag - from 1 to 100 pcs.
5. Mini USB cable - from 1 to 100 pcs.
6. Cuff - from 1 to 100 pcs.

30. Medical pulse oximeter AORTIS® F01LT consisting of:

- Medical pulse oximeter AORTIS® F01LT -1 pc;
- Operation manual -1pc.

Accessories:

1. Rechargeable battery - from 1 to 10 pcs.
2. Lace - from 1 to 100 pcs.
3. Cover - from 1 to 100 pcs.
4. Bag - from 1 to 100 pcs.
5. Mini USB cable - from 1 to 100 pcs.
6. Cuff - from 1 to 100 pcs.

31. Medical pulse oximeter AORTIS® W01LT consisting of:

- Medical pulse oximeter AORTIS® W01LT -1 pc;
- Operation manual -1pc.

Accessories:

1. Rechargeable battery - from 1 to 10 pcs.
2. Lace - from 1 to 100 pcs.
3. Cover - from 1 to 100 pcs.
4. Bag - from 1 to 100 pcs.
5. Mini USB cable - from 1 to 100 pcs.
6. Cuff - from 1 to 100 pcs.

(hereinafter referred to as - product, pulseoxymetry device or device).

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1.2 Purpose: they are intended for percutaneous
measurement and display of hemoglobin saturation with
oxygen (SpO₂) and pulse rate.

1.3 Domain of usage: They should be used in healthcare
settings, medical preventive institutions and domiciliary

1.4 Operational conditions:

at the temperatures from 5 to 40°C and relative humidity
from 30 to 80%, at +25 °C and at lower temperatures
without moisture condensation; atmospheric pressure: 700-
1060 hPa (525-795 mm Hg).

1.5. Indications for use:

The device is suitable for the following patients:

- ✓ patients with vascular diseases, such as: coronary
heart disease, hypertension, hyperlipidemia, brain
thrombosis;
- ✓ patients with respiratory diseases such as: asthma,
trachitis, chronic bronchitis, chronic obstructive
pulmonary disease;
- ✓ senior citizens that are over 60 years old;
- ✓ people, who work more than 12 hours a day;
- ✓ people, who work with extreme physical loads or
under the conditions of alpine hypoxic
environment;
- ✓ people, who drink alcohol for a long time.

1.6 Contraindications:

There are no contraindications to pulse oximetry.

But do not use the device as follows:

1 Do not use pulseoxymetry device in nuclear magnetic
resonance environment or in a computerized tomography
environment.

2. It is forbidden to use pulseoxymetry device in situations,
when the alarm is required. This device is not equipped
with an alarm.

3. Explosive hazard: it is forbidden to use a pulseoxymetry
device in an explosive atmosphere.

4 The pulseoxymetry device only should be used as a
means to evaluate the patient's condition. The one should
be used in combination with other methods for evaluation
of both clinical signs and symptoms.

5. Time and again check the pulseoxymetry device sensor
application area in order to determine the position of the
sensor, the sensitivity of blood circulatory system and the
patient's skin.

6. The pulseoxymetry device sensor should not be fixed
with adhesive plaster. It can lead to device inaccurate
readings or skin blisters.

7. Before use, carefully read the Operations manual.

8. The pulseoxymetry device does not have SpO₂ alarms;
the one is not intended for continuous monitoring.

9. Prolonged use or patient's condition may require a
periodic change of sensor application place. Please, change
the place of sensor application, taking into account the skin

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integrity, the state of blood circulatory system and the alignment correctness, at least, every 2 hours.

10. Measurement inaccuracy may come out from autoclaving, ethylene oxide sterilization or sensors immersion into liquid.

11. SpO₂ measurements can be harmfully affected by the presence of bright ambient light; (protect the sensor area with a screen - a surgical towel, for example, from direct sunlight, if any)

12. Excessive patient movement.

13. Interference of high-frequency electrosurgical devices and defibrillators.

14. Venous pulsations.

15. Placement of a sensor on an extremity with a cuff for blood pressure measuring, an arterial catheter, or an intravascular line.

16. In the case when the patient is addicted to hypertension, severe vasoconstriction, severe anemia or hypothermia.

17. The patient is in the state of cardiopulmonary arrest or in shock.

18. Nail polish or false nails can cause inaccurate readings of SpO₂.

1.7 Measures of precaution:

Keep the device away from the sources of fire and / or heat.

- Avoid strong exposure to the pulseoxymetry device.
- Avoid intense direct sunlight exposure to the pulseoxymetry device.
- DO NOT leave the pulseoxymetry device in a closed vehicle.
- DO NOT open or disassemble the device. Opening and disassembling can damage the pulseoxymetry device.
- Use ONLY the charger that is supplied with the device. Other chargers SHOULD NOT be used. If this recommendation is not followed, there is a risk of pulseoxymetry device damage or danger to human health.
- DO NOT charge the battery for more than 5 hours, otherwise there is a risk of its damage.
- DO NOT charge or store the battery in extremely hot or cold environmental conditions. The most suitable temperature for charging is 5~35°C.
- Charge the pulseoxymetry device in due time, when the low-charge indicator appears. Otherwise, it may affect the pulseoxymetry device settings.
- Avoid contact with water. Store the pulseoxymetry device under the specified conditions. The accumulator storage battery will respond to liquid exposure, when it enters the pulseoxymetry device, and to exposure to high temperatures.
- Actions if snow, rain, sweat, etc. substances get on a pulseoxymetry device. You should stop using the device, otherwise it may catch fire. Then contact the technical support center and get the necessary

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assistance.

- It is not recommended to use pulseoxymetry device, while charging it.
- The pulse oximeter takes measurements at the level of the peripheral circulation (in the finger). Errors are possible in the assessment of oxygenation (it may be normal in the central bloodstream) if blood circulation is disturbed or slowed down in the limb on which the measurement is taken (blood loss, shock, cooling of the limb, vascular clamping in an uncomfortable position, tight clothing or blood pressure cuff.)

CAUTION: During long-term storage of the pulseoxymetry device, it is assumed that its battery will be discharged and require charging once every three months. The charge level of battery should be checked in both modes: storage and use.

1.8 Safety instructions

1 This device is not intended for treatment.

2 Before use of this pulseoxymetry device, carefully read safety, operation, and specifications manual.

3 Do not open the device case. Otherwise, you may damage the device. All maintenance and subsequent upgrades should be carried out only by the personnel that is trained and authorized by our company.

4 Do not use this device with a defibrillator or other high-frequency equipment.

5 Explosive hazard - Do not use this device in explosive atmosphere.

6. Prolonged use or patient's condition may require a periodic change of sensor application place. Please, change the place of sensor application, taking into account the skin integrity, the state of blood circulatory system and the alignment correctness, at least, every 2 hours.

7. The measuring data that are displayed on the device are for reference only and they cannot be directly used for diagnostic interpretation.

8. This device is not intended for use by infant and newborn.

9. In some cases, the device may interpret incorrect data when moving the finger. Minimize patient movement as much as possible.

10. The presence of bright light can lead to inaccurate SpO₂ measurements.

11. Only charge the device by connecting it to a specific device that meets the requirements of IEC 60601-1 for electrical safety, and make sure that the voltage and current of the specified device meet the requirements of the standard.

12. Do not connect it to any external device or use it for measurement, when this device is charging.

13. Do not use this device in situations, when an alarm is

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required. The pulseoxymetry device has no alarms.

14. The users, who have allergic reactions to rubber cannot use this product.

1.9 Risk class and applicable classification rules

Class of potential risk application 2b.

1.10 Medical device designer

Jiangsu Konsung Bio-Medical Science and Technology Co., Ltd

Address: NO.8, Shengchang West Road, Danyang Development Zone, Jiangsu Province, 212300, P.R.China

1.11 Manufacturer of a medical device

Jiangsu Konsung Bio-Medical Science and Technology Co., Ltd

production facility address: NO.8, Shengchang West Road, Danyang Development Zone, Jiangsu Province, 212300, P.R.China

1.12. General construction of medical device

Pulseoximeter is a portable device for measuring the saturation index. It consists of a sensor and a monitor. The sensor is attached to the finger or around the wrist. The sensor has a source of infrared and red light, these sources supply two beams that pass through the fabric.

Depending on how well the hemoglobin is saturated with O₂, the length of the light beam that it absorbs changes. The detector registers the light that has not been absorbed.

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 Full capacity of battery, the middle part shows capacity.

 Battery low, charge it immediately (or replace the batteries for model F03).

2. Power On / Off Button

Initiation of pulseoxymetry device: Press on / off button to turn on the device.

Switch off the pulseoxymetry device: Press this button for 2 seconds to switch off the device.

Press this button to switch on the display, when it is in low brightness mode.

3. Display screen

4. SpO2 value

5. PR value: heart rate per minute.


6. Bluetooth icon (not available in F03 model)

ON: means the connecting to Bluetooth device

Off: means disconnecting from Bluetooth device

7. Heart symbol

8. Heart Rate Histogram

1-Connector via mini USB cable for versions F01LT and W01LT

2-Finger compartment

The scope of delivery includes:

1. Medical pulseoxymetry device in the option - 1 pc.

2. Operations manual - 1 pc.

Accessories – if need

1.13 List of accessories

Battery, Rechargeable battery, Lace, Case, Bag, Mini USB cable, Cuff.

1.14 General description of key functional elements, parts, composition, functional capabilities

Pulseoxymetry device is a portable, battery-powered photoelectric device. Signals that are generated by light emitting diodes (LEDs) and received by the detector in the sensor or directly by the built-in detector are used for measurements, using the spectrophotometry principle. The

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Oxygen saturation is a term that specifies the proportion of oxygenated hemoglobin (HbO₂) as to total hemoglobin (Hb) in blood, i.e., oxygen saturation of the blood, which is an important physiological parameter of the respiratory and blood circulatory systems. Many respiratory diseases can cause oxygen saturation decreasing. Automatic regulatory dysfunction of the body, caused by anesthesia, serious surgical injury and injury that is caused by some physical examination, etc., will likely lead to a decrease of oxygen saturation that will lead to dizziness, weakness, vomiting and other symptoms in patients; severe cases will threaten life. SpO₂ is important in the scope of clinical medicine. And timely understanding of the patient's saturation of blood with oxygen will help the doctor find the problems. The device is worn on the tip of the finger and around the wrist (for versions W01T, W01LT).

The F01LT and W01LT versions are battery-powered, while all other pulse oximeter versions are powered by conventional AAA batteries.

Product Specifications:

SpO₂ display range: 0 ~ 100%
Measurement accuracy: 70-100% ±2%
Measurement discreteness 1%
Pulse rate display mode: Pulse Amplitude graph
Operating mode: continuous
Refresh Rate: 1 sec

The range of wavelengths in the maximum radiation and the maximum value of the light intensity of the pulse oximeter sensor emitters: Wavelength: red light: 660±3nm; infrared light: 940±5nm
Maximum power consumption, W: <30
Software version: 01 from 01.2019
Software security class: But
Water and dust protection type: IP21
Safety class: Class II type BF
Power cable length: 800mm
Noise level, Db: not applicable. No noise
Pulse measurement accuracy: 25 ~ 250, ±3 beats per minute
Measurement discreteness 1 beat per minute

Battery:

2pcs type AAA, 1.5 V alkaline battery
Working time: 30h

Rechargeable battery:

Typical Li-ion battery: 300 cycles
Power off delay: Max. 5 minutes after the low battery message appears.
Charging time: 2.5 hours (charge up to 90% capacity)
Battery: 300mAh, 3.7 V
For charging, use a mini USB 2.0 cable

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Any known interference that affects the operation or accuracy of the instrument: Temperature, humidity, physical excitation. Whether a hazard signal is given: yes, with a SpO₂ value of 90%.

The external surfaces of the product should not have violations of protective decorative coatings. On the surface of the products there should be no cracks, bubbles, chips, shells and other defects that impair their presentation and technical characteristics. Metal and non-metal parts of products must be made of a material that is resistant to corrosion or reliably protected from it.

The maximum permissible time for set of operating mode of pulseoxymetry device from the moment of switching on is 1 min.

The device should be mechanically strong and maintain operability in process and (or) after exposure to mechanical loads in accordance with ISO 80601-2-61.

The device should ensure safe operation in process and(or) after exposure to climatic factors:

- air temperature from +5°C to + 40°C;

- relative humidity from 30 to 80% at +25 °C and at lower temperatures without moisture condensation;

- atmospheric pressure: 700-1060 hPa (525-795 mm Hg)

- Recommended temperature for charging: from +5°C to + 35°C;

The pulseoxymetry device should comply with IEC 60601-1-2 requirements (group 1, class B)..

The cuff for the W01T, W01LT versions is a strap made of durable fabric. The cuff length is 20 cm.

The cord is used to wear a pulse oximeter around your neck. Length 35cm.

Case-designed to protect against scratches, splashes, bumps. size: corresponds to the size of the pulse oximeter.

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Bag-designed for storing, protecting and transporting the pulse oximeter.

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1.15 Description of the materials, contained in key functional elements / blocks / parts of a medical device, including materials, contacting with the human body

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2. Technical description of medical device

2.1 Safe use conditions

- 2.1.1 In order to avoid undesirable consequences, the pulseoxymetry device use is possible, when observing the safety procedures and excluding possible contraindications.
- 2.1.2 Before starting work, check the presence of physical damages of the devices and all connections of the parts to the device (batteries, cuff).
- 2.1.3 Operating conditions at the temperature from 5°C to + 40°C.
- 2.1.4 The procedure should be carried out by the personnel, who has familiarized with the operations manual.
- 2.1.5 The operational procedure with the pulseoxymetry device is described in the operations manual

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
2.2 Working procedure

Before unpacking, examine the packing case carefully for signs of damage. If any damage is detected, contact the carrier or us.


If the packing case is intact, open the package and remove the device and accessories carefully. Check whether there is any damage to the surface of the oximeter such as notches, dents, abrasions and so on.

Do not use it if the oximeter has mechanical damage.

Do not store it in a place where it is easy to splash water, direct sunlight, high temperature, humidity, dust, and corrosive gas, storage of chemicals or corrosive gases.

2.2.1 Install battery/ Rechargeable battery

Install two AAA batteries into battery compartment in correct polarities and put cover back on. (see Figure) or install the battery pack on the same principle as the batteries.

WARNING! Do not attempt to recharge normal alkaline batteries, they may leak and may cause fire or explosion.

2.1.2 Turn the Pulseoximeter on Put one of your fingers into rubber hole of the oximeter (as much area as possible) with nail surface facing up(as Figure), then release the clamp.

Press the ON/OFF button to turn the pulseoximeter on. Press the button and pulseoximeter will go into the working state.

SpO2 and pulse rate value displays on screen.

The oximeter will automatic fall into standby or sleep mode after 10 seconds without finger in it.

2.1.3 Read correspondent data from display screen (see

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
p.1.12).

NOTE:

1. The pulseoximeter cannot be measured while charging.
2. Finger can not apply nail polish and other cosmetics.
3. Do not shake the finger, hand or body during the measuring.
4. Jewelry or other items around the finger when measured can affect the accuracy of SpO2 and heart rate.

Install Lace

Let the thin end of the rope go through the lace hole, next let the big point of lace go through the hole, then tighten the lace.

WARNING: The pulseoximeter is NOT intended for diagnosis and the is only for warning, if there is any abnormal about the readings please go to the hospital for further examination. Indication indicates the severity of hypoxia during sleep, if this number is greater than 25, please go to the hospital for further examination.

2.3 Maintenance

For safe operation, it is recommended to carry out maintenance procedures. For safe operation of the pulseoxymetry device, the regular maintenance is required. Description and frequency of regular maintenance are given in the manual that came with the device. The frequency of inspections is specified, based on the calculation that the pulseoxymetry device is used 16-20 hours a day.

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+7 499 161 77 54

For damage that is not related to the user's actions and compliance with the operations manual of the pulseoxymetry device, contact the authorized representative of the manufacturer.

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2.4 Cleaning and disinfection

To avoid damage of the pulseoxymetry device, follow these rules:

Always dilute in accordance with the manufacturer's instructions or use as low as practicable concentration.

Do not immerse the pulseoxymetry device into liquid.

Do not pour liquid onto the pulseoxymetry device c accessories.

Do not allow liquid to enter the device case.

Never use abrasive materials (such as steel wool or silver polish) or erosive cleaners (such as acetone or acetone-based cleaners).

Warning

1 Before cleaning, turn off the pulseoxymetry device and disconnect the charging cable (or remove the batteries).

2 Do not disinfect the pulseoxymetry device with high temperature, high pressure, gas vapor, or liquid immersion. Clean and disinfect the oxygen according to manufacturer requirements.

3. If you spilled liquid on the device or accessories, wipe the

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immediately. If the pulseoxymetry device does not work
properly, contact the manufacturer or your service personnel.

2.4.1 Cleaning

Carrying out the following steps to clean the pulseoxymetry
device:

1. Without fail, turn off the pulseoxymetry device.
2. Clean outside surface of the pulseoxymetry device with soft
cloth, soaked in detergent ("Lotus" type).

3. If necessary, wipe out the entire cleaning solution with dry
cloth after cleaning.

4. Dry the pulseoxymetry device in a ventilated and cool place.

2.4.2 Disinfection

The manufacturer recommends to disinfect the pulseoxymetry
device only in case of need. Before disinfection, clean the
pulseoxymetry device.

The recommended detergent is ethanol (75%); the
recommended disinfection time is 300 minutes.

Note: never use ethylene oxide or formaldehyde for
disinfection.

Sterilization:

The device is not sterile.

2.5 Environmental requirements, including the requirements of safe destruction and disposal. Environmental protection

All production work should be carried out in the premises that
are equipped with local supply and exhaust ventilation. The use
of open flame in these premises is prohibited. Environmental
protection in accordance with applicable sanitary and
epidemiological rules and standards. The devices do not emit
toxic substances into environment and do not have irritating
sensitizing effect on humans.

Environmental protection in accordance with applicable
sanitary and epidemiological rules and standards of
manufacturer. Unorganized burning of used cushions is
prohibited. Disposal of device should be carried out in
accordance with applicable legislation.

Used devices should be disposed as Class A wastes, with
particular attention to the disposal of electrical components and
accumulator battery. The pulseoxymetry devices should be
operated in accordance with operations manual. Safety
requirements for the use of electrical components of the device
should be comply with IEC 60601-1-6, IEC 60601-1-2, IEC
60601-1.

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2.6 Operating life data

3 years

2.7 Transportation and storage requirements

The devices should be transported by all means of transport providing:

Temperature: from - 20°C to + 60°C.

Humidity: 10% - 95%.

Pressure: 700 - 1060 hPa.

In a well-ventilated premise.

Be careful, when transporting: avoid humidity and shaking.

It is necessary to store the pulseoxymetry device in a dry premise at the temperature from 5 to 40 °C. Shelf life is 12 months from the date of manufacture.

2.8 Manufacturer's Warranty

The manufacturer gives 12-months warranty for pulseoxymetry device, 6-month warranty for accumulator storage battery from the date of sale, which should be indicated on warranty card.

Correct use of device is warranty of its long service life.

Warranty service does not cover the repair of a damaged surface that is damaged due to improper use of the device.

Warranty is not valid in the following cases:

a) damage, caused by improper handling during transportation.

b) subsequent damage, caused by improper use or maintenance.

c) damage, caused by a change or repair by someone, who does not authorized by the manufacturer.

d) damage, caused by accidents.

e) replacement or removal of serial number label and production label.

If the product that is covered by this warranty is found to be defective due to defective materials, components or workmanship, and the warranty claim is submitted within the warranty period, the Manufacturer or the authorized representative of the manufacturer will repair or replace the defective part (- s) free of charge, at its own convenience.

2.9 Responsibility

Failure to comply with the paragraphs of the operation manual, as well as performing of repair work by persons without relevant qualifications, in particular, technical changes and refitting, without the consent of the manufacturer to do it may void the warranty, as well as responsibility for the products.

3. Product labeling and packaging

3.1. MI marking description

The marking is done by typographic or printing method on paper tag, with indication of:

- name of the product;

- product model designation;

- name and (or) trademark of the manufacturer;

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- catalog number (article);
- operating and storage conditions;
- date of manufacture;
- the symbols: "Refer to the usage instruction"; "Protect against
sunlight";
- registration certificate number;
Marking includes the marking according to IEC 60601-1.

If the size of the pulseoxymetry device does not allow to carry
out complete marking of the product, then, at least, the
following information should be applied to the pulseoxymetry
device:

- name of the manufacturer;
- serial number (or symbol 3.16 according to ISO 15223) or
batch or series identification number (or symbol 3.14 according
to ISO 15223) and
- the words: "Attention, take a look at accompanying
documents" or ISO 7000-0434 designation;
- the serial number (or symbol 3.16 according to ISO 15223) or
the identification number of a batch or series (or symbol 3.14
according to ISO 15223);
- detachable sensors should be marked in the form of a type
number and a batch number (or symbol 3.14 according to ISO
15223) or a serial number (or symbol 3.16 according to ISO
15223) that are affixed to them or to the packaging, if any;

- postprocessed pulseoxymetry device sensors should have the
same marking.

Displayed values

The pulseoxymetry device monitor should display the values of
functional oxygen saturation SpO_2 in percentas and have
marking: % SpO_2 or SpO_2 . The display image of the heart rate
should be expressed in inverse minutes units (1/min), for
example, heart rate/min. All other measured values should be
indicated with the appropriate units.

Do not reuse

If the sensor is not intended for reuse, then packaging or sensor
itself should have marking with an indication that the
pulseoxymetry device sensor should not be reused, or the one
should have a symbol 3.2 according to ISO 15223-1.

If the pulseoxymetry device sensor is intended for one patient
only, package or sensor should have marking, indicating that
the sensor is intended for one patient.

Where required, the safe operating deadline should be indicated
with the symbol 3.12 according to ISO 15223-1. The date
should include four digits to indicate the year and two digits to
indicate the month, as well as, if necessary, additional two
digits in order to indicate the day. The date should be indicated
next to the symbol.

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The pulseoxymetry device monitor and its parts should be
marked for their proper location.

Note - It is allowed to put on the marking both graphic symbols
and pictures, explaining the storage conditions of products,
their disposal, etc. according to ISO 15223-1.

The shipping marking includes the following:

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- address of the manufacturer;

- date of manufacture;

- "Refer to the usage instruction" symbols according to ISO
15223-1;

Marking of transportation packing should be carried out with
indication of the following handling marks: "Keep away from
moisture", "Fragile. Caution".

Each consumer package should include operations manual.

Other symbols are allowed to use according to ISO 15223-1.

3.2 Description of symbols that are put on the marking

Date of manufacture

Manufacturer

Caution! Refer to the instructions.

Serial number

BF applied part

Not intended for continuous monitoring.

IP 21 Degree of protection provided by the shells (IP code)

The symbol indicates that the device should be sent to
special facilities, according to the local regulations for separate
collection at the end of its useful life.

Fragile

Keep dry

It complies with EU Directive

3.3 Product packaging process

The pulseoxymetry devices are packed in polyethylene bag,
then in a cardboard box, ensuring products safety.

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Weight: from 6 to 16 kg

In a transport package from 20 to 150 pieces.

4. Production Information


4.1 Methods of control and testing

4.1.1 Control of suppliers and raw materials

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4.3 Manufacturing process of pulseoxymetry devices production

Prepare all the material; inspect the material; arrange the production; producing inspection ; inspection before delivery; packing.

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5. Significant principles correspondence table

5.1. The list of applicable national and international standards for the purpose to ensure safety, effectiveness and quality of a medical device - (minimal requirement)

- Quality management system EN ISO 13485
- Quality System Annex II - excluding Section 4 of Council Directive 93/42 / EEC
- Application of risk management to medical devices GOST ISO 14971: 2019

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5.2. International standards

ISO 14971:2019 Medical devices. Application of risk management to medical devices. Amendment 1. Validation of requirements

EN ISO 13485: 2016 Medical devices - Quality Management Systems. Normative Goals Requirements

IEC 60529: 2013 Degrees of protection, provided by enclosures (IP code)

IEC 60601-1: 2012 Medical electrical equipment Part 1. General safety requirements

IEC 60601-1-2: 2014 Medical electrical equipment. Part 1-2. General requirements to both basic safety and basic functional characteristics. Parallel standard. Electromagnetic

compatibility Requirements and Tests.

IEC 60601-1-6: 2013 Medical electrical equipment Part 1-6.

General safety requirements, taking into account the basic functional characteristics. Supplementary standard. Serviceability

ISO 10993-1: 2018 Medical devices. Assessment of the biological effects of medical devices. Part 1. Assessment and

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ISO 80601-2-61: 2017 Medical electrical equipment - part 2-61: special requirements for basic safety and basic characteristics of pulse oximeter equipment.

IEC 60068-2-31: 2008 Basic environmental testing procedures - Part 2: Tests - Test Ed: Free fall (Tests for environmental effect. Part 2. Tests. Test Ed: Free fall).

IEC 60068-2-64: 2008 Environmental testing - Part 2: Test methods - Test Fh: Vibration, broad-band random (digital control) and guidance (Tests for environmental effect). Part 2. Tests. The Fh: Broadband random vibration (digital control) and manual).

IEC 61000-4-2: 2008 Electromagnetic compatibility (EMC) - Part 4-2: Testing and measuring techniques - Electrostatic discharge immunity test GOST ISO 9919-2011 Medical electrical equipment. Particular safety requirements and basic characteristics of pulse oximeters (Electromagnetic compatibility (EMC). Part 4-2. Methods of measurements and tests. Electrostatic discharge immunity test).

5.4. The list of certificates, declarations of conformity with national, international standards

- Quality Management System EN ISO 13485: 2016
- Quality System Appendix II - excluding Section 4 of Council Directive 93/42 / EEC

6. There are no medicinal products, that are used in construction of pulseoxymetry devices.

7. Authorized representative of the manufacturer in the Russian Federation (name, address, phone, mail)

Limited Liability Company "Samed"
registered address: RF, 107497, Moscow, Ext. ter. g.
Golyanovo Municipal District, Shchelkovskoe sh., 79, room 1,
sq. 444, room 1
Phone: +7 499 161 77 54

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TECHNOLOGY CO., LTD.)*3211815015196/