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М.П. / Stamp

Instruction for Use
for the Medical Device

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97000230

Symbol description


- WARNING:** Failure to observe may result in injury or even death.
- CAUTION:** Failure to observe may result in damage to or even destruction of the product.
- NOTE:** Special information on the operation of the instrument.
- Consult instructions for use
- Manufacturer

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- ① Videoscope connector
- ② Sheath
- ③ Handle
- ④ Deflection lever
- ⑤ Distal tip
- ⑥ Serial port
- ⑦ E-Box connector
- ⑧ Serial port 1
- ⑨ Serial port 2
- ⑩ Mains connection
- ⑪ On/Off button
- ⑫ Charging indicator lamp (orange)

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U.S. Instruction

1 Intended use

Flexible CMOS Video-Endoscopes SSU are intended for visualization purposes during ENT procedures.

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2 Application area

Ears, nose and throat medicine.

NOTE: The product is sterilized using EtO.

NOTE: The product is packaged individually in protective packaging.

WARNING: The product is intended for a single use on a single patient. Otherwise, there is a risk of cross-contamination.

3 Indications

The CMOS Video Rhino-Laryngoscope SSU (091330-01) is intended to provide visualization of nasal lumens and airway anatomy (including nasopharynx and trachea) during diagnostic procedures. The E-Box serves as an adaptor for operating the flexible single-use videoscope on the compatible CCU.

The CMOS Video Esophagoscope SSU (091370-01) is intended to provide visualization of nasal sinuses, larynx, esophagus and gastro-esophageal junction during diagnostic procedures.

The E-Box serves as an adaptor for operating the flexible single-use videoscope on the compatible CCU.

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4 Contraindications

Rhinolaryngoscopy (091330-01):

The use of the instrument set for paranasal sinus surgery and for transnasal access for external skull base surgery is contraindicated if, in the opinion of the attending physician, the endoscopic surgical method as such is contraindicated, or if the patient is not able to undergo surgery or anesthesia due to his or her general condition. Instruments for paranasal sinus surgery and for transnasal access for external skull base surgery must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

The following contraindications also apply:

- A tendency to bleed and coagulation disorders which can neither be reversed nor treated.

Esophagoscopy (091370-01):

No contraindications relating to the product are known.

Use of the rigid and flexible esophagoscopy is contraindicated if, in the opinion of the responsible physician, the health of the patient is endangered by such use, due to the general condition of the patient or if the endoscopic method itself is contraindicated.

Esophagoscope and instrument sets must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

The following contraindications also apply:

- Contraindications for general anesthesia
- Coagulation disorders
- ASA > 3
- Multiple diverticula
- Esophageal stricture

5 Patient population

Men and women. Use is not limited to a certain patient profile (sex, weight, etc.), except for age.

NOTE: The attending physician decides which patients to use the product on.

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6 User qualification

KARL STORZ flexible videoscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the flexible laryngoscopy and esophagoscopy (including gastric cardia) technique. The instructions given in this manual are merely a guide for correct usage of the flexible videoscope. They are not intended as an introduction to the techniques of flexible laryngoscopy and esophagoscopy.

7 Compatibility

E-Box TP 012
C-MAC® Monitor 8403ZX
C-HUB® II camera control unit 20290320
TELE PACK+ TP101

8 Technical Data

091330-01:

Direction of view: 0°
Angle of view: 90°
Field of view: 110°
Working length: 30 cm
External diameter: 3.5 mm
Deflection up/down: 140°/140°

091370-01:

Direction of view: 0°
Angle of view: 90°
Field of view: 110°
Working length: 75 cm
External diameter: 3.5 mm
Deflection up/down: 210°/140°

9 Safety instructions

WARNING: Risk of injury and damage to the products: Failure to observe this instruction manual and the instructions for use of products used in combination can result in injury to patients, users and third parties, as well as to damage to the product. Please read all relevant instruction manuals and always follow the instructions for use precisely. Check the functioning of the products used in combination. Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended use and interface specifications of the products used in combination permit this.

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WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend, or malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original position.

WARNING: Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

WARNING: Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, adjacent components and tissue in front of the light emission window to heat up. This can cause burns to the patient, user and surgical accessories.

- Avoid direct tissue contact with the distal end.
- Do not rest the endoscope on the patient or where it is in direct contact with surgical accessories.
- The tip of the device can be safely focused in one stationary position when in direct contact with tissue for a maximum duration of 5 seconds.
- If the tip of the device is more than 5 mm away from tissue, there is no limit to the length of time it can stay focused in one stationary position.

WARNING: Risk of infection: Reprocessing a disposable product can jeopardize the functionality and safety of the product and lead to infections. You must not reprocess disposable products. Dispose of disposable products in accordance with the applicable regulations.

WARNING: This product is only intended for the specified purpose.

WARNING: Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope with a Type CF applied part is used. In this case, an endotherapy device with a Type CF applied part must be used in order to minimize the total patient leakage current.

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10 Safety precautions when using flexible videoscopes

The flexible videoscopes must only be used according to the medical rules and procedures of endoscopy recognized for this method.

NOTE: The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts. In this case, the reduction of the patient leakage current should be effected using only CF application parts as far as this is possible.

WARNING: Danger of explosion if flammable gases are used in the immediate proximity of the flexible videoscope. The flexible videoscope must not be used in the presence of flammable or combustible media.

CAUTION: Do not allow any live components or instruments to come into contact with the flexible videoscope.

CAUTION: Do not place the videoscope on devices that become very hot. Do not expose the device to direct sunlight or excessive heat, to avoid damage to the device.

CAUTION: Do not hit the distal tip of the flexible videoscope against hard objects, and never use a sharp object to remove dirt.

CAUTION: Do not squeeze the distal tip of the flexible videoscope too tightly.

CAUTION: Never place heavy objects on the flexible videoscope.

CAUTION: Never clamp or fold the sheath of the flexible videoscope. This can result in damage to the videoscope.

11 Electromagnetic Compatibility

This device has been tested and is found to comply with the EMC limit values as per EN/IEC 60601-1-2: CISPR 11 Class B.

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12 Handling

It is recommended that the suitability of the devices for the intended procedure be checked prior to use.

WARNING: Before use, ensure that the packaging is undamaged. Products with damaged packaging may not be used.

WARNING: Before use, check the expiration date. Products which have expired must not be used.

WARNING: Risk of infection: The use of non-sterile products poses a risk of infection for patients. Do not remove the products from the sterile packaging until immediately before use.

NOTE: When removing the product from the primary packaging, the rules of aseptic handling must be followed.

WARNING: The single-use product must only be assembled in a sterile area.

WARNING: A functional test must be performed before use on the patient.

12.1 General information

CAUTION: Do not carry the flexible videoscope by the sheath. Do not pull, clamp or twist the sheath of the flexible videoscope.

NOTE: Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

WARNING: If the image fails or there is interference during use, put the distal tip into the neutral position and remove from the patient's body slowly and carefully. If the flexible videoscope continues to be used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

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
12.2 Connection to the E-Box TP012

Insert the connector of the flexible videoscope ① into the serial port ⑥ of the E-Box.

CAUTION: Note the arrows on the connector; the elevation fits exactly in the recess in the serial port.

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
12.3 Connection to the C-MAC® monitor

The flexible videoscope can be used with the KARL STORZ C-MAC® monitor 8403ZX. The instruction manual for the monitor must be observed (see manual 96076008D).

Older monitor models (8402ZX/8401ZX) are not compatible.

Insert the E-Box connector (7) into one of the monitor's serial ports (8) or (9).


The connectors (1) and (7) can be connected or disconnected during operation. If necessary, connect the cord of the power supply up to the corresponding serial port (10) of the C-MAC® monitor and plug it into the power supply. The power cord can be fixed to the serial port using the union nut.

The battery is charged automatically as soon as the C-MAC® monitor is connected to the power supply. This is indicated by the lighting up of the orange charging indicator lamp (12).

12.4 Functions of the C-MAC® monitor 8403ZX

A more detailed description of the operation and functions of the monitor 8403 ZX can be found in the corresponding instruction manual supplied with every monitor (see manual 96076008D).

12.5 Connection to the C-HUB® II

The C-HUB® II (20290320) is a plug & play device for processing and transmitting image and control signals to external devices in medical video endoscopy. The flexible videoscope can be used in combination with the C-HUB® II and an external monitor.

To do so, insert the connector of the E-Box (7) into the serial port of the C-HUB® II.

A more detailed description of the operation and functions of the C-HUB® II can be found in the corresponding manual (see manual 96206529D).

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12.6 Deflection mechanism

The deflection lever can be used to angle the distal tip of the flexible videoscope. To do so, the deflection lever needs to be moved upwards (tip angled downwards) or downwards (tip angled upwards) using the thumb.

NOTE: The neutral position of the distal tip is straight, without any deflection whatsoever.

WARNING: The deflection lever must never be operated with force or abruptly to avoid injuring the patient or damaging the device!

WARNING: Before removing the flexible videoscope from the patient, the distal tip must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.

WARNING: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the flexible endoscope must not be used.

WARNING: If the flexible videoscope malfunctions during use on the patient, discontinue use immediately. Place the distal tip of the flexible videoscope in the neutral position and remove the flexible videoscope slowly and carefully from the patient.

CAUTION: The resistance that can be felt when deflecting the distal tip using the deflection lever increases as the angle of deflection increases. Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop

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12.7 Switching the flexible videoscope On/Off

Pressing the On/Off button ① on the C-MAC* monitor switches on the flexible videoscope with the integrated LED light source. When using the C-HUB* II, the flexible videoscope switches on automatically once the external monitor is activated.

WARNING: Risk of burns and glare! Never look directly into the light emitted by the flexible videoscope. The LED light emitted through the flexible videoscope may cause glare or high temperatures, and thereby burns, at the distal light outlet.

WARNING: The integrated LED light source will warm up during operation of the flexible videoscope. Switch off the flexible videoscope before putting it down, and do not place it in the vicinity of the patient (as per IEC 60601-1).

NOTE: The C-MAC* monitor can be switched on or off by pressing the On/Off button ①. The flexible videoscope does not have its own On/Off button.

Pressing and holding the On/Off button switches off the monitor together with the connected endoscope. If the softkey is actuated only briefly during operation, the monitor will switch to standby mode. This is indicated by the flashing LED on the monitor.

NOTE: The image focus is adjusted automatically. Manual adjustment is not necessary.

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12.8 Connection to the TELE PACK +

12.8.1 Switching on


To switch on the device, press the ON/OFF key on the front of the device and hold down for approx. 2 s.

12.8.2 Connecting the camera head / light cable/videoscopes

Insert the camera connector into the corresponding socket. Insert the light cable into the designated connection.

NOTE: Please note that 2 cameras can be connected to the device at the same time; however, only one camera can be displayed. If you have connected 2 cameras, you can select the relevant image with the "Swap Cameras" button in the menu.

NOTE: The use of original KARL STORZ light cables only is recommended. Light cables from other manufacturers may not be optimally adjusted for light transmission.

12.8.3 White balance

NOTE: Perform a white balance before every use.

The white balance can be performed directly after connecting the camera head with the TELE PACK + by pressing the White balance button in the menu.

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13 Reprocessing compatible devices

In order to reprocess compatible devices, please observe the following instructions for use below:

- C-MAC® Monitor 8403ZX: IFU 96076008D
- C-HUB® II: IFU 96206529D
- E-Box TP012: IFU 96286006D
- TELE PACK + TP101: IFU 96286007D

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14 Operation and storage conditions

Store in the original packaging. Protect against moisture and direct and indirect sunlight or any other form of UV radiation.

CAUTION: Do not place any objects on the product or the packaging.

CAUTION: Do not store in the proximity of chemicals, disinfectants or radioactive emissions.

NOTE: The shelf life and storage period can be found on the label.

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15 Disposal

Dispose of the product according to instructions after use. National laws and regulations must be observed.

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16 Complaint shipments

In the interests of our employees' health, only sterilized or disinfected products are accepted as return shipments. These should be marked as such; otherwise no further processing may be performed.

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17 Standard compliance

Provides BF protection when used with KARL STORZ E-Box TP012 and the C-MAC® Monitor 8403ZX with BF symbol on the connector.

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18 Directive compliance

This medical device bears a CE mark in accordance with the Medical Devices Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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19 Symbol description

Symbol	Meaning
	Manufacturer
	Date of manufacture
	Use-by date
LOT	Batch code
REF	Article number
STERILE EO	Sterilized with ethylene oxide
	Do not re-sterilize
	Do not re-use
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Consult the instructions for use
RxOnly	In accordance with US federal law (21 CFR 801.109), this product must only be sold on prescription from a licensed physician.

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A handwritten signature in blue ink, likely belonging to the acting notary.

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Addendum to Instruction for Use for the Medical Device

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
The information presented in this appendix to the operating instructions should be considered priority on the territory of the Russian Federation

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AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER:

LLC "KШТЭ ВОСТОК"
 115114, Moscow, Derbenevskaya embankment, 7, p. 4
 Phone: 8 (495) 983-02-40
 Email: info@karlstorz.com

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The medical device is intended for use in medical institutions.


Potential consumer

KARL STORZ flexible videoscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the flexible laryngoscopy technique.


Patient population

Men and women as well as children and adolescents. Use is not limited to a certain patient profile (sex, age, weight, etc.).

Main technical and functional characteristics

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Medical Devices and Non-Medical Devices, to Be Used in Combination with Medical Device Submitted for Registration

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Name	RC
Adapter E-Box TP012	In the process of registration
Endoscopic video camera C-HUB II	FSZ 2012/12443 from 19.12.2017

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	Manufacturer
	Date of manufacture
	Use-by date
LOT	Batch code
REF	Catalogue number
QTY	Number of products in the product packaging
STERILE EO	Sterilized using ethylene oxide
	Do not resterilize
	Do not reuse
	Do not use if packaging is damaged
	Validated sterile barrier symbol
	Validated sterile barrier and protective layer symbol

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ISO 14644-1	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness by particle concentration
ISO 14698-1	Cleanrooms and associated controlled environments - Biocontamination control - Part 1: General principles and methods
ISO 14698-2	Cleanrooms and associated controlled environments - Biocontamination control - Part 2: Evaluation and interpretation of biocontamination data
ISO 14971	Medical devices - Application of risk management to medical devices
ISO 15223-1	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

2 List of applied standards (Product Design Specification)

Identification	Designation
ANSI AAMI ST 72	Bacterial endotoxin - Test methods, routine monitoring and alternatives to batch testing
ASTM D 4169	Standard Practice for Performance Testing of Shipping Containers and Systems
ASTM D 4332	Standard Practice for Conditioning Containers, Packages, or Packaging Components for Testing
ASTM F 1886/ F1886M	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection
ASTM F 1929	Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration
ASTM F1980	Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices
ASTM F 2096	Standard Test Method for Detecting Gross Leaks in Packaging by Internal Pressurization (Bubble Test)
ASTM F 88/F88M	Standard Test Method for Seal Strength of Flexible Barrier Materials

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ISO 11135	Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices
ISO 11737-1	Sterilization of health care products - Microbiological methods - Part 1: Determination of a population of microorganisms on products
ISO 11737-2	Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
ISO 11138-1	Sterilization of health care products - Biological indicators - Part 1: General requirements
ISO 11138-2	Sterilization of health care products - Biological indicators - Part 2: Biological indicators for ethylene oxide sterilization processes
ISO 11607-1	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
ISO 11607-2	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes

3 List of partially applied standards

Identification	Designation/ Applied parts
IEC 60601-1-9	Medical electrical equipment - Part 1-9: General requirements for basic safety and essential performance - Collateral Standard: Requirements for environmentally conscious design

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The procedure for submitting complaints and responses to them is governed by civil law. A complaint may be brought only on such issues that were not the subject of acceptance of the goods made in accordance with the terms of the contract.

For all issues related to product maintenance, please contact the Authorized Representative of the manufacturer.

In case of complaint, contact the manufacturer and / or authorized representative of the manufacturer.

Manufacturer:

Karl Storz SE & Co.KG, Germany
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen, Germany
Tel. +49 7461 708-0, +49 7461 708-105
info@karlstorz.com

Authorized representative of the manufacturer:

Limited Liability Company “KARL STORZ - Endoscopy WOSTOK” (LLC “KSHTTE WOSTOK”)
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 (495) 983-02-40, Fax: +7 (495) 983-02-41
e-mail: info-ru@karlstorz.com

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