BE IT KNOWN that I Laura Anne Pawley, Mezzanine Floor, Epsom Square, Epsom, Surrey, KT19 8AG, a duly authorised Notary Public, duly admitted, sworn and practice throughout England and Wales.

CERTIFY ONLY that Brian Michael Howes, a United Kingdom Citizen, who is well known to me, who, is duly authorised by Mark Ainsworth of Reckitt Benckiser 'the Company' to represent them in this matter, has today caused the annexed Technical Documentation and Maintenance Documentation for the Medical Device Durex® Naturals Natural Rubber Latex Condoms for use in the Russian Federation to be produced to me and that he has represented to me on behalf of the Company that the said document is a true electronic copy of the original document.

SIGNED and sealed at Mezzanine Floor, Epsom Square, Epsom, Surrey, KT19 8AG on the 10th of December 2021.

Protocol No: 2021-1049

Laura Anne Pawley

Laura Anne Pawley
Notary Public

Dear Sir / Madam

I Mark Ainsworth of Reckitt Benckiser CERTIFY that the attached Technical and Maintenance Documentation for the Medical Device Durex ® Naturals condoms for use in the Russian Federation is a true electronic copy of the original document.

Name: Mark Ainsworth

Signature: *Mark Ainsworth*

Electronically signed by
Mark Ainsworth
Position: Technical File
Manager
July 2nd 2021 16:45
GMT

Date: 09-Dec-2021

Reckitt Group plc (Global Headquarters)
Turner House, 103-105 Bath Road, Slough,
Berkshire, SL1 3UH, UK

     17800 W: Reckitt.com
+44 (0)1753 717899

Reckitt refers to entities in the Reckitt Group plc group of companies.

Approved by:

Andrew Brownhill
Global Regulatory Affairs Manager

Date: Dec 5, 2021

Signature: 
Andrew Brownhill, Global Regulatory Affairs Manager

TECHNICAL DOCUMENTATION OF THE MEDICAL DEVICE

Durex® Naturals natural rubber latex condoms

Reckitt Benckiser Healthcare (UK) Ltd., UK

NAME OF THE MEDICAL DEVICE

Durex® Naturals natural rubber latex condoms

PURPOSE OF THE MEDICAL DEVICE

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Medical Device is intended to be used by laypersons.

MEDICAL DEVICE POTENTIAL RISK CLASS: 2b**STERILITY**

The device is non-sterile.

MEDICAL DEVICE DESCRIPTION

The condoms are made of natural latex, lubricated with water-based lubricant, smooth, colourless, easy-on shaped, teat-ended.

TECHNICAL CHARACTERISTICS

Technical parameter	Specification	Test Method used by the manufacturer	Frequency
Condom length, mm	195 ±10	ISO 4074:2015 clause 9.3.1 Annex D	Every batch
Condom width, mm (measured at 30 mm from the bead)	56 ±2	ISO 4074:2015 clause 9.3.2 Annex E	Every batch
Teat length, mm	Min 11	Appendix 1	Once a month
Thickness, mm	0.055 ± 0.005	ISO 4074:2015 clause 9.3.3 Annex F	Every batch
Bead thickness, mm	1-2	Appendix 2	Every batch
Condom weight (after removed from the primary package), g	2.00 ± 0.5	Appendix 3	Once a month
Elongation, %	Median 700 Min	ISO 4074:2015 Annex J	Once a month
Elongation, % (aged)	Median 600 Min	ISO 4074:2015 Annex J	Once a month
Freedom from holes	Impermeable	ISO 4074:2015 clause 12 Annex M	Every batch
Tensile strength, mPa	Median 20 Min	ISO 4074:2015 Annex J	Once a month
Tensile strength, mPa (aged)	Median 17 Min	ISO 4074:2015 Annex J	Once a month

RAW MATERIALS (% material is based on theoretical value)

Material	Nominal content of the material, % $\pm$ 10%	Trade name
Basic material:		
- Compounded Natural Rubber latex	68.97	Exceltex D
Dusting material:		
- Calcium Carbonate	2.31	Calcium carbonate light Bluestar
- Magnesium Carbonate	0.15	Magnesium Carbonate Extra Light
Lubricant:		
- Water	25.37	N/A
- Refine Glycerine	1.43	Refined Glycerine 99.5%' MIN EP/USP
- Phenethyl Alcohol	0.18	Conarom P2
- 1,2 Octanediol	0.06	
- Propanediol	0.89	Zemea Propanediol
- Xanthan Gum	0.20	Xanthan Gum FNCS-PC
- Sodium hyaluronate	0.029	Cristalhyal
- Lactic acid	0.023	Purac BF P/41
- Potassium lactate	0.054	

Raw Material	Additional Information
Compounded Natural Rubber latex	Exceltex D is a trade name of the natural latex used in compounded NRL
Phenethyl Alcohol 1,2 Octanediol	These two substances are the key antimicrobial aromatic blended substances of Conarom P2.
Lactic acid Potassium lactate	These two substances are the key pH adjusters substances Purac BF P/41

PACKAGING

The condom is packaged in an individual package of multi-layer material (foil).

Package dimensions:

Individual foil: 60 $\pm$ 2 mm x 60 $\pm$ 2 mm.

Primary package material:

Frontal foil layer: PET/AL/LDPE

Back foil layer:

PET/AL/LDPE

From the outer to the inner foil layer:

PET = Polyethylene terephthalate (or Polyethyl benzene-1,4-dicarboxylate)

AL = Aluminium

LDPE = low Density Polyethylene

PRESENTATION OF THE MEDICAL DEVICE

Secondary packaging is a carton pack, wrapped with the transparent tape where applicable. The quantity of condoms is indicated in the package.

PACKAGE CONTENTS OF THE MEDICAL DEVICE

The medical device is packed in a secondary / consumer package of 3 pcs. a pack or 12 pcs. a pack. It is accompanied by the instruction for consumers that can be enclosed in the package or printed in the inside and / or outside of the package.

LABELLING INFORMATION

Each individual package is labelled with:

- trademark;
- lot number (LOT);
- expiry date.

Symbols used for labelling and what they mean:

	Use by/Expiry date
	Lot number

Additional information can be printed on the individual package including reference to international standards.

The following information is applied to the outer surface of the secondary/ consumer package:

- condom description;
- quantity of condoms in consumer package;
- nominal condom length;
- nominal condom width;
- indication that the condom is lubricated;
- indication that the lubricant does not contain spermicide;
- indication that the condom contains natural latex;
- storage guidelines;
- name and address of the manufacturer;
- country of origin (if applicable);
- lot number;
- expiry date.

Symbols used for labelling and what they mean:

	Use by/Expiry date
---	--------------------

	Lot number
	Contains natural latex

Additional information can be printed on the outer surface of the secondary/ consumer package including reference to international standards.

Outer or inner surface of consumer packages and/or package insert contains prescribing information, which includes:

- reminder to handle the condom carefully, including when removing it out of the package;
- indications, contraindications, precautions for use, possible side effects;
- method of application;
- necessity to remove condom right after sexual intercourse;
- necessity to exclude petroleum- and oil-based lubricants;
- necessity to use safe additional lubricant if needed;
- necessity to consult a specialist on compatibility with topical medications;
- disposal instructions for used condoms;
- condoms are disposable products.

INDICATIONS

- prevention of unwanted pregnancy;
- protection from HIV and other sexually transmitted diseases.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

PRECAUTIONS FOR USE

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections. Check the expiry date on the condom wrapper before you use it. Do not use the condom after the expiration date. If the individual container is obviously damaged or deterioration in the properties of the condom is detected discard that condom and use a new one from an undamaged package.

Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil), some pessaries and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist:

- about the compatibility of condoms with other contraceptives;
- if you are using topical medicines applied to penis or vagina;
- if irritation/discomfort persists after use;
- if you are unsure if you are allergic to latex;
- if a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

For anal sex or if it is necessary use additional lubrication on the outside of the condom.

Condoms may present a choking hazard. Keep out of reach of children.

If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows.

Use a condom only once. Reuse may increase risk of failure or infection.

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which, like any component of the condom, might cause allergic reactions including anaphylactic shock if user is allergic to latex.

APPLICATION

Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.

Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

Dispose in a bin. Do not flush the condom down the toilet.

TRANSPORTATION CONDITIONS

Transportation temperature range guidance: -10°C up to +30°C (max. humidity of 65%RH for 30°C)

USAGE CONDITONS

There are no specific conditions required when the medical device is used for its purpose, indications, contraindications, usage application, taking into account possible side effects and precautions for use.

STORAGE CONDITIONS AND SHELF LIFE

Store the product in a dry place at 0-25 °C. Avoid mechanical impact and direct sunlight.

The condom must not come into contact with petroleum products, phenols and their derivatives, oils of petroleum origin, gasoline, kerosene and other organic substances.

Shelf life - 5 years.

RISK MANAGEMENT OVERVIEW

Risk assessment relating to the use of the medical device has been conducted in accordance with ISO 14971, Medical devices - Application of risk management to medical devices.

SUMMARY OF TESTS and RESULTS

BIOCOMPATABILITY TESTS

Test	International standard	Result
Cytotoxicity	ISO 10993-5	Cytotoxic found at: 100 (Undiluted), 50 and 25% dilution*. Non-cytotoxic: at dilution 12.5%, 6.25%, and 3.125%
Sensitization	ISO 10993-10	Non-sensitizer
Vaginal irritation	ISO 10993-10	Non-irritant
Penile irritation	ISO 10993-10	Non-irritant
Rectal irritation	ISO 10993-10	Non-irritant
Systemic toxicity	ISO 10993-11	No mortality

*ISO 10993-5 suggests that a response grade of moderate or severe should be considered a cytotoxic effect. The typical condom result has been found not to be associated with any problem in the regular use of condoms by the male or female partner. If the results of other testing (typically sensitization and irritation tests) and associated evaluations are seen as acceptable then the cytotoxicity result should be considered in context and the condom considered to meet the requirement of ISO 10993. The results indicate that Durex Naturals natural rubber latex condoms can be considered to be potentially biologically safe.

MANUFACTURING PROCESS

See Appendix 4

MEDICINAL OR PHARMACEUTICAL SUBSTANCE CONTAINED IN THE MEDICAL DEVICE

Not applicable.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

Compounded natural rubber latex contains casein (0.02% $\pm$ 10% of the contents of it, functions as a stabilizer) that is considered as a derivative from tissue/materials of animal origin as it is extracted from milk.

Casein (milk protein) is a natural by-product of milk. It has no known toxic effects although it may cause a reaction in individuals with an allergy to milk. The exposure from the condoms is considered not to pose any additional risk to users of the condoms.

The device does not contain any human blood derivatives.

STERILIZATION

The medical device is presented as a non-sterile device and is not intended to be sterilized thereafter.

SOFTWARE DESIGN, DEVELOPMENT AND VALIDATION

Not applicable.

MAINTENANCE AND REPAIR

Not applicable, as the product should be disposed of after a single use by a user.

GUIDELIENSS FOR DESTRUCTION AND DISPOSAL

The devices should be used in accordance with the instructions for use.

In household use, the product is disposed of as household waste. Do not put it down the drain.

In all other cases, the condoms are disposed of in accordance with local requirements.

ENVIRONMENT SAFETY REQUIREMENTS

With proper usage, storage and disposal Durex® Naturals natural rubber latex condoms have no negative impact on the environment.

Test Method used by the manufacturer

Appendix 1: Test method for condom teat length determination (D8291875)

- Sample preparation:
 - Dipped condom
Dry dipped condom before rolling will be used for this test. Stretch the condom slightly twice to smooth out the wrinkles.
 - Condom in container
Move the condom inside the individual container such that it is away from the area where the container is to be torn. Tear the container and remove the condom. Do not use scissors or other sharp instruments to open the package.
- Unroll the condom, stretch it slightly twice to smooth out the wrinkles. Lubricants may be removed by washing with a suitable solvent. If a solvent is used to remove the lubricant an adequate period of drying should be allowed before testing.
- Put the condom over the mandrel, having a rounded end with a diameter of (25 ± 2) mm with a scale divided into millimeters and dimensions with the zero beginning at the rounded end, and let it hang freely, stretched only by its own mass.
 - Mark a point at the junction of condom teat and the shoulder of the condom. See Figure 1 and Picture 1.



Figure 1. Length of condom teat

Picture 1. A point marked at the junction of condom teat and the shoulder.

- Remove the marked condom out from the mandrel. Lay it flat and put the ruler on top and measure the distance between the marked point and teat end by using a calibrated ruler and record the result.

Appendix 2: Test method for measurement of condom bead thickness (LAB/LT-ET/9, D0007024)

- Take a condom out of the pack or use a dried condom and spread it out gently.
- Measure the thickness of the bead by thickness gauge with the nearest reading of 0.01 mm.
- Repeat the measurement another 2 times.
- Average condom bead thickness for individual condom will be calculated.

Appendix 3: Test method for weight of dried condom + powder + lubricant (remove from foil) and Condom weight in primary package (D8314944)

- Weight individual foiled condom and record the weight (W1).
- Weigh a dry and clean petri dish and tare zero.
- Move the condom inside the individual package such that it is away from the area where the package is to be torn. Tear the package. Do not use scissors or other sharp instruments to open the package.
- Remove the condom carefully out of the opened foil package by using a dry and clean forceps and place it on the petri dish. Record the weight of condom (W2).
- Ensure that petri dish and forceps are dry and clean before testing next sample.
- Weight of foil package and remaining lubricant and powder after condom removal ($W3 = W1 - W2$)

PRODUCT MANUFACTURING FLOW DIAGRAM

Durex Naturals_TF_RTQ_Jan'2022


Final Audit Report

2021-12-05

Created:	2021-12-03
By:	Danyelle Jordan (Danyelle.Jordan@rb.com)
Status:	Signed
Transaction ID:	CBJCHBCAABAAhRKJunA-AN61gkwLCXfVgHprWuCgn5s9

"Durex Naturals_TF_RTQ_Jan'2022" History

-  Document created by Danyelle Jordan (Danyelle.Jordan@rb.com)
2021-12-03 - 11:27:22 AM GMT
-  Document emailed to Andrew Brownhill (andrew.brownhill@rb.com) for signature
2021-12-03 - 11:27:57 AM GMT
-  Document e-signed by Andrew Brownhill (andrew.brownhill@rb.com)
Signature Date: 2021-12-05 - 4:31:41 PM GMT - Time Source: server
-  Agreement completed.
2021-12-05 - 4:31:41 PM GMT

approved by:

Andrew Brownhill
Global Regulatory Affairs Manager

Date: Dec 5, 2021

Signature:

Andrew Brownhill (Dec 5, 2021 16:30 GMT)

**MAINTENANCE DOCUMENTATION
(Instruction for use)
OF THE MEDICAL DEVICE**

Durex® Naturals natural rubber latex condoms

Reckitt Benckiser Healthcare (UK) Ltd., UK

NAME OF THE MEDICAL DEVICE

Durex® Naturals natural rubber latex condoms

MANUFACTURER

Reckitt Benckiser Healthcare (UK) Ltd., 103-105 Bath Road, Slough, Berkshire, SL1 3 UH, UK.

MANUFACTURING SITE*

SSL Manufacturing (Thailand) Ltd., Wellgrow Industrial Estate 96 & 100, Moo 5, Bangna-Trad Road KM.36, Bangsamak, Bangpakong, Chachoengsao, 24180, Thailand.

**information about the place of manufacture is presented in the outer side of the secondary/consumer packaging in the format "The country of origin: [country]"*

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER

Reckitt Benckiser Healthcare LLC, Floor 3, 4 Shlyuzavaya emb., Moscow, 115114, Russia.
Phone: +7 (495) 961 25 65; fax: +7 (495) 961 25 66

PURPOSE OF THE MEDICAL DEVICE

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Medical Device is intended to be used by laypersons.

MEDICAL DEVICE COMPOSITION AND DESCRIPTION

The condoms are made of natural latex, lubricated with water-based lubricant, smooth, colourless, easy-on shaped, teat-ended.

TECHNICAL CHARACTERISTICS**

Technical parameter	Specification
Condom length, mm	195 ±10
Condom width, mm (measured at 30 mm from the bead)	56 ±2
Teat length, mm	Min 11
Thickness, mm	0.055 ± 0.005
Bead thickness, mm	1-2
Condom weight (after removed from the primary package), g	2.00 ± 0.5
Elongation, %	Median 700 Min
Elongation, % (aged)	Median 600 Min
Freedom from holes	Impermeable
Tensile strength, mPa	Median 20 Min
Tensile strength, mPa (aged)	Median 17 Min

***main technical characteristics are presented in the outer side of the secondary/consumer packaging as "nominal length – 195 mm; nominal width – 56 mm"*

PACKAGE CONTENTS OF THE MEDICAL DEVICE***

The medical device is packed in a secondary / consumer package of 3 pcs. a pack or 12 pcs. a pack. It is accompanied by the instruction for consumers that can be enclosed in the package or printed in the inside and / or outside of the package.

****quantity of condoms is presented in the outer side of the secondary/consumer packaging, the instruction for consumers is printed in the inner side of the secondary/consumer packaging*

LABELLING INFORMATION

Each individual package is labelled with:

- trademark;
- lot number (LOT);
- expiry date.

Symbols used for labelling and what they mean:

	Use by/Expiry date
	Lot number

Additional information can be printed on the individual package including reference to international standards.

The following information is applied to the outer surface of the secondary/ consumer package:

- condom description;
- quantity of condoms in consumer package;
- nominal condom length;
- nominal condom width;
- indication that the condom is lubricated;
- indication that the lubricant does not contain spermicide;
- indication that the condom contains natural latex;
- storage guidelines;
- name and address of the manufacturer;
- country of origin (if applicable);
- lot number;
- expiry date.

Symbols used for labelling and what they mean:

	Use by/Expiry date
	Lot number
	Contains natural latex

Additional information can be printed on the outer surface of the secondary/ consumer package including reference to international standards.

Outer or inner surface of consumer packages and/or package insert contains prescribing information, which includes:

- reminder to handle the condom carefully, including when removing it out of the package;
- indications, contraindications, precautions for use, possible side effects;
- method of application;
- necessity to remove condom right after sexual intercourse;
- necessity to exclude petroleum- and oil-based lubricants;
- necessity to use safe additional lubricant if needed;
- necessity to consult a specialist on compatibility with topical medications;
- disposal instructions for used condoms;
- condoms are disposable products.

INDICATIONS

- prevention of unwanted pregnancy;
- protection from HIV and other sexually transmitted diseases.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

PRECAUTIONS FOR USE

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections. Check the expiry date on the condom wrapper before you use it. Do not use the condom after the expiration date. If the individual container is obviously damaged or deterioration in the properties of the condom is detected discard that condom and use a new one from an undamaged package.

Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil), some pessaries and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist:

- about the compatibility of condoms with other contraceptives;
- if you are using topical medicines applied to penis or vagina;
- if irritation/discomfort persists after use;
- if you are unsure if you are allergic to latex;
- if a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

For anal sex or if it is necessary use additional lubrication on the outside of the condom.

Condoms may present a choking hazard. Keep out of reach of children.

If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows.

Use a condom only once. Reuse may increase risk of failure or infection.

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which, like any component of the condom, might cause allergic reactions including anaphylactic shock if user is allergic to latex.

APPLICATION

Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.

Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a

new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

Dispose in a bin. Do not flush the condom down the toilet.

TRANSPORTATION CONDITIONS

Transportation temperature range guidance: -10°C up to +30°C (max. humidity of 65%RH for 30°C)

STORAGE CONDITIONS AND SHELF LIFE

Store the product in a dry place at 0-25 °C. Avoid mechanical impact and direct sunlight.

The condom must not come into contact with petroleum products, phenols and their derivatives, oils of petroleum origin, gasoline, kerosene and other organic substances.

Shelf life - 5 years.

MEDICINAL OR PHARMACEUTICAL SUBSTANCE CONTAINED IN THE MEDICAL DEVICE

Not applicable.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN****

Compounded natural rubber latex contains casein (0.02% ± 10% of the contents of it, functions as a stabilizer) that is considered as a derivative from tissue/materials of animal origin as it is extracted from milk.

Casein (milk protein) is a natural by-product of milk. It has no known toxic effects although it may cause a reaction in individuals with an allergy to milk. The exposure from the condoms is considered not to pose any additional risk to users of the condoms.

The device does not contain any human blood derivatives.

*****warning about possible side effect in the form of the allergic reaction because of the natural latex contained in the medical device composition is presented in the inner side of the secondary/consumer package*

STERILIZATION

The medical device is presented as a non-sterile device and is not intended to be sterilized thereafter.

USAGE WITH OTHER MEDICAL DEVICES*****

Condoms can be used on their own to prevent from STIs and pregnancy and aren't intended to be used in conjunction with other medical devices.

However, consumers can use condoms with condom compatible lubricants if needed.

******recommendation to use compatible with condoms lubricant when it is necessary is presented in the inner side of the secondary/consumer packaging*

INSTALLATION, MAINTENANCE, CALIBRATION

Not applicable.

CLEANING, DISINFECTION, STERILISATION OF THE MEDICAL DEVICE

Not applicable.

INTERNATIONAL STANDARDS

ISO 13485

ISO 14971

ISO 4074

PRECAUTIONS TO BE TAKEN BY USER RELATED TO THE DISPOSAL OF THE DEVICE

Durex Naturals natural rubber latex condoms should be disposed in a bin. The condoms should not be flushed down the toilet.

DATE OF THE ISSUE / THE LATEST REVISION OF THE MAINTENANCE DOCUMENTATION

The medical device includes an indication of the latest revision within the instruction for use.

GUIDELIENSS FOR DESTRUCTION AND DISPOSAL

The devices should be used in accordance with the instructions for use.

In household use, the product is disposed of as household waste. Do not put it down the drain.

In all other cases, the condoms are disposed of in accordance with local requirements.

ENVIRONMENT SAFETY REQUIREMENTS

With proper usage, storage and disposal Durex® Naturals natural rubber latex condoms have no negative impact on the environment.

*Information
to be printed in the inner side
of the secondary/ consumer package*

Read the instruction carefully before use, especially if condoms are used for oral or anal sex.

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections.

ATTENTION: Check the expiry date on the condom wrapper before you use it. Do not use the condom after the expiration date. If the individual container is obviously damaged or deterioration in the properties of the condom is detected discard that condom and use a new one from an undamaged package.

PURPOSE/ INDICATIONS

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases. The medical device is intended to be used by laypersons.

APPLICATION

1. Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.
2. Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis.
3. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

4. Dispose in a bin. Do not flush the condom down the toilet.

For anal sex or if it is necessary use additional lubrication on the outside of the condom.

Condoms may present a choking hazard. Keep out of reach of children. If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows.

Use a condom only once. Reuse may increase risk of failure or infection.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which, like any component of the condom, might cause allergic reactions including anaphylactic shock if user is allergic to latex.

PRECAUTIONS FOR USE

Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil), some pessaries and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist:

- about the compatibility of condoms with other contraceptives;
- if you are using topical medicines applied to penis or vagina;

- if irritation/discomfort persists after use;
- if you are unsure if you are allergic to latex;
- if a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

Date of the latest revision of the instruction:

Durex Naturals_MD (IFU)_RTQ_Jan'2022

Final Audit Report

2021-12-05

Created:	2021-12-03
By:	Danyelle Jordan (Danyelle.Jordan@rb.com)
Status:	Signed
Transaction ID:	CBJCHBCAABAABkWGVL4tBP12BjEHH1YscE5zRKzRDb5u

"Durex Naturals_MD (IFU)_RTQ_Jan'2022" History

-  Document created by Danyelle Jordan (Danyelle.Jordan@rb.com)
2021-12-03 - 11:28:17 AM GMT
-  Document emailed to Andrew Brownhill (andrew.brownhill@rb.com) for signature
2021-12-03 - 11:28:42 AM GMT
-  Email viewed by Andrew Brownhill (andrew.brownhill@rb.com)
2021-12-05 - 4:27:50 PM GMT
-  Document e-signed by Andrew Brownhill (andrew.brownhill@rb.com)
Signature Date: 2021-12-05 - 4:30:11 PM GMT - Time Source: server
-  Agreement completed.
2021-12-05 - 4:30:11 PM GMT

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