Medtronic

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
tel: 001-763-514-4000
fax: 001-763-514-4879

To Whom It May Concern

It is hereby certified that the document attached is the 'User Manual' (M979336A006) for the Symplicity G3 Renal Denervation RF Generator.

Yours sincerely,

Anna Venegoni Dyer

Signed on behalf of Medtronic Inc.
Anna Venegoni Dyer
Sr. Regulatory Affairs Manager

10 May 2021

Date

This record was signed or attested to by use of communication technology on **10 May 2021** by **Varina V Philavong** who is located in Anoka County of Minnesota and that this record is to be filed with or relates to a matter before a court, government entity, public official, or other entity located in the territorial jurisdiction of the United States.

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
© 2018 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. *** Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

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Table of Contents

Table of Contents

1 Introduction	2
Overview	2
Contraindications	2
Warnings and Precautions	2
2 Product Overview	4
Product Description	4
System Components	4
System Notes	4
3 Installation, Setup, and Operation	6
Symplcity G3™ Generator Installation	6
Symplcity G3™ Generator Setup	6
Attaching the Dispersive Electrode	7
Treatment Procedure	7
4 Accessories	9
Cart (Optional)	9
Remote	9
5 Product Details	11
Front Panel	11
Touchscreen Interface	11
Side Panel Connections	23
Back Panel	24
Detection of previously used catheter	25
6 System Software	26
Catheter Compatibility	26
7 Audio and Display Information	27
8 Check Status Conditions	28
9 Fault Conditions	31
10 Troubleshooting	32
11 Cleaning	33
Symplcity G3™ Renal Denervation RF Generator	33
Foot Switch	33
Remote Control	33
Cart	33
12 Maintenance	34
Recommended Annual Electrical Safety Inspection: IEC82353	34
End of Life Disposition	34
13 Product Specifications	35
14 Technical Information	36
Essential Performance Characteristics	36
Power Output Diagrams	36
Electromagnetic Compatibility (EMC)	36
Electromagnetic Emissions per IEC 60601-1-2	37
Electromagnetic Immunity per IEC 60601-1-2	37
Limited Warranty	38

Chapter 1

Introduction

1

Overview

The Symplicity G3™ renal denervation RF generator is designed to deliver low-level radio frequency (RF) energy, via the Symplicity Spyral™ catheter through the wall of the renal artery to achieve renal denervation. Components include the disposable multichannel Symplicity Spyral™ catheter, the Symplicity G3™ generator, remote control, and an optional foot switch. The Symplicity Spyral™ catheter is used to deliver RF energy to a target site. The Symplicity G3™ generator provides an automatically controlled treatment of RF energy delivered by the catheter. The foot switch provides simple, hands-free RF activation of the system, while the remote provides a method for channel selection and RF activation from the sterile field.

The single-use Symplicity Spyral™ catheter is 6 Fr guide-catheter compatible. Using standard interventional techniques, the Symplicity Spyral™ catheter is placed in the renal artery, so RF energy can be delivered through the catheter's electrodes to the arterial wall and surrounding tissue where sympathetic nerves reside. Refer to the Symplicity Spyral™ catheter's Instructions for Use for further information about using the catheter.

The unit is specifically configured for use with electrodes conforming to the IEC 60601-2-2 standard. The Symplicity G3™ generator has been tested and found to comply with the electromagnetic emissions limits for medical devices given in IEC 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

The date of manufacture of the unit can be found by checking the unit's back panel.

Contraindications

There are no known absolute contraindications to the use of radio frequency surgery. Please refer to the Warnings and Precautions section for more details.

Warnings and Precautions

Related to the Use of Radiofrequency in Catheterization Laboratories

- Carefully read all instructions prior to use.
- Radio frequency surgery uses high-frequency output. Do not perform procedures if flammable or explosive media are present, such as flammable anesthetics or skin preparation agents.
- Interference produced by the operation of the Symplicity G3™ generator may adversely influence the operation of other electronic equipment.
- The Symplicity G3™ generator may produce hazardous electrical output and is designed to be used only by qualified medical personnel trained in the use of this equipment.
- Safe use of monopolar radio frequency energy demands proper attachment of a dispersive electrode to the patient. Follow all of the manufacturer's directions for skin preparation and the placement of the dispersive electrode and for proper insulation between the patient and any metal surfaces.
- Failure to achieve good skin contact by the entire adhesive surface of the dispersive electrode may result in a burn or measured high impedance. The dispersive electrode should be placed on the thigh or other nonbony area out of the angiogram field (if necessary, shave the placement area for good contact). Do not apply the dispersive electrode where fluid may pool.
- The entire area of the dispersive electrode should be reliably attached to a suitably prepared and appropriate area of the patient's body per the above instructions.
- Skin-to-skin contact (for example, between the arms and body of the patient) should be avoided. Use dry gauze to prevent skin-to-skin contact.
- The patient should not come into contact with metal parts that are grounded or have an appreciable capacitance to ground (for example, operating table supports). The use of antistatic sheeting is recommended.
- Monitoring electrodes of any physiological monitoring equipment used on the patient at the same time as RF energy delivery should be placed as far away as possible from the lower back region. Needle monitoring electrodes are not recommended. Physiological monitoring equipment incorporating high-frequency current-limiting devices is recommended. In all cases, monitoring systems incorporating high-frequency current-limiting devices are recommended.

Introduction

Related to RF Treatment

- Remove any guidewires that are not contained within the Symplicity Spyral™ catheter from the treatment site before activating the RF output.
- Implantable pacemakers (IPGs) and implantable cardioverter defibrillators (ICDs) or other active implants may be adversely affected by RF ablation. Consider deactivating IPGs and ICDs during ablation, having temporary external sources of pacing and defibrillation available during ablation, and performing a complete analysis of the implanted device function after ablation.
- Do not allow the Symplicity Spyral™ catheter electrode to come into contact with a metal instrument or surface during energy delivery. Superficial skin burns could occur.
- The Symplicity Spyral™ catheter and the Symplicity G3™ generator have not been evaluated in patients with prior renal angioplasty, indwelling renal stents, or aortic stent grafts.
- Do not touch the Symplicity Spyral™ catheter electrodes and the dispersive electrode at the same time during energy delivery. Superficial skin burns could occur.
- Initiate RF energy delivery only when the system is displaying the READY screen, as described in Chapter 5.

Related to Symplicity G3™ Generator and Ablation Catheter

- The Symplicity G3™ generator should be powered on and allowed to complete the system self-tests prior to introducing the Symplicity Spyral™ catheter.
- The Symplicity G3™ generator should be placed at least 30 cm (12 in) away from walls or corners.
- The Symplicity G3™ generator should sit for a minimum of 10 minutes with power ON before initiating a treatment.
- Use the Symplicity G3™ generator in a well-ventilated area. Do not drape the Symplicity G3™ generator because proper ventilation is required throughout the procedure.
- Check all external surfaces of the Symplicity G3™ generator, as well as components such as the power cord, foot switch, remote, etc. for mechanical damage prior to each use. The equipment should not be used if any damage is found. If inspection reveals a defect, the Symplicity G3™ generator should not be used until it has been properly repaired/serviced. The operator must immediately notify Medtronic of the defect.
- A failure of the Symplicity G3™ generator could result in an unintended increase of output power.
- To avoid the risk of electric shock, this equipment must only be connected to AC supply mains with protective earth.
- If the Symplicity G3™ generator malfunctions, stop using it and contact a Medtronic representative.
- The Symplicity G3™ generator is intended to be used only with the Symplicity Spyral™ catheter.
- Use the Symplicity G3™ generator only with the accessories in Chapter 2. Other active accessories or equipment may not have the appropriate ratings for the maximum output voltage of the Symplicity G3™ generator.
- Do not connect any cables or equipment to the Symplicity G3™ generator that are not specified by Medtronic. Doing so could potentially damage the system or cause injury.
- Do not remove the cover of the Symplicity G3™ generator as electric shock may occur. Refer servicing to qualified personnel. There are no user-serviceable parts inside the Symplicity G3™ generator.
- Do not disassemble the Symplicity G3™ generator. No modification of this equipment is allowed. This will void the warranty and may cause serious injury to the operator or damage to the unit. If any problems are not resolved by the directions in Chapter 10, contact Medtronic for further instructions.
- The covered connectors are for Medtronic use only.
- When connecting to an auxiliary secondary monitor, that is, DVI-D, the Symplicity G3™ generator should be connected to a monitor that is certified to IEC60950.

Note: Consult the catheter's Instructions for Use for additional warning and precautions.

Chapter 2

Product Overview

2

Product Description

The Symplicity G3™ generator is designed to deliver precisely controlled RF energy treatments using the Symplicity Spyral™ catheter. The Symplicity G3™ generator is an automatically controlled electrosurgical generator that produces RF energy according to a programmed algorithm. During energy delivery, power is continuously monitored and controlled based on temperature and impedance measured at the electrode-tissue interface to ensure proper operation. The Symplicity G3™ generator is compact for ease of use within the interventional suite. It is activated by using an optional foot switch, hand-held remote control, or through the interface on the Symplicity G3™ generator. The remote is non-sterile and needs to be draped if operated within the sterile field.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Symplicity G3 generator, including cables specified by Medtronic. Otherwise, degradation of the performance of this equipment could result.

Figure 1. Symplicity G3™ generator

System Components

Components provided with the Symplicity G3™ generator include:

- Hospital-grade AC power cord
- User manual
- Remote control
- DVI-D cable

Components provided separately:

- Symplicity Spyral™ catheter
- Optional foot switch

Additional accessories:

- Cart

System Notes

Use only Symplicity Spyral™ catheter with the Symplicity G3™ generator.

Product Overview

The Symplicity Spyral™ catheter and dispersive electrodes are intended for single use only. Do not attempt to sterilize and reuse them. Once used, catheters and dispersive electrodes should be treated as biohazards and disposed of in accordance with hospital policies, local governing ordinances, and recycling plans.

Chapter 3

Installation, Setup, and Operation

3

Note: Review all warnings and precautions prior to using the Symplicity G3™ generator.

Symplicity G3™ Generator Installation

Prior to initial use of the Symplicity G3™ generator, complete the performance tests listed in this section. This ensures that the Symplicity G3™ generator is operating properly when a new Symplicity G3™ generator is initially put into service, or when the Symplicity G3™ generator has been moved between facilities. There is no special equipment required for the performance tests outlined in this section.

Caution: Do not use the Symplicity G3™ generator if it does not meet the inspection criteria detailed below. Contact your local Medtronic representative if the system does not behave as indicated for the following performance tests:

1. Visual Inspection

Complete a visual inspection according the following instructions:

- Verify that the labels on the Symplicity G3™ generator are legible.
- Inspect the Symplicity G3™ generator, connectors, AC power cord, remote control, and optional foot switch for any evidence of damage.

2. System Self-test

Power on the Symplicity G3™ generator and verify the following:

- Verify that the display illuminates with the appropriate start-up screen.
- Verify that all indicator lights temporarily illuminate.
- Verify that the Speaker Test audio alert sound indicating that the speaker is functioning.
- Verify that the system self-test was successful and that no faults occur. See Chapter 9 for more information on faults.

3. Cooling Fan Operation

While the Symplicity G3™ generator is powered on, verify the following:

- Verify that the fan at the back of the Symplicity G3™ generator is functioning.
- Verify that the fan is free of dust or any foreign particles. If not, take a clean brush and dust the fan to provide air flow.

Symplicity G3™ Generator Setup

Note: Perform the following steps prior to accessing the patient's vasculature or opening/connecting a catheter:

1. If applicable and desired, attach the Symplicity G3™ generator to the cart by placing it into the latching mechanism on the top of the cart. Ensure the Symplicity G3™ generator latches into place and is secure within the cart.
2. If applicable and desired, connect the remote control to the back of the Symplicity G3™ generator. While in use, ensure the remote is properly draped in order to be accessed within the sterile field.
3. If applicable and desired, insert the foot switch connector into the receptacle on the back panel. Keep the foot switch in an area where it is easy to recognize and operate without any confusion between the foot switch for the Symplicity G3™ generator and any other foot pedal or foot switch present in the operating area.
4. Ensure the Symplicity G3™ generator is at least 30 cm (12 in) away from any corner or wall so that it can properly ventilate during a procedure.
5. Power on the Symplicity G3™ generator:
 - Verify the system self-test was successful (see the Symplicity G3™ Generator Installation section).
 - If the Symplicity G3™ generator does not complete the system self-test, do not use the Symplicity G3™ generator, (see Chapter 10). If the problem cannot be addressed, contact your Medtronic representative.
 - Following system self-test, the system is in STANDBY. No measurements or power controls are enabled.
6. Ensure the appropriate date and time settings are correct on the Symplicity G3™ generator screen. If necessary, adjust the date and time in the Settings menu.

Note: The Symplicity G3™ generator will not automatically update for local and annual time changes.

Installation, Setup, and Operation

Attaching the Dispersive Electrode

Prepare the patient using standard techniques for electrosurgery. Ensure that the patient's entire body, including extremities, is insulated from contact with grounded metal parts. A disposable dispersive electrode with a molded connector specified for adult use that complies with IEC 60601-2-2 is required, such as a Covidien™ Valleylab™ non-REM (return electrode monitoring) Polyhesive™ patient return electrode (REF E7506). Closely follow the instructions provided by the manufacturer of the dispersive electrode. Failure to achieve good skin contact by the entire adhesive surface of the dispersive electrode may result in a burn or measured high impedance. The dispersive electrode should be placed on the thigh or another nonbony area out of the angiogram field (if necessary, shave the placement area for good contact). Do not apply the dispersive electrode where fluid may pool.

Treatment Procedure

Prior to activating the RF output, verify that the following steps are taken:

1. Connect the dispersive electrode to the appropriate receptacle on the Symplicity G3™ generator.
2. Prepare the catheter per its Instructions for Use.
3. Connect the catheter to the Symplicity G3™ generator.
4. Select New Patient or Same Patient, as appropriate, on the touchscreen interface on the Symplicity G3™ generator.
5. Position the catheter to the targeted renal anatomy per the appropriate interventional procedure (see the catheter's Instructions for Use).
6. During catheter positioning, instructions to guide the workflow of the procedure will appear on the screen of the Symplicity G3™ generator.
7. When the electrodes of the catheter are not fully exposed to the vasculature, that is, if any electrodes are present within the guide catheter, the messaging screens will persist.
8. Once all catheter electrodes are positioned in the vasculature, distal to the guide catheter, screen messages will be replaced with impedance measurements. This is part of the software's "READY" state.
 - a. The impedance measurements should be used along with fluoroscopy to assess catheter positioning. The impedance plot for each channel provides the last 20 seconds of measurements for the respective channel.
 - b. Electrode stability is portrayed by the stability of the plot; electrode movement may be reflected by fluctuations of the plot.

Note: Electrode stability does not necessarily reflect adequate electrode contact against the arterial wall. Fluoroscopic imaging is required to ensure proper electrode contact.
 - c. If the screen continues to show messages and does not show impedance graphs and values after the electrodes have been positioned in the vasculature, then ensure that the catheter and the dispersive electrode are properly connected to the patient and generator. Proper and full contact of the entire adhesive surface of the dispersive electrode is required for operation of the generator and catheter; repositioning of the dispersive electrode may be required.
9. Electrodes may be deselected and selected using the respective channel numbers on the touch screen or remote control. Channel 1 reflects the most distal electrode of the catheter, while channel 4 reflects the most proximal electrode. When an electrode is deselected, RF energy is not supplied to the deselected electrode.
10. Electrode contact to the renal artery needs to be established so that the therapy is delivered to the site of the electrodes. Refer to the catheter's Instructions for Use for techniques to improve electrode contact with arterial wall.
11. Prior to initiation of RF energy, ensure that the RF channels and kidney demarcation are appropriately selected.
12. When it is appropriate to apply RF energy, press the foot switch, remote control RF button, or the RF button on the front of the Symplicity G3™ generator to begin RF energy delivery.
13. During the procedure, the display shows temperature, percent impedance drop, the starting impedance, and the remaining treatment time. The RF generator delivers power using an automated algorithm. When the treatment is complete, power delivery automatically stops. At any point in the procedure, treatment on all active channels may be stopped by pressing the foot switch, remote control RF button, or the RF button on the front of the Symplicity G3™ generator.

Note: If the impedance rises above the appropriate limits or if the temperature exceeds threshold levels, the Symplicity G3™ generator stops delivery of RF energy on the affected channels (see Chapter 8).

Chapter 3

14. From the Report screen, the treatment summary reports can be downloaded to an external storage device (such as SanDisk 2GB USB 2.0 flash drive, FAT32 formatted with at least 100MB of free space) via the USB port. In addition, from the Report screen treatment summary data can be viewed.
15. See the Symplicity Spyral™ catheter IFU for specific instructions regarding catheter use.
16. At the end of the procedure, disconnect the catheter, power off the Symplicity G3™ generator, disconnect the power cord from the power source, and fasten the power cord on the straps located on the back of the generator. Store the Symplicity G3™ generator in the appropriate storage location.

Accessories

4

Cart (Optional)

An optional wheeled cart is shown in Figure 2. If available, the cart can be used for keeping the Symplcity G3™ generator stable during therapy delivery and also to move the Symplcity G3™ generator from one location to another. When kept on the cart, the Symplcity G3™ generator automatically latches to the cart using its own weight. If the Symplcity G3™ generator has to be detached from the cart, a latch lever (located on the back of the cart) releases the Symplcity G3™ generator.

Figure 2. Symplcity G3™ renal denervation RF generator on a wheeled cart

Remote

As shown in Figure 3, a remote is available and can be used by connecting it on the back panel port provided for the remote. Refer to the Back Panel section for details. The remote can be stored in a detachable holster, which can be connected to the side panel of Symplcity G3™ generator as shown in Figure 4.

Figure 3. Remote Control

Chapter 4

Figure 4. Remote Control stored in holster on the side panel of the Symplicity G3™ renal denervation RF generator

The remote duplicates some of the control functionalities of the touch screen as mentioned.

1. The RF button on the remote control provides the same functionality as the RF button on the front panel of the Symplicity G3™ generator.
2. The buttons numbered 1, 2, 3, and 4 provide the same functionality as electrode selection buttons on the touch screen.
3. The L and R buttons (with kidney shadows) allow the user to toggle between the L and R buttons on the touch screen. This functionality is available only during the READY state, and provides an annotation in the ablation tab header indicating whether the ablation was performed in a renal artery leading to the left or right kidney. Labeling tabs is an optional step, and can be completed only prior to performing an ablation (it is not possible to label a tab once an ablation has been completed).
4. The arrow keys duplicate the function of the left and right tab buttons on the touch screen, which can be used to move between different treatment screens. Pressing the left arrow displays the previous treatment and pressing the right arrow displays the latest treatment.

Note: The remote is not a sterile device. If the remote is being used in the sterile area, then draping is required, as shown in Figure 5.

Figure 5. Remote control in draping for use in the sterile field

Product Details

5

Front Panel

Figure 6 shows the front panel display, RF button, and indicator lights.





Figure 6

1. **Display:** The display is the primary means of providing information to the user. During start-up and use, the status of the unit and energy delivery is displayed on these screens. The display has a touchscreen function that enables the user to interact with the system.
2. **RF Button:** Can be used to activate RF energy output while the system is in READY or deactivate RF energy output when the system is in RF ON (see Chapter 6). The RF Button functions the same as the remote or foot switch (see Foot Switch connection under the Back Panel description).
3. **Indicator Lights:** blue=RF ON; red=FAULT. See Chapter 7 for a functional description of indicator lights.




Touchscreen Interface

The display of the Symplicity G3™ generator contains a touchscreen interface that enables the user to interact with the system throughout the various states of the software, with the exception of the RF ON state. The touchscreen can be used to set the date and time, select a new patient or the same patient when a catheter is connected, deselect and select channels for RF activation, scroll through the ablation history plots of a given patient, select which kidney is being treated, access and download treatment reports and summaries, and other activities (see Chapter 6 for a description of the referenced states).

Chapter 6

Icon	Description
STANDBY STATE	
Connect Catheter screen	
	
The Connect Catheter screen is displayed in the STANDBY state of the system when no catheter is connected to the system.	
	The label points to the side panel where the Symplicity Spyrat™ catheter is attached. Only one catheter may be attached to the Symplicity G3™ generator at any given time.
	The Medtronic logo is displayed on the Connect Catheter screen.
	From the Connect Catheter screen, making this selection will enable the user to adjust the date, time, and language. (English is the only language available.)
	From the Connect Catheter screen, making this selection will enable the user to select, view, and download treatment summaries.



Product Details

Icon	Description
STANDBY STATE Settings screen 	
The Settings screen can be accessed from the Connect Catheter screen when no catheter is connected. This screen is accessible in the STANDBY state.	
	In the Settings screen, making selections in this field will modify the default language. Note: In order to change the setting, the SAVE selection must be made.
	In the Settings screen, making selections in these fields will modify the date. Note: In order to change the setting, the SAVE selection must be made.
	In the Settings screen, making selections in these fields will modify the time. Note: In order to change the setting, the SAVE selection must be made.
	In the Settings screen, these icons are used to scroll up or down through the values of date, time, and language selections. (English is the only language available).


Chapter 5

Icon	Description
STANDBY STATE	
	In the Settings screen, pressing the SAVE icon will save the currently selected date, time, and language. (English is the only language available).
	In the Settings screen, pressing the EXIT icon exits out of the Settings screen.
Patient Selection screen	
	
The Patient Selection screen is displayed once a catheter is connected in the STANDBY state.	
	In the Patient Selection screen, making this selection associates the connected catheter with a new patient which resets the ablation counter.
	In the Patient Selection screen, making this selection associates the connected catheter with the previously connected catheter's patient history and maintains the ablation counter from the last connected catheter.

Product Details

Icon	Description
READY STATE	
Message screen	
	
In the READY state, the Message screen is used to display troubleshooting messages.	
	In the Message screen, this symbol is followed by troubleshooting and guidance messages.
	In the READY state, the left and right arrows can be used to move between different treatment screens.
	In the READY state, the tab is used to view previously completed ablations. The indicator L or R within the tab indicates that Left or Right selection was made. The number represents the ablation count. If the tab is marked by an orange box, it represents that one or more channels ended in a treatment out of range or check status condition.
	From the READY state, making this selection will enable the user to select, view, and download treatment summaries.


Chapter 5

Icon	Description
READY STATE	
Ready screen	
	
In the READY state, the Ready screen is used to monitor the impedance of all active channels.	
 OR 	In the READY state, the icon indicates the most distal channel on the four-channel catheter. The icon is shaded green when selected or shaded gray when deselected.
 OR 	In the READY state, the icon indicates the second most distal channel on the four-channel catheter. The icon is shaded green when selected or shaded gray when deselected.

Product Details

Icon	Description
READY STATE	
 OR 	In the READY state, the icon indicates the second most proximal channel on the four-channel catheter. The icon is shaded green when selected or shaded gray when deselected.
 OR 	In the READY state, the icon indicates the most proximal channel on the four-channel catheter. The icon is shaded green when selected or shaded gray when deselected.
	In the READY state, the icon indicates the active channel on the single-electrode catheter.
	In the Ready screen, making this selection associates the current ablation counter with the renal artery leading to the left or right kidney.
	In the READY state, the left and right arrows appearing on the screen can be used to move between different treatment screens. The left arrow is pressed to display the previous treatment and right arrow is pressed to display the latest treatment.
	In the READY state, the text displays the current status of the RF delivery.
	This is the symbol of ohm associated with the impedance value displayed.
	In the READY state, the tab is used to view previously completed ablations. The indicator L or R within the tab indicates that Left or Right selection was made. The number represents the ablation count. If the tab is marked by an orange box, it represents that one or more channels ended in a treatment out of range or check status condition.
	From the READY state, making this selection will enable the user to select, view, and download treatment summaries.

Chapter 6

Icon	Description
RF ON state	
RF ON screen	
	
In the RF ON state, the RF ON screen provides feedback on the progress of the ablation.	
	In the RF ON screen, this icon indicates the most distal channel on the four-channel catheter.
	In the RF ON screen, this icon indicates the second most distal channel on the four-channel catheter.
	In the RF ON screen, this icon indicates the second most proximal channel on the four-channel catheter.
	In the RF ON screen, this icon indicates the most proximal channel on the four-channel catheter.

Product Details


Icon	Description
	In the RF ON screen, this icon indicates the single-channel catheter in RF ON state.
	In the RF ON state, the text displays the current status of the RF delivery.
	This icon appears next to the temperature value displayed in degrees Celsius.
	This is the percent icon associated with impedance drop from the baseline impedance value in RF ON condition.
	In the RF ON screen, these icons show if the impedance has decreased, stayed the same or increased from the baseline impedance.

Icon	Description
RF OFF state Check Status screen 	
	This icon is displayed to show that RF delivery has been shut off due to a check status condition.


Chapter 5

Icon	Description
	The channel that had the check status condition is displayed in orange color.
	These icons, which show if the impedance decreased, stayed the same, or increased from the baseline impedance, are grayed out to illustrate that the impedance will not change as RF delivery has stopped.
	The Stop icon is displayed on channels that are shut off due to a check status condition on a different channel.
	From the Check Status screen, making this selection will enable the user to select, view, and download treatment summaries.
User-Ended screen	
	
The User-Ended screen is displayed when a treatment is stopped by the user.	
	The Stop icon is displayed on channels due to a user-ended treatment.
	From the User-Ended screen, making this selection will enable the user to select, view, and download treatment summaries.


Product Details

Icon	Description
Report state	
Report screen	
	
The Report screen is used to display, select, and download treatment summaries.	
	In the Report screen, the left and right arrows can be used to move between different pages with sets of up to five reports per page.
	In the Report screen, the box icon indicates if the report associated with that data set is selected for download. When selected, the box icon will be filled with a check symbol.
	In the Report screen, making this selection will select all available treatment summary reports for download.
	In the Report screen, making this selection will download all selected treatment summaries to the USB memory stick connected to the USB-A port.
	From the Report screen, making this selection will exit out of the Report state.
	In the Report screen, touching the Summary icon will bring up the Treatment Summary screen.

Chapter 5

Icon	Description
Report state	
Treatment Summary screen	
	
The Treatment Summary screen is used to view treatment summaries of each ablation of the selected patient.	
	In the Treatment Summary screen, the left and right arrows can be used to move between different treatment summary pages within the same selected patient.
Electrode #	In the Treatment Summary screen, the label Electrode # refers to the electrodes in the four-channel catheter with 1 being the most distal and 4 being the most proximal. In the single-channel catheter, the fields under Electrode 1 are populated while others will be blank.
Baseline Imp	In the Treatment Summary screen, the label Baseline Imp refers to the impedance value of the channel when the treatment was initiated.
Max Imp Drop	In the Treatment Summary screen, the label Max Imp Drop refers to the maximum percent reduction in impedance seen during the treatment.
Max Temperature	In the Treatment Summary screen, the label Max Temperature refers to the maximum temperature reached during the treatment.
Max Power	In the Treatment Summary screen, the label Max Power refers to the maximum power delivered during the treatment.
End Power	In the Treatment Summary screen, the label End Power refers to the power level at the end of the treatment.
RF Time	In the Treatment Summary screen, the label RF Time refers to the time (indicated in seconds) that the therapy was delivered.
Kidney selection	In the Treatment Summary screen, the label Kidney selection displays R or L if the kidney selection button was pressed in the Ready screen for the associated treatment.
Code	In the Treatment Summary screen, the label Code will display Norm for a normal treatment or an associated Check Status Code per Chapter 8.

Product Details

Icon	Description
FAULT state Fault screen 	
The Fault screen is displayed when a Fault condition develops. 	This icon indicates that a fault has taken place (see Chapter 9 for more information).

Side Panel Connections

The side panel connections are shown in Figure 7. The Symplicity Spyral™ catheter and dispersive electrode connect to the Symplicity G3™ generator on the side panel. The USB-A port is enabled or disabled as per the software configuration.

Chapter 5

Figure 7. Side Panel Layout

1. **USB-A**: This port is for downloading treatment summary data to a memory stick.
2. **Type CF Applied Part (Defibrillator Proof)**: The Symplicity Spyral™ catheter connects to the Symplicity G3™ generator via the connector adjacent to each respective catheter symbol. The catheter is type CF.
3. **Refer to User Manual Symbol**: This symbol is for referring the user to the User Manual.
4. **Four-Channel Catheter Connector**: The four-channel catheter is connected on this port.
5. **Dispersive Electrode (Floating)**: The dispersive electrode connects to the Symplicity G3™ generator via this port. The dispersive electrode provides a path for the electrical current through the patient and back to the Symplicity G3™ generator. It is important to properly attach the dispersive electrode to the patient (see Attaching the Dispersive Electrode in Chapter 3).

Back Panel


Figure 8. Back Panel Layout

The back panel is represented in Figure 8. The power switch, AC power connector, Potential Equalization Terminal, and volume adjust are located on the back panel of the Symplicity G3™ generator.

1. Allows for projecting the display of the Symplicity G3™ generator onto an external monitor using a DVI-D cable.

Product Details

2.  **Remote connection:** Allows for the connection of the remote.
3.  **Volume Adjust:** Allows for adjustment of the output level of the volume-adjustable audio alerts.
4.  **Potential Equalization Terminal:** Provides common electrical ground for other electronic devices used during the procedure.
5.  **Foot Switch connection:** Allows for connection of the optional foot switch. The pneumatic foot switch is connected via the connector adjacent to this symbol. Depress the foot switch to begin delivering RF energy. To stop the energy delivery during operation, depress the foot switch again.
6. **AC Power Connector:** Connector for the AC power cable.
7. **Power Switch:** Turns the Symplicity G3™ generator on and off. Position the Symplicity G3™ generator so that the power switch is easy to operate to disconnect power from the Symplicity G3™ generator.
8. **Fuse Door:** For use by trained Medtronic service personnel only.
9. **Straps:** Secures any cord slack, or is used for storage of the cord on the back of the generator when the generator is not in use.

Detection of previously used catheter

If the feature is enabled, then the system will recognize a previously used catheter and display a message to the user to replace the catheter with a new catheter. Refer to Chapter 7.

Chapter 6

System Software

6

The Symplicity G3™ generator has the following system states that are controlled by software and user inputs. See Chapter 7 for audio alerts and information provided to the user in each state.

State	Description
STANDBY	System state after the Symplicity G3™ generator is turned on and has completed a successful system self-test. Entered only after self-test. Medtronic logo and software version are displayed.
READY	System State following STANDBY. Short pulses of nonharmful low-level RF power are delivered to allow for impedance measurements, which may aid in positioning the electrode. This state is accompanied by impedance plots showing on the display when the Symplicity G3™ generator is ready to initiate power delivery. An informational message may be displayed if certain conditions are detected. See Chapter 7 for more information.
RF ON	System state during RF power delivery that executes the programmed control algorithm. RF ON will be displayed along with the RF progress bar. If the programmed algorithm detects a measured treatment condition, for example, temperature or impedance out of its acceptable range, RF output power is automatically discontinued to the respective electrodes within this state, or the system transitions to the RF OFF state.
RF OFF	System state when RF power delivery is automatically discontinued after treatment time, or delivery is manually stopped by the user (through use of the RF button, foot switch, or remote), or stopped by the control algorithm. This system state is also reached when a combination of low temperature rise and impedance decrease is detected at the end of the treatment period (see Chapter 8).
REPORT	System state when treatment summary data is being accessed and viewed on the display or downloaded onto an external storage device.
FAULT	System state when the Symplicity G3™ generator detects an incorrect connection or setting, a self-test failure, or a failure in internal circuitry (see Chapter 9).

Catheter Compatibility

The Symplicity G3™ generator is designed to work specifically with the Symplicity Spyral™ catheter, and should only be used with Medtronic manufactured or recommended devices. Closely follow the Instructions for Use provided with the Symplicity Spyral™ catheter.

Audio and Display Information

Audio and Display Information

7

The Symplcity G3™ generator uses audio alerts, indicator lights, and the Display screen to communicate information to the user. The information conveyed to the user during a procedure by these methods is summarized below

Audio Alert and Display Information			
Display Message	Indicator Light	Audio Alert	Event
N/A [READY state]	"RF ON" pulsed	None	System in READY Low-level RF output pulses for impedance measurements
RF ON	"RF ON" lit (steady)	A repeating tone	System in RF ON and RF output power is being delivered.
RF ON (treatment nearing completion)	"RF ON" lit (steady)	A repeating tone faster than the standard RF ON tone to indicate that the treatment is about to complete	System in RF ON and RF output power is being delivered.
RF ON (with check status code during treatment)	"RF ON" lit (steady)	Two tones in short succession played over the standard RF ON tone	System in RF ON and RF output power is being delivered, where one or more channels have had RF discontinued
RF OFF (standard)	"RF ON" not lit	Three short tones	System transition to RF OFF after treatment has completed
RF OFF with check status code	"RF ON" pulsed	Two tones in short succession	System transition to RF OFF due to Check Status condition See Chapter 8.
FAULT ## Turn off power and Contact Medtronic, Inc	"FAULT" lit	One continuous beep (lasting ~5 sec.)	Fault condition exists See Chapter 9
CHECK CATHETER POSITION INSERT CATHETER	"RF ON" pulsed	None	Normally indicates that the catheter is partially in the guide catheter.
HIGH IMPEDANCE CHECK CATHETER POSITION OR DISPERSIVE ELECTRODE REPLACE CATHETER	"RF ON" pulsed	None	Indicates that the catheter has been detected but has not been inserted, or has been partially inserted, into the patient.
PREVIOUSLY USED CATHETER CONNECTED PLEASE INSERT A NEW CATHETER	"RF ON" pulsed	None	Indicates that the catheter is inserted into the patient but could still be in the guide catheter, or the dispersive electrode is not properly connected
MEDTRONIC CATHETER NOT RECOGNIZED PLEASE UNPLUG AND FULLY RE- INSERT	"RF ON" pulsed	None	Indicates a potentially faulty connection in either the catheter or the cable. The catheter should be replaced.
REPORT DOWNLOAD IN PROGRESS SELECTED REPORTS HAVE BEEN SUCCESSFULLY EXPORTED	N/A	None	If enabled in software this message indicates that a catheter used in a previous procedure has been attached to the Symplcity G3™ generator
MEMORY STICK NOT USABLE PLEASE INSERT A DIFFERENT MEMORY STICK	N/A	None	Indicates that the catheter has not been fully inserted into the receptacle. Remove and reinsert the catheter into the Symplcity G3™ generator
	N/A	None	Indicates that the selected reports are being downloaded to an external storage device
	N/A	One tone	Indicates that the selected reports have finished downloading onto the external storage device
	N/A	None	Indicates that the selected reports cannot be downloaded onto the external storage device. Consider using a different external storage device

User Manual English 27

global:20180610, altDocId:20170618, altDocType:20180728, master_pageSet:20180728, language:lookup:20181110, title:20140828 TOC:20180225, chapter:20170404, 3of9:20110211, 3of9_009:20110211

Chapter 8

Check Status Conditions

8

A Check Status condition occurs when the Symplicity G3™ generator detects a treatment condition outside of its acceptable range. When a Check Status condition is detected, the Symplicity G3™ generator will go into the RF OFF state or remain in the RF ON state, depending on the Check Status Condition, and:

- Prevent or stop RF energy delivery on the respective electrodes
- Sound the Check Status audio alert
- Display a Check Status Code and message

The following table provides the list of Check Status Codes and conditions that could occur during a treatment. Suggested actions for addressing the condition are provided. The Check Status or user-ended treatment message can be cleared by pressing the '>' key on the display or remote, or by pressing the foot switch two times within 3 seconds. If the condition persists after performing the suggested action, contact your Medtronic representative.

Symplicity G3™ Generator Check Status Conditions			
Check Status Code	Description	Action	Corresponding Condition on the Other Channels
20	Treatment impedance out of range: High	If Check Status code 20 occurs within 15 seconds of starting treatment: Ensure no excessive tissue/catheter electrode contact, that is, vessel distension and that catheter is not positioned in guide, a branch, or small artery. Reposition the catheter electrode if necessary and attempt to restart treatment. If Check Status code 20 occurs after 15 seconds of starting treatment: Check the catheter/dispersive electrode connections and/or replace catheter/dispersive electrode. Proceed as appropriate.	Other channels stop power delivery, and system transitions to the check status state.
21	Treatment impedance out of range: Low	Remove or disconnect the catheter. Turn unit off, then on again. Ensure that other devices are not interfering with the Symplicity G3™ generator. Proceed as appropriate.	Other channels stop power delivery, and system transitions to the check status state.
22a	Temperature too high	Image treatment site; wait approximately 1 minute and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Remove and inspect the catheter electrode. Proceed to the next treatment site as appropriate.	Other channels stop power delivery, and system transitions to the check status state.
22b	Temperature too high with indications of electrode movement	Image treatment site; wait approximately 1 minute and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Electrodes may have moved. Remove and inspect the catheter electrode. Proceed to the next treatment site as appropriate.	Other channels stop power delivery, and system transitions to the check status state.
23	Insufficient temperature rise	Check the catheter electrode position and reposition as necessary to improve contact with vessel wall.	Other channels continue power delivery. After RF activation is completed, check status message for the appropriate channel is shown.
24a	Unexpected impedance rise with indications of electrode movement	Image treatment site; wait approximately 1 minute and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Electrodes may have moved. Remove and inspect the catheter electrode. Proceed to the next treatment site as appropriate.	Other channels stop power delivery, and system transitions to the check status state.

Check Status Conditions

Simplicity G3™ Generator Check Status Conditions			
Check Status Code	Description	Action	Corresponding Condition on the Other Channels
24b	Insufficient impedance decrease	Electrodes may have moved. If not, reposition as necessary to improve contact with vessel wall. Proceed to the next treatment site as appropriate.	Other channels continue power delivery. After RF activation is completed, check status message for the appropriate channel is shown.
24c	Insufficient impedance decrease	Check the catheter electrode position and reposition as necessary to improve contact with vessel wall. Proceed to the next treatment site as appropriate.	N/A
26	Initial impedance out of range: High	Check that the catheter electrode is not within the guide catheter and at the desired treatment location. Ensure no excessive tissue/catheter electrode contact, ie, vessel distension and catheter is not positioned in a small branch. Reposition the catheter electrode to achieve lower starting impedance. Check catheter/dispersive electrode connections. Remove and inspect catheter electrode. Replace catheter/dispersive electrode if necessary. Proceed as appropriate.	Other channels stop power delivery, and system transitions to the check status state.
27	Initial impedance out of range: Low	Remove or disconnect catheter. Turn unit off, then on again. Ensure that other devices are not interfering with the Simplicity G3™ generator. Proceed as appropriate.	Other channels stop power delivery, and system transitions to the check status state.
28	Starting temperature too high	Ensure that other devices are not interfering with the Simplicity G3™ generator. Replace catheter.	Other channels stop power delivery, and system transitions to the check status state.
29	Starting temperature too low	Ensure that the catheter electrodes are not within the guide catheter and are within the desired treatment location. Also ensure that other devices are not interfering with the Simplicity G3™ generator.	Other channels stop power delivery, and system transitions to the check status state.
30	Power shutoff due to continual temperature rise	Image treatment site; wait approximately 1 minute and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Remove and inspect catheter electrode. Proceed to the next treatment site as appropriate.	Other channels continue power delivery. After RF activation is completed, check status message for the appropriate channel is shown.
32	Unexpected temperature decrease	Electrodes may have moved. Reposition as necessary.	Other channels continue power delivery. After RF activation is completed, check status message for the appropriate channel is shown.
33	Low temp during RF ON	Replace catheter.	Other channels stop power delivery, and system transitions to the check status state.
34	Unexpected temperature rise on an electrode which has been turned OFF (not delivering RF)	Image treatment site; wait approximately one minute and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Remove and inspect catheter electrode. Proceed to the next treatment site as appropriate.	Other channels stop power delivery, and system transitions to the check status state.
35	Catheter movement/torquing	Entire catheter may have moved. Image treatment site to determine current position of electrodes. Proceed to next treatment site.	Other channels stop power delivery.

Chapter 8

Symplicity G3™ Generator Check Status Conditions			
Check Status Code	Description	Action	Corresponding Condition on the Other Channels
36	Shorting of the electrodes (in the Symplicity Spyral™ catheter) during RF ON	The catheter may have moved. Image treatment site, wait approximately 1 minute and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Remove and inspect catheter electrode. Proceed to the next treatment site as appropriate.	Other channels stop power delivery, and system transitions to the check status state.
50	Suboptimal treatment detected while using the catheter, at the end of a normally Symplicity G3™ generator-terminated treatment	Move to next treatment site. At the next selected treatment site, ensure proper, stable wall contact between electrode and tissue prior to RF initiation. Use the starting impedance of the suboptimal treatment as a guide or reference to try to achieve a higher starting impedance at the new site, if possible.	No effect on RF delivery. Evaluation made at the end of a normally ending treatment.
90	Malfunctioning catheter	Replace the catheter.	All channels stop power delivery, and system transitions to the check status state.

Fault Conditions

Fault Conditions

9

A Fault condition occurs when the Symplicity G3™ generator detects an incorrect connection or setting, a self-test failure, or a failure in internal circuitry. When a Fault condition is detected, the Symplicity G3™ generator will go into the FAULT state and:

- Stop RF energy delivery
- Turn on the Fault Indicator light
- Attempt to display a Fault code
- Sound the Fault audio alert
- No longer allow user input

The following table provides Fault conditions and codes that could occur during a treatment. Suggested actions for addressing the condition are provided. Fault conditions can only be reset by turning off power. Prior to turning off power, record the Fault code number.

Symplicity G3™ Generator Fault Conditions		
Fault Code	Description	Action
Undefined	Internal software or hardware fault	Turn unit off and contact your Medtronic representative.
11	Button failure	A button on the front panel or the foot switch may have been pressed while the unit was powered on. Turn unit off, check buttons/foot switch, and then turn unit on. If condition persists, contact your Medtronic representative.
46	Symplicity G3™ generator internal temperature fault (onboard sensor)	If the internal Symplicity G3™ generator temperature is outside of the allowed operating temperature range, this fault may occur. The unit must be operated in an environment where the temperature is within the bounds noted in the Operating Environmental Conditions under Chapter 13. If the unit is transferred from an environment where the temperature is outside these bounds, it will take some time for the internal temperature to come to the temperature of its surroundings. If this fault occurs, let the unit sit for 5 to 10 minutes with the power ON (the hardware fault beep will stop after 10 seconds). Power the unit off and then back on. If the fault occurs again, let the unit sit for several more minutes and power the unit off and then back on again. If the fault occurs yet again, turn unit off and contact your Medtronic representative.
48	Symplicity G3™ generator internal temperature fault (CJC sensor)	If the internal Symplicity G3™ generator temperature is outside of the allowed operating temperature range, this fault may occur. Follow the same course of action given for fault code #46.
76	Catheter connected at start-up	If a catheter is connected when the Symplicity G3™ generator is powered on, the unit should be powered off. The catheter needs to be disconnected before the unit is powered back on. If condition persists, contact your Medtronic representative.
All other codes	Internal software or hardware fault	Turn unit off and contact your Medtronic representative.

Chapter 10

Troubleshooting

10

The following table provides a list of probable symptoms, which may occur during routine operation. If a problem is encountered that is not listed here or cannot be addressed with the suggested actions, contact your Medtronic representative.

Troubleshooting	
Symptom	Action
The display or indicator lights do not illuminate when the power switch is turned on	Ensure that the Symplicity G3™ generator is plugged into a working electrical outlet. Check the power switch on the rear of the unit. If problem persists, contact your Medtronic representative.
Fault indicator light comes on and a Fault code appears in message window	See Chapter 9.
Symplicity G3™ generator appears frozen and will not accept user input	In certain fault conditions the Symplicity G3™ generator may not be able to turn the Fault indicator on, beep, or display a Fault code. With the Symplicity G3™ generator in this state, RF output cannot be turned ON and user input is disabled. Disconnect the catheter from the Symplicity G3™ generator (if connected), turn unit off, and then back on again. If problem persists, contact your Medtronic representative.
RF power does not turn on when foot switch is depressed	Check that the Symplicity G3™ generator is displaying the READY screen, as defined in Chapter 5. Check the foot switch—unplug the foot switch from the Symplicity G3™ generator, place your thumb over the foot switch hose and depress the foot switch. If the foot switch is functioning properly, you should feel air pressure. Alternatively, RF power can be activated from the front panel of the generator.
RF interferes with ultrasound and other equipment	If the Symplicity Spyrat™ catheter along with the Symplicity G3™ generator, causes interference with other equipment, which can be isolated by turning the Symplicity G3™ generator off and on, try to correct the interference by one or more of the following measures: - Reorient or relocate the receiving equipment. - Increase the separation between the equipment. - Connect the Symplicity G3™ generator to an outlet on a circuit different from that to which the other equipment is connected. - Contact your Medtronic representative for assistance.
Impedance is less than 100 ohms	Remove catheter and inspect electrode and shaft; replace if necessary. Turn unit off, then on again. Ensure that other devices are not interfering with the Symplicity G3™ generator. If the problem persists, contact your Medtronic representative.
Unable to clear informational message, "INSERT CATHETER"	Confirm that the catheter is inserted into the guide catheter and into the patient's vasculature. Replace the catheter if message persists. If the problem persists, contact your Medtronic representative.
Unable to clear informational message, "CHECK CATHETER POSITION OR DISPERSIVE ELECTRODE"	Indicates that the catheter is inserted into the patient but could still be in the guide catheter, that the dispersive electrode is not properly connected, or that the system has encountered a high impedance condition on one or more of the channels. Check that all catheter electrodes are within the patient's vasculature, outside the guide catheter. Check that the dispersive electrode is properly connected to the Symplicity G3™ generator and properly coupled to the patient. Consider changing the location of the dispersive electrode pad on the patient. If the problem persists, contact your Medtronic representative.
"REPLACE CATHETER"	Indicates a potentially faulty connection in the catheter. Prior to replacing the catheter, consider checking the dispersive electrode connection and changing its location on the patient. Replace the catheter if message persists. If the problem persists, contact your Medtronic representative.
"CHECK CATHETER POSITION"	Indicates that the catheter is inserted into the patient but some electrodes may still be in the guide catheter. Verify that all catheter electrodes are within the patient's vasculature, outside the guide catheter. Check that the dispersive electrode is properly connected to the Symplicity G3™ generator and properly coupled to the patient. Consider changing the location of the dispersive electrode pad on the patient. If the problem persists, contact your Medtronic representative.

Cleaning

Symplicity G3™ Renal Denervation RF Generator

- Do not use steam or heat sterilization. Do not soak in any disinfectants or fluids. Do not allow liquid to enter into any of the electrical connections or the interior of the unit.
- Disconnect the unit from either the AC or power inlet before cleaning.
- To clean, wipe the exterior surface of the unit with 70 to 90% isopropyl alcohol only.
- Allow all surfaces and connections to dry before reconnecting the system.

Foot Switch

- To clean, wipe the exterior surface with 70 to 90% isopropyl alcohol only.
- Care should be taken to prevent fluid from entering the foot switch hose.

Remote Control

- Do not use steam or heat sterilization. Do not soak in any disinfectants or fluids. Do not allow liquid to enter into any of the electrical connections or the interior of the remote.
- Disconnect the remote from the Symplicity G3™ generator before cleaning.
- To clean, wipe the exterior surface of the remote and cable with 70 to 90% isopropyl alcohol only.
- Allow all surfaces and connections to dry before reconnecting the remote to the generator.

Cart

- Do not use steam or heat sterilization. Do not soak in any disinfectants or fluids. Do not allow liquid to enter the interior of the cart.
- Disconnect the Symplicity G3™ generator from the cart before cleaning.
- To clean, wipe the exterior surface of the cart with 70 to 90% isopropyl alcohol only.
- Allow all surfaces and connections to dry before reattaching the Symplicity G3™ generator.

Chapter 12

Maintenance

12

Recommended Annual Electrical Safety Inspection: IEC62353

Follow test instructions in the appropriate area of the standard. The acceptable limits (normal condition) for the Symplicity G3™ generator are as follows:

- Protective earth resistance (resistance between Mains Plug protective earth connector and protectively earthed accessible conductive parts): 0.3 ohm
- Earth (or Equipment) leakage current (all other AC operating voltages): 0.5 mA (500 µA)
- Enclosure (or Touch) leakage current: 0.1 mA (100 µA)
- Patient (or Applied Part) leakage current: 0.01 mA (10 µA)

End of Life Disposition

Do not dispose of the Symplicity G3™ generator or its accessories (cables, remote, cart, and foot switch) in the unsorted municipal waste stream. Follow local regulations for proper disposal.

Please refer to Medtronic website <http://recycling.medtronic.com> for proper disposal guidance.

Product Specifications

Product Specifications

13

RF Output Power	6.5 W maximum per electrode, ± 1.5 W accuracy
Impedance Measurement	175 Ω - 200 Ω , 28% 200 Ω - 1200 Ω , 15%
Fuse	2 Slow-Blow (T3.15A T, 250V), 20X5 mm
Temperature Range	Measurement 37°C - 65°C, $\pm 3^\circ\text{C}$
Input Power	100 to 240 V ~, 50/60 Hz universal power supply, 200 VA input power rating Disconnecting the AC power cord removes AC mains power from the unit
Dimensions	29.9 cm x 36.7 cm x 18.9 cm (11.8 in x 14.5 in x 7.5 in) (height x width x depth)
Weight	9.07 kg (20 lb)
Front Controls	Touchscreen buttons, RF button
Rear Controls	Power On/Off, Volume adjust
Display	1 touchscreen display on front panel, 1024x768 XGA monitor displaying graphics, messages, impedance, temperature, and time
Connections	AC power, isolated USB, foot switch, four-channel catheter connector, remote control connector, DVI-D connector, and dispersive electrode connector
Protection	Class I, Defibrillator Proof—Type CF, intermittent operation with a duty cycle of 995 seconds on, 60 seconds off. The Symplixity Spiral™ catheter and Symplixity G3™ generator are not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide. IEC 60601-1 Class I, IEC 60601-2-2 Safety of high-frequency surgical equipment, IEC 60601-1-2 Electromagnetic compatibility (EMC)

Operating Environmental Conditions

Temperature	15°C to 40°C (59°F to 104°F)
Humidity	30% to 70% relative humidity, noncondensing
Pressure	70 to 106 kPa [-0.7 to 1.05 ATM]

Storage/Transit Conditions

Temperature	-35°C to +57°C (-31°F to 135°F)
Humidity	30 to 75%
Pressure	58.5 to 106 kPa [-0.585 to 1.05 ATM]

Chapter 14

Technical Information

14

Essential Performance Characteristics

The Symplicity G3™ generator functions safely, properly, and maintains default conditions. The essential performance characteristics related to the safety of the Symplicity G3™ generator, as tested in compliance with IEC 60601-1 and IEC 60601-1-2, include the proper function of the following features:

- The operator should be able to stop delivery of power.
- The system shall be able to accurately measure temperature and impedance per system specification.
- Accurate delivery of RF energy per the defined power delivery algorithm.
- Safety features associated with minimizing the risk of electric shock to the patient and user.
- Hardware related Power On Self-Test (POST) activities and fault recognition.
- Internal and external watchdog timers.
- The system must not allow the user to alter the RF energy delivery algorithm, that is, by adjusting the clinical and system operating parameters through the front panel interface or via the communication ports with an external device, or by any other means.

Note: The Symplicity G3™ generator display is not included as part of Essential Performance. The display is for informational purposes only and is not used for treatment determination. Transient loss of display readability or the display of unintended characters, for example, due to interference, would not impact the safety of the procedure. Permanent loss of display could result in a terminated procedure, but would not result in an unacceptable risk to the patient.

Power Output Diagrams

Figure 9 shows the maximum power from the catheter as a sum of wattage from the active channels. Each individual electrode of the Symplicity Spyral™ catheter can deliver a maximum of 6.5W. The maximum output voltage is 150Vp.


Figure 9. Symplicity Spyral™ Catheter Total Power Delivered per Active Channel

Note: Total Power delivered by the Symplicity Spyral™ catheter is presented as a function of the number of active electrodes.

Electromagnetic Compatibility (EMC)

Interference produced by the operation of high-frequency surgical equipment, such as the Symplicity G3™ generator, may adversely influence the operation of other electronic medical equipment such as monitors and imaging systems. Special precautions should be taken. The Symplicity Spyral™ catheter and Symplicity G3™ generator should be installed and put into service according to the EMC information provided in this section.

Warning: The Symplicity Spyral™ catheter and Symplicity G3™ generator are intended for use by healthcare professionals only. The Symplicity G3™ generator may cause radio interference and disrupt the operation of nearby

Technical Information

equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the Symplicity G3™ generator or shielding the location. Interference with electronic medical equipment, such as monitors and imaging systems, is usually resolved or minimized by rearranging the cables such that the unit cables do not overlap with the cables from the monitoring equipment.

Warning: The Symplicity G3™ generator should not be located immediately adjacent to or stacked with other equipment. If adjacent or stack use is necessary, the Symplicity G3™ generator should be observed to verify normal operation in the configuration in which it will be used. Portable and mobile RF communications equipment can affect the Symplicity G3™ generator.

Warning: The AC power cord set (AC power cable) used with the Symplicity G3™ generator should be non-shielded, have a maximum length of less than 3 meters, be at least 18 AWG (0.823 mm²) in USA/Canada, or 0.75 mm² in other countries and be approved for use in the country for which it is intended.

Warning: Use of accessories, transducers, and cables, other than those specified or provided by Medtronic, could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Note: The Emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Electromagnetic Emissions per IEC 60601-1-2

The Symplicity Spyrat™ catheter and Symplicity G3™ generator are suitable for use in the electromagnetic environment specified below. The user should ensure that it is used in such an environment.			
Emissions test	Compliance	Electromagnetic environment—guidance	
RF emissions CISPR 11 (EN55011)	Group 1	When the Symplicity G3™ generator is in the STANDBY state, RF energy is only used for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.	
RF emissions CISPR 11 (EN55011)	Class A	During the RF ON state, the Symplicity G3™ generator must emit electromagnetic energy to perform its intended function. Nearby electronic equipment may be affected.	
Harmonic emissions IEC 61000-3-2	Class A	The Symplicity Spyrat™ catheter and Symplicity G3™ generator are suitable for use in all establishments other than domestic, that is, residential and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded.	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the Symplicity G3™ generator or shielding the location.	

Electromagnetic Immunity per IEC 60601-1-2

The Symplicity Spyrat™ catheter and Symplicity G3™ generator are suitable for use in the electromagnetic environment specified below. The user should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment—guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	AC mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV, ±1 kV line to line & ±0.5 kV, ±1 kV, ±2 kV line to ground	±0.5 kV, ±1 kV line to line & ±0.5 kV, ±1 kV, ±2 kV line to ground	AC mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11 UT=230 V _{an}	0% UT (100% dip in UT) for 0.5 cycle 0% UT (100% dip in UT) for 1 cycle 70% UT (30% dip in UT) for 25/30 cycles 0% UT (100% dip in UT) for 5 s	0% UT (100% dip in UT) for 0.5 cycle 0% UT (100% dip in UT) for 1 cycle 70% UT (30% dip in UT) for 25/30 cycles 0% UT (100% dip in UT) for 5 s	AC mains power quality should be that of a typical commercial or hospital environment. If the user of the Symplicity G3™ generator requires continued operation during AC power mains interruptions, it is recommended that the Symplicity G3™ generator be powered from an uninterruptible power supply or a battery.

Chapter 14

The Symplicity Spyr™ catheter and Symplicity G3™ generator are suitable for use in the electromagnetic environment specified below. The user should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment—guidance
Power frequency (50 Hz) magnetic field IEC 61000-4-8 NOTE UT is the a.c. mains voltage prior to application of the test level.	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms ISM/Amateur Radio bands inside 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz 6 Vrms ISM/Amateur Radio bands inside 150 kHz to 80 MHz	WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Symplicity G3™ generator, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz RF communication equipment inside 80 MHz to 6 GHz	3 V/m 80 MHz to 2.7 GHz RF communication equipment inside 80 MHz to 6 GHz	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b .
NOTE These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Symplicity G3™ generator is used exceeds the applicable RF compliance level above, the Symplicity G3™ generator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Symplicity G3™ generator.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Limited Warranty

MEDTRONIC, INC. WARRANTS THAT THE SYMPPLICITY G3™ RF GENERATOR, HEREAFTER REFERRED TO AS "GENERATOR," SHALL BE FREE FROM MATERIAL DEFECTS AND DEFECTS IN WORKMANSHIP FOR A PERIOD OF 12 MONTHS COMMENCING ON THE DATE OF DELIVERY TO THE CUSTOMER, PROVIDED THAT THE CUSTOMER PROVIDES MEDTRONIC WITH WRITTEN NOTICE OF THE DEFECT PRIOR TO THE EXPIRATION OF THE WARRANTY PERIOD AND RETURNS THE GENERATOR TO MEDTRONIC NO LATER THAN 10 DAYS FOLLOWING THE EXPIRATION OF THE WARRANTY PERIOD.

MEDTRONIC, INC. WILL, AT ITS DISCRETION, EITHER:

- WITHOUT CHARGE TO THE CUSTOMER, REPAIR THE GENERATOR AND BRING IT INTO FULL CONFORMANCE WITH ITS SPECIFICATIONS, INCLUDING RE-TEST/RE-INSPECTION BY MEDTRONIC, INC.; OR,
- WITHOUT CHARGE TO THE CUSTOMER, REPLACE THE GENERATOR WITH A NEW DEVICE THAT IS FUNCTIONALLY COMPARABLE TO THE GENERATOR.

TO QUALIFY FOR THIS LIMITED WARRANTY, ALL OF THESE CONDITIONS MUST BE MET:

- THE GENERATOR MUST BE USED IN ACCORDANCE WITH THE LABELING, INCLUDING BUT NOT LIMITED TO THE INSTRUCTIONS FOR USE, RECOMMENDATIONS AND MANUALS;
- THE GENERATOR MUST NOT BE SUBJECTED TO UNAUTHORIZED ALTERATION, MODIFICATION OR REPAIR;
- THE GENERATOR MUST NOT BE SUBJECTED TO MISUSE, ABUSE, ACCIDENT, OR IMPROPER HANDLING;
- THE GENERATOR MUST BE PROPERLY MAINTAINED IN ACCORDANCE WITH THE MAINTENANCE SCHEDULE AND USED SOLELY IN CONJUNCTION WITH PRODUCTS APPROVED BY MEDTRONIC, INC. FOR USE WITH THE GENERATOR;
- THE GENERATOR MUST BE RETURNED TO MEDTRONIC IN ACCORDANCE WITH THE RETURN GUIDELINES STIPULATED IN THE USER MANUAL; THE GENERATOR MUST NOT BE MOVED OUTSIDE THE CUSTOMER'S SITE WITHOUT PRIOR NOTICE AND APPROVAL FROM MEDTRONIC, INC.

THIS LIMITED WARRANTY IS LIMITED TO ITS EXPRESS TERMS. IN PARTICULAR:

- EXCEPT AS EXPRESSLY PROVIDED BY THIS LIMITED WARRANTY, MEDTRONIC, INC. IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT,

Technical Information

FAILURE, OR MALFUNCTION OF THE GENERATOR, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

- B. THIS LIMITED WARRANTY IS MADE ONLY TO THE CUSTOMER TO WHOM THE GENERATOR WAS ORIGINALLY SOLD. AS TO ALL OTHERS, MEDTRONIC, INC. MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOM, OR OTHERWISE. NO SUCH EXPRESS OR IMPLIED WARRANTY TO THE CUSTOMER SHALL EXTEND BEYOND THE TWELVE (12) MONTH PERIOD COMMENCING ON THE DATE OF DELIVERY TO THE CUSTOMER. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
- C. THE EXCLUSIONS AND LIMITATIONS SET OUT ABOVE ARE NOT INTENDED TO, AND SHOULD NOT BE CONSTRUED SO AS TO, CONTRAVENE ANY MANDATORY PROVISIONS OF APPLICABLE LAW. IF ANY PART OR TERM OF THIS LIMITED WARRANTY IS HELD TO BE ILLEGAL, UNENFORCEABLE, OR IN CONFLICT WITH APPLICABLE LAW BY A COURT OF COMPETENT JURISDICTION, THE VALIDITY OF THE REMAINING PORTIONS OF THE LIMITED WARRANTY SHALL NOT BE AFFECTED, AND ALL RIGHTS AND OBLIGATIONS SHALL BE ENFORCED AS IF THIS LIMITED WARRANTY DID NOT CONTAIN THE PARTICULAR PART OR TERM HELD TO BE INVALID.
- D. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC, INC. TO ANY REPRESENTATION, CONDITION, OR WARRANTY WITH RESPECT TO THE GENERATOR THAT DEVIATES IN ANY RESPECT FROM THIS LIMITED WARRANTY.
- E. THIS LIMITED WARRANTY IS NOT APPLICABLE TO CATHETERS AND ACCESSORIES USED WITH THE GENERATOR.

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2018/APR/24 at 4:57 p.m. Doc number: M979336A006 Rev. 1A [Global]
Printspec full-b - 8.5 x 11 inches

global:20150610; alldocMod:20170518; alldocStyle:20150728; master_pageset:20150728; languageLookup:20161110; table:20140620; TOC:20160225; characterJoin:20170404; 3of9:20110211;
3of9_svg:20110211

2018/APR/24 at 4:57 p.m. Doc number: M87B336A006 Rev. 1A [Global]
Paper size: 11.6 x 11 inches

global: 20150810, altDocMod: 20170818, altDocStyle: 20160726, master_pageset: 20160726, languageLookup: 20181110, table: 20140820, TOC: 20160226, characterJoin: 20170404, 3d/9: 20110211, 3d/9_svg: 20110211

2018/APR/24 at 4:57 p.m. Doc number: M979336A006 Rev. 1A [Global]
Printspec full-b - 11.5 x 11 inches

global.20150810, alldocMed:20170618, alldocStyle:20150728, master_pageant:20150728, language1.com:up:20181110, title:20140820, TOC:20180225, characterJoin:20170404, Join:20180211, Join_svg:20110211

2016/APR/24 at 4:57 p.m. Doc number M979336A006 Rev 1A [Global]
Printspec: 9x14in - 8.5 x 11 inches

global:20160810, sldocMod:20170518, sldocStyle:20180728, master_pageSet:20180728, languageLookup:20161110, table:20140820, TOC:20160226, characterSet:20170404, sub:20110211
3c@_rvg:20110211

2018/APR/24 at 4:57 p.m. Doc number: M879336A008 Rev 1A [Global]
Printspec full-b - 8.5 x 11 inches

global:20150610; attachMod:20170618; attachStyle:20150728; master_pageset:20150728; languageLookup:20151110; table:20140820; TOC:20150225; characterJen:20170404; JoPP:20110211; Soft_avg:20110211

DocId:20150810, attachment:20170518, address:20150728, mail: jigsaw.20150728, language:colony 20181110, web: 20140820 TC:20140228, email:20170404, info:20110211
info:20110211

2016/APR/24 at 4:57 p.m. Doc Number 44879338A006 Ray 1A (Global)
Printed: 48x48 - 5 x 1 inches

global20150610, adddoc20170618, adddoc20170728, master_jmgmt20150728, mngmntal_backup_20181110, weba_20160820, TOC_20190225, chmascwjohn_20170404, 3082010211, 308_vw_2010211

2016/APR/24 at 4:57 pm Doc number: M979335A006 Rev. 1A [Global]
Printed: Sub-B - 8 5 x 11 inches

Global 20150810, addoklad 20170518, addokladya 20150728, mesnyy_pogovorok 20150728, myagkoye_kolobay 20151110, kabin 20150520 TOC 20150325, chupchakdun 20170504, sub 20150215
2015_gvg 20150215

2015/APR24 at 4:57 p.m. Doc number M979336A006 Rev 1A [Global]
Printed: 5x6 - 5 x 11 inches

2018/APR/24 at 4:57 p.m. Doc number: MB79336A006 Rev 1A [Global]
Printspec full-o - 8.5 x 11 inches

global.20150610; sldocMod:20170518; sldocStyle:20150728; master_pageSet:20150728; languageLookup:20161110; table:20140820; TOC:20160225; characterJoin:20170404; 3of9:20110211;
3of9_avg:20110211

2018/APR/24 at 4:57 p.m. Doc number M979338A006 Rev 1A (Global)
Printspec full-b - 8.5 x 11 inches

global:20190910, altDocMod:20170510, altDocStyle:20150728, master_pageset:20160728, languageLookup:20181110, table:20140620, TDC:20160225, characterJam:20170404, soft:20110211
soft_svg:20110211

2018/APR/24 at 4:57 p.m. Doc number: M979336A006 Rev. 1A [Global]
Paper size: A4-B - 8.5 x 11 inches

global.20150610, eldocMod.20170618, eldocStyle:20150728, master_pageset:20150728, language:lookup:20151110, table:20140820, TOC:20160228, characterJoin:20170404, Soft:20110211,
Soft_sug:20110211

bioRxiv preprint doi: <https://doi.org/10.1101/102111>; this version posted October 10, 2016. The copyright holder for this preprint (which was not certified by peer review) is the author/funder, who has granted bioRxiv a license to display the preprint in perpetuity. It is made available under aCC-BY-NC-ND 4.0 International license.

2018/APR/24 04:57 p.m. Doc Number M979338A006 Rev. 1A [Global]
Printed on: 4/24/2018 4:57 PM

global 20150610 addokMod 20170618, addokStyle 20150728, master_pageout 20180728, language 20181110, nls 20140620, TOC 20150225, characterSet 20170404, job 201 10211, job_wg 20110211

2018APV24 at 4:57 p.m. Doc number: M979336A008 Rev 1A (Global)
Printed: 11/10/2018 8:57 p.m.

:

2018/APR/24 at 4:57 p.m. Doc number: M979336A006 Rev. 1A [Global]
Printspec full-b - 8.5 x 11 inches

global:20150610; alldoc/mtd:20170518; alldoc/style:20150728; master_page/mtd:20150728; language/lookup:20161110; table:20140520; TOC:20160225; character/rev:20170404; spm:20110211; 3rd_eng:20110211

2018/APR/24 at 4:57 p.m. Doc number: M979336A008 Rev. 1A [Global]
Printspec full-b - 8.5 x 11 inches

global.20190610, afdocMed.20170818, afdocStyle.20190728, master_pageenr.20150728, languageatoolup.20161119, table.20140620, TOC.20160225, characterlen.20170404, toB.20110211, 3x10_sug.20110211

2018/APR/24 at 4:57 p.m. Doc number: M879336A006 Rev: 1A [Global]
Printspec full-s - 8.5 x 11 inches

global:20150610, adddocMod:20170518, adddocStyle:20150728, master_page.mst:20150728, language.mst:20181110, table:20140620, TOC:20180226, characterJoin:20170404, SetB:20110211
SetB_avg:20110211

2018/APR/24 at 4:57 p.m. Doc number: M979336A006 Rev 1A {Global}
Pantapac full-b - 8.5 x 11 inches

global:20180410, allDocMod:20170518, allDocStyle:20180728, master_page:20150738, languageLookup:20161110, table:20140820, TOC:20180226, characterJam:20170404, 3d9:20180211, 3d9_svg:20180211

2018/APR/24 at 4:57 p.m. Doc number: M979336A006 Rev. 1A [Global]
Printspec: full-b - 8.5 x 11 inches

global: 20150810, globalMod: 20170518, globalStyle: 20150728, master_page: 20150728, languageLookup: 20161110, table: 20140620, TOC: 20160225, characterJoin: 20170404, font: 20110211, font_size: 20110211

2018/APR/24 at 4:57 p.m. Doc number: M978336A006 Rev. 1A [Global]
Pentapac full-b - 8.5 x 11 inches

global:20160610, sldocMod:20170818, sldocStyle:20160728, master_pageSet:20160728, languageLookup:20161110, table:20140620 TOC:20160226, characterJoin:20170404, Job:20110211
30/8_avg:20110211

2018/APR/24 at 4:57 p.m. Doc number: MS79336A006 Rev. 1A [Global]
Printed: full-b - 8.5 x 11 inches

global:20150610; aldocMod:20170518; aldocStyle: 20150728; master_page: 20150728; languageLookup:20161110; table:20160620; TOC:20160225; character: 20170404; Set:20110211; Set2: 20110211

2018/APR/24 at 4:57 p.m. Doc number: MS79336A006 Rev. 1A [Global]
Printspec: full-size - 8.5 x 11 inches

global:20150610, altDocMod:20170615, altDocStyle:20160726, master_pageset:20160726, languageLookup:20181110, table:20140620 TOC:20160226, characterJoin:20170404, 3of9:20110211, 3of9_sug:20110211

2018/APR/24 at 4 57 p.m. Doc number: MS79336A006 Rev. 1A [Global]
Printsize A4-B - 8.6 x 11 inches

| global 20150610, aIdocId=20170518, aIdocStyle=20150728, master_page=20150728, language=Latvian, table 20140820, TOC=20180225, characterJoin 20170404, 3of9,20110211,
| 3of9_eng 20110211

2018/APR/24 at 4:57 p.m. Doc number: MS793364006 Rev: 1A [Global]
Pentapac full-b - 8.5 x 11 inches

global:20150610, effdocMod:20170518, nldocStyle:20160726, master_pageset:20160726, language:en, setup:20161110, table:20140620, TOC:20160228, character:Jan:20170404, 3of9:20110211
3of9_arg:20110211

2018:APR/24 at 4:37 p.m. Doc number: M979335A006 Rev. 1A (Global)
Printsize full-b - 8.5 x 11 inches

global:20180510; altDocMod:20170518; altDocStyle:20160726; master_pageest:20160728; language:dotup:20181110; table:20140620; TOC:20160225; characterJoin:20170404; Soft:20110211; Soft_svg:20110211

2018/APR/24 at 4:57 p.m. Doc number M979336A006 Rev 1A [Global]
Printarea full-size - 6.9 x 11 inches

picture:20150610, alldocMod 20170618, alldocStyle 20150728 master_pageest:20180728 languageLookup:20181110 table:20140629, TOC:20160226, characterJoin 20170404, Soft:20110211
Soft_avg 20110211

2018/APR/24 at 4:57 p.m. Doc number: M879336AD06 Rev 1A [Global]
Pentapac full-b - 8.5 x 11 inches

global:20150810, afdocMod:20170518, afdocStyle:20150728, master_pageSet:20150728, language:ockup:20161110, table:20140630, VOC:20160225, character:Jan-20170404, 3cB:20110211, 3cB_avg:20110211

global:20150610, modified:20170518, modifiedBy:20150720, master, page:20150728, image:global:20161110, image:20160820 TOC 20160226, image:20170404, 369/20110211, 369/20110211

2016APR24 at 4:57 p.m. Doc number 16978336A006 Rev 1A [Okobal]
Printed: 369 - 8.5 x 11 inches

2018/APR/24 at 4:57 p.m. Doc number: M979336A006 Rev. 1A [Global]
Pentapac A4-b - 8.5 x 11 inches

global:20160610, atDocMod:20170519, atDocStyle:20160728, master_pagecat:20160728, languageLookup:20161110, table:20140620, TOC:20160225, characterJoin:20170404, SetB:20110211, SetB_avg:20110211

2018/APR/24 at 4 57 p.m. Doc number: M979336A006 Rev. 1A [Global]
Printspec A0-b - 8.6 x 11 inches

global.20150610. allDocMod:20170518. allDocStyle:20150728. master_pageSet:20150728. language:colup.20181118. info:20140820. TOC:20160228. characterJoin:20170404. info:20110211. ,
3of9_vvg.20110211

gethui:2015010, attached: 20170510, addocdy: 20160728 master_jpgset20160728, kimgage:okkup:20161110, kbbk: 20160820, TCC: 20160726, chmckckckck:20170404, kpb:20110211, kpb: 20110211

2018APR/24 at 4:57 p.m. Doc number M879336A006 Rev: 1A [Jobel]
Printed kbbk - 8 x 11 inches

2018/APR/24 at 4:57 p.m. Doc number: M879338A006 Rev. 1A [Global]
Printspec full-b - 8.5 x 11 inches

global:20190810. sldocMed.20170518. sldocStyle:20190728. master_pagewet.20190728. languageLookup:20181110. title:20140820 TOC:20180225. characterJoin:20170404. 3of9:20110211.
3of9_svg 20110211

[github.com:50150610, method:20170518, requestId:20190728, name:pagefor:20160728, language:java, time:20170404, size:28110211, sha512sum:20110211](#)

2018/APP/24 at 4, 57 p.m. Doc number 1807235A006 Rev. 1A [Global]
 Printout File: 0 - 8.5 x 11 inches

2018/APR/24 at 4:57 p.m. Doc number: M879336A006 Rev. 1A (Global)
Pentapac 34-b - 8.8 x 11 inches

global:20150810, alldoc:20170518, alldocStyle:20180728, master_page:20150728, language:lookup:20161110, table:20140820, TOC:20180228, character:join:20170404, Soft:20180211, Soft_avg:20110211

Medtronic

Manufacturer:
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com

Authorized Representative in the European Community

Medtronic Ireland
Parkmore Business Park West
Galway, Ireland
+353-91708000

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
tel: 001-763-514-4000
fax: 001-763-514-4879

To Whom It May Concern


Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, as the legal manufacturer of 'Symplicity G3 Renal Denervation RF Generator' hereby confirm that the attached documentation is a true and original document.

Yours sincerely,

Signed on behalf of Medtronic Inc.
Leisa Martinez
Regulatory Affairs Director

18 Nov 2021
Date:

Saadat M. Ali 18 Nov 2021

As a registered Notary Public for the State of Minnesota,
I hereby certify that this is a true and accurate copy of
the original document. The original is filed at
Medtronic, Inc.

Appendix to the User Manual

Simplicity G3 Renal Denervation RF Generator

PRESCRIPTION OF MEDICAL DEVICE

Designed to deliver low-level radio frequency (RF) energy, via the Simplicity Spyral catheter through the wall of the renal artery to achieve renal denervation

INFORMATION ABOUT STICKER

The final packaging for the Russian market is accompanied with additional sticker.

Devices supplied in packaging marked with the following possible information:

- Name of medical device in Russian language;
- Name of the manufacturer;
- Country of origin;
- Code UPN (unique product number);
- Code CFN (customer facing number);
- Date and number of Registration Certificate;
- Warning and informing symbols and inscriptions (if applicable)
- EAC

Sample of the sticker:

	Name of medical device	
	Manufacturer:	Code UPN
	Country of origin:	Code CFN
	Manufacturer's authorized representative::	
	Registration Certificate № _____ or _____	
	Date of manuf./Lot/UBD:	See User Manual
	See on pack	See packaging symbols

Available upon request

MAIN TECHNICAL DESCRIPTION

RF Output Power	6.5 W maximum per electrode, ± 1.5 W accuracy
Impedance Measurement	175 Ω - 200 Ω , 28% 200 Ω - 1200 Ω , 15%
Fuse	2 Slow-Blow (T3.15A T, 250V), 20X5 mm
Temperature Range	Measurement 37°C - 65°C, $\pm 3^\circ\text{C}$
Input Power	100 to 240 V ~, 50/60 Hz universal power supply, 200 VA input power rating Disconnecting the AC power cord removes AC mains power from the unit.
Dimensions	29.9 cm $\times$ 36.7 cm $\times$ 18.9 cm (11.8 in $\times$ 14.5 in $\times$ 7.5 in) (height $\times$ width $\times$ depth)
Weight	9.07 kg (20 lb)
Front Controls	Touchscreen buttons, RF button
Rear Controls	Power On/Off, Volume adjusts
Display	1 touchscreen display on front panel, 1024 $\times$ 768 XGA monitor displaying graphics, messages, impedance, temperature, and time
Connections	AC power, isolated USB, foot switch, four-channel catheter connector, remote control connector, DVI-D connector, and dispersive electrode connector
Protection	Class I, Defibrillator Proof—Type CF, intermittent operation with a duty cycle of 995 seconds on, 60 seconds off. The Simplicity Spyral™ catheter and Simplicity G3™ generator are not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide. IEC 60601-1 Class I, IEC 60601-2-2 Safety of high-frequency surgical equipment; IEC 60601-1-2 Electromagnetic compatibility (EMC).

CONDITIONS OF TRANSPORTATION AND STORAGE

Temperature: from -35°C to $+57^\circ\text{C}$

Relative Humidity: from 30% to 75%

Pressure: from 59.5 to 106 kPa

PROCEDURE AND CONDITIONS FOR THE DISPOSAL OR DESTRUCTION OF A MEDICAL DEVICE

Do not dispose of the Simplicity G3 generator or its accessories (cables, remote, cart, and foot switch) in the unsorted municipal waste stream. Follow local regulations for proper disposal (hazard class A).

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OPERATING CONDITIONS

Temperature 15°C to 40°C (59°F to 104°F)

Humidity 30% to 70% relative humidity, noncondensing

Pressure 700 to 1060 hPa [~ 0.7 to 1.05 ATM]

Available upon request

CUSTOMER SERVICE

Contact information:

Contact for all questions, including complaints, repair and service of medical devices the Authorized Representative in Russia:

Medtronic, LLC.

123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41


Tel.: (495) 580-73-77, Fax: (495) 580-73-78,

www.medtronic.ru

Signed on behalf of Medtronic Inc.

Name: Leisa Martinez

Job Title: Regulatory Affairs Director

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© 2018 Medtronic. All rights reserved. Medtronic и Medtronic logo являются товарными знаками Medtronic. *** Третьи бренды являются товарными знаками их соответствующих владельцев. Все другие бренды являются товарными знаками их соответствующих владельцев.