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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Jiangxi Yikang Medical Instrument Group Co., Ltd. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Zou Ruirui

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(Date: Feb. 25, 2022)

证明网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

C CP27
(2)
APPROVED

Sales manager
Jiangxi Yikang Medical
Instrument Group Co., Ltd.

Guo Hualiang

江西益康医疗器械集团有限公司
JIANGXI YIKANG MEDICAL INSTRUMENT GROUP CO., LTD.

14.02.2022

USER'S MANUAL

on medical item

STERILE DISPOSABLE URETHRAL NELATON CATHETER, CERTUS

manufactured by

JIANGXI YIKANG MEDICAL INSTRUMENT GROUP CO., LTD.,

No.188 Liduaihua Ave, Jinxian Country, 331725, Nanchang, Jiangxi, China

TABLE OF CONTENTS

Product Description.....	4
1.1 Product Name	4
1.2 Intended Purpose of the Product	4
1.3 General Characteristics of the Product.....	5
1.4 Product Features	5
1.5 Product Sterile disposable urethral Nelaton catheter, versions	6
2 Product Configuration	8
2.1 Design Features of the Product	8
Figure 2.1 – Nelaton catheter with cylindrical distal tip	8
Figure 2.2 – Nelaton catheter with Tiemann distal tip.....	8
Figure 2.3 – Section of the Nelaton catheter Body	9
Table 2.1 – Catheter Components according to Figures 2.1-2.3.....	9
Table 2.2 – Color Marking of the Product depending on the Catheter Size	9
Table 2.3 – Nelaton catheter Physical Dimensions (according to Figures 2.1-2.3).....	10
2.2 Description of the Product Components and Materials Used	11
2.2.1 Drainage Connector.....	11
Figure 2.4 – Drainage Connector	11
Table 2.4 – Proportions of Dye Mixing to Produce Finite Colors of Product Connectors	11
2.2.2 Catheter Body.....	12
2.2.3 Drainage Holes.....	12
Figure 2.5 – Distal tip type – cylindrical tip	13
Figure 2.6 – Distal tip type – Tiemann tip	13
Table 2.5 – Dimensions of drainage holes, depending on catheter size.....	13
Figure 2.7 – Dimension of drain holes (values are given in Table 2.5).....	13
2.2.4 Catheter Distal Tip	14
3 General Indications, Contra-Indications and Adverse Effects	15
3.1 Indications:	15
3.2 Contra-Indications:	15
3.3 Potential Complications:	15
3.3.1 False Passage.....	15
3.3.2 Urinary Tract Infections (UTI).....	15
3.3.3 Catheter Obstruction	16
3.3.4 Urethrostenosis.....	16
4 Safety Requirements	17
5 Application Notes.....	18
6 Packaging	19
Table 6.1 – Properties of materials of individual product packaging	19

Marking	21
7.1 Individual Package Marking:	21
Figure 7.1 – Sample of Product Individual Package Marking	21
7.2 Consumer Package Marking	21
Figure 7.2 – Sample of Product Consumer Package Marking	22
7.3 Transport Package Marking	23
7.4 Explanation of Symbols used in Product Marking.....	23
8 Product Maintenance Requirements.....	25
9 Manufacturer Warranty, Storage and Transportation Conditions.....	26
10 Disposal	27

1 PRODUCT DESCRIPTION

1.1 Product Name

Sterile disposable urethral Nelaton catheter, CERTUS (hereinafter referred to as «Product», «Medical Product», «Catheter»)

Manufacturer: Jiangxi Yikang Medical Instrument Group Co., Ltd

Place of Manufacture: No.188 Liduohua Ave, Jinxian Country, 331725, Nanchang, Jiangxi, China

The Product is designed, manufactured according to the Quality Management System and complies with the procedure relating to the EC Declaration of Conformity, in accordance with Appendix II of Directive 93/42/EEC. This Declaration covers all medical devices and is valid for all CE-marked products. The requirements of EN ISO 13485 standards are met and confirmed by ISO Certificate. The application of risk management for this medical device is made in accordance with ISO 14971

Product code according to the All-Russian Classifier of Products by Type of Economic Activity OKPD2: 32.50.13.110.

On issues concerning purchase and operation of the Product, as well as other clarifying issues, please contact the authorized representative of the manufacturer in the Russian Federation – MBT LLC (Medical Business Technologies Limited Liability Company), registered at the address: 25/2, apt.81, Tikhoretsky prospekt, St. Petersburg, Russia. 195427. Tel 8 (911) 928-53-51.

1.2 Intended Purpose of the Product

Sterile disposable urethral Nelaton catheter, CERTUS is intended for intermittent (periodic) bladder catheterization in case of impossibility of natural urination.

Product class, depending on the degree of potential risk of medical use in accordance with the nomenclature classifier of medical devices is 2a. Type of human body contact – short-term use (no more than 60 minutes).

Product application conditions are the use in medical and medical and preventive treatment institutions by qualified and trained personnel specializing in urology.

Product classification according to the Nomenclature Classifier of Medical Devices by Type – 209970.

1.3 General Characteristics of the Product

Sterile disposable urethral Nelaton catheter, CERTUS is a polymeric (PVC material or polyurethane) catheter with a rounded tip and draining side holes at the distal tip. The distal tip has the form of two types: a straight cylindrical (cylindrical distal tip) and a curved, called the Tiemann distal tip.

The size range of the Product is expressed in units (Ch/Fr) according to the Charriere French Gauge and in millimeters (mm): 6 FR (2.0 mm), 8 FR (2.7 mm), 10 FR (3.3 mm), 12 FR (4.0 mm), 14 FR (4.7 mm), 16 FR (5.3 mm), 18 FR (6.0 mm), 20 FR (6.7 mm).

Sterile disposable urethral Nelaton catheter, CERTUS is a single use product and is supplied sterile.

Product sterilization is carried out in accordance with the requirements of EN ISO 11135-1. The sterilization method is gas. The sterilizing agent used is ethylene oxide, the gas mixture is 30 % ETO. Maximum gas residue: 10 ppm (10 parts per million)

Shelf life (sterility) of the product in individual package: 5 years.

Sterile disposable urethral Nelaton catheter, CERTUS does not contain materials of animal and/or human origin.

1.4 Product Features

The Product is produced in the following versions:

- Sterile disposable urethral Nelaton catheter, PVC, cylindrical distal tip.
- Sterile disposable urethral Nelaton catheter, PVC, Tiemann distal tip.
- Sterile disposable urethral Nelaton catheter, PVC, hydrophilic coated, cylindrical distal tip.
- Sterile disposable urethral Nelaton catheter, PVC, hydrophilic coated, Tiemann distal tip.
- Sterile disposable urethral Nelaton catheter, Polyurethane, cylindrical distal tip.
- Sterile disposable urethral Nelaton catheter, Polyurethane, Tiemann distal tip.

Each version of Product is represented by its size range, expressed in units (Ch/Fr) according to the Charriere French Gauge and in millimeters, indicating the length of the catheter (mm), indicating the length of the catheter (mm).

The Product size range has a color marking of the drainage connector (item 1 in Figures 2.1 and 2.2), in accordance with Table 2.2.

1.5 Product Sterile disposable urethral Nelaton catheter, versions

- Sterile disposable urethral Nelaton catheter, PVC, cylindrical distal tip

Typical sizes:

- 6 FR (2.0 mm), length 200 mm;
- 8 FR (2.7 mm), length 200 mm;
- 10 FR (3.3 mm), length 200 mm;
- 12 FR (4.0 mm), length 200 mm;
- 14 FR (4.7 mm), length 200 mm;
- 16 FR (5.3 mm), length 200 mm;
- 18 FR (6.0 mm), length 200 mm;
- 20 FR (6.7 mm), length 200 mm;
- 6 FR (2.0 mm), length 400 mm;
- 8 FR (2.7 mm), length 400 mm;
- 10 FR (3.3 mm), length 400 mm;
- 12 FR (4.0 mm), length 400 mm;
- 14 FR (4.7 mm), length 400 mm;
- 16 FR (5.3 mm), length 400 mm;
- 18 FR (6.0 mm), length 400 mm;
- 20 FR (6.7 mm), length 400 mm.

- Sterile disposable urethral Nelaton catheter, PVC, Tiemann distal tip

Typical sizes:

- 10 FR (3.3 mm), length 400 mm;
- 12 FR (4.0 mm), length 400 mm;
- 14 FR (4.7 mm), length 400 mm;
- 16 FR (5.3 mm), length 400 mm;
- 18 FR (6.0 mm), length 400 mm;
- 20 FR (6.7 mm), length 400 mm.

- Sterile disposable urethral Nelaton catheter, PVC, hydrophilic coated, cylindrical distal tip

Typical sizes:

- 6 FR (2.0 mm), length 200 mm;
- 8 FR (2.7 mm), length 200 mm;
- 10 FR (3.3 mm), length 200 mm;
- 12 FR (4.0 mm), length 200 mm;
- 14 FR (4.7 mm), length 200 mm;
- 16 FR (5.3 mm), length 200 mm;
- 18 FR (6.0 mm), length 200 mm;
- 20 FR (6.7 mm), length 200 mm;
- 6 FR (2.0 mm), length 400 mm;
- 8 FR (2.7 mm), length 400 mm;
- 10 FR (3.3 mm), length 400 mm;

- 12 FR (4.0 mm), length 400 mm;
- 14 FR (4.7 mm), length 400 mm;
- 16 FR (5.3 mm), length 400 mm;
- 18 FR (6.0 mm), length 400 mm;
- 20 FR (6.7 mm), length 400 mm.

- Sterile disposable urethral Nelaton catheter, PVC, hydrophilic coated, Tiemann distal tip

Typical sizes:

- 10 FR (3.3 mm), length 400 mm;
- 12 FR (4.0 mm), length 400 mm;
- 14 FR (4.7 mm), length 400 mm;
- 16 FR (5.3 mm), length 400 mm;
- 18 FR (6.0 mm), length 400 mm;
- 20 FR (6.7 mm), length 400 mm.

- Sterile disposable urethral Nelaton catheter, Polyurethane, cylindrical distal tip

Typical sizes:

- 6 FR (2.0 mm), length 200 mm;
- 8 FR (2.7 mm), length 200 mm;
- 10 FR (3.3 mm), length 200 mm;
- 12 FR (4.0 mm), length 200 mm;
- 14 FR (4.7 mm), length 200 mm;
- 16 FR (5.3 mm), length 200 mm;
- 18 FR (6.0 mm), length 200 mm;
- 20 FR (6.7 mm), length 200 mm;
- 6 FR (2.0 mm), length 400 mm;
- 8 FR (2.7 mm), length 400 mm;
- 10 FR (3.3 mm), length 400 mm;
- 12 FR (4.0 mm), length 400 mm;
- 14 FR (4.7 mm), length 400 mm;
- 16 FR (5.3 mm), length 400 mm;
- 18 FR (6.0 mm), length 400 mm;
- 20 FR (6.7 mm), length 400 mm.

- Sterile disposable urethral Nelaton catheter, Polyurethane, Tiemann distal tip

Typical sizes:

- 10 FR (3.3 mm), length 400 mm;
- 12 FR (4.0 mm), length 400 mm;
- 14 FR (4.7 mm), length 400 mm;
- 16 FR (5.3 mm), length 400 mm;
- 18 FR (6.0 mm), length 400 mm;
- 20 FR (6.7 mm), length 400 mm.

2 PRODUCT CONFIGURATION

2.1 Design Features of the Product

The Product is produced in six versions (see item 1.5).

Figures 2.1-2.3 show dimensions and design features of the Product versions, regardless of the material and the presence of a hydrophilic coating: with a cylindrical distal tip – in Figure 2.1, with a Tiemann distal tip – in Figure 2.2.

Table 2.1 shows the names of the components, and Table 2.2 shows the physical values of the size range of the considered catheters.


Figure 2.1 – Nelaton catheter with cylindrical distal tip
(items are specified in Table 2.1, A-A Section – in Figure 2.3)


Figure 2.2 – Nelaton catheter with Tiemann distal tip
(items are specified in Table 2.1, A-A Section – in Figure 2.3)

A-A (enlarged image)


Figure 2.3 – Section of the Nelaton catheter Body

Table 2.1 – Catheter Components according to Figures 2.1-2.3

Item No.	Name
1	Drainage Connector
2	Catheter Body
3	Drainage Hole
4	Distal Tip

Table 2.2 – Color Marking of the Product depending on the Catheter Size

Catheter size in French Units, Ch/Fr (in millimeters)	Color
6 (2.0 mm)	light green
8 (2.7 mm)	blue
10 (3.3 mm)	black
12 (4.0 mm)	white
14 (4.7 mm)	green
16 (5.3 mm)	orange
18 (6.0 mm)	red
20 (6.7 mm)	yellow

Table 2.3 – Nelaton catheter Physical Dimensions (according to Figures 2.1-2.3)

Distal Tip Type	Length L, (± 20) mm	Diameter O.D, mm (± 10 %)	Size (Fr/Ch)	Flow rate (fluid discharge) through the drainage port, mL/min, not less
cylindrical	200	2.0	6	10
cylindrical	200	2.7	8	15
cylindrical	200	3.3	10	30
cylindrical	200	4.0	12	50
cylindrical	200	4.7	14	70
cylindrical	200	5.3	16	100
cylindrical	200	6.0	18	100
cylindrical	200	6.7	20	100
cylindrical	400	2.0	6	10
cylindrical	400	2.7	8	15
cylindrical	400	3.3	10	30
cylindrical	400	4.0	12	50
cylindrical	400	4.7	14	70
cylindrical	400	5.3	16	100
cylindrical	400	6.0	18	100
cylindrical	400	6.7	20	100
Tiemann	400	3.3	10	30
Tiemann	400	4.0	12	50
Tiemann	400	4.7	14	70
Tiemann	400	5.3	16	100
Tiemann	400	6.0	18	100
Tiemann	400	6.7	20	100

Catheter body force at break:

not less than 15 N

Product weight without package:

no more than 20 g

2.2 Description of the Product Components and Materials Used

Following are descriptions of the catheter components according to Table 2.1 and Figures 2.1-2.3 of this document. The color marking of the catheter type color coding connector shown in Table 2.2 is obtained by mixing various dyes (in Table 2.4. proportions are presented).

2.2.1 Drainage Connector

Designed to connect the article to the urinary manifold and has the shape of a funnel (Fig.2.4). The total length of the drain connector L2 not less than 30 mm.


Figure 2.4 – Drainage Connector

Table 2.4 – Proportions of Dye Mixing to Produce Finite Colors of Product Connectors

Catheter size in French Units, Ch/Fr (in millimeters)	Color	Dye fraction (%)			
		White (Titanium Oxide)	Yellow (Pigment Yellow 151)	Red (Pigment Red 23)	Blue (Pigment Blue 29)
6 (2.0 mm)	light Green	20 %	50%		30 %
8 (2.7 mm)	Blue	70 %			30 %
10 (3.3 mm)	Black	5%	25 %	30 %	40 %
12 (4.0 mm)	White	100 %			
14 (4.7 mm)	Green	10%	40 %		50 %
16 (5.3 mm)	Orange	10 %	70 %	20 %	
18 (6.0 mm)	Red	10%		90%	
20 (6.7 mm)	Yellow	10 %	90 %		

It is made of PVC Compound, Type MT-3075D B1020 manufactured by Shenzhen Hopefinder Polymer Sci.&Tech. Co., Ltd. (China), components: polyvinyl chloride (PVC) CAS No.9002-86-2 85-35% ; di-(2-ethylhexylphthalate) CAS No.117-81-7 15-65%

The drainage connector is color-coded depending on the size (see Table 2.2). The following materials (dyes) are used to obtain the final colors of the connectors:

- CAS No.1317-80-2 – White (*Titanium Oxide*) manufactured by Clearsynth Labs Pvt. Ltd. (India);
- CAS No.31837-42-0 – Yellow (*Pigment Yellow 151*) manufactured by Clearsynth Labs Pvt. Ltd. (India);
- CAS No.6471-49-4 – Red (*Pigment Red 23*) manufactured by Clearsynth Labs Pvt. Ltd. (India);
- CAS No.57455-37-5 – Blue (*Pigment Blue 29*) manufactured by Clearsynth Labs Pvt. Ltd. (India).

2.2.2 Catheter Body

It is intended for urine diversion from the bladder into the urine-collecting bag.

The catheter body, depending on the version, is made of PVC or polyurethane. On the surface of the catheter body there are no foreign inclusions, shells, dents and other mechanical damage.

Hydrophilic coating of catheters made of PVC provide an even easier penetration of the catheter into the urethra, because when the surface of the catheter is wetted with water, the hydrophilic coating makes the surface of the catheter more slippery and smooth. No hydrophilic coating is used for polyurethane. Thickness of hydrophilic coating layer is not more than 0.5 mm.

Materials used:

- For the PVC Product: PVC Compound, Type MT-3075D B1020 manufactured by Shen-zhen Hopefinder Polymer Sci.&Tech. Co., Ltd. (China), components: poly-vinyl chloride (PVC) CAS No.9002-86-2 85-35% ; di-(2-ethylhexylphthalate) CAS № 117-81-7 15-65%, is used.
- For the hydrophilic coated PVC Product: PVC Compound, Type MT-3075D B1020 manufactured by Shen-zhen Hopefinder Polymer Sci.&Tech. Co., Ltd. (China), components: poly-vinyl chloride (PVC) CAS No.9002-86-2 85-35% ; di-(2-ethylhexylphthalate) CAS № 117-81-7 15-65%, is used and the HYDROMER Coating Solution No.2204-199 hydrophilic coating manufactured by Hydromer, Inc. (USA) are used in the following composition: methyl ethyl ketone (2-Butanone)CAS No.78-93-3, tetrahydrofuran CAS No.109-99-9, ethyl lactate CAS No.97-64-3.
- For the Polyurethane Product, the ESTANE®58315 NAT035 polyurethane manufactured by The Lubrizol Corporation (USA) is used.

2.2.3 Drainage Holes

Designed to divert urine through the catheter out from the bladder. Drainage holes are located on the catheter body in front of the distal tip (see Figures 2.5, 2.6).

NOTE – The location of the drainage holes in the figures is conditional.


Figure 2.5 – Distal tip type – cylindrical tip


Figure 2.6 – Distal tip type – Tiemann tip

Dimensions of drain holes are given in Table 2.5 as per Figure 2.7.

Table 2.5 – Dimensions of drainage holes, depending on catheter size

Foley catheter size according to French Gauge, Ch/Fr	Drain hole size mm, not less than	
	Width W1	Length L1
6	0,5	1,0
8		
10		
12	1,0	2,0
14		
16	2,0	4,0
18		
20		


Figure 2.7 – Dimension of drain holes (values are given in Table 2.5)

2.2.4 Catheter Distal Tip

Designed to facilitate the passage of the catheter inside the urethra and is made of the same materials as the catheter body (see item 2.2.2).

Characteristics of the distal tip: closed atraumatic, two types available (see Figures 2.5, 2.6):

Tiemann tip catheter is recommended for catheterization of the bladder in patients with urinary tract obstructive disorder.

Materials used:

For the following Product versions:

- Sterile disposable urethral Nelaton catheter, PVC, cylindrical distal tip.
- Sterile disposable urethral Nelaton catheter, PVC, Tiemann distal tip.

The PVC Compound, Type MT-3075D B1020 manufactured by Shen-zhen Hopefinder Polymer Sci.&Tech. Co., Ltd. (China), components: poly-vinyl chloride (PVC) CAS No.9002-86-2 85-35% ; di-(2-ethylhexylphthalate) CAS № 117-81-7 15-65% is used;

For the following Product versions:

- Sterile disposable urethral Nelaton catheter, PVC, hydrophilic coated, cylindrical distal tip.
- Sterile disposable urethral Nelaton catheter, PVC, hydrophilic coated, Tiemann distal tip.

The PVC Compound, Type MT-3075D B1020 manufactured by Shen-zhen Hopefinder Polymer Sci.&Tech. Co., Ltd. (China), components: poly-vinyl chloride (PVC) CAS No.9002-86-2 85-35% ; di-(2-ethylhexylphthalate) CAS № 117-81-7 15-65% is used; HYDROMER Coating Solution No.2204-199 hydrophilic coating manufactured by Hydromer, Inc. (USA) in the following composition: methyl ethyl ketone (2-Butanone) CAS No.78-93-3, tetrahydrofuran CAS No.109-99-9, ethyl lactate CAS No.97-64-3.

For the following Product versions:

- Sterile disposable urethral Nelaton catheter, Polyurethane, cylindrical distal tip.
- Sterile disposable urethral Nelaton catheter, Tiemann distal tip.

ESTANE® 58315 NAT035 polyurethane manufactured by The Lubrizol Corporation (USA) is used;

3 GENERAL INDICATIONS, CONTRA-INDICATIONS AND ADVERSE EFFECTS

3.1 Indications:

- Urinary incontinence.
- Urinary retention.
- Measurement of residual volume and urine removal monitoring.
- Urethra edema due to inflammatory processes.
- Neoplasms of the prostate, bladder and (or) ureter.

NOTE – Instructions for standard use of the described article are given in item 5, taking into account paragraphs 3 and 4. of this document. The decision to apply in each case should be made directly at the clinic by the doctor.

3.2 Contra-Indications:

- Inflammation of the prostate gland.
- Pelvic organs/urinary tract injuries.
- Severe urethritis.
- Severe urethral stricture, urinary tract blunt or acute injuries.

3.3 Potential Complications:

3.3.1 False Passage

Etiology. The catheter is too small or too rigid; unjustified use of force; urinary tract stricture.

Prevention. It is necessary to use a catheter of the appropriate size and apply only a slight force to advance the catheter.

Actions taken. You should not try to advance a catheter that does not move using light effort; if there is a suspected false passage (acute pain, urethremorrhagia), you must immediately remove the instrument and not resume attempts to introduce it within 5 days; urological consultations.

3.3.2 Urinary Tract Infections (UTI)

Etiology. Long-term use of the catheter (catheters are foreign bodies that eventually become colonized).

Prevention. Minimization of catheterization duration.

Actions taken. Systemic antibiotics (oral or intravenous administration). Catheter removal as soon as possible.

3.3.3 Catheter Obstruction

Etiology. Blocking the catheter lumen with blood or heterogeneous particles of another substance.

Prevention. Increased oral intake of liquids; Methenaminum; urinary acidification with vitamin C.

Actions taken. Catheter removal.

3.3.4 Urethrostenosis

Etiology. Injury during catheterization. Urethritis.

Prevention. Careful advancement of the catheter.

Actions taken. Urological consultation regarding the urinary tract dilatation if necessary.

4 SAFETY REQUIREMENTS

Precautions during using the Product:

- The Product is intended only for operational use by medical personnel.
- The product is for single use only. Do not reuse the Product.
- Do not use the expired Product or Product in damaged individual packaging. Please read the instructions on application carefully before using the Product.
- The Product is suitable only within the time limits indicated on individual and consumer packaging and must not be re-sterilized.
- You must select the appropriate catheter size; a small size should be chosen for children or patients with severe urethral stricture.
- To avoid damage to the urethral mucosa, insertion of the catheter should be careful. It is not allowed to pull and push the catheter for several times.
- If the bladder is too full and the patient cannot empty it on his own, drainage should be slow and the first drainage volume should not exceed 1000 ml, in order to avoid bleeding or fainting due to rapid decompression of the bladder.
- Before measuring the urine residual volume, you should ask the patient to urinate independently, and then perform a catheterization. The normal urine residual volume should be 5-10 ml, if it exceeds 100 ml, the use of permanent catheters is required.
- Always check the fixation during catheterization; if the catheter slips out, the catheter should be replaced in accordance with the current guidelines in clinical settings. It is intended to take a break for several hours before reintroducing the catheter into the urethral orifice.
- Do not use force to remove the catheter in the presence of obstruction, because this may injure the urethra.
- After use, the Product and packaging materials must be disposed in accordance with the current medical standards of the hospital, or in accordance with environmental protection requirements.
- If an employee, while handling the Product, receives an injury that is potentially dangerous in terms of infection (damage to skin and/or mucosa integrity), emergency preventive measures must be taken. The workplace should have a first aid kit for injuries.

During operation, the product is resistant to temperature influences from 10 °C to 42 °C and relative humidity of 100 %.

5 APPLICATION NOTES

- Explain to the patient the essence and course of the forthcoming procedure and get his/her consent or consent of his/her family members; mention the discomfort during the procedure and how to reduce it, possible complications and care after catheterization.
- Prepare the required catheter by type and size.
- Help the patient to take the necessary position: lying on his/her back with legs apart, bent at the knee joints; place an absorbent diaper (or drawsheet and bed mat) under the patient's pelvis.
- Prepare a container for the used material near the patient and carry out hygienic treatment of the external genital organs, disinfect the site of the Product introduction.
- Open the bottle with sterile lubricant.
- Check the integrity of the individual packaging; check the packaging for damage. Open the catheter pack.
- Lubricate the catheter with sterile grease or moisten the catheter with water if the catheter is hydrophilic coated.
- Carefully insert the catheter into the urethra and gently, intercepting it, push the catheter deeper into the urethra until urine appears (direct urine into the tray).
- Advance the catheter 3-6 cm further.
- Connect the catheter through the drain connector to the urine-collecting bag.
- Empty the bladder.
- Carefully remove the catheter.
- Remove the waterproof diaper (drawsheet and bed mat) and help the patient lie comfortably.
- Place used material in a container with a disinfectant solution in preparation for disposal.

IMPORTANT! Perform catheterization in clean sterile medical gloves.

6 PACKAGING

The following types of packaging are used for the Product: individual, consumer and transport packaging.

Individual packaging provides protection against mechanical damage to the Product itself, and also maintains its sterility throughout its whole shelf life.

Individual packaging of the Product consists of A-2 medical paper manufactured by Nanchang Meiguang Packing Co., LTD (China) or medical paper manufactured by Quzhou Z-SHINE Medical Packing Co., Ltd (China); PE/PET polymer film – composition: polyethylene CAS No.9002-88-4, polyethylene terephthalate CAS No.25038-59-9, manufactured by Quzhou Z-SHINE Medical Packing Co., Ltd (China) or film XPP ("XPP medical package films") manufactured by "Fujian green sail advanced medical material Co., Ltd.", (China).

Table 6.1 – Properties of materials of individual product packaging

Packagin materil	Properties	Value
Paper +PE/PET film:	- seam strength is not less than	1.2 N/15 mm
Paper:	- paper thickness	0.08 (± 0.02) mm
	- density of paper	60 g/m ² ± 20 %
	- forcing strength, more than	70 kPa
	- tensile strength, more than	0,5 kN/m
	- elongation, not less than	0.7 %
	- air permeability, more than	3.4 um/pa.s
PE/PET film:	- film thickness	0.06 (± 0.02) mm
	- film density	55 g/m ² (± 20 %)
	- tensile strength in longitudinal/transverse direction, not less than	30 MPa
	- elongation, not less than	30 %

Dimensions of an individual package depends on the size of the Product enclosed therein:

- Length 200 mm – individual package dimensions: (246 ± 20) mm x (33 ± 5) mm;
- Length 400 mm – individual package dimensions: (470 ± 20) mm x (33 ± 5) mm.

If necessary, the inner packaging materials may be replaced with materials of other manufacturers, while retaining the characteristics and verified requirements for medical packaging materials specified in Table 6.1.

Individual packaging is marked in accordance with item 7.1 of this document.

Consumer packaging of the Product is made of cardboard and is intended for packaging and storage of 100 pieces of catheters of a certain version and size.

Dimensions of consumer package depends on the size of the individually packaged Products enclosed therein. If the length of the catheter is:

- 200 mm – consumer package dimensions (± 20 mm): 280 mm x 175 mm x 120 mm;
- 400 mm – consumer package dimensions (± 20 mm): 500 mm x 175 mm x 120 mm.

Consumer packaging is marked in accordance with item 7.2 of this document.

Consumer packaging is covered with polyethylene tape with adhesive layer.

Transport packaging of the Product is made of cardboard and is intended for transportation and storage of consumer packaging of Nelaton Catheters packaged therein: from 5 consumer packages (500 pieces of catheters) and more.

Instructions on the transport packaging marking are given in item 7.3 of this document.

7 MARKING

7.1 Individual Package Marking:

- Product name.
- Product size.
- Lot (batch) number.
- Date of manufacture.
- Information on sterility and sterilization method.
- Sterility expiry date.
- Name, address, of the manufacturer, country of goods origin.
- Trademark or other additional information that attracts consumer attention, including advertising.
- Name, address, telephone number of the authorized representative of the manufacturer in the territory of the Russian Federation.
- Marking (indication) of opening location
- Inscription: «Non-toxic».
- Marking with the symbol «Do not re-sterilize».
- Marking with the symbol «For single use only, do not reuse».
- Marking with the symbol «Do not use in case of packaging damage».
- Marking with the symbol «See Operating Manual».
- Marking with the Declaration of Conformity symbol (PCTconformity mark).
- Marking with the symbol «Do not dispose as normal waste! Hand over to disposal facilities!».
- Marking with the symbol «Do not expose to sunlight».
- Marking with the symbol «Protect from moisture».
- Number and date of Certificate of Registration.
- Marking (indication) of the opening place.


Figure 7.1 – Sample of Product Individual Package Marking

7.2 Consumer Package Marking:

- Product name.
- Product size.
- Number of Products per package.
- Lot (batch) number.
- Date of manufacture.
- Information on sterility and sterilization method.
- Sterility expiry date.
- Name, address, of the manufacturer, country of goods origin.

- Logo.
- Name, address, telephone number of the authorized representative of the manufacturer in the territory of the Russian Federation.
- Application notes.
 - Safety precautions.
 - Storage conditions.
 - Inscription: «Non-toxic».
 - Marking with the symbol «Do not re-sterilize».
 - Marking with the symbol «For single use only, do not reuse».
 - Marking with the symbol «Do not use in case of packaging damage».
 - Marking with the symbol «See Operating Manual».
 - Marking with the Declaration of Conformity symbol (PCTconformity mark).
 - Marking with the symbol «Do not dispose as normal waste! Hand over to disposal facilities!».
 - Marking with the symbol «Do not expose to sunlight».
 - Marking with the symbol «Protect from moisture».
 - Number and date of Certificate of Registration.


Figure 7.2 – Sample of Product Consumer Package Marking

7.3 Transport Package Marking

- Product name.
- Number and date of Certificate of Registration.
- Lot (batch) number.
- Date of manufacture
- Sterility expiry date.
- Number of Products per package.
- Name, address, of the manufacturer, country of goods origin.
- Name, address, telephone number of the authorized representative of the manufacturer in the territory of the Russian Federation.
- Transportation requirements (in symbols, including marking with the symbol of the temperature range for transportation).
- Weight (kg).
- Package dimensions (cm).
- Marking with the symbol «Do not use in case of packaging damage».

7.4 Explanation of Symbols used in Product Marking

See Operating Manual

Declaration of Conformity symbol (PCT conformity mark).

Do not use in case of packaging damage.


Disposal prohibition: "Do not dispose as normal waste! Pass to disposal facilities!"

Do not exposure to sunlight


Protect from moisture

Temperature conditions (during transportation)

For single use only, do not reuse

ETO sterilized

Use before ____ (expiration date)

Do not re-sterilize

Lot/Batch number

Date of manufacture

Name, address, telephone number of the manufacturer, country of goods origin

WEIGHT: _____ kg Weight**SIZE (cm):** _____ Transport package dimensions

8 PRODUCT MAINTENANCE REQUIREMENTS

The Product is disposable.

There are no maintenance requirements.

When using the Product please pay special attention to the instructions in item 3 and 4 of this document.

9 MANUFACTURER WARRANTY, STORAGE AND TRANSPORTATION CONDITIONS

The warranty shelf life of the Product, the sterilization validity is 5 years.

Keep the Products in dry and ventilated rooms with the humidity not more than 80 %, at the temperature of $(25 \pm 10) ^\circ\text{C}$.

Range of permissible temperatures during transportation is from minus 30 to plus 49 $^\circ\text{C}$.

During transportation and storage, the packaging integrity violation, direct sunlight, moisture and exposure to high temperatures should be avoided. During storage and transportation, Products should be located at a distance of not less than 1 meter away from heat-radiating devices.

10 DISPOSAL

The Product disposal should be carried out in accordance with the guidelines for class B wastes (epidemiologically hazardous wastes).

The Product after use are subject to mandatory purifying (disinfection). Selection of disinfection method is determined by the capabilities of the organization engaged in medical activities, and is carried out when developing a scheme for the treatment of medical waste.

In any method of class B medical waste disinfection, the disinfectants and equipment registered in the Russian Federation shall be used in accordance with the instructions for their use.

The use of recycling technologies, including waste sorting, is possible only after preliminary hardware disinfection of class B wastes by physical methods. It is forbidden to use the secondary raw materials obtained from medical waste for the manufacture of children's goods, for materials and products in contact with drinking water and food products, for medical devices.

Class B waste is collected in disposable soft (bags) or solid (non-punctured) package (containers) of yellow color or having yellow markings. Selection of packaging depends on the morphological composition of the waste. Soft package (disposable bags) used to collect Class B waste should ensure that they can safely collect no more than 10 kg of waste.

After filling the bag by no more than 3/4, the person responsible for the collection of waste in this medical unit fastens the bag or closes it using tie-wraps or other devices that prevent the spillage of Class B waste. Solid (non-punctured) containers are closed with lids. It is not permitted to move Class B waste outside the unit in open containers.

At the final packaging of class B waste in order to remove it from the unit (organization), disposable containers (bags, bins) with class B waste are marked with the inscription: "Waste. Class B" with the name of the organization, unit, date and surname of the person responsible for the waste collection.

Accumulation and temporary storage of non-disinfected Class B waste are carried out separately from waste of other classes in special rooms, which exclude access by unauthorized persons.

In small medical organizations (health centers, offices, medical and obstetric stations, etc.), temporary storage and accumulation of Class B waste in containers located in utility rooms is allowed (refrigeration equipment is used for storage for more than 24 hours). The use of refrigeration equipment for other purposes than for the accumulation of waste is not allowed.

Reusable containers for transportation of Class B waste products are subject to washing and disinfection after each emptying, in a special area for washing, disinfection and desinsection of containers and vehicles. The requirements for personnel involved in the transport of medical waste, as well as the sanitary and epidemiological requirements for vehicles intended for the transport of non-disinfected Class B wastes, must comply with the relevant legislation.

Disposal of neutralized Class B wastes at the landfill is allowed only when their presentation (grinding, sintering, pressing, etc.) is changed and their reuse is impossible.

When collecting medical waste, it is prohibited to:

- destroy manually Class B waste;
- pour out (reload) unpacked Class B waste from one container to another;
- compact Class B waste;
- carry out any operations with waste without gloves or the necessary personal protective equipment and work clothing;
- use soft disposable packaging to collect sharp medical instruments and other sharp objects;
- install disposable and reusable containers for collecting waste at the distance less than 1 m from heating devices.

Industrial control of collection, temporary storage and disposal of medical waste includes:

A) Visual and documentary inspection (at least once a month):

- amount of consumables (stock of packages, containers, etc.), small-scale mechanization means, disinfectants;
- provision of personnel with personal protective equipment, organization of a centralized washing of work clothing and its regular change;
- sanitary condition and disinfection regime for temporary storage facilities and/or medical waste treatment areas, garbage chutes, container sites;
- compliance with disinfection/decontamination regimes;
- regularity of waste removal.

B) Laboratory and instrumental verification:

- microbiological monitoring of the effectiveness of disinfection/neutralization of waste at plants according to approved methods (at least once a year);
- control of microclimate parameters (at least once a year);
- control of the working zone air in the areas of disinfection/neutralization of waste for the content of volatile toxic substances (it should be carried out in accordance with the technological regulations of the equipment).

Emergency prevention measures must be taken if an employee receives an injury that is potentially dangerous in terms of infection (damage to the integrity of the skin and/or mucous membranes) when handling medical waste. The workplace should have a first aid kit for injuries.

Notification, recording and investigation of cases of personnel infection by infectious agents associated with professional activities are carried out in accordance with established requirements.

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