RanD S.p.A.
Via Statale 12, n. 62
Medolla (Modena)
Tel. +39 0535 49283 Fax +39 0535 660636
E-mail: info@rand-biotech.com Web: www.rand-biotech.com

RAND

RanD S.p.A.
Via Statale 12, n. 62 - 41036 MEDOLLA (MO)
Tel. +39 0535 49283 Fax +39 0535 660636
E-mail: info@rand-biotech.com Web: www.rand-biotech.com

INNOVATION IN MEDICAL TECHNOLOGY

INSTRUCTIONS FOR USE

1. Name of the product

HANG&GO HT BASIC

Catalogue number R9900088

2. Manufacturer address, phone, email

RanD S.p.A.
via Statale 12, 62
41036 Medolla (MO)
Italy
phone: +39 0535 442983
fax: +39 0535 660636
email: info@rand-biotech.com
WebSite: www.rand-biotech.com

3. Authorized representative in Russia

Company name:	Biocommerce Ltd
Legal and manufacturing site:	Kulakova street 20, building 1A Moscow 123592 Russia
Telephone/Fax No.	Tel. +7(495) 781-17-87
Email:	info@biocommerce.ru
Website:	www.biocommerce.ru

4. Designation of the set with indication of the potential consumer

Disposable device manufactured/assembled by RAND. Not made with natural rubber latex. Made of non-pyrogenic, non-toxic and biocompatible materials according to ISO 10993 standards and in relation with the intended use.

Specifically designed for use in oncology for intraperitoneal and intrapleural hyperthermic perfusion.

Each "HANG&GO HT BASIC" (code R9900088) contains the following devices:

- 1 "Hang&Go HT Equipment" line set,
- 1 "Hang&Go Table Pack" line set.

The device is intended for professional use. The procedure must be carefully and constantly controlled by qualified personnel trained in use of the device.

5. Risks of the use, possible side effects, their elimination

There is no evidence of side effects specifically related to the use of the product. The device is contraindicated for procedures outside of the intended use.

The device contains phthalates. Although it has been demonstrated on experimental animal models that DEHP and other phthalates contained in PVC may cause hepatic and endocrine alterations, especially in the reproductive system, there is no proof of a similar effect on man. It is however reasonable to exercise due prudence having a physician assess the cost/benefit ratio of a treatment during which there is extended exposure by direct contact of PVC with blood or plasma in subjects considered at risk, such as pregnant or breast-feeding women and children.

PVC medical devices must not be used for infusion of lipid-based drugs (e.g. Paclitaxel) capable of extracting DEHP or otherwise incompatible with phthalates. It is recommended to carefully read the information leaflet on the drug to be infused before using the device in order to exclude any incompatibility.

INNOVATION IN MEDICAL TECHNOLOGY

6. Specifications

Recommended perfusion flow: up to 2000 ml/min

Operating temperature: 37 to 46°C

Overall priming volume of the extracorporeal circuit (including the table pack lines): 500 ml

Overall priming volume: 1000 ml (including a minimum 500 ml level within the filtered reservoir).


7. Information on the main characteristics of the operation, the conditions of transportation and storage of the set

For the main characteristics of the operation, see additional information.

Storage conditions are indicated on label in form of symbols. They are represented here below. All symbols are explained in local language in the instructions for use.

Keep away
from heat

Temperature
limitation

Keep dry

The storage temperature must be within +5°C e +35°C. Keep the product into its original packaging out of light and humidity. Transport and/or storage conditions other than those prescribed may cause damages to the device. Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken.

8. Information on the sterile condition of the product and the method of its sterilization

The product is sold sterile and for single-use/single-patient only. It cannot be reused or resterilized.

The product is sterilized by ethylene oxide, and it has a shelf life of 3 years from the date of sterilization, as indicated on the label.

9. Information of known restriction on the joint use of medical devices

The HANG&GO HT BASIC has been specifically designed for use in combination with the Rand PERFORMER HT.

10. Information on precautions for the use of a medical device

All the precautions are given in the instructions for use provided with each set. Please refer to annex 1 IFU.

11. Manufacturer's warranty

Rand gives a warranty of 3 years from the date of sterilization, which corresponds to the product shelf life.

12. Disposal

As described in the instructions for use, disposal of used device should be in compliance with the applicable regulations in force in the country of use for contaminated products.

INNOVATION IN MEDICAL TECHNOLOGY

ADDITIONAL INFORMATION

13. Functional characteristics

HANG&GO HT BASIC is a disposable set of extracorporeal circulation lines and bags, specifically designed, realized and assembled by Rand to be used in conjunction with Rand PERFORMER HT equipment in HIPEC (hyperthermic intra-peritoneal/intra-pleural perfusion) treatment. The set is composed of components made of plastic polymers. The device is sterilized by ethylene oxide, non-pyrogenic and intended for single-use only, biocompatible for its intended use according to ISO 10993 standards.

The HANG&GO HT BASIC is composed of the following main parts:

Equipment lines set
(reservoir + tubes)



Table pack lines set
(infusion and withdrawal lines)


Heating bag and waste bag

The picture below represents the fluidic scheme of an HIPEC treatment:

Fluidic scheme of Hyperthermic Intraperitoneal Perfusion

NOVATION IN MEDICAL TECHNOLOGY

Fluidic scheme of Hyperhermic Intrapleural Perfusion

14. Information on the main characteristics of the operation

The HIPEC procedure consists of perfusing the peritoneal or pleural cavity with a defined volume of clear solution, also enriched by the physician with chemotherapeutic agents, at a controlled temperature typically ranging from 41 to 43°C, for a duration according to a protocol selected by the physician. The treatment can be carried out on open or closed peritoneal or pleural cavity. Flow typically ranges from 500 to 2000 mL/min. The procedure consists of three phases:

CIRCUIT PRIMING / SOLUTION WARM-UP: before connecting the Patient, the system automatically carries out the HANG&GO HT tubing / reservoir filling with a volume according to the therapeutic scheme in use. By means of the PERFORMER HT PMI pump movement, the fresh physiological solution is taken from the bag, poured in the reservoir up to the desired volume, and recirculated from the reservoir to the heater bag and then back to the reservoir to be warmed up to the desired temperature.

PATIENT FILLING / SOLUTION RECIRCULATION: once the surgeon has placed the inflow and outflow catheters (up to 6) in the peritoneal / pleural cavity via surgical incision in the skin, and the solution has reached the desired temperature, the surgeon connects the Patient's catheters to the HANG&GO HT Table Pack lines. Once the fluid in the Patient's cavity has reached the set volume, the solution is pumped back from the Patient into the reservoir, then again through the heater and to the Patient and so on, keeping a continuous circulation.

RINSING / EMPTYING: at the end of the procedure, the system gives the possibility to rinse the cavity, then the peristaltic pump transfers the exhausted solution to the reservoir and waste bag.

15. List of applicable national standards

General standards

EN ISO 14971:2012	Application of risk management to medical devices (ISO 14971:2007)
CEI EN 60812:2006	Analysis techniques for system reliability - Procedure for failure mode and effects analysis (FMEA)
UNI ISO 2859-1:2007	Sampling procedures for inspection by attributes - Part 1: Sampling schemes indexed by acceptance quality limit (AQL) for lot-by-lot inspection

Packaging and labeling

EN 980:2008	Symbols for use in the labeling of medical devices
-------------	--

INNOVATION IN MEDICAL TECHNOLOGY

EN ISO 15223-1:2016	Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: general requirements (ISO 15223-1:2016)
EN 1041:2008	Information supplied by the manufacturer of medical devices
EN ISO 11607-1:2009	Packaging for terminally sterilized medical devices – Part 1: requirements for materials, sterile barrier system and packing system (ISO 11607-1:2006)
EN ISO 11607-2:2006	Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes (ISO 11607-2:2006)
UNI EN 868-5:2009	Packaging for terminally sterilized medical devices. Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods
ASTM D4169-09	Standard practice for Performance testing of shipping Containers and Systems
ASTM F1980-07	Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices
ASTM F1886-98 (2004)	Standard Test method for determining integrity of seals of medical packaging by visual inspection
ASTM F1929-98 (2012)	Standard Test method for detecting seal leaks in porous medical packaging by dye penetration

Biocompatibility

EN ISO 10993-1:09/AC:10	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process (ISO 10993-1:2009)
EN ISO 10993-4:2009	Biological evaluation of medical devices - Part 4: Selection of tests for interactions with blood (ISO 10993-4:2002, including Amd 1:2006)
EN ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity (ISO 10993-5:2009)
EN ISO 10993-7:08/AC:09	Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals
EN ISO 10993-10:2013	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization (ISO 10993-10:2010)
EN ISO 10993-11:2009	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity (ISO 10993-11:2006)
EN ISO 10993-17:2009	Biological evaluation of medical devices - Part 17: Establishment of allowable limits for leachable substances (ISO 10993-17:2002)
Notice n°36	"Guidelines for Basic Biological tests of Medical Materials and devices" in "Guide to Medical Device Registration in Japan" published by Yakiji Nippo Ltd, 1995.

Sterilization and shelf-life

EN 556-1:2001/AC:06	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
EN ISO 11135-1:2007	Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices (ISO 11135-1:2007)
EN ISO 11138-1:2017	Sterilization Of Health Care Products - Biological Indicators - Part 1: General Requirements
EN ISO 11138-2:2009	Sterilization Of Health Care Products - Biological Indicators - Part 2: Biological Indicators For Ethylene Oxide Sterilization Processes (ISO 11138-2:2006)

INNOVATION IN MEDICAL TECHNOLOGY

EN ISO 11737-1:06/AC:09	Sterilization of medical devices - Microbiological methods - Part 1: Determination of a population of microorganisms on products (ISO 11737-1:2006)
EN ISO 11737-2:2009	Sterilization Of Medical Devices - Microbiological Methods - Part 2 - Tests Of Sterility Performed In The Validation Of A Sterilization Process (ISO 11737-2:2009)
ASTM F1980-07(2011)	Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices
Controlled environment	
EN ISO 14644-1:2015	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness (ISO 14644-1:2015).
EN ISO 14698-1:2003	Cleanrooms and associated controlled environments - Biocontamination control. Part 1: General principles and methods (ISO 14698-1:2003)
EN ISO 14698-2:03/AC:06	Cleanrooms and associated controlled environments - Biocontamination control - Part 2: Evaluation and interpretation of biocontamination data (ISO 14698-2:2003/Cor.1:2004)
Design	
EN 1707:1996	Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment - Lock fittings
EN 62366:2008	Medical devices - Application of usability engineering to medical devices (IEC 62366:2007)

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ENGLISH

DESCRIPTION

Disposable device manufactured assembled by RAND. Not made with natural rubber latex. Made of non-pyrogenic, non-toxic and biocompatible materials according to ISO 10993 standards and in relation with the intended use. Specifically designed for use in oncology for intraperitoneal and intrapleural hyperthermic perfusion.

Each "HANG&GO HT" (code R9900067) kit contains the following devices:

- 1 "Hang&Go HT Equipment" line set (code R9900074).
- 1 "Hang&Go Table Pack" line set (code R9900075).
- 6 round silicone Catheters Ch 24 (code R9900071).
- 3 "Mon-a-therm" temperature probes (code 90050).

Each "HANG&GO HT FR" (code R9900110) kit contains the following devices:

- 1 "Hang&Go HT Equipment" line set (code R9900074).
- 1 "Hang&Go Table Pack" line set (code R9900075).
- 6 round silicone Catheters Ch 24 (code R9900071).
- 3 "Mon-a-therm" temperature probes (code 90050).
- 1 "Flow Reverse" (code R9900101).

Each "HANG&GO HT BASIC FR" (code R9900127) kit contains the following devices:

- 1 "Hang&Go HT Equipment" line set (code R9900074).
- 1 "Hang&Go Table Pack" line set (code R9900075).
- 1 "Mon-a-therm" temperature probe (code 90050).
- 1 "Flow Reverse" (code R9900101).

Each "HANG&GO HT BASIC" (code R9900088) contains the following devices:

- 1 "Hang&Go HT Equipment" line set (code R9900074).
- 1 "Hang&Go Table Pack" line set (code R9900075).

TECHNICAL FEATURES

Recommended perfusion flow: up to 2000 ml/min

Operating temperature: 37 to 45°C

Overall priming volume of the extracorporeal circuit (including the table pack lines): 500 ml

Overall priming volume: 1000 ml (including a minimum 500 ml level within the filtered reservoir).

INTENDED USE

Disposable devices intended to recirculate, filtrate and perfuse clear solutions, also enriched with chemotherapeutic agents, in the thoracic or abdominal cavity. The devices are intended for use up to 6 hours, in regional cancer hyperthermic therapies following surgical reduction. The devices can be used safely in conjunction with the medical devices listed in the paragraph "Medical devices useable in combination".

SIDE EFFECTS AND CONTRAINDICATIONS

There is no evidence of side effects specifically related to the use of the product. The device is contraindicated for procedures outside of the intended use.

SAFETY INFORMATION

The information to draw the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device are indicated in the text as follows:

WARNING

Indicates serious adverse reactions and potential safety hazards for the user and/or the patient that can occur both during proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such case.

PRECAUTION

Indicates any special precaution than the user must exercise for safe and effective use of the device.

Below follows some general safety information for proper use of the device. Furthermore, specific safety-related information is given in the instructions for use where they may condition the operation to be performed.

For further information and/or in case of complaints, contact Rand directly or the authorised local representative.

Hang&Go HT-1

ENGLISH

WARNING

- The user must read and understand all user information prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- The devices are intended for professional use only. The procedure must be carefully and constantly controlled by qualified personnel trained in use of the device.
- The device contains phthalates. Although it has been demonstrated on experimental animal models that DEHP and other phthalates contained in PVC may cause hepatic and endocrine alterations, especially in the reproductive system, there is no proof of a similar effect on man. It is however reasonable to exercise due prudence having a physician assess the cost/benefit ratio of a treatment during which there is extended exposure by direct contact of PVC with blood or plasma in subjects considered at risk, such as pregnant or breast-feeding women and children.
- Visually inspect and carefully check the devices before use. Transport and/or storage conditions not conform to the specifications may have caused damage to the devices.
- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the devices if sterility cannot be guaranteed.
- Check the expiry date shown on the label attached. Do not use the devices after this date.
- The devices must be used immediately after opening the sterile packaging.
- For single use only. Reuse of the devices brings with it the risk of cross-infection regardless of the cleaning or sterilization method used.
- Do not resterilise, do not undergo any further processing. Resterilization may not be effective and may compromise the integrity of the device thus resulting in unintended injury.
- Tubing should be connected in such a manner as to prevent kinks or narrowings.
- Avoid alcohol, anaesthetic fluids, or corrosive solvents from coming into contact with the devices as they may jeopardize the structural integrity.

PRECAUTION

- Never unpack or prepare the table pack lines close to the patient's surgical wound. Detachable parts not well secured could fall into the patient's cavity without being noticed.
- After use, dispose of the devices in accordance with the applicable regulations in force in the country where the device is used.

INSTALLATION

Refer to "Tubing set-up" in the Performer HT user's manual.

WARNING

Perform the procedure using an aseptic technique.

CIRCUIT PRIMING PROCEDURE

Carefully follow the step-by-step instructions provided in the PERFORMER HT user's manual for a correct priming of the disposable circuit.

PRECAUTION

- 1) For priming and rinsing procedures it is recommended to use at least 2000 ml of sterile saline solution.
- 2) During these phases, it is advisable to lightly tap the outside of the tubing in order to facilitate bleeding all the air from the circuit, especially from low flow segment.

THERAPY

PREPARATION, PATIENT FILLING PHASE

Refer to "Preparation to Patient filling phase" in the Performer HT user's manual.

PRECAUTION

Select the appropriate volume and temperature on the equipment and load the volume of solution to perfuse the patient, in this phase carefully ensure that the air coming from the withdrawal table pack is adequately held back and/or eliminated by the filtration system integrated in the reservoir hung on the device scale.

WARNING

The surgical team is directly responsible for the correct catheters placement, as well as for selecting catheters different from those included in the kits R9900067 and R9900110.

Hang&Go HT-2

ENGLISH

WARNING

The above mentioned devices have been designed for use up to a maximum pressure of 500 mmHg; Rand is not responsible for improper setting of the parameters or accidental over-pressurization of the devices exceeding this safety limit.

TREATMENT, GENERAL NOTES:

- For a correct execution of this phase, follow the conditions and the indications provided with the PERFORMER HT user's manual.
- During perfusion it is advisable to check the state of saturation of the filters integrated in the semi-rigid reservoir hung on the equipment scale. Irrespectively on the patient's pathology, treatment duration and surgical approaches (closed or open abdomen perfusion), make sure that the filters are pervious throughout the treatment phase, by monitoring the visualized pressures on the PERFORMER HT screen and by visually cross-checking.

PRECAUTION

RAND advises to reach the right treatment parameters gradually, to ensure a proper therapy start as well as a correct treatment execution.

TREATMENT, DRUG INJECTION (OPTIONAL):

- Rand suggests injecting medication or active molecules through the Clave connectors (blue) integrated at different points in the circuit.
- In order to obtain the best therapeutic result, it is advisable to inject drug only when the perfused district has already reached the required treatment temperature.

WARNING

Every decision about the temperature level, duration of the perfusion and catheter positioning is remitted to the physician responsible of the therapy, as well as the dosage and drug injection method, consistent with the established usage limits of Rand devices.

The compatibility of use of the quoted devices with some commonly used medications - at the discretion of the physician - in intraperitoneal and/or intrapleural hyperthermal perfusion treatments has been appropriately verified by careful review of scientific literature and analysis of clinical test data relating to use of commercial medications based on the following active ingredients of consolidated use: CARBOPLATIN, CISPLATIN, MITOMYCIN C, ADRIAMYCIN/DOXORUBICIN, IRINOTECAN, MELPHALAN and OXALIPLATIN. Rand does not assume any responsibility for use of the devices in combination with medications or preparations not previously evaluated.

PRECAUTION

PVC medical devices must not be used for infusion of lipid-based drugs (e.g. Paclitaxel) capable of extracting DEHP or otherwise incompatible with phthalates. It is recommended to carefully read the information leaflet on the drug to be infused before using the device in order to exclude any incompatibility.

END OF TREATMENT LIQUID RINSING AND RECOVERY

Intra-peritoneal and/or intra-pleural hyperthermal perfusion treatments include a recovery and a rinsing phase. At the end of the therapy phase a simultaneous recovery and rinsing phase is started for the removal of the solution used for hyperthermal perfusion and of possible waste products. These steps are accomplished by perfusing the anatomical district with a predetermined volume of fresh washing solution. Due to potential high residual toxicity of the perfusion solutions and to allow the medical staff to perform safe perfusion and rinsing liquids recovery and disposal procedures, follow the indications below:

RINSING PHASE:

Refer to "Preparation to rinsing phase" in the Performer HT user's manual.

COMPONENTS REPLACEMENT

As a precaution a spare device should always be available during priming and therapy phases. If you believe that a particular event (such as insufficient performances, leaks, filter saturation, etc.) may affect the patient/user safety, interrupt the treatment as explained in the PERFORMER HT equipment operating manual (Section "Forced Interruption of the Treatment") and proceed to the line replacement.

MEDICAL DEVICES TO BE USED IN COMBINATION

The disposable devices manufactured by RAND have been specifically designed for use in combination with the RAND PERFORMER HT equipment.

Hang&Go HT-3

INNOVATION IN MEDICAL TECHNOLOGY

ENGLISH

The compatibility of use of the Mon-a-therm temperature probes (code 90050) with the PERFORMER HT machine has been adequately established.

The Ch 24 round silicone Catheters (code R9900071) have been specifically designed for use in combination with the "Infusion Table Pack" (code R2500158) and "Withdrawal Table Pack" (code R2500157) lines.

The Flow Reverse (code R9900101) is an accessory intended to be connected between the Hang&Go equipment lines and the table pack lines. For set-up and use of this accessory, please refer to the instructions for use provided in the Flow Reverse package.

WARNING

The compatibility of use of the devices listed in this instruction leaflet has been appropriately proved and verified through in vitro and in vivo tests. Rand will not be liable for any consequences, problems or malfunctions deriving from the use of the system in conjunction with different devices not previously evaluated.

LIMITED WARRANTY

This limited warranty is in addition to any statutory rights of the purchaser pursuant to applicable law.

RAND warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. RAND warrants that the medical device is capable of functioning as indicated in the current instruction for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, RAND cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performances and effectiveness of the device with damaging consequences for the patient, even though the specified instruction for use have been respected. RAND, whilst emphasizing the need to adhere strictly to the instruction for use and to adopt all the precautions necessary for correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. RAND undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by RAND up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of RAND or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device excepted as expressly stated herein. RAND disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this limited warranty and in particular agrees, in the event of a dispute or litigation with RAND, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

Hang&Go HT-4

ENGLISH

CAUTION

- 1) Before starting, check once more that the devices connected to the circuit and/or used in combination with it are correctly filled.
- 2) For an adequate therapy monitoring connect the temperature probes cabling to the appropriate jacks on the HTS module of the PERFORMER HT system (from T1 to T8):
 - By using the appropriate cable, connect the temperature probes located on the A/V line to channels T1 and T3 on the HTS module;
 - By using the appropriate cable, connect the temperature probe located on the A/V line to channel T2 of the HTS module;
 - Connect any other temperature probe to the corresponding channels on the HTS module, if applicable;
- 3) Open the external package, release the R2500183 A/V line table pack and clamp both the venous cannulae return line and withdrawal cannulae line. Remove the "U" silicone connector which joins the withdrawal and infusion segment line and connect them respectively to the venous and arterial cannulae previously positioned in the perfusion district vascular accesses.
- 4) When the blood begins to flow into the circuit, start the oxygen flow towards the oxygenator gas inlet. To avoid gas emboli in the blood circuit, make sure that the pressure inside the blood compartment is always higher than the pressure in the gas compartment.
- 5) RAND advises to reach the right treatment parameters gradually, to ensure a correct therapy start.

WARNING

The ILP devices have been designed for use up to a maximum pressure of 500 mmHg. Rand accepts no responsibility whatsoever for improper setting of the parameters or accidental pressurisation of the devices to levels incompatible with the stated limit.

THERAPY

GENERAL INSTRUCTIONS:

For a correct execution of this phase, follow the indications provided with the PERFORMER HT equipment (user's manual and equipment display).

- During the extracorporeal circulation phase it is recommended to continuously monitor the level of the blood drained from the limb / organ venous oxygenator reservoir, and properly adjust the blood flow to avoid the total depletion of the venous cardiomy reservoir and the consequent possible access of air in the blood circuit.
- Few minutes after the therapy start it is recommended to check the blood gas composition (by means of a hemogasalyzer); if necessary correct the oxygen concentration of the ventilation gas (FiO2) by turning the mixing flowmeter adjustment knob.
- Constantly monitor the correct dose and efficacy of the adopted anticoagulant treatment, before, during and after the therapy.

DRUG INJECTION (OPTIONAL)

- Rand recommends to inject the drug directly into the oxygenator venous reservoir or through the appropriate clave-connector fitted to the oxygenator inlet line, code R2500184.
- In order to obtain the best therapeutic result it is recommended to inject the drug only when the perfused limb/organ has already reached the therapy established temperature.

ILP HT-4

ENGLISH

WARNING

Every decision about the temperature level, perfusion duration and cannulae positioning is remitted to the physician responsible of the therapy, as well as the dosage and drug injection method, in accordance with the established usage limits of Rand devices.

The compatibility of use of the quoted devices with some commonly used medications - at the discretion of the physician - in isolated hyperthermic perfusion treatments has been appropriately verified by careful review of scientific literature and analysis of clinical test data relating to use of commercial compounds based on the following active ingredients of consolidated use: TNF- α , ADRIAMYCIN/DOXORUBICIN, OXALIPLATIN, and MELPHALAN. Rand does not assume any responsibility for use of the devices in combination with medications or preparations not previously evaluated.

CAUTION

PVC medical devices must not be used for infusion of lipid-based drugs (e.g. Paclitaxel) capable of extracting DEHP or otherwise incompatible with phthalates. It is recommended to carefully read the information leaflet on the drug to be infused before using the device in order to exclude any incompatibility.

THERAPY COMPLETION; RECOVERY OF FLUIDS

"LRT" therapies contemplate the recovery of the fluids used during the hyperthermic perfusion, which may be used in combination with drugs. Due to their high residual toxicity level, the perfusion fluids must be physically isolated and collected separately. To this purpose the "ILP" circuits provide a special system for the recovery and management of the fluids used during the therapy. The entire volume of such fluids, in the worst-case scenario, can be collected in suitable waste bags (code R2500188) that allow a safe management of the emptying phase.

RINSING//REINFUSION PHASE:

Before activating this phase:

- Hang the solution bag to be used for the rinsing procedure on the stand;
- Connect the additional 5 litre waste bag (code R2500188, included in the set) to the connector fitted to the A/V line, code R2500183.

PARTS REPLACEMENT

As a precaution, a spare device must always be available during the priming phase and in the course of the therapy.

If you feel that a particular event (such as insufficient performances, leaks, abnormal blood parameters, et c.) might affect the patient/user safety, proceed to the line replacement as explained in the PERFORMER HT equipment operating manual.

MEDICAL DEVICES THAT CAN BE USED IN COMBINATION

The ILP HT line set, code R9900086, has been specifically designed for use in combination with the RAND PERFORMER HT apparatus.

The line set has been optimised for connection to oxygenators whose specifications are compatible with those set out in the paragraph "technical specifications" and which are fitted with the following connectors: 1/4" connector for blood inlet and outlet, 3/8" connector for rigid venous reservoir and Hansen-type connector for the heat exchanger.

WARNING

The compatibility of the devices listed in the present instructions has been proved and verified with significant in vitro and in vivo tests. Rand will not be liable for any consequences, problems or malfunctions deriving from the use of the system in conjunction with different devices not previously evaluated.

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



(Modena)

Tel. +39 0535/660636

Design of the HANG&GO HT BASIC label (Rand S.p.A.)

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The label will feature the following handling and precaution symbols:

	Catalogue number		Date of manufacture
	Use by		Sterile - Ethylene oxide sterilized
	Non-pyrogenic		Do not reuse
	The device contains phthalates		Kit in conformity with article 12 of MDD 93/42/EEC

INNOVATION IN MEDICAL TECHNOLOGY

Fax 0535/660636
0535/5850360

	Do not use if package is damaged		CE Mark. This medical device fully complies with the European Council directive MDD 93/42 EEC
	Temperature limitation		Keep dry
	Keep away from heat		Caution, consult accompanying documents

INNOVATION IN MEDICAL TECHNOLOGY

Rand S.p.A.
Via Statale 12, n. 62
41036 Medolla (Modena)
Tel. +39 0535 49283 Fax +39 0535 660636
P.I. 02278000369

statement and signature

I certify that, in my capacity as Executive Vice President of Rand S.p.A., I believe to the best of my knowledge that all data and information contained in this document are truthful and accurate at time of publication.

This file has been fully revised and approved by:

Corrado Bellini
Executive Vice President
Rand S.p.A.

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CHAMBER OF COMMERCE, INDUSTRY, HANDICRAFT AND AGRICULTURE MODENA

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We certify that the signer is a duly authorized signatory on behalf of the Company for
any documents and acts to be used in trade with foreign countries.

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sign of the HANG&GO HT BASIC label (RanD S.p.A.)

- Company logo
- Name of the medical device
- CE mark
- Catalogue number
- Batch code
- Date of manufacture
- Use by...
- Sterile, Ethylene oxide sterilized
- Non-pyrogenic
- The product does not contain latex
- Do not reuse
- Temperature limitation: from 5° to 35°C
- Keep away from heat
- Keep dry
- The product cannot be re-sterilized
- The device contains phthalates
- Do not use if package is damaged
- Caution, consult accompanying documents
- Number and date of issue of the Medical Product Registration Certificate
- For questions related to the medical device, contact the Authorized Representative: ООО «Биокоммерц», 123592, Москва, ул. Кулакова, д. 20, стр. 1А. Тел: (495) 781-17-87; www.biocommerce.ru
- Barcode of the country of origin
- Name and address of the manufacturer company

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	Catalogue number		Date of manufacture
	Use by		Sterile - Ethylene oxide sterilized
	Non-pyrogenic		Do not reuse
	The device contains phthalates	§12 - 93/42EEC	Kit in conformity with article 12 of MDD 93/42/EEC

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