Medtronic

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As a responsible official of Medtronic, Inc., I hereby certify that the attached Duran AnCore Annuloplasty ring / band Information for Use (M059950T001) is a true copy of the original.

Joven Almazan
Regulatory Affairs
Medtronic, Inc.

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.


State of California
County of Orange

On July 30, 2018 before me, T. Gray-Ginger, Notary Public, personally appeared Joven Almazan, who proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that he executed the same in his authorized capacity, and that by his signature on the instrument the person, or the entity upon behalf of which the person acted, executed the instrument.


I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

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 **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

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Annuloplasty Ring and Band

1. Device description

The annuloplasty ring and band are available in the following models:

- Duran AnCore™ annuloplasty ring, Model 620R
- Duran AnCore™ annuloplasty band, Model 620B

The ring and band each consist of a flexible ring or band of polyester fabric with a 3 mm cross-section and a holder feature to facilitate replacement of the mitral valve chordae tendineae. The ring and band may be implanted in either the mitral or tricuspid position. The ring and band are marked at 3 points by colored suture. Two markers correspond to the trigones (mitral valve) or commissures (tricuspid valve), and a central marker identifies the midpoint of the ring or band. The individual ring or band size is identified by the circumferential distance between the trigone markers, as measured along the inside of the ring or band. A silicone marker, impregnated with barium sulfate, enables radiographic visualization.

2. Indications for use

The Duran AnCore annuloplasty ring and the Duran AnCore annuloplasty band are indicated for the reconstruction or remodeling of pathological mitral and tricuspid valves. Appropriate repair and annular remodeling may correct combined valvular insufficiency and stenosis.

Note: The chordal guide feature of the ring and band is indicated only for chordal replacement surgery of pathological mitral valves.

3. Contraindications

- Heavily calcified valves
- Valvular retraction with severely reduced mobility
- Active bacterial endocarditis

4. Warnings and Precautions

4.1. Warnings

- This device was designed for single use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Only surgeons who have received adequate training to determine whether incompetent, stenotic, or diseased heart valves are capable of being repaired or replaced should use this device.
- Only surgeons who have received appropriate training in valve repair, including ring and band implant and sizing techniques, should use this device.
- Do not use rings and bands that have been damaged or contaminated from patient contact.
- Correct annuloplasty ring and band sizing is an important element of a successful valve repair. Using a ring or band that is too small can result in dehiscence or valve stenosis. Using a ring or band that is too large can result in valve regurgitation.
- During implantation, place horizontal, interrupted, or continuous sutures into the annulus following the line of leaflet insertion.
- Ensure that sutures are not placed in the atrial tissue, as this may result in impairment of the cardiac conduction system.
- Secure the ring or band into fibrous trigone tissue to maintain permanent attachment.
- Securely tie all suture knots. Loose knots and long suture tails may be a source for hemolysis, thrombosis, or thromboembolism.
- Do not cut the ring or band. Loose threads can be a source of hemolysis, thrombosis, or thromboembolism.
- Correct placement and anchoring of artificial chords are important elements of successful valve repair. Improper placement and anchoring of artificial chords can result in leaflet tethering or regurgitation.
- Determining the correct length of artificial chords is an important element of successful valve repair. Artificial chords that are too long can result in mitral valve prolapse. Artificial chords that are too short can result in restricted leaflet motion.
- During implantation, the holder can be used to assist in placement of the ring or band, tying artificial chords, and testing for valve competency. The holder must be removed from the ring or band at the end of the procedure for proper function.
- Use intraoperative or postoperative echocardiography to evaluate the effectiveness of the valve repair. Minimizing regurgitation and preventing systolic anterior motion (SAM) are important elements of an effective repair.
- Surgeons who use annuloplasty rings and bands should be current on all anticoagulation regimens.
- Carefully monitor the patient's anticoagulant status when using postoperative anticoagulant therapy.
- Patients with intra-atrial thrombi or a giant left atrium may benefit from long-term anticoagulation therapy.
- Patients in atrial fibrillation may remain on anticoagulation therapy until sinus rhythm is established.

4.2. Precautions

- Do not use cutting edge needles. Cutting edge needles may damage the ring or band and potentially lead to dehiscence and possible mitral or tricuspid regurgitation.

5. Potential adverse events

While infrequent, certain complications have been reported with the use of annuloplasty ring and bands. Reported complications include:

- Uncorrected or recurrent regurgitation
- Stenosis
- Ring or band dehiscence
- Hemolysis (even with mild regurgitation)
- Low cardiac output
- Heart block
- Systolic anterior motion (SAM) and left ventricular outflow tract obstruction (LVOTO)
- Damage to coronary arteries
- Damage to coronary sinus
- Endocarditis
- Thrombosis
- Thromboembolism
- Anticoagulant-related hemorrhage

Consider the potential for these complications when selecting the most beneficial surgical procedure for each patient.

Use the sizing and implantation methods described in these Instructions for Use to avoid or minimize occurrence of these adverse events. Ensure that the physician has appropriate training and experience in valve repair and chordal replacement.

6. Instructions for use

6.1. Sizing

Select the appropriate ring or band size to help restore proper function of the valve. Use the Duran AnCore™ sizer set, Model 7620S, and the Medtronic annuloplasty handle, Model 7686, to select the proper ring or band size for repair of the mitral or tricuspid valve. The reusable sizers and handle are provided nonsterile and must be thoroughly cleaned and sterilized prior to use. Refer to the appropriate instructions for use for further information.

Align and insert the handle into the snap cavity of the sizer (Figure 2). The thin section of the handle is malleable to align the sizers with the annulus.

Warning: Do not use other manufacturers' annuloplasty sizers or sizers from other Medtronic annuloplasty products to size the Duran AnCore ring or band. Other sizers may not indicate the appropriate Duran AnCore ring or band size.

Note: The sizers can be used to gauge either mitral or tricuspid ring or band size. Rotate the sizer until the letter "M" (Mitral) or "T" (Tricuspid) is upright for proper orientation, assuming a median sternotomy has been performed.

Mitral valve

To determine the proper mitral ring or band size, both the distance between the annular trigones and the area of the anterior leaflet must be measured.

1. Lower the sizer onto the valve annulus and align the sizer notches with the annular trigones.
2. Gently extend the anterior (A) leaflet and cover its surface with the selected sizer.

Note: The sizer that has a notch spacing most nearly matching the intertrigonal distance and a surface area most nearly matching that of the anterior leaflet corresponds to the appropriate size of the ring or band (Figure 3).

Tricuspid valve

To determine the proper tricuspid ring or band size, both the distance between the commissures of the septal (S) leaflet and the area of the anterior leaflet must be measured (Figure 4).

1. Lower the sizer onto the valve annulus and align the sizer notches with the commissures of the septal leaflet.
2. Gently extend the anterior leaflet and cover its surface with the selected sizer.

The sizer that has a notch spacing most nearly matching the intercommissural distance of the septal leaflet and a surface area most nearly matching that of the anterior leaflet corresponds to the ring or band size that should be selected.

6.2. Handling and preparation instructions

1. Open the box. Remove the product literature and Patient Registration Form.

2. Remove the double-aseptic transfer pouch containing the ring or band mounted to the flange and holder.
3. Inspect the outer pouch. If the outer pouch is damaged, the exterior surface of the inner pouch may not be sterile.
4. Open the outer transfer pouch, and while still holding the bottom of the outer pouch, pass the inner pouch into the sterile field.
5. Inspect the inner pouch. The ring and band assemblies are sterile as long as the inner pouch has not been compromised. If the inner pouch is damaged, do not implant the ring or band.
6. Open the inner pouch only in the sterile field.

6.3. Device implantation

The ring and band assembly may be used with or without the Medtronic annuloplasty handle, Model 7686, provided separately.

1. Thoroughly clean and sterilize the handle prior to each use.
2. With the printed side of the holder facing up, align and engage the sterile handle into the snap fit cavity of the flange (Figure 5). The thin section of the handle is malleable to align the device with the annulus.
3. Remove the serial number identification tag (Figure 6). Record the serial number in the patient's record. Verify that the serial number matches the serial number on the Patient Registration Form.

Warning: The serial tag must be removed from the ring or band for proper function. Do not cut or tear the fabric of the ring or band during removal of the serial tag.

Mitral annulus suture placement

1. Place sutures in each trigone, approximately 4 mm in width.
2. Place additional interrupted sutures, approximately 4 mm in width, in the anterior and posterior portions of the mitral annulus (Figure 7). Place a total of 11 to 13 sutures for the ring. Place a total of 8 to 10 sutures for the band.

Warning: Do not place sutures in the circumflex coronary artery.

Tricuspid annulus suture placement

1. Place septal sutures at each commissure (C), approximately 4 mm in width, through the annulus at the base of the leaflet.
2. Place additional interrupted sutures, approximately 4 mm in width, along the tricuspid annulus (Figure 8). Place a total of 11 to 13 sutures for the ring. Place a total of 8 to 10 sutures for the band.

Warning: Do not place sutures in the coronary sinus, right coronary artery, atrioventricular (AV) node, or other conduction tissue (Figure 8).

Ring and band suture placement

1. Orient the ring or band assembly with the "M" or "T" on the holder for the mitral (M) annulus or the tricuspid (T) annulus (Figure 9 and Figure 10).
2. Pass sutures through the ring or band, approximately 2 mm in width, entering at the bottom of the ring or band and exiting through the periphery of the ring or band (Figure 11).

Flange and handle removal

1. Push the ring or band down onto the valve annulus while pulling back on the sutures.
2. Use a scalpel to cut the 1 flange retention suture (Figure 12).
3. Remove the flange and handle from the valve annulus (Figure 13). Dispose of the flange.

Knots and Testing

1. Tie all knots around the device securely.
2. Trim all excess sutures and test valvular competency with the holder still in place.
3. Cut the holder retention sutures at the 4 locations (Figure 14).
4. Grasp the central bar of the holder and pull the holder upward (Figure 15). Dispose of the holder.

Warning: Do not cut the fabric of the ring or band while cutting the flange or holder retention sutures.

Warning: The holder and flange must be removed from the ring or band at the end of the procedure for proper function. Never leave the holder attached to the ring or band.

Caution: The ring and band are designed to facilitate the implantation by interrupted suturing techniques only. If continuous suturing technique is desired, the ring or band must be removed from the holder assembly as previously described.

Note: The holder consists of 2 pieces that remain captive to each other. Do not attempt to separate the 2 holder pieces.

Tying mitral valve artificial chords

Note: Use conventional surgical chordal replacement techniques to place artificial, mitral valve chords.

1. Size the annuloplasty ring or band.
2. Place sutures in the mitral valve annulus and the ring or band. (See "Mitral suture placement".)
3. Pull the artificial chord(s) through the appropriate half of the holder before pushing the ring or band down onto the mitral valve annulus (Figure 16).
4. Push the ring or band down onto the annulus.
5. Remove the flange. (See "Flange and handle removal".)
6. Tie and trim all ring or band sutures. For each artificial chord, tie 2 to 3 loose knots using the holder chordal guide to set the length (Figure 17).
Caution: Do not push the initial knots for the artificial chords against the holder chordal guide. Tie the initial knots loosely to avoid displacing the papillary muscles and creating artificial chords that are too short.
7. Grasp the initial knots with forceps. Tie 7 to 8 tight square knots against the forceps to secure the artificial chord(s) (Figure 18).
8. Trim excess suture from the artificial chord(s).
9. Test valvular competency with the ring or band and holder in place.
10. Remove the holder and chordal guide by carefully cutting the holder retention sutures at the 4 locations (Figure 19).
11. Grasp the central bar of the holder and pull the holder upward (Figure 20). Dispose of the holder.

7. Individualization of treatment

Consider postoperative anticoagulation therapy for at least 6 weeks to allow for healing and incorporation of the annuloplasty ring or band by host tissue, regardless of cardiac rhythm.

8. Patient counseling information

Patients with annuloplasty rings or bands are at increased risk for endocarditis due to bacteremia (for example, when undergoing dental procedures). Advise patients about prophylactic antibiotic therapy.

Encourage patients to carry the Implanted Device Identification Card, provided by Medtronic, with them at all times.

9. How Supplied

9.1. Packaging

A single annuloplasty ring or band assembly (device, holder, and flange) is packaged in sterile, double-aseptic transfer pouches (Figure 1). The packaging system is designed to ease placement of the assembly into the sterile field. The assemblies are sterile if the pouches are undamaged and unopened. The outer surfaces of the outer pouch are nonsterile and must not be placed in the sterile field.

9.2. Available sizes

The Duran AnCore annuloplasty ring, Model 620R, and Duran AnCore annuloplasty band, Model 620B, are available in the following sizes for the mitral or tricuspid position: 23 mm, 25 mm, 27 mm, 29 mm, 31 mm, 33 mm, and 35 mm.

9.3. Storage

Store the product in its original packaging, including the outer shelf carton, in a clean, cool, and dry area to protect the product and minimize the potential for contamination.

The sterility and nonpyrogenicity of the Duran AnCore annuloplasty ring and band is validated to remain unaffected until the Use-by date identified on the shelf carton, provided the pouches are not opened or damaged.

10. Patient Information

10.1. Registration Information

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information, including countries from the European Union.

A patient registration form is included in each device package. After implantation, please complete all requested information. The serial number may be found on the package and on the identification tag attached to the device. Return the original form to the Medtronic address indicated on the form and provide the temporary identification card to the patient prior to discharge.

An Implanted Device Identification Card is provided to the patient. The card contains the name and telephone number of the patient's physician, as well as information that medical personnel would require in the event of an emergency.

11. Postoperative Information

11.1. MRI safety information

The Duran AnCore annuloplasty ring and band are magnetic resonance (MR) Safe. The device will not cause any harm to the patient when exposed to MR scanning immediately after implantation.

12. Accessories

Use the Duran AnCore sizer set, Model 7620S, to determine the appropriate ring or band size.

If using a handle, only use the Medtronic annuloplasty handle, Model 7686, to interface with the flange and sizers.

Caution: Do not use other manufacturers' annuloplasty sizers, or sizers for another Medtronic annuloplasty product, to size the Duran AnCore annuloplasty ring or band.

Caution: Do not use the accessories until they have been thoroughly cleaned and sterilized. Refer to the appropriate instructions for use for further instructions.

13. Disclaimer of warranty

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO UNITED STATES CUSTOMERS ONLY:

ALTHOUGH THE DURAN ANCORE ANNULOPLASTY RING, MODEL 620R AND DURAN ANCORE ANNULOPLASTY BAND, MODEL 620B, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

14. Disclaimer of warranty

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES:

ALTHOUGH THE DURAN ANCORE ANNULOPLASTY RING, MODEL 620R AND DURAN ANCORE ANNULOPLASTY BAND, MODEL 620B, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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

[Handwritten signature]

Medtronic

As a responsible official of Medtronic, Inc., I hereby certify that the attached Technical Details (IFU) is a true original.

Vivian Hirata
Regulatory Affairs
Medtronic, Inc.

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On October 9, 2018 before me, T. Gray-Ginger, Notary Public, personally appeared Vivian Hirata, who proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that he executed the same in his authorized capacity, and that by his signature on the instrument the person, or the entity upon behalf of which the person acted, executed the instrument.

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Duran AnCore Annuloplasty Ring and Band:

Description	Material (Trade Name)
Fabric Tube	Polyester (Polyethylene Terephthalate)
Radiopaque Strip	Silicone Medical Adhesive Barium Sulfate
Suture (Green)	Polyester (Braided Polytetrafluoroethylene-Impregnated) Colorant green
Suture (Green)	Polyester (Braided) Colorant green

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Suture (Colourless)	Polyester (Braided)
Holder Assembly/ Flange	Acetal Homopolymer
Holder Assembly Pins	Stainless Steel
Holder Assembly Sutures: - Colourless - Black	Polyester (Braided) Polypropylene Monofilament Colorant black
Identification Tag	Polytetrafluoroethylene

Duran AnCore Annuloplasty Ring and Band Accessories:

Description	Material (Trade Name)
Sizer Set Duran AnCore (7620S)	Polysulfone
Accessory Tray Duran AnCore (T7660): Base and Lid	Polyphenylsulfone
Annuloplasty Handle Duran AnCore (7620H): - stem - head	Stainless steel

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Devices supplied in packaging marked with the following information:

1. device name and reference number;
2. LOT of the product;
3. trademark and address of the manufacturer;
4. Single-Use (for band and rings);
5. Sterile, EO (for band and rings);
6. label complies with European Council 93/42 EEC;
7. date of production and expiry date;
8. storage conditions;
9. the characteristics of the product;
10. warning and informing signs.

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Medtronic
Duran AnCore® Ring with Chordal Guide 620R

Keep Away From Heat
Keep Dry
STERILE EO
Sterilized Using Ethylene Oxide
Do Not Reuse
Attention: See Instructions for Use
USA Rx Only
For US Audiences Only

Size
Serial Number
Use By

CE
0123

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
© 2015 Medtronic
140674002 Rev. 1.0

1409074002

Open Here

Medtronic
Duran AnCore® Band with Chordal Guide 620B
Annuloplasty Band - Mitral or Tricuspid / Bande pour annuloplastie - mitral ou tricuspidal / Annuloplasty Band - mitral oder tricuspidal / Banda para anulopectia - Mitral o tricuspidal / Annuloplastiering - mitral o tricuspidal / Аннулопластическая лента - митральная или трикуспидальная / Annuloplastyring - Mitral eller tricuspidal / Banda de anuloplastie - mitral ou tricuspidal

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
Manufactured at:
Earth City, CA USA
Attention: US Distributor
Ochsner, 2014
12117201 Rev 1.0

STERILE EO
Sterilized Using Ethylene Oxide

ISO 13485
Certification

CE
0123

USA Rx Only
For US Audiences Only

Size
Serial Number
Use By

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
© 2015 Medtronic
140674002 Rev. 1.0

1409076002

Open Here

Medtronic
Duran AnCore® Band with Chordal Guide 620B

Keep Away From Heat
Keep Dry
STERILE EO
Sterilized Using Ethylene Oxide
Do Not Reuse
Attention: See Instructions for Use
USA Rx Only
For US Audiences Only

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

REF ▶ T7660 LOT ▶

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-9004 USA

Manufactured at:
Minneapolis, MN USA
© Medtronic, Inc. 2007
MCE0000002 Rev. 1.0

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- code CFN (customer facing number);
- date and number of Registration Certificate;
- sign of conformity in Russia ;
- warning and informing symbols and inscriptions.
- The authorized official of manufacturer in Russian Federation

Sample of the sticker

	Product name	
	Name of the manufacturer	Code UPN
	Country of origin	Code CFN
	See IFU, Date and Number of RC	
NON-TOXIC		
The authorized official of manufacturer in Russian Federation:		

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Trademarks may be registered and are the property of their respective owners.