
APPROVED BY
REPROLIFE, Inc., Japan

Signature, seal/stamp, date
2021 8. 24 

Instructions for use for the medical device “Warming Kit for Vitri-fied Oocytes and Embryos, in variants”

Prepared in accordance with the requirements of the Recommendations of International Medical Devices Regulators Forum (IMDRF) and the Order of the Ministry of Health of the Russian Federation No. 11n dated January 19th, 2017.

CONTENT:

Instructions for use for medical professionals for variant 1: Warming Kit 102

Instructions for use for medical professionals for variant 2: Warming Kit 205

INSTRUCTIONS FOR USE FOR MEDICAL PROFESSIONALS

NAME OF THE MEDICAL DEVICE

Warming Kit for Vitrified Oocytes and Embryos, in variants:

Warming Kit 102 for Vitrified Oocytes and Embryos, consisting of:

1. Warming solution TS, 1.8 ml;
2. Diluent solution DS, 0.5 ml;
3. Washing solution WS, 1.0 ml;
4. Warming Plate;
5. Instructions for use.

RESPONSIBLE MANUFACTURER OF THE MEDICAL DEVICE

REPROLIFE, Inc., Japan
2-5-3-9F Shinjuku, Shinjuku-ku Tokyo,
160-0022,
Japan
Tel.: +813 59258931
info@reprolife.jp

MEDICAL DEVICE MANUFACTURING ADDRESS

REPROLIFE, Inc., Japan
2-5-3-9F Shinjuku, Shinjuku-ku Tokyo,
160-0022,
Japan
Tel.: +813 59258931
info@reprolife.jp

INTENDED USE OF THE MEDICAL DEVICE

The kit is intended for warming of the patient's oocytes and embryos in order to restore the activity of human cells after vitrification (freezing) as part of assisted reproductive technologies (ART) implementation.

FIELD OF APPLICATION OF THE MEDICAL DEVICE

The Kit is used in the field of reproductive technologies (IVF) for treating infertility associated with any disorders of the reproductive function, including uterine tubes patency disorder and commissures of female pelvic organs, as well as with disorders related to male infertility factor. The Medical Device is also used for warming of the oocytes of cancer patients, previously frozen before the beginning of treatment (chemotherapy).

Medical Device can also be used in scientific studies that require warming of biological objects.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Warming Kit for Vitrified Oocytes and Embryos, in variants, is intended for use in medical and health facilities that are certified to provide help in the field of assisted reproductive technologies (ART), such as clinics, laboratories, and IVF centers.

The Medical Device is to be handled by qualified embryologists with the necessary specialization and with relevant professional knowledge and skills, as well as experience in the field of assisted reproductive technologies (ART).

Specialized training in vitrification and warming method is recommended since the method has many technical nuances which play an important role in minimizing cell damage.

OPERATING PRINCIPLE OF THE MEDICAL DEVICE

Warming is a reverse process to vitrification. During warming the cryoprotector is removed from biological objects and is gradually being substituted by water, with the equalization of osmotic pressure in a controlled environment.

Kit component	Component role	Principle of operation	Contact duration
Warming solution TS, 1.8 ml	Warms the biological object	Hypo-osmotic solution.	1 min
Diluent solution DS, 0.5 ml	Removes cryoprotectors from the object	Hypo-osmotic solution. Removes cryoprotectors from the object	3 min
Washing solution WS, 1.0 ml	Saturates cells with water	Washes the object and equalizes the osmotic pressure	5 min
Warming plate	Solution vessel	Solutions are transferred to the plate wells, where manipulations with the objects are performed	10 min
Instructions for use	Medical Device Consumer Guidelines	-	-
The Kit is intended for 1 application.			

INDICATIONS FOR USE OF THE MEDICAL DEVICE

- After hormonal stimulation of the patient, vitrification (freezing) of embryos is performed to facilitate the 'natural cycle' transfer.
- It is recommended for creation of donor egg cell bank.
- It is recommended for creating the bank of donor oocytes from patients in childbearing age who were diagnosed with cancer and are planning to have children, before the start of chemotherapy.

CONTRAINDICATIONS FOR USE OF THE MEDICAL DEVICE

There is no special contraindications.

MEDICAL DEVICE APPLICATION CONDITIONS

Warming Kit for Vitrified Oocytes and Embryos, in variants, can only be used in the specialized IVF laboratory with an incubator and cryopreservation facility that is equipped for working with liquid nitrogen. The medical device is not to be used in any other conditions.

Operating temperature 25–27 °C, RH — 45–80 %, atmospheric pressure 84.0–106.7 kPa (630–800 mm Hg).

Since the shelf-life of the medical device is limited (1 year), it is necessary to examine the label and expiry date before use. Do not use the medical device if it is expired.

Warming Kit for Vitrified Oocytes and Embryos, in variants, should be stored and used in accordance with the instructions provided on the label and in the Instructions for use. Do not use the medical devices stored improperly.

Before use, inspect the package for damage, ruptures in the perforation lines, or control stickers. Do not use a package exhibiting external signs of damage, loss of the perforation or stickers integrity. Inspect the vials with solutions for damage after opening the package. Do not use the vials with signs of violated integrity, lacking or violated labeling, or solution streaks.

STERILIZATION

Sterile

Sterilization type – solutions are sterilized by filtration.

Warming plate and vials – by radiation sterilization.

Do not resterilize.

PRECAUTIONS

In vitro fertilization (IVF) and other assisted reproductive technologies departments should be in accordance with the 'Sanitary and Epidemiological Requirements for Health-Care Institutions'.

The staff of the institution should adhere to the security and safety measures generally established in laboratories, and use personal protection equipment — special clothing, gloves, and safety goggles.

The staff should be appropriately instructed and take precautions when working with liquid nitrogen and dewars.

Precautions for handling of the solutions

The solutions (TS, DS, WS) are not hazardous or toxic, and do not contain any substances classified as 'hazardous' or 'requiring special precautions' in handling and storage. Solutions are incombustible and non-flammable, stable under normal conditions.

Along with the precautions provided in Section 'Precautions' hereof, it shall be considered that upon a contact with the skin, solutions may cause irritation — redness, inflammatory reaction or eczema, depending on the individual sensitivity. Avoid spilling solutions, contact of the solutions with the eyes or unprotected areas of the body. Avoid swallowing the solutions.

Do not use a package exhibiting external signs of damage.

Do not re-use or re-sterilize.

Shelf life: 1 year. Examine the label and expiry date before use. Do not use the medical device if it is expired.

The use of these solutions in combination with solutions from other manufacturers for vitrification of oocytes and embryos, as well as warming after the vitrification, is strictly prohibited.

Access of unauthorized persons and unqualified staff to the laboratory should be restricted.

Access of unauthorized persons and unqualified staff to the Medical Device should be restricted.


MEDICAL DEVICE OPERATION PROCEDURE

Before use, inspect the package for damage, ruptures in the perforation lines, or control stickers. Do not use a package exhibiting external signs of damage, loss of the perforation or stickers integrity. Inspect the vials with solutions for damage after opening the package. Do not use the vials with signs of violated integrity, lacking or violated labeling, or solution streaks.

Check «Warming Kit 102 for Vitrified Oocytes and Embryos» contents, it should contain the following:

1. Warming Solution (TS), 1.8 ml (1 vial);
2. Diluent Solution (DS), 0.5 ml (1 vial);
3. Washing Solution (WS), 1.0 ml (1 vial);
4. Warming plate, 1 pc. (with 4 wells);
5. Instructions for use, 1 pc.

Warming Kit for Vitrified Oocytes and Embryos, in variants, can only be used in the specialized IVF laboratory with an incubator and cryopreservation facility that is equipped for working with liquid nitrogen. The medical device is not to be used in any other conditions.

Preparation

- The whole process should be performed at room temperature (25–27 °C)
Important! Use Pasteur pipettes with the inner diameter appropriate for oocyte and embryo (140–150 µm) or blastocyst (160–220 µm).
- Place the Warming Plate (closed with a lid) and sealed TS vial in the incubator at 37 °C for NLT 3 hours before warming procedure (overnight storage is recommended).
- Bring DS and WS vials to room temperature (25–27 °C) at least 1 hour before warming procedure.

- Retrieve the Cryotec vitrification carrier (hereinafter referred to as 'Cryotec') from the liquid nitrogen reservoir, remove the cap and bring Cryotec to the inner wall of the cooling rack filled with liquid nitrogen, leaning it up against the wall.

Warming and dilution

- Retrieve the Warming Plate from the incubator, open it and fill the second well of the Plate with 300 µl of DS.
- Take the TS out of the incubator and empty the contents into the first rectangular well of the Plate.
- Quickly (within 1 s) put the Cryotec into rectangular well with TS and wait for 1 min.
- Aspirate the oocyte/embryo and 3 mm of TS into the tip of the pipette and expel slowly at the bottom of the second well (with DS). Leave it for 3 min.
- While waiting, fill the third (WS1) and fourth (WS2) wells with 300 µl of WS each.
- Aspirate the oocyte/embryo and 3 mm of DS into the tip of the pipette and expel slowly at the bottom of the third well (WS1). Leave it for 5 min.
- Estimate the oocyte/embryo quality at this step depending on its recovery from the shrunken state.
- Place the oocyte/embryo on the surface of the fourth well (WS2). When oocyte/embryo sinks to bottom, aspirate it into the tip of the Pasteur pipette, and the once again put on the surface of the same well (WS2). Repeat twice.
- Put the oocyte/embryo in the droplet of the culture media for the recovery and subsequent micro-procedures.

Note: the time in the culture media is 2 to 4 hours for oocytes and 3 hours for embryos.

MEDICAL DEVICE STORAGE CONDITIONS

Store at 2–8 °C, keep away from direct sunlight.

Avoid freezing of the solutions. Proceed with caution during transfer.

Do not use after expiry date.

TRANSPORT CONDITIONS

Transportation at a temperature between 2 and 8 °C. Avoid freezing of the solutions.

Transportation is possible by any means of transport in covered vehicles, in compliance with the transportation conditions and in accordance with the shipping requirements and regulations for particular types of vehicles.

MEDICAL DEVICE DISPOSAL CONDITIONS

After being used, the Medical Device shall be disposed of as potentially infectious waste in accordance with the public standards and regulations adopted by the health-care institution.

Hazard class of medical device waste, considering biological hazard: epidemiologically hazardous waste - class B.

An unused expired Medical Device is to be disposed of as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution.

Medical Device is not to be utilized along with household waste.

Instructions for use for medical professionals
Warming Kit for Vitrified Oocytes and Embryos, in variants:
Warming Kit 102
Draft the Russian Federation market

REQUIREMENTS TO MAINTENANCE AND REPAIR

The Medical Device does not contain any components requiring maintenance or repair.
All products are single-use.

MEDICAL DEVICE COMPOSITION

Solutions composition

Solutions	Warming solution TS	Diluent solution DS	Washing solution WS
Components			
HEPES (4-(2-hydroxyethyl)-1-piperazineethanesulfonic acid) in the basal culture medium M-199)	contains	contains	contains
Hydroxypropyl cellulose	contains	contains	contains
Trehalose, free from endotoxins	contains	contains	contains
Phenol red	contains	contains	contains

Warming plate composition

Polystyrene, liquid paraffin, other additives.

Medical device does not contain pharmaceuticals as well as materials of animal and/or human origin.

CHARACTERISTICS OF THE SOLUTIONS

Characteristic	Value
Appearance	Clear pink solutions
Odor	Odourless
pH	7.2 (range between 7.1 and 7.25)
Stability	Solutions are stable

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE RUSSIAN FEDERATION

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Marketing Authorization: **** dated ****

SYMBOLS USED FOR LABELING

Symbol	Decryption
	Please refer to the instructions for use
	Caution! Please refer to the instructions for use
	Date of manufacture (DD-MM-YYYY)
	Use by date (DD-MM-YYYY)
	Batch code
	Catalogue number
	Sterile
	Sterilized using irradiation
	Do not reuse
	Do not resterilize
	Do not use if the package is damaged
	Storage conditions: temperature range from 2 °C to 8 °C
	Keep away from sunlight
	Keep dry
	Separate collection for waste
	Manufacturer
	Fragile, handle with care

The lot of solutions is indicated on the vials, the lot of the plate is on its packaging.

INSTRUCTIONS FOR USE FOR MEDICAL PROFESSIONALS

NAME OF THE MEDICAL DEVICE

Warming Kit for Vitrified Oocytes and Embryos, in variants:

Warming Kit 205 for Vitrified Oocytes and Embryos, consisting of:

1. Warming solution TS, 1.8 ml, 5 pcs.;
2. Diluent solution DS, 1.8 ml;
3. Washing solution WS, 1.8 ml, 2 pcs.;
4. Instructions for use.

RESPONSIBLE MANUFACTURER OF THE MEDICAL DEVICE

REPROLIFE, Inc., Japan
2-5-3-9F Shinjuku, Shinjuku-ku Tokyo,
160-0022,
Japan
Tel.: +813 59258931
info@reprolife.jp

MEDICAL DEVICE MANUFACTURING ADDRESS

REPROLIFE, Inc., Japan
2-5-3-9F Shinjuku, Shinjuku-ku Tokyo,
160-0022,
Japan
Tel.: +813 59258931
info@reprolife.jp

INTENDED USE OF THE MEDICAL DEVICE

The kit is intended for warming of the patient's oocytes and embryos in order to restore the activity of human cells after vitrification (freezing) as part of assisted reproductive technologies (ART) implementation.

FIELD OF APPLICATION OF THE MEDICAL DEVICE

The Kit is used in the field of reproductive technologies (IVF) for treating infertility associated with any disorders of the reproductive function, including uterine tubes patency disorder and commissures of female pelvic organs, as well as with disorders related to male infertility factor. The Medical Device is also used for warming of the oocytes of cancer patients, previously frozen before the beginning of treatment (chemotherapy).

Medical Device can also be used in scientific studies that require warming of biological objects.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Warming Kit for Vitrified Oocytes and Embryos, in variants, is intended for use in medical and health facilities that are certified to provide help in the field of assisted reproductive technologies (ART), such as clinics, laboratories, and IVF centers.

The Medical Device is to be handled by qualified embryologists with the necessary specialization and with relevant professional knowledge and skills, as well as experience in the field of assisted reproductive technologies (ART).

Specialized training in vitrification and warming method is recommended since the method has many technical nuances which play an important role in minimizing cell damage.

OPERATING PRINCIPLE OF THE MEDICAL DEVICE

Warming is a reverse process to vitrification. During warming the cryoprotector is removed from biological objects and is gradually being substituted by water, with the equalization of osmotic pressure in a controlled environment.

Kit component	Component role	Principle of operation	Contact duration
Warming solution TS, 1.8 ml	Warms the biological object	Hypo-osmotic solution.	1 min
Diluent solution DS, 1.8 ml	Removes cryoprotectors from the object	Hypo-osmotic solution. Removes cryoprotectors from the object	3 min
Washing solution WS, 1.8 ml	Saturates cells with water	Washes the object and equalizes the osmotic pressure	5 min
Instructions for use	Medical Device Consumer Guidelines	-	-
The Kit is intended for 5 applications.			

INDICATIONS FOR USE OF THE MEDICAL DEVICE

- After hormonal stimulation of the patient, vitrification (freezing) of embryos is performed to facilitate the 'natural cycle' transfer.
- It is recommended for creation of donor egg cell bank.
- It is recommended for creating the bank of donor oocytes from patients in childbearing age who were diagnosed with cancer and are planning to have children, before the start of chemotherapy.

CONTRAINDICATIONS FOR USE OF THE MEDICAL DEVICE

There is no special contraindications.

MEDICAL DEVICE APPLICATION CONDITIONS

Warming Kit for Vitrified Oocytes and Embryos, in variants, can only be used in the specialized IVF laboratory with an incubator and cryopreservation facility that is equipped for working with liquid nitrogen. The medical device is not to be used in any other conditions.

Operating temperature 25–27 °C, RH — 45–80 %, atmospheric pressure 84.0–106.7 kPa (630–800 mm Hg).

Since the shelf-life of the medical device is limited (1 year), it is necessary to examine the label and expiry date before use. Do not use the medical device if it is expired.

Warming Kit for Vitrified Oocytes and Embryos, in variants, should be stored and used in accordance with the instructions provided on the label and in the Instructions for use. Do not use the medical devices stored improperly.

Before use, inspect the package for damage, ruptures in the perforation lines, or control stickers. Do not use a package exhibiting external signs of damage, loss of the perforation or stickers integrity. Inspect the vials with solutions for damage after opening the package. Do not use the vials with signs of violated integrity, lacking or violated labeling, or solution streaks.

STERILIZATION

Sterile

Sterilization type – solutions are sterilized by filtration.

Warming plate and vials – by radiation sterilization.

Do not resterilize.

PRECAUTIONS

In vitro fertilization (IVF) and other assisted reproductive technologies departments should be in accordance with the ‘Sanitary and Epidemiological Requirements for Health-Care Institutions’.

The staff of the institution should adhere to the security and safety measures generally established in laboratories, and use personal protection equipment — special clothing, gloves, and safety goggles.

The staff should be appropriately instructed and take precautions when working with liquid nitrogen and dewars.

Precautions for handling of the solutions

The solutions (TS, DS, WS) are not hazardous or toxic, and do not contain any substances classified as ‘hazardous’ or ‘requiring special precautions’ in handling and storage. Solutions are incombustible and non-flammable, stable under normal conditions.

Along with the precautions provided in Section ‘Precautions’ hereof, it shall be considered that upon a contact with the skin, solutions may cause irritation — redness, inflammatory reaction or eczema, depending on the individual sensitivity. Avoid spilling solutions, contact of the solutions with the eyes or unprotected areas of the body. Avoid swallowing the solutions.

Do not use a package exhibiting external signs of damage.

Do not re-sterilize.

Shelf life: 1 year. Examine the label and expiry date before use. Do not use the medical device if it is expired.

The use of these solutions in combination with solutions from other manufacturers for vitrification of oocytes and embryos, as well as warming after the vitrification, is strictly prohibited.

Access of unauthorized persons and unqualified staff to the laboratory should be restricted.

Access of unauthorized persons and unqualified staff to the Medical Device should be restricted.


MEDICAL DEVICE OPERATION PROCEDURE

Before use, inspect the package for damage, ruptures in the perforation lines, or control stickers. Do not use a package exhibiting external signs of damage, loss of the perforation or stickers integrity. Inspect the vials with solutions for damage after opening the package. Do not use the vials with signs of violated integrity, lacking or violated labeling, or solution streaks.

Check «Warming Kit 205 for Vitrified Oocytes and Embryos» contents, it should contain the following:

1. Warming Solution (TS), 1.8 ml (5 vials);
2. Diluent Solution (DS), 1.8 ml (1 vial);
3. Washing Solution (WS), 1.8 ml (2 vials);
4. Instructions for use, 1 pc.

Warming Kit for Vitrified Oocytes and Embryos, in variants, can only be used in the specialized IVF laboratory with an incubator and cryopreservation facility that is equipped for working with liquid nitrogen. The medical device is not to be used in any other conditions.

Preparation

- The whole process should be performed at room temperature (25–27 °C)
Important! Use Pasteur pipettes with the inner diameter appropriate for oocyte and embryo (140–150 µm) or blastocyst (160–220 µm).
- Place the Warming Plate (closed with a lid) and sealed TS vial in the incubator at 37 °C for NLT 3 hours before warming procedure (overnight storage is recommended).
- Bring DS and WS vials to room temperature (25–27 °C) at least 1 hour before warming procedure.

- Retrieve the Cryotec vitrification carrier (hereinafter referred to as 'Cryotec') from the liquid nitrogen reservoir, remove the cap and bring Cryotec to the inner wall of the cooling rack filled with liquid nitrogen, leaning it up against the wall.

Warming and dilution

- Retrieve the Warming Plate from the incubator, open it and fill the second well of the Plate with 300 µl of DS.
- Take the TS out of the incubator and empty the contents into the first rectangular well of the Plate.
- Quickly (within 1 s) put the Cryotec into rectangular well with TS and wait for 1 min.
- Aspirate the oocyte/embryo and 3 mm of TS into the tip of the pipette and expel slowly at the bottom of the second well (with DS). Leave it for 3 min.
- While waiting, fill the third (WS1) and fourth (WS2) wells with 300 µl of WS each.
- Aspirate the oocyte/embryo and 3 mm of DS into the tip of the pipette and expel slowly at the bottom of the third well (WS1). Leave it for 5 min.
- Estimate the oocyte/embryo quality at this step depending on its recovery from the shrunken state.
- Place the oocyte/embryo on the surface of the fourth well (WS2). When oocyte/embryo sinks to bottom, aspirate it into the tip of the Pasteur pipette, and the once again put on the surface of the same well (WS2). Repeat twice.
- Put the oocyte/embryo in the droplet of the culture media for the recovery and subsequent micro-procedures.

Note: the time in the culture media is 2 to 4 hours for oocytes and 3 hours for embryos.

MEDICAL DEVICE STORAGE CONDITIONS

Store at 2–8 °C, keep away from direct sunlight.

Avoid freezing of the solutions. Proceed with caution during transfer.

Do not use after expiry date.

TRANSPORT CONDITIONS

Transportation at a temperature between 2 and 8 °C. Avoid freezing of the solutions.

Transportation is possible by any means of transport in covered vehicles, in compliance with the transportation conditions and in accordance with the shipping requirements and regulations for particular types of vehicles.

MEDICAL DEVICE DISPOSAL CONDITIONS

After being used, the Medical Device shall be disposed of as potentially infectious waste in accordance with the public standards and regulations adopted by the health-care institution.

Hazard class of medical device waste, considering biological hazard: epidemiologically hazardous waste - class B.

An unused expired Medical Device is to be disposed of as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution.

Medical Device is not to be utilized along with household waste.

REQUIREMENTS TO MAINTENANCE AND REPAIR

The Medical Device does not contain any components requiring maintenance or repair.

MEDICAL DEVICE COMPOSITION

Solutions composition

Solutions Components	Warming solution TS	Diluent solution DS	Washing solution WS
HEPES (4-(2-hydroxyethyl)-1-piperazineethanesulfonic acid) in the basal culture medium M-199)	contains	contains	contains
Hydroxypropyl cellulose	contains	contains	contains
Trehalose, free from endotoxins	contains	contains	contains
Phenol red	contains	contains	contains

Warming plate composition

Polystyrene, liquid paraffin, other additives.

Medical device does not contain pharmaceuticals as well as materials of animal and/or human origin.

CHARACTERISTICS OF THE SOLUTIONS

Characteristic	Value
Appearance	Clear pink solutions
Odor	Odourless
pH	7.2 (range between 7.1 and 7.25)
Stability	Solutions are stable

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE RUSSIAN FEDERATION

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Marketing Authorization: **** dated ****

SYMBOLS USED FOR LABELING

Symbol	Decryption
	Please refer to the instructions for use
	Caution! Please refer to the instructions for use
	Date of manufacture (DD-MM-YYYY)
	Use by date (DD-MM-YYYY)
	Batch code
	Catalogue number
	Sterile
	Sterilized using irradiation
	Do not reuse
	Do not re sterilize
	Do not use if the package is damaged
	Storage conditions: temperature range from 2 °C to 8 °C
	Keep away from sunlight
	Keep dry
	Separate collection for waste
	Manufacturer
	Fragile, handle with care

The lot of solutions is indicated on the vials.

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
Registration No. 772 of 2021

NOTARIAL CERTIFICATE

This is to certify that Ms. Harumi Sato, an attorney-in-fact of Mr. Masahiro Mukai, CEO of REPROLIFE, INC. , declared in my very presence that said Mr. Masahiro Mukai had signed and sealed to the attached document and so the signature and seal were true and genuine.

Dated this 30th day of August, 2021

Yamashiki Yaichiro

YAMASHIKI Yaichiro
Notary, attached to
The Tokyo Legal Affairs Bureau
No. 1-10, Nihombashi, Kabuto-cho,
Chuo-ku, Tokyo, JAPAN

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APOSTILLE

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1. Country: JAPAN
This public document
2. has been signed by **YAMASHIKI Yaichiro**
3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
4. bears the seal/stamp of **YAMASHIKI Yaichiro, Notary**
Certified
5. at Tokyo
6. **AUG. 30. 2021**
7. by the Ministry of Foreign Affairs
8. 21- **№ 030079**
9. Seal/stamp:
10. Signature

HAMAMOTO Hiroki

For the Minister for Foreign Affairs

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