5th of March, 2021

Approved by

(signature)

Yoo Jae Hong

Director

Dentium Co., Ltd.(ICT Branch)

**INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE**

Intraoral dental digital sensor, versions: IOX1, IOX2

Edition 1

Medical device name**Intraoral dental digital sensor, versions: IOX1, IOX2****1. Version IOX1 consisting of:**

- IOX1 sensor with 2-m cable – 1 pc.
- Electronic media with the software – 1 pc.
- Wall mounted sensor cradle – 1 pc.
- Silicone sensor case – 1 pc.
- USB extension cable - 1 pc.
- Instruction for use – 1 pc.

1. Version IOX2 consisting of:

- IOX2 sensor with 2-m cable – 1 pc.
- Electronic media with the software – 1 pc.
- Wall mounted sensor cradle – 1 pc.
- Silicone sensor case – 1 pc.
- USB extension cable - 1 pc.
- Instructions for use – 1 pc.

Hereinafter – intraoral sensor, IOX1 sensor, IOX2 sensor, medical device, device.

2. Information about the manufacturer, developer, and authorized representative of the manufacturer**Manufacturer:**

Dentium Co., Ltd (ICT Branch)

76, Changnyong-daero 256beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do 16229, (KOREA)

Phone: + 82-70-7098-4027

Fax: + 82-31-281-5561

www.dentium.co.kr

**Developer:**

Dentium Co., Ltd (ICT Branch)

76, Changnyong-daero 256beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do 16229, (KOREA)

Phone: + 82-07-7098-4027

Fax: + 82-31-281-5561

www.dentium.co.kr

Authorized representative of the manufacturer in the Russian Federation

Limited Liability Company IMPLANT.RU (IMPLANT.RU, LLC)

Bolshoy Gnezdnikovskii alley, 3, Moscow, 125009

Phone/Fax: +7(495)118-76-20, +7(495)941-02-61

Classification

Risk class IIa according to the Directive MDD/93/42/EEC+2007/47/EC

Type of contact with the human body: intraoral sensor does not contact with the body; no contact. Medical personnel should perform all manipulations with the device using medical gloves. When using the intraoral sensor for its intended purpose, a disposable plastic hygienic cover should be put on it.

INSTRUCTION FOR USE

FOR THE MEDICAL DEVICE

Intraoral dental digital sensor, versions: IOX1, IOX2

4. Intended purpose of the device and potential consumer

4.1 Intended use of the device

The device "Intraoral dental digital sensor, versions: IOX1, IOX2" is intended to acquire digital intraoral X-ray images of teeth for the following purposes: diagnosis and treatment of diseases in the jaw and the oral cavity.

4.2 Potential consumers

The intraoral sensor is intended for use by the qualified dentists or radiologists.

Device design


Table 1 Device design

No.	Components	Quantity	Intended use
①②③	Sensor with cable	1 pc.	Basic device (sensor module) ①, cable ②, USB connector ③
④	Electronic media with the software (USB memory)	1 pc.	Drive for storing the device installation file, USB driver, calibration card file and user manual
⑤	Wall mounted sensor cradle	1 pc.	Sensor holder and protection device
⑥	Silicone sensor case	1 pc.	Protective case for sensor storage
⑦	USB extension cable	1 pc.	Sensor cable extension


5.1. Device diagrams*

* - All dimensions are in mm.


Sensor with cable


IOX1 sensor module


IOX2 sensor module

Electronic media with the software (USB memory)




Wall mounted sensor cradle


5


General tolerances (mm)			
Size	Tolerance	Size	Tolerance
0.5 ~ 3.0	± 0.05	120 ~ 315	± 0.5
3.0 ~ 6.0	± 0.05	315 ~ 1000	± 0.8
6.0 ~ 30	± 0.1	1000 ~ 2000	± 1.2
30 ~ 120	± 0.3	More than 2000	± 1.5

IOX1 sensor silicone case


IOX2 sensor silicone case

General tolerances (mm)			
Size	Tolerance	Size	Tolerance
0.5 ~ 3.0	± 0.05	120 ~ 315	± 0.5
3.0 ~ 6.0	± 0.05	315 ~ 1000	± 0.8
6.0 ~ 30	± 0.1	1000 ~ 2000	± 1.2
30	± 0.2	More than 2000	± 1.5

USB extension cable

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE
Intraoral dental digital sensor, versions: IOX1, IOX2

5.2. Intraoral sensor configuration


Table 2 Intraoral sensor configuration

Components	Description
Cover (top, bottom)	Internal module protection
Gasket plate	Support for control board, sensor and FOS
Control board	Sensor control
Sensor	Image acquisition
FOS (Fiber-optic plate with Csl scintillator)	Conversion of X-rays into visible light
Microporous polyvinyl chloride	Microcellular polyurethane foam
Carbon plate	Windowed cover and scintillator plate

6. Indications, contraindications and possible side effects

In all types of diagnostic dental procedures and treatment of diseases in the jaw and oral cavity, when images acquired with an X-ray machine and an intraoral sensor are required.

Contraindications: not detected.

Possible side effects: not detected.

7. Precautions

7.1. Precautions before use

- 1) The intraoral sensor can be used by duly trained dentists and radiologists.
- 2) The medical device should be used strictly in accordance with the operation manual.
- 3) Do not place identical medical electrical products near the device. This may lead to deterioration of the device performance. Keep the device away from the electromagnetic sources to avoid electromagnetic interference.
- 4) Damage to the device due to disassembly by the user is not covered by the warranty.
- 5) Pulling or twisting the cable connected to the device may cause deterioration and malfunction.

7.2. Precautions during operation

- 1) The intraoral sensor shall be kept in the sensor cradle after use.
- 2) Falling, vibration, impact or careless handling can cause damage that affects the image quality.
- 3) Do not use the device without a disposable hygienic cover.
- 4) Take care so that the sensor is not bitten by the patient during placement in the oral cavity.
- 5) Do not switch off the PC while the device is in operation.
- 6) Do not apply excessive force to the product, as this may lead to equipment malfunction or deterioration of its performance.
- 7) Failure to use device in accordance with the intended purpose may lead to deterioration of its performance and malfunction.
- 8) When continuously acquiring images for a long time, the sensor can generate heat. The temperature of the sensor module can exceed 41 °C (up to 41.8 °C).
- 9) Use the intraoral sensor only with a disposable hygienic cover registered in the territory of the Russian Federation. After each use, the cover is disposed of in accordance with the established procedure and replaced with a new one for each patient. (see clause 20.).
- 10) The silicone sensor cover is intended only to protect the device. Before using the sensor, remove the silicone cover.
- 11) The X-ray dose influences the image quality, so it is necessary to adjust the X-ray dose and the exposure time according to the recommendations.
- 12) To ensure high-quality operation of the intraoral sensor, use a PC only with the recommended characteristics (specifications).
- 13) When performing scanning, the patient should put on protective clothing and thyroid protection prior to scanning.
- 14) Pregnant women can undergo scanning only after consultation with a doctor.
- 15) If any malfunctions occur, please contact the authorized representative of the manufacturer on the territory of the Russian Federation.

16) Always keep the device clean to avoid foreign particles penetration.

7.3. Electrical precautions

- 1) Do not connect the USB connector of the intraoral sensor to any device other than the PC.
- 2) Do not install the device in the places where the USB connector may come into contact with water.
- 3) If the device does not function properly, it is necessary to contact an authorized representative of the manufacturer on the territory of the Russian Federation.

7.4. Warnings and cautions

- 1) Do not spill or spray liquid substances, such as water or aerosols, on the device.
- 2) Do not keep flammable materials close to the device. This may cause a malfunction.
- 3) In case of inflammation of the electrical wiring, use a fire extinguisher intended for this purpose. Use of water or liquids in case of wiring inflammation can lead to serious damage of the device and injury to people. Prior to fire extinguishing, disconnect the USB connector to avoid the risk of electric shock.
- 4) Neither sterilize nor disinfect the device with sprays.

 Warning	Neither disassemble nor modify the device!
--	--

 CAUTION	The intraoral sensor does not guarantee normal operation when being operated with an unregistered Microsoft Windows operating system. For this reason, it is required to use a registered original version of Microsoft Windows. Operating systems other than Windows are not supported.
---	--

8. Application procedure and operation principle

The sensor uses an analog-to-digital signal conversion technique to record X-rays absorption pattern for the diagnosis and treatment of diseases in the jaw and oral cavity. The device is used in combination with an X-ray source (X-ray tube), which is placed outside the patient's mouth during the X-ray imaging procedure; the resulting digital image is transmitted to the PC display.

Place the sensor into the oral cavity and adjust the X-ray tube so that it is perpendicular to the sensor placement position. The resulting image from the sensor is reconstructed on the basis of the primary image in the FPGA, transmitted to the USB controller and then to the PC via USB interface.

The transmitted image is reconfigured on the PC, after which the image is checked in the viewer.

Do not use the sensor without a disposable protective hygienic cover, which is registered in the Russian Federation in accordance with the established procedure.

8.1. Software driver installation


The intraoral sensor includes an integrated installer that installs USB and TWAIN drivers.

 WARNING	To use the intraoral sensor, it is required to install the device driver. The USB drive shall be connected to this or any other device meeting the requirements of IEC60601-1 or IEC 60950-1 standards.
--	---

8.1.1. Preparing for installation

Connect USB memory to the PC and check the folder.

Folder configuration includes "Installation", "Calibration" and the user manual as follows. The installation program required to acquire an image is presented in the form of two files – Install and USB Driver.

8.1.2. Installation of the USB driver


Go to the corresponding file on the USB memory.

Path: Installation\USB Driver

- 1) Select the folder suitable for the operating system for installation on the PC (starting from Vista, 7, 8, 8.1, 10, XP).
- 2) Select the folder matching the operating system bitness (32-bit: x86, 64-bit: x64)
- 3) Run the installation file (32-bit: DPInst32.exe, 64-bit: DPInst64.exe).
- 4) Click the Next button to complete the installation.

 Warning	Do not connect the intraoral sensor to the USB interface on the PC until the installation is completed successfully.
--	--

8.1.3. Description of the USB driver installation screen


- (1) When the Run screen appears, click the Next button.
- (2) A window will appear to confirm the driver installation, click "Install this driver software" to continue the installation.
- (3) As soon as the installation is complete, click Finish to exit.

8.1.4. TWAIN driver installation

The intraoral sensor can acquire images with the use of TWAIN-compatible applications. To use this option, both the images program and the patient program shall support the TWAIN interface. Path: Run the installation \ Install \ Dentium-IntraoralSensor-TWAINx.x.x. msi


Click the Next button to complete the installation.

8.1.5. Cable connections

As soon as the driver installation is complete, connect the intraoral sensor to the USB port on the PC.

- When connecting to the USB connector, we recommend to use the connector on the back of the PC.
- USB extension cable can also be used depending on the installation conditions and needs.

* USB extension cable: USB 2.0 limits the data transmission distance to 5 m, so it is recommended to use repeaters that amplify the signals for stable data transmission when using a USB extension cable.

 Warning	When using a USB extension cable, only use the power adapter connection type. In case of any problems in device operation after connection of the USB extension cable, it is necessary to contact the manufacturer's authorized representative in the Russian Federation.
--	--



8.1.6. Device driver verification

Launch the Device Manager on the Windows control panel. If you see "Dentium IntraOralSensor" in the "Universal Serial Bus Controller", the driver operates properly.

Verification method:

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE
Intraoral dental digital sensor, versions: IOX1, IOX2


- 1) Windows 7 or 10: Control Panel → System and Security → System → Device Manager
- 2) Windows XP: Settings → Control Panel → System → Hardware → Device Manager


If the device is not installed properly, it is necessary to contact the manufacturer's authorized representative on the territory of the Russian Federation.

8.2.1. System configuration

To acquire an image with the use of the intraoral sensor (①), connect the USB data connector (②) to the PC or laptop (④) which can monitor the sensor and be used to view the images. There are two types of generators that can generate X-radiation: portable and standard. When X-ray radiation is transmitted to the sensor with the use of the X-ray generator, the image is acquired automatically and can be viewed on the PC or laptop.

8.2.2. Image acquisition process


The intraoral sensor enables to acquire images through three simple steps. When the *sensor* is ready for imaging, it emits X-rays to produce an image.

8.2.3. Intraoral Sensor Imaging Program

The imaging method uses the TWAIN driver. TWAIN is one of the industrial standards for the common image acquisition interface, and it is possible to acquire images using an intraoral sensor with the help of an auxiliary program. This operation manual describes how to acquire images using an intraoral sensor with the help of Dentium image processing program. When using a third-party program, see the operation manual for it.

Image acquisition process is described below. When X-rays are generated into the intraoral sensor, the X-ray signal is converted to visible light by FOS, and the image is detected by the sensor. The produced image signal is transmitted to the PC by USB data transfer and can be viewed with the use of the image processing program.


8.3. Switching on the PC and checking the USB connection

As a rule, if the USB cable of the intraoral sensor is connected to the PC, image acquisition is possible. See the procedure below when the device is connected.


8.3.1. Launching “rainbow Image Viewer”

After starting the Image Viewer, it is necessary to select the patient. Click the Capture tab and then double-click on the Standard button.

8.3.2. Preparing for imaging

Click the Capture button to select the desired device in the Image Viewer. Click Select Source and then select Dentium IOX DS 1.0.

8.3.3. Preparation of the X-ray generator

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE
Intraoral dental digital sensor, versions: IOX1, IOX2


Prepare a portable X-ray generator or a standard X-ray generator for exposure.

8.3.4. Determining the sensor position for imaging

Place the intraoral sensor (after placing it in a disposable hygienic cover registered in accordance with the established procedure on the territory of the Russian Federation) in the patient's oral cavity so that the lower plane of the sensor is directed towards the X-ray generator.

8.3.5. Starting the imaging program

When the sensor is in place, press the Capture button to start the imaging program. At this stage, the device is ready to receive X-rays.

Warning

The settings button (⚙️) in the upper-right corner of the imaging program can only be used by the administrator. Users are not allowed to use this button.

8.3.6. Status description: the device information status is divided into 6 levels

No device: Indicates that device cable is disconnected.

Initialize: Initialization is performed after the device is connected to the PC. When the cable is disconnected, the connection is attempted only once.

Connected: Device initialization is performed properly.

Ready for operation: By clicking the Capture button when the device is connected, an X-ray image will be acquired.


In progress: Indicates the image acquisition stage upon completion of X-ray exposure.

Capture end: The image is acquired normally. When you click the Accept button, the acquired image will be transmitted to the image processing software (rainbow Image Viewer).

8.3.7. X-rays exposure

When the device is properly positioned in the patient's image area and the Capture button is clicked, the device is ready to acquire images. Exposure to X-rays is carried out after setting the appropriate parameters depending on the patient's condition.

8.3.8. Image capturing

The image acquired from the device is displayed on the screen as shown below. When the device is properly positioned in the patient's image area and the Capture button is clicked, the device is ready to acquire images. Exposure to X-rays is carried out after setting the appropriate parameters depending on the patient's condition.

Warning

The sanitary plastic cover is designed for single use and is replaced with a new one for each patient. Use a sanitary and hygienic plastic cover which meets the requirements of ISO 10993-1 and is registered in accordance with the established procedure in the Russian Federation.

Always keep the device clean to avoid foreign particles penetration.
Sensor's cover is intended for the protection and shall be removed prior to use.


Warning

Position the sensor on the X-ray generator so that it is directed exactly to the device bottom.

Take care to avoid patient movement during image acquisition.

When the intraoral sensor is located in the patient's mouth, care must be taken to avoid biting or hitting the device.

If you click the Capture button, the device will be ready to capture for 3 minutes . After receiving an image, the device will be ready to work again. The captured image can be previewed by clicking the left/right buttons. The resulting image can be saved by clicking the Accept button (in the continuous imaging mode, all captured images are saved).

8.3.9. Adjustment of the preset parameters for image processing

The acquired images can be selected from five levels of image filters according to the user's preferences. The default value is set to Default, and you can change sharpness by switching between Smooth + (Blur +) (left) to Sharp + (Sharp+) (right).

It is possible to select a filter level prior to image capturing or change the image filter after imaging. In the case of continuous capture mode, the filter can be changed individually at the preview stage, and the image can be saved as is.

Examples of images with processing filters are shown below.

Smooth+

Smooth


Default

Sharp


Sharp+

8.3.10. Check of the final image

It is possible to view the resulting image with the help of the patient management program. The saved image can be edited with the use of the options supported by the patient management program. See the program's user guide for more information.

Actual screen of the resulting image after imaging completion

8.3.11. Guide on exposure to X-rays

X-ray radiation dose required for the best image quality is determined as follows:

- ◆ X-ray source: tube voltage, tube current, tube module, etc.
- ◆ SID (distance from the source to the receiver of the image): the distance is measured between the focal point on the X-ray tube and the image receiver of the X-ray cassette
- ◆ Target teeth and tissues
- ◆ Patient's age and bone density

The X-ray dose affects the image quality. An insufficient dose increases the image noise, making it impossible to distinguish the details. Conversely, too high doses reduce distinguishing in the dark areas.

The X-ray dose is determined on the basis of the tube voltage (kVp), tube current (mA), and exposure time. The exposure time parameter for the dose may vary depending on the X-ray generator. See the table below for the recommended X-ray exposure time for different areas.

Table 3 X-ray exposure parameters for different areas

Exposure conditions	70 kV, 2 mA
SID	200 mm
Front teeth/incisors	600 μ Gy
	0.06~0.10 s
Molars	750 μ Gy
	0.10~0.13 s

SID - distance from the source to the receiver of the image.

Warning

X-ray exposure conditions can be changed on the basis of an expert opinion.

- For a large body-build: increase the source current by 25 %

- For a small body-build: reduce the source current by 20 %

- For edentulous patients: reduce the source current by 20 %

The X-ray dose required to produce an image may vary depending on the X-ray source and the environmental conditions.

9. Consumer properties

CMOS-based high-resolution images

Fully moisture-proof/dust-proof design with IP68 protection rating

Support of high-speed data transfer with the use of USB 2.0

Flexible cable made of silicone

Easy to use with PC

TWAIN interface (standard for image acquisition devices)

10. Specifications

Table 4 Intraoral sensor specification

Components	Specification	
	IOX1	IOX2
Pixel pitch	20 × 20 μm	
Dimensions (W x L x H)	24 × 36.7 × 12.1 mm	29.1 × 42.3 × 12.3 mm
Number of pixels	1000 × 1500	1280 × 1801
Analog-to-digital converter	12 bits	
Frame rate	1 frame/s	
Dynamic range	70 dB	
Data transfer interface	Parallel data bus (USB 2.0)	
USB cable length	2.0 m	
Rated electrical performance	5 V DC, 1 VA	
Ambient temperature	from 10 °C to 30 °C (operation) from - 20 °C to 60 °C (transport and storage)	
Relative humidity	30-95 % (operation) 10-95 % (transportation and storage)	
Atmospheric pressure	700-1060 HPa	
Electric shock protection	Type B applied part	
Scintillator type	CsI	
Sensor technology	CMOS	
High-contrast resolution, pl/mm, min.	6.3	
Electrical device class	IIa	

Protection against penetration of foreign substances/water	IP68	
Sensitive area (width × height)	20 x 30 mm	25.6 x 36 mm
Modulation transfer function (MTF) at 6.0 lp/mm, min.	30 %	
Quantum efficiency (DQE) at 2.5 lp/mm and the dose of 0.53 μGy, min.	50 %	
Rated electrical parameters		
Input voltage and frequency	5 V DC, 50/60 Hz	
Electric shock protection and applied part	Type 1, IP B	
Power consumption	1 ± 3% VA	
Voltage (from the USB port)	5 ± 5% V	

Table 5 Recommended PC specification

Component	Specification	
	Minimum	Recommended
Operating system	Microsoft Windows XP SP3 (32 bit)	Microsoft Windows 7 (32 bit/64 bit)
CPU	Core i5-3470 3.20 GHz or higher	
RAM	2 Gb	4 Gb
HDD	10 GB of free space on the hard disk	
Interface	USB 2.0	
VGA	Video Memory 256 Mb	NVIDIA GeForce GTX660 D5/1 Gb
Resolution	1920080 × 32 bit/pixel	

11. Electromagnetic compatibility data

Table 6 Electromagnetic emission

Manufacturer's recommendations and statement – electromagnetic radiation		
The intraoral sensor is intended for use in the electromagnetic environment described below. The customer or device user must ensure that the device is used in such environment.		
Emission tests	Conformity	Electromagnetic environment – recommendations
Radiated emission CISPR 11	Group 1 Class A	The intraoral sensor uses radio frequency energy only for the internal functioning. For this reason, its RF emission is extremely low and can hardly cause interference with the nearby equipment.
Disturbance voltage of the mains terminals CISPR 11	Group 1 Class A	
Harmonics IEC 61000-3-2	Class A	
Voltage fluctuations	Conforms	

NOTE: EMISSION characteristics of the device make it suitable for use in the industrial zones and hospitals (CISPR 11, Class A). If it is used in the residential environment (which generally requires CISPR 11, Class B), this equipment may not provide adequate protection against emission from the radio frequency communications. The user may need to take some mitigation measures, such as change of the device location or orientation.

Table 7 Resistance to electromagnetic interferences

Manufacturer's recommendations and statement – resistance to electromagnetic interferences			
The <i>intraoral sensor</i> is intended for use in the electromagnetic environment described below. The customer or device user must ensure that the device is used in such environment.			
Immunity tests	Test level according to IEC 60601	Compliance level	Electromagnetic environment – recommendations
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV direct contact discharge ± 15 kV air gap discharge	± 8 kV direct contact discharge ± 15 kV air gap discharge	The floor should be wooden, concrete or covered with ceramic tiles. If the floor is covered with a synthetic material, relative humidity should be at least 30 %.
Resistance to emissive radio-frequency electromagnetic fields IEC 61000-4-3	3 V/m 80 MHz-2.7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz-2.7 GHz 80 % AM at 1 kHz	The intraoral sensor is suitable for use in the professional environment (general commercial environment, clinics and hospitals)
Resistance to fields from the nearby wireless communication equipment IEC 61000-4-3	28 V/m max. 385-5785 MHz according to Table 6	28 V/m max. 385-5785 MHz according to Table 6	Portable and transportable RF communication equipment (including peripheral devices such as antenna cables and external antennas) shall be used at least 30 cm (12 inches) away from any part of the <i>intraoral sensor</i> , including the cables specified by Dentium. Otherwise, this may lead to deterioration of the equipment performance
Electrical fast transient/bursts 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines Recurrence rate 100 kHz	± 2 kV for power supply lines ± 1 kV for input/output lines Recurrence rate 100 kHz	The mains quality should be consistent with the standard commercial or hospital environment.
Voltage surges IEC 61000-4-5	± 1 kV line-to-line ± 2 kV line-to-ground	± 1 kV line-to-line ± 2 kV line-to-ground	The mains quality should be consistent with the standard commercial or hospital environment.
Immunity to	3 V	3 V	Intensity of the RF field in the

conducted interference caused by RF fields IEC 61000-4-6	0.15-80 MHz 6 V in the frequency ranges for the industrial, scientific and medical equipment (ISM) between 0.15 and 80 MHz 80 % AM at 1 kHz Power supply lines and input/output lines	0.15-80 MHz 6 V in the ISM ranges between 0.15 and 80 MHz 80 % AM at 1 kHz Power supply lines and input/output lines	frequency range 150 kHz~80 MHz should be below 3 V. In case of abnormal functioning, additional measures may be required, such as change of the intraoral sensor's orientation or location.
Voltage slumps, short interruptions and voltage fluctuations in the incoming power supply lines IEC 61000-4-11	0 % UT: 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles	0 % UT: 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles	The mains quality should be consistent with the standard commercial or hospital environment. If the user needs continuous operation during power interruption, it is recommended to use an uninterruptible power source (UPS) or a battery for device power supply.

12. Explanation of the symbols used in the labeling of consumer and transport package:

Table 8 Symbols and meanings

Symbol	Description
	Read user's instructions, manual/booklet
	Type B applied part in accordance with IEC 60601-1.

	Read instructions for use for important information, namely warnings and cautions.
	Caution! Read the instructions for use
IP68	Dust-proof housing. Protection against the consequences of temporary immersion into water. Use of the customer's specification and specific tests may be required.
	Date of manufacture of the medical device.
	Medical device manufacturer
EC REP	EU Authorized Representative
REF	Catalog number
SN	Serial number
CE 0068	Compliance with European standards
R_x Only	Use only as prescribed
	Correct position of the package with this side up
	Fragile. Handle with care
	Handle the package with care
	Keep away from moisture
	Temperature range
	Do not dispose of with household wastes

13. Commissioning, maintenance and repair

The intraoral sensor requires not only proper use, but also visual inspection and regular periodic monitoring before use. These precautions enable the device to function properly, safely and efficiently. Prior to using the device, it is necessary to check it for physical damages and defects. If any malfunction is suspected, contact the manufacturer's authorized representative in the Russian Federation.

- 1) Check the cable for detachments, keep it twisted into a ring.
- 2) In case of any malfunction, such as heat generation or an error of device recognition during use, immediately stop device operation and contact the manufacturer's authorized representative in the Russian Federation.
- 3) Use the device in a place not exposed to high temperatures, ventilation, light, dust, salts, etc.
- 4) Do not keep the device at an angle and do not expose it to shocks.
- 5) Do not use in places where chemicals or gases are stored.
- 6) When cleaning the device, lightly wipe its surface with a dry disposable towel.

In case of disorders contact the authorized representative in the Russian Federation.

Table 9 Control prior to use

Period	Checklist
Daily	• Are the sensor and the PC connected appropriately? • Are there any damages on the surface of the sensor module? • Are the sensor surface and cover clean?
Monthly	• Are there any damages to cable insulation? • Are there any damages to the cable connector (bushing)? • Are there any damages to the USB connector which is connected to the PC?
Annually	• Are there permanent defects in the image? • Is there any difference in the image brightness at identical exposure to X-rays?

14. Transportation, storage and operation conditions

Storage. The devices shall be stored in well-ventilated areas away from heat sources, sunlight, direct exposure to light and identical electrical equipment.

Transportation. Transportation of the device packaged into a transport container shall be carried out by all types of covered vehicles in compliance with the precautions indicated on the transport container and in accordance with the rules of cargo transportation applicable to that type of transport.

Operation. The devices are used in medical institutions, in the conditions requiring protection of the hands of a medical professional, protective clothing and thyroid protection for the patient when using the device according to its intended purpose.

Ambient temperature	From 10 °C to 30 °C (operation) From -20 °C to 60 °C (transportation and storage)
Relative humidity	30-95 % (operation) 10-95 % (transportation and storage)
Atmospheric pressure	700-1060Pa (operation, transportation and storage)

15. Device disinfection and cleaning

 Warning	Always clean the electronic parts after power shutdown. Do not use compressed air for cleaning.
--	---

The intraoral sensor is used with a disposable case, and therefore does not undergo disinfection, only cleaning.

Cleaning:

1. Ethyl or isopropyl alcohol (70~90 %) should be used as a detergent.
2. Clean the surface thoroughly with a clean, lint-free cloth and detergent. Clean until visible contaminations are removed. Wipe with a clean, lint-free, disposable cloth until all detergent residues is removed. Wipe it thoroughly.
3. Dry in the air or manually with a disposable lint-free cloth (do not use mechanical drying.)
4. Do not reuse a disposable cloth.
5. Change the solution for each cleaning or if there is visible contamination.

The intraoral sensor meets IP68 protection requirements, but if it is exposed to alcohol for a long time, this may lead to device malfunction and damage. For this reason, it is necessary to dry the device as quickly as possible after wiping.

 Warning	Do not clean the silicone case as described above. Do not use hard brushes and scrapers, concentrated acids and alkalis.
--	---

16. Disposal

When disposing of the device and its components, it is necessary to comply with the requirements of the instructions applicable to the Russian Federation.

Medical devices shall be disposed of in accordance with the medical procedures adopted by the medical institutions, and the local, state, or federal regulations.

Medical device "Intraoral dental digital sensor, versions: IOX1, IOX2" cannot be disposed together with solid household wastes and shall be disposed separately as industrial wastes in accordance with the federal and regional legislation. Medical device that needs to be disposed of is transferred to a specialized company holding a license for such type of activity.

Disposable hygienic covers, depending on the degree of their epidemiological, toxicological and radiation hazards, are classified as epidemiologically hazardous medical wastes: materials and tools, items contaminated with blood and/or other biological fluids.

After use, the covers are placed into special bags for medical wastes, which differ in color depending on the waste hazard class (e.g., yellow). After filling of no more than three-quarters of the volume, the bags are tied with braces and transported to a place of temporary storage. Transportation of wastes from medical institutions to the place of disposal is carried out by the specialized companies.

Consumer package, defective devices and unused devices with the expired shelf life are classified as epidemiologically safe wastes, similar to solid household waste (SHW) in composition.

The devices shall be disposed of in accordance with the requirements of the legislation and regulatory documents on the disposal of the devices adopted in the country of use of the medical device.

17. Warranty period and service life

The shelf life for this type of device is not stipulated.

Warranty period of the medical device is 2 years after installation.

18. Warranty statement

The manufacturer warrants to a customer that the product is free from defects in materials and workmanship under normal conditions of use, and also guarantees service for one year after the purchase of the intraoral transducer and accessories. The product warranty is valid only if the product is used correctly by a customer. During the warranty period, the manufacturer's authorized representative in the Russian Federation undertakes to repair or replace defective parts of the product free of charge. However, if it is determined that the warranty does not apply, a customer should bear the costs for all parts, shipping and repair.

When applying for a guarantee, it is necessary to provide a guarantee obligation to an authorized representative of the manufacturer in the territory of the Russian Federation.

Free warranty service is not provided in the following cases:


- 1) Defects or damage due to accident, misuse, abnormal use, abnormal conditions, improper storage due to user neglect.
- 2) Damage and harm caused by external causes, such as collision with objects, fire, flood, sandstorm, earth, storm, lightning, earthquake, weather.
- 3) Damage and harm caused by improper testing, operation, maintenance, installation, service or adjustment not provided or approved by the manufacturer.

4) Defects or harm caused by ignoring warnings in the user manual (instructions for use).
When applying for a guarantee it is necessary to complete the following application:

Product name	Intraoral sensor	Customer information (please complete an application)	
Date of purchase		Customer	
Date of production		Address	
Serial number			
Warranty period		Contact information	

19. List of international standards that the medical device complies with.

EN ISO 13485-2016, EN ISO 14971-2012, EN 60601-1-2006, EN 60601-1-6-2010, IEC 60529-2013, EN 60601-1-2-2015, EN 62304-2006, EN IEC 62366-1-2015, EN ISO 15223-1-2016, EN 1041-2008, EN ISO 14155-2011, MEDDEV 2.4/1 Ред. 9, 2010; MEDDEV 2.7/1, Ред. 4, 2016; MEDDEV 2.12/1, Ред. 8, 2013; MEDDEV 2.12/2, Ред. 2, 2012; Directive 1993.

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This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
7, Feb. 2020 Under Law No.138.

APOSTILLE

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1. Country : Republic of Korea

This public document

2. has been signed by SEONG BAE NAMKUNG

3. acting in the capacity of Notary Public

4. bears the seal/stamp of THE EAST SUWON LAW CORPORATION

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