BE IT KNOWN that I, GURNAM SINGH SHERGILL, of the County of Middlesex, England a duly authorised Notary Public.

DO HEREBY CERTIFY:


THAT Brian Michael Howes, a United Kingdom Citizen, of 28 Woodthorpe Road, Ashford, Middlesex TW15 2RH, holder of United Kingdom Passport No. 530973145, who is well known to me, has represented to me that he is duly authorised by Reckitt Benckiser Healthcare (UK) Limited ("the Company"), to represent them in this matter has today caused the annexed:

- Maintenance Documentation of the Medical Device regarding Context® Romantic Love natural rubber latex condoms
- Instruction for Consumers document regarding Context® Romantic Love natural latex condoms
- Technical Documentation of the Medical Device regarding Context® Romantic Love natural rubber latex condoms

to be produced to me and that he has represented to me on behalf of the Company that the annexed documents are electronic copies for use in Russia.

SIGNED & SEALED at 22 Bath Road, Hounslow, Middlesex TW3 3EB this 10th August 2021.

Protocol Number: 5507/21

Gurnam Singh Shergill

Notary Public

Gurnam Singh Shergill, N.A. (HONS)
Notary Public
SHERGILL & Co.
22 Bath Road, Hounslow, Middx
TW3 3EB
United Kingdom Tel: 020 8570 2323

approved by:

Andrew Brownhill

Global Regulatory Affairs Manager

Date: 21-Jul-2021 | 15:55 BST

Signature:

DocuSigned by:

Andrew Brownhill

Signer Name: Andrew Brownhill

Signing Reason: I approve this document

Signing Time: 21-Jul-2021 | 15:54 BST

AE0016585B70414C98AB6582059DBFF4

MAINTENANCE DOCUMENTATION OF THE MEDICAL DEVICE

Contex® Romantic Love natural rubber latex condoms

Reckitt Benckiser Healthcare (UK) Ltd., UK

NAME OF THE MEDICAL DEVICE

Contex® Romantic Love natural rubber latex condoms

MANUFACTURER

Reckitt Benckiser Healthcare (UK) Ltd., 103-105 Bath Road, Slough, Berkshire, SL1 3 UH, UK.

MANUFACTURING SITES

SSL Manufacturing (Thailand) Ltd., Wellgrow Industrial Estate 96 & 100, Moo 5, Bangna-Trad Road KM.36, Bangsamak, Bangpakong, Chachoengsao, 24180, Thailand.

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER

Reckitt Benckiser Healthcare LLC, Floor 3, 4 Shlyuzavaya emb., Moscow, 115114, Russia.

Phone: +7 (495) 961 25 65; fax: +7 (495) 961 25 66

PURPOSE OF THE MEDICAL DEVICE

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Medical Device is intended to be used by laypersons.

MEDICAL DEVICE DESCRIPTION

Condom Contex® Romantic Love - smooth, straight walled, flavored, colored, teat ended, in silicone lubricant. Made from natural rubber latex.

TECHNICAL CHARACTERISTICS

Technical parameter	Specification
Condom length, mm	180 ± 10
Condom width, mm (measured at 30 mm from the bead)	52 ± 2
Teat length, mm	Min 11
Thickness, mm	0.070+0.004/-0.008
Bead thickness, mm	1-2
Color	Red
Condom weight (after removed from the primary package), g	1.6 ± 0.5
Elongation, %	Median 700 Min
Elongation, % (aged)	Median 600 Min
Freedom from holes	Impermeable
Tensile strength, mPa	Median 20 Min
Tensile strength, mPa (aged)	Median 17 Min
Color retention	No signs of staining on absorbent paper

INDICATIONS

- prevention of unwanted pregnancy;
- protection from HIV and other sexually transmitted diseases.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

PRECAUTIONS FOR USE

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections. Check the expiry date on the condom wrapper before you use it. If the individual container is obviously damaged, discard that condom and use a new one from an undamaged package.

Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil), some pessaries and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist:

- about the compatibility of condoms with other contraceptives
- if you are using topical medicines applied to penis or vagina;
- if irritation/discomfort persists after use;
- if you are unsure if you are allergic to latex;
- if a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

For anal sex, or if necessary, use only additional, recommended lubricants on the outside of the condom.

Always read the instructions for use on the pack of the lubricant before use.

Condoms may present a choking hazard. Keep out of reach of children. If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows.

Use a condom only once. Reuse may increase risk of failure or infection.

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which might cause allergic reactions including anaphylactic shock if user is allergic to latex.

APPLICATION

Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.

Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

Dispose in a bin. Do not flush the condom down the toilet.

MEDICINAL OR PHARMACEUTICAL SUBSTANCE CONTAINED IN THE MEDICAL DEVICE

Not applicable.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

Compounded natural rubber latex contains casein (0.02% ± 10% of the contents of it, functions as a stabilizer) that is considered as a derivative from tissue/materials of animal origin as it is extracted from milk.

Casein (milk protein) is a natural by-product of milk. It has no known toxic effects although it may cause a reaction in individuals with an allergy to milk. The exposure from the condoms is considered not to pose any additional risk to users of the condoms.

The device does not contain any human blood derivatives.

STERILIZATION

The medical device is presented as a non-sterile device and is not intended to be sterilized thereafter.

USAGE WITH OTHER MEDICAL DEVICES

Condoms can be used on their own to prevent from STIs and pregnancy and aren't intended to be used in conjunction with other medical devices.

However, consumers can use condoms with condom compatible lubricants if needed.

INSTALLATION, MAINTENANCE, CALIBRATION

Not applicable.

CLEANING, DISINFECTION, STERILISATION OF THE MEDICAL DEVICE

Not applicable.

INTERNATIONAL STANDARDS

ISO 13485

ISO 14971

ISO 4074

PRECAUTIONS TO BE TAKEN BY USER RELATED TO THE DISPOSAL OF THE DEVICE

Contex® Romantic Love natural rubber latex condoms should be disposed in a bin. The condoms should not be flushed the condom down the toilet.

DATE OF THE ISSUE / THE LATEST REVISION OF THE MAINTENANCE DOCUMENTATION

The medical device includes an indication of the latest revision within the instruction for use.

GUIDELIENSS FOR DESTRUCTION AND DISPOSAL

The devices should be used in accordance with the instructions for use.

In household use, the product is disposed of as household waste. Do not put it down the drain.

In all other cases, the condoms are disposed of in accordance with local requirements.

ENVIRONMENT SAFETY REQUIREMENTS

With proper usage, storage and disposal Contex® Romantic Love natural rubber latex condoms have no negative impact on the environment.

Approved by:

Andrew Brownhill
Global Regulatory Affairs Manager

Date: 21-Jul-2021 | 15:55 BST

Signature:

DocuSigned by:

Andrew Brownhill

Signer Name: Andrew Brownhill

Signing Reason: I approve this document

Signing Time: 21-Jul-2021 | 15:55 BST

AE0016585B70414C98AB6582059DBFF4

INSTRUCTION for CONSUMERS

Contex® Romantic Love natural rubber latex condoms

Reckitt Benckiser Healthcare (UK) Ltd., UK

2.

3.

4.

Read the instruction carefully before use, especially if condoms are used for oral or anal sex.

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections.

ATTENTION: Check the expiry date on the condom wrapper before you use it. If the individual container is obviously damaged, discard that condom and use a new one from an undamaged package.

PURPOSE

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases.

APPLICATION

1. Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.

2. Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis.

3. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

4. Dispose in a bin. Do not flush the condom down the toilet.

Only use lubricants recommended for use with condoms.

For anal sex, or if necessary, use only additional, recommended lubricants on the outside of the condom. Always read the instructions for use on the pack of the lubricant before use.

Condoms may present a choking hazard. Keep out of reach of children.

If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows.

Use a condom only once. Reuse may increase risk of failure or infection.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which might cause allergic reactions including anaphylactic shock if user is allergic to latex.

PRECAUTIONS FOR USE

Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil), some pessaries and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist:

about the compatibility of condoms with other contraceptives

- if you are using topical medicines applied to penis or vagina;
- if irritation/discomfort persists after use;
- if you are unsure if you are allergic to latex;
- if a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

Date of the latest revision of the instruction:

Approved by:

Andrew Brownhill

Global Regulatory Affairs Manager

Date: 21-Jul-2021 | 15:55 BST

Signature:

DocuSigned by:

Andrew Brownhill

Signer Name: Andrew Brownhill

Signing Reason: I approve this document

Signing Time: 21-Jul-2021 | 15:55 BST

AE0016585B70414C98AB6582059DBFF4

TECHNICAL DOCUMENTATION OF THE MEDICAL DEVICE

Contex® Romantic Love natural rubber latex condoms

Reckitt Benckiser Healthcare (UK) Ltd., UK

NAME OF THE MEDICAL DEVICE

Contex® Romantic Love natural rubber latex condoms

PURPOSE OF THE MEDICAL DEVICE

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Medical Device is intended to be used by laypersons.

MEDICAL DEVICE POTENTIAL RISK CLASS: 2b**STERILITY**

The device is non-sterile.

MEDICAL DEVICE DESCRIPTION

Contex® Romantic Love natural rubber latex condoms - smooth, straight walled, flavored, colored, teat ended, in silicone lubricant. Made from natural latex.

TECHNICAL CHARACTERISTICS

Technical parameter	Specification	Test Method used by the manufacturer	Frequency
Condom length, mm	180 ± 10	ISO 4074:2015 clause 9.3.1 Annex D	Every batch
Condom width, mm (measured at 30 mm from the bead)	52 ± 2	ISO 4074:2015 clause 9.3.2 Annex E	Every batch
Teat length, mm	Min. 11	Appendix 1	Once a month
Thickness, mm	0.070+0.004/-0.008	ISO 4074:2015 clause 9.3.3 Annex F	Every batch
Bead thickness, mm	1-2	Appendix 2	Every batch
Color	Red	Visual	Every batch
Condom weight (after removed from the primary package, W2), g	1.6 ± 0.5	Appendix 3	Once a month
Elongation, %	Median 700 Min	ISO 4074:2015 Annex J	Once a month
Elongation, % (aged)	Median 600 Min	ISO 4074:2015 Annex J	Once a month
Tensile strength, mPa	Median 20 Min	ISO 4074:2015 Annex J	Once a month
Tensile strength, mPa (aged)	Median 17 Min	ISO 4074:2015 Annex J	Once a month
Freedom from holes	Impermeable	ISO 4074:2015 clause 12 Annex M	Every batch
Color retention	No signs of staining on absorbent paper	Appendix 4	Once a month

RAW MATERIALS (% material is based on theoretical value)

Materials	Nominal content of the material, %+/-10%	Trade name
Basic material:		
Compounded Natural Rubber latex	76.41	NR Latex Concentrate (CL60 HA)
Red pigment	0.52	KROMACRYL 4393 RED FDA
Dusting material:		
Calcium Carbonate	2.41	Calcium carbonate light Bluestar
Magnesium Carbonate	0.15	Magnesium Carbonate Extra Light
Lubricant:		
Silicone oil	20.10	BRB Silicone Oil DM 350 EP*/ Angeläy® 350CS**
Sodium saccharin	0.02	Saccharin, sodium salt
Flavoring agent strawberry	0.39	Strawberry Flavor

Raw Material	Additional information
Compounded Natural Rubber latex	NR Latex Concentrate (CL60 HA) is a trade name of natural latex raw material used in compounded NRL
Red pigment	KROMACRYL 4393 RED FDA is a trade name of red pigment
Silicone	*BRB Silicone Oil DM 350 EP and **Angeläy® 350CS are trade names of the silicone
Sodium saccharin	Saccharin, sodium salt - is a trade name of sodium saccharin
Flavouring agent strawberry	Strawberry Flavor is a trade name of flavouring agent strawberry

PACKAGING

The condom is packaged in an individual package of multi-layer material (foil).

PACKAGE DIMENSIONS:

Individual foil: 60±2 mm x 60± 2 mm.

Primary package material:

Front foil layer:
PET/AL/LDPE

Back foil layer:
PET/AL/LDPE

From the outer to the inner foil layer:

PET = Polyethylene terephthalate (Polyethyl benzene-1,4-dicarboxylate)

AL = Aluminium

LDPE = Low Density Polyethylene

PRESENTATION OF THE MEDICAL DEVICE

Secondary packaging is a carton pack, wrapped with the transparent tape where applicable. The number of condoms is indicated on the package.

LABELLING INFORMATION

Each individual package is labelled with:

- trademark;
- lot number (LOT);
- Expiry date (EXP).

Symbols used for labelling and what they mean:

	Expiry date
	Lot number

Additional information can be printed on the individual package including reference to international standards.

The following information is applied to the outer surface of the secondary/consumer package:

- condom description;
- quantity of condoms in consumer package;
- nominal condom length;
- nominal condom width;
- indication that condom contains lubricant;
- indication that the lubricant does not contain spermicide;
- storage guidelines;
- name and address of the manufacturer;
- country of origin (if applicable);
- shelf life formatted as "Expiry date".

Additional information can be printed on the outer surface of the secondary/ consumer package including reference to international standards.

Outer or inner surface of consumer packages and/or package insert contains prescribing information, which includes:

- reminder to handle the condoms carefully, including when removing them out of the package;
- indications, contraindications, precautions for use, possible side effects;
- method of application;
- necessity to remove condom right after sexual intercourse;
- necessity to exclude petroleum- and oil-based lubricants;
- necessity to use safe additional lubricant if needed;
- necessity to consult a specialist on compatibility with topical medications;

- disposal instructions for used condoms;
- condoms are disposable products.

INDICATIONS

- prevention of unwanted pregnancy;
- protection from HIV and other sexually transmitted diseases.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

PRECAUTIONS FOR USE

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections. Check the expiry date on the condom wrapper before you use it. If the individual container is obviously damaged, discard that condom and use a new one from an undamaged package.

Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil and some pessaries) and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist:

- about the compatibility of condoms with other contraceptives
- if you are using topical medicines applied to penis or vagina;
- if irritation/discomfort persists after use;
- if you are unsure if you are allergic to latex;
- if a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

Only use lubricants recommended for use with condoms.

For anal sex, or if it is necessary, use additional lubrication on the outside of the condom. Always read the instructions for use on the pack of the lubricant before use.

Condoms may present a choking hazard. Keep out of reach of children. If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows.

Use a condom only once. Reuse may increase risk of failure or infection.

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which might cause allergic reactions including anaphylactic shock if user is allergic to latex.

APPLICATION

Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.

Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

Dispose in a bin. Do not flush the condom down the toilet.

TRANSPORTATION

There are no special conditions required for transportation.

EXPLOITATION CONDITIONS

There are no special conditions required when the medical device is used for its purpose and as indicated in the labelling.

STORAGE AND SHELF LIFE

Store the product in a dry place at 0-25 °C. Avoid direct sunlight.

Shelf life - 5 years.

RISK MANAGEMENT SUMMARY

In accordance with the BS EN ISO 14971, Medical devices - Application of risk management to medical devices, the international standard for medical device risk management an assessment of risk relating to the use of medical device has been conducted.

INFORMATION ABOUT BIOCOMPATIBILITY

Test	Standard	Biological effect
In vitro cytotoxicity testing	ISO 10993-5	Toxic at 1:4 dilution
Sensitization	ISO 10993-10	The test articles showed no evidence of causing delayed dermal contact sensitization
Vaginal Irritation	ISO 10993-10	Test articles was considered a non-irritant to vaginal tissue
Penile Irritation	ISO 10993-10	Test article extract was considered a non-irritant to penile tissue
Rectal Irritation	ISO 10993-10	Test article extract was considered a non-irritant to rectal tissue
Acute systemic toxicity	ISO 10993-11	There was no mortality or evidence of systemic toxicity from test article

*ISO 10993-5 suggests that a response grade of moderate or severe should be considered a cytotoxic effect. The typical condom result has been found not to be associated with any problem in the regular use of condoms by male or the female partner. If the results of other testing (typically sensitization and irritation tests) and associated evaluations are seen as acceptable then the cytotoxicity result should be considered in context and the condom considered to meet the requirements of ISO 10993. The results indicate that Context® Romantic Love natural rubber latex can be considered to be potentially biologically safe.

MANUFACTURING PROCESS

See Appendix 5

MEDICINAL OR PHARMACEUTICAL SUBSTANCE CONTAINED IN THE MEDICAL DEVICE
Not applicable.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

Compounded natural rubber latex contains casein (0.02% $\pm$ 10%) of the contents of it, functions as a stabilizer) that is considered as a derivative from tissue/materials of animal origin as it is extracted from milk. Casein (milk protein) is a natural by-product of milk. It has no known toxic effects although it may cause a reaction in individuals with an allergy to milk. The exposure from the condoms is considered not to pose any additional risk to users of the condoms.

The device does not contain any human blood derivatives.

STERILIZATION

The medical device is presented as a non-sterile device and is not intended to be sterilized thereafter.

SOFTWARE DESIGN, DEVELOPMENT AND VALIDATION

Not applicable.

MAINTENANCE AND REPAIR

Not applicable, as the product should be disposed of after a single use by a user.

GUIDELIENSS FOR DESTRUCTION AND DISPOSAL

The devices should be used in accordance with the instructions for use.

In household use, the product is disposed of as household waste. Do not put it down the drain.

In all other cases, the condoms are disposed of in accordance with local requirements.

ENVIRONMENT SAFETY REQUIREMENTS

With proper usage, storage and disposal Natural latex condoms Contex® Romantic Love have no negative impact on the environment.

Test Method used by the manufacturer

Appendix 1: Test method for condom teat length determination (D8291875)

- Sample preparation:
 - Dipped condom
Dry dipped condom before rolling will be used for this test. Stretch the condom slightly twice to smooth out the wrinkles.
 - Condom in container
Move the condom inside the individual container such that it is away from the area where the container is to be torn. Tear the container and remove the condom. Do not use scissors or other sharp instruments to open the package.

Unroll the condom, stretch it slightly twice to smooth out the wrinkles. Lubricants may be removed by washing with a suitable solvent. If a solvent is used to remove the lubricant an adequate period of drying should be allowed before testing.
- Put the condom over the mandrel, having a rounded end with a diameter of (25 $\pm$ 2) mm with a scale divided into millimetres and dimensions with the zero beginning at the rounded end, and let it hang freely, stretched only by its own mass.
- Mark a point at the junction of condom teat and the shoulder of the condom. See Figure 1 and Picture 1.



Figure 1. Length of condom teat

Picture 1. A point marked at the junction of condom teat and the shoulder.

- Remove the marked condom out from the mandrel. Lay it flat and put the ruler on top and measure the distance between the marked point and teat end by using a calibrated ruler and record the result.

Appendix 2: Test method for measurement of condom bead thickness (LAB/LT-DMS/16, D8291143)

- Take a condom out of the pack or use a dried condom and spread it out gently.
- Measure the thickness of the bead by thickness gauge with the nearest reading of 0.01 mm.
- Repeat the measurement another 2 times.
- Average condom bead thickness for individual condom will be calculated.

Appendix 3: Test method for condom weight (after removed from the primary package) (D8291956)

- Weigh individual foiled condom and record the weight (W1).
- Weigh a dry and clean petri dish and tare zero.
- Move the condom inside the individual package such that it is away from the area where the package is to be torn. Tear the package. Do not use scissors or other sharp instruments to open the package.
- Remove the condom carefully out of the opened foil package by using dry and clean forceps and place it on the petri dish. Record the weight of condom (W2).
- Ensure that petri dish and forceps are dry and clean before testing next sample.
- Weight of foil package and remaining lubricant and powder after condom removal (W3) = W1 - W2


Appendix 4: Test method for colour retention (QSWI 8.2-5.0117)

- Completely wet the inside and outside of the condom with deionized water.
- Wrap the wet condom in the white absorbent paper so that the largest possible surface area of the condom is in contact with the paper and wet the ends of the wrapper and bend over ends. Seal the whole in a container to prevent loss of moisture.

- Allow the container and its content to stand for 16 - 24 hours at room temp.
- After removing the absorbent paper form the container, examine it visually for any indicator of staining.
- Repeat the test for other 12 condoms (N = 13). Accept/ Reject criteria is 0:1.

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