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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of YANTAI ADDCARE BIO-TECH LIMITED COMPANY on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized

Signature:

Lv Yan

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(Date: Aug. 14, 2020)

证明书 网站 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

11
Yantai Addcare Bio-Tech Limited Company, China

Ms. Carolyn Wang 王利华

Supervisor International Department

Disposable tips MilaLab-TP

Disposable non-sterile tips MilaLab-TP
for use in a fully automated liquid handling system MilaLab-ASP/8

Instructions for use

Document	Language variants:	Version	Dated
1.3.04.0070-IU	Russian (only for the Russian market)	0	2020-06-24
-	English translation		

Instructions for use

Disposable non-sterile tips MilaLab-TP for use in a fully automated liquid handling system MilaLab-ASP/8

1. Medical device information

Medical device name: Disposable non-sterile tips MilaLab-TP for use in a fully automated liquid handling system MilaLab-ASP/8.

The medical device is manufactured in two various:

MilaLab-TP, 300 µl and MilaLab-TP, 800 µl.

The instructions for use covers to the following versions of a medical device for in vitro diagnostics:

MilaLab-TP, 300 µl and MilaLab-TP, 800 µl.

Medical device codes are listed in Table 1.

Table 1

Medical device code	Title	Tips per tray, pcs	Tips per group package, pcs	Tips per transport package, pcs
1.3.04.0071	MilaLab-TP, 300 µl	96	480	9600
1.3.04.0072	MilaLab-TP, 800 µl	96	480	7200

2. Intended use

For sample and reagent dosage in a fully automated liquid handling system MilaLab-ASP/8 by the trained health care professionals in clinical laboratory conditions.

3. Indications/Contraindications

For sample and reagent dosage in a fully automated liquid handling system MilaLab-ASP/8.

For in vitro diagnostic use.

For professional use.

For use by the trained health care professionals in clinical laboratory conditions.

There are no specific contraindications for use.

4. Warnings and precautions

Non-sterile IVD medical device.

Disposable device. Do not re-use.

Do not use the medical device without labelling or in case labelling/package is damaged..

Do not use after expiration date.

Do not use with other robotic stations.

Follow general laboratory precautions, wear protective clothing and protective equipment (gloves).

Low level of personal and environmental risk.

5. Principal technical characteristics of medical device

Table 2

Potential risk class: 1	MilaLab-TP, 300 µl	MilaLab-TP, 800 µl
Colour	white	white
Sterility	non-sterile	non-sterile
Maximum volume, µl	300	800
External diameter, mm	7.80 (±5%)	7.80 (±5%)
Outlet, mm	0.65 (±5%)	0.65 (±5%)
Length, mm	56.0 (±5%)	78.0 (±5%)
Weight, g	0.49 (±10%)	0.66 (±10%)

6. Disposal medical device

The medical device must be disposed of as potentially hazardous medical waste. In accordance with Sanitary Rules and Regulations (SanPiN) 2.1.7.2790-10 «Sanitary and epidemiological requirements to treatment of medical wastes», this type of waste is classified as class B (epidemiologically hazardous waste).

The system of collection, temporary storage and transportation of medical wastes must be organized in accordance with SanPiN 2.1.7.2790-10 «Sanitary and epidemiological requirements to treatment of medical wastes».

Collection, temporary storage and disposal of wastes must be carried out in accordance with the procedure of medical waste management adopted by the organization carrying out medical and/or pharmaceutical activities.

Warning: single-use medical device! Do not re-use! Dispose of in a medical waste container! Do not dispose of used medical device together with household waste.

After the expiration date, unused devices can be disposed of as non-hazardous medical waste of class A (epidemiological safe waste, close in composition to solid domestic waste).

SanPiN 2.1.7.2790-10 as well as local/national/federal requirements to the disposal of medical waste must be followed and general precautions must be taken.

7. The order of application medical device

Open the transport package and take the group package with tips. Open the group package with the tips. Put the required quantity of trays with the tips into the specialized compartment of a fully automated liquid handling system MilaLab-ASP/8.

Further actions are carried out inside used fully automated liquid handling system MilaLab-ASP/8 and described in appropriate Operating manual (User manual) of it.

Conditions of Use: at 5-30 °C, humidity of 40-60 %. Atmosphere pressure is not specified.

8. Interaction with other medical devices

There are no adverse effects when the device is used in combination with fully automated liquid handling system MilaLab-ASP/8, provided it is used as intended in accordance with the instruction for use/user manual.

Not intended for use with other medical devices.

9. Shelf life. Storage and transportation medical devices

Shelf life is 3 years. Please check the manufacturing date before use.

To be stored at 5-30°C, humidity of 40-60 % rH.

Atmosphere pressure is not specified.

Can be transported by all types of enclosed transport at 5-30°C, humidity of 40-60 %.

Atmosphere pressure is not specified.

10. Requirements for maintenance and repair medical devices

There are no parts, which need service, in the product and it is not subject to maintenance and repair.

11. Guarantee

This product is designed, manufactured, tested and packaged in accordance with all applicable requirements. Manufacturer guarantees conformity of the product to the requirements of regulatory and technical documentation. Quality and safety of the product is guaranteed within established shelf life in accordance with storage/transportation conditions. Do not use the product without labelling or in case labeling/package is damaged. Complaints should be sent to the authorized representative of the manufacturer in the Russian Federation.


Responsible manufacturer:

Yantai Addcare Bio-Tech Limited Company
No.7, Hanjiang Road, Development Zone, Yantai
264006, Shandong Province, China
Phone: +86 535 6398536
Fax: +86 535 6398506
Email: infor@adcare.cn
Website: www.addcare.cn/en/

Authorized representative of the manufacturer in the the Russian Federation:

Research and production company "Diagnostic Systems",
Ltd. (RPC "Diagnostic Systems", Ltd.)
Legal address: 47 Komintern Street, 603014, Nizhny
Novgorod, Russia
The address to send complaints: 22 Yablonevaya Street,
603093, Nizhny Novgorod, Russia
Phone: +7 (831) 434-86-83, 8 (800) 555-03-00
Email: info@npods.ru, sos@npods.ru
Website: www.npods.ru

12. Key to symbols used

	Manufacturer
	Date of manufacture
	Expiry date...
	Batch code
	Catalog number
	Non-sterile
	Do not use if package is damaged
	Storage temperature limitation
	Can not be reused
	Consult Instructions for use
	For In vitro diagnostic use
	Quantity (in a package)
	Fragile, handle with care
	Up
	Keep dry
	No stepping

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