BE IT KNOWN that I Laura Anne Pawley, Mezzanine Floor, Epsom Square, Epsom, Surrey, KT19 8AG, a duly authorised Notary Public, duly admitted, sworn and practice throughout England and Wales.

CERTIFY ONLY that Katie Pearce, a United Kingdom Citizen, who is well known to me, who, is duly authorised by Samantha Neilson of Reckitt Benckiser Healthcare (UK) Limited 'the Company' to represent them in this matter, has today caused the annexed Technical Documentation of the Medical Device Contex Textured Natural Rubber Latex Condoms for use in the Russian Federation to be produced to me and that they have represented to me on behalf of the Company that the said document is a true electronic copy of the original.

SIGNED and sealed at Mezzanine Floor, Epsom Square, Epsom, Surrey, KT19 8AG on the 27th January 2022.

Protocol No: 2022-537

Laura Anne Pawley

Laura Anne Pawley
Notary Public

Dear Sir / Madam

I Samantha Neilson of Reckitt Benckiser Healthcare (UK) Limited CERTIFY that the attached
**Technical Documentation of the Medical Device for - Context Textured Natural rubber Latex
Condoms -.**

Name: Samantha Neilson

Signature: *Samantha Neilson*

Electronically signed by:
Samantha Neilson
Reason: I approve this
document
Date: Jan 25, 2022 14:54
GMT

Date: 25-Jan-2022

Registered office at 103 - 105 Bath Road, Slough, Berkshire, SL1 3UH. Registered in England & Wales, No 6270876

Approved by:

Andrew Brownhill
Global Regulatory Affairs Manager
Date:
Signature:

Andrew Brownhill

Electronically
signed by:
Andrew
Brownhill
Reason: I
approve this
document.
Date: Jan 21,
2022 09:17
GMT

TECHNICAL DOCUMENTATION OF THE MEDICAL DEVICE

Contex® textured natural rubber latex condoms

Reckitt Benckiser Healthcare (UK) Ltd., UK

NAME OF THE MEDICAL DEVICE

Contex® textured natural rubber latex condoms

VARIANTS

- 1) Dotted – with dots
- 2) Ribbed– with ribs.

PURPOSE OF THE MEDICAL DEVICE

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Medical Device is intended to be used by laypersons.

MEDICAL DEVICE POTENTIAL RISK CLASS: 2b

STERILITY

The device is non-sterile.

MEDICAL DEVICE VARIANTS AND DESCRIPTIONS

Condoms Contex Ribbed, with ribs, straight walled, teat ended, in silicone lubricant. Made from natural latex

Condoms Contex Dotted, with dots, straight walled, teat ended, in silicone lubricant. Made from natural latex

TECHNICAL CHARACTERISTICS

Contex Ribbed			
Technical parameter	Specification	Test Method used by the manufacturer	Frequency
Condom length, mm	180 ± 10	ISO 4074:2015 clause 9.3.1 Annex D	Every batch
Condom width, mm (measured at 30 mm from the bead)	53 ± 2	ISO 4074:2015 clause 9.3.2 Annex E	Every batch
Teat length, mm	Min 11	Appendix 1	Once a month
Thickness, mm	Not applicable to textured condoms	ISO 4074:2015 clause 9.3.3 Annex F	Every batch
Bead thickness, mm	1-2	Appendix 2	Every batch
Color	Coral	Visual	Every batch
Condom weight (after removed from the primary package, W2), g	2.0 ± 0.5	Appendix 3	Once a month
Elongation, %	Median 700 Min	ISO 4074:2015 Annex J	Once a month
Elongation, % (aged)	Median 600 Min	ISO 4074:2015 Annex J	Once a month

Tensile strength, mPa	Median 20 Min	ISO 4074:2015 Annex J	Once a month
Tensile strength, mPa (aged)	Median 17 Min	ISO 4074:2015 Annex J	Once a month
Freedom from holes	Impermeable	ISO 4074:2015 clause 12 Annex M	Every batch
Color retention	No signs of staining on absorbent paper	Appendix 4	Once a month

Contex Dotted			
Technical parameter	Specification	Test Method used by the manufacturer	Frequency
Condom length, mm	180 ± 10	ISO 4074:2015 clause 9.3.1 Annex D	Every batch
Condom width, mm (measured at 30 mm from the bead)	52 ± 2	ISO 4074:2015 clause 9.3.2 Annex E	Every batch
Teat length, mm	Min 11	Appendix 1	Once a month
Thickness, mm	Not applicable to textured condoms	ISO 4074:2015 clause 9.3.3 Annex F	Every batch
Bead thickness, mm	1-2	Appendix 2	Every batch
Condom weight (after removed from the primary package, W2), g	1.8 ± 0.5	Appendix 3	Once a month
Elongation, %	Median 700 Min	ISO 4074:2015 Annex J	Once a month
Elongation, % (aged)	Median 600 Min	ISO 4074:2015 Annex J	Once a month
Tensile strength, mPa	Median 20 Min	ISO 4074:2015 Annex J	Once a month
Tensile strength, mPa (aged)	Median 17 Min	ISO 4074:2015 Annex J	Once a month
Freedom from holes	Impermeable	ISO 4074:2015 clause 12 Annex M	Every batch

RAW MATERIALS (% material is based on theoretical value)

Item Name	Raw Materials	Trade name	Nominal content of the material, % +/- 10%
Contex Ribbed	Basic material: - Compounded Natural Rubber latex - red pigment	- NR Latex Concentrate (CL60 HA) - KROMACRYL 4393 RED FDA	77.71 0.02
	Dusting material: - Calcium Carbonate - Magnesium Carbonate	- Calcium carbonate light Bluestar - Magnesium Carbonate Extra Light	1.97 0.13

Item Name	Raw Materials	Trade name	Nominal content of the material, % $\pm$ 10%
	Lubricant: - Silicone oil	*BRB Silicone Oil DM 350 EP/ **Angeläy® 350CS	20.12-20.15
	- Odor masker	*Zeus Odour masking agent FLA70281/ **Transatak TAK-181256	0.02-0.05

Item Name	Raw Materials	Trade name	Nominal content of the material, % $\pm$ 10%
Context Dotted	Basic material: Compounded Natural Rubber latex	NR Latex Concentrate (CL60 HA)	78.57
	Dusting material: - Calcium Carbonate - Magnesium Carbonate	- Calcium carbonate light Bluestar - Magnesium Carbonate Extra Light	2.24 0.14
	Lubricant: - Silicone oil	*BRB Silicone Oil DM 350 EP/ **Angeläy® 350CS	19.00-19.03
	- Odor masker	*Zeus Odour masking agent FLA70281/ **Transatak TAK-181256	0.02-0.05

Raw Material	Additional information
Compounded Natural Rubber latex	NR Latex Concentrate (CL60 HA) is a trade name of natural latex raw material used in compounded NRL
Silicone	*BRB Silicone Oil DM 350 EP and **Angeläy® 350CS are trade names of the silicone
Odor masker	*Zeus Odour masking agent FLA70281 and **Transatak TAK-181256 are trade names of the odour masker

PACKAGING

The condom is packaged in an individual package of multi-layer material (foil).

PACKAGE DIMENSIONS:

Individual foil: 60 $\pm$ 2 mm x 60 $\pm$ 2 mm.

Primary package material:

Context® Ribbed	Frontal foil layer:
-----------------	---------------------

	PET/AL/LDPE Back foil layer: PET/AL/LDPE
Contex® Dotted	Frontal foil layer: PET/AL/LDPE Back foil layer: PET/AL/LDPE

From the outer to the inner foil layer:

PET = Polyethylene terephthalate (Polyethyl benzene-1,4-dicarboxylate)

AL = Aluminium

LDPE = Low Density Polyethylene

PRESENTATION OF THE MEDICAL DEVICE

Secondary packaging is a carton pack, wrapped with the transparent tape where applicable. The number of condoms is indicated on the package.

PACKAGE CONTENTS OF THE MEDICAL DEVICE

Each separate medical device variant is packed in a secondary / consumer package of 3 pcs. a pack. It is accompanied by the instruction for consumers that can be enclosed in the package or printed in the inside and / or outside of the package.

LABELLING INFORMATION

Each individual package is labelled with:

- trademark;
- lot number (LOT);
- Expiry date (EXP).

Symbols used for labelling and what they mean:

	Expiry date
	Lot number

Additional information can be printed on the individual package including reference to international standards.

The following information is applied to the outer surface of the secondary/consumer package:

- condom description;
- quantity of condoms in consumer package;
- nominal condom length;
- nominal condom width;
- indication that condom contains lubricant;
- indication that the lubricant contains spermicide;
- indication that the condom contains natural latex;
- storage guidelines;
- lot number;
- name and address of the manufacturer;
- country of origin (if applicable);
- shelf life formatted as "Expiry date".

Symbols used for labelling and what they mean:

	Use by/Expiry date
	Lot number
	Contains natural latex
	For single use only
	Caution!

Additional information can be printed on the outer surface of the secondary/ consumer package including reference to international standards.

Outer or inner surface of consumer packages and/or package insert contains prescribing information, which includes:

- reminder to handle the condoms carefully, including when removing them out of the package;
- indications, contraindications, precautions for use, possible side effects;
- method of application;
- necessity to remove condom right after sexual intercourse;
- necessity to exclude petroleum- and oil-based lubricants;
- necessity to use safe additional lubricant if needed;
- necessity to consult a specialist on compatibility with topical medications;
- disposal instructions for used condoms.
- the condom is a disposable product

INDICATIONS

- prevention of unwanted pregnancy;
- protection from HIV and other sexually transmitted diseases.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

PRECAUTIONS FOR USE

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections. Check the expiry date on the condom wrapper before you use it. If the individual container is obviously damaged, discard that condom and use a new one from an undamaged package.

Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil and some pessaries) and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist:

- about the compatibility of condoms with other contraceptives
- if you are using topical medicines applied to penis or vagina;
- if irritation/discomfort persists after use;
- if you are unsure if you are allergic to latex;
- if a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

For anal sex, if it is necessary use additional lubrication on the outside of the condom.

Condoms may present a choking hazard. Keep out of reach of children.

If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows. Use a condom only once. Reuse may increase risk of failure or infection.

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which might cause allergic reactions including anaphylactic shock if user is allergic to latex.

APPLICATION

Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.

Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

Dispose in a bin. Do not flush the condom down the toilet.

TRANSPORTATION

There are no special conditions required for transportation.

EXPLOITATION CONDITONS

There are no special conditions required when the medical device is used for its purpose and as indicated in the labelling.

STORAGE AND SHELF LIFE

Store the product in a dry place at 0-25 °C. Avoid direct sunlight.

Shelf life - 5 years.

RISK MANAGEMENT SUMMARY

In accordance with the ISO 14971, Medical devices - Application of risk management to medical devices, the international standard for medical device risk management an assessment of risk relating to the use of the medical device has been conducted.

INFORMATION ABOUT BIOCOMPATIBILITY

Test	Standard	Biological effect
In vitro cytotoxicity testing	ISO 10993-5	Toxic at 1:4 dilution*
Sensitization	ISO 10993-10	The test articles showed no evidence of causing delayed dermal contact sensitization
Vaginal Irritation	ISO 10993-10	Test articles was considered a non-irritant to vaginal tissue
Penile Irritation	ISO 10993-10	Test article extract was considered a non-irritant to penile tissue
Rectal Irritation	ISO 10993-10	Test article extract was considered a non-irritant to rectal tissue
Acute systemic toxicity	ISO 10993-11	There was no mortality or evidence of systemic toxicity from test article

*ISO 10993-5 suggests that a response grade of moderate or severe should be considered a cytotoxic effect. The typical condom result has been found not to be associated with any problem in the regular use of condoms by male or the female partner. If the results of other testing (typically sensitization and irritation tests) and associated evaluations are seen as acceptable then the cytotoxicity result should be considered in context and the condom considered to meet the requirements of ISO 10993. The results indicate that Contex® textured natural rubber latex condoms can be considered to be potentially biologically safe.

MANUFACTURING PROCESS

See Appendix 5

MEDICINAL OR PHARMACEUTICAL SUBSTANCE CONTAINED IN THE MEDICAL DEVICE
Not applicable.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

Compounded natural rubber latex contains casein (0.02% ± 10%) of the contents of it, functions as a stabilizer) that is considered as a derivative from tissue/materials of animal origin as it is extracted from milk. Casein (milk protein) is a natural by-product of milk. It has no known toxic effects although it may cause a reaction in individuals with an allergy to milk. The exposure from the condoms is considered not to pose any additional risk to users of the condoms.

The device does not contain any human blood derivatives.

STERILIZATION

The medical device is presented as a non-sterile device and is not intended to be sterilized thereafter.

SOFTWARE DESIGN, DEVELOPMENT AND VALIDATION

Not applicable.

MAINTENANCE AND REPAIR

Not applicable, as the product should be disposed of after a single use by a user.

GUIDELIENSS FOR DESTRUCTION AND DISPOSAL

The devices should be used in accordance with the instructions for use.

In household use, the product is disposed of as household waste. Do not put it down the drain.

In all other cases, the condoms are disposed of in accordance with local requirements.

ENVIRONMENT SAFETY REQUIREMENTS

With proper usage, storage and disposal Contextex® textured natural rubber latex condoms have no negative impact on the environment.

Test Method used by the manufacturer

Appendix 1: Test method for condom teat length determination (D8291875)

- Sample preparation:
 - Dipped condom
Dry dipped condom before rolling will be used for this test. Stretch the condom slightly twice to smooth out the wrinkles.
 - Condom in container
Move the condom inside the individual container such that it is away from the area where the container is to be torn. Tear the container and remove the condom. Do not use scissors or other sharp instruments to open the package.

Unroll the condom, stretch it slightly twice to smooth out the wrinkles. Lubricants may be removed by washing with a suitable solvent. If a solvent is used to remove the lubricant an adequate period of drying should be allowed before testing.
- Put the condom over the mandrel, having a rounded end with a diameter of (25 ± 2) mm with a scale divided into millimetres and dimensions with the zero beginning at the rounded end, and let it hang freely, stretched only by its own mass.
- Mark a point at the junction of condom teat and the shoulder of the condom. See Figure 1 and Picture 1.



Figure 1. Length of condom test

Picture 1. A point marked at the junction of condom test and the shoulder.

- Remove the marked condom out from the mandrel. Lay it flat and put the ruler on top and measure the distance between the marked point and test end by using a calibrated ruler and record the result.

Appendix 2: Test method for measurement of condom bead thickness (LAB/LT-DMS/16, D8387639)

- Take a condom out of the pack or use a dried condom and spread it out gently.
- Measure the thickness of the bead by thickness gauge with the nearest reading of 0.01 mm.
- Repeat the measurement another 2 times.
- Average condom bead thickness for individual condom will be calculated.

Appendix 3: Test method for condom weight (after removed from the primary package) (D8291956)

- Weigh individual foiled condom and record the weight (W1).
- Weigh a dry and clean petri dish and tare zero.
- Move the condom inside the individual package such that it is away from the area where the package is to be torn. Tear the package. Do not use scissors or other sharp instruments to open the package.
- Remove the condom carefully out of the opened foil package by using dry and clean forceps and place it on the petri dish. Record the weight of condom (W2).
- Ensure that petri dish and forceps are dry and clean before testing next sample.
- Weight of foil package and remaining lubricant and powder after condom removal (W3) = W1 - W2

Appendix 4: Test method for colour retention (D8305772, QSWI 8.2-5.0117)


- Completely wet the inside and outside of the condom with deionized water.
- Wrap the wet condom in the white absorbent paper so that the largest possible surface area of the condom is in contact with the paper and wet the ends of the wrapper and bend over ends. Seal the whole in a container to prevent loss of moisture.

- Allow the container and its content to stand for 16 - 24 hours and record sample stored time in form number D8390807
- After removing the absorbent paper from the container, examine it visually for any indicator of staining.
- Repeat the test for other 12 condoms (N = 13). Accept/ Reject criteria is 0:1.

Appendix 5

OVERVIEW PROCESS CONTROL FLOWCHART

PRODP-1 Version 1.0

proved by:

Andrew Brownhill
Global Regulatory Affairs Manager

Date:

Signature:

Andrew Brownhill

Electronically
signed by:
Andrew
Brownhill
Reason: I
approve this
document.
Date: Jan
21, 2022
09:17 GMT

MAINTENANCE DOCUMENTATION
(Instruction for use)
OF THE MEDICAL DEVICE
Contex® textured natural rubber latex condoms

Reckitt Benckiser Healthcare (UK) Ltd., UK

NAME OF THE MEDICAL DEVICE

Contex® textured natural rubber latex condoms

VARIANTS

- 1) Dotted – with dots
- 2) Ribbed – with ribs.

MANUFACTURER

Reckitt Benckiser Healthcare (UK) Ltd., 103-105 Bath Road, Slough, Berkshire, SL1 3 UH, UK.

MANUFACTURING SITE*

SSL Manufacturing (Thailand) Ltd., Wellgrow Industrial Estate 96 & 100, Moo 5, Bangna-Trad Road KM.36, Bangsamak, Bangpakong, Chachoengsao, 24180, Thailand;

**information about the place of manufacture is presented in the outer side of the secondary/consumer packaging in the format "The country of origin: [country]"*

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER

Reckitt Benckiser Healthcare LLC, Floor 3, 4 Shlyuzavaya emb., Moscow, 115114, Russia.

Phone: +7 (495) 961 25 65; fax: +7 (495) 961 25 66

PURPOSE OF THE MEDICAL DEVICE

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Medical Device is intended to be used by laypersons.

MEDICAL DEVICE COMPOSITION AND DESCRIPTION

Condoms Contex Ribbed, with ribs, straight walled, teat ended, in silicone lubricant. Made from natural latex

Condoms Contex Dotted, with dots, straight walled, teat ended, in silicone lubricant. Made from natural latex

TECHNICAL CHARACTERISTICS**

Contex Ribbed	
Technical parameter	Specification
Condom length, mm	180 ± 10
Condom width, mm (measured at 30 mm from the bead)	53 ± 2
Teat length, mm	Min 11
Thickness, mm	Not applicable to textured condoms
Bead thickness, mm	1-2

Contex Ribbed	
Technical parameter	Specification
Color	Coral
Condom weight (after removed from the primary package, W2), g	2.0 ± 0.5
Elongation, %	Median 700 Min
Elongation, % (aged)	Median 600 Min
Tensile strength, mPa	Median 20 Min
Tensile strength, mPa (aged)	Median 17 Min
Freedom from holes	Impermeable
Color retention	No signs of staining on absorbent paper

Contex Dotted	
Technical parameter	Specification
Condom length, mm	180 ± 10
Condom width, mm (measured at 30 mm from the bead)	52 ± 2
Teat length, mm	Min 11
Thickness, mm	Not applicable to textured condoms
Bead thickness, mm	1-2
Condom weight (after removed from the primary package, W2), g	1.8 ± 0.5
Elongation, %	Median 700 Min
Elongation, % (aged)	Median 600 Min
Tensile strength, mPa	Median 20 Min
Tensile strength, mPa (aged)	Median 17 Min
Freedom from holes	Impermeable

****main technical characteristics are presented in the outer side of the secondary/consumer packaging as "Contex® textured natural rubber latex condoms Dotted – with dots: nominal length – 180 mm; nominal width – 52 mm"; Contex® textured natural rubber latex condoms Ribbed – with ribs: nominal length – 180 mm; nominal width – 53 mm".**

PACKAGE CONTENTS OF THE MEDICAL DEVICE***

Each separate medical device variant is packed in a secondary / consumer package of 3 pcs. a pack. It is accompanied by the instruction for consumers that can be enclosed in the package or printed in the inside and / or outside of the package.

****quantity of condoms is presented in the outer side of the secondary/consumer packaging, the instruction for consumers is printed in the inner side of the secondary/consumer packaging*

LABELLING INFORMATION

Each individual package is labelled with:

- trademark;
- lot number (LOT);
- expiry date.

Symbols used for labelling and what they mean:


	Use by/Expiry date
	Lot number

Additional information can be printed on the individual package including reference to international standards.

The following information is applied to the outer surface of the secondary/ consumer package:

- condom description;
- quantity of condoms in consumer package;
- nominal condom length;
- nominal condom width;
- indication that the condom is lubricated;
- indication that the lubricant does not contain spermicide;
- indication that the condom contains natural latex;
- storage guidelines;
- name and address of the manufacturer;
- country of origin (if applicable);
- lot number (LOT);
- expiry date.

Symbols used for labelling and what they mean:

	Use by/Expiry date
	Lot number
	Contains natural latex

	For single use only
	Caution!

Additional information can be printed on the outer surface of the secondary/ consumer package including reference to international standards.

Outer or inner surface of consumer packages and/or package insert contains prescribing information, which includes:

- reminder to handle the condom carefully, including when removing it out of the package;
- indications, contraindications, precautions for use, possible side effects;
- method of application;
- necessity to remove condom right after sexual intercourse;
- necessity to exclude petroleum- and oil-based lubricants;
- necessity to use safe additional lubricant if needed;
- necessity to consult a specialist on compatibility with topical medications;
- disposal instructions for used condoms;
- condoms are disposable products.

INDICATIONS

- prevention of unwanted pregnancy;
- protection from HIV and other sexually transmitted diseases.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

PRECAUTIONS FOR USE

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections. Check the expiry date on the condom wrapper before you use it. If the individual container is obviously damaged, discard that condom and use a new one from an undamaged package.

Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil), some pessaries and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist:

- about the compatibility of condoms with other contraceptives
- if you are using topical medicines applied to penis or vagina;
- if irritation/discomfort persists after use;
- if you are unsure if you are allergic to latex;
- if a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

For anal sex or if it is necessary use additional lubrication on the outside of the condom.

Condoms may present a choking hazard. Keep out of reach of children. If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows.

Use a condom only once. Reuse may increase risk of failure or infection.

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which might cause allergic reactions including anaphylactic shock if user is allergic to latex.

APPLICATION

Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.

Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

Dispose in a bin. Do not flush the condom down the toilet.

STORAGE AND SHELF LIFE

Store the product in a dry place at 0-25 °C. Avoid direct sunlight.

Shelf life - 5 years.

MEDICINAL OR PHARMACEUTICAL SUBSTANCE CONTAINED IN THE MEDICAL DEVICE

Not applicable.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN****

Compounded natural rubber latex contains casein (0.02% ± 10% of the contents of it, functions as a stabilizer) that is considered as a derivative from tissue/materials of animal origin as it is extracted from milk.

Casein (milk protein) is a natural by-product of milk. It has no known toxic effects although it may cause a reaction in individuals with an allergy to milk. The exposure from the condoms is considered not to pose any additional risk to users of the condoms.

The device does not contain any human blood derivatives.

*****warning about possible side effect in the form of the allergic reaction because of the natural latex contained in the medical device composition is presented in the inner side of the secondary/consumer package*

STERILIZATION

The medical device is presented as a non-sterile device and is not intended to be sterilized thereafter.

USAGE WITH OTHER MEDICAL DEVICES*****

Condoms can be used on their own to prevent from STIs and pregnancy and aren't intended to be used in conjunction with other medical devices.

However, consumers can use condoms with condom compatible lubricants if needed.

******recommendation to use compatible with condoms lubricant when it is necessary is presented in the inner side of the secondary/consumer packaging*

INSTALLATION, MAINTENANCE, CALIBRATION

Not applicable.

CLEANING, DISINFECTION, STERILISATION OF THE MEDICAL DEVICE

Not applicable.

INTERNATIONAL STANDARDS

ISO 13485

ISO 14971

ISO 4074

PRECAUTIONS TO BE TAKEN BY USER RELATED TO THE DISPOSAL OF THE DEVICE

Contex® textured natural rubber latex condoms should be disposed in a bin. The condoms should not be flushed the condom down the toilet.

DATE OF THE ISSUE / THE LATEST REVISION OF THE MAINTENANCE DOCUMENTATION

The medical device includes an indication of the latest revision within the instruction for use.

GUIDELIENSS FOR DESTRUCTION AND DISPOSAL

The devices should be used in accordance with the instructions for use.

In household use, the product is disposed of as household waste. Do not put it down the drain.


In all other cases, the condoms are disposed of in accordance with local requirements.

ENVIRONMENT SAFETY REQUIREMENTS

With proper usage, storage and disposal Contex® textured natural rubber latex condoms have no negative impact on the environment.

**Appendix to the Maintenance
documentation of the medical device
«Contex® textured natural rubber latex condoms »**

*Information
to be printed in the inner side
of the secondary/ consumer package*

Read the instruction carefully before use, especially if condoms are used for oral or anal sex.

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections.

ATTENTION: Check the expiry date on the condom wrapper before you use it. If the individual container is obviously damaged, discard that condom and use a new one from an undamaged package.

PURPOSE

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases.

APPLICATION

1. Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.

2. Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis.

3. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

4. Dispose in a bin. Do not flush the condom down the toilet.

For anal sex or if it is necessary use additional lubrication on the outside of the condom.

Condoms may present a choking hazard. Keep out of reach of children.

If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows.

Use a condom only once. Reuse may increase risk of failure or infection.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which might cause allergic reactions including anaphylactic shock if user is allergic to latex.

PRECAUTIONS FOR USE

Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil), some pessaries and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist:

- about the compatibility of condoms with other contraceptives
- if you are using topical medicines applied to penis or vagina;
- if irritation/discomfort persists after use;
- if you are unsure if you are allergic to latex;
- if a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

Date of the latest revision of the instruction:

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