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INSTRUCTIONS FOR USE OF THE MEDICAL DEVICE

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CUTTING INSTRUMENTS FOR BONE TISSUE HOLES HANDLING WITH ACCESSORIES

Official Certification:

The undersigned Notary certifies, that the preceding signature of

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is authentic and that the signature currently act in the position of Ass. Director RA EU, Synthes GmbH, a company currently existing an organized under the laws of Switzerland, that the adress is Eimattstrasse 3, 4436 Oberdorf.

The Notary Public is aware and attests to the above facts upon examination of the corresponding documents.

Given and signed in 4500 Solothurn, Switzerland

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- mm
16. Spare Reamer Tube, for No. 309.065
 17. Spare Reamer Tube, for No. 309.150
 18. Hollow Reamer, anticlockwise cutting, for Screws Ø 1.5mm, complete
 19. Spare Reamer Tube, for No. 309.200
 20. Hollow Reamer, anticlockwise cut complete
 21. Hollow Reamer, anticlockwise cutting, complete
 22. Spare Reamer Tube, for No. 309.250
 23. Hollow Reamer, anticlockwise cutting, for Screws Ø 4.5mm, complete
 24. Spare Reamer Tube, for No. 309.450
 25. SynReam Medullary Reamer Head Ø 8.5 mm
 26. SynReam Medullary Reamer Head Ø 9.0 mm
 27. SynReam Medullary Reamer Head Ø 9.5 mm
 28. SynReam Medullary Reamer Head Ø 10.0 mm
 29. SynReam Medullary Reamer Head Ø 10.5 mm
 30. SynReam Medullary Reamer Head Ø 11.0 mm
 31. SynReam Medullary Reamer Head Ø 11
 32. SynReam Medullary Reamer Head Ø 12.0 mm
 33. SynReam Medullary Reamer He
 34. SynReam Medullary Reamer Head Ø 13.0 mm
 35. SynReam Medullary Reamer Head Ø 13.5 mm
 36. SynReam Medullary Reamer Head Ø 14.0 mm
 37. SynReam Medullary Reamer Head Ø 14.5 mm
 38. SynReam Medullary Reamer Head Ø 15.0 m
 39. SynReam Medullary Reamer Head Ø 15.5 mm
 40. SynReam Medullary Reamer Head Ø 16.0 mm

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Further in the text the following names of the medical device can be used: cutting instruments, device, product, MD, instruments

The names of medical devices can be mentioned in the text form for the convenience of perceiving in

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

Functional and technical characteristics of the medical device

Note: The following images are for informational purposes only and intended to demonstrate the General appearance and design of the MD, and are not to scale.

Note: The parameters specified in the MD name are part of the name and presented for the convenience of the user's choice

The table of technical characteristics contains the accurate parameters of MD.

Ref	Name	
03.111.032	Reamer for VA-Locking Intercarpal Fusion Plate Ø 15 mm	
03.111.033	Reamer for VA-Locking Intercarpal Fusion Plate Ø 17 mm	
03.211.016	Cannulated Distal Reamer Ø 16 mm	
03.211.018	Cannulated Distal Reamer Ø 18 mm	
03.211.020	Cannulated Distal Reamer Ø 20 mm	
03.211.022	Cannulated Distal Reamer Ø 22 mm	
03.211.024	Cannulated Distal Reamer Ø 24 mm	

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Table 2. Transportation conditions:

Designation	Range
Temperature	-
Relative humidity	30
Barometric pressure	No special

Table 3. Storage conditions:

Designation	
Temperature	
Relative humidity	
Barometric pressure	

MAIN STANDARDS USED BY MANUFACTURER

The medical devices meet the requirements of the following
 Complete list of applicable standards is available upon request.

Standard	Title
EN ISO 13485	Medical device – Requirements for
EN ISO 14971	Medical devices management to medical devices
EN ISO 10993-1	Biological evaluation of medical devices Part 1: Evaluation and testing within a risk management system

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Description of sterilization method

Synthes products supplied in a non-sterile condition must be cleaned and steam sterilized prior to surgical use. Prior to cleaning, remove and dispose all original disposable packaging (e.g. protection before first and every use, and before returning for maintenance and repair. Prior to steam sterilization, place the product in an approved sterilization wrap or container.

The first and most important step in reprocessing all re thorough (manual and/or mechanical) cleaning and rinsing. Thorough cleaning is a complex process whose success depends on various interrelated factors: Water quality, quantity and type of (manual, ultrasonic bath, washer/disinfector), thorough rinsing and drying, proper product preparation, time, temperature, and thoroughness of the individual responsible for cleaning.

Residual organic matter and/or a large reduce the effectiveness of the sterilization process.

Cautions

All devices must be thoroughly cleaned and inspected prior to sterilization. Blind holes, moving and intricate parts require particular attention during cleaning and inspection. During c detergents that are labelled for use on medical devices and in accordance with the manufacturer's instructions (e.g. temperature, contact time, and rinse time). Cleaning agents with a used dilution pH of within 7 recommended. Highly alkaline conditions (pH >11) can damage components/devices. Do not use saline, environmental disinfection (including chlorine solutions) or surgical antiseptics (such as iodine

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chlorhexidine-containing products). Do not use a damage the surface of instruments such as steel wool, abrasive cleaners or wire brushes.

- Only place Synthes devices with items of similar metallic composition together in an ultrasonic cleaner.
- Soiled or used Synthes devices should not be cleaned in a mechanical washer. Soiled Synthes devices must be processed separate from trays and cases.
- The sterilization parameters are only valid for devices that are adequately cleaned.
- The parameters listed are only valid if the reprocessing equipment is maintained, calibrated and compliant reprocessing equipment in accordance with standards such as the ISO 15883 and ISO 17665 series.
- Surgical patients identified as at risk for prion disease (CJD) and related infections. Dispose of instruments used or suspected of use on a patient with CJD after surgery and/or follow current national recommendations.
- Consult national regulations and guidelines for additional information. Compliance is additionally required with hospital policies and procedures and recommendations of manufacturers of detergents, disinfectants, and any clinical processing equipment.

Limits on reprocessing

- Repeated processing cycles that include ultrasonic washing and sterilization have minimal effects on Synthes surgical instrumentation.
- End of life of a device is normally determined by wear and damage

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due to use. Evidence of damage and wear on a device may include but is not limited to corrosion (e.g. rust, pitting, excessive scratches, flaking, wear and cracks. Improperly functioning devices, devices with unrecognizable markings, missing or removed (buffed off) part numbers, damaged and excessively worn devices should not be used.

- Further details regarding end of life indicators are available from Authorized representative in Russia.

Number of reprocessing cycles

Synthes surgical instruments are delivered non clinical reprocessing cycles and a limited shelf on functional control requirements. The instruments is multifactorial, which is dependent upon but not limited to following factors: the design of the instrument, the cleaning and handling of the instrument and the intensity of use of the instrument.

Point of Use Care

Wipe blood and/or debris from device throughout surgical procedure to prevent it from drying onto the surface.

- Soiled devices should be separated from non avoid contamination of personnel or surroundings,
- Devices should be covered with a towel dampened with sterile or purified water to prevent blood and/or debris from drying.

Containment and Transportation:

- Soiled devices should be transported separate from devices to avoid contamination.

Preparation for Cleaning (for all cleaning methods):

- It is recommended that devices should be reprocessed as soon as

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reasonably practical following use.

- Disassemble device, if device is able to be disassembled for reprocessing. Further detailed instrument dismantling instructions are available from Authorized representative in Russia .
- Open devices with hinges.
- Care should be taken in the handling and cleaning of sharp devices. These are recommended to be cleaned to prevent injury.
- Soak and/or rinse heavily soiled devices prior to cleaning to loosen any dried soil or debris. Use a neutral or mild alkaline detergent solution. Follow the detergent manufacturer's correct exposure time, temperature, water quality and concentration. Use cold tap water to rinse devices.

Synthes devices must be cleaned separately from instrument trays and cases. Lids should be removed from cases for the cleaning process.

Cleaning – Manual Method

1. Rinse soiled device under running cold water for a minimum of ten minutes. Use a soft-bristled brush to assist in the removal of gross soil and debris.
2. Soak device in a neutral or mild alkaline detergent solution for ten minutes. Follow the detergent manufacturer's instructions for correct exposure time, temperature, water quality and concentration.
3. Rinse device with cold water for a minimum of two minutes. Use a syringe, pipette, or water jet to flush hard to reach areas.
4. Manually clean device for a minimum of five minutes in a prepared neutral or mild alkaline detergent solution. Use a soft brush to remove soil and debris. Actuate joints, handles and

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features to expose all areas to the detergent solution, if applicable. Cle device under water to prevent aerosolization of contaminants. Note: fresh solution is a newly-made, clean solution.

5. Rinse device thoroughly using cold o two minutes.. Actuate joints, handles and other moveable device order to rinse thoroughly under running water, if applicable.
6. Visually inspect device. Repeat steps 2 device.

Ultrasonic process: Pre-cleaning steps 1

7. Prepare a fresh detergent solution for the ultrasonic bath using a neutral or mild alkaline detergent solution. Follow the detergent manufacturer's instructions for use for correct exposure time, te concentration. Note: fresh solution is a newl
8. Clean Synthes device ultrasonically for a minimum of 15 minutes, using a minimum frequency of 38 kHz.
9. Rinse device thoroughly with deionized (DI) or pur for a minimum of two minutes. Actuate joints, handles a device features in order to rinse thoroughly under running water, if applicable.
10. Visually inspect device. Repeat steps 2 on device.
11. Perform a final rinse on device using DI or PURW water for a mini of 15 seconds.
12. Dry devices using a clean, soft, lint grade compressed air. Ensure that all lumens and articulated are using compressed air.

Cleaning-Disinfection Automated Method

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1. Rinse soiled device under running cold tap water for a minimum of one minute. Remove gross soil using a soft-bristled brush or soft, lint

2. Manually clean device for a minimum of two minute prepared neutral or mild alkaline detergent solution. Follow manufacturer's instructions for the correct dilution, temperature, water quality and exposure time. Use a soft debris. Actuate joints, handles and other moveable device features to expose all areas to detergent solution, if applicable. Clean device under water to prevent aerosolization of contaminants, Note: fresh solution is a newly made, clean solution.

3. Rinse device using cold to lukewarm runni one minute. Actuate joints, handles and other moveable device features in order to rinse thoroughly under running water, if applicable.

4. Prepare a fresh detergent solution for the ultra mild alkaline detergent. Follow the detergent manufacturer's instruc the correct dilution, temperature, water quality and exposure time. Note: fresh solution is a newly-made, clean solution.

5. Clean Synthes devices ultrasonically for a mini a minimum frequency of 38 kHz.

6. Rinse device using DI or PURW water for a minimum of two minutes. DI or PURW water must be used for final rinse.

7. Visually inspect device. Repeat steps 2 device.

8. Automated washing shall be conducted in a validated washer in compliance to ISO 15883-1 and -2, or to an equivalent standard. Load the device components in the washer manufacturer's instructions, ensuring tha Automated washing can be included as part disinfection, and/or drying cycle in accordance to manufacturer's instructions.

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An example of a validated cycle used for cleaning va

Cycle	Minimum Time (minutes)	Minimum Temperature/Water
Pre-wash	2	Cold tap water
Wash I	2	Cold tap water (<40°C)
Wash II	5	Warm tap water (>40°C)
Rinse	2	Warm DI or PURW (>40°C)
Thermal Disinfection	5	> 93°C
Dry	40	> 90°C

* see Additional information

Thermal disinfection

For automated cleaning-disinfection, thermally disinfect at a minimum of 93 °C for a minimum of 5 minutes (see Cleaning Method).

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If a dry cycle is not included in the mechanical washer:

- Dry each device thoroughly inside and out to prevent rusting and malfunction.
- Use a clean, soft, lint-free single surface.
- Pay special attention to threads and hinges or areas accumulate. Open and close devices so that all areas are hollow parts using the air jet with medical grade compressed air.

Inspection

Synthes instruments should be inspected after processing, prior to sterilization, for end of life indicators such as:

- Cleanliness, if any residual soil is discovered repeat the cleaning steps on those devices until all visible soil is removed from the device.
 - Lack of moisture, carefully inspect device is detected, manual drying should be performed.
 - Damage, including but not limited to corrosion (e.g. rust, pitting), discoloration, excessive scratches, flaking, cracks and wear.
 - Proper function, including but not limited to movement hinges/joints and moveable features
- Damaged or worn devices should not be used.

Further detailed function control instructions and end of life indicators are available from your Authorized representative in Russia.

Lubricate instruments with moving parts, threaded parts. It is recommended to lubricate and maintain Synthes instruments with Synthes Special Oil only.

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Packaging

Care should be taken to protect pointed and sharp instruments with other objects that may damage the surface.

Sterilization

Steam (moist heat) sterilization shall be performed in a locally approved, pre-vacuum (forced air removal) cycle. The steam sterilizer should be validated to the requirements of any local standards and 285 or AAMI/ANSI ST8. The steam sterilizer should be maintained in compliance to manufacturer's instructions and local requirements.

The following steam sterilization cycles are examples of vali

Cycle Type	Minimum Sterilization Exposure Time (minutes)	Minimum Sterilization Expo Temperature
Prevacuum		
Saturated steam-	4	132°C

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forced air
removal
(pre-
vacuum) 3 134°C
(minimum
three
pulses)

* The dry time is most often influenced by the presence of polymer based (plastic) materials; therefore changes such as elimination of silicone mats and/or change in sterile barrier system (i.e. heavy grade to light grade wrap) can reduce the necessary dry time. The user should use methods (e.g. visual inspections) to confirm adequate dry time. Dry times generally range from 20 to 60 minutes due to differences in packaging materials (Sterile Barrier System, e.g. wraps or reprocessor systems), steam quality, device materials, total mass, sterilizer performance, and varying cool-down time.

Storage

Packaged products should be stored in a dry, clean environment, protected from direct sunlight, pests, and extremes of temperature and humidity.

Additional Information

Further information regarding the use of specific cleaning agents, washers, washer-disinfector, packaging materials or sterilizers during validation studies are available on request. Synthes used the following validation of these reprocessing recommendations.

- Manual cleaning: Manual Pre-Detergent 8 mL/L at 16–17 °C and Ultrasonic Cleaning with Enzol Enzymatic Detergent 8 mL/L at 18 °C

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- Automated Cleaning: Manual Pre Detergent 8 mL/L at 16–17 °C and Ultrasonic Cleaning with Enzol Enzymatic Detergent 8 mL/L at 18 with (Wash 1) Prolystica 2X Concentrate
- Enzymatic Cleaner 1 mL/L at 23–
- Neutral Detergent 1 mL/L at 44–46
- Lint-Free Cloth: Berkshire Durx 670.

The cleaning and sterilization information is provided in accordance with ISO 17664.

The recommendations provided above have been validated by the medical device manufacturer as being capable of prepa medical device. It remains the responsibility o the processing is performed using equipment, materials and personnel in the reprocessing facility, and achieves the desired res and routine monitoring of the process. Likewise, any deviation by th processor from the recommendations provided should be properly evaluated for effectiveness and potential adverse consequences.

The chemical quality of the water used during reproce device safety. Facilities should use the recommended wate requirements for device reprocessing in accordance with local guidance (such as AAMI TIR 34, Water for the reprocessing of medical devices) and these instructions for use.

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