21

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CERTIFICATE

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China Council for the Promotion of International Trade

China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of BLUE SAIL MEDICAL CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion of International Trade

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Authorized Signature:

Wei Shanshan

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(Date: Aug. 04, 2021)

INSTRUCTION FOR USE FOR MEDICAL DEVICE:

Medical diagnostic (examination) disposable non-sterile gloves "A.D.M."
vinyl smooth powder free with inner polyurethane coating

APPROVED BY:

Sales Manager Blue Sail Co., Ltd.

Mr. Bin Li /

*Initial issue of operational documentation (instruction for use for medical device:
«Medical diagnostic (examination) disposable non-sterile gloves "A.D.M."
vinyl smooth powder free with inner polyurethane coating»): 10/07/2020.*

2020

TABLE OF CONTENTS

1 Data on the medical device	3
1.1 Name of the medical device.....	3
1.2 Manufacturer (producer) of the medical device.....	3
1.3 Purpose of the medical device	3
1.4 Scope of the medical device.....	3
1.5 Intended users of the medical device	3
1.6 Application (operating) conditions of the medical device	3
1.7 Operating principle of the medical device	4
1.8 Indications for use of the medical device.....	4
1.9 Contradictions for use of the medical device.....	4
1.10 Shelf life of the medical device.....	4
2 Warning and precautions	5
3 Instruction for use of the medical device	6
4 Transportation, storage and disposal conditions of the medical device.....	7
4.1 Transportation conditions of the medical device	7
4.2 Storage conditions of the medical device	7
4.3 Disposal conditions of the medical device.....	7
5 Cleaning and disinfection methods of the medical device.....	7
6 Specification of the medical device	8
6.1 Description of the medical device.....	8
6.2 Physical specification of the medical device	10
6.2.1 Physical and mechanical performance of the medical device.....	10
6.2.2 Inspection levels and acceptable quality level of the medical device.....	10
7 Adverse reactions	11
8 Circumstances, in which the consumer must consult healthcare professional	11
9 Manufacturer's warranties	11
10 Data on Authorized representative of the manufacturer	11
11 Scope of delivery of the medical device	11
12 List of national standards, applied by manufacturer (producer) of the medical device.....	12

1 Data on the medical device

1.1 Name of the medical device

Medical diagnostic (examination) disposable non-sterile gloves “A.D.M.” vinyl smooth powder free with inner polyurethane coating, in hardware versions: (beaded cuff), color: clear, white, black, pale-yellow, blue, size: XS, S, M, L, XL.

1.2 Manufacturer (producer) of the medical device

Manufacturer (producer) of the medical device: Blue Sail Medical Co., Ltd. (No.21 Qingtian Road, Qilu Chemical Industrial Area, 255414, Zibo, Shandong Province, P.R. China).

Address of the medical device manufacture (production): Blue Sail Medical Co., Ltd. (No.21 Qingtian Road, Qilu Chemical Industrial Area, 255414 Zibo, Shandong Province, PEOPLE'S REPUBLIC OF CHINA)..

1.3 Purpose of the medical device

The gloves are intended for single use as protective barrier in the process of wearing on healthcare professional's hands during examination/treatment of patient.

1.4 Scope of the medical device

The gloves are used during medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces, providing hand hygiene, as also protection of patient and healthcare professional from cross-contamination.

1.5 Intended users of the medical device

Medical personnel, who carry out medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces and other activities, which require compliance with the requirements for protection of patient/healthcare professional from cross-infection.

1.6 Application (operating) conditions of the medical device

Medical examination disposable gloves are used in medical or healthcare facilities indoors in a dry place, away from direct sunlight, at temperature from 10°C to 30°C and relative humidity up to 80% in accordance with the instruction for use on medical device.

Gloves for external use only (the device is non-invasive).

1.7 Operating principle of the medical device

The basis of examination disposable vinyl gloves limits the movement of infectious substances, improves the efficiency of antibacterial protection in the healthcare facility, protecting the healthcare professional against possible infection during examination/treatment of the patient and also contributes to the reliable capture of medical instruments and increases the safety of the work.

1.8 Indications for use of the medical device

The indications for use of the medical device are the necessity to carry out medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces.

1.9 Contradictions for use of the medical device

The contradiction of the medical device is the patient's or healthcare professional's allergic reaction to residual process chemicals in the medical device.

1.10 Shelf life of the medical device

The established shelf life of examination disposable non-sterile vinyl gloves is 3 years.

2 Warning and precautions

- The device is intended for single use;
- For individual use only;
- Every time it is necessary to perform hygienic treatment of your hands before putting on or changing gloves, as well as after removing the gloves;
- It is necessary to choose the gloves of the right size – they should be easily put on and should not impede movement or squeeze the hands. Do not use in case of size mismatch;
- Touching is recommended to the glove surface corresponding to the wrist (upper edge of the cuff);
- It is necessary to change the gloves in case of any defect, puncture, damage, contamination, fluid contact with skin under the glove, as well as stickiness;
- If one glove is damaged, change both ones;
- Depending on the type of the performed manipulation, to keep an optimal level of protection, it is recommended to change gloves every 60-120 minutes;
- Do not wet and keep away from moisture;
- It is necessary to put the gloves on dry hands only;
- Avoid exposure to direct sunlight;
- The device contains plasticizers (DINP, DOTP), which can be unsafe for user;
- Gloves contain residual process chemicals that may cause an allergic reaction;
- Gloves do not guarantee 100% hand protection against punctures, cuts and damage caused by exposure to chemicals;
- In case of skin allergic reactions when using the device it is recommended to consult a healthcare professional.

3 Instruction for use of the medical device

Preliminary hygienic treatment of your hands is required before use of the gloves (treat hands with skin antiseptic or wash with soap and water). Gloves are put on dry hands only.

How to put the gloves on: Take a glove from the box. It is recommended to touch the bounded surface of the glove corresponding to the wrist (upper edge of the cuff). Pull on the first glove. Take the second glove with your bare hand, touching the bounded surface of the glove corresponding to the wrist. To avoid touching the skin of the forearm with a gloved hand, rotate the outer surface of the glove to be put on with folded fingers in the glove so that you could pull the glove over the other hand. Having put the gloves on, it is required to touch only those objects provided for by the indications for use and conditions of the use of the gloves. To avoid contamination, it is required to change gloves before proceeding with the new procedure.

The used gloves should be put off carefully, avoiding splashing from the glove's surfaces where microbial contamination of hands and environment is possible.

How to take the gloves off: To take a glove off, it is required to pick it up at the level of the wrist without touching the skin of the forearm and pull it off your hand, turning it inside out. Holding the removed glove with the gloved hand, it is necessary to place the fingers of your bare hand between the glove and wrist. Remove the second glove by rolling it down the hand and put into the first glove. Put the gloves in the container designated for them. Then perform hygienic treatment of your hands.

4 Transportation, storage and disposal conditions of the medical device

4.1 Transportation conditions of the medical device

The transportation of the device may be made by all means of transport (sea, land, railway or air transport) according to the rules of transportation, which are applied to such means of transport. During transportation, the storage conditions should not be violated.

4.2 Storage conditions of the medical device

Devices should be stored in a dry place protected from direct sunlight, heat, ignition and ultraviolet sources, as well as ozone, at temperature from 10°C to 30°C and relative humidity of air – up to 80%.

4.3 Disposal conditions of the medical device

Disposal of the medical device in the Russian Federation is carried out by the medical organization as class B waste in accordance with SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of urban and rural settlements, for water facilities, drinking water and drinking water supply, atmospheric air, soils, living quarters, the operation of industrial, public premises, the organization and implementation of sanitary and anti-epidemic (preventive) measures" and approved guideline of the organization.

5 Cleaning and disinfection methods of the medical device

Examination gloves are disposable medical device. Reuse of the medical device is a gross violation of application methods.

6 Specification of the medical device

6.1 Description of the medical device

Form of the medical device. Medical examination gloves have straight fingers and are ambidextrous (non-anatomic). The glove's cuff has a standard length. Rolled cuff.

Surface of the medical device. Gloves have a powder free surface – with inner polyurethane coating.

Texture of the outer surface of the medical device. Gloves have a smooth surface.

Sizes. Depending on the size, table 1 gives the values of width, length and thickness of the medical device.

Table 1 – Sizes of the medical device:

Value name	Size				
	XS	S	M	L	XL
Width W (figure 1), mm	≤ 80	80±10	95±10	110±10	≥110
Length l (figure 1), mm, min	220		230		
Thickness (in points, shown on figure 2), mm, min	0.08				
Thickness of rolled cuff, mm, max	2.50				


Figure 1 – Points for width and length measurement of the glove


Figure 2 – Points for wall thickness measurement of the glove

6.2 Physical specification of the medical device

6.2.1 Physical and mechanical performance of the medical device

Physical and mechanical parameters of the medical device are represented in table 2.

Table 2 – Values of physical and mechanical parameters of the gloves:

Description	Specification
<u>Before/ after accelerated aging:</u>	
- elongation at break, %, min	350
- force at break, N, min	7.0

6.2.2 Inspection levels and acceptable quality level of the medical device

For the purposes of quality control medical gloves are selected and examined in accordance with standard ISO 2859-1. Inspection levels and acceptable quality level (AQL) of the gloves are represented in table 3.

Table 3 – Inspection levels and acceptable quality level (AQL):

Performance	Inspection level	Acceptable quality level (AQL)
Physical dimension (width, length, thickness)	S-2	4.0
Hermeticity	G-1	2.5
Force and elongation at break	S-2	4.0

7 Adverse reactions

Manifestation of dermal hypersensitivity reactions, when using the device is possible.

8 Circumstances, in which the consumer must consult healthcare professional

The patient is recommended to consult healthcare professional in the event of dermal hypersensitivity reactions, when using the device.

9 Manufacturer's warranties

The manufacturer (producer) guarantees preservation of all the parameters and characteristics of the medical device taking into account acceptable quality level (AQL), as well as the quality of the medical gloves within all shelf life from the date of manufacture (3 years), upon the condition of package integrity, observance of transportation, storage and operation conditions in accordance with the instruction.

10 Data on Authorized representative of the manufacturer

Please contact the Authorized representative of the manufacturer for any and all purposes, related to circulation of the medical device on the territory of the Russian Federation:

GL Limited Liability Company (GL LLC)

Russia, 194044, St. Petersburg, 9, Mendelevskaya st., litera V, office 10-N, room 7

Tel.: +7 (812) 309-92-66

e-mail: GL3099266@gmail.com

11 Scope of delivery of the medical device

Gloves can be delivered in the following scopes of supply:

- bulk packaging – 2000 gloves, 1000 pairs in consumer package;
- 400 gloves, i.e. 200 pairs in consumer package;
- 200 gloves, i.e. 100 pairs in consumer package;
- 100 gloves, i.e. 50 pairs in consumer package;
- 50 gloves, i.e. 25 pairs in consumer package;
- 40 gloves, i.e. 20 pairs in consumer package;
- 30 gloves, i.e. 15 pairs in consumer package;
- 20 gloves, i.e. 10 pairs in consumer package;
- 10 gloves, i.e. 5 pairs in consumer package;
- 4 gloves, i.e. 2 pairs in consumer package;
- 2 gloves, i.e. 1 pair in consumer package.

12 List of national standards, applied by manufacturer (producer) of the medical device

Table 4 – List of standards, applied by manufacturer:

Designation	Name
ISO 13485:2016	Medical devices. Quality management systems. Requirements for regulatory purposes.
ISO 14971:2007	Medical devices – Application of risk management to medical devices.
ISO 2859-1	Statistical methods. Sampling procedures for inspection by attributes. Part 1. Sampling schemes indexed by acceptance quality limit for lot-by-lot inspection.
ISO 15223-1:2012	Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied. Part 1. General requirements.
ISO 10993-1:2003	Biological evaluation of medical devices. Part 1: Evaluation and testing.
ISO 10993-10:2002	Biological evaluation of medical devices. Part 10: Tests for irritation and delayed-type hypersensitivity.
ISO 10993-11:2006	Biological evaluation of medical devices. Part 11: Tests for systemic toxicity.
ISO 11193-2:2006	Single-use examination gloves. Part 2: Specification for gloves made from poly (vinyl chloride)

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