APOSTILLE**(Convention de La Haye du 5 octobre 1961)****1. Country:** Sweden**This public document****2. has been signed by** Gabriel Rudbeck**3. acting in the capacity of** Notary Public**4. bears the seal/stamp of** Notary Public in
Stockholm**Certified****5. at** Stockholm**6. the** 2020-06-29**7. by** Adrienne Bonde

Deputy Notary Public

8. No. 5078**9. Seal/stamp:****10. Signature:**

APPROVED

General Director, Boule Medical LLC

INSTRUCTION
FOR USE OF MEDICAL DEVICE
"REAGENTS IN VITRO FOR HEMATOLOGICAL TESTS"

Order Information:

Art. No.	Art. Name	Russian Description
1504124	Swelab AlfaDiluent 20L	Isotonic Diluent Swelab AlfaDiluent, 20 L
1504462	Swelab AlfaDiluent, RFID	Isotonic Diluent Swelab AlfaDiluent, RFID, 20L
1504125	Swelab AlfaLyse 5L	Hemolytic Reagent Swelab AlfaLyse, 5 L
1504463	Swelab AlfaLyse, RFID	Hemolytic Reagent Swelab AlfaLyse RFID, 5L
1504127	Swelab Alfa ComboPack 200	Pack Swelab Alfa ComboPack 200, (5+1) L
1504464	Swelab Alfa ComboPack, RFID	Pack Swelab Alfa ComboPack, RFID, (5+1) L
1504112	Boule Enzymatic Cleaner 100mL	Cleaning Solution Boule Enzymatic Cleaner 100 mL
1504113	Boule Hypochlorite 2% Cleaner 500mL	Cleaning Solution Boule Hypochlorite 2% Cleaner 500 mL
1504111	Boule Cleaning Kit, 3x450mL	Cleaning Solution Boule Cleaning Kit, 3x450 mL

1. INTENDED USE

1.1. "Reagents in vitro for hematological tests" are intended for use in blood tests and maintenance of the Swelab Alfa and Swelab Alfa Plus instruments in clinical diagnostic laboratories.

1.2. Field of use – clinical diagnostic laboratories.

1.3. Intended user – staff of the clinical diagnostic laboratories who have undergone training in working with hematology analyzers.

Not intended for use near patient.

1.4. Complete blood count (CBC) is prescribed for screening and dispensary examinations, differential diagnosis of blood diseases and monitoring of ongoing therapy.

2. SPECIFICATION OF THE DEVICE

2.1. Components and content of the reagents:

Swelab AlfaDiluent¹⁾, pH 6.70 – 6.90:

Inorganic salts to support isotonic stability	< 1.5 %
Bacteriostatic agents (cationic)	< 0.1 %
Buffer salts	< 0.3 %
Antimicrobials (kathon) ²⁾	< 0.0015 %

Swelab AlfaLyse¹⁾, pH 7.80 – 8.30:

Quaternary ammonium salts	< 1.0 %
Buffer salts	< 1.5 %

Boule Enzymatic Cleaner³⁾, pH 7.6 – 8.4

Subtilisin for proteolytic action	< 0.1 %
Surfactants	< 0.2 %
Stabilizers incl. salts and buffers	< 2 %
Coloring substances	<< 0.1 %
Antimicrobials (kathon)	< 0.0015 %

Boule Hypochlorite 2% Cleaner⁴⁾, pH 11.5 – 12.1

Sodium hypochlorite for lipid/protein degradation	2.0 – 2.4 % active chlorine
Sodium hydroxide for stability	< 0.05 %
Surfactants	< 0.05%

Boule Detergent Cleaner⁵⁾, pH 7.6 – 8.4

Sodium dodecyl sulfate for protein dissociation	< 0.1 %
Surfactants	< 0.2 %
Stabilizers incl. salts and buffers	< 1 %
Coloring substances	<< 0.1 %
Antimicrobials (kathon)	< 0.0015 %

Notes:

- 1) Also applicable to Swelab Alfa ComboPack combination packaging and reagents with RFID tags.
- 2) Stabilizing substances with antimicrobial activity against bacteria and yeast are biodegradable.
- 3) Plastic bottle, content 100 mL or 450 mL as part in Boule Cleaning Kit.
- 4) Plastic bottle, content 500 mL or 450 mL as part in Boule Cleaning Kit.
- 5) Plastic bottle, content 450 mL part in Boule Cleaning Kit.

2.2. The number of reagent tests performed depends on the Swelab hematology analyzer model and the number of tests performed per day. Approximately 22 ml of diluent and 4.5 ml of hemolytic reagent are consumed per analysis cycle.

2.3. Principal of the method.

The conductometric (Coulter) method of sizing and counting cells is based on changes in electrical impedance produced by particles suspended in an electrolyte as they pass through a small opening (aperture) between electrodes. A continuous current runs between the external and internal electrodes, and as each cell passes through the aperture it generates specific changes in impedance in the aperture channel. These changes are detected in the form of an increase in voltage between the electrodes. The number of pulses is proportionate to the number of particles and the height of each pulse is proportional to the volume of the particle. Particle size distribution is recorded as histograms of leucocytes (WBC), erythrocytes (RBC) and platelets (PLT) on the display screen.

Measuring principle for hemoglobin (HGB): the reagent lyses the blood cells, releasing hemoglobin from the erythrocytes. The hemoglobin in the lysed sample is measured photometrically. Methemoglobin is chemically stabilized and then also measured photometrically.

Swelab AlfaDiluent is an isotonic buffer solution with defined pH, conductivity and osmolality intended for pre-dilution of blood samples in automated hematology analyzers before the direct differential count. Its main purpose is to maintain the osmotic pressure required to retain blood cell volume and viability throughout the entire measurement cycle, create a conductive medium for an impedance-based differential blood cell count, and prevent fibrin clot formation and platelet aggregation.

Swelab AlfaLyse is a cyanide-free reagent for lysing erythrocytes and converting hemoglobin into a measurable complex in automated hematology analyzers for differentiation of leukocytes into three subpopulations. When combined with the diluted blood suspension, the lytic solution quickly destroys the cytoplasm membrane of the red blood cells, shrinking cell volumes and leaving the WBC subpopulation count intact. The surface-active hemolytic compounds affect the size of the different forms of WBC at different rates, and these differences are used as the basis for dividing the WBCs into subpopulations.

Swelab Alfa ComboPack is a combination pack with hematology diluent and a cyanide-free lytic reagent used for cell counting and sizing in the Swelab Alfa Series, Automated, hematology analyzers.

Boule Enzymatic Cleaner and Boule Hypochlorite 2% Cleaner are used for maintenance and cleaning of Boule hematology analyzers. Active components of the reagents provide the breakdown of protein and lipid deposits on the dosing and hydraulic system of the analyzer.

Boule Cleaning Kit is a combination pack with Enzymatic Cleaner, Hypochlorite 2% Cleaner and Detergent Cleaner and intended for use with Boule Hematology Systems. The kit contains 3 different cleaning solutions as indicated on their labels. Using this cleaning kit at the stated regular intervals will ensure:

- Longer lifetime of instrument tubing
- Lower risk of bacterial growth
- Lower background counts
- Increased mean time between failure (higher MTBF).

3. ANALYTICAL PARAMETERS

3.1. Using the reagents in the system with Swelab hematology analyzers with differentiation of leukocytes onto 3 subpopulations provides following parameters of measurements of count of the blood cells and hemoglobin concentration.

Table 3.1.1 Analytical parameters; Swelab Alfa

Parameter		RBC	MCV	HGB	PLT	WBC
Linear Range ¹		0.30 – 7.00 x 10 ¹² /L	N/A	2.0 – 24.0 g/dL	20 – 1800 x 10 ⁹ /L	0.5 – 99.9 x 10 ⁹ /L
Coefficient of variation, (%)	Specification OT, CP, AS	≤ 1.8	≤ 1.5	≤ 1.5	≤ 5.2	≤ 3.5
	Specification MCI, PD	≤ 2.7	≤ 1.5	≤ 2.4	≤ 5.2	≤ 4.0
Sensitivity ⁴		0.30 x 10 ¹² /L	N/A	2.0 g/dL	20 x 10 ⁹ /L	0.5 x 10 ⁹ /L

OT: Open Tube mode; MCI: Micro Capillary Inlet mode; CP: Cap Piercer mode; AS: Auto sampler mode; PD: Pre-dilute mode

¹ Linearity measured based on Standard EP6-A. Control blood material.

² Results presented in CV% tot. max, three instruments as a set.

³ Results presented as CV% tot. of three sample series

⁴ For circulation in the Russian Federation. Sensitivity is the lowest limit of the linearity range.

Table 3.1.2. System Performance Parameters; Swelab Alfa Plus

Parameter		RBC	MCV	HGB	PLT	WBC
Linear Range ¹		0.30-7.00 x 10 ¹² /L	N/A	2.0-24.0 g/dL	10 ⁽²²⁰⁾ -1800 x 10 ⁹ /L	0.20-130.0 x 10 ⁹ /L
Coefficient of variation, (%)	Specification OT, CP, AS	≤ 1.8	≤ 1.5	≤ 1.5	≤ 5.2	≤ 3.5
	Specification MPA, PD	≤ 2.7	≤ 1.5	≤ 2.4	≤ 5.2	≤ 4.0
Sensitivity ⁴		0.30 x 10 ¹² /L	N/A	2.0 g/dL	10 ⁽²²⁰⁾ x 10 ⁹ /L	0.20 x 10 ⁹ /L

OT: Open Tube mode; MPA: Micro Pipette Adapter mode; CP: Cap Piercer mode; AS: Auto sampler mode; PD: Pre-dilute mode

¹ Linearity measured based on Standard EP6-A. Control blood material.

² If your analyzer was manufactured before December 2016.

³ Results presented in CV% tot. max, three instruments as a set.

⁴ For circulation in the Russian Federation. Sensitivity is the lowest limit of the linearity range.

Table 3.1.3. Accuracy and permissible deviation (Swelab Alfa & Swelab Alfa Plus)

Parameter	RBC	MCV	HGB	PLT	WBC
The accuracy of measurement, % of deviation not higher than	5	5	5	5	5
Permissible deviation of results in parallel definitions of one sample by different sets of one lot, % of deviation not higher than	10	10	10	10	10

3.2 It is recommended to specify the range of normal values in the population of patients in each laboratory when using reagents. Refer to the Table 13-2 for reference information on reference intervals of indicators for adult patients.

The obtained results of a single blood test can not be used as a single indicator for the diagnosis of diseases but is to be used in conjunction with other laboratory indicators and other diagnostic information.

4. PRECAUTIONS

4.1. Potential risk of use – class 2a.

4.2. Professional staff of the clinical diagnostic laboratories must undergo training in working with hematology analyzers.

4.3. When performing a laboratory test, the precautions taken when handling potentially infectious material should be observed. The main rules are stated in «Instructions on measures to prevent the spread of infectious diseases when working in clinical diagnostic laboratories of medical institutions», approved by the Ministry of health of the USSR 17 Jan 1991 r. And in “Methodology Instructions” MU 287-113.

4.4. When working, use personal protective equipment, at the risk of splashing reagent-safety glasses. To prevent any leakage through the lid, always store the bottles in an upright position. In case of accidental spill, flush the surface with water.

4.5 Possible toxic effects on the human body in case of non-compliance with safety measures and first aid measures:

	If product is inhaled	If product comes into contact with skin	If product comes into contact with eyes	If product is swallowed
Boule Hypochlorite 2% Cleaner First aid measures	If the product is inhaled, and symptoms like shortness of breath or other symptoms of illness occur, fresh air and rest is recommended. If simple first aid does not produce a quick recovery, call a POISON CENTER or doctor/physician.	Wash with soap and water. In contact with chemic	Rinse immediately with a soft or low pressure waterjet or eye wash for at least 5 minutes. Remove contact lenses and keep eyes wide open. If symptoms persist (intense stinging, pain, light sensitivity, poor vision) continue rinsing and call a POISON CENTER or doctor/physician.	Drink a couple of glasses of water. If more than a small quantity has been ingested seek medical advice. If symptoms arise (intense stinging or pain), or if more than a small quantity has been ingested, call a POISON CENTER or doctor/physician.
Boule Hypochlorite 2% Cleaner Information on toxicological effects	Inhalation of vapours may cause irritation. Heating of the product or in case of contact with acids irritating substances, e.g. chlorine, may be released which can cause cough, sore nose and throat, nausea, chest pressure, shortness of breath.	Dry skin, stain may bleach skin.	Eye contact may cause tears, red eyes, pain, blurred vision, impaired but reversible vision.	Ingestion may cause irritation, nausea and vomiting.

	If product is inhaled	If product comes into contact with skin	If product comes into contact with eyes	If product is swallowed	
Boule Enzymatic Cleaner First aid measures	If the product is inhaled, and symptoms like shortness of breath or other symptoms of illness occur, fresh air and rest is recommended. If simple first aid does not produce a quick recovery, call the emergency number.	Wash with soap and water. In contact with chemical substances exposed clothes and shoes should normally be removed. The product does normally not possess any hazard to the exposed person or to first-aiders.	To prevent eye irritation, rinse immediately with a tempered, soft or low pressure water jet or eye wash for at least 5 minutes. If symptoms persist (intense stinging, pain, light sensitivity, poor vision) continue rinsing and seek medical assistance.	Drink a couple of glasses of water. If more than a small quantity has been ingested seek medical advice.	Swelab Alfa Lyse First aid measures
Boule Enzymatic Cleaner Information on toxicological effects	Inhalation of vapors may cause irritation, cough, sore nose and throat.	Not expected	Eye contact may cause tears red eyes, pain, blurred vision, impaired but reversible vision.	Irritation, nausea, vomiting,	Swelab Alfa Lyse Information on toxicological effects
Boule Detergent Cleaner First aid measures	If the product is inhaled, and symptoms like shortness of breath or other symptoms of illness occur, fresh air and rest is recommended. If simple first aid does not produce a quick recovery, call the emergency number.	Wash with soap and water. In contact with chemical substances exposed clothes and shoes should normally be removed. The product does normally not possess any hazard to the exposed person or to first-aiders.	To prevent eye irritation, rinse immediately with a tempered, soft or low pressure water jet or eye wash for at least 5 minutes. If symptoms persist (intense stinging, pain, light sensitivity, poor vision) continue rinsing and seek medical assistance.	Drink a couple of glasses of water. If more than a small quantity has been ingested seek medical advice.	4.6. Instr reagents few 4.7. environmental and hypind on 4.8
Boule Detergent Cleaner Information on toxicological effects	Not expected	Not expected	Tears, red eyes, pain, blurred vision, impaired but reversible vision.	Irritation, nausea, vomiting	Fc
Swelab Alfa Diluent First aid measures	If the product is inhaled, and symptoms like shortness of breath or other symptoms of illness occur, fresh air and rest is recommended. If simple first aid does not produce a quick recovery, call the emergency number.	Wash with soap and water. In contact with chemical substances exposed clothes and shoes should normally be removed. The product does normally not possess any hazard to the exposed person or to first-aiders.	To prevent eye irritation, rinse immediately with a tempered, soft or low pressure water jet or eye wash for at least 5 minutes. If symptoms persist (intense stinging, pain, light sensitivity, poor vision) continue rinsing and seek medical assistance.	Drink a couple of glasses of water. If more than a small quantity has been ingested seek medical advice.	- 1. F
Swelab Alfa Diluent Information on toxicological effects	Not expected.	Not expected.	Tears, red eyes, pain, blurred vision, impaired but reversible vision.	Irritation, nausea, vomiting	

	If product is inhaled	If product comes into contact with skin	If product comes into contact with eyes	If product is swallowed
Swelab Alfa Lyse First aid measures	If the product is inhaled, and symptoms like shortness of breath or other symptoms of illness occur, fresh air and rest is recommended. If simple first aid does not produce a quick recovery, call the emergency number.	Wash with soap and water. In contact with chemical substances exposed clothes and shoes should normally be removed. The product does normally not possess any hazard to the exposed person or to first-aiders.	To prevent eye irritation, rinse immediately with a tempered, soft or low pressure water jet or eye wash for at least 5 minutes. If symptoms persist (intense stinging, pain, light sensitivity, poor vision) continue rinsing and seek medical assistance.	Drink a couple of glasses of water. If more than a small quantity has been ingested seek medical advice.
Swelab Alfa Lyse Information on toxicological effects	Not expected	Not expected.	Tears, red eyes, pain, blurred vision, impaired but reversible vision.	Irritation, nausea, vomiting

4.6. When working with the medical device, follow local procedures SP 1.3.2322, GOST R 52905. Instruments and equipment, as well as the surfaces on which the work was carried out using the reagents, should be treated with a solution of ethyl alcohol with a volume fraction of 70 % within a few seconds.


4.7. The Hypochlorite 2% Cleaner contains sodium hypochlorite in a concentration that is health and environmentally hazardous. The bottle is marked with a hazard pictogram. Sodium Hypochlorite 2% is skin irritating and causes serious eye irritation. The solution shall be handled with care, and skin and eyes shall be protected if there is a risk of splashing. For more information, read SDS. Since hypochlorite solutions are oxidizing with a basic pH, avoid contact with materials other than those indicated in the instructions (spills can leave a remaining smell, e.g. on skin, and stains (bleaching) on clothes etc.)

4.8. Close containers with cleaning reagents immediately after use.

5. EQUIPMENT AND MATERIALS

For performing CBC analysis of the whole blood, it's recommended to use:

- 1) Automated hematology analyzers Swelab Alfa, Swelab Alfa Plus;
- 2) Tubes, microtubes and micropipettes with anticoagulants which are listed in User's Manual for analyzers and Section 6.1;
- 3) Use personal protective equipment;
- 4) Control materials certified for 16 parameters of whole blood for hematology analyzers Swelab - Boule control materials for hematology analyzers (hematological control blood 16 parameters, levels low, normal and high) and hematological calibrator produced by Boule Medical AB, Sweden, RU № FSZ 2008/02770).
- 5) For the cleaning procedure: distilled water, vial or empty tube 10 ml, absorbent soft cloth/material.

6. SAMPLES ANALYSED

6.1. Samples used are whole human blood (capillary or venous), stabilized by anticoagulant K2 EDTA or K3-EDTA.

In order to obtain accurate results and prevent blockage of the hydraulic system of the analyzer samples used must be collected and prepared in accordance with the respective Swelab Alfa and Swelab Alfa Plus hematology analyzer User Manual (sections 3-4).

6.2. Samples in the EDTA tubes can be stored less than 4 h at 18 - 25 °C or less than 24 h from 2 to 8° C according to [1].

6.3. The results:
Blood e
be stirr
6.4. Re
compa
the ide

7.1 A
User's
Dilue
7.2. P
to use
7.3. '
analy
wher
7.4.
cont

The
Cal
list
leu

9
9.
Cl
Sv
9.
9.
re
1
2
3
9

6.3. The following cases were detected in which the type of anticoagulant affected the analytical results:

Blood anti coagulated with EDTA: Samples of venous blood collected in test tubes of EDTA, must be stirred cautiously and retained for 15 minutes before testing.

6.4. Results of the clinical tests performed using Swelab hematology analyzers and reagents, confirm comparability of the results obtained with the use of venous and capillary blood samples, as well as the identity of the results with the samples of anticoagulants K2EDTA and K3EDTA.

7. PREPARATION OF COMPONENTS FOR TESTING

7.1 A detailed description of reagent installation / reagent change can be found in the instrument User's Manual under section 2 – Installation and Reagent Setup. It is recommended that both the Diluent and the Lyse reagents are placed at the instrument level or below.

7.2. Place the tag reader or barcode reader on the marked position on the reagent container you want to use. When the reagent tag / barcode has been read a screen will display that the it has been accepted.

7.3. Transfer the reagent level sensor from the used container to the new reagent container. The analyzer is now ready to resume operation or analyze samples. No priming or fill cycle is necessary when putting on a new reagent container, if indicator and warning messages are followed.

7.4. IMPORTANT! Remember to verify calibration prior to running patient samples using Boule control material. If it's necessary, start calibration.

8. SAMPLE ANALYSIS AND CALCULATIONS

The analysis of whole blood should be carried out according to the instruments User's Manual. Calibration and calculations are performed automatically using the analyzer software. The complete list of blood test parameters determined by a set of parameters (erythrocyte indices, platelet indices, leukocyte formula and hemoglobin concentration) depends on the analyzer model used.

9. CLEANING AND MAINTAINANCE PROCEDURES USING CLEANING REAGENTS

9.1. A detailed description of use of the Boule Cleaners (Clot prevention, Monthly cleaning, Boule Cleaning Kit procedure) can be found in Section 8 of Swelab Alfa User Manuals and Section 10 of Swelab Alfa Plus the instrument User Manual and in the Quick Reference Guides.

9.2. Boule Cleaning Kit procedure.

9.2.1. It is very important that the cleaning procedure is performed in the correct order (or it may result in high PLT background counts). The following order must be followed:

- 1) Enzymatic Cleaner (enzyme solution - blue cap)
- 2) Hypochlorite (2% solution - white cap)
- 3) Detergent Cleaner (detergent solution - red cap)

9.2.2. Follow the procedure described in

- Section 8; Swelab Alfa User Manuals
- Section 10; Swelab Alfa Plus User Manual
- Quick Reference Guides; Swelab Alfa and Swelab Alfa Plus

Notes - As the volume in the cleaning bottles will decrease during the cleaning cycles, be sure to secure the bottles to avoid tilting. - During the described cleaning process above, some foaming might occur on the aspiration needle. Remove foam with a soft tissue. - In case of heavy contaminated systems, it is recommended to warm the detergent to 45-55 °C prior to cleaning the system. - Only the Hypochlorite Cleaner has the ability to act as a disinfectant, significantly reducing the number of viable microorganisms. - It is recommended to check calibration by running a control after cleaning

10. TERMS OF USE OF THE REAGENTS

10.1. Reagents should be used according to the Instruction for Use approved in the prescribed manner.

10.2. Swelab AlfaDiluent and Swelab AlfaLyse can be used within 90 days after opening. By the addition of antimicrobial agent, microbial growth is prohibited during storage as well as during use. Avoid the ingress of dust into the container and tubes. Accessory caps for reagent containers are designed with holes fitting for the reagent probes. This keeps the reagents sealed and preventing particles from entering the container while being used.

The open cleaner container is stable at least up to the expiration date as shown on the label when stored at a temperature of 4-30 °C.

10.3. In case of contamination of reagents, signs of deterioration, as well as when the background value exceeds the limits specified in the analyzer user Manual, dispose of reagents as indicated below.

10.4. Quality control shall be performed with Boule control materials (see section 5 of current IFU) in accordance with Russian GOST R 53133.2-2008 [2] and instrument User's Manual (section 5).

10.5. If the results of quality control that exceed the permissible limits of the control plot of the analyzer or the parameter limits of the control material, it is necessary to check the analyzer and make a new calibration with the hematological calibrator.

10.6. Adulteration by dilution or addition of any materials to the Swelab reagents and cleaners invalidates any diagnostic use of the product.

10.7. Reagents are used at 18 to 32 °C; Avoid direct sunlight.

11. TRANSPORTATION AND STORAGE OF THE REAGENTS

11.1. Reagents must be transported by all types of indoor transport in compliance with the conditions and requirements set for this type of transport, at a temperature of 4–35 °C for Swelab AlfaDiluent, Swelab AlfaLyse and 4–30 °C for cleaning reagents. Avoid overturning containers. Freezing of kit components is not allowed.

11.2. Storage in the original cardboard packs at 4–35 °C for Swelab AlfaDiluent, Swelab AlfaLyse and 4–30 °C for Cleaners until expiration date. Do not freeze.

12. DISPOSAL OF THE REAGENTS

12.1. The product should be disposed of as Class B medical waste in accordance with the applicable laws of the Russian Federation (Sanitary Regulation and Standard (SanPiN) N2.1.7.2790-10 dated 12 December 2010 "Sanitary and Epidemiological Requirements on the Handling of Medical Waste").

12.2 Before disposing of used reagent, 2 ml of hypochlorite solution per liter of waste product must be added and the solution retained for 1 hour. It should then be disposed of as Class B medical waste.

12.3 Any unused reagent (including any that is out-of-date) should be disposed of as Class B medical waste. Liquid components (including expired) are destroyed by dilution with 1:100 water, draining the solution into the sewer and removing the remainder of the packages as industrial or household garbage [3].

13. WARRANTIES OF THE MANUFACTURER

13.1. Warranties extend to all characteristics indicated on the label and in the Instructions for Use throughout the entire shelf life, provided that the conditions for use, storage and shipping are complied with.

Shelflife for AlfaLyse and AlfaDiluent is 36 months, for Cleaners – 24 months.

13.2. Manufacturer information: Boule Medical AB, Domnarvsgatan 4, SE-163 53 Spånga, Sweden
E-mail: info@boule.com Website: www.boule.com.

Manufacturing site: Boule Medical AB, Fagerstagatan 7, SE-163 53 Spånga, Sweden.

13.3. **Claims on the quality of the product in the territory of the Russian Federation should be sent to:** LLC "Boule Medical", Russian Federation, 125047, Moscow, 4th Lesnoy Lane, b.4, office 519. E-mail: boule.registration@mail.ru, Website: www.boulemedical.ru.

This IFU is also available on <http://www.boulemedical.ru>.

The medical device is allowed for circulation in the territory of the Russian Federation.

RC № FSZ 2008/01097 of

Table 13-1. Symbols on the labels and package:

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I, the undersigned, Gabriel Rudbeck, Notary Public in
Stockholm, Sweden, hereby certify that

Jesper Filip Gustaf Busck Söderqvist---

has/have issued and signed the foregoing document.

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