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ATTESTATION STATEMENT

As a responsible member of Medtronic MiniMed, I hereby attest that I believe the Appendix to IFU, attached herein, to be a true and accurate copy.

Katherine Rojas
Associate RA Specialist - International

Katherine Rojas
Signature

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document


State of California
County of Los Angeles

On 08/07/2020 before me, Yvette Lee: Notary public
(date) (insert name and title of the officer)

personally appeared Katherine Rojas, who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument. I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature *[Signature]* (Seal)


Sensor Guardian Sensor (3) in variants

Appendix to IFU

Name of medical device	Guardian™ Sensor (3) percutaneous glucose monitoring sensor in variants: 1. Guardian™ Sensor (3) percutaneous glucose monitoring sensor, catalog number MMT-7020D3 including: - Guardian™ Sensor (3) percutaneous glucose monitoring sensor – 1 p. - Adhesive patch – no more 2 p. 2. Guardian™ Sensor (3) percutaneous glucose monitoring sensor, catalog number MMT-7020C3 including: - Guardian™ Sensor (3) percutaneous glucose monitoring sensor – 5 p. - Adhesive patch – no more 10 p.	
Prescription of medical device	A sterile product for individual use, designed to continuously measure the concentration of glucose in the extracellular fluid in order to determine the trend and monitor the glucose level of a patient with diabetes mellitus.	
Indications and Contraindications for medical device	Indications: the sensor is intended for use with Guardian Link (3) and Guardian Connect transmitters for continuous monitoring of glucose levels in people with diabetes all ages. Contraindications: none.	
Information on potential consumers of medical device, use conditions	Intended for use by patient or physician. Intended for use by adult and children both in inpatient and home conditions. The indication for children under the age of 18 is limited to those who are supervised by a caregiver. The caregiver must be at least 18 years of age. The performance of the sensor has not been evaluated for pregnant women or persons on dialysis.	
Description of the main functional parts of medical device	Permissible tolerances ±10%, unless otherwise indicated	
	Guardian Sensor (3) glucose monitoring sensor	
		
	Mass, g	14,75
	Volume, cm³	5
	Volume (in package), cm³	86

Dimensions (LxWxH), mm Length and width are of the sensor pedestal only, height is of the sensor pedestal and the needle hub (not including the needle)	43.3 mm x 63.9 mm x 61.9 mm
Cutting angle of needle, °	45
Effective length of needle, mm, not more	13
Needle outer diameter, in (needle tube), mm	0.03 ± 0.001
Needle inner diameter, in (needle tube), mm	0.0155 ± 0.0005
Needle Cross Section Shape	
Needle head color	Uncolored (steel)
Puncture depth	No more than 13 mm
Surface roughness	0.19 and 0.21 microns
Rockwell Hardness, HRC	Not less than 70
The force required to introduce the sensor by a person, N (not to exceed)	5
Durability of connection of a needle and a sensor	An applied force exceeding 10LbF shall be required to separate the sensor needle from the plastic hub assembly. An applied force exceeding 0.13 lbf shall be required to break the sensor from the sensor base.
Needle condition (defects, scratches, cracks, foreign particles)	Needle shall be sharp with no round, flat or bent edges at 10X magnification. Parts inspected at 5X magnification shall be free of adhered or embedded particles, dried surface film or visible oil. Parts shall be burr-free on needle tip facets and on all edges with the effective length. A 0.0003 inches maximum burr is allowable on all other areas.
Invasiveness	Inserted into the user's subcutaneous tissue. The profile (thickness) of the sensor assembly as worn on the body (measured from the bottom of the patch to the highest point on the base) shall be no larger than 0.32 inches.
Multiplicity of application	Single-use
Sensor water resistance	Watertight (IPX8)
Sensor Care	Only store sensors at room temperature between 2 °C to 27 °C (36 °F to 80 °F). Discard sensor after the "Use by date"

		indicated on the label, if the package is damaged, or the seal is broken.											
Sensitivity and specificity of the sensor in relation to blood glucose level, accuracy of sensor measurement, coefficient of variation	Sensors were used in both abdominal and arm sites. The overall mean and median absolute relative difference (MARD) was calculated for sensors located in both the abdomen and arm areas compared to a Yellow Springs Instrument (YSI) reference. <table><tr><th colspan="3">Absolute relative difference (%) compared to YSI</th></tr><tr><th></th><th>Mean</th><th>Median</th></tr><tr><td>Abdominal</td><td>10.5</td><td>7.7</td></tr><tr><td>Arm</td><td>9.1</td><td>6.9</td></tr></table>	Absolute relative difference (%) compared to YSI				Mean	Median	Abdominal	10.5	7.7	Arm	9.1	6.9
Absolute relative difference (%) compared to YSI													
	Mean	Median											
Abdominal	10.5	7.7											
Arm	9.1	6.9											
Sterilization mode	Sterility assurance level (SAL) of 10 ⁻⁶												
Sterilization method	Radiation Sterilization												
Adhesive patch													
													
Dimensions (LxW), inches (mm)	38.1 mm x 12.7 mm (from flat edge to oval opening)												
Mass, g	0,7												
Patch type	fixing												
Patch design and shape	Tape with 2 pieces of liner (paper) and an adhesive, oval shape												
Functional properties of the patch	Non-sterile, single-use, disposable adhesive tape designed to aid in adhering a CGM sensor or transmitter/recorder onto the patient. Does not contain medicines or biologically active compounds.												
Type of adhesive	Nylon Reinforced Polyurethane Film with medical grade acrylic adhesive (synthetic adhesive)												
It does not have a pungent smell of other chemical compounds, except for the adhesive layer	Match												

	The sticky layer is colorless or white	Match colorless
	It does not have the following external defects: unevenly cut edges, warp folds, winding defects, holes, foreign matter, lumps, traces of uneven impregnation with a plaster mass and places uncovered with a plaster mass	Match Tape must be supplied clean and free of loose particulates. Tape is inspected for general workmanship and proper construction, and is verified for contamination.
	Protected by release coating: The protective coating is removed without effort; patches of plaster do not remain on the surface of the protective coating.	Match
	Non-toxic and does not cause local irritating and allergic effects when the sticky layer (adhesive) comes into contact with intact skin and a functional pad with a wound surface	Match The tape has been evaluated and complies with EN ISO 10993-1 and FDA Guidance for biocompatibility. The tape passed cytotoxicity (ISO 10993-5:2009) and sensitization/irritation (ISO 10993-10:2010) testing.
	Resistance to peeling of the sticky layer, N / m (not less than 25 - not more than 200)	The peel force of the adhesive tape shall be between 0.155 lbf and 0.770 lbf when tested per AP2380.
	The sorption capacity of the functional pad, the number of drops / cm (at least 1) (only for patch dressings with sorption pad, including atraumatic, antimicrobial, bactericidal, hemostatic)	Not applicable (not a sorption pad)
	The wetting time of the functional pad, s, not more than 5 (for dressings of a plaster type with a sorption pad, including with atraumatic, antimicrobial, bactericidal, hemostatic)	Not applicable (not a sorption pad)
	The vapor permeability of the fixing layer, mg / cm / h, not less than 1.5 (vapor-permeable for dressings of a plaster type)	Not applicable (not a plaster-type)
	The water resistance of the fixing layer, mm water column, at least 300 (for dressings of the plaster type "waterproof" (with a waterproof adhesive layer))	Not applicable (not waterproof)
	The degree of adhesion to the model of the wound surface, % to medical gauze, not more than 30 (for atraumatic dressings of a plaster type)	Not applicable (not atraumatic)
	Sterility (For sterile patch dressings in sealed packaging)	Non-sterile
	Tightness (For sterile adhesive dressings in sealed packaging)	Not applicable (non-sterile)

Medtronic

Utilization	Discard the sensor according to local regulations for disposal, or contact your healthcare professional for disposal information. Medical waste class: B	
Adverse events (Potential complications)	General risks with sensor use include: • Skin irritation or other mild adverse skin reactions • Bruising • Discomfort • Redness • Bleeding • Pain • Rash • Infection • Raised bump • Appearance of a small "freckle-like" dot where needle was inserted • Allergic reaction • Fainting secondary to anxiety or fear of needle insertion • Soreness or tenderness • Swelling at insertion site • Sensor fracture, breakage or damage • Minimal blood splatter associated with sensor needle removal • Residual redness associated with adhesive or tapes or both • Scarring	
All-Russian classification of products by type of economic activity (OKPD2)	32.50.50.190 "Medical devices, including surgical, other, not included in other groups".	
Warranty	Manufacturer guarantees compliance of the device with standards and product requirements subject to instructions of use, transport, storage and utilization regulations Shelf life (expiration date) is 12 months.	
The list of national and international regulatory documents/standards, which matches medical device	Standard	Title
	(EC) 1272/2008	THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006
	2011/65/EU	Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment (ROHS)
	(EC) 1907/2006	Regulation (EC) No 1907/2006 of the European Parliament and of the Council of 18 December 2006 concerning the Registration, Evaluation, Authorization, and Restriction of Chemicals (REACH), establishing a European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No 793/93 and Commission Regulation (EC) No 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC, and 2000/21/EC.

EN ISO 13485:2012	Medical Devices – Quality Management Systems – Requirements for Regulatory Purposes (ISO 13485:2003)
EN ISO 14971:2012	Application of Risk Management to Medical Devices
EN 1041:2008	Information Supplied by the manufacturer with Medical Devices
EN 556-1:2001/AC: 2006	Sterilization of Medical Devices – Requirements for Medical Devices to be Designated ‘STERILE’ – part 1: Requirements for Terminally Sterilized Medical Devices
EN ISO 15223-1:2016	Medical Devices – Symbols to be used with Medical Device Labels, Labelling, and Information to be Supplied – Part 1: General Requirements (ISO 15223-1:2016)
EN 10993-1:2009/AC: 2010	Biological Evaluation of Medical Devices. Part 1: Evaluation and Testing with a Risk Management Process (ISO 10933-1:2009) CORRIGENDUM July 31, 2010
EN ISO 10993-11:2009	Biological Evaluation of Medical Devices Part 11: Tests for Systemic Toxicity (ISO 10993-11:2006)
EN ISO 10993-12:2012	Biological Evaluation of Medical Devices Part 12: Samples Preparation and Reference Materials (ISO 10993-12:2012)
EN ISO 10993-17:2009	Biological Evaluation of Medical Devices Part 17: Establishment of Allowable Limits for Leachable Substances (ISO 10993-17:2002)/
EN ISO 10993-18: 2009	Biological Evaluation of Medical Devices Part 18: Chemical Characterization of Materials (ISO 10993-18: 2005)
EN ISO 10993-3: 2014	Biological Evaluation of Medical Devices, Part 3: Tests for Genotoxicity, Carcinogenicity, and Reproductive Toxicity.
EN ISO 10993-5: 2009	Biological Evaluation of Medical Devices, Part 5: Tests for in vitro Cytotoxicity (ISO 10993-5:2009)
EN ISO 10993-10:2014	Biological Evaluation of Medical Devices Part 10: Tests for Irritation and Skin Sensitization
EN ISO 10993-6: 2016	Biological Evaluation of Medical Devices Part 6: Tests for Local Effects After Implantation.
EN ISO 11137-1:2015	Sterilization of Health Care Products - Radiation Part 1: Requirements for Development, Validation, and Routine Control of a Sterilization Process for Medical Devices (ISO 11137-1: 2006 Amendment 1:2013)
EN ISO 11137-2: 2015	Sterilization of Health Care Products – Radiation – Part 2: Establishing the Sterilization dose (OSO 11137-2:2013)
EN ISO 11737-1: 2006/AC: 2009	Sterilization of Medical Devices – Microbiological Methods – Part 1: Determination of a Population of Microorganisms on Products
EN ISO 11737-2:2009	Sterilization of Medical Devices – Microbiological Methods – Part 2 Tests of Sterility Performed in the Definition, Validation, and Maintenance of a Sterilization Process.
EN ISO 14155:2011	Clinical Investigation of Medical Devices for Human Subjects. Good Clinical Practice (ISO 14155:2011)
EN ISO 14644-1:2015	Cleanrooms and Associated Controlled Environments – Part 1: Classification of Air Cleanliness

Medtronic

	EN ISO 14644-:2015	Cleanrooms and Associated Controlled Environments – Part 2: Specifications for Testing and Monitoring
	EN ISO 14644-3:2005	Cleanrooms and Associated Controlled Environments – Part 3: Test Methods
	EN ISO 14644-4:2001	Cleanrooms and Associated Controlled Environments – Part 4: Design Construction and Start-up
	EN ISO 11607-1: 2009 + A1:20014	Packaging for Terminally Sterilized Medical Devices Part 1: Requirements for Materials, Sterile Barrier Systems and Packaging Systems
Customer Service	Contact for all questions, including complaints, repair and service of medical devices the Authorized Representative in Russia: Medtronic, LLC, 123112, Russia, Moscow, Presnenskaya Naberezhnaya., h.10, pom. III, ap. 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru	

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© 2018 Medtronic Medinol, Inc. Все права защищены.