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NIALL SHEEHY

NOTARY PUBLIC

6/7, Glentworth St.,
Limerick
Ireland

Telephone +353 61 414355


Facsimile + 353 61 414738

Email nsheehy@sellors.ie

TO ALL TO WHOM THESE PRESENTS SHALL COME, I, Niall Sheehy, Notary Public, duly authorised, admitted, and sworn, practicing at 6/7 Glentworth St., Limerick, Republic of Ireland, DO HEREBY CERTIFY that I am familiar with documents as issued from time to time by Cook Ireland Ltd and that the original document hereunto annexed namely, **Instructions for Use-EVOLUTION COATED ESOPHAGAL NITINOL STENT-ENDOPROTHESIS** upon which I have signed my name and affixed my Seal of Office is a true and faithful original and that the original appears to me to be genuine.

AND I HAVE CERTIFIED as herein at the request of Cook Ireland Ltd having its place of business at O'Halloran Road, National Technology Park, in the city of Limerick, Ireland, and for the benefit of the proper authorities in **Russian Federation** and for all to whom it shall come.

IN FAITH AND TESTIMONY whereof I have hereunto set my hand and affixed my Seal of Office.

Niall Sheehy
Notary Public

Date

NIALL SHEEHY
Notary Public
6/7 Glentworth St., Limerick.
Commissioned For Life

APOSTILLE (Convention de La Haye du 5 octobre 1961)			
1. Country: Pays/País:		IRELAND	
This public document Le présent acte public / El presente documento público			
2. has been signed by a été signé par ha sido firmado por		Niall Sheehy	
3. acting in the capacity of agissant en qualité de quien actúa en calidad de		Notary Public	
4. bears the seal / stamp of est revêtu du sceau / timbre de y está revestido del sello / timbre de			
Certified Attesté / Certificado			
5. at à / en	Cork	6. the le / el día	16/08/2021
7. by par / por	Department of Foreign Affairs		
8. No sous no bajo el número	0909742021		
9. Seal / stamp: Sceau / timbre: Sello / timbre:	10. Signature: Signature: Firma: Mary Kelleher		
This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. To verify the issuance of this Apostille, see www.authentication.dfat.ie			327419

INSTRUCTIONS FOR USE

EVOLUTION COATED ESOPHAGAL NITINOL STENT-ENDOPROTHESIS

PURPOSE:

This device is used to maintain patency of malignant esophageal strictures and/or to seal tracheoesophageal fistulas.

The Evolution Esophageal Stent System is represented by the versions mentioned below:

-Evolution Esophageal Stent System - Partially Covered (EVO-20-25-8-E; EVO-20-25-10-E; EVO-20-25-12.5-E; EVO-20-25-15-E)

-Evolution Esophageal Stent System - Fully Covered (EVO-FC-R-18-23-8-E; EVO-FC-R-18-23-10-E; EVO-FC-R-18-23-12-E; EVO-FC-R-20-25-8-E; EVO-FC-R-20-25-10-E; EVO-FC-R-20-25-12-E)

The Evolution Esophageal Stent System consists of the following two components:

- Self Expanding Metal Stent
- Stent Delivery System (made up of the handle and introducer)

Stent Delivery System

The Stent Delivery System is common to the aforementioned devices - it is very similar to traditional 'Push-Pull' stent delivery systems. It consists of two catheters that operate co-axially to one another (inner and outer catheter) and a handle. The handle allows for desheathing to deploy the stent and resheathing recapturing the stent during stent deployment.

Differences in the construction of the Self Expanding Metal Stents in the aforementioned versions are described below:

Evolution Esophageal Stent System - Partially Covered (EVO-20-25-8-E; EVO-20-25-10-E; EVO-20-25-12.5-E; EVO-20-25-15-E)

The stent is woven from a single strand of nitinol wire. The stent foreshortens due to its design. It is formed with a flange at either end. The increased diameter of the stent at the flanged ends provides resistance to migration. The proximal and distal ends of the stent are designed with rounded "crowns", to potentially reduce tissue trauma. A silicone membrane applied to the stent covers the complete body of the stent along with part of both flanges. The covering is intended to reduce the risk of tissue ingrowth and provides a seal for tracheoesophageal fistulas. To aid in visibility under fluoroscopy there are four tantalum bands at either end of the stent. The device has a lasso loop woven through the crowns on the proximal end of the stent (first lasso loop). A second lasso loop is woven through the loops formed by the first lasso loop. This lasso loop assembly can be used to reposition the stent during the initial placement procedure if desired. (Fig. 1)

Niall Sheehy
NIAL SHEEHY
Notary Public
6/7 Glentworth St., Limerick.
Commissioned For Life

09-08-2021


Figure 1: Evolution Esophageal Stent- Partially Covered (EVO-20-25-8-E; EVO-20-25-10-E; EVO-20-25-12.5-E; EVO-20-25-15-E)

Evolution Esophageal Stent System - Fully Covered (EVO-FC-R-18-23-8-E; EVO-FC-R-18-23-10-E; EVO-FC-R-18-23-12-E; EVO-FC-R-20-25-8-E; EVO-FC-R-20-25-10-E; EVO-FC-R-20-25-12-E)

The stent is woven from a single strand of nitinol wire. The stent has flanges at both ends of the stent. The increased diameter of the stent at the flanged ends is intended to provide resistance to migration. The stents have a silicone membrane covering. For this version of the stent - a silicone membrane applied to the stent covers the complete stent except for sections of every crown (at both ends of the stent) where a hole is left open to allow the lasso loop to pass through. The covering is intended to reduce the risk of tissue ingrowth and to provide a seal for tracheoesophageal fistulas. To aid in visibility under fluoroscopy, there are four tantalum marker bands at either end of the stent. The device at both stent ends has a lasso loop woven through the crowns on the end of the stent (first lasso loop). A second lasso loop (grasping lasso loop) is woven through the loops formed by the first lasso loop and contains a green grasping feature. The green grasping feature is gripped to allow repositioning of the stent directly after placement or is used for stent removal in benign indications. (reference Fig. 2).


Figure 2: Evolution Esophageal Stent-Fully Covered (EVO-FC-R-18-23-8-E; EVO-FC-R-18-23-10-E; EVO-FC-R-18-23-12-E; EVO-FC-R-20-25-8-E; EVO-FC-R-20-25-10-E; EVO-FC-R-20-25-12-E).

Construction

1 Fr=0.33

Part Number	Introduction System Diameter	Introduction System Length	Stent Length	Stent Diameter
Evolution Esophageal Stent System - Partially Covered				
EVO-20-25-8-E	24 Fr/ 7.92 mm (±10%)	78cm ± 3cm	8cm (±10%)	Body – 20mm (±10%), Flange – 25mm (±10%)
EVO-20-25-10-E	2 4Fr/ 7.92 mm (±10%)	78cm ± 3cm	10cm (±10%)	Body – 20mm (±10%), Flange – 25mm (±10%)

Part Number	Introduction System Diameter	Introduction System Length	Stent Length	Stent Diameter
EVO-20-25-12.5-E	24 Fr/ 7.92 mm (±10%)	78cm ± 3 cm	12.5cm (±10%)	Body – 20mm (±10%), Flange – 25mm (±10%)
EVO-20-25-15-E	24 Fr/ 7.92 mm (±10%)	78cm ± 3 cm	15cm (±10%)	Body – 20mm (±10%), Flange – 25mm (±10%)
Evolution Esophageal Stent System - Fully Covered				
EVO-FC-R-18-23-8-E	24 Fr/ 7.92 mm (±10%)	78cm ± 3 cm	8cm (±10%)	Body – 18mm (±10%), Flange – 23mm (±10%)
EVO-FC-R-18-23-10-E	24 Fr/ 7.92 mm (±10%)	78cm ± 3 cm	10cm (±10%)	Body – 18mm (±10%), Flange – 23mm (±10%)
EVO-FC-R-18-23-12-E	24 Fr/ 7.92 mm (±10%)	78cm ± 3 cm	12 cm (±10%)	Body – 18mm (±10%), Flange – 23mm (±10%)
EVO-FC-R-20-25-8-E	24 Fr/ 7.92 mm (±10%)	78cm ± 3 cm	8 cm (±10%)	Body – 20mm (±10%), Flange – 25mm (±10%)
EVO-FC-R-20-25-10-E	24 Fr/ 7.92 mm (±10%)	78cm ± 3 cm	10 cm (±10%)	Body – 20mm (±10%), Flange – 25mm (±10%)
EVO-FC-R-20-25-12-E	24 Fr/ 7.92 mm (±10%)	78cm ± 3 cm	12 cm (±10%)	Body – 20mm (±10%), Flange – 25mm (±10%)

Allowance / tolerances of the dimensional parameters

Unit: mm Tolerance: ±10%

Stent Length *	Silicone Membrane length	Stent Diameter	
		Stent Body Diameter	Stent Flange Diameter
80	50	20	25
100	70		
125	95		
150	120		

* The length indicates the actual length of the stent at nominal stent diameter.

Unit: mm Tolerance: ± 10%

Stent Length * \ Silicone Membrane Length	Stent Diameter	
	Stent Body Diameter	Stent Flange Diameter
80	18	23
100		
120		
80	20	25
100		
120		

* The Length indicates the actual stent length at nominal stent diameter.

Unit: mm Tolerance: ± 10%

Introducer System Tip Length	Introducer System Tip ID	
F	G1	G2
51	1.83	1.0


Figure 3: Introducer System Tip (Inner Catheter)

Unit: mm Tolerance: $\pm 10\%$

The length of the outer catheter
758

Change of length

The length of the stent must be $+15\%$ / -5% of its nominal dimension post deployment.

All test articles met the predetermined acceptance criteria, demonstrating an accurate stent deployment.

Radial Resistance at 37° C

Radial Force Testing

Acceptance criteria – the stent radial force must be $> 20\text{N}$ at a stent diameter of 10mm (pass: min 29.22N, max 46.52N).

Radial force testing was conducted on 8cm stents. The testing was carried out on the machine solutions radial expansion force gauge tester. The temperature within the tester was maintained at $37^\circ\text{C} \pm 2^\circ\text{C}$. Radial force at a diameter of 10mm was measured. Each stent was compressed from an initial, fully expanded diameter to a final diameter of 9mm at a rate of 0.5mm/s. The test article was allowed to expand back to the initial diameter at the same rate. Three compression expansion cycles were performed. Forces were established from the 3rd cycle.

Deployment & Recapturing Force Testing

To allow the testing to be completed the handle was removed from the test articles. The position of the “point-of-no-return” position and the 75% deployment indicator were marked on the test articles. The testing took place in a water bath set to water bath to $37^\circ\text{C} \pm 2^\circ\text{C}$ (this simulated body conditions). A 0.035” wire guide was introduced through the delivery system. The test article was attached to the force gauge using the luer lock.

Deployment Force:

The force gauge was set to record maximum compressive load. Using a slow uniform motion, the shuttle was moved towards the force gauge until 50% of the stent was deployed. This position corresponds to the “point-of-no-return” position. The maximum force was recorded. Note: This corresponds to stent deployment.

Recapturing Force:

The force gauge was set to record maximum tensile load. Using a slow uniform motion, the shuttle was moved away from the force gauge from the “point-of-no-return” position until the stent is fully recaptured. The maximum force was recorded. Note: This test corresponds to stent recapturing.

As part of this test each test article was recaptured a total of 5 times.

Results and Acceptance Criteria:

All the test articles met the acceptance criteria outlined in the following table.

Table: Deployment and recapturing forces: Acceptance criteria and results.

Test Type	Acceptance Criteria	Results	Pass/ Fail
Deployment Force Testing	The stent deployment force must be below 40N when deployed from a full length delivery system into a water bath heated to $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$ after the stent has been recaptured up to 5 times.	<i>Deployment Forces:</i> Min: 19.15N Max: 31.35N Note: all deployment values for each of the 6 cycles were <40N.	Pass
Recapturing Force Testing	The stent recapture force must be below 40N when recaptured up to 5 times into a full length delivery system using water heated to $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$.	<i>Recapture Forces:</i> Min: 24.10N Max: 31.70N Note: all recapture values for each of the 5 cycles were <40N.	Pass

Connection Strength

- 1. Distal Tip Joint Strength**
 - Mean -3.981 SD. > 15N
- 2. Marker Band Joint Strength**
 - Mean -3.981 SD. > 5N
- 3. Bilumen Tubing Joint Strength**
 - Mean-3.064 x standard deviation > 45N
- 4. Outer Tubing Flare to connector cap Joint strength**
 - Mean – 3.064 x standard deviation > 5kg
- 5. PEEK Tubing Flare Joint Strength**
 - Force at failure for all samples > 45N
- 6. PERT Tubing to Handle Cannula Joint Strength**
 - Force at failure for all samples > 45N
- 7. Cannula to Nitinol Wire Joint strength**
 - Mean – 3.064 x standard deviation > 15N

Fatigue Testing

Objective: To evaluate the fatigue life of the Evolution Esophageal Self Expanding Metal Stent by subjecting it to time-accelerated, physiologically-modelled, diameter-controlled radial compression for a specified number of cycles.

Number of Samples: A total number of ten (10) stents were used in this study (Five 18 mm body x 23 mm flange diameter and five 22 mm body x 27 mm flange diameter). The 2 sizes tested bracket the 20mm body x 25mm flange size option.

Acceptance Criteria and Results:

Table: Acceptance Criteria and Results: Stent 18mm Body, 23mm Flange.

Test Article	Acceptance Criteria	Result
Evolution Esophageal Fully Covered Stent (18 mm body x 23 mm flange diameter) n= 5	The stents must remain intact.	Pass
	The stents must keep their cylindrical geometry.	Pass
	The stent wires must not break during the test.	Pass
	The marker bands must remain on the stents and not become detached.	Pass
	The lasso loop material shall not break during the test.	Pass
	The silicone covering must remain on the stent and not become detached from the body of the stent. Some localized delamination or holes in the silicone are acceptable.	Pass

Table: Acceptance Criteria and Results: Stent 22mm Body, 27mm Flange.

Test Article	Acceptance Criteria	Result
Evolution Esophageal Fully Covered Stent (22 mm body x 27 mm flange diameter) n= 5	The stents must remain intact.	Pass
	The stents must keep their cylindrical geometry.	Pass
	The stent wires must not break during the test.	Pass
	The marker bands must remain on the stents and not become detached.	Pass
	The lasso loop material shall not break during the test.	Pass
	The silicone covering must remain on the stent and not become detached from the body of the stent. Some localized delamination or holes in the silicone are acceptable.	Pass

Analysis: The test samples met all of the predetermined acceptance criteria demonstrating acceptable radial sling fatigue.

All of the above results apply to all stent designs: Evolution esophagus stent system - partially coated (EVO-20-25-8-E; EVO-20-25-10-E; EVO-20-25-12.5-E; EVO-20-25-15-E); The Evolution esophagus stent system is fully coated (EVO-FC-R-18-23-8-E; EVO-FC-R-18-23-10-E; EVO-FC-R-18-23-12-E; EVO-FC-R-20-25-8-E; EVO-FC-R-20-25-10-E; EVO-FC-R-20-25-12-E).

INSTRUCTIONS FOR USE

Evolution Esophageal Stent System - Partially Covered (EVO-20-25-8-E; EVO-20-25-10-E; EVO-20-25-12.5-E; EVO-20-25-15-E)

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

IMPORTANT INFORMATION

Please review prior to use.

INTENDED USE

This device is used to maintain patency of malignant esophageal strictures and/or to seal tracheoesophageal fistulas.

STENT DESCRIPTION

This flexible, self-expanding stent is constructed of a single, woven, nitinol wire. The stent foreshortens due to its design. The stent has a flange at both ends to aid in preventing migration. A silicone membrane applied to the stent covers the complete body of the stent along with part of both flanges in the partially

covered stent. The covering is intended to reduce the risk of tissue ingrowth and to provide a seal for tracheoesophageal fistulas. The total length of the stent is indicated by radiopaque markers on the inner catheter, indicating the actual length of the stent at nominal stent diameter. To aid in visibility under fluoroscopy there are four tantalum bands at either end of the stent. There is a lasso at the proximal end of the stent, the purpose of which is to reposition the stent in the event of incorrect placement. The proximal and distal ends of the stent are designed with numerous rounded "crowns" to potentially reduce tissue trauma.

DELIVERY SYSTEM DESCRIPTION

Stent is mounted on an inner catheter, which accepts a .035 inch wire guide, and is constrained by an outer catheter. A pistol-grip delivery handle allows controlled stent deployment or recapture. The position of the red marker on top of the handle in relation to the "point-of-no-return" position on the handle label indicates the point after which it is not possible to recapture the deployed stent.

CONTRAINDICATIONS

Those specific to upper GI endoscopy and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: total esophageal obstruction, strictures that cannot be dilated a minimum size as outlined in the precautions section, placement requiring positioning of stent within 2 cm of the cricopharynx, surgical resection candidates, hiatal hernia and gastric prolapse in the esophageal area, patients with a perforated esophagus, placement in actively bleeding tumors, benign diseases.

Relative contraindications include, but are not limited to: uncooperative patient, coagulopathy, tracheal compression, recent myocardial infarction, cervical arthritis with fixed cervical spine, large tumor mass occupying the mediastinum, nonobstructive tumor, gastric outlet obstruction, necrotic esophageal mucosa, acutely angled stenosis.

POTENTIAL ADVERSE EVENTS

Those associated with upper GI endoscopy include, but are not limited to: allergic reaction to contrast or medication, aspiration, cardiac arrhythmia or arrest, fever, hemorrhage, hypotension, infection, perforation, reflux, respiratory depression or arrest, vomiting.

Additional adverse events include, but are not limited to: airway compression; allergic reaction to nickel; aortoesophageal fistula and arterioesophageal fistula; chest or retrosternal pain; death (other than due to normal disease progression); dysphagia; edema; erosion or perforation of stent into adjacent vascular structures; esophageal ulceration and erosion; esophagitis; fistula involving trachea, bronchi or pleural space; food bolus impaction; foreign body sensation or reaction; gas bloat; inadequate stent expansion; intestinal obstruction secondary to migration; mediastinitis or peritonitis; nausea; pain/discomfort; reocclusion; sensitivity to metal components; sepsis; stent misplacement and/or migration; tracheal obstruction; tumor ingrowth or overgrowth; wire entrapment.

PRECAUTIONS

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

Stent should be placed using fluoroscopic monitoring.

Stent should only be placed with the Cook delivery system, which is provided with each stent.

Note: Prior to advancing system, area to be stented should be dilated to:

For (18mm x 23mm) stent – a minimum of 9 mm and a maximum of 11 mm. If area is dilated greater than 11 mm, stent may migrate.

For (20mm x 25mm) stent – a minimum of 10 mm and a maximum of 14 mm. If area is dilated greater than 14 mm, stent may migrate.

This device is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

Long-term patency of this device has not been established. Periodic evaluation is advised.

This device shortens upon deployment. With proximal strictures near the upper esophageal sphincter, deployment should be performed under fluoroscopic visualization as this may enhance placement accuracy.

Do not use this device for any purpose other than stated intended use.

Stent cannot be retrieved after the deployment threshold has been passed. The position of the red marker on top of the handle in relation to the point-of-no-return position on the handle label indicates when the threshold has been passed.

Use of this device is restricted to a trained healthcare professional.

WARNINGS

The **stent is not intended to be removed** and is considered a permanent implant. Attempts to remove stent after placement may cause damage to esophageal mucosa.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect product with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working conditions, do not use. Please notify Cook for return authorization.

This stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

Do not introduce the endoscope into the stent as displacement may occur.

After stent placement, alternative methods of treatment such as chemotherapy and radiation should not be administered as this may increase risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

The risk of perforation and erosion into adjacent vascular structures or aorto-esophageal and arterio-esophageal fistulas may be increased with pre- or post-operative chemotherapy and radiation, longer implantation times, aberrant anatomy, and/or mediastinal contamination or inflammation.

MR INFORMATION

This symbol means the device is MR Conditional

Non-clinical testing has demonstrated that Evolution Esophageal Stents are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

Static Magnetic Field

Static magnetic field of 3.0 Tesla or less

Maximum spatial field gradient of 1900 gauss/cm (19 T/m)

Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

MRI-Related Heating

1.5 and 3.0 Tesla Systems

Non-clinical testing was conducted under the following conditions and produced a maximum temperature rise of 4.4°C:

- A maximum whole body average specific absorption rate (SAR) of 3.8 W/kg (corresponding to a calorimetry measured value of 3.1 W/kg) for 15 minutes of MR scanning in a 1.5 Tesla Magnetom (Siemens Medical Solutions, Malvern, PA, Software Numaris/4) MR scanner.
- A maximum whole body average specific absorption rate (SAR) of 3.0 W/kg (corresponding to a calorimetry measured value of 2.8 W/kg) for 15 minutes of MR scanning in a 3.0 Tesla Excite (GE Electric Healthcare, Milwaukee, WI, Software G3.0-052B) MR scanner.

Image Artifacts

In nonclinical testing, the image artifact caused by the device extends approximately 7 mm from the Evolution Esophageal Stent when imaged with a spin echo pulse sequence and a 3 T MR system.

The safety of MR imaging using a local transmit/receive RF coil in areas immediately adjacent or directly involving the implant has not been established and is therefore not recommended.

The safety of performing MR imaging in patients with overlapping stents or other device(s) in the esophagus has not been established and is therefore not recommended.

EQUIPMENT REQUIRED

Stent and delivery system

Esophageal Balloon Dilators or Bougie Dilators

.035 inch wire guide

PATIENT PREPARATION

1. Intubate the patient using an endoscope per standard procedure. Assess the stricture location and dimensions upon direct visualization.
2. If required, dilate stricture to a minimum size as outlined in the precautions section using Balloon Dilators or Bougie Dilators.
3. Endoscopically locate and mark upper and lower margins of lesion with radiopaque marker. **Caution:** Accurate marking of lesion borders is essential for proper stent placement.
4. The esophageal stent should be at least 4 cm longer than the stricture, allowing for approximately 2 cm longer at either end of the stricture.
5. Place a .035-inch wire guide, floppy tip first, through accessory channel, through stricture, until it is fluoroscopically visualized in stomach. Leave wire guide in place and remove endoscope.

INSTRUCTIONS FOR USE

Stent Placement

1. Remove stent delivery system from the package and protective tubing.
2. Using fluoroscopic guidance, introduce and advance the delivery system in short increments over the pre-positioned wire guide into the esophagus.
3. To ensure that stent will bridge stricture after deployment, fluoroscopically position radiopaque markers on inner catheter beyond extremities of stricture to be crossed.
4. Confirm desired stent position fluoroscopically and deploy stent by removing red safety guard from the handle. (Fig. 1)
5. Continue deploying stent by squeezing trigger. **Note:** Each trigger squeeze will deploy stent by an equal amount.
6. If stent repositioning is required during deployment, it is possible to recapture stent. **Note:** It is not possible to recapture stent after passing point-of-no-return, indicated when red marker on top of the handle has passed the point-of-no-return position on the handle label. (Fig. 2)
7. To recapture stent, push direction button on side of delivery system to opposite side. (Fig. 3) **Note:** Hold thumb on button when squeezing trigger for first time to recapture.
8. Continue squeezing trigger as required to recapture stent by desired amount.

9. To resume deployment, push button to opposite side again and hold thumb on button for first stroke while squeezing trigger. (Fig. 4)
10. When stent point-of-no-return has been passed, pull safety wire out of delivery handle near wire guide port. (Fig. 5)
11. Continue deploying stent by squeezing trigger.
12. After deployment, fluoroscopically confirm full stent expansion. Once full expansion is confirmed, delivery system can be safely removed. Note: You may introduce endoscope and advance to the upper margin of the stent to endoscopically confirm position and patency of the stent.
13. To reposition stent directly after placement, grip the white lasso at proximal-end of stent with forceps and reposition stent to desired area.



POST PROCEDURE






- Observe the patient for development of any complications of endoscopy, esophageal dilation or stent placement. During the first 24 hours, vital signs should be monitored, and clear liquid given in an upright position. Patients being treated for fistula should have no fluids or food by mouth until successful seal of fistula has been confirmed.
- Subsequent follow-ups (such as endoscopy and/or chest x-rays) with patient should be performed at the physician's direction.

PATIENT COUNSELLING INFORMATION

- After stent placement, post 24 hours, the patient should be instructed to chew food well or eat soft or pureed food and to drink fluids before and after meals.
- Patients with stents placed in the distal esophagus or across the GE junction may be instructed to elevate the head of the bed and may be prescribed acid suppression therapy to minimize gastric reflux into the stent.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

INSTRUCTION FOR USE

Evolution Esophageal Stent System - Fully Covered (EVO-FC-R-18-23-8-E; EVO-FC-R-18-23-10-E; EVO-FC-R-18-23-12-E; EVO-FC-R-20-25-8-E; EVO-FC-R-20-25-10-E; EVO-FC-R-20-25-12-E)

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

IMPORTANT INFORMATION

Please review prior to use.

INTENDED USE

This device is used to maintain luminal patency of the esophagus in cases of: obstruction caused by intrinsic or extrinsic malignancies, refractory benign strictures, or to seal tracheoesophageal fistulas.

MALIGNANT AND BENIGN

The above symbol means the following: Malignant and Benign.

STENT DESCRIPTION

This flexible, self-expanding stent is constructed of a single, woven, nitinol wire. The stent foreshortens due to its design. The stent has a flange at both ends to aid in preventing migration. The stent has a silicone membrane covering. The covering is intended to reduce the risk of tissue ingrowth and to provide a seal for tracheoesophageal fistulas. The total length of the stent is indicated by radiopaque markers on the inner catheter, indicating the actual length of the stent at nominal stent diameter. To aid in visibility under fluoroscopy there are four tantalum bands at either end of the stent. There is a lasso at the proximal and distal ends of the stent that allows repositioning of the stent directly after placement or is used for stent removal in benign indications. There are radiopaque markers on the proximal and distal ends of the stent. The proximal and distal ends of the stent are designed with numerous rounded "crowns" to potentially reduce tissue trauma.

DELIVERY SYSTEM DESCRIPTION

The stent is mounted on an inner catheter, which accepts a .035-inch wire guide, and is constrained by an outer catheter. A pistol-grip delivery handle allows controlled stent deployment or recapture. The position of the red marker on top of the handle in relation to the "point-of-no-return" position on the handle label indicates the point after which it is not possible to recapture the deployed stent.

CONTRAINDICATIONS

Contraindications include those specific to upper gastrointestinal (GI) endoscopy and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to:

- total esophageal obstruction
- strictures that cannot be dilated to a minimum size as outlined in the precautions section
- placement requiring positioning of stent within 2 cm of the cricopharynx
- surgical resection patients
- hiatal hernia and gastric prolapse in the esophageal area
- patients with a perforated esophagus
- placement in actively bleeding tumors
- placement in polypoid lesions

Relative contraindications include, but are not limited to:

- uncooperative patient
- coagulopathy
- tracheal compression
- recent myocardial infarction
- cervical arthritis with fixed cervical spine
- large tumor mass occupying the mediastinum
- nonobstructive tumor

- gastric outlet obstruction
- necrotic esophageal mucosa
- acutely angled stenosis
- placement in an esophago-jejunostomy (following gastrectomy)

POTENTIAL ADVERSE EVENTS

Those associated with upper GI endoscopy include, but are not limited to: allergic reaction to contrast or medication, aspiration, cardiac arrhythmia or arrest, fever, hemorrhage, hypotension, infection, perforation, reflux, respiratory depression or arrest, vomiting.

Additional adverse events include, but are not limited to: airway compression; allergic reaction to nickel; aortoesophageal fistula and arterioesophageal fistula; chest or retrosternal pain; death (other than due to normal disease progression); dysphagia; edema; erosion or perforation of stent into adjacent vascular structures; esophageal ulceration and erosion; esophagitis; fistula involving trachea; bronchi or pleural space; food bolus impaction; foreign body sensation or reaction; gas bloat; inadequate stent expansion; intestinal obstruction secondary to migration; mediastinitis or peritonitis; nausea; pain/discomfort; reocclusion; sensitivity to metal components; sepsis; stent misplacement and/or migration; tracheal obstruction; tumor overgrowth; wire entrapment.

PRECAUTIONS

- A complete diagnostic evaluation must be performed prior to use to determine proper stent size.
- The Evolution Esophageal Stent – Fully Covered should be used with caution and only after careful consideration in patients with:
 - Significant preexisting pulmonary or cardiac disease.
 - Conditions where esophageal stenting may cause a secondary tracheal stenosis.
 - Prior gastric pull up operations resulting in altered anatomy.
- Stent should be placed using fluoroscopic monitoring.
- Prior to advancing the system, physician discretion should be used to determine if dilation is required. If required, dilate area to be stented to a diameter that is at least equal to the diameter of the delivery system and a maximum of 11 mm.
- Long-term patency of this device has not been established. Periodic evaluation is advised.
- This device shortens upon deployment. With proximal strictures near the upper esophageal sphincter, deployment should be performed under fluoroscopic visualization as this may enhance placement accuracy.
- For Malignant use, alternate methods of therapy should be investigated prior to placement.
- The stent should only be placed with the Cook delivery system, which is provided with each stent.
- Do not use this device for any purpose other than the stated intended use.
- Use of this device is restricted to a trained healthcare professional.
- Stent cannot be retrieved after the deployment threshold has been passed. The position of the red marker on the top of the handle in relation to point-of-no-return position on the handle label indicates when threshold has been passed.

WARNINGS

- The safety and effectiveness of leaving this stent in place, or removing the Evolution Esophageal Stent – Fully Covered from a benign lesion beyond 8 weeks has not been established.
- Exercise care when removing a stent that has significant overgrowth around the stent ends.
- It is not intended to remove a stent from a malignant lesion once original placement has been completed and confirmed. Attempts to remove the stent after placement in a malignant lesion may cause damage to the esophageal mucosa.
- Stenting across the gastroesophageal junction may increase the risk of migration.

- This stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- The effectiveness of the delivery system cannot be guaranteed if the stent has been recaptured more than 5 times.
- Testing of overlapping stents has not been completed and is not recommended.
- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.
- Do not introduce the endoscope into the stent, as displacement may occur.
- After stent placement, alternative methods of treatment such as chemotherapy and radiation should not be administered as this may increase risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.
- The migration rate may be increased by absence of a tumor abutment in benign esophageal stenoses, esophageal fistula without stenoses, incorrect sizing of the stent, or in cases of malignant stenoses where radiotherapy and/or chemotherapy is employed.
- The risk of perforation and erosion into adjacent vascular structures or aorto-esophageal and arterio-esophageal fistulas may be increased with pre- or post-operative chemotherapy and radiation, longer implantation times, aberrant anatomy, and/or mediastinal contamination or inflammation.

MR INFORMATION

This symbol means the device is MR Conditional

Non-clinical testing has demonstrated the Evolution Esophageal Stents are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

Static Magnetic Field:

- Static magnetic field of 3.0 Tesla or less
- Maximum spatial field gradient of 1900 gauss/cm (19 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/Kg (Normal Operating Mode).

MRI-Related Heating

1.5 and 3.0 Tesla Systems:

Non-clinical testing was conducted under the following conditions and produced a maximum temperature rise of 4.4 °C:

- a maximum whole body average specific absorption rate (SAR) of 3.8 W/kg (corresponding to a calorimetry measured value of 3.1 W/kg) for 15 minutes of MR scanning in a 1.5 Tesla Magnetom (Siemens Medical Solutions, Malvern, PA, Software Numaris/4) MR scanner.
- a maximum whole body average specific absorption rate (SAR) of 3.0 W/kg (corresponding to a calorimetry measured value of 2.8 W/kg) for 15 minutes of MR scanning in a 3.0 Tesla Excite (GE Electric Healthcare, Milwaukee, WI, Software G3.0-052B) MR scanner.

Image Artifacts

In non-clinical testing, the image artifact caused by the device extends approximately 7mm from the Evolution Esophageal Stent when imaged with a spin echo pulse sequence and 3T MR system.

The safety of MR imaging using a local transmit/receive RF coil in areas immediately adjacent or directly involving the implant has not been established and is therefore not recommended.

The safety of performing MR imaging in patients with overlapping stents or other device(s) in the esophagus has not been established and is therefore not recommended.

EQUIPMENT REQUIRED

- Stent and delivery system
- Esophageal Balloon Dilators or Bougie Dilators
- .035 inch wire guide

PATIENT PREPARATION

1. Intubate the patient using an endoscope per standard procedure. Assess the stricture location and dimensions upon direct visualization.
2. **If required**, dilate stricture to a minimum size as outlined in the Precautions section using Balloon Dilators or Bougie Dilators.
3. Endoscopically locate and mark upper and lower margins of lesion with radiopaque marker, internal or external, or use an injectable radiopaque substance. **Caution:** Accurate marking of lesion borders is essential for proper stent selection and placement.
4. The esophageal stent should be at least 4 cm longer than the stricture, allowing for approximately 2 cm longer at either end of the stricture.
5. Place a .035 inch wire guide, floppy tip first, through accessory channel, through stricture, until it is fluoroscopically visualized in stomach. Leave wire guide in place and remove endoscope.

INSTRUCTIONS FOR USE

Stent Placement

1. Remove stent delivery system from the package and protective tubing.
 2. Using fluoroscopic guidance, introduce and advance the delivery system in short increments over the pre-positioned wire guide into the esophagus.
 3. To ensure that the stent will bridge stricture after deployment, fluoroscopically position radiopaque markers on inner catheter beyond extremities of stricture to be crossed.
 4. Confirm desired stent position fluoroscopically and to enable stent deployment remove the red safety guard from the handle. **(Fig.1)**
 5. Continue deploying stent by squeezing the trigger. **NOTE:** Each trigger squeeze will deploy the stent by an equal amount.
 6. If stent repositioning is required during deployment, it is possible to recapture the stent. **NOTE:** It is not possible to recapture the stent after passing point-of-no-return, indicated when the red marker on top of the handle has passed the point-of-no-return position on the handle label. **(Fig.2)**
 7. To recapture the stent, push the direction button on the side of the delivery system to the opposite side. **(Fig.3)**
- NOTE:** Hold thumb on button when squeezing the trigger for the first time to recapture.
8. Continue squeezing the trigger as required, to recapture the stent by desired amount.
 9. To resume deployment, push the button to the opposite side again and hold thumb on the button for first stroke while squeezing the trigger. **(Fig.4)**
 10. When stent point-of-no-return has been passed, pull and remove the safety wire completely out of the delivery handle near the wire guide port to allow full stent deployment. **(Fig.5)**
 11. Continue deploying the stent by squeezing the trigger.
 12. After deployment, fluoroscopically confirm full stent expansion. Once full expansion is confirmed, the delivery system and wire guide can be safely removed.

NOTE: You may introduce endoscope and advance to the upper margin of the stent to endoscopically confirm position and patency of the stent.

Stent Repositioning Technique

To reposition the stent directly after placement, grip either the proximal or distal-end green lasso with forceps and reposition the stent to the desired area. (Fig.6)

Stent Removal Technique (Benign Indications Only)

1. Locate the proximal and distal-end green lasso of the stent by direct endoscopic visualization.
2. With forceps, gently apply tension to each green lasso to mobilize the distal and then the proximal end of the stent from the esophageal wall. (Fig.7)
3. Once the stent has been mobilized, grasp the green lasso at proximal end of stent securely and pull gently to collapse proximal end of stent and gently remove the stent. (Fig.8)


POST PROCEDURE







- Observe the patient for development of any complications of endoscopy, esophageal dilation or stent placement. During the first 24 hours, vital signs should be monitored and clear liquid given in an upright position. Patients being treated for fistula should have no fluids or foods by mouth until successful sealing of fistula has been confirmed.
- Subsequent follow-ups (such as endoscopy and/or chest x-rays) with patient should be performed at the physician's discretion.

PATIENT COUNSELLING INFORMATION

- After stent placement, post 24 hours, the patient should be instructed to chew food well or eat soft or pureed food and to drink fluids before and after meals.
- Patients with stents placed in the distal esophagus or across the GE junction may be instructed to elevate the head of the bed, and may be prescribed acid suppression therapy to minimize gastric reflux into the stent.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

MATERIALS

Evolution Esophageal Stent System-Partially Covered (EVO-20-25-8-E; EVO-20-25-10-E; EVO-20-25-12.5-E; EVO-20-25-15-E)

All materials mentioned in the table below are common to both the Evolution Esophageal Stent System-Partially Covered and Evolution Esophageal Stent System-Fully Covered unless indicated as follows:

Note A: Used in the Evolution Esophageal Stent System-Partially Covered.

Note C: Used in the Evolution Esophageal Stent System-Fully Covered EVO-FC-R-18-23-8-E; EVO-FC-R-18-23-10-E; EVO-FC-R-18-23-12-E; EVO-FC-R-20-25-8-E; EVO-FC-R-20-25-10-E; EVO-FC-R-20-25-12-E

Component	Material	Patient
Stent		
Stent wire	Nitinol wire Nitinol alloy (nickel-titanium) Nickel 54.5 -57.0 Carbon 0.05 Max Oxygen 0.05 Max Hydrogen 0.005 Max Titan Decisive	Yes
Silicone membrane	Nusil silicone grade MED-4755	Yes
Marker Bands	Tantalum NM-0010	Yes
Lasso loop	Force Fibre Suture Ultra high molecular weight Polyethylene Tecafine	Yes
Lasso loop ^{Note A}	Polytetrafluorethylene (PTFE) Impregnated Polyester 1113-T TevdekSuture	Yes
Grasping Feature ^{Note C}	Polytetrafluorethylene (PTFE) Impregnated Polyester 1113-T TevdekSuture	Yes
Introducer System		
Introducer system tip (inner catheter)	Pebax with 40% barium	Yes
Tubing (inner catheter)	Polyetheretherketone PEEK 381G	Yes
Marker bands	Platinum with 10% iridium	Yes
Pusher (inner catheter)	Polyethylene radiopaque tubing 30413	Yes
Cannula	Stainless steel 304	No
Support cannula (inner catheter)	Stainless steel 304	Indirect
Outer catheter ^{Note A}	Pebax, braided with stainless steel, fluorinated ethylene propylene lined. Bismuth subcarbonate in Pebax tip	Yes
Outer catheter ^{Note C}	Pebax, braided with stainless steel, PTFE lined. Bismuth subcarbonate in Pebax tip	Yes
Bilumen tube (inner catheter)	Vestamid	Yes

Retention wire	Nitinol Nitinol alloy (nickel-titanium) Nickel 54.5 -57.0 Carbon 0.05 Max Oxygen 0.05 Max Hydrogen 0.005 Max Titan Decisive	Yes
Retention wire cannula loop	Stainless steel 304	No
Introducer support cannula (inner catheter)	Stainless steel 304	No
Retention wire protection	Polytetrafluorethylene (PTFE) Impregnated Polyester 1113-T TevdekSuture	Yes
PERT Distal Collar ^{Note C}	Stainless Steel 304	Yes
Handle		
Various Components	Acetal / Polycarbonate grade Lexan / Stainless Steel 304 with PTFE coating / Polychloropren Neoprene / glasscord Kevlar-49/ polyamide 6/ Stainless Steel 304 / Aluminium / 10% glass-filled polycarbonate (compound 80020) / polycarbonate (compound 80130) / polycarbonate (compound 80130) / Acrylonitrile butadiene styrene blend (ABS POLIMAXX GA800)	No
Consumables		
Lubricante	Silicone liquid MED-400	Yes
Adhesive	Loctite 4061	Yes
Adhesive	Loctite 20268/26073/ 4013	Yes
Primer	Loctite 7701	Yes
Solvent	Toluene	No
Nusil Primer	MED6-161	Yes
Suture ^{Note A} ^{Note C}	Tevdek Polytetrafluorethylene (PTFE) 1113-T (TevdekSuture) Impregnated Polyester 6211	No

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Cook Ireland Ltd manufacture all medical products in accordance with Medical Devices Directive 93/42/EEC.

All of Cook facilities are also certified to ISO 13485 (2016) *Quality Management System*.

DISPUTE RESOLUTION PROCEDURE

Please notify Cook for return authorization.

Cook Medical Europe Ltd

O'HALLORAN ROAD

NATIONAL TECHNOLOGY PARK

LIMERICK, IRELAND

TEL: +353 61 334440 FAX: +353 61 334441
WWW.COOKMEDICAL.COM

WARRANTY INFORMATION

Authorized representative on the territory of the Russian Federation - Company "SCHAG" (ZAO "SCHAG"). Russia, 119002, Moscow, Karmanitzki per. 9, floor 5, office 501

Tel. (495) 956-13-09;

Fax (495) 956-1310,

e-mail: info@schag.ru

The shelf life of products is 2 (two) years. Manufacturer's warranty extends to the quality of products.

SHELF LIFE

The "use by" date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Evolution Esophageal Stent System - Partially Covered (EVO-20-25-8-E; EVO-20-25-10-E; EVO-20-25-12.5-E; EVO-20-25-15-E), Evolution Esophageal Stent System - Fully Covered (EVO-FC-R-18-23-8-E; EVO-FC-R-18-23-10-E; EVO-FC-R-18-23-12-E; EVO-FC-R-20-25-8-E; EVO-FC-R-20-25-10-E; EVO-FC-R-20-25-12-E) – have a two (2) year shelf life.

STERILIZATION DETAILS

These devices are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} . These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile.

CONTENTS OF DELIVERY

EVOLUTION coated esophageal stents with nitinol are delivered individually in a sterile package, which is hermetically sealed, complete with a delivery system and instructions for use.

RPN	Weight, g ($\pm 5\%$)	Product weight without packaging, g ($\pm 5\%$)
EVO-20-25-8-E	780	726.5
EVO-20-25-10-E	775.4	721.9
EVO-20-25-12.5-E	775.9	722.4
EVO-20-25-15-E	788.2	734.7
EVO-FC-R-18-23-8-E	1,182.9	778.0
EVO-FC-R-18-23-10-E	1,206.5	801.6
EVO-FC-R-18-23-12-E	1,181.1	776.2
EVO-FC-R-20-25-8-E	1,191.5	786.6
EVO-FC-R-20-25-10-E	1,190.3	785.4
EVO-FC-R-20-25-12-E	1,185.5	780.6

CONDITIONS OF USE

These products are used in medical institutions.

STORAGE RULES

These devices should be stored in a dry location, away from temperature extremes / these devices should be stored in a dark, dry, cool place, with avoidance of extended exposure to light.

TRANSPORTATION RULES


Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling required during transportation.

PACKAGING AND LABELLING

Supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.

The label of the device shows the following device parameters:

- the name of the medical device;
- trademark of the manufacturing company;
- catalog number;
- device diagram with dimensions;
- link to the glossary of symbols;
- method of sterilization;
- number of the series (lot);
- symbol "use before";
- symbol "date of manufacture";
- symbol "prohibition of reuse";
- the symbol "refer to the instructions for use";
- symbol "be careful!";
- symbol "to sale by or on the order of a physician";
- site of the manufacturing company;
- barcodes;
- the logo of the manufacturing company;
- name of the manufacturing company;
- address of the manufacturing company;
- symbol of the presence of the CE certificate.

REPAIRS

Devices are supplied sterile and for single-use only. Do not attempt to reuse or carry out maintenance.

DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste. For EU countries - in accordance with Council Directive 93/42/EEC/ For Russian Federation – in accordance with national requirements, regulated by sanitary rules and norms 2.1.3684-21 of January 28, 2021 'Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures'. Epidemiologically hazardous waste - class B.

Maryna Zhukava
Regulatory Affairs Specialist

29-July-2021
Date

Cook Ireland Ltd.

O'Halaran Road
National Technology Park
Limerick, Ireland

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