
20

INSTRUCTIONS FOR USE

Blood tubing lines for hemodialysis in the set (arterial tubing line and venous tubing line), models: AV18AFA-P, AV18AFA-HFP, AV18AFD-P, AV18AFA-SNP

Approved

Hiroki Moriyama

General Manager

Quality Assurance Department

Medical Division


NIKKISO CO., LTD.

NIKKISO

- (en) Instructions for use – BLOOD TUBING LINES FOR HEMODIALYSIS
- (de) Gebrauchsanweisung – BLUTSCHLAUCHSYSTEM FÜR DIE HÄMODIALYSE
- (fr) Instructions d'utilisation – TUBULURES À SANG POUR HÉMODIALYSE
- (it) Istruzioni per l'uso – LINEE EMATICHE PER EMODIALISI
- (es) Instrucciones de uso – LÍNEAS DE SANGRE PARA HEMODIÁLISIS
- (nl) Instructies voor gebruik – BLOEDLIJNEN VOOR HEMODIALYSE
- (pl) Instrukcja użytkowania – LINIE KRWI DO HEMODIALIZY
- (sv) Bruksanvisning – BLODSLÄNGAR FÖR HEMODIALYS
- (da) Brugsanvisning – BLODSLÄNGER TIL HÆMODIALYSE
- (no) Brukerveiledning – BLODSLÄNGER FOR HEMODIALYSE
- (fi) Käyttöohjeet – VERILETKUSTO HEMODIALYYSIA VARTEN
- (ru) Инструкция по эксплуатации – МАГИСТРАЛИ КРОВИ ДЛЯ ГЕМОДИАЛИЗА

- (el) Οδηγίες χρήσης – ΓΡΑΜΜΕΣ ΑΙΜΑΤΟΣ ΓΙΑ ΑΙΜΟΚΑΘΑΡΣΗ
- (tr) Kullanım talimatları – HEMODİYALİZ İÇİN KAN HORTUM HATLARI
- (cs) NÁVOD K POUŽITÍ – KREVŇÍ LINKY PRO HEMODIALÝZU
- (pt) INSTRUÇÕES DE UTILIZAÇÃO – LINHAS DE SANGUE PARA HEMODIÁLISE
- (hr) UPUTE ZA UPOTREBU – KRVNE LINIJE ZA HEMODIJALIZU
- (lt) NAUDOJIMO INSTRUKCIJA – HEMODIALIZĖS KRAUJO ŽARNELĖS
- (sk) NÁVOD NA POUŽITIE – HADIČKY VEDÚCE KRV PRE HEMODIALÝZU
- (ro) INSTRUȚIUNI DE UTILIZARE – LINII DE SÂNGE PENTRU HEMODIALIZĂ
- (sr) UPUTSTVO ZA UPOTREBU – KRVNE LINIJE ZA HEMODIJALIZU
- (uk) ІНСТРУКЦІЯ З ВИКОРИСТАННЯ – КРОВОПРОВІДНІ МАГІСТРАЛІ ДЛЯ ГЕМОДІАЛІЗУ
- (bg) ИНСТРУКЦИИ ЗА УПОТРЕБА – КРЪВНИ ЛИНИИ ЗА ХЕМОДИАЛИЗА
- (lv) NORĀDĪJUMI PAR LIETOŠANU – HEMODIĀLĪZEI PAREDZĒTAS ASINS CĪRKULĀCIJAS CAURULĪTES

NIKKISO Europe GmbH

Desbrocksriede 1, D-30855, Langenhagen, Germany

NIKKISO CO., LTD.

20-3, Ebisu 4-Chome, Shibuya-ku, Tokyo 150-6022, Japan

Revised Date: 2020-02-18

97-9565-004G - 1 -

LOT	Lot number / Chargennummer / Número de lot / Número di lotto / Número de lote / Batchnummer / Numer serii / Lotnummer / Varennummer / Varenpartinummer / Erännumero / Номер партии / Αριθμός παρτίδας / Lot numbers / Căso șarže / Número de lote / Broj serije / Partijos numeris / Číslo šarže / Număr lot / Broj serije / Номер партии / Партиден номер / Partijos numurs
	Date of manufacture / Herstellungsdatum / Date de fabrication / Data di produzione / Fecha de fabricación / Productiedatum / Data produkcji / Tilveeringsdatum / Fremstillingsdato / Produktionsdato / Valmistuspäivä / Дата изготовления / Ημερομηνία παραγωγής / Üretim tarihi / Datum výroby / Data de fabrico / Datum proizvodnje / Pagineiminen data / Datum výroby / Data fabricației / Datum proizvodnje / Дата изготовления / Дата на произвеждане / Izzatvošanas datums
	Expiry date / Verfalldatum / Date de péremption / Data di scadenza / Fecha de caducidad / Vervaldatum / Data ważności / Utgångsdatum / Utlöbsdato / Vimeinen käyttöpäivä / Срок годности / Ημερομηνία λήξης / Son kullanma tarihi / Datum expirace / Data de validade / Datum isteka valjanosti / Tinkamumo laikas / Datum expirácie / Data expirării / Rok upotrebe / Термін придатності / Срок на годност / Derīguma termiņš
STERILE	Sterile. Sterilized by steam / Steril. Mit Dampf sterilisiert / Stérile. Stérilisé à la vapeur / Sterile. Sterilizzato a vapore / Estéril. Estéril por vapor / Sterile. Gesteriliseerd met stoom / Sterilne. Sterilizacija para / Steril. Steriliserad med ånga / Steril. Steriliseret med damp / Steril. Steriliseret med damp / Sterill. Steriloitu höyryllä / Стерильно. Стерилизовано паром / Αποστειρωμένο / Αποστειρωμένο με ατμό / Steril. Buharla sterilize edilmistir / Sterilni. Sterilizováno parou / Esterilizado. Esterilizado por vapor / Sterilino. Sterilizirano parom / Sterilius. Sterilizuoti garais / Sterilné. Sterilizované parou / Produs steril. Produs sterilizat cu abur / Sterilino. Sterilisanu parom / Стерильный. Стерилизовано паром / Стерилно. Стерилизирано с пара / Steris. Sterilizēts ar tvaiku
	Do not re-use / Nicht wiederverwenden / Ne pas réutiliser / Non riutilizzare / No reutilizar / Niet hergebruiken / Do jednorazowego użycia / Får ej återanvändas / Må ikke genbruges / Må ikke brukes om igjen. / Ålä kåytä uudelleen / Повторное использование запрещено / Μη επαναχρησιμοποιείτε / Yeniden kullanmayın / Neponožujte opakovaně / Nāo reutilize / Nije za ponovnu upotrebu / Nensudoti pakaratolnā / Neponožujte opakovaně / A nu se reutilize / Ne upotrebljavati ponovo / Повторно не використовувати / Не използвайте повторно / Nelietot atkārtoti
	Refer to instructions for use / Gebrauchsanleitung beachten / Reportez-vous aux instructions d'utilisation / Rispettare le istruzioni per l'uso / Consulte las instrucciones de uso / Raadpleeg de gebruiksaanwijzing / Zapoznać się z instrukcją użycia / Se användningsanvisningarna / Lees brugsvejledningen for anvendelse / Se bruksanvisningen / Lue käyttöohjeet / См. инструкцию по эксплуатации / Ανατρέξτε στις οδηγίες χρήσης / Kullannin talimatlarina bakin / Viz návod k použití / Consulte as instruções de utilização / Pogledajte upute za upotrebu / Zr. naudojimo instrukcija / Pozrite si návod na použitie / A se vedea instrucțiunile de utilizare / Pogledajte uputstvo za upotrebu / Δια. инструкции за використання / Виекте инструкциите за употреба / Skatīt norādījumus par lietošanu
EC REP	Authorised Representative in the European Community / Zugelassene Vertretung in der EU / Représentant agréé dans la Communauté européenne / Rappresentante autorizzato nella Comunità Europea / Representante autorizado en la Comunidad Europea / Gemachtigde in de Europese Gemeenschap / Autoryzowany przedstawiciel w Unii Europejskiej / Auktorisatad EU-representant / Godkendt repræsentant i EU / Autoriseret repræsentant i EU-land / Valtuutettu edustaja Euroopan Yhteisössä / Авторизованный представитель в странах ЕС / Εξουσιοδοτημένος Αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Avrupa Birliğinde yetkili temsilcisi / Autorizovaný zástupce v Evropském společenství / Representante autorizado na Comunidade Europeia / Ovlašteni predstavnik u Evropskoj zajednici / Igaltais atstāvs Eiropas Bendrijā / Autorizovaný zástupca v Európe / Repræsentant autorizat în Comunitatea Europeană / Ovlašćeni zastupnik u Evropskoj zajednici / Уповноважений представник у Європейській спільноті / Упълномощен представител в Европейската общност / Pilnvarotais pārstāvis Eiropas Kopienā
	Manufacturer / Hersteller / Fabricant / Produttore / Fabricante / Fabrikant / Wytwórca / Tilveikare / Producent / Produsent / Valmistaja / Производител / Κατασκευαστής / Üretici / Výrobce / Fabricante / Proizvođač / Gamintojas / Výrobca / Producător / Proizvođač / Виробник / Производител / Ražotājs
	Recyclable package / Verpackungsrecycling / Emballage recyclable / Materiale Riciclabile / Envase reciclable / Recyclebare verpakking / Opakowanie nadaje się do recyklingu / Återvinningsbar förpackning / Genbrugsemballage / Resirkulerbar innpakning / Kierrätettävä pakkaus / Вторичная переработка упаковки / Ανακυκλώσιμη συσκευασία / Geri dönüştürülebilir paket / Recyklovateľný obal / Embalagem reciclável / Pakiranje namijenjeno recikliranju / Pardirbama pakuošė / Recyklovateľný obal / Ambalaj riciclabil / Ambalaža pogodna za reciklažu / Упаковка, що підлягає вторинній переробці / Рециркулируема опаковка / Pārstrādājams iepakojums
	Black universal recycling symbol / Schwarzes internationales Recyclingsymbol / Symbole universel noir de recyclage / Simbolo di riciclaggio universale nero / símbolo negro de reciclaje universal / Zwart universeel recycling symbool / Czarny uniwersalny symbol recyklingu / svart universell återvinningsymbol / Sort. universelt genbrugssymbol / Svart universellt resirkuleringsymbol / Musta kansainvälinen kierrätysymbol / Черный универсальный символ переработки / Προτόν καθόληλη για ανακύκλωση / Evrensel siyah geri dönüşüm sembolü / Černý univerzální symbol recyklace / Símbolo de reciclagem preto universal / Crni univerzalni simbol za recikliranje / Juoda bendrasis perdirbimo simbolis / Čierny univerzálny symbol recyklácie / Simbol universal de reciclare, de culoare neagră / Crni univerzalni simbol reciklaže / Універсальний символ вторинної переробки чорного кольору / Черен универсален символ за рециклиране / Melns universālais pārstrādes simbols
	Fragile. Handle with care / Vorsicht, zerbrechlich / Fragile, manipuler avec précaution / Fragile, maneggiare con cura / Frágil, manipular con precaución / Voorzichtig, breekbaar / Ostrożnie, produkt delikatny / Ömtålig, hanteras varsamt / Skrbelig, forsigtig / Ömtålig gods, må behandlas försiktigt / Särkyvä, käsiteltävä varoen / Хрупкое. Осторожно / Ευθραusto, χειρίζετο με προσοχή / Krikkabil, dikkatli taşym / Křehké, opatrně zacházejte / Frágil, manuseie cu cuidado / Lomljivo, rukujte pažljivo / Neatsparus, elgtis atsargiai / Křehké, manipulujte opatrně / Produs fragil. A se manipula cu grijă. / Lomljivo, pažljivo rukovati / Κρυσκός, ποσодитисεσ оберажно / Чупливо, работете внимателно / Trausis, rikoties uzmanigi
	Keep dry / Vor Nässe schützen / Protéger de l'humidité / Conservare in luogo asciutto / Mantener seco / Droog houden / Chronić przed wilgocią / Māste hālies torr / Opbevaras tørt / Oppbevaras tørt / Pidettävä kuivana / Бережть от влаги / Διατηρείτε το στεγνό / Kuru tutun / Chrániť před vlhkem / Manterha seco / Čuvajte na suhom / Laikyti sausal / Udržite suhé / A se păstra uscat / Čuvati na suvom mestu / Зберігати сухим / Поддержайте в сухо состоянии / Glatbāt sausā vietā
	Temperature limitation / Temperaturbegrenzung / Température de stockage / Limite de temperatura / Limite de temperatura / Temperaturlimit / Ograniczenie temperatury / Temperaturgränser / Temperaturbegrensning / Temperaturbegrensning / Lämpötilarajat / Ograniczenie temperatury / Περιορισμός θερμοκρασίας / Sicaklık limiti / Omezení teploty / Limite de temperatura / Ograničenje temperature / Temperatūros ribojimas / Obmedzenie teploty / Limită de temperature / Ograničenje temperature / Температурні обмеження / Ograničenje na temperaturata / Temperatūras ierobežojums
	This way up / Hier oben / Haut / Verso lalto / Esta parte hacia arriba / Deze kant boven / Tā šajā gāji / Denne side upp / Denne side op / Denne siden opp / Tämä puoli ylös / Вверх / Αυτή η πλευρά προς τα πάνω / Bu taraf yukarı / Touto stranou nahoru / Este a pusi uz augšu

	Stacking limit by mass / Stapelbegrenzung nach Menge / Poids de gerbage maximum (kg) / Limite di stoccaggio in base al peso / Limite de apilamiento por masa / Stapellimiet naar massa / Uklādiešanas pietvorne ogrācināšana pēc masas / Stäplingsgränse efter mängd / Stäplingsbegrensning etter vekt / Painorajoitus pinttaessa / Штабелирование ограничено / Определили количество штабелей по массе / Күтпеге göre istif limiti / Hmotnostní limit stohování / Limite de empihlamento por peso / Ograničenje za okomito slaganje prema masi / Krovimo vienas ant kito svorio riba / Obmedzenie stohovanie podľa hmotnosti / Limită de depozitare a produselor unul peste altul în funcție de greutate / Ograničenje slaganje prema masi / Обмеження на укладку за масою / Ограничение по тегло за подреждане едно върху друго / Sakraušanas ierobežojums pēc masas
	Do not use a sharp tool when opening the box / Kiste nicht mit einem scharfen Gegenstand öffnen / Ne pas utiliser d'outil tranchant pour ouvrir le carton / Non utilizzare uno strumento tagliente per l'apertura della scatola / No utilizar una herramienta afilada si abrir la caja / Geen scherp voorwerp gebruiken om de doos te openen / Nie używać ostrych narzędzi podczas otwierania pudełka / Öppna inte lådan med hjälp av ett vasst verktyg / Använd ikke et skarpt instrument ved åbning af kassen / Ikke bruk skarpe redskaper ved åpning av kassen / Ålä avas laatikkoa terävällä työkalulla / He использовать острый инструмент при открывании коробки / Μη χρησιμοποιείτε σκληρό εργαλείο για το άνοιγμα της συσκευασίας / Kutuyu açarken keskin alet kullanmayın / Ph otvırání bedny nepoužívat ostré nástroje / Nāo utilize uma ferramenta afiada para abrir a caixa / Prilikom otvaranja kutije ne koristite oštar alat / Dēžs neatdarinēti astri ierakni / Ne otvāranie krabīce nepoužijajte ostré nástroje / Nu utilizati un instrument ascuțit la deschiderea cutiei / Kutiju ne otvarati oštrom alatom / Ne використовувати гострий інструмент при відкриванні коробки / He използвайте остър инструмент при отваряне на кутията / Alverot kasti, nelietot asu instrumentu
	Non-pyrogenic / Nicht pyrogen / Non pyrogène / Non pirogenico / Apirogēno / Niet-pyrogēen / Niepirogēny / Ikke-pyrogent / Ikke-pyrogenisk / Pyrogeenitön / Niepirogēno / Μη πυρογόνος / Pirogenik değıdir / Niepyrogēni / Nāo pirogénico / Nije pirogēno / Niepirogēns / Nepyrogēny / Niepirogenic / Nije pirogēno / Niepirogēno / Niepirogēny / Niepirogēns
	Do not use if package is damaged / Nicht verwenden, wenn die Verpackung beschädigt ist / Ne pas utiliser si l'emballage est endommagé / Non utilizzare se la confezione è danneggiata / No usar si el paquete está dañado / Niet gebruiken als de verpakking beschadigd is / Nie używać, jeśli opakowanie jest uszkodzone / Använd inte om paketet är skadat / Må ikke anvendes, hvis emballagen er beskadiget / Må ikke brukes hvis emballasjen er skadet / Ålä kåytä, jos pakkaus on vahingoittunut / He использовать, если упаковка повреждена / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά / Ambalaj hasarlıysa kullanmayın / Neponožujte, je-li obal poškozen / Nāo utilizar se a embalagem estiver danificada / Ne upotrebljavati ako je pakiranje oštećeno / Nenaudoti, jei pakuotė pažeista / Neponožujte, ak je obal poškodený / A nu se utilizeza dacă ambalajul este deteriorat / Nemojte koristiti ako je pakovanje oštećeno / He використовувати, якщо упаковка пошкоджена / He използвайте, ако опаковката е повредена / He използвайте, ако опаковката е повредена

INSTRUCTIONS FOR USE


BLOOD TUBING LINES FOR HEMODIALYSIS

[GENERAL INFORMATION]

- Read all the instructions for use carefully before using this device.
- Read all instructions and safety messages in the DBB-EXA operator's manual before using this device.
- Refer to product or carton label for Type / Code of the blood tubing lines.
- The materials in direct or indirect contact with blood are: Polycarbonate(PC), Polypropylene(PP), Polyvinylchloride(PVC), Silicone rubber(SR).
- The blood volume inside the circuit during treatment is stated on the individual packages.

[INDICATIONS OF USE]

- This device should be used under the supervision and direction of a physician.
- This device consists of an arterial line and a venous line. It is designed for use with the DBB-EXA hemodialysis instrument. Do not use with any other hemodialysis instrument. The following figure shows a typical layout of this device.

[CONTRAINDICATIONS]

- This device is single use only. Re-use or reprocessing of a single use product may lead to contamination and compromised function or structural integrity of the product, possibly resulting in infection, injury, or death of the patient.

[WARNINGS AND CAUTIONS]

- Storage**
- Store in a clean, dry location between +5°C (41°F) and +30°C (86°F) temperature. Avoid direct sunlight.
- Do not store at extreme temperature or humidity.

Prior to Treatment

- The bloodlines are equipped with Luer-Lock connectors to provide safe connections to dialyzers, vascular accesses and various perfusion lines.
- This device does not contain any natural latex.
- This device may contain polycarbonate, which can be cracked by administration of fat emulsion preparations, drugs solubilized with surface activating agents ethanol, etc. Carefully check if this device is securely connected when using above mentioned medications during treatment. Do not fasten excessively and repeatedly so as to reduce the probability of crack.
- Locking connectors may separate if either the male or female part is exposed to a lubricant, for example, by transfer from a lubricated needleless valve.
- The air detector will not detect any air introduced by a syringe through an access port distal to the air detector.
- Prior to use, check device expiration date. Do not use if past expiration date.
- Inspect device packaging prior to use. Do not use if the package is opened or damaged.

- Do not use if the device is damaged.
- Avoid damage by hemostats or other instruments.
- This device is designed to provide a sterile and non-pyrogenic fluid pathway. To prevent contamination, aseptic technique is required throughout device set-up and prime/rinse.
- Ensure the hemodialysis instrument blood pump is calibrated for the appropriate device pump segment diameter according to the blood pump calibration instructions of the hemodialysis instrument manufacturer.
- If transducer protectors are not included with this device, use transducer protectors together with this device for treatment.

- To prevent contamination, use this device immediately after breaking the package seal.
- Monitor carefully to ensure all device connections are secure.

- To avoid potentially serious injury, ensure visibility of bloodline connection to fistula needle or catheter at all times during the treatment. Do not cover the access or bloodlines with any opaque material prior to or during the treatment to prevent blood leakage not being detected.

- Ensure that all not used female Luer-Lock connectors are secured by closed clamps to prevent leakage.
- Please tighten all Luer-Lock caps to prevent contamination.

During Treatment

- To decrease the risk of air embolization or leakage of blood, use this device at a blood flow rate of ≤ 500 ml/min. and with a pressure resistance level between +750 mmHg (100 kPa) and -500 mmHg (66,7 kPa).
- Do not continue treatment longer than the time recommended. Continuing treatment beyond that time could cause the actual blood flow to be slower than indicated by the machine, which might decrease the efficacy of treatment.
- Air entering the device may cause fatal air embolism to the patient. Use of an air/bubble detectors is required. Continually observe the drip chambers for adequate blood levels.

- Leaks at joints, at connections or at any point in the device may cause patient blood loss, air embolism or substitution fluid leaks. Monitor carefully for leaks and, if detected, take immediate corrective actions.

- Avoid device kinking or occlusion during the treatment. Excessive pressure may damage the extracorporeal circuit or patient blood access site, or may cause a disconnection and/or patient blood loss.

- Use a Luer-Lock type syringe for infusion lines. Keep infusion lines clamped when not in use.
- If transducer protectors become wetted with blood/saline or any other fluid.

- 000 ЦРММ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

- When taking blood samples from or injecting solutions through the access ports, do not use needles larger than 21 gauge. Insert needle at a 90° angle into the center of the sample/injection port and withdraw/insert to prevent leaks or patient injury.
- Disinfect the access port before use. Ethylalcohol and isopropyl alcohol are recommended for disinfection. If other disinfectants are used, check their compatibility with access port before use.
- Do not apply excessive force when connecting Luer lock connectors. It may cause deformation or breakage of connectors which leads to connectors stuck, improper tightening, air entry, and leakage of blood or substitution fluid.
- Check the accuracy of the machines as the actual blood flow rate may differ from the blood flow rate indicated by the hemodialysis machine and may change over time.

- Please be aware of using any drugs, chemicals and other medical devices according to their instruction for use.
- During the infusion of any drugs or chemicals be aware of air and blood clots.

After Treatment

- Discard device using universal precautions for potentially contaminated equipment.
- Comply with all applicable laws when discarding this device.

[INSTRUCTIONS FOR USE]

- Aspecially prepare and set up the bloodlines to the dialyzer and hemodialysis instrument.
- Ensure that all connectors are secure and there are no leaks in the blood tubing lines and substitution lines.
- Prime the extracorporeal circuit and dialyzer according to the dialyzer manufacturer's directions/instructions for use.
- If there is no pre-pump arterial drip chamber in this device, ensure that the pre-pump arterial pressure monitor line is filled with saline. When priming the blood tubing lines with normal saline or other priming solution, clamp the arterial patient line to trigger a maximum negative pressure alarm. This should cause the pre-pump arterial pressure monitor line to fill with approximately 30cm (12 inches) of priming solution.

- Adjust drip chamber levels as indicated in operator's manual.
- During treatment check fluid levels in drip chambers and adjust as required.
- Initiate treatment according to the dialyzer and hemodialysis instrument manufacturer's instructions for use.

- Discontinue patient treatment according to the dialyzer and hemodialysis instrument manufacturer's instructions for use.

[WARRANTY AND LIMITED LIABILITY]

- The manufacturer warrants that the device has been manufactured in accordance with its specifications and in compliance with ISO13485, applicable industry standards, and regulatory requirements.
- The manufacturer will not be liable for any misuse, improper handling, non-compliance with warnings, cautions and instructions, or for any damage incurred subsequent to the manufacturer's delivery of the device.
- If this device is deemed to contain manufacturing defects, it will be replaced if the manufacturer is notified of the device defect along with the device lot number.

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© 2000 ЦРМИ. INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU