POLYTECH

Instructions for Use

Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it

Manufacturer:

POLYTECH Health & Aesthetics GmbH

Alzheimer Strasse 32

64807 Dieburg

Germany

The device name:

Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it

in design versions:

I. The Mème designs

1. Mème Mesmo MS moderate profile (Article number: width (mm) x height (mm) x maximum projection (mm) x distance between the highest point of projection and lower edge of prosthesis (mm) x volume (ml)):

-15725-70: 70 x 70 x 28 x 47 x 70;

-15725-90: 75 x 75 x 30 x 50 x 90;

-15725-105: 80 x 80 x 31 x 53 x 105;

-15725-120: 85 x 85 x 32 x 56 x 120;

-15725-145: 90 x 90 x 33 x 59 x 145;

-15725-165: 95 x 95 x 35 x 62 x 165;

-15725-195: 100 x 100 x 36 x 65 x 195;

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	1/21

POLYTECH

-15725-220: 105 x 105 x 37 x 68 x 220;
-15725-255: 110 x 110 x 38 x 71 x 255;
-15725-270: 113 x 113 x 39 x 72 x 270;
-15725-285: 115 x 115 x 40 x 74 x 285;
-15725-300: 118 x 118 x 40 x 75 x 300;
-15725-320: 120 x 120 x 41 x 77 x 320;
-15725-360: 125 x 125 x 43 x 80 x 360;
-15725-400: 130 x 130 x 44 x 83 x 400;
-15725-445: 135 x 135 x 45 x 86 x 445;
-15725-495: 140 x 140 x 47 x 89 x 495;
-15725-545: 145 x 145 x 48 x 92 x 545;
-15725-600: 150 x 150 x 49 x 95 x 600;
-15725-660: 155 x 155 x 51 x 98 x 660.

2. Même Mesmo HS high profile (Article number: width (mm) x height (mm) x maximum projection (mm) x distance between the highest point of projection and lower edge of prosthesis (mm) x volume (ml)):

-15726-090: 70 x 70 x 34 x 51 x 90;
-15726-110: 75 x 75 x 36 x 54 x 110;
-15726-125: 80 x 80 x 37 x 58 x 125;
-15726-155: 85 x 85 x 39 x 61 x 155;
-15726-180: 90 x 90 x 40 x 64 x 180;
-15726-210: 95 x 95 x 42 x 67 x 210;
-15726-235: 100 x 100 x 44 x 71 x 235;
-15726-255: 103 x 103 x 45 x 72 x 255;
-15726-275: 105 x 105 x 45 x 74 x 275;
-15726-295: 108 x 108 x 46 x 75 x 295;
-15726-315: 110 x 110 x 47 x 77 x 315;
-15726-335: 113 x 113 x 48 x 79 x 335;
-15726-360: 115 x 115 x 49 x 80 x 360;
-15726-400: 120 x 120 x 51 x 83 x 400;

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary protheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	2/21

POLYTECH

-15726-450: 125 x 125 x 53 x 87 x 450;
-15726-505: 130 x 130 x 55 x 90 x 505;
-15726-560: 135 x 135 x 57 x 93 x 560;
-15726-620: 140 x 140 x 59 x 96 x 620;
-15726-690: 145 x 145 x 61 x 100 x 690;
-15726-760: 150 x 150 x 62 x 103 x 760;
-15726-835: 155 x 155 x 64 x 106 x 835.

3. Même Mesmo XS extra-high profile (Article number: width (mm) x height (mm) x maximum projection (mm) x distance between the highest point of projection and lower edge of prosthesis (mm) x volume (ml)):

-15727-105: 70 x 70 x 40 x 55 x 105;
-15727-125: 75 x 75 x 43 x 59 x 125;
-15727-155: 80 x 80 x 45 x 62 x 155;
-15727-180: 85 x 85 x 47 x 66 x 180;
-15727-215: 90 x 90 x 49 x 69 x 215;
-15727-250: 95 x 95 x 51 x 73 x 250;
-15727-285: 100 x 100 x 53 x 77 x 285;
-15727-305: 103 x 103 x 54 x 78 x 305;
-15727-330: 105 x 105 x 55 x 80 x 330;
-15727-355: 108 x 108 x 56 x 82 x 355;
-15727-380: 110 x 110 x 57 x 84 x 380;
-15727-430: 115 x 115 x 59 x 87 x 430;
-15727-485: 120 x 120 x 61 x 91 x 485;
-15727-545: 125 x 125 x 63 x 94 x 545;
-15727-610: 130 x 130 x 65 x 98 x 610;
-15727-675: 135 x 135 x 67 x 101 x 675;
-15727-755: 140 x 140 x 69 x 105 x 755;
-15727-830: 145 x 145 x 71 x 108 x 830;
-15727-920: 150 x 150 x 73 x 112 x 920;
-15727-1010: 155 x 155 x 75 x 115 x 1010.

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary protheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	3/21

POLYTECH

II. The REPLICON designs

1. RPL Mesmo MS moderate profile (Article number: width (mm) x height (mm) x maximum projection (mm) x distance between the highest point of projection and lower edge of prosthesis (mm) x volume (ml)):

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2. RPL Mesmo HS high profile (Article number: width (mm) x height (mm) x maximum projection (mm) x distance between the highest point of projection and lower edge of prosthesis (mm) x volume (ml)):

-15736-75: 70 x 70 x 34 x 42 x 75;

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary protheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	4/21

POLYTECH

-15736-95: 75 x 75 x 36 x 44 x 95;
-15736-115: 80 x 80 x 37 x 47 x 115;
-15736-130: 85 x 85 x 39 x 49 x 130;
-15736-160: 90 x 90 x 40 x 52 x 160;
-15736-180: 95 x 95 x 42 x 54 x 180;
-15736-210: 100 x 100 x 44 x 57 x 210;
-15736-240: 105 x 105 x 45 x 59 x 240;
-15736-260: 108 x 108 x 46 x 60 x 260;
-15736-275: 110 x 110 x 47 x 62 x 275;
-15736-300: 113 x 113 x 48 x 63 x 300;
-15736-315: 115 x 115 x 49 x 64 x 315;
-15736-330: 118 x 118 x 51 x 65 x 330;
-15736-350: 120 x 120 x 51 x 67 x 350;
-15736-395: 125 x 125 x 53 x 69 x 395;
-15736-445: 130 x 130 x 55 x 72 x 445;
-15736-495: 135 x 135 x 57 x 74 x 495;
-15736-550: 140 x 140 x 59 x 77 x 550;
-15736-605: 145 x 145 x 61 x 79 x 605;
-15736-670: 150 x 150 x 62 x 82 x 670;
-15736-735: 155 x 155 x 64 x 84 x 735.

3. RPL Mesmo XS extra-high profile (Article number: width (mm) x height (mm) x maximum projection (mm) x distance between the highest point of projection and lower edge of prosthesis (mm) x volume (ml)):

-15737-95: x 70 x 70 x 40 x 46 x 95;
-15737-110: 75 x 75 x 43 x 49 x 110;
-15737-130: 80 x 80 x 45 x 52 x 130;
-15737-160: 85 x 85 x 47 x 55 x 160;
-15737-185: 90 x 90 x 49 x 58 x 185;
-15737-220: 95 x 95 x 51 x 61 x 220;
-15737-255: 100 x 100 x 53 x 63 x 255;

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary protheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	5/21

POLYTECH

-15737-290: 105 x 105 x 55 x 66 x 290;
-15737-330: 110 x 110 x 57 x 69 x 330;
-15737-380: 115 x 115 x 59 x 72 x 380;
-15737-425: 120 x 120 x 61 x 75 x 425;
-15737-480: 125 x 125 x 63 x 78 x 480;
-15737-535: 130 x 130 x 65 x 80 x 535;
-15737-600: 135 x 135 x 67 x 83 x 600;
-15737-665: 140 x 140 x 69 x 86 x 665;
-15737-730: 145 x 145 x 71 x 89 x 730;
-15737-810: 150 x 150 x 73 x 92 x 810;
-15737-890: 155 x 155 x 75 x 95 x 890.

III. The OPTICON designs

I. Opticon Mesmo MS moderate profile (Article number: width (mm) x height (mm) x maximum projection (mm) x distance between the highest point of projection and lower edge of prosthesis (mm) x volume (ml)): Модели OPTICON

-15745-055: 70 x 59 x 28 x 35 x 55;
-15745-065: 75 x 63 x 30 x 37 x 65;
-15745-075: 80 x 67 x 31 x 39 x 75;
-15745-095: 85 x 71 x 32 x 41 x 95;
-15745-105: 90 x 76 x 33 x 43 x 105;
-15745-125: 95 x 80 x 35 x 45 x 125;
-15745-145: 100 x 84 x 36 x 47 x 145;
-15745-165: 105 x 88 x 37 x 49 x 165;
-15745-185: 110 x 92 x 38 x 51 x 185;
-15745-210: 115 x 97 x 40 x 53 x 210;
-15745-235: 120 x 101 x 41 x 55 x 235;
-15745-265: 125 x 105 x 43 x 57 x 265;
-15745-295: 130 x 109 x 44 x 59 x 295;
-15745-330: 135 x 113 x 45 x 61 x 330;
-15745-365: 140 x 118 x 47 x 63 x 365;

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	6/21

POLYTECH

- 15745-400: 145 x 122 x 48 x 65 x 400;
- 15745-445: 150 x 126 x 49 x 67 x 445;
- 15745-490: 155 x 130 x 51 x 69 x 490.

2. Opticon Mesmo HS high profile (Article number: width (mm) x height (mm) x maximum projection (mm) x distance between the highest point of projection and lower edge of prosthesis (mm) x volume (ml)):

- 15746-065: 70 x 59 x 34 x 40 x 65;
- 15746-075: 75 x 63 x 36 x 42 x 75;
- 15746-095: 80 x 67 x 37 x 44 x 95;
- 15746-115: 85 x 71 x 39 x 47 x 115;
- 15746-130: 90 x 76 x 40 x 49 x 130;
- 15746-155: 95 x 80 x 42 x 51 x 155;
- 15746-175: 100 x 84 x 44 x 54 x 175;
- 15746-205: 105 x 88 x 45 x 56 x 205;
- 15746-230: 110 x 92 x 47 x 58 x 230;
- 15746-265: 115 x 97 x 49 x 61 x 265;
- 15746-295: 120 x 101 x 51 x 63 x 295;
- 15746-335: 125 x 105 x 53 x 65 x 335;
- 15746-375: 130 x 109 x 55 x 68 x 375;
- 15746-420: 135 x 113 x 57 x 70 x 420;
- 15746-460: 140 x 118 x 59 x 72 x 460;
- 15746-510: 145 x 122 x 61 x 75 x 510;
- 15746-560: 150 x 126 x 62 x 77 x 560;
- 15746-615: 155 x 130 x 64 x 80 x 615.

IV. The OPTIMAM designs

1. Optimam Mesmo MS moderate profile (Article number: width (mm) x height (mm) x maximum projection (mm) x distance between the highest point of projection and lower edge of prosthesis (mm) x volume (ml)):

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary protheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	7/21

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2. Optimam Mesmo HS high profile (Article number: width (mm) x height (mm) x maximum projection (mm) x distance between the highest point of projection and lower edge of prosthesis (mm) x volume (ml)):

- 15776-095: 70 x 81 x 34 x 44 x 95;
- 15776-110: 75 x 87 x 36 x 46 x 110;
- 15776-130: 80 x 93 x 37 x 49 x 130;
- 15776-155: 85 x 99 x 39 x 52 x 155;
- 15776-180: 90 x 104 x 40 x 54 x 180;
- 15776-215: 95 x 110 x 42 x 57 x 215;
- 15776-250: 100 x 116 x 44 x 60 x 250;
- 15776-280: 105 x 122 x 45 x 62 x 280;

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	8/21

POLYTECH

- 15776-320: 110 x 128 x 47 x 65 x 320;
- 15776-365: 115 x 133 x 49 x 68 x 365;
- 15776-415: 120 x 139 x 51 x 70 x 415;
- 15776-460: 125 x 145 x 53 x 73 x 460;
- 15776-515: 130 x 151 x 55 x 76 x 515;
- 15776-570: 135 x 157 x 57 x 78 x 570;
- 15776-640: 140 x 162 x 59 x 81 x 640;
- 15776-705: 145 x 168 x 61 x 84 x 705;
- 15776-775: 150 x 174 x 62 x 86 x 775;
- 15776-855: 155 x 180 x 64 x 89 x 855.

Accessories:

- Additional polyethylene sterile delivery sleeve for mammary prostheses insertion

In the sets:

- 1) Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO outfitted with a sleeve:
 - Prosthesis in double blister sterile package
 - Instructions for use
 - Additional polyethylene sterile delivery sleeve for mammary prostheses insertion
 - Prosthesis passport (patient identification card)
 - Information label (8 pcs)
- 2) Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO:
 - Prosthesis in double blister sterile package
 - Instructions for use
 - Prosthesis passport (patient identification card)
 - Information label (8 pcs)

The peculiar feature of physical structure of the POLYTECH Health & Aesthetics prostheses is to closely resemble the main properties of human body tissue, thus allowing prostheses to be widely used in aesthetic and reconstructive surgery. Every type of prosthesis is available in different sizes and shapes.

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	9/21

POLYTECH

The present instructions for use applies to permanent mammary prostheses. Breast prostheses of the Sublime Line series are single-lumen prostheses filled with one type of silicone gel. Prostheses have a multi-layered silicone shell with microtextured surface MESMO. For more information on the shell and the surface, please refer to the chapter "Information on materials".

Design as well as production, delivery and customer service are subject to our quality management control according to EN ISO 13485. The products are certified and comply with the requirements of the European Medical Device Directive (MDD) and other national regulations.

All materials used are tested for biocompatibility and conform to the medical quality standards. All products are tested for the materials' mechanical properties such as elongation, tear resistance and fatigue have been investigated.

INFORMATION ON MATERIALS

Prosthesis shell

POLYTECH prostheses have a chemically and mechanically resistant, soft silicone elastomer shell which consists of several polydimethylsiloxane layers. The shell of gel-filled prostheses is additionally equipped with a special barrier layer. The prosthesis shell surface is microtextured MESMO. Prostheses with a textured silicone shell reduce, due to their surface structure, the length and orientation of fibres in the tissue capsule, thus diminishing the capsular contracture rate. Wall thickness of the shell depends on the size of the product: for microtextured MESMO prostheses it is roughly 0.35–1.0 mm. The average surface roughness of textured prostheses is approx. 150 µm.

Filler material

Prostheses are filled with a definite amount of clear, highly elastic silicone gel. Silicone gel, filler, consist of polydimethylsiloxane. This silicone gel shows a "memory effect", which means that after moderate compression it returns to its original form.

Accessories

For textured prostheses is offered one polyethylene, opened on both ends sleeve

INDICATIONS AND POSSIBLE CONTRA-INDICATIONS

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	10/21

POLYTECH

Indications

- contour augmentation (e.g. aplasia, hypoplasia) and contour correction (e.g. asymmetry, ptosis, atrophy)
- reconstruction (e.g. after mastectomy)
- prosthesis replacement
- combined breast and thorax deformities (e.g. Poland Syndrome, pectus excavatum, pectus carinatum)

Possible contra-indications

General clinical contra-indications in addition to general invasive measures can include:

- patients with mental instabilities
- repeated attempts or failures of contour correction
- persistent infections or systemic diseases
- allergies or extraordinary immune response to prostheses
- abscesses, cysts or tumours in the area of prostheses, especially cancer and recurrent metastases (e.g. persistent or intermittent mammary cancer), systemic connective tissue diseases
- severe radiation injuries in the prospective area of prostheses
- heavy burn scars in the prospective area of prostheses
- insufficient tissue covering in the prospective area of implantation (e.g. after preceding breast reduction) or reduced vascularization
- wound healing impairments
- costal injuries (concerning implantation in the area of mammae).
- pregnancy, breast feeding (concerning implantation in the area of mammae).

Custom-made prostheses

Customer-specific prostheses (normally mould-based) are subject to the same quality standards as applied to the above mentioned products. Any information given in this brochure also refers to custom-made prostheses.

TRANSPORT AND PACKAGING

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	11/21

POLYTECH

Products are delivered individually packed using sterile double packaging systems. Care should be taken that the inner package is neither opened nor damaged before the prosthesis to be definitely implanted.

Sealed double packages are placed into individual cardboard box.



The package of each product has personal labelling showing manufacturer, volume or size, order and serial number. Such labels should be affixed to the patient's records and to the prosthesis passport.

The following symbols are used on the labelling:

	Order number
	Serial number
	„Use-by date“ (month/year)
	Sterilisation with steam or dry heat
	Do not re-use
	Consult Instructions for use

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary protheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	12/21

POLYTECH

	Humidity limitation (0% - 70% rH)
	Upper limit of temperature: 35°C maximum for storage
	CE marking
	A: width in cm or mm B: height in cm or mm C: max. projection in cm or mm D: distance between the highest point of projection and the lower edge of the prosthesis. This distance shows the ratio between the position of the areola and the crease of the breast in mm

STERILIZATION, STORAGE, SHELF LIFE

WARNING: The prostheses are delivered in a sterile double-packaging. In the operating theatre, product sterility is ensured as long as the sterile inner package has neither been damaged nor opened. A replacement prosthesis should be available during implantation surgery. Product sterility is not guaranteed when the sterile packaging is damaged. We do not recommend to use such a prosthesis. When in doubt, please contact POLYTECH Health & Aesthetics or your local authorized distributor.

Beauty Matrix LLC, 119002, Moscow, Kaloshin pereulok, 4, building 1, 3rd Floor, office 3-K3; tel. +7 (499) 918-77-70; Info@bmatrix.ru

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	13/21

POLYTECH

Products are sterilized using dry heat. Sterility is guaranteed for 5 years from sterilization date (see label for maximum storage period ☒ YYYY-MM).

Product sterility is ensured as long as the sterilized inner package has neither been damaged nor opened. Sterility of products is not guaranteed with damaged packaging.

Products must be stored under dry conditions at room temperature so as not to exceed a relative humidity of 70% and a temperature of 35°C. Since the product packaging is gas-permeable, do not allow exposure of the product package to steam or gas. Prostheses should be stored in horizontal position so that labels can be easily read.

Products are made for single use only and must neither be reused nor resterilized. This is due to the following:

- The surface quality does not permit safe decontamination and no adequate decontamination procedure can be recommended.
- The products are low pressure / vacuum-sensitive and may not be subjected to the usual sterilization procedures.
- The product integrity especially of the shell is subject to a maximum of mechanical stress during implantation and explantation and cannot be sufficiently rated by the user.
- The mechanical safety limits as required by international standards for new prostheses cannot be guaranteed after implantation / explantation of the prostheses because of the mechanical forces applied during the procedure. Therefore, the required product safety cannot be guaranteed at reuse.

PATIENT INFORMATION

Every patient who has to undergo surgery, or any person authorized to act in his / her name on the occasion of surgery, must be informed by the attending physician about any operative or postoperative risks in general as well as those related to the prostheses. Patients should have realistic expectations of possible cosmetic result. Patients should be advised that prostheses are not intended for life and may have to be replaced after some time.

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	14/21

POLYTECH

Patients should be invited to consult their physician at regular intervals (every 6 to 12 months) to have the prostheses including surrounding tissues checked. Also, they should be informed that abnormal strain, excessive movement or traumatization of surrounding tissue can result in prosthesis rupture under the skin and, therefore, in a premature reoperation and explantation.

So that all information on the prosthesis in uses to be available for specialists in any time it is recommended to give to patients a passport of the corresponding prosthesis.

The physician in charge should advise patients on how to distinguish the prosthesis from the tissues after surgery and how to examine for tumours themselves.

Prosthesis passport

Prosthesis passports are attached to each POLYTECH Health & Aesthetics GmbH product. This passport should be completed and given to patients to ensure availability of implantation date and prosthesis type in case of need.

Screening diagnostics

Silicone prostheses can interfere with screening techniques and have a falsifying effect on results. It is, therefore, necessary that patients inform the medical staff involved about the existence of prostheses.

For a satisfactory mammography, we recommend to use the Eklund Technique as prostheses prove, to a certain extent, to be X-ray opaque. Ultrasonic or MRI examination may be an alternative to obtain more precise results.

POSSIBLE SIDE EFFECTS (COMPLICATIONS)

Possible complications associated with prostheses may be due to individual anamnesis, medication, operation technique, and postoperative care. Intra- and postoperative risks related to prostheses are specified below:

- Skin necroses, sloughings and extrusions may occur due to excessive tension of the skin overlaying the prosthesis (e.g. oversizing the prosthesis in relation to the prosthesis space or at quick tissue expansion by expanders) or by interrupted or reduced skin vascularization (e.g.

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	15/21

POLYTECH

traumatization of skin flap during the surgical intervention). Treatment with steroids may also provoke necroses. Smokers are high-risk patients as blood circulation in the overlaying skin may be impaired. Under such circumstance revisions may be necessary.

- Haematomas may occur prior to an infection, increasing the risk of capsular fibrosis, and therefore shall be avoided or treated immediately. Large postoperative oedemas or haematomas may, if untreated, provoke extrusion or displacement of prostheses. Non-resorbable haematomas may require an immediate explantation. Meticulous haemostasis during a surgery is the best prevention.
- Seromas, due to traumatization or excessive postoperative mobility, may cause swellings and pain. Possible therapies: immobilization, compression, drainage or explantation, if necessary.
- Swellings may be caused by serous fluid accumulation (refer to “seromas”).
- **Displacement or continuous pain** may be due to improperly sized or placed prostheses. Over-sized prostheses, capsular contracture as well as irritations due to excessive prosthesis movement may provoke such symptoms. Pain may occur in the operated area as well as in pectoralis and shoulder-arm areas after mammary surgery.
- **Folding** (also referred to as “wrinkling” in medical literature): palpable prosthesis edges or distortions may be due to a discrepancy between prosthesis size and cavity, insufficient tissue covering.
- **Infections** in connection with breast prostheses are very scarce (0.114%, Brand, 1993). Infections with unclear etiology after implantation should be treated immediately. The use of antibiotics, drainage or explantation can be necessary. A culture of the wound secretion should be made. Typical germs are: Staphylococcus aureus, Staphylococcus epidermidis or Streptococcus. Not all infections can be treated with the product remaining implanted. The “toxic shock” syndrome has been reported in single cases in connection with breast prostheses.
- **Reddening of the skin**, in the literature described as “rash”, can be observed with specific indications in different frequencies, depending on the surface structure (approximately in 1-5% of cases with Microthane prostheses) This reddening of the skin should not be confused with an infection. It differs from an infection by itching and the absence of systemic infection symptoms. It usually occurs 7 to 10 days after the implantation and can last 2 to 3 weeks. The use of steroids is recommended in literature for the case that a treatment is absolutely necessary.

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	16/21

POLYTECH

- **Capsular contracture** / fibrosis: formation of a connective tissue capsule around an prosthesis is a normal host response to the foreign body . Form and position of the prosthesis can change due to a capsular contracture. A clinically relevant (Baker III, IV) contracture may incur aesthetic impairment and pain. An explantation may become necessary. A closed capsulotomy is not recommended since this may damage the product. If such a capsulotomy cannot be avoided, a screening procedure should be employed to investigate possible prosthesis damage.
- **Desensibilization or hypersensibilization** may occur transiently or permanently (e.g. totally or partially in the breast area). There may be changes in sensibility in the presence of mammary prostheses, especially when a peri- or transareolar incision was made.
- **Prosthesis rupture** or prosthesis leakage are absolute indications for an explantation. Damaging of the prosthesis (intraoperatively, e.g. by too short incisions, by surgical instruments, or postoperatively, in case of punctures, biopsies, closed capsulotomies etc.), traumata or material fatigue are, in our opinion, the most frequent reasons. After rupture of a gel-filled prosthesis the consistency of the silicone gel prevents migration. However, it is not guaranteed that the gel remains a complete unity. After a rupture the prosthesis should be exchanged.
- **Permeation of silicone**: permeation of low molecular silicone components cannot be completely excluded. Normally, the minute amounts of gel remain within the capsule. Granulomas of unknown genesis are an indication for a biopsy or explantation.
- **Cancer risk**: epidemiologic studies of women having prostheses reveal that breast cancer risk for them is not higher if compared with control groups. By comparison: breast cancer risk is 1:9 in Western industrialized nations.
- **Possible association of autoimmune and connective tissue diseases** is discussed in connection with mammary prostheses. However, there is no proof of a causal coincidence.
- **The ability to breast feed** may, irrespective of breast prostheses, be affected by breast surgery. There is no evidence that breast prostheses could be detrimental to the **health of infants**.
- **Anaplastic large cell lymphoma (ALCL)**: FDA researchers report on a possible link between breast prostheses and the development of ALCL. FDA representatives are confident that women **with breast prostheses has a small increased risk of formation of such a disease. Specialists are recommended to inform their patients on the possibility of development this disease.**

SURGICAL TECHNIQUES

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	17/21

POLYTECH

Products are exclusively intended for implantation by skilled surgeons. Before implantation, each product should be checked for possible damage. In no case should prostheses be modified prior to or during implantation. Before implantation, each product should be checked for possible damage. Damaged or modified products must not be implanted!

A replacement prosthesis should be available during the implantation!

A direct combination of the prosthesis with a medicinal preparation has to be avoided. Implantation may be performed by a variety of surgical techniques. It is only up to the surgeon to decide on the technique which should be in line with the medical requirements as well as the individual needs of patients.

Patients should be examined prior to surgery in order to guarantee sufficient, tension-free tissue covering above the prosthesis and the incision. Pressure, traction or other stress factors, to which prosthesis and incision area may react, must be avoided.

Implantation must be performed carefully. Prostheses can be damaged by excessive manual strain, manipulation or surgical instruments. Should the elastomer shell be damaged, filler material might leak from the prosthesis.

We recommend to use a polyethylene sleeve for implantation of textured prostheses (refer to "Dolsky Technique") that may go with the prosthesis.

Following recommendations should be observed:

- Incisions should be made so as to correspond with the prosthesis size and surface.
- Haemostasis should be carried out meticulously. Cauterizers may damage the prosthesis shell.
- Prosthesis size should be selected carefully. It is recommended to store different sizes of prostheses in order to have the appropriate one available. We recommend using sizers to determine the appropriate size during surgery (see current price list).
- Prosthesis cavity should correspond to prosthesis size. Prostheses should be placed flatly and without wrinkles.

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	18/21

POLYTECH

- Appropriate tissue covering must be ensured or, if necessary, preference be given to submuscular placement of prostheses.

Prevention of contamination

The prosthesis surface must not be contaminated with talcum, dust or skin oils. Cotton fibres, talcum or other surface contaminants can be attracted by electrostatic charging and, thus, provoke an increased foreign body reaction after implantation.

Extraordinary care has to be taken to prevent contamination. Medical staff involved should, therefore, clean sterile gloves with sterile saline solution before touching the prosthesis.

Recommended passing procedure


1. Open outer package under aseptic conditions. Remove outer sterile packaging. Attach self-adhesive product label to patient's records.
2. Open outer sterile packaging. Offer inner sterile package to the nurse for putting it on the sterile instrument table under sterile conditions. The same handling procedure is followed for separately packed accessories.
3. Open inner sterile package under sterile conditions shortly before implantation. Prepare the sterile prosthesis according to the requested surgical technique and hand it to the surgeon.

Dolsky Technique

- Moisten prosthesis and delivered sleeve with sterile saline solution.
- Introduce the prosthesis into the sleeve (1).
- Fold the proximal end of the sleeve over the top of the prosthesis.
- Twist the distal end of the sleeve until the sleeve tightly encloses the prosthesis.
- Introduce the sleeve and the prosthesis by the proximal end first with the open sleeve folded over the top (2).
- Bring the prosthesis into its correct position and keep it there (3) while removing the sleeve (4).

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary prostheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	19/21

POLYTECH

POST-OPERATIVE CARE

The patient's condition should be controlled after surgery during an appropriate period of time in order to follow up as quickly as possible and to promptly provide qualified medical care (refer to Complications)

WARRANTY

Products are produced and controlled in conformity with the Medical Devices Act. Individual components have been selected thoroughly and comply with the latest standards for medical excellence, as does manufacturing.

Any information laid down in the "Instructions for use" must be followed. Replacement of products by POLYTECH Health & Aesthetics GmbH is only possible if manufacturing faults are definitely proven.

There is no warranty whatever, expressed or implied, including any warranty under commercial law and any warranty for modifications the user may apply to products.

Title	Status	Version	Date	Identification	Page
Instructions for Use Silicone gel filled sterile mammary protheses of Sublime line series with microtextured surface MESMO, outfitted with a sleeve and without it	Approved	3	2021-07-09	AB	20/21

POLYTECH

No representative or distributor of the manufacturer is entitled to change the above-listed warnings or precautions. Neither is any natural person nor legal party allowed to assume any liability or responsibility except the manufacturer.

Return of products

On all matters relating to the circulation of this medical product in Russia please apply to our Authorized representative in Russia: Beauty Matrix LLC; Russian Federation 119002, Moscow Kaloshin Pereulok, 4, building 1, 3rd Floor, office 3-K3; Tel. 7 499 918-77-70; info@bmatrix.ru

Disposal


Disposal of contaminated or used prostheses and prostheses of end shelf-life must be made according to standard practice, adopted in the medical institutions and procedure and guidelines established by local, federal and state regulations.

The regulations for the disposal of surgical waste also apply to used prostheses. Decontaminated and sterilized prostheses can be disposed of as domestic waste in a procedure established by local, federal and state regulations and rules of law.

Our business activities are based on the DIN EN ISO quality standards guide. Our production effort, manufacturing and quality management is governed by the Medical Devices Act and ISO 9001 / EN ISO 13485, ranging from the choice of raw material, product design, production, sterilization, storage, final distribution to customer service.

Herewith I declare, that the present Instruction for Use is applicable to mammary prostheses with textured surface (MESMO). It corresponds to the latest version of this document at time of this declaration.

Dieburg, 15.09.2021

Dr. Oliver Bögershausen
Senior Director Operations

Title
Instructions for Use
Silicone gel filled sterile mammary
prostheses of Sublime line series with
microtextured surface MESMO, outfitted
with a sleeve and without it

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I hereby certify the above signature
of Mr.

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