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
NGƯỜI DỊCH

Sam

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Protective coverall disposable non-sterile
(Russian edition)

HEN EDWARDS PEIKANG
Chairman
Environstar Limited Company
April 28, 2021

1. Name of the medical device

Protective coverall disposable non-sterile

Execution options:

1. Protective coverall with hood STANDART, sizes: S, M, L, XL, 2XL, 3XL, 4XL, 5XL.
2. Protective coverall with hood OPTIMAL, sizes: S, M, L, XL, 2XL, 3XL, 4XL, 5XL.
3. Protective coverall with PREMIUM hood, sizes: S, M, L, XL, 2XL, 3XL, 4XL, 5XL.

2. Scope of the medical device

2.1 Purpose

Coveralls are intended to be used as a barrier to protect against harmful biological factors in the conditions of medical institutions.

2.2. Potential customers

Employees of institutions related to public health services

3. Use of a medical device

3.1 Indications for use

Coveralls are used to protect medical personnel from pathogens of infectious diseases when providing medical and preventive care to patients in medical institutions.

3.2 Contraindications:

Individual intolerance of raw material components of the product. Do not use during surgical operations and in other medical premises with similar requirements. Do not sterilize before use.

3.3 Side effects:

No.

3.4 Appearance and application

Coveralls are an effective means of protecting medical personnel and reducing the transmission of potentially pathogenic microorganisms.

It is necessary to use the correct technique for removing and putting on coveralls.

If the integrity of the coveralls is compromised (for example, if they are punctured), they should be replaced as soon as possible.

The exterior of комбинезонов, the insulating coveralls must be dry, clean, free from dust and mildew, and the surface must not have scratches, cracks or other defects.

The connecting parts of overalls are processed with the help of seams. The seam is sewn with a needle to seal the tape, the distance between the seams every 3 cm should be from 8 to 14 stitches. The stitch should be smooth and straight. The treated part of the tape must be smooth and airtight, free of air bubbles.

The zipper cannot be opened, and the pull rod must be self-locking.

The structure of the coverall should be loose, comfortable to wear and remove, tight-fitting parts.

The coverall has elastic cuffs, ankles, front side of the hood and waist.

ATTENTION!

- Always follow the instructions for putting on coveralls to avoid contamination.
- Protect coveralls and packaging from mechanical damage.
- Do not use coveralls if the packaging is damaged.
- Do not use coveralls after the expiration date.
- Never wash or reuse disposable coveralls after use.

3.5 Putting on and taking off coveralls

3.5.1 How to put on a coverall:

1 Wear work clothes

2 Take the coverall as follows: with your right hand, hold the right sleeve by the cuff and the right leg by the bottom, with your left hand, hold the left sleeve by the cuff and the left leg by the bottom

3 Start your leg, gently releasing the leg so that it does not touch the floor, while not letting go of the sleeve cuffs

4 Take the zipper with your hand by the lower part and pull it down, spreading it out, your right hand to fasten the zipper to the end

3.5.2 How to remove the coverall:

1 Unzip the zipper. Remove the coverall, turning it inside out. Submerge in a decontamination container

2 Gather the edges of the cape at the back of your head in one hand and remove the hood.

3 Remove your work clothes

4. Description of the medical device

4.1 Color-white.

4.2 There are eight standard sizes of coveralls: S, M, L, XL, 2XL, 3XL, 4XL, 5XL.

4.3 Non-sterile products.

4.4 The products are intended for one-time use only.

4.10 Coveralls are classified as Class I according to the classification for potential use risk.

5. Technical Specifications

5.1 Main technical characteristics

5.1.1 Linear dimensions of coveralls depending on their size are shown in Table 1.

Table 1. Linear dimensions of coveralls

Size	Growth	Chest circumference	Sleeve Length	Cuffs
S	165	120	84	18
M	169	125	86	18
L	173	130	90	18
XL	178	135	93	18
2XL	181	140	96	18
3XL	188	145	99	18
4XL	195	150	102	18
5XL	202	155	105	18
Deviation	±2	±2	±2	±2

A schematic representation of the coveralls is shown in Figure 1.


Fig. 1: Schematic representation of the coverall

5.1.2 Strength characteristics.

The tensile strength of the material in key areas of coveralls should be at least 45 N.

The elongation at break of the material of the key parts of coveralls should be at least 15%.

5.1.3 Performance indicators

Coveralls must withstand a hydrostatic pressure of 1.67 kPa. Coveralls should not leak.

Filtration efficiency of non-oily material particles in key areas of coveralls should be at least 70%.

The moisture permeability of the material of coveralls should be at least 2500 g / m².

The charge of protective clothing should not exceed 0.6 μS / pc.

5.2 Materials in contact with human skin

5.2.1 STANDART protective coverall with hood is made of the following materials:

- nonwoven fabric (SPP, white, 35 g / m²).

5.2.2 OPTIMAL protective coverall with hood is made of the following materials:

- nonwoven fabric (SMS, white, 55 g / m²).

5.2.3 PREMIUM protective coverall with hood is made of the following materials:

- nonwoven fabric (Microporous, white, 60 g / m²).

6. Storage and transportation conditions

6.1 Products can be transported by all modes of transport in accordance with the transport regulations applicable to a particular mode of transport in the original factory transport packaging (cartons). The mode of transport used must be protected from direct sunlight.

The coveralls are expected to maintain their performance characteristics when used in the following environmental conditions: temperature and relative humidity:

temperature - from +10 to +35 °C;

relative humidity — up to 80% at a temperature of 25 °C.

Coveralls in transport packaging must be resistant to the following climatic influences:

Temperature - from -30 to +50 °C;

relative humidity — up to 98% at a temperature of 35 °C without condensation.

6.2. Coveralls should be stored in a well-ventilated area away from heat sources, direct sunlight, light sources and electrical equipment. Storage temperatures below 10 °C may cause condensation. If the storage temperature drops below 10 °C, warm up the coveralls to 25-30 °C before use, otherwise they may be difficult to put on due to loss of elasticity.

6.3 Coveralls should be stored in a dry place, away from direct sunlight and high humidity at temperatures between 10 °C and 30 °C.

6.4 Coveralls should be kept away from sources of ignition.

7. Terms of use of the medical device

Coveralls are used in the following environmental conditions:

temperature - from +10 to +35 °C; relative humidity - up to 80% at +25 °C.

8. Packaging

8.1 Packaging should provide protection against mechanical and climatic factors during transportation and storage, as well as ensure full loading (capacity) of vehicles and ease of loading and unloading operations.

8.2 Consumer packaging

Coveralls are packed in plastic bags of 1 pc. Consumer packaging protects coveralls from external mechanical damage and allows you to visually determine the size of coveralls.

8.3 Transport package

Coveralls in transport packaging are packed in cardboard boxes. One transport package includes 30 consumer packages of coveralls. Cardboard boxes must be hard enough to ensure the safety of products during transportation and storage.

8.4 Package dimensions and quantity of products

Number of coveralls in consumer packaging	Packaging Parameters						Number of consumer packages in one transport package
	Consumer packaging ±2mm			Transport package ±2mm			
	L	W	In	HL	W	In	
H I	460	280	120	410	255	320 360 400	30

9. Marking

9.1. Symbols used in marking:

- "Batch number"	LOT	
- "Use once"	⊗	

- "Production date"		(day/year/Month)
- "Manufacturer"		
- "Use before ..."		(day/year/month)
- "Avoid exposure to sunlight"		
- "Refer to the instructions for use"		

9.2. Marking design

Protective coverall disposable non-sterile Version: Protective coverall with hood STANDART Size: L	
 LOT LF-202024-O	 Eaviron Star Limited Company
 15.04.2020	Address: Lot L 1.5 and 1/3 L1. 6 Do Son Industrial Zone, Ngoc Sien Ward, DO Son District, Hai Phong City, Vietnam
 14.04.2025	Authorized Representative: Avangard Safety Group LLC, Address: 109052, Moscow, Ryazansky Prospekt, 2, bld 49
  	Roszdraznadzor registration Certificate No. _____

10. Scope of delivery

No	Full name	Quantity
I	Protective coverall with hood	1 pc
II	Packing	1 pc.

11. Warranty and expiration date

The manufacturer guarantees that the product will have all the characteristics specified in this document for the entire shelf life, provided that the integrity of the packaging is preserved and the conditions of storage, transportation and use are observed in accordance with the instructions for use.
Shelf life – 5 years.

12. Procedure for disposing and destroying a medical device

The system for collecting, temporarily storing, and transporting medical waste should include the following steps:

- waste collection by organizations engaged in medical and / or pharmaceutical activities;
- waste removal and temporary storage on the territory of the waste producing organization;
- disinfection / decontamination of waste;
- waste removal from the territory of the waste producing organization;
- disposal of medical waste.
- The method of chemical disinfection of waste at the place of its formation is used as a mandatory temporary measure in the absence of a site intended for processing medical waste in medical

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