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
Siraj Syeda Taj
Manager Quality
Control & Legal
Affairs

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MUHAMMAD ASIF
B.A.L.L.B. Advocate
NOTARY PUBLIC
Karachi-Pakistan

Medical device name

SANIPLAST® First Aid Bandage in different design options:

1. SANIPLAST® water-resistant First Aid Bandage with Acrinol pad, medium-sized, 19 mm x 72 mm, 10; 20; 100 plasters
2. SANIPLAST® water-resistant First Aid Bandage with Acrinol spots pad, round, diameter 25 mm, 20; 100 plasters
3. SANIPLAST® water-resistant First Aid Bandage with Acrinol pad, large, 25 mm x 72 mm, 20; 100 plasters
4. SANIPLAST® water-resistant First Aid Bandage with Acrinol pad, square, 38 mm x 38 mm, 20 plasters
5. SANIPLAST® water-resistant First Aid Bandage with Acrinol pad, mixed package, 20 plasters: SANIPLAST® medium — 8 pcs; SANIPLAST® round — 4 pcs; SANIPLAST® large — 4 pcs; SANIPLAST® square — 4 pcs

Developed by

UNIFEROZ (PVT) LTD

32/8, Sector-15, Korangi Industrial Area Karachi, Pakistan

Manufacturer

UNIFEROZ (PVT) LTD

32/8, Sector-15, Korangi Industrial Area Karachi, Pakistan

Place of manufacture of a medical device

UNIFEROZ (PVT) LTD

32/8, Sector-15, Korangi Industrial Area Karachi, Pakistan

Authorized representative of the manufacturer in the Russian Federation:

REGICOM Limited Liability Company (REGICOM, LLC)

123290 Moscow, 39 Shmitovsky Driveway, building 2, apt. 10
Tel.+7 926 385 5611

Intended use/indications of the medical device

Intended use: Protection and healing of small wounds, cuts, blisters and insect bites.

Potential consumer: Personal care products

Indications: SANIPLAST® is used to protect small wounds when a full dressing is impractical.

Technical specifications of the medical device

Description and appearance

SANIPLAST® First Aid Bandage protects wounds from germs, dust and water, thus protecting the wound from infection. It consists of a porous, waterproof material that enables skin breathability. The antiseptic pad protects the wound and aids in healing response.

Table. Components and description of SANIPLAST® First Aid Bandage

Component	Configuration
SANIPLAST® First Aid Bandage	Siliconized paper reinforced self-adhesive, porous, light brown tape with a yellow acrinol pad in the center. The plaster is packed between flat surface paper sheets and typographic printing by cold sealing.

Key specifications

Plaster Specifications

No.	Name	Specifications
1	SANIPLAST® water-resistant First Aid Bandage with Acrinol pad, medium-sized, 19 mm x 72 mm, 10, 20, 100 plasters	The middle stripe is 19 x 72 mm +/- 2 mm, pad is 110 to 115 g/m ² , thickness without pad is 0.13 mm +/- 0.1 mm, with pad 0.55 mm +/- 0.1 mm. – pad size is 11 mm x 25 mm +/- 1 mm – color is yellow
2	SANIPLAST® water-resistant First Aid Bandage with Acrinol spots pad, round, diameter 25 mm, 20; 100 plasters	The spot stripe is 25 x 25 mm +/- 2 mm, pad is 110 to 115 g/m ² , thickness without pad is 0.13 mm +/- 0.1 mm, with pad 0.58 mm +/- 0.1 mm. – pad size is 11 mm x 11 mm +/- 1 mm – color is yellow
3	SANIPLAST® water-resistant First Aid Bandage with Acrinol pad, large, 25 mm x 72 mm, 20; 100 plasters	The large stripe is 25 x 72 mm +/- 2 mm, the pad is 110 to 115 g/m ² , thickness without pad is 0.13 mm +/- 0.1 mm, with pad 0.52 mm +/- 0.1 mm. – pad size is 14mm x 25mm +/- 1mm – color is yellow
4	SANIPLAST® water-resistant First Aid Bandage with Acrinol pad, square, 38 mm x 38 mm, 20 plasters	The square stripe is 38 x 38 mm +/- 2 mm, the pad is 110 to 115 g/m ² , thickness without pad is 0.13 mm +/- 0.1 mm, with pad 0.55 mm +/- 0.1 mm. – pad size is 20 mm x 20 mm +/- 1 mm – color is yellow
5	SANIPLAST® water-resistant First Aid Bandage with Acrinol pad, mixed package, 20 plasters: SANIPLAST® medium — 8 pcs; SANIPLAST® round — 4 pcs; SANIPLAST® large — 4 pcs; SANIPLAST® square — 4 pcs	This package contains all the above parameters

Pad Specifications

No.	Parameter	Value
1	Functional pad loading capacity	0.076 cm ³ /cm ² - 0.084 cm ³ /cm ²
2	Functional pad wetting time	4 to 5 sec.
3	Bacterial purity (total count of aerobic bacteria and yeasts and fungi (cumulatively))	200 CFU/cm ²
4	Antibacterial activity (bacteriostatic) measures	1 to 2 mm
5	Noninvasive functional pad (adhesive rate for wound surface model)	Functional pad without noninvasive layer

Package Specifications

Packaging name	Material	Manufacturer
Device surface protective layer	Siliconized paper White	NICHIBAN CO., LTD Japan
Primary package	Coated paper stock White	NICHIBAN CO., LTD Japan
Secondary package	Cardboard	

Performance:

Adhesion is strong and the patch is easily peelable.

Peel force is established in the range from 140 N/m to 600 N/m

Physical properties:

Acrinol is dissolvable in isopropyl alcohol and hot water.

Name	Material, used in Manufacturing Content range	MD Composition	Manufacturer
Polyvinyl chloride film	Polyvinyl chloride Terpolymer, 99.7% CAS: 9002-86-2 Light brown textured PVC film Terpolymer has a molecular weight of 60,000–15,000 g/mol. Auto ignites at 435 °C, but stable at normal temperature and pressure	70.2 %	Hanwa Chemical Corporation Korea
Acrinol pad	Ethacridine lactate, 0.45-0.55% CAS: 6402-23-9 Non-woven fabric, 100 % cotton with perforated padding	10.4 %	Young Chemical Company Limited Korea
Adhesive; Acrylate adhesive	Solvent Ethyl acetate, Solvent content from 50% to 54%, polymer acrylate content from 46 % to 50 % CAS: 141-78-6 2-ethylhexyl acrylate, CAS: 103-11-7 Clear viscous liquid	19.4%	Woori Chemicals South Korea

Description of application method

Clean and dry the wound and skin before use.

Replace SANIPLAST® daily for better wound care.

In case of an allergic skin reaction or deep cut, consult your physician.

Store at room temperature, protect from sunlight.

Remarks

Prolonged skin contact, thus the adhesive used may, in rare cases, cause skin allergies.

Plaster removal may cause a painful sensation if the adhesive properties missed the proper condition.

Elderly and pediatric patients may experience skin stretching.

Storage and transportation conditions for the medical device

Store in a cool, dry place. Optimum conditions include temperature not exceeding + 25 °C.
Keep away from direct sunlight.

The product can be shipped by any means of transport in covered vehicles according to the applicable transportation rules.

The product remains safe and stable at a temperature of + 25 °C to + 40 °C degrees, with a humidity of 60 % to 75 %,.

Compliance with Russian national standards

GOST P 53498-2019 — Russian national standard. Dressings of adhesive type;
GOST R ISO 15223-1-2014 — Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied.

Part 1. General requirements;

GOST ISO 10993-1-2011 — Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and testing;

GOST ISO 10993-2-2009 — Medical devices. Biological evaluation of medical devices. Part 2. Animal welfare requirements;

GOST ISO 10993-5-2011. — Medical devices. Biological evaluation of medical devices. Part 5. Tests for *In Vitro* Cytotoxicity;

GOST ISO 10993-10-2011. — Medical devices. Biological evaluation of medical devices. Part 10. Tests for irritation and delayed-type hypersensitivity;

GOST ISO 10993-11-2011. — Medical devices. Biological evaluation of medical devices. Part 11. Tests for systemic toxicity;

GOST ISO 10993-12-2015. — Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and reference materials;

GOST R 52770-2016 — Medical devices. Safety requirements. Methods of sanitary and chemical and toxicological studies;

GOST 31214-2016. — Medical devices. Requirements for samples and documentation presented for toxicological tests, sanitary and chemical analyses, tests for sterility and pyrogenicity.

Information on the sterility of the medical device, sterilization method and procedures, if the sterile package compromised or, if the medical device supplied non-sterile, indication to sterilize it prior to use.

The SANIPLAST® First Aid Bandage is supplied non-sterile and does not require sterilization before use.

Information on proper processing the medical device for reuse, including cleaning, disinfection, packaging and, if necessary, a re-sterilization method (if the medical device intended for multiple use), as well as non-usability criteria for the medical device;

The SANIPLAST® First Aid Bandage is not reusable and does not require re-sterilization.

Information on restrictions on the joint use of medical devices (for medical devices intended for use with other medical devices and (or) with general-purpose devices)

SANIPLAST® First Aid Bandage is not intended for mixed use with other medical devices.

Counter-indications, warnings and precautions

Hypersensitivity, dermatitis.

Side effects

No side effects are known to date. Hypersensitivity to components cannot be ruled out in each specific case.

Residual risks

Not detected.

Warnings and/or precautions taken by the consumer when disposing of the medical device, accessories and consumables used therewith (if applicable), including information on infection, microbial, ecological or physical dangers posed by the medical device;

When using a medical device (waste disposal), disposal and/or recycling of materials should be performed in accordance with the current legislation.

The product should be disposed as Class A waste according to SanPiN 2.1.3684-21 (epidemiological safe waste, similar to solid waste composition).

Information about the circumstances where the user should consult a medical professional (for medical devices intended for the persons without medical education)

In case of an allergic skin reaction or deep cut, consult your physician.

Procedure and conditions of disposal or destruction of the medical device

When using a medical device (waste disposal), disposal and/or recycling of materials should be performed in accordance with the current legislation.

The product should be disposed as Class A waste according to SanPiN 2.1.3684-21 (epidemiological safe waste, similar to solid waste composition).

Information about non-conforming event that have the signs of an adverse event (incident) is accepted by an authorized representative of the manufacturer

REGICOM, LLC

The Russian Federation, 123290 Moscow, 39 Shmitovsky Driveway, building 2, apt. 10

Tel.+7 926 385 5611

Warranty statement

The manufacturer guarantees the manufactured product compliance with the requirements of technical documentation, subject to the proper conditions of transportation and storage.

The shelf life is four (4) years.

Information on the issue or the latest revision of the operating instructions

April, 2022

List of symbols of the product packaging

Ser. No.	Designation	Explanation
1		Non-sterile
2		Do not reuse
3		The device complies with the essential requirements of EU Directive

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