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20

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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of ZHEJIANG BANGLI MEDICAL PRODUCTS CO., LTD. on the annexed INSTRUCTION FOR USE is genuine.

The attached Russian text of the said DOCUMENT is true to the English text.

China Council for the Promotion
of International Trade

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
Name and description of the medical product

Surgical adhesive strip-tape bandage Reneplast® is a non-woven rectangular with adhesive ends and adsorbing pad in the center. The adhesive layer of the dressing is protected by a thin anti-adhesive siliconized paper, which is easily removed before use.

The dressing consists of several layers, shown below in Figure 1.

Figure 1.

The top layer consists of soft non-woven base material with an adhesive layer. The base material with the adhesive layer of the dressing does not come into contact with the wound surface.

The lower layer (skin-contacting) is a pad of absorbent material, covered with a microscopic film, due to which the bandage does not stick to the wound.

Anti-adhesive siliconized paper adheres to the adhesive layer with a smooth surface that makes its easy removal before use.

Once the protective liner is removed, the bandage is ready for use.

Versions

The self-adhesive postoperative bandage of the plaster type Reneplast® is presented in the following versions:

- 1 Surgical adhesive strip-tape bandage Reneplast®, 5x7 cm;
- 2 Surgical adhesive strip-tape bandage Reneplast®, 9x10 cm;
- 3 Surgical adhesive strip-tape bandage

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- Reneplast®, 10x15 cm;
 4 Surgical adhesive strip-tape bandage Reneplast®, 10x20 cm;
 5 Surgical adhesive strip-tape bandage Reneplast®, 10x25 cm.

Composition of the bandage

Description	Name of material	Purpose of material
Self-adhesive postoperative dressing of Reneplast® plaster type	Non-woven material (viscose 90 %), dacron 10 % (polyethylene terephthalate) with polyethylene coating (100 % polyethylene)	Pad with a microscopic film
	Non-woven material (viscose 70 %), Dacron 30 % (polyethylene terephthalate)	Base material
	White oil 75 %, rosin 15 %, rubber 10 %	Adhesive layer
	Anti-adhesion siliconed paper	Protective liner

Purpose and principle of the medical product

The self-adhesive postoperative dressing of the Reneplast® plaster type is sterile and intended for the protection of wounds after surgery and traumatic skin lesions with mild or moderately severe exudation.

The absorbent pad, made in the center of the product is designed atraumatic (does not irritate the wound and does not stick to it), removes excess exudate and prevents the risk of wound infection when it is located directly on postoperative wounds or minor skin lesions. It removes excess exudate and prevents the risk of wound infection when it placed directly upon the postoperative wounds or minor skin lesions. Non-woven base reliably fixes the bandage and provides high breathability. The

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bandage is painlessly removed from the skin, isolates the wound from the environment and protects against mechanical irritations.

Indications

Protection and care of postoperative surgical wounds, traumatic skin lesions.

Counterindications

Individual sensitivity to the components of the bandage (the presence of predisposition to allergy).

Method of use

Open the package, remove the protective coating of the dressing from its surface and, put it upon the wound without touching the pad with fingers, stick the dressing to dry clean skin. In this case, the pad must completely cover the wound.

The dressing should be removed with a gentle movement over hair growth.

The bandage is recommended to be on the skin no more than 24 hours, after which it must be removed and, if necessary, replaced with a new one.

Use and potential consumer

The product is intended for one-time individual use at home, as well as in inpatient and outpatient medical organizations.

No special knowledge is required to use a self-adhesive postoperative Reneplast® plaster type bandage.

Potential undesired effects

- Skin redness in the application area;
- skin itching.

In case of such situations, immediately stop using the bandage and consult your doctor.

Precautions

1. Do not apply the bandage to the areas with skin diseases, open wounds, ulcers.

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2. Use immediately after pack opening.

3. Do not use if the pack is opened or damaged.

4. The bandage should be in contact with the skin no more than 24 hours.

5. Do not reuse it.

Pack size

Option 1

1 bandage 5x7 cm or 9x10 cm, or 10x15, or 10x20 cm, or 10x25 is packed in the individual strip paper pack (consumer package).

Option 2

5 or 10 bandages 5x7 cm or 10 bandages 9x10 cm, or 10x15, or 10x20 cm, or 10x25 are packed in the individual strip paper pack and inserted in the carton package (consumer package).

Storage conditions

Store in the original packaging at a temperature of + 5 °C to + 25 °C and a relative humidity of not more than 80 %.

Keep out of the reach of children.

Shelf life

3 years.

Do not use after expiration date.

Environmental requirements

When using, transporting, storing and disposing of the bandage according to the operational documentation, negative effects on humans and the environment are excluded.

Information about disinfection and cleaning methods

The bandage is individual and single used. Re-use is not provided, so no cleaning or disinfection is required.

Sterilization

The bandages are supplied sterilized.

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Method of sterilization: sterilization with ethylene oxide.

Disposal and destruction procedure

Dispose of the used bandage in accordance with SanPiN 2.1.3684-21 as class A waste (epidemiologically safe waste, similar in composition to solid municipal waste).

Maintenance and repair requirements

Since the bandage is an individual and single-use product, re-use is not provided, thus, it is not subject maintenance and repair.

Warranty

The manufacturer, Zhejiang Bangli Medical Products Co., Ltd., warrants the compliance of "Surgical adhesive strip-tape bandage Reneplast®" with this document provided that the use, transportation and storage conditions are complied with.

The manufacturing company Jiangsu Nanfang Medical Co., Ltd. is not responsible for damage or accidents caused by improper use of the product, due to a fault of customers or any external unforeseen situations.

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