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CERTIFICATE

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International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of ZHEJIANG BANGLI MEDICAL PRODUCTS CO., LTD. on the annexed INSTRUCTION FOR USE is genuine. The attached Russian text of the said DOCUMENT is true to the English text.

China Council for the Promotion
of International Trade

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Name and description of the medical product

Reneplast® medical adhesive plaster (hereinafter referred to as adhesive plaster) is a fabric base with an adhesive layer and a pad of absorbent pad in the center. The adhesive base of the adhesive plaster is protected with siliconed paper, which is easily removed before use.

Designs

Reneplast® medical band-aid, versions:

1. Medical band-Aid Reneplast® size 72x19 mm.

Formulation: fabric base (60 % – polyester, 40 % – viscose), adhesive mass (75 % – white oil, 15% – rosin, 10 % – rubber), absorbent pad (90 % – viscose, 10 % – dacron), anti-adhesive siliconed paper.

Purpose and mode of action of medical product

Adhesive plaster protects minor skin damage: (scratches, abrasions, cuts) from external infection and dirt. An absorbent pad ensures that the wound surface fluid (exudate) is absorbed, and also creates favorable conditions for natural healing of the damaged skin area. Using an adhesive layer, the adhesive plaster adheres tightly to skin and is securely fixed on the body surface. The plaster is single-use.

Indications

Protection of minor skin damages: abrasions, cuts, scratches and other skin injuries.

Contradictions

Hypersensitivity to the components of the product (predisposition to allergy).

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Method of use

Clean and dry the skin in the damage area. It is advisable to additionally treat the damaged skin area with an antiseptic solution. Do not remove damaged skin!

Before use, open the package on the edge and remove the adhesive plaster strip. Remove the protective paper from the adhesive plaster surface and, without touching the pad with your fingers, place it on the wound, and apply the adhesive base to dry clean skin.

Avoid contact of the adhesive base with the damage area, and the absorbent pad must completely cover the wound.

If necessary, remove the adhesive plaster by gently moving in the direction of hair growth. Adhesive plaster is recommended to be changed at least twice a day.

Apply adhesive plaster regularly until, the damaged skin is completely repaired.

Adhesive plaster is not waterproof, therefore, when performing hygienic procedures it requires additional protection from moisture.

Conditions of use and potential consumer

Adhesive plaster can be used at home, as well as in inpatient and outpatient medical institutions. No special knowledge is required to use Reneplast® adhesive plaster.

Potential side effects

- allergic reactions;
- contact dermatitis.

In case of side effects, immediately stop using the adhesive plaster and consult your doctor.

Interaction with other medical products.

Do not use with other external products on the same skin areas.

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Precautions

1. Adhesive plaster is for external use only.
2. Do not reuse.
3. Avoid contact of the adhesive base with a damage surface.
4. Do not use adhesive plaster with a damaged protective coating;
5. Do not allow adhesive plaster to be applied or removed.
6. If the damage persists or application of adhesive plaster brings discomfort, stop using the product and consult your doctor.

Pack size

Reneplast® medical adhesive plaster is put into individual non-cell does-tray paper packaging (consumer packaging).

The number of adhesive plasters in individual packaging:

- 72x19 mm – 1 pc.

Medical adhesive plaster Reneplast® in individual packaging is packed in a cardboard box (consumer packaging):

- 72x19 mm – 10 pc.

- 72x19 mm – 20 pc.

Storage conditions

Store in the original packaging at a temperature of + 5 °C to + 25 °C and a relative humidity of not more than 80%. Keep out of the reach of children.

Shelf life

3 years.

Do not use after expiration date.

Environmental requirements

When using, transporting, storing and disposing of the plaster according to the terms of the operational documentation, negative effects on humans and the environment are excluded.

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Disposal and destruction procedure

Dispose of used adhesive plaster in accordance with SanPiN 2.1.3684-21 as class A waste (epidemiologically safe waste, similar in composition to solid municipal waste).

Sterilization and Reuse

The plasters are supplied non-sterile and are not to be sterilized. Adhesive plasters are for single use, and shall not be reused.

Maintenance and repair requirements

Since the plaster is individual and single-use product, re-use is not provided, and it cannot be maintained or repaired.

Warranty

The manufacturing company Jiangsu Nanfang Medical Co., Ltd., guarantees the Reneplast® Medical adhesive plaster compliance with the requirements of this document, subject to the conditions of use, transportation and storage. The manufacturing company Jiangsu Nanfang Medical Co., Ltd. is not responsible for damage or accidents caused by improper use of the product, due to a fault of customers or any external unforeseen situations.

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