
Subscribed and sworn to before me
this date 19 October 2021

Notary Public

STEPHEN G. NESBITT
Attorney at Law
Notary Public, State of Ohio
My Commission Has No Expiration
Date. Section 147.03 O.R.C.

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
М.П. / Stamp

INSTRUCTIONS FOR USE OF THE MEDICAL DEVICE

Circular stapler with staples and adjustable height of their closure
ECHELON CIRCULAR®, electric, autonomous in versions

manufactured by
Ethicon Endo-Surgery, LLC,
475 Calle C, Guaynabo, Puerto Rico 00969, USA

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2021 year / г.

1. MEDICAL DEVICE (MD) NAME AND CLASSIFICATION

1.1 MD name, models or variations of MD

Circular stapler with staples and adjustable height of their closure ECHELON CIRCULAR®, electric, autonomous in versions:

1. ECHELON CIRCULAR® Powered Stapler, 23 mm (CDH23P), contents:
 - 1.1. Circular stapler with staples and adjustable height of their closure ECHELON CIRCULAR®, electric, autonomous, model CDH23P – 3 pc.
 - 1.2. Battery Pack (model D34924, manufacturer Ethicon Endo-Surgery LLC) – 3 pc.
 - 1.3. Instruction for use – 1 pc.
2. ECHELON CIRCULAR® Powered Stapler, 25 mm (CDH25P), contents:
 - 2.1. Circular stapler with staples and adjustable height of their closure ECHELON CIRCULAR®, electric, autonomous, model CDH25P – 3 pc.
 - 2.2. Battery Pack (model D34924, manufacturer Ethicon Endo-Surgery LLC) – 3 pc.
 - 2.3. Instruction for use – 1 pc.
3. ECHELON CIRCULAR® Powered Stapler, 29 mm (CDH29P), contents:
 - 3.1. Circular stapler with staples and adjustable height of their closure ECHELON CIRCULAR®, electric, autonomous, model CDH29P – 3 pc.
 - 3.2. Battery Pack (model D34924, manufacturer Ethicon Endo-Surgery LLC) – 3 pc.
 - 3.3. Instruction for use – 1 pc.
4. ECHELON CIRCULAR® Powered Stapler, 31 mm (CDH31P), contents:
 - 4.1. Circular stapler with staples and adjustable height of their closure ECHELON CIRCULAR®, electric, autonomous, model CDH31P – 3 pc.
 - 4.2. Battery Pack (model D34924, manufacturer Ethicon Endo-Surgery LLC) – 3 pc.
 - 4.3. Instruction for use – 1 pc.

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(Hereinafter, the devices may be referred to as “MD”, “Device”, “ECHELON CIRCULAR® Powered Staplers”, “circular stapler, model CDH23P” (etc.), “stapler model CDH23P” (etc.), “model CDH23P” (etc.), “CDH23P” (etc.), “Vegas” and etc.

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1.2 MD classification

Table 1.

Classification	Class
Sterility	Sterile (ethylene oxide sterilization)
Scope of use	Surgery

2 INTENDED USE AND PRINCIPLE OF OPERATION

MD is intended for use during open or laparoscopic surgical procedures on the digestive tract to create anastomoses (for example, end-to-end, end-to-side or side-to-side).

Device Description

The ECHELON CIRCULAR® Powered Staplers are anastomotic staplers available in four sizes to permit proper matching of device to diameter of the lumen.

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Each device has a detachable anvil that allows the surgeon to place the anvil in the desired location. The devices allow the surgeon to control tissue compression by varying the height of the closed staple. The devices have been designed to facilitate insertion, operation, and removal. The device is designed to prevent the device from being fired with an unsecure anvil and/or outside the device's firing range.

Each device is packaged with a primary lithium Battery Pack that must be installed prior to use. There are specific requirements for disposing of the Battery Pack. Refer to Battery Pack Disposal for Battery Pack disposal instructions.

A Staple Retaining Cap on each device's head protects the staple leg points during shipping and transporting.

ECHELON CIRCULAR® Powered Stapler Product Codes Table

Product Code	Open Staple Leg Length	Adjustable Closed Staple Height	Staple Housing Outside Diameter	Knife Outside Diameter	Number of Staples	Staple Rows
CDH23P	5,2 mm	1,5-2,2 mm	23 mm	14,6 mm	20	2
CDH25P	5,2 mm	1,5-2,2 mm	25 mm	16,5 mm	22	2
CDH29P	5,2 mm	1,5-2,2 mm	29 mm	20,4 mm	26	2
CDH31P	5,2 mm	1,5-2,2 mm	31 mm	22,5 mm	30	2

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4. INFORMATION ON POTENTIAL USERS OF MD AND DIRECTIONS FOR USE

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
Illustration and Nomenclature (Illustration 1)


1. Orange Band
2. Device Trocar
3. Staple Housing
4. Knife Outside Diameter
5. Tissue Compression Scale
6. Adjusting Knob
7. Battery Pack

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








Instructions for Use
Prior to Using Device

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1. Verify compatibility of all devices and accessories prior to using the device.

Preparing the Device for Use

2. Using sterile technique, remove the device and battery from the package. To avoid damage, do not flip or drop the device or battery into the sterile field.

3. Install the Battery Pack. The Battery Pack must be installed prior to use. Insert the Battery Pack by lining up the release tabs on the Battery Pack with the slots on the rear of the device. (Illustration 2). Ensure Battery Pack is fully inserted into the device. An audible click will be heard and the Tissue Compression Scale (Illustration 1-5) backlight will be illuminated when the Battery Pack is fully inserted.

CAUTION: The device must be used within 12 hours of inserting the Battery Pack. Failure to use within this time frame may cause the battery to be depleted. **Refer to Battery Pack Disposal** for Battery Pack disposal instructions.

4. Turn the Adjusting Knob counter-clockwise to remove the Staple Retaining Cap (Illustration 1-10) which contains the Ancillary Trocar (Illustration 1-12) and is attached to the anvil (Illustration 3). If use of the Ancillary Trocar is anticipated, do not discard the Staple Retaining Cap. Keep the retainer in the sterile field for later use. **Refer to Using the Ancillary Trocar.**

Using the Device

CAUTION: Do not inadvertently release the Red Safety or Firing Trigger until ready to fire the device. Releasing the Red Safety removes protection from premature firing which will compromise the ability to select the desired staple height.

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5. The device can be used with either of the following two techniques: A) the Open Lumen Purse String Technique, or B) the Closed Lumen Stapling Technique (double or triple stapling technique). (Illustration 4)

6. Open the device by turning the Adjusting Knob counter clockwise (Illustration 1-6) until the Anvil Shaft (Illustration 1-15) is fully exposed. Remove the anvil to expose the Device Trocar (Illustration 1-2). Retract the Device Trocar by rotating the Adjusting Knob clockwise until a stop is reached. Check the Device Trocar to verify that it is fully retracted before proceeding. (Illustration 5)

The device may also be inserted without removing the anvil if the preferred application is a sutured purse-string technique. In this case, however, prior to insertion, ensure the anvil is properly attached and the device is fully closed by rotating the Adjusting Knob clockwise.

7. Insert the anvil into the lumen using either technique ensuring the tissue is located at the Suture Tying Area. (Illustration 1-14, 6)

8. With the anvil removed and the Device Trocar fully retracted, insert the device up to the closed lumen. (Illustration 7A)

9. Fully extend the Device Trocar and pierce tissue by holding the device while gently rotating the Adjusting Knob counter-clockwise. Continue to push the tissue down until the Orange Band is visible. (Illustration 7B)

CAUTION: Keep the Device Trocar visible at all times to prevent injury or inadvertent trauma to adjacent structures. DO NOT use the Anvil Shaft to assist in piercing the tissue by placing it over the unexposed Device Trocar to avoid inclusion of tissue within the Anvil Shaft. DO NOT use it for piercing.

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10. Grasp the Anvil Grasping Area and reattach the anvil by sliding the Anvil Shaft over the Device Trocar and pushing until the anvil snaps into its fully seated position. The Orange Band on the Device Trocar will be covered when the anvil is fully attached. (Illustration 8). If necessary, hold the adjusting knob to prevent the Device Trocar from retracting during laparoscopic anvil attachment (Illustration 19).

CAUTION: DO NOT clamp across or grip on the Locking Springs (Illustration 1-16) when attempting to reattach the anvil as this will prevent anvil attachment.

CAUTION: Ensure the Anvil is attached by verifying the Orange Band is covered and the anvil is secure. The device will not fire with an unsecured anvil.

11. While closing the device, keep the organ segments in proper orientation. Inspect to ensure extraneous tissue is excluded. To close the device, turn the Adjusting Knob clockwise.

CAUTION: Ensure tissue thickness lays flat and is evenly distributed within the device. Excess tissue on one side may result in unacceptable staple formation and can result in staple line leakage.

12. As the final adjusting revolution is approached, the orange staple height indicator moves into the green range of the Tissue Compression Scale. (Illustration 9). (A: Green Range) (B: Orange Staple Height Indicator) Adjust the device until appropriate tissue resistance is felt for a secure anastomosis. Wait 15 seconds to allow for adequate tissue compression, and adjust if needed to maintain appropriate tissue resistance. Providing the orange staple height indicator falls fully within the green range of the Tissue Compression Scale upon firing, the device will form staples at the height indicated for the selected tissue compression.

Markings are provided in the Tissue Compression Scale to serve as reference

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points only. **Refer to the Product Codes Table** for adjustable closed staple height range.

13. Prefire checklist:

- Orange staple height indicator is fully within green range.
- Anvil is properly attached. (Illustration 8)

CAUTION: The device will not fire if the anvil is not properly attached or if the orange staple height indicator is not fully within the green range of the Tissue Compression Scale.

14. To fire the device, draw the Red Safety (Illustrations 1-8 and 10) back toward the handle. To ensure the device remains in the safe firing range, once the Red Safety has been released, do not turn the Adjusting Knob. Activate the firing sequence by completely depressing the Firing Trigger (Illustrations 1-9 and 11). Remain still until the full firing sequence is completed which will be indicated by the illuminated green check mark on the indicator window (Illustration 12).

It is not necessary to hold the trigger once the firing sequence has initiated. The user may notice audible feedback during the firing sequence when cutting through the breakaway washer. Once fired, the device cannot be fired again.

 **CAUTION:** Do not attempt to manipulate the Adjusting Knob during the firing sequence. Doing so may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.

CAUTION: If excessive force is required to close the device, this may indicate there is too much tissue or thickened tissue in the device. Attempting to fire the device in this condition may result in malformed staples,

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incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.

WARNING: If the firing cycle is not complete as indicated by the illuminated green check mark, the knife may remain exposed and staples may not be fully formed. Use caution when removing a device with an exposed knife. Since staples may not be fully formed, consider further reinforcement of staple lines.

WARNING: The device is designed to stop automatically at the completion of the stroke. If device fails to automatically stop, remove the battery. Use caution when removing the device. Failure to automatically stop may compromise the integrity of the anastomosis. Use caution when inspecting the anastomosis and consider further reinforcements of staple lines.

15. To safely release the device from the newly formed anastomosis, turn the Adjusting Knob counter-clockwise for two complete revolutions or four half turns. (Illustration 13)

16. To withdraw the open device, rotate it 90° in both directions taking care to minimize movement of the distal tip. This ensures the tissue is released. Gently pull out the device while simultaneously rotating. (Illustration 14).

CAUTION: If the device does not freely release from the anastomosis or withdraw easily after rotation, turn the adjusting knob counterclockwise one additional complete revolution (360°) and attempt removal again by rotating the device 90° in both directions, taking care to minimize movement of the distal tip. Gently pull out the device while simultaneously rotating.

17. To inspect the donuts, remove the anvil, breakaway washer (if present), and donuts from within the circular knife (Illustration 15). Examine the integrity of the donuts and white breakaway washer. Donuts should be intact and include all tissue layers. Ensure that donuts are complete, and the white

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breakaway washer is completely transected into inner and outer concentric rings. If donuts or breakaway washer are not complete, the anastomosis should be carefully checked for leakage and appropriate repairs made.

18. Dispose of all opened devices whether used or unused. This device is packaged and sterilized for single use only.

19. Dispose of the Battery Pack. See **Battery Pack Disposal**.

Using the Ancillary Trocar

The radiopaque Ancillary Trocar is designed to facilitate controlled penetration through tissue selected by the surgeon. For example, a surgeon may elect to use the Ancillary Trocar to perform a triple stapling technique or end-to-side anastomosis.

CAUTION: Injury hazard. Keep the Ancillary Trocar visible at all times to prevent personal injury or inadvertent trauma to adjacent structures.

Attaching the Ancillary Trocar to the Anvil

1. The Ancillary Trocar is shipped in the Staple Retaining Cap of the ECHELON CIRCULAR® Powered Stapler (Illustration 16).

2. To remove the Ancillary Trocar from the Staple Retaining Cap, grasp the Finger Notches (Illustration 1-11) and pull straight out, parallel to the retainer body.

3. Place the blunt end of the Ancillary Trocar into the Anvil Shaft. (Illustration 17) Rotate the Ancillary Trocar with respect to the Anvil Shaft approximately 45 degrees to ensure secure attachment (Illustration 18).

CAUTION: If Ancillary Trocar is not properly seated in the anvil, the Ancillary Trocar could detach from the anvil during manipulation.

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

designed to correspond with the heads of the Intraluminal Staplers (ILS), both Mechanical and Powered. Class I devices.

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8 LIST AND DESCRIPTION OF MD MATERIALS

Material information for permanent patient contact materials

Part Name	Generic Material Name
Materials in the device classified as implant device, tissue, permanent contact (greater than 30 days) as defined by ISO 10993-1	
Staple	Titanium Alloy

Material information for limited patient contact materials

Part Name	Generic Material Name
Materials in the device classified as external communicating device, tissue, limited contact (less than 24 hours) as defined by ISO 10993-1	
Distal End of Device Shaft with End-Effector Components	Styrene Ethylene Butylene Styrene Block Copolymer Acrylonitrile-Butadiene-Styrene (ABS) Stainless Steel Polyether Block Polyamide Copolymer Polycarbonate Polyetherimide, 30% Glass Fiber Nylon, 60% Glass Fiber Aluminum Polyurethane Oligomer Mixture Epoxy-Based Printing Ink Hardener Sodium Stearate

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


trauma to adjacent structures. DO NOT use the Anvil Shaft to assist in piercing the tissue by placing it over the unexposed Device Trocar to avoid inclusion of tissue within the Anvil Shaft, DO NOT use it for piercing.

- DO NOT clamp across or grip on the Locking Springs when attempting to reattach the anvil as this will prevent anvil attachment.
- Ensure the Anvil is attached by verifying the Orange Band is covered and the anvil is secure. The device will not fire with an unsecure anvil.

• Ensure tissue thickness lays flat and is evenly distributed within the device. Excess tissue on one side may result in unacceptable staple formation and can result in staple line leakage.

• The device will not fire if the anvil is not properly attached or if the orange staple height indicator is not fully within the green range of the Tissue Compression Scale.

• Do not attempt to manipulate the Adjusting Knob during the firing sequence. Doing so may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.

• If excessive force is required to close the device, this may indicate there is too much tissue or thickened tissue in the device. Attempting to fire the device in this condition may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.

• If the firing cycle is not complete as indicated by the illuminated green check mark, the knife may remain exposed and staples may not be fully formed. Use caution when removing a device with an exposed knife. Since staples may not be fully formed, consider further reinforcement of staple lines.

• The device is designed to stop automatically at the completion of the stroke. If device fails to automatically stop, remove the battery. Use caution when removing the device. Failure to automatically stop may compromise the integrity of the anastomosis. Use caution when inspecting the anastomosis and consider further reinforcements of staple lines.

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WARNING: Portable radio and mobile RF communications equipment can affect Medical Electrical Equipment. Follow the directions of the following Electromagnetic Emissions and Electromagnetic Immunity Guidance tables when using the ECHELON CIRCULAR® Powered Stapler.

WARNING: Avoid use of the ECHELON CIRCULAR® Powered Stapler adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, monitor the ECHELON CIRCULAR® Powered Stapler and the other equipment to assure normal operation.

Table of Toxic and Hazardous Substances / Elements and Their Content

As required by China's Management Methods for Controlling Pollution

By Electronic Information Products

Product Codes: ECHELON CIRCULAR® Powered Stapler series

Toxic and Hazardous Substances and Elements						
Part Name	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent Chromium (Cr(VI))	Polybrominated biphenyls (PBB)	Polybrominated diphenyl ethers (PBDE)
Handle	X	O	O	O	O	O
Shaft	O	O	O	O	O	O

This table is prepared in accordance with the regulations of SJ/T11364.

O: indicates that this toxic or hazardous substance contained in all of the

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homogenous materials for this part is below the limit requirement in GB/T 26572.

X: indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement in GB/T 26572.

• All parts named in this table are in compliance with the European union's RoHS legislation "Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of hazardous substances in electrical and electronic equipment."

The referenced Environmental Protection Use Period Marking was determined according to normal operating conditions of the product such as temperature and humidity.

Electromagnetic Emissions Guidance

Guidance and manufacturer's declaration—electromagnetic emissions		
The assembled ECHELON CIRCULAR Powered Stapler is intended for use in the electromagnetic environment specified below. The customer or the user of the assembled ECHELON CIRCULAR Powered Stapler should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
Harmonic emissions IEC 61000-3-2	Not applicable	The assembled ECHELON CIRCULAR Powered Stapler is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	
RF emissions CISPR 14-1	Complies	The assembled ECHELON CIRCULAR Powered Staplers are not suitable for interconnection with other equipment.

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Electromagnetic Immunity Guidance

For electromagnetic immunity, essential performance is no knife movement without activation of the Firing Trigger.

Guidance and manufacturer's declaration—electromagnetic immunity			
The assembled ECHELON CIRCULAR Powered Staplers are intended for use in the electromagnetic environment specified below. The customer or the user of assembled ECHELON CIRCULAR Powered Staplers should ensure that they are used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth		
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycles 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s		
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m		
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

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Guidance and manufacturer's declaration—electromagnetic immunity			
The ECHELON CIRCULAR Powered Stapler is intended for use in the electromagnetic environment specified below. The customer or the user of the ECHELON CIRCULAR Powered Stapler should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the ECHELON CIRCULAR Powered Stapler, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ECHELON CIRCULAR Powered Stapler is used exceeds the applicable RF compliance level above, the ECHELON CIRCULAR Powered Stapler should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ECHELON CIRCULAR Powered Stapler. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Recommended separation distances between portable and mobile RF communications equipment and the ECHELON CIRCULAR Powered Stapler			
The ECHELON CIRCULAR Powered Stapler is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ECHELON CIRCULAR Powered Stapler can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ECHELON CIRCULAR Powered Stapler as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

12 INFORMATION ABOUT MEDICINES AND (OR) PHARMACEUTICAL SUBSTANCES INCLUDED IN MD
N/A

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13. DESCRIPTION OF STERILIZATION METHOD

The ECHELON CIRCULAR® Powered Staplers are intended for single patient use and are not intended to be reused or resterilized. The ECHELON CIRCULAR® Powered Staplers are disposable devices intended for single patient use and shall be sterilized by Ethylene Oxide (EO) to a Sterility Assurance Level (SAL) of 10^{-6} .

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14 REQUIREMENTS FOR TECHNICAL SUPPORT AND REPAIR OF MD

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14.1 Applicable Standards**HARMONIZED STANDARDS LIST (EU)**

EN ISO 13485:2016
EN ISO 14971:2012
EN ISO 15223-1:2016
EN 1041: 2008
EN ISO 16061:2015
EN ISO 14630:2012
EN 62366:2008
EN 60601-1:2006/A1:2013/A12:2014
EN 60601-1-2:2015
EN 60601-1-6:2010/A1:2015
EN 60601-1-9:2008
EN 60086-4:2015
EN 50581:2012
EN 62321-1:2013
EN ISO 10993-1:2009
EN ISO 10993-2:2006
EN ISO 10993-3:2014
EN ISO 10993-5:2009

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IDENTIFICATION OF APPLICABLE STANDARDS

EN 62366:2015 (IEC 62366:2007 1st Ed) Project: IEC 62366 2nd Ed.
EN ISO 11737- 1:2006/AC:2009/Cor:2011) (ISO 11737-1:2006/Cor 1:2007/Cor:2011)
EN ISO 11737-2:2009 (ISO 11737-2:2009 2nd Ed)
EN ISO 14644-1:1999 (ISO 14644-1:1999)
EN ISO 14644-2:2000 (ISO 14644-2:2000)
EN ISO 14644-5:2004 (ISO 14644-5:2004)
EN 556-1:2001/AC:2006
EN ISO 11607-1:2009 (ISO 11607-1:2006) Work Item: ISO 11607-1:2006/NP Amd 1 EN ISO 11607-1:2009/prA1
EN ISO 11607-2:2006 (ISO 11607-2:2006) Work Item: EN ISO 11607-2:2006/prA1 ISO 11607-2:2006/NP Amd 1
EN ISO 11135-1:2007 (ISO 11135-1:2007 1st Ed)

EN ISO 10993-6:2016
EN ISO 10993-10: 2013
EN ISO 10993-11:2018
EN ISO 10993-12:2012
EN ISO 10993-18:2009
EN 556-1:2001
EN ISO 11737-1:2018
EN ISO 11737-2:2009
EN ISO 11135-1:2014
EN ISO 10993-7:2008
EN ISO 11607-1: 2017
EN ISO 11607-2: 2017
EN ISO 14644-1:2015
EN ISO 14644-2:2015
EN ISO 14644-5:2004

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15 INFORMATION ABOUT IMPLANT REMOVAL AND EXTRACTION

Once the staple has completed its functional lifetime, the device remains in the patient for the lifetime of the patient or may be expelled after healing during normal cellular renewal or, if the physician deems it necessary the implant may be removed either intra-operatively or post-operatively. The clinical safety results of the literature assessed over the period of 1 January 1991 to 1 June 2018 showed no events of either staple loosening or staple being expelled. Permanent surgical titanium alloy staples are well-established technologies with extensive history of safety over the lifetime of the patient.

Based on product testing using accelerated aging, the ECHELON CIRCULAR® Powered Stapler has a shelf life of 3 years.

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16 CONDITIONS OF TRANSPORTATION AND STORAGE

Table 5. Transportation conditions:

Designation	Range
Temperature	-22 °C – 60 °C
Relative humidity	10% – 80%
Barometric pressure	500 – 1060 hPa

Keep Dry

Keep Away from Heat

Table 6. Storage conditions:

Designation	Range
Temperature	-20 °C – 30 °C
Relative humidity	10% – 65%
Barometric pressure	N/A

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17 MD DISPOSAL RULES AND CONDITIONS

Dispose of the device in accordance with applicable local/regional/national/international regulations

- At the end of product life, dispose of all components in accordance with your governmental regulations.
- Dispose of batteries properly in accordance with your local, state, and/or government requirements.

Medical devices in Russian Federation must be disposed of in the local public entities in accordance with SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste".

After the use the product and packaging should be disposed of in accordance with hospital, administrative and / or approved local national standards (regulations) on disposal of infected or potentially infected materials.

Unused product is utilized in accordance with the handling of medical waste of class A.

The used product that was in contacted with biological fluids could represent epidemiological risk and must be disposed of in accordance with the rules of handling medical waste of class B.

Battery Pack Disposal

The Battery Pack contains a built-in battery drain and must be installed into the device to trigger the automatic battery drain feature. The battery does not need to remain in the device – it will continue to drain after it has been removed from the device. Once removed, depending on your local regulations, the battery can be put directly into the battery recycling bin or normal waste stream.

Batteries used in this device contain 1,2-dimethoxyethane.

Prior to Installation into the Device

If the Battery Pack needs to be disposed of prior to installation into device (e.g. product is beyond the Use by Date indicated on the package, or the

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Battery Pack is dropped), first install the Battery Pack into the device, then remove, to trigger built-in battery drain.

After Device Use

The Battery Pack must be removed from the device prior to disposal.

If the Battery Pack requires decontamination prior to disposal, follow hospital protocol or the Battery Pack Cleaning and Disinfection instructions.

To remove the Battery Pack, squeeze the release tabs and pull the battery pack straight back. Do not disassemble the Battery Pack.

Battery Pack Cleaning and Disinfection

WARNING: Do not use autoclaves, ethylene oxide, or radiation to sterilize or disinfect the Battery Pack.

Manual Cleaning

1 Remove the Battery Pack from the device before cleaning. The Battery Pack should not be submerged in water or cleaning solutions.

2 Clean the Battery Pack surfaces with a neutral pH detergent or neutral pH enzymatic detergent, prepared according to the manufacturer's instructions.

3 Use soft bristle brush to manually clean the Battery Pack with the cleaning solution.

4 Ensure areas containing crevices are scrubbed thoroughly.

5 Wipe off detergent thoroughly with lukewarm tap water.

6 Perform visual inspection to determine if debris is removed.

7 Repeat cleaning as necessary to obtain a visually clean Battery Pack.

8 Ensure excess water has drained from the battery before final disposal.

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