
K 11

INSTRUCTION FOR USE
THE NAVIGATION CATHETER FOR CARDIAC ABLATION:
THERMOCOOL SF NAV UNI-DIRECTIONAL AND BI-DIRECTIONAL CATHETERS

Issued on behalf of
Biosense Webster, Inc.

Rey Jacinto
Senior Regulatory Affairs Specialist

6/13/17
Date

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy or validity of that document.


State of California)
County of Orange)

On June 13th, 2017 before me, Richard G. Perez, NOTARY PUBLIC, personally appeared Rey Jacinto, who proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that she executed the same in her authorized capacity and that her signature on the instrument the person, or the entity upon behalf of which the person acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature: 
Richard G. Perez
Notary Public

ThermoCool SF NAV Uni-Directional and Bi-Directional Catheters

- **STERILE.** Sterilized using ethylene oxide.
- For single use only.
- Do not resterilize.
- Do not use if the package is open.
- Do not use if the package is damaged.

Device Description

THERMOCOOL SF NAV Uni-directional Catheter

The THERMOCOOL SF NAV Uni-directional Catheter is a steerable multi electrode (6), luminal catheter with a deflectable tip via moveable thumb knob, designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) current to the catheter tip dome electrode for ablation purposes. The catheter shaft measures 7.5 F with 7.5 F ring electrodes. The catheter is introduced through a guiding sheath with a minimum size of 8 F. For ablation, the catheter is used in conjunction with an RF generator and a dispersive pad (indifferent electrode).

The catheter has a high-torque shaft with a braided deflectable tip section containing an array of ring electrodes. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site. All ring electrodes are manufactured from noble metals and, unless otherwise noted, may be used for recording and stimulation purposes.* The tip dome electrode serves to deliver RF current from the RF generator to the desired ablation site. The connector type of the catheter is 34-pin for connection to the PIU of the CARTO 3 Navigation System (version 2.3 and higher) and 25-pin for connection to the PIU of the CARTO XP Navigation System.

Two proximal electrodes (electrode 5 and 6) on the catheter cannot be used for stimulation and recording purpose.

The catheter incorporates a temperature sensor that is embedded in the 3.5 mm tip dome electrode. The catheter is available in B, D, F, J curves.

At the proximal end of the catheter, a saline input port with a standard Luer fitting terminates from the open lumen. This saline port serves to permit the injection of normal saline to irrigate the tip dome electrode. During ablation, heparinized normal saline is passed through the internal lumen of the catheter and through the multi-holed, irrigated, tip dome electrode. The saline irrigates and cools the ablation site as well as the electrode tip. A compatible irrigation pump is used to control the saline irrigation.

The catheter features a magnetic location sensor embedded in the tip dome electrode that transmits location information to a compatible CARTO EP Navigation System. An appropriate reference device is required for location reference position purposes. The catheter interfaces with standard recording equipment and a compatible RF generator via accessory extension cables with the appropriate connectors.

For use in mapping procedures, for information on appropriate reference devices, and for further description of the operation of the CARTO EP Navigation System, refer to the User Manual for your CARTO EP Navigation System.

For further description of the operation of the irrigation pump and the RF generator, refer to the operating instructions for these instruments.

THERMOCOOL SF NAV Diagnostic/Ablation Deflectable Tip Catheter (Bi-directional Catheter)

The THERMOCOOL SF NAV Diagnostic/Ablation Deflectable Tip Catheter is a steerable multi-electrode (6), luminal catheter with a deflectable tip via rocker lever, designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) current to the catheter tip dome electrode for ablation purposes. The catheter shaft measures 7.5 F with 7.5 F ring electrodes. The catheter is introduced through a guiding sheath with a minimum size of 8 F. For ablation, the catheter is used in conjunction with an RF generator and a dispersive pad (indifferent electrode). The catheter has a high-torque shaft with a braided deflectable tip section containing an array of ring electrodes. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site. All ring electrodes are manufactured from noble metals and, unless otherwise noted, may be used for recording and stimulation purposes.* The tip dome electrode serves to deliver RF current from the RF generator to the desired ablation site. The connector type of the catheter is 34-pin for connection to the PIU of the CARTO 3 Navigation System (version 2.3 and higher) and 25-pin for connection to the PIU of the CARTO XP Navigation System.

Two proximal electrodes (electrode 5 and 6) on the catheter cannot be used for stimulation and recording purpose. The catheter incorporates a temperature sensor that is embedded in the 3.5mm tip dome electrode. The catheter is available in BB, DD, FF, JJ, BD, FJ, DJ, BF, DF curves.

At the proximal end of the catheter, a saline input port with a standard luer fitting terminates from the open lumen. This saline port serves to permit the injection of normal saline to irrigate the tip dome electrode. During ablation, heparinized normal saline is passed through the internal lumen of the catheter and through the multi-holed, irrigated, tip dome electrode. The saline irrigates and cools the ablation site as well as the electrode tip. A compatible irrigation pump is used to control the saline irrigation.

The catheter features a magnetic location sensor embedded in the tip dome electrode that transmits location information to a compatible CARTO EP Navigation System. An appropriate reference device is required for location reference position purposes. The catheter interfaces with standard recording equipment and a compatible RF generator via accessory extension cables with the appropriate connectors.

For use in mapping procedures, for information on appropriate reference devices, and for further description of the operation of the CARTO EP Navigation System, refer to the User Manual for your CARTO EP Navigation System.

For further description of the operation of the irrigation pump and the RF generator, refer to the operating instructions for these instruments.

Indications and Use

The THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter and related accessories are indicated for catheter-based cardiac electrophysiological mapping (stimulating and recording) and, when used in conjunction with a compatible radiofrequency generator, for cardiac ablation.

The THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter provides location information when used with CARTO EP Navigation Systems. The THERMOCOOL SF NAV Catheters is not compatible with CARTO 3 EP Navigation Systems prior to Version 2.3.

Contraindications

Do not use this catheters:

1. If the patient has had a ventriculotomy or atriotomy within the preceding eight weeks because the recent surgery may increase the risk of perforation.
2. In the patient with a myxoma or an intracardiac thrombus as the catheter could precipitate an embolus.
3. In patients with prosthetic valves as the catheter may damage the prosthesis.
4. In the coronary vasculature due to risk of damage to the coronary arteries.
5. In patients with an active systemic infection because this may increase the risk of cardiac infection.
6. Via the transseptal approach in a patient with an interatrial baffle or patch because the opening could persist and produce an iatrogenic atrial shunt.
7. Via the retrograde trans-aortic approach in patients who have had aortic valve replacement.

Warnings and Precautions

1. Do not use the temperature sensor to monitor tissue temperature. The temperature sensor located within the tip section of the catheter does not reflect either electrode-tissue interface or tissue temperature due to the cooling effects of the saline irrigation of the electrode. The temperature displayed on the RF generator is the temperature of the cooled electrode, not tissue temperature. The temperature sensor is used to verify that the irrigation flow rate is adequate. Before initiating the application of RF current, a decrease in electrode temperature confirms the onset of saline irrigation of the ablation electrode. Monitoring the temperature from the electrode during the application of RF current ensures that the irrigation flow rate is being maintained.
2. Do not rely on electrode temperature rise to determine if tissue heating is occurring during RF energy delivery as bench and animal studies showed no significant electrode temperature rise during RF ablation.
3. It is important to carefully follow the power titration procedure as specified in the instructions for use. Too rapid an increase in power during ablation may lead to perforation caused by steam pop.
4. This catheter may damage the prosthetic tricuspid valve of a patient if the catheter is accidentally advanced through the valve.
5. The patient who has had a prior atrial flutter ablation procedure may be at greater risk for perforation and/or pericardial effusion with the use of this catheter system.
6. In accordance with your hospital's protocol, monitor the patient's fluid balance throughout the procedure to avoid fluid volume overload. Some patients may have factors that reduce their ability to handle the volume overload, making them susceptible to developing pulmonary edema or heart failure during or after the procedure. Patients with congestive heart failure or renal insufficiency, and the elderly are particularly susceptible. Prior to the procedure, always identify the patient's risk of volume overload.
7. The safety of discontinuing anticoagulation therapy following catheter ablation of atrial fibrillation has not been established; anticoagulation therapy in such patients should be administered in accordance with the ACC/AHA/ESC Guidelines for the Management of Patients in Atrial Fibrillation.
8. To avoid thromboemboli, intravenous heparin should be used when entering the left heart during ablation. Follow the clinical guideline and common practice for optimal post-procedural anticoagulation.
9. The safety and effectiveness of radiofrequency ablation for the treatment of atrial fibrillation in patients with significant left ventricular dysfunction, advanced heart failure, substantial left atrial enlargement, and structural heart disease have not been established.
10. The catheter has not been shown to be safe at electrode temperatures above 40° C. Verify that the CATHETER SELECTION KNOB on the compatible RF generator is on the "Tcool SF" or THERMOCOOL SF similar option and ensure that the maximum temperature is set at 40°C.
11. Implantable pacemakers and implantable cardioverter/defibrillator (ICDs) may be adversely affected by RF current. It is important to have temporary external sources of pacing and

defibrillation available during ablation and to temporarily reprogram the pacing system to minimum output or OFF mode to minimize the risk of inappropriate pacing. Exercise extreme caution during ablation when in close proximity to atrial or ventricular permanent leads; program the ICD to the OFF mode during the ablation procedure; and, perform complete implantable device analysis on all patients after ablation.

12. Patients undergoing septal accessory pathway ablation are at risk for complete AV block which requires the implantation of a permanent pacemaker. Patients who experience inadvertent complete AV block as a result of RF ablation may also require permanent pacing.

13. During the trans-aortic approach, adequate fluoroscopic visualization is necessary to avoid placement of the catheter in the coronary vasculature. Intracoronary placement of the ablation catheter, RF energy application, or both have been associated with myocardial infarction.

14. When ablating near adjacent anatomical structures, take precautions to minimize collateral damage to the adjacent structures.

15. When ablating near the esophagus (along the posterior wall of the left atrium), take precautions to avoid injuring the esophagus, including appropriately reducing RF power.

16. When ablating near the phrenic nerve, take precautions to avoid injuring the phrenic nerve, including appropriately reducing RF power and pacing to identify the proximity of the nerve.

17. Minimize X-ray exposure during the procedure. Catheter ablation procedures present the potential for significant X-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the X-ray beam intensity and duration of the fluoroscopic imaging. Catheter ablation should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps have been taken to minimize this exposure. Careful consideration must therefore be given for the use of the device in pregnant women.

18. Do not expose catheter to organic solvents such as alcohol.

19. Do not autoclave the catheter.

20. Do not immerse proximal handle or cable connector in fluids; electrical performance could be affected.

21. Do not scrub or twist the distal tip electrode during cleaning.

22. Inspect irrigation saline for air bubbles prior to its use in the procedure. Air bubbles in the irrigation saline may cause emboli.

23. Purge catheter and irrigation tubing with heparinized normal saline.

24. Electrophysiology catheters and systems are intended for use only in X-ray shielded rooms due to electromagnetic compatibility requirements and other hospital safety guidelines.

25. Do not attempt to operate The Biosense Webster THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter or the RF generator prior to completely reading and understanding the applicable instructions for use.

26. Cardiac ablation procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory. Appropriate clinical instruction in the use of the THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter should also be completed.

27. The long-term risks of protracted fluoroscopy and creation of RF induced lesions have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children. Furthermore, the risk/benefit in asymptomatic patients has not been studied.

28. When using the THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter with conventional systems (using fluoroscopy to determine catheter tip location), or with the CARTO EP Navigation System, careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered. The firmness of the braided tip dictates that care must be taken to prevent perforation of the heart.

29. Always place the Rocker Lever of THERMOCOOL SF NAV Bi-directional Catheter in the neutral position to straighten the catheter tip before insertion or withdrawal of the catheter. Always pull the thumbknob of THERMOCOOL SF NAV Uni-directional Catheter back to straighten the catheter tip before insertion or withdrawal of the catheter.
30. Always maintain a constant heparinized normal saline infusion to prevent coagulation within the lumen of the catheter. Do not use this catheter without irrigation flow.
31. When RF current is interrupted for either a temperature or an impedance rise (the set limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum, if present. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode. Make sure the irrigation holes are not plugged prior to re-insertion.
32. Apparent low power output, high impedance reading, or failure of the equipment to function correctly at normal settings may indicate faulty application of the indifferent electrode(s) or failure of an electrical lead. Do not increase power before checking for obvious defects or misapplication of the indifferent electrode or other electrical leads.
33. Read and follow the indifferent electrode manufacturer's instructions for use; the use of indifferent electrodes that meet or exceed ANSI/AAMI requirements (AAMI IEC 60601-2-2), is recommended.
34. The THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter is intended for use with a compatible RF generator, a compatible irrigation pump, CARTO EP Navigation System and Biosense Webster cables, and other appropriate interface cables and connectors. Use of a compatible irrigation pump is recommended to assure proper irrigation flow rate.
35. The THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter has been shown to create larger lesions than standard non-irrigated RF ablation catheters. Care should be taken when ablating near structures such as the sino-atrial and atrioventricular nodes.
36. The sterile packaging and catheter should be inspected prior to use. Do not use if the packaging or catheter appears damaged.
37. The catheter is sterilized with ethylene oxide gas and should be used by the "Use By" date on the device package. Do not use the catheter if past the "Use By" date.
38. The THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter is intended for single patient use only.
39. Do not resterilize and reuse.
40. Do not use near MRI equipment since movement or heating of the catheter may occur and the image on the display may become distorted.
41. Use both fluoroscopy and electrogram data to monitor catheter advancement and reduce risk of tissue injury.
42. The THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter used in conjunction with a compatible RF generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the catheter and indifferent electrode, particularly when operating the catheter. During energy delivery, the patient should not be allowed to come in contact with grounded metal surfaces. If during ablation, temperature does not rise, discontinue delivery of energy and check set-up.
43. The risk of igniting flammable gases or other materials is inherent in electrosurgery. Precautions must be taken to restrict flammable materials from the electrosurgical suite.
44. Electromagnetic interference (EMI) produced by the THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter, when used in conjunction with a compatible RF generator during normal operation, may adversely affect the performance of other equipment.
45. Electrodes and probes for monitoring and stimulating devices can provide paths for high frequency current. The risk of burns can be reduced but not eliminated by placing the electrodes and probes as far away as possible from the ablation site and the indifferent electrode. Protective impedances may reduce the risk of burns, and permit continuous monitoring of the electrocardiogram during energy delivery.

46. The temperature sensor measures electrode tip temperature, not tissue temperature. The temperature displayed on the RF generator is for the cooled electrode only and does not represent tissue temperature. If the RF generator does not display temperature, verify that the appropriate cable is plugged into the RF generator. If temperature still is not displayed, there may be a malfunction in the temperature sensing system that must be corrected prior to applying RF power.
47. The temperature measurement accuracy of the THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter, as with any temperature measurement electrophysiology catheter, is largely determined by the temperature accuracy specification of the RF generator used. Please consult the user manual for the RF generator to be used for the temperature accuracy specification.
48. Before use, check irrigation ports are patent by infusing heparinized normal saline through the catheter and tubing.
49. Regularly inspect and test reusable cables and accessories.

RF Ablation

For RF ablation, the catheter must be connected to the PIU of the CARTO EP Navigation System. The PIU is then connected to the appropriate input connector(s) on the RF generator. For proper interface of the catheter to PIU and PIU to RF generator, use only Biosense Webster or compatible interface cables.

To complete the electrical circuit, an indifferent electrode must be connected to the indifferent electrode input on the RF generator. Verify that circuit impedance prior to RF ablation is within expected parameters. Verify that the RF generator displays a temperature not above 37°C after the catheter is inserted into the patient and before applying RF power.

RF Generator Operation

Refer to the applicable RF generator manual for proper operation of the RF generator.

RF ablation application parameters will vary depending on the ablation site, the specific conditions present in each procedure and the RF generator control circuitry. Based on data obtained from prior animal and clinical studies, recommended RF application parameters are provided below in the “Directions for Use” and in Table 1. Always monitor temperature and impedance rise when using the THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter.

How Supplied

- The THERMOCOOL SF NAV Uni-directional and Bi-directional Catheter is supplied sterile (EtO) unless the package is damaged or opened.
- Refer to the Device Description section of this Instructions for Use for available curve types for the various catheter configurations.
- Additional catheter accessory devices are provided separately.

Storage

Stored in a cool, dry, dark place. Storage temperature should be between 5 and 25°C (41 and 77°F).

Sterilization/“Use By” Date

This catheter has been sterilized with ethylene oxide gas. Product and package testing have been conducted to support the “Use By” date printed on the product labels. **DO NOT USE** after the “Use By” date.

This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or

lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

Disposal

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

Directions for Use

THERMOCOOL SF NAV Uni-directional

Please refer to the User Manuals for the CARTO EP Navigation System, irrigation pump, RF generator and the irrigation tubing for instructions on connecting and operating these systems in conjunction with the THERMOCOOL SF NAV Uni-directional Catheter. Use appropriate Biosense Webster accessory cables to connect the THERMOCOOL SF NAV Uni-directional Catheter to the appropriate accessory equipment.

1. Using aseptic technique, remove the catheter from the package and place in a sterile work area. Inspect the catheter carefully for electrode integrity and overall condition.
2. Create a vascular access in a large central vessel using aseptic techniques.
3. To verify compatibility between the sheath and catheter, advance the catheter through sheath prior to insertion.
4. Connect the catheter to the Patient Interface Unit (PIU) via the appropriate Biosense Webster cable. Connect the PIU to the RF generator via the appropriate Biosense Webster cable. Connect the PIU to the appropriate recording and mapping systems, including a compatible CARTO EP Navigation System, with appropriate interface cables. Use only Biosense Webster interface cables. To complete the electrical circuit, connect an indifferent electrode to the indifferent electrode input on the generator.
5. Connect the irrigation pump tubing to a room temperature, heparinized (1 IU heparin/ml) normal saline bag using standard safe hospital practices. Open the stopcock on the end of the irrigation tubing and fill the irrigation tubing as slowly as possible. Remove any trapped air and then close the stopcock.
6. Load the irrigation tubing into the irrigation pump. Open the stopcock and flush the irrigation tubing per the irrigation pump instructions until the air is expelled through the open end of the tubing.
7. Connect the stopcock on the end of the irrigation tubing to the Luer fitting of the THERMOCOOL SF NAV Uni-directional Catheter.
8. Flush the catheter and tubing per standard technique to ensure purging of trapped air bubbles and to verify that the irrigation holes are patent.
9. Start continuous irrigation at a flow rate of 2 ml/min.
10. Insert the THERMOCOOL SF NAV Uni-directional Catheter via the entrance site, using a guiding sheath with a minimum size of 8 F. Advance the catheter to the area under investigation. Use both fluoroscopy and electrograms to aid in proper positioning.
11. Deflection of the Catheters:
The catheter tip can be deflected to facilitate positioning by using the thumbknob to vary tip curvature. Tip deflection is controlled at the proximal end of the catheter by a handpiece in which a piston slides; a thumbknob on the piston controls piston travel. When the thumbknob is pushed forward, the tip is deflected (curved). When the thumbknob is pulled back, the tip straightens. The shape of the curve depends on the deflectable tip length.
12. For RF application and flow settings when using Stockert's RF generator, refer to Table 1.

Table 1: RF and Flow Settings when using Stockert's RF Generators

RECOMMENDED RF APPLICATION PARAMETERS		
PARAMETER	ATRIAL ABLATION	VENTRICULAR ABLATION
Power Range	15 W to 30 W*	31 W to 50 W
Temperature Monitoring	<40°C**	<40°C**
Irrigation Flow Rate During RF Application	8 ml/min	15 ml/min
Application Time	30 to 120 seconds	60 to 120 seconds
* Power levels exceeding 30 Watts may be used when transmural lesions cannot be achieved at lower energy levels. For power settings > 30 Watts, the recommended irrigation flow rate is 15 ml/min. ** The temperature displayed on the RF generator does not represent tissue temperature or electrode tissue interface temperature. Additional Recommendation: For isthmus dependent flutter ablation, power applications exceeding 30 watts and up to 50 watts should only be used if conduction block cannot be achieved at lower power levels.		

13. Recommendation for irrigation when using Stockert's RF generator:

Increase the irrigation to high flow rate starting up to 5 seconds before the onset of RF energy delivery and maintaining this higher flow rate up to 5 seconds after termination of the energy application. For power levels up to 30 W, a high flow rate of 8 ml/min should be used. For power levels between 31-50 W, a high flow rate of 15 ml/min should be used. Do not use the catheter without irrigation flow.

14. Monitor the catheter tip temperature throughout the procedure to ensure adequate irrigation. If the temperature increases to 40°C during RF energy delivery, power delivery should be interrupted. The irrigation system must be rechecked prior to restarting RF delivery. Note: the displayed temperature represents the temperature of the electrode only, not the temperature of the tissue.

15. Start a procedure at 15-20 W. After 15 seconds, power may be increased by 5-10 W increments as needed, until a transmural lesion is achieved, defined by > 80% reduction in unipolar atrial electrogram amplitude, or emergence of double potentials of equal and low amplitude. It is recommended that power not exceed 50 W when the catheter is parallel to the tissue and 35 W if the catheter is perpendicular to the tissue. The duration of each RF ablation should not exceed 120 seconds. Dragging the catheter to the next location is permissible during the 120-second energy application. RF current may be reapplied to the same or alternate sites using the same catheter.

16. In the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip dome electrode inspected for coagulum before RF current is reapplied. To remove any coagulum, if present, a sterile gauze pad dampened with sterile saline may be used to gently wipe the tip section clean; do not scrub or twist the tip dome electrode as damage to the tip dome

electrode bond may occur and loosen the tip dome electrode. Prior to reinsertion, ensure that the irrigation holes are not plugged as follows:

- a. Fill a 1 or 2 ml syringe* with sterile saline and attach to the stopcock on the end of the irrigation tubing.
- b. Carefully inject the saline from the syringe into the catheter. Uniform streams of fluid should be visible from the tip of the catheter.
- c. Repeat steps a and b, if necessary.
- d. Flush catheter and tubing per standard technique to ensure purging of trapped air bubbles and to verify that the irrigation holes are patent.
- e. The catheter can now be introduced into the patient.

WARNING: Do not continue use of the catheter if still occluded or if it is not functioning properly.

***NOTE:** A small syringe provides sufficient pressure to produce a visible stream of fluid.

THERMOCOOL SF NAV Bi-directional Catheter

Please refer to the User Manuals for the CARTO EP Navigation System, irrigation pump, RF generator and the irrigation tubing for instructions on connecting and operating these systems in conjunction with the THERMOCOOL SF NAV Bi-directional Catheter. Use appropriate Biosense Webster accessory cables to connect the THERMOCOOL SF NAV Bi-directional Catheter to the appropriate accessory equipment.

1. Using aseptic technique, remove the catheter from the package and place in a sterile work area. Inspect the catheter carefully for electrode integrity and overall condition.
2. Create a vascular access in a large central vessel using aseptic techniques.
3. To verify compatibility between the sheath and catheter, advance the catheter through sheath prior to insertion.
4. Connect the catheter to the Patient Interface Unit (PIU) via the appropriate Biosense Webster cable. Connect the PIU to the RF generator via the appropriate Biosense Webster cable. Connect the PIU to the appropriate recording and mapping systems, including a compatible CARTO EP Navigation System, with appropriate interface cables. Use only Biosense Webster interface cables. To complete the electrical circuit, connect an indifferent electrode to the indifferent electrode input on the generator.
5. Connect the irrigation pump tubing to a room temperature, heparinized (1 IU heparin/ml) normal saline bag using standard safe hospital practices. Open the stopcock on the end of the irrigation tubing and fill the irrigation tubing as slowly as possible. Remove any trapped air and then close the stopcock.
6. Load the irrigation tubing into the irrigation pump. Open the stopcock and flush the irrigation tubing per the irrigation pump instructions until the air is expelled through the open end of the tubing.
7. Connect the stopcock on the end of the irrigation tubing to the luer fitting of the THERMOCOOL SF NAV Bi-directional Catheter.
8. Flush the catheter and tubing per standard technique to ensure purging of trapped air bubbles and to verify that the irrigation holes are patent.
9. Start continuous irrigation at a flow rate of 2 ml/min.
10. Insert the THERMOCOOL SF NAV Bi-directional Catheter via the entrance site, using a guiding sheath with a minimum size of 8 F. Advance the catheter to the area under investigation. Use both fluoroscopy and electrograms to aid in proper positioning.
11. Deflection of the Catheters:

Use the Rocker Lever to deflect the catheter tip (Figure 1). When the lever is pulled back from neutral, the tip will deflect relative to the direction of rotation. The amount of deflection is relative to the amount of lever rotation.

When the lever is pushed forward, the tip will deflect in the opposite direction. To straighten the tip, return the Rocker Lever to neutral position.



Figure 1

The handle has an adjustable friction control that allows the operator to use the Rocker Lever and deflecting tip in a “free” state or adjust the friction to where the Rocker Lever and tip curve are “locked” in place (Figure 2). This knob is located on the opposite side of the Rocker Lever. Out of the package, the knob will be in the “off” position, which allows the freest movement for the lever and deflecting tip. The amount of friction increases as the Friction Control Knob is rotated clockwise until it reaches the fully “on” position. Clockwise rotation from the “off” position increases the friction within the Deflection Mechanism. Counter-Clockwise rotation from the “on” position decreases the friction within the Deflection Mechanism.

Clockwise rotation from the “off” position increases the friction within the Deflection Mechanism.

Counter-Clockwise rotation from the “on” position decreases the friction within the Deflection Mechanism.

Figure 2

12. For RF application and flow settings when using Stockert’s RF generator, refer to Table 1.

13. Recommendation for irrigation when using Stockert’s RF generator:

Increase the irrigation to high flow rate starting up to 5 seconds before the onset of RF energy delivery and maintaining this higher flow rate up to 5 seconds after termination of the energy

application. For power levels up to 30 watts, a high flow rate of 8 ml/min should be used. For power levels between 31-50 watts, a high flow rate of 15 ml/min should be used.

14. The delivery of RF energy must not be initiated until the increase in irrigation flow rate is confirmed by a minimum of 2°C decrease in tip dome electrode temperature.

15. Monitor the catheter tip temperature throughout the procedure to ensure adequate irrigation. Note: the displayed temperature represents the temperature of the electrode only, not the temperature of the tissue.

16. Start a procedure at 15-20 Watts. After 15 seconds, power may be increased by 5-10 W increments as needed, until a transmural lesion is achieved, defined by > 80% reduction in unipolar atrial electrogram amplitude, or emergence of double potentials of equal and low amplitude. It is recommended that power not exceed 50 W when the catheter is parallel to the tissue and 35 W if the catheter is perpendicular to the tissue. The duration of each RF ablation should not exceed 120 seconds. Dragging the catheter to the next location is permissible during the 120-second energy application. RF current may be reapplied to the same or alternate sites using the same catheter.

17. In the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip dome electrode inspected for coagulum before RF current is reapplied. To remove any coagulum, if present, a sterile gauze pad dampened with sterile saline may be used to gently wipe the tip section clean; do not scrub or twist the tip dome electrode as damage to the tip dome electrode bond may occur and loosen the tip dome electrode. Prior to reinsertion, ensure that the irrigation holes are not plugged as follows:

- a. Fill a 1 or 2 ml syringe* with sterile saline and attach to the stopcock on the end of the irrigation tubing.
- b. Carefully inject the saline from the syringe into the catheter. Uniform streams of fluid should be visible from the tip of the catheter.
- c. Repeat steps a and b, if necessary.
- d. Flush catheter and tubing per standard technique to ensure purging of trapped air bubbles and to verify that the irrigation holes are patent.
- e. The catheter can now be introduced into the patient.

WARNING: Do not continue use of the catheter if still occluded or if it is not functioning properly.

***NOTE:** A small syringe provides sufficient pressure to produce a visible stream of fluid.

Adverse Reactions

A number of serious adverse reactions have been documented for catheter ablation procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death.

The following complications were also noted to have occurred during prior studies or have been reported in the literature:

1. Catheterization/catheter procedure related: vascular bleeding/local hematomas, thrombosis, arteriovenous (AV) fistula, pseudoaneurysm, thromboembolism and vasovagal reactions, cardiac perforation, pericardial effusion/tamponade, thrombi, air embolism, arrhythmias and valvular damage, pneumothorax and hemothorax, pulmonary edema, hypoxia, pleural effusion, acute respiratory distress syndrome (ARDS), congestive heart failure, aspiration pneumonia, pneumonia, asthmatic attack, hypotension, implantable cardioverter-defibrillator (ICD) lead malfunction, anemia, thrombocytopenia, disseminated intravascular coagulation, epistaxis, systemic infection, urinary tract infection, sedation induced apnea, sedation induced CO₂ retention with lethargy and cholecystitis.
2. RF related: chest pain/discomfort, ventricular tachyarrhythmia, transient ischemic attack (TIA), cerebrovascular accident (CVA), complete heart block, coronary artery spasm, coronary

artery thrombosis, coronary artery dissection, cardiac thromboembolism, pericarditis, cardiac perforation/tamponade, valvular damage and increased phosphokinase level.

3. Unrelated to device or procedure: urinary retention, transient extremity numbness, Parkinson’s disease and gastrointestinal diverticulosis.

DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW.

WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product at the time of manufacture and are not made or given as a warranty of the prescribed product in any way.

Appendix 1 (Addition to Instructions for Use)

1.1 Operation conditions

Intended for use in the internal environment of man.

This medical device must be used only by qualified medical specialists or medical personal who have been training sessions on working with similar products, under the guidance of an experienced specialist.

Table 1 — Environmental specifications

Environment Settings	Values
Temperature, °C	(+32 - +42) ±2 °C
Air humidity (relative), non-condensing, %	(80-100%) ±2 %
Atmosphere pressure, hPa	101 hPa

2 Information of manufacturers of the medical device (name, location address, mailing address, phone number, address of place of production)

DEVELOPER, MANUFACTURER:

Name: Biosense Webster, Inc.

Address: 33 Technology Drive Irvine, California 92618 USA

MANUFACTURER

1. Biosense Webster, Inc.

33 Technology Drive Irvine, California 92618 USA

2. Biosense Webster, Inc.

Circuito Interior Norte #1820, Parque Industrial Salvarcar, Juarez, Chihuahua, Mexico, 32574

3. Biosense Webster, Inc.

15715 Arrow Highway, Irwindale, California, 91706 USA

AUTHORIZED REPRESENTATIVE OF MANUFACTURER

LLC «Johnson & Johnson»

121614, Moscow, st. Krylatskaya 17/2, Russian Federation

3 Classification of medical device

Table 2 — Classification

Classification	Class
Class risk	3
Invasiveness	Surgically invasive MD
Sterility	Sterile MD
Type of contact with the body	Circulating Blood
Duration of contact	Limited duration of contact – less than 24 hours

4 Technical specifications of the medical device

4.1 Photo of MD

Figure 1 — Photo of THERMOCOOL SF NAV UNI-DIRECTIONAL CATHETER


Figure 2 — Photo of THERMOCOOL SF NAV BI-DIRECTIONAL CATHETER

4.2 General description

The THERMOCOOL SF NAV Bi-directional Catheter is a steerable multi electrode (6), luminal catheter with a deflectable tip via moveable thumb knob, designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) current to the catheter tip dome electrode for ablation purposes. The catheter shaft measures 7.5 F with 7.5 F ring electrodes. The catheter is introduced through a guiding sheath with a minimum size of 8 F. For ablation, the catheter is used in conjunction with an RF generator and a dispersive pad (indifferent electrode).

The catheters have a high-torque shaft with a bi-directional deflectable tip section containing an array of electrodes, which includes a 3.5mm tip dome.

All ring electrodes are manufactured from noble metals and, unless otherwise noted, may be used for recording and stimulation purposes. The tip dome electrode serves to deliver RF current from the RF generator to the desired ablation site.

The THERMOCOOL SF NAV Bi-directional Catheter incorporate a temperature sensor, which is embedded in the 3.5mm tip electrode. A Rocker Lever is used to deflect the tip.

The Rocker Lever pulls one of the two puller wires to deflect the catheter in one direction then can alternatively select the second (opposite) puller wire to deflect the catheter in a second direction. Selective deflection of the catheter is accomplished by rotating the control knob in the direction of the desired deflection.

The handle also has an adjustable Friction Control Knob that allows the operator to use the Rocker Lever and deflecting tip in a “free state” or adjust the friction to where the Rocker Lever and tip curve are “locked” in place.

The amount of friction increases as the Friction Control Knob is rotated clockwise until it reaches the fully “locked” (“on”) position. The Friction Control Knob is located on the opposite side of the Rocker Lever. The high-torque shaft also allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site.

Additionally, a variety of curve types are available in symmetric or asymmetric combinations, providing two 180° opposed, single planed curves.

Currently, the available curves for the THERMOCOOL SF NAV Bi-directional Catheter

include BB, DD, FF, JJ, BD, FJ, DJ, BF, DF. Refer to Figure 3 for some of the available curve types for the THERMOCOOL SF NAV Bi-directional Catheter.


Figure 3 —

At the proximal end of the catheter, a saline input port with a standard luer fitting terminates from the open lumen. This saline port serves to permit the injection of normal saline to irrigate the tip electrode. During ablation, heparinized normal saline is passed through the internal lumen of the catheter and through the tip electrode, to irrigate and cool the ablation site as well as the electrode tip. An irrigation pump is used to control the saline irrigation. The catheter interfaces with standard recording equipment and a compatible RF generator via accessory extension cables with the appropriate connectors.

The THERMOCOOL SF NAV Bi-directional Catheter feature a magnetic location sensor embedded in the tip section which transmits location information to the Carto EP Navigation System. An appropriate reference device is required for location reference position purposes.


Figure 4 — THERMOCOOL SF NAV UNI-DIRECTIONAL CATHETER


Figure 5 — THERMOCOOL SF NAV BI-DIRECTIONAL CATHETER

4.3 The handle of the catheter is laser etched for easy identification of catheter and direction of curve deflection. Versions / Options of MD

1. ThermoCool SF NAV Uni-Directional Catheter:
 - Curve type B; electrodes configuration 2-5-2 mm;
 - Curve type D; electrodes configuration 2-5-2 mm;
 - Curve type F; electrodes configuration 2-5-2 mm;
 - Curve type J; electrodes configuration 2-5-2 mm.

2. ThermoCool SF NAV Bi-Directional Catheter:
 - Curve type B-B; electrodes configuration 2-5-2 mm;
 - Curve type B-D; electrodes configuration 2-5-2 mm;
 - Curve type B-F; electrodes configuration 2-5-2 mm;
 - Curve type D-D; electrodes configuration 2-5-2 mm;
 - Curve type D-F; electrodes configuration 2-5-2 mm;
 - Curve type D-J; electrodes configuration 2-5-2 mm;
 - Curve type F-F; electrodes configuration 2-5-2 mm;
 - Curve type F-J; electrodes configuration 2-5-2 mm;
 - Curve type J-J; electrodes configuration 2-5-2 mm.
3. ThermoCool SF NAV Uni-Directional Catheter with curve visualization:
 - Curve type B; electrodes configuration 2-5-2 mm;
 - Curve type D; electrodes configuration 2-5-2 mm;
 - Curve type F; electrodes configuration 2-5-2 mm;
 - Curve type J; electrodes configuration 2-5-2 mm.
4. ThermoCool SF NAV Bi-Directional Catheter with curve visualization:
 - Curve type B-D; electrodes configuration 2-5-2 mm;
 - Curve type B-F; electrodes configuration 2-5-2 mm;
 - Curve type J-J; electrodes configuration 2-5-2 mm.

Table 3 — Options of MD

Catheter	Product Code	Configuration	Manufacturing Part Number
ThermoCool SF NAV Uni-Directional Catheter	D131801	Curve type B; electrodes configuration 2-5-2 mm	D-1318-XX-S
	D131802	Curve type D; electrodes configuration 2-5-2 mm	
	D131803	Curve type F; electrodes configuration 2-5-2 mm	
	D131804	Curve type J; electrodes configuration 2-5-2 mm	
ThermoCool SF NAV Bi-Directional Catheter	BNI35BBH	Curve type B-B; electrodes configuration 2-5-2 mm	D-1317-XX-S
	BNI35BDH	Curve type B-D; electrodes	

		configuration 2-5-2 mm	
	BNI35BFH	Curve type B-F; electrodes configuration 2-5-2 mm	
	BNI35DDH	Curve type D-D; electrodes configuration 2-5-2 mm	
	BNI35DFH	Curve type D-F; electrodes configuration 2-5-2 mm	
	BNI35DJH	Curve type D-J; electrodes configuration 2-5-2 mm	
	BNI35FFH	Curve type F-F; electrodes configuration 2-5-2 mm	
	BNI35FJH	Curve type F-J; electrodes configuration 2-5-2 mm	
	BNI35JJH	Curve type J-J; electrodes configuration 2-5-2 mm	
ThermoCool SF NAV Uni-Directional Catheter with curve visualization	D131501	Curve type B; electrodes configuration 2-5-2 mm	D-1315-XX-S
	D131502	Curve type D; electrodes configuration 2-5-2 mm	
	D131503	Curve type F; electrodes configuration 2-5-2 mm	
	D131504	Curve type J; electrodes	

		configuration 2-5-2 mm	
ThermoCool SF NAV Bi-Directional Catheter with curve visualization	BNI35BDCT	Curve type B-D; electrodes configuration 2-5-2 mm	D-1313-XX-S
	BNI35BFCT	Curve type B-F; electrodes configuration 2-5-2 mm	
	BNI35JJCT	Curve type J-J; electrodes configuration 2-5-2 mm	

4.4 The possibility and ways of integrating with other MD

The following devices are required to conduct an Electrophysiology procedure with the ThermoCool SF Catheters:

- Catheter interface cable
- External reference device
- Compatible RF generator
- Irrigation pump and tubing set
- Electrophysiology (EP) recording equipment (not marketed by Biosense Webster)
- Grounding pad (not marketed by Biosense Webster)
- Pacing stimulator (not marketed by Biosense Webster)
- Carto XP or Carto 3 EP Navigation System (navigational catheters only)

The catheter interface cable, irrigation pump, external reference device, RFGenerator (currently only Stockert RF Generators) and the Carto XP / Carto 3 EP Navigation systems are CE marked devices marketed by Biosense Webster and they have technical files of their own. The grounding pad, pacing stimulator and EP recording equipment are CE marked devices marketed by other companies. There are no accessories included with the ThermoCool SF Catheters.

4.5 The materials from which it is made or its medical device basic parts, biocompatibility

Table 4 — The materials from which it is made or its medical device basic parts

Description	Material
Materials of following catheters: THERMOCOOL SF NAV UNI-DIRECTIONAL AND BI-DIRECTIONAL CATHETERS	
Tip/ Shaft Adhesive	Polyurethane (PU) P-9231
Luer Hub	Polycarbonate Plastic Polycarbonate Resin Blue or White - HP3NREU-8H9D273
Irrigation Pathway	Polycarbonate (Plastic Polycarbonate Resin Blue or White), HP3NREU-8H9D273

		Polyimide 75 D, Palladium/Platinum Victrex 450G PEEK Black, Polyurethane adhesive P-9231	
		Tip	
Outer Coat Tip		Polyimide 75 D, Teflon Tubing 125ID	
Core		Pallethane 55D/80A	
Tip Electrode Insert		Palladium 80%/ Platinum 20% Victrex 450G PEEK Black	
Ring Electrode		90% Platinum / 10% Iridium ASTM: B684-97 P/TR 8007 contain P/TR 8007-01	
Soft tip		Pellethane 2363	
Tip of Quad Lumen Junction		Pellethane 80% 55D 20% 80A	
Braiding		Stainless Steel 304v	
Temperature Sensor (Thermocouple)		Polyimide Sleeve 75D, Epoxy Encapsulant Constantan/ Cooper (THERMAL CLASSIFICATION 155C)	
Connector/Leadwire Joints		Solder SN63/PB37	
		Shaft	
Shaft		Pebax 2553-SA-01 Vestamid VESTAMID, CMPND W/20% BS-80% 75D	
Braiding		Stainless Steel 304v	
Insulation		SPN (SINGLE POLY NYLON) (THERMAL CLASSIFICATION 155C)	
Puller Wire		Stainless steel 304v	
Compression Coil		Stainless Steel Flat Wire 304v	
Compression Coil Protection		Polyimide Tube 75D	
Shaft Stiffener		Polyimide 75D	
Tip Conductor/Lead Wire Joints		Solder SN96/AG4-Weld	
		Handle	
THERMOCOOL SF NAV UNI DIRECTIONAL		THERMOCOOL SF NAV Bi-Directional Catheter	
Description	Material	Description	Material
Deflection mechanism	Valox 325	Deflection mechanism	Nylon Zytel® 42A NC010 w/ 30% glass and 15% PTFE filler Nylon

Outer Barrel	Delrin/ Acetal ASTM-D6778 POM 0211	Outer Barrel	Polycarbonate LEXAN™ Resin HP3NREU w/30% Glass and 5% PTFE filler
--------------	--	--------------	---

Animal Materials

Pellethane material used in both the Bi-directional and Uni-directional ThermoCool SF Catheter families contains animal tallow. Tallow material used in the device is compliant with EMEA/410/01 Rev. 2.

Tallow material contains materials derived from animal sources. These materials are primarily derived from bovine (beef) sources, although porcine (hog) fat can also be used.

Biocompatibility

Testing was conducted in accordance with the testing derived from “Table 1 – Guidance for Initial Evaluation Tests” of the ISO 10993-1 Standard. Per the ISO Standard, the catheters were classified as Externally Communicating Device; Limited Contact Duration (<24 hours); Circulating Blood. The testing was conducted on the distal tip shaft / tip section of the Celsius DS catheter (D-1194-22) which is incorporated with a 80% Palladium / 20% Platinum tip electrode. The testing conducted on the Celsius DS catheter (D- 1194-22) is applicable to all BWI catheters incorporating an 80% Palladium /20% Platinum tip electrode. The following testing was conducted following GLP guidelines:

- Hemolysis Test ASTM (F 756) Method – Direct Contact Method
- MEM Elut ion Using L-929 Mouse Fibroblast Cells (ISO) (Cytotoxicity)
- Intracutaneous Irritation Test (ISO) – 2 extract
- Acute Systemic Injection Test (ISO) – 2 extract
- Materials Mediated Rabbit Pyrogen
- Guinea Pig Maximization Sensitization Test – 2 extract
- Complement Activation C3a Assay and SC5b-9 Assay
- Partial Thromboplastin Time (PTT)
- Platelet and Leukocyte Counts
- Thrombosis (in vivo) – 2 dogs

Results of the biocompatibility testing performed are contained in Table 4.

Table 5 — Biocompatibility Testing Results

Test	Result
Hemolysis Test (ASTM (F756) Method – Direct Contact	Pass
MEM Elution Using L-929 Mouse Fibroblast Cells (ISO) (Cytotoxicity)	Pass
Intracutaneous Irritation Test (ISO)	Pass
Acute Systemic Injection Test (ISO)	Pass
Materials Mediated Rabbit Pyrogen	Pass
Guinea Pig Maximization Sensitization	Pass
Complement Activation C3a Assay and SC5b-9 Assay	Pass
Partial Thromboplastin Time (PTT)	Pass
Platelet and Leukocyte Counts	Pass
Thrombosis (In vivo)	Pass

The results of the testing demonstrate the biocompatibility of the 80% Palladium / 20% Platinum tip electrode material.

5 Packaging

The Uni-directional ThermoCool SF Nav Catheters packaged in a two-piece thermoform tray with a handle and shaft component. The two-piece thermoform tray is snapped together to form a catheter support tray that is not part of the sterile barrier. The catheter is loaded into the tray. The loaded tray is placed in a Tyvek/ nylon pouch, heat-sealed forming a single sterile barrier. The sealed pouch is then placed into a unit carton with an IFU insert.

Figure 6 — Primary packaging of The Uni-directional ThermoCool SF Nav Catheter

Figure 7 — Secondary packaging of The Uni-directional ThermoCool SF Nav Catheter

The ThermoCool SF NAV Bi-Directional Catheter will be packaged in a two-piece thermoform tray with a handle and shaft component. The two-piece thermoform tray is snapped together to form a catheter support tray that is not part of the sterile barrier. The catheter is loaded into the tray. The loaded tray is placed in a Tyvek/ nylon pouch, heat-sealed forming a single sterile barrier. The sealed pouch is then placed into a unit carton with an IFU insert.

Figure 8 — Primary packaging of The Bi-directional ThermoCool SF Nav Catheter

Figure 9 — Secondary packaging of The Bi-directional ThermoCool SF Nav Catheter

Table 6 — Packaging components for ThermoCool SF Nav Catheters

The Uni-directional ThermoCool SF Nav Catheters	
Description	Material
Shaft Protector, Straight Catheter	Dural 714 FD PVC,Natural
Mounting Card, NaviStar ThermoCool	Solid Bleached Sulfate,White 609.6 (0.24")
Pouch	Nylon, 0,15 mm (6 mil)
Carton Unit	Carton 109,54mm x 31,75mm x 1524mm (4 5/16" x 1 1/4" x 60")
Patent insert	Opaque paper

IFU	1.1.1.1	Printed with black ink on white paper (40# Opaque Smooth).
	1.1.1.2	Booklet form, printed double-sided (on both front and back sides), and having stitched binding.
	1.1.1.3	Size and format of finished IFU booklet will be 95.25mm x 323.85mm (3.75" x 12.75") and will be portrait orientation.
Transport packaging (materials, dimensions, qty.)	Paper, Carton	
The ThermoCool SF NAV Bi-Directional Catheter		
Chevron Pouch	Tyvek/Nylon(NYLON POUCH, 0,15 mm (6 mil))	
Directional Handle Tray	PETG(350 mmx250 mmx150 mm)	
Catheter Shaft Tray	PETG(BOX, 304.8mm X 304.8mm X 22.225mm (12"X12"X7/8"),SQUARE)	
Unit Carton	Corrugated Paper Board (CARTON UNIT 4.20mm X 1.40mm X59.75mm)	
Transport packaging (materials, dimensions, qty.)	Paper, Carton HEADER BAG, 381mm X 330.2mm (15" X 13")	

6 Information of labeling


Figure 10 — Label on primary packaging


Figure 11 — Label on secondary packaging

Table 7 — List of symbols

Symbol	Definition	Symbol	Definition
	Sterilized Using Ethylene Oxide		Keep Dry
	Do Not Re-Use		Use-By Date
	Caution		Batch Code
	Consult Instructions for Use		Catalog Number
	Do not use if package is damaged		Contents: 1
	Do not use if package is opened.		Manufacturer
	Keep Away from Sunlight		Authorized Representative in the European Community
	Date of Manufacture		Temperature Limit
	Electrodes		Spacing
	Curve Type.		Pin Connector
	Navigational Catheter		Compatible with Carto 3 EP Navigation System
	Compatible with Carto XP EP Navigation System		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	THERMOCOOL SF Catheter		

7 Conditions of transportation and storage

Table 8 — Storage

Environment Settings	Values
Air temperature, °C	(5-25) ± 2 °C
Air humidity (relative), non-condensing, %	(20-60%) $\pm 2\%$
Atmosphere pressure, hPa	101 ± 1 hPa

Table 9 — Transportation

Environment Settings	Values
Air temperature, °C	(+5 - +25) ± 2 °C
Air humidity (relative), non-condensing, %	(20-60%) $\pm 2\%$
Atmosphere pressure, hPa	101 ± 1 hPa

8 Main parameters and characteristics of the medical device

Table 10 — Main parameters

1. ThermoCool SF NAV Uni-Directional Catheter				
Characteristics	D131801	D131802	D131803	D131804
General				
Usable Length, cm	115 ± 3	115 ± 3	115 ± 3	115 ± 3
Puller Wire material Length mm	1524 $\pm 76,2$	1524 $\pm 76,2$	1524 $\pm 76,2$	1524 $\pm 76,2$
Outer diameter, mm	0 to 2,6416	0 to 2,6416	0 to 2,6416	0 to 2,6416
Application Time, sec	not exceed 120	not exceed 120	not exceed 120	not exceed 120

Power Settings During RF Application (atrial Ablation), W	15 to 30	15 to 30	15 to 30	15 to 30
Power Settings During RF Application (ventricular Ablation), W	31 to 50	31 to 50	31 to 50	31 to 50
Flow Rate During ablation, ml/min	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)
	15 ml/min-5% +15% (> 30 W)	15 ml/min-5% +15% (> 30 W)	15 ml/min-5% +15% (> 30 W)	15 ml/min-5% +15% (> 30 W)
Pre-curve °	0 to 15°	0 to 15°	0 to 15°	0 to 15°
Full length, cm*	146 ± 3	146 ± 3	146 ± 3	146 ± 3
Mass, g*	90±10%	90±10%	90±10%	90±10%
Luer fitting				
Inner diameter, mm	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286
Outer diameter, mm	7.747 +0.000 - 0.0762	7.747 +0.000 - 0.0762	7.747 +0.000 - 0.0762	7.747 +0.000 -0.0762
Length, mm	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254
Side arm Cord length,	215,9 +/- 3,175	215,9 +/- 3,175	215,9 +/- 3,175	215,9 +/- 3,175

mm				
Tip Specifications				
Tip electrode diameter, mm (Fr)	2,3 (7,5)	2,3 (7,5)	2,3 (7,5)	2,3 (7,5)
Number of Irrigation Holes	56	56	56	56
Nominal irrigation hole diameter, mm	0,0889 +/- 0,508	0,0889 +/- 0,508	0,0889 +/- 0,508	0,0889 +/- 0,508
Number of Lumens in the Deflectable Tip	4	4	4	4
Deflectable Tip Lumen Dimensions, mm	2 x 0,4826 +/- 0.0254 2 x 0,762 +/- 0.0254	2 x 0,4826 +/- 0.0254 2 x 0,762 +/- 0.0254	2 x 0,4826 +/- 0.0254 2 x 0,762 +/- 0.0254	2 x 0,4826 +/- 0.0254 2 x 0,762 +/- 0.0254
Braiding	Braided: 16 strands, 45 +/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45 +/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45 +/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45 +/-1 picks per 2,54 cm (1")
Tip Electrode Length, mm	3,6703 +/- 0,5083,5	3,6703 +/- 0,5083,5	3,6703 +/- 0,5083,5	3,6703 +/- 0,5083,5
Tip Electrode Construction	Two piece with stem	Two piece with stem	Two piece with stem	Two piece with stem
Number of Ring Electrodes	3	3	3	3

Additional Rings, mm	2	2	2	2
Ring Electrode Dimensions, mm	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length= 1,27 +/- 0.127	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length= 1,27 +/- 0.127	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length= 1,27 +/- 0.127	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length= 1,27 +/- 0.127
	Shaft Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3622 +/- 0.0254 Length= 0,7112 +/- 0.127	Shaft Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3622 +/- 0.0254 Length= 0,7112 +/- 0.127	Shaft Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3622 +/- 0.0254 Length= 0,7112 +/- 0.127	Shaft Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3622 +/- 0.0254 Length= 0,7112 +/- 0.127
Ring Electrode Spacing	2-5-2	2-5-2	2-5-2	2-5-2
RF Impedance ohms when 5 Killohertz	0 to less than 5	0 to less than 5	0 to less than 5	0 to less than 5
RF Impedance ohms when 500 Killohertz	0 to less than 25	0 to less than 25	0 to less than 25	0 to less than 25
Isolation and Phase Angle Between Circuits	At 5 kHz: : greater than 100k ohms to infinity. At 500 kHz: : greater than 1k ohms to infinity.. The phase angle must be between -86° to -90°.	At 5 kHz: : greater than 100k ohms to infinity. At 500 kHz: : greater than 1k ohms to infinity.. The phase angle must be between -86°	At 5 kHz: : greater than 100k ohms to infinity. At 500 kHz: : greater than 1k ohms to infinity.. The phase angle must be	At 5 kHz: : greater than 100k ohms to infinity. At 500 kHz: : greater than 1k ohms to infinity.. The phase angle must be between -86° to -90°.

		to -90°.	between -86° to -90°.	
RF Leakage Current, mA	0 to 252	0 to 252	0 to 252	0 to 252
Temperature Sensors	Thermo-couple: 14004 Ohms 37°+/- 2°C	Thermo-couple: 14004 Ohms 37° +/- 2°C	Thermo-couple: 14004 Ohms 37°+/- 2°C	Thermo-couple: 14004 Ohms 37° +/- 2°C
Compression Coil Location	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.
Curve Types	B	D	F	J
Deflection Angle Range	0-150°	0-150°	0-150°	0-150°
Tip Stiffness, N	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647
Buckle Force, N	0 to less than 0,5589	0 to less than 0,5589	0 to less than 0,5589	0 to less than 0,5589
Side Force, N	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity
Pull Entire Catheter, N	Greater than 17,793 to infinity	Greater than 17,793 to infinity	Greater than 17,793 to infinity	Greater than 17,793 to infinity
Irrigation Pathway Specifications				
Irrigation Tubing Dimensions, mm	ID = 0,6358+/- 0.0254	ID = 0,6358+/- 0.0254	ID = 0,6358+/- 0.0254	ID = 0,6358+/- 0.0254
Shaft Specifications				
Shaft Size, mm (Fr)	2,3 (7,5 r)	2,3 (7,5)	2,3 (7,5)	2,3 (7,5)

Shaft Waviness, mm	0 to 20	0 to 20	0 to 20	0 to 20
Internal Lead Conductors	38 AWG Copper Wire	38 AWG Copper Wire	38 AWG Copper Wire	38 AWG Copper Wire
Temperature Sensor Conductor	Thermo-couple: • 38 AWG Constantan Wire • 38 AWG Copper Wire	Thermo-couple: • 38 AWG Constantan Wire • 38 AWG Copper Wire	Thermo-couple: • 38 AWG Constantan Wire • 38 AWG Copper Wire	Thermo-couple: • 38 AWG Constantan Wire • 38 AWG Copper Wire
location sensor conductor	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper
Puller Wire material length, mm	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2
Compression Coil	0,1016 x 0,1524	0,1016 x 0,1524	0,1016 x 0,1524	0,1016 x 0,1524
Handle Specifications				
Construction	Handle with adjustable deflection knob, tension knob and luer fitting	Handle with adjustable deflection knob, tension knob and luer fitting	Handle with adjustable deflection knob, tension knob and luer fitting	Handle with adjustable deflection knob, tension knob and luer fitting
Length, cm	21.03+/- 0.05	21.03+/- 0.05	21.03+/- 0.05	21.03+/- 0.05
Length of Piston, cm	6.985 +/- 0.0762	6.985 +/- 0.0762	6.985 +/- 0.0762	6.985 +/- 0.0762
Connector	25 pin	25 pin	25 pin	25 pin
Handle Temperature, °C	0 to 48	0 to 48	0 to 48	0 to 48

2. ThermoCool SF NAV Bi-Directional Catheter									
Charact eristics	BNI35 BBH	BNI35 BDH	BNI35 BFH	BNI35 DDH	BNI35D FH	BNI35 DJH	BNI35FF H	BNI35 FJH	BNI35 JJH
General									
Usable Length, cm	115 ± 3	115 ± 3	115 ± 3	115 ± 3	115 ± 3	115 ± 3	115 ± 3	115 ± 3	115 ± 3
Puller Wire material Length, mm	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2
Outer diameter, mm	0 to 2.6416	0 to 2.6416	0 to 2.6416	0 to 2.6416	0 to 2.6416	0 to 2.6416	0 to 2.6416	0 to 2.6416	0 to 2.6416
Appli cation Time, sec	not exceed 120	not exceed 120	not exceed 120	not exceed 120	not exceed 120	not exceed 120	not exceed 120	not exceed 120	not exceed 120
Power Settings Durin g RF Appli cation (Atria l Ablati on) W	15 to 30	15 to 30	15 to 30	15 to 30	15 to 30	15 to 30	15 to 30	15 to 30	15 to 30
Power Settings Durin g RF Appli cation (Vent ricula r Ablati on) W	31 to 50	31 to 50	31 to 50	31 to 50	31 to 50	31 to 50	31 to 50	31 to 50	31 to 50

Flow Rate, ml/min	8 ml/min 5% +15% (≤ 30 W)	8 ml/min 5% +15% (≤ 30 W)	8 ml/min 5% +15% (≤ 30 W)	8 ml/min 5% +15% (≤ 30 W)	8 ml/min 5% +15% (≤ 30 W)	8 ml/min 5% +15% (≤ 30 W)	8 ml/min 5% +15% (≤ 30 W)	8 ml/min 5% +15% (≤ 30 W)	8 ml/min 5% +15% (≤ 30 W)
	15 ml/min 5% +15% (> 30 W)	15 ml/min 5% +15% (> 30 W)	15 ml/min 5% +15% (> 30 W)	15 ml/min 5% +15% (> 30 W)	15 ml/min 5% +15% (> 30 W)	15 ml/min 5% +15% (> 30 W)	15 ml/min 5% +15% (> 30 W)	15 ml/min 5% +15% (> 30 W)	15 ml/min 5% +15% (> 30 W)
Pre-curve °	0 to 15°	0 to 15°	0 to 15°	0 to 15°	0 to 15°	0 to 15°	0 to 15°	0 to 15°	0 to 15°
Full length, cm*	146 ± 3	146 ± 3	146 ± 3	146 ± 3	146 ± 3	146 ± 3	146 ± 3	146 ± 3	146 ± 3
Mass, g*	0±10%	0±10%	0±10%	0±10%	0±10%	0±10%	0±10%	0±10%	0±10%
Luer fitting									
Inner diameter, mm	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286
Major Outer diameter, mm	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762
Length, mm	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254
Cord length, mm	2,71 - 3,5	2,71 - 3,5	2,71 - 3,5	2,71 - 3,5	2,71 - 3,5	2,71 - 3,5	2,71 - 3,5	2,71 - 3,5	2,71 - 3,5
Side arm Cord length, mm	76.2 +/- 2.54	76.2 +/- 2.54	76.2 +/- 2.54	76.2 +/- 2.54	76.2 +/- 2.54	76.2 +/- 2.54	76.2 +/- 2.54	76.2 +/- 2.54	76.2 +/- 2.54
Tip Specifications									

Tip electrode size, mm (Fr)	2, 3 (7,5)	2,3 (7,5)	2,3 (7,5)	2,3 (7,5)	2,3 (7,5)	2,3 (7,5)	2,3 (7,5)	2,3 (7,5)	2,3 (7,5)
Number of Irrigation Holes	56	56	56	56	56	56	56	56	56
Irrigation Hole Diameter, mm	0,0889 +/- 0,508	0,0889 +/- 0,508	0,0889 +/- 0,508	0,0889 +/- 0,508	0,0889 +/- 0,508	0,0889 +/- 0,508	0,0889 +/- 0,508	0,0889 +/- 0,508	0,0889 +/- 0,508
Number of Lumens in the Deflectable Tip	4	4	4	4	4	4	4	4	4
Deflectable Tip Lumen Dimensions, mm	2 x 0,4826 +/- 0,0254 2 x 0,889 +/- 0,0254 2 x 0,762	2 x 0,4826 +/- 0,0254 2 x 0,889 +/- 0,0254 2 x 0,762	2 x 0,4826 +/- 0,0254 2 x 0,889 +/- 0,0254 2 x 0,762	2 x 0,4826 +/- 0,0254 2 x 0,889 +/- 0,0254 2 x 0,762	2 x 0,4826 +/- 0,0254 2 x 0,889 +/- 0,0254 2 x 0,762	2 x 0,4826 +/- 0,0254 2 x 0,889 +/- 0,0254 2 x 0,762	2 x 0,4826 +/- 0,0254 2 x 0,889 +/- 0,0254 2 x 0,762	2 x 0,4826 +/- 0,0254 2 x 0,889 +/- 0,0254 2 x 0,762	2 x 0,4826 +/- 0,0254 2 x 0,889 +/- 0,0254 2 x 0,762
Braiding	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")
Tip	3,6703	3,6703	3,6703	3,6703	3,6703	3,670	3,6703	3,670	3,670

Electro de Length, mm	+/- 0,5083, 5	+/- 0,5083, 5	+/- 0,5083, 5	+/- 0,5083, 5	+/- 0,5083, 5	3 +/- 0,5083, 5	+/- 0,5083, 5	3 +/- 0,5083, 5	3 +/- 0,5083, 5
Tip Electro de Constru ction	Two piece with stem	Two piece with stem	Two piece with stem	Two piece with stem	Two piece with stem	Two piece with stem	Two piece with stem	Two piece with stem	Two piece with stem
Numbe r of Ring Electro des	3	3	3	3	3	3	3	3	3
Additio nal Ring Electro des, mm	2	2	2	2	2	2	2	2	2
Ring Electro de Dimens ions, mm	Tip Ring Electro des OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length = 1,27 +/- 0.127	Tip Ring Electro des OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length = 1,27 +/- 0.127	Tip Ring Electro des OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length = 1,27 +/- 0.127	Tip Ring Electro des OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length = 1,27 +/- 0.127	Ring Electrode Dimensions, mm Tip Ring Electrodes OD = 2,4638 ID = 2,3876 Length = 1,27	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length = 1,27	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length = 1,27	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length = 1,27	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876 +/- 0.0254 Length = 1,27
	Shaft Ring Electro des OD = 2,4638 +/-	Shaft Ring Electro des OD = 2,4638 +/-	Shaft Ring Electro des OD = 2,4638 +/-	Shaft Ring Electro des OD = 2,4638 +/-	Shaft Ring Electrodes OD = 2,4638 ID =	Shaft Ring Electrodes OD = 2,4638 +/-	Shaft Ring Electrodes OD = 2,4638 +/-	Shaft Ring Electrodes OD = 2,4638 +/-	Shaft Ring Electrodes OD = 2,4638 +/-

	0.0254 ID = 2,3622+ /- 0.0254 Length = 0,7112+ /- 0.127	0.0254 ID = 2,3622 +/- 0.0254 Length = 0,7112 +/- 0.127	0.0254 ID = 2,3622+ /- 0.0254 Length = 0,7112+ /- 0.127	0.0254 ID = 2,3622 +/- 0.0254 Length = 0,7112 +/- 0.127	2,3622 Length= 0,7112	0.025 4 ID = 2,362 2+/- 0.025 4 Lengt h= 0,711 2+/- 0.127	0.0254 ID = 2,3622+/ - 0.0254 Length= 0,7112+/ - 0.127	0.025 4 ID = 2,362 2+/- 0.025 4 Lengt h= 0,711 2+/- 0.127	0.025 4 ID = 2,362 2+/- 0.025 4 Lengt h= 0,711 2+/- 0.127
Ring Electrode Spacing	2-5-2	2-5-2	2-5-2	2-5-2	2-5-2	2-5-2	2-5-2	2-5-2	2-5-2
RF Impedance, ohms when 5 Killohertz	0 to less than 5	0 to less than 5	0 to less than 5	0 to less than 5	0 to less than 5	0 to less than 5	0 to less than 5	0 to less than 5	0 to less than 5
RF Impedance ohms when 500 Killohertz	0 to less than 25	0 to less than 25	0 to less than 25	0 to less than 25	0 to less than 25	0 to less than 25	0 to less than 25	0 to less than 25	0 to less than 25
Isolation and Phase Angle Between Circuits	At 5 kHz: greater than 100k ohms. At 500 kHz: >1k ohms to infinity. The phase	At 5 kHz: greater than 100k ohms to infinity. At 500 kHz: greater than 1k ohms	At 5 kHz: greater than 100k ohms to infinity. At 500 kHz: greater than 1k ohms	At 5 kHz: greater than 100k ohms to infinity. At 500 kHz: greater than 1k ohms	At 5 kHz: greater than 100k ohms to infinity. At 500 kHz: greater than 1k ohms	At 5 kHz: greater than 100k ohms to infinity. At 500 kHz: greater than 1k ohms	At 5 kHz: greater than 100k ohm to infinity. At 500 kHz: greater than 1k ohms to infinity	At 5 kHz: greater than 100k ohm to infinity. At 500 kHz: greater than 1k ohms to infinity	At 5 kHz: greater than 100k ohms. At 500 kHz: greater

	angle must be between -86° to -90° .	than 1k ohms to infinity ohms. The phase angle must be between -86° . The phase angle must be between -86° to -90° .	to infinity ohms. The phase angle must be between -86° . The phase angle must be between -86° to -90° .	1k ohms to infinity ohms. The phase angle must be between -86° . The phase angle must be between -86° to -90° .	to infinity. The phase angle must be between -86° to -90° .	: greater than 1k ohms to infinity. The phase angle must be between -86° to -90° .	. The phase angle must be between -86° to -90° .	: greater than 1k ohms to infinity. The phase angle must be between -86° to -90° .	than 1k ohms to infinity. The phase angle must be between -86° to -90° .
RF Leakage Current, mA	0 to 252	0 to 252	0 to 252	0 to 252	0 to 252	0 to 252	0 to 252	0 to 252	0 to 252
Temperature Sensors (Thermistor and Thermocouple)	Thermocouple: 14004 Ohms $37^{\circ}\pm 2^{\circ}\text{C}$	Thermocouple: 14004 Ohms $37^{\circ}\pm 2^{\circ}\text{C}$	Thermocouple: 14004 Ohms $37^{\circ}\pm 2^{\circ}\text{C}$	Thermocouple: 14004 Ohms $37^{\circ}\pm 2^{\circ}\text{C}$	Thermocouple: 14004 Ohms $37^{\circ}\pm 2^{\circ}\text{C}$	Thermocouple: 14004 Ohms $37^{\circ}\pm 2^{\circ}\text{C}$	Thermocouple: 14004 Ohms $37^{\circ}\pm 2^{\circ}\text{C}$	Thermocouple: 14004 Ohms $37^{\circ}\pm 2^{\circ}\text{C}$	Thermocouple: 14004 Ohms $37^{\circ}\pm 2^{\circ}\text{C}$
Compression Coil Location	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.

Curve Types	BB	BD	BF	DD	DF	DJ	FF	FJ	JJ
Deflection Angle Range	0-180°	0-180°	0-180°	0-180°	0-180°	0-180°	0-180°	0-180°	0-180°
Tip Stiffness, N	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647
Buckle Force, N	0 to less than 0,5589	0 to less than 0,5589	0 to less than 0,5589	0 to less than 0,5589	0 to less than 0,5589	0 to less than 0,5589	0 to less than 0,5589	0 to less than 0,5589	0 to less than 0,5589
Side Force, N	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity
Pull Entire Catheter, N	Greater than 17,793 to infinity	Greater than 17,793 to infinity	Greater than 17,793 to infinity	Greater than 17,793 to infinity	Greater than 17,793 to infinity	Greater than 17,793 to infinity	Greater than 17,793 to infinity	Greater than 17,793 to infinity	Greater than 17,793 to infinity
Irrigation Pathway Specifications									
Irrigation Tubing Dimensions, mm	ID = 0,635+/- 0.0127	ID = 0,635+/- 0.0127	ID = 0,635+/- 0.0127	ID = 0,635+/- 0.0127	ID = 0,635+/- 0.0127	ID = 0,635+/- 0.0127	ID = 0,635+/- 0.0127	ID = 0,635+/- 0.0127	ID = 0,635+/- 0.0127
Shaft Specifications									
Shaft Size, mm (Fr)	2,3 (7.5)	2,3 (7.5)	2,3 (7.5)	2,3 (7.5)	2,3 (7.5)	2,3 (7.5)	2,3 (7.5)	2,3 (7.5)	2,3 (7.5)

Shaft Waviness, mm	0 to 20	0 to 20	0 to 20	0 to 20	0 to 20	0 to 20	0 to 20	0 to 20	0 to 20
Internal Lead Conductors	38 AWG Copper Wire	38 AWG Copper Wire	38 AWG Copper Wire	38 AWG Copper Wire	38 AWG Copper Wire	38 AWG Copper Wire	38 AWG Copper Wire	38 AWG Copper Wire	38 AWG Copper Wire
Temperature Sensor Conductor	Therm o- couple : • 38 AWG Const antan Wire • 38 AWG Copp er Wire	Ther mo- couple: • 38 AW G Const anta n Wire • 38 AW G Cop per Wire	Ther mo- couple: • 38 AWG Const antan Wire • 38 AWG Copp er Wire	Ther mo- couple: • 38 AW G Const anta n Wire • 38 AW G Cop per Wire	Therm o- couple : • 38 AWG Const antan Wire • 38 AWG Copp er Wire	The rmo - couple: • 38 AW G Co nsta nta n Wire • 38 AW G Co ppe r Wire	Therm o- couple: • 38 AWG Const antan Wire • 38 AWG Copp er Wire	The rmo - couple: • 38 AW G Co nsta nta n Wire • 38 AW G Co ppe r Wire	The rmo - couple: • 38 AW G Co nsta nta n Wire • 38 AW G Co ppe r Wire
Tip Conduc tor/ Lead Wire Joints	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper
Puller Wire Material Length, mm	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2
Compre	0,1016	0,1016	0,1016	0,1016	0,1016 x	,1016	,1016 ⁰ x	,1016	,1016

ssion Coil	x 0,1524	x 0,1524	x 0,1524	x 0,1524	0,1524	x 0,152 4	0,1524	x 0,152 4	x 0,152 4
Handle Specifications									
Constru ction	Handl e with adjust able deflec tion knob, tensio n knob and luer fitting	Hand le with adju stable defle ction knob, tensi on knob and luer fittin g	Handl e with adjust able deflec tion knob, tensio n knob and luer fitting	Hand le with adju stable defle ction knob, tensi on knob and luer fittin g	Handl e with adjust able deflect ion knob, tensio n knob and luer fitting	Han dle with adju stable defl ecti on knob, tens ion knob and luer fitti ng	Handle with adjusta ble deflecti on knob, tension knob and luer fitting	Han dle with adju stable defl ecti on knob, tens ion knob and luer fitti ng	Han dle with adju stable defl ecti on knob, tens ion knob and luer fitti ng
Length, cm	21.03 +/- 0.05	21.03 +/- 0.05	21.03 +/- 0.05	21.03 +/- 0.05	21.03 +/- 0.05	21.0 3+/- 0.05	21.03+ /- 0.05	21.0 3+/- 0.05	21.0 3+/- 0.05
Connec tor	25 pin	25 pin	25 pin	25 pin	25 pin	25 pin	25 pin	25 pin	25 pin
andle Temper ature, °C	to 48	to 48	to 48	to 48	to 48	to 48	to 48	to 48	to 48
3. ThermoCool SF NAV Uni-Directional Catheter UNI with curve visualization									
Characteristi cs	D131501			D131502		D131503		D131504	
General									
Usable Length, cm	115 ± 3			115 ± 3		115 ± 3		115 ± 3	
Puller wire material length, mm	1524 ± 76,2			1524 ± 76,2		1524 ± 76,2		1524 ± 76,2	

Outer diameter, mm	0 to 2.6416	0 to 2.6416	0 to 2.6416	0 to 2.6416
Application Time, sec	not exceed 120	not exceed 120	not exceed 120	not exceed 120
Power Settings During RF Application (atrial Ablation), W	15 to 30	15 to 30	15 to 30	15 to 30
Power Settings During RF Application (ventricular Ablation), W	31 to 50	31 to 50	31 to 50	31 to 50
Flow Rate, ml/min	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)
	15 ml/min-5% +15% (> 30 W)	15 ml/min-5% +15% (> 30 W)	15 ml/min-5% +15% (> 30 W)	15 ml/min-5% +15% (> 30 W)
Pre-curve, °	0 to 15°	0 to 15°	0 to 15°	0 to 15°
Full length, cm*	146 ± 3	146 ± 3	146 ± 3	146 ± 3
Mass, g*	90±10%	90±10%	90±10%	90±10%
Luer fitting				
Inner diameter, mm	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286
Length, mm	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254
Side arm Cord length, mm	215,9 +/- 3,175	215,9 +/- 3,175	215,9 +/- 3,175	215,9 +/- 3,175
Tip Specifications				
Tip electrode	2,3 (7.5)	2,3 (7.5)	2,3 (7.5)	2,3 (7.5)

size, mm (Fr)				
Number of Irrigation Holes	56	56	56	56
Irrigation Hole Diameter, mm	0,0889+/- 0,508	0,0889+/- 0,508	0,0889+/- 0,508	0,0889+/- 0,508
Number of Lumens in the Deflectable Tip	4	4	4	4
Deflectable Tip Lumen Dimensions, mm	2 x 0,4826+/- 0.0254 2 x 0,762+/- 0.0254	2 x 0,4826+/- 0.0254 2 x 0,762+/- 0.0254	2 x 0,4826+/- 0.0254 2 x 0,762+/- 0.0254	2 x 0,4826+/- 0.0254 2 x 0,762+/- 0.0254
Braiding	Braided: 16 strands, 45 +/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")
Tip Electrode Length, mm	3,6703 +/- 0,5083,5	3,6703 +/- 0,5083,5	3,6703 +/- 0,5083,5	3,6703 +/- 0,5083,5
Tip Electrode Constructio n	Two piece with stem	Two piece with stem	Two piece with stem	Two piece with stem
Number of Ring Electrodes	3	3	3	3
Additional Ring Electrodes, mm	2	2	2	2
Ring Electrode Dimensions, mm	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876+/- 0.0254	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876+/-	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876+/-	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876+/-

	Length= 1,27 +/- 0.127Tip Ring Electrodes OD = 2,4638 ID = 2,3876 Length= 1,27	0.0254 Length= 1,27+/- 0.127Tip Ring Electrodes OD = 2,4638 ID = 2,3876 Length= 1,27	0.0254 Length= 1,27+/- 0.127Tip Ring Electrodes OD = 2,4638 ID = 2,3876 Length= 1,27	0.0254 Length= 1,27+/- 0.127Tip Ring Electrodes OD = 2,4638 ID = 2,3876 Length= 1,27
	Shaft Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3622+/- 0.0254 Length= 0,7112+/- 0.127Shaft Ring Electrodes OD = 2,4638 ID = 2,3622 Length= 0,7112	Shaft Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3622+/- 0.0254 Length= 0,7112+/- 0.127Shaft Ring Electrodes OD = 2,4638 ID = 2,3622 Length= 0,7112	Shaft Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3622+/- 0.0254 Length= 0,7112+/- 0.127Shaft Ring Electrodes OD = 2,4638 ID = 2,3622 Length= 0,7112	Shaft Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3622+/- 0.0254 Length= 0,7112+/- 0.127Shaft Ring Electrodes OD = 2,4638 ID = 2,3622.028 Length= 0,7112
Ring Electrode Spacing	2-5-2	2-5-2	2-5-2	2-5-2
RF Impedance ohms when 5 Killohertz	0 to less than 5	0 to less than 5	0 to less than 5	0 to less than 5
RF Impedance ohms when 500 Killohertz	0 to less than 25	0 to less than 25	0 to less than 25	0 to less than 25
Isolation and Phase Angle Between Circuits	At 5 kHz: greater than 100k to infinity ohms. At 500 kHz: greater than 1k to infinity ohms. The phase angle must be between -86° to -90°.	At 5 kHz: greater than 100k to infinity ohms. At 500 kHz: greater than 1k to infinity ohms. The phase angle must be	At 5 kHz: greater than 100k to infinity ohms. At 500 kHz: greater than 1k to infinity ohms. The phase angle must be between -86°	At 5 kHz: greater than 100k to infinity ohms. At 500 kHz: greater than 1k to infinity ohms. The phase angle must be

		between -86° to -90°.	to -90°.	between -86° to -90°.
RF Leakage Current, mA	0 to 252	0 to 252	0 to 252	0 to 252
Temperature Sensors (Thermistor and Thermocouple)	Thermo-couple: 14004 Ohms 37°+/- 2°C	Thermo-couple: 14004 Ohms 37°+/- 2°C	Thermo-couple: 14004 Ohms 37°+/- 2°C	Thermo-couple: 14004 Ohms 37°+/- 2°C
Compression Coil Location	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.
Curve Types	B	D	F	J
Deflection Angle Range	0-150°	0-150°	0-150°	0-150°
Tip Stiffness, N	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647
Buckle Force, N	0 to less than 0,5589	0 to less than 0,5589	0 to less than 0,5589	0 to less than 0,5589
Side Force, N	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity
Pull Entire Catheter, N	Greater than 17,793 to infinity	Greater than 17, 793 to infinity	Greater than 17, 793 to infinity	Greater than 17, 793 to infinity
Irrigation Pathway Specifications				
Irrigation Tubing Dimensions, mm	ID = 0,635+/-+/- 0.0127	ID = 0,635+/-+/- 0.0127	ID = 0,635+/-+/- 0.0127	ID = 0,635+/-+/- 0.0127
Shaft Specifications				
Shaft Size, mm (Fr)	2,3 (7,5)	2,3 (7,5)	2,3 (7,5)	2,3 (7,5)
Shaft Waviness, mm	0 to 20	0 to 20	0 to 20	0 to 20
Internal	38 AWG Copper	38 AWG Copper	38 AWG Copper	38 AWG Copper

Lead Conductors	Wire	Wire	Wire	Wire
Temperature Sensor Conductor	Thermo-couple: • 38 AWG Constantan Wire • 38 AWG Copper Wire	Thermo-couple: • 38 AWG Constantan Wire • 38 AWG Copper Wire	Thermo-couple: • 38 AWG Constantan Wire • 38 AWG Copper Wire	Thermo-couple: • 38 AWG Constantan Wire • 38 AWG Copper Wire
Tip Conductor/ Lead Wire Joints	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper
Puller Wire Material length, mm	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2
Compression Coil	0,1016 x 0,1524	0,1016 x 0,1524	0,1016 x 0,1524	0,1016 x 0,1524
Handle Specifications				
Construction	Handle with adjustable deflection knob, tension knob and luer fitting	Handle with adjustable deflection knob, tension knob and luer fitting	Handle with adjustable deflection knob, tension knob and luer fitting	Handle with adjustable deflection knob, tension knob and luer fitting
Length, cm	21.03+/- 0.05	21.03+/- 0.05	21.03+/- 0.05	21.03+/- 0.05
Length of the Piston, cm	6.985 +/- 0.0762	6.985 +/- 0.0762	6.985 +/- 0.0762	6.985 +/- 0.0762
Connector	34 pin	34 pin	34 pin	34 pin
Handle Temperature, °C	0 to 48	0 to 48	0 to 48	0 to 48
4. ThermoCool SF NAV Bi-Directional Catheter with curve visualization				
Characteristics	BNI35BDCT	BNI35BFCT		BNI35JJCT
General				

Usable Length, cm	115 ± 3	115 ± 3	115 ± 3
Puller wire material length, mm	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2
Outer diameter, mm	0 to 2.6416	0 to 2.6416	0 to 2.6416
Application Time, sec	not exceed 120	not exceed 120	not exceed 120
Power Settings During RF Application, W	15 to 30	15 to 30	15 to 30
Flow Rate, ml/min	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)
	15 ml/min-5% +15% (> 30 W)	15 ml/min-5% +15% (> 30 W)	15 ml/min-5% +15% (> 30 W)
Pre-curve, °	0 to 15°	0 to 15°	0 to 15°
Full length, cm*	146 ± 3	146 ± 3	146 ± 3
Mass, g*	90±10 %	90±10%	90±10%
Luer fitting			
Inner diameter, mm	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286
Major Outer diameter, mm	7.747 +0.000 - 0.0762	7.747 +0.000 -0.0762	7.747 +0.000 - 0.0762
Length, mm	20.066 +/- 0.254	20.066 +/- 0.254	20.066 +/- 0.254
Side arm Cord length, mm	76.2 +/- 2.54	76.2 +/- 2.54	76.2 +/- 2.54
Tip Specifications			
Tip electrode size, mm (Fr)	2,3 (7.5)	2, 3 (7.5)	2, 3 (7.5)
Number of Irrigation Holes	56	56	56
Irrigation Hole Diameter, mm	0,0889+/- 0,508	0,0889+/- 0,508	0,0889+/- 0,508

Number of Lumens in the Deflectable Tip	4	4	4
Deflectable Tip Lumen Dimensions, mm	2 x 0,4826 2 x 0,762	2 x 0,4826 2 x 0,762	2 x 0,4826 2 x 0,762
Braiding	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45+/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45 +/-1 picks per 2,54 cm (1")
Tip Electrode Length, mm	3,6703 +/- 0,5083,5	3,6703 +/- 0,5083,5	3,6703 +/- 0,5083,5
Tip Electrode Construction	Two piece with stem	Two piece with stem	Two piece with stem
Number of Ring Electrodes	3	3	3
Additional Ring Electrodes, mm	2	2	2
Ring Electrode Dimensions, mm	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876+/- 0.0254 Length= 1,27+/- 0.127	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876+/- 0.0254 Length= 1,27+/- 0.127	Tip Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3876+/- 0.0254 Length= 1,27+/- 0.127
	Shaft Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3622+/- 0.0254 Length= 0,7112+/- 0.127	Shaft Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3622+/- 0.0254 Length= 0,7112+/- 0.127	Shaft Ring Electrodes OD = 2,4638 +/- 0.0254 ID = 2,3622+/- 0.0254 Length= 0,7112+/- 0.127
Ring Electrode Spacing	2-5-2	2-5-2	2-5-2
RF Impedance, ohms	0 to less than 5	0 to less than 5	0 to less than 5
	0 to less than 25	0 to less than 25	0 to less than 25
Isolation and Phase	At 5 kHz:	At 5 kHz: greater than 100k	At 5 kHz: greater

Angle Between Circuits	greater than 100k ohms to infinity. At 500 kHz: greater than 1k ohms to infinity. The phase angle must be between -86° to -90°.	ohms to infinity. At 500 kHz: greater than 1k ohms to infinity. The phase angle must be between -86° to -90°.	than 100k ohms to infinity. At 500 kHz: greater than 1k ohms to infinity. The phase angle must be between -86° to -90°.
RF Leakage Current, mA	0 to 252	0 to 252	0 to 252
Temperature Sensors (	Thermo-couple: 14004 Ohms 37°+/- 2°C	Thermo-couple: 14004 Ohms 37°+/- 2°C	Thermo-couple: 14004 Ohms 37°+/- 2°C
Compression Coil Location	Glued in stiffener.	Glued in stiffener.	Glued in stiffener.
Curve Types	BD	BF	JJ
Deflection Angle Range	0-180°	0-180°	0-180°
Tip Stiffness, N	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647
Buckle Force, N	0 to less than 0,589	0 to less than 0,5589	0 to less than 0,5589
Side Force, N	Greater than 0,04 to infinity	Greater than 0,04 to infinity	Greater than 0,04 to infinity
Pull Entire Catheter, N	Greater than 17,793 to infinity	Greater than 17, 793 to infinity	Greater than 17, 793 to infinity
Irrigation Pathway Specifications			
Irrigation Tubing Dimensions, mm	ID = 0,635+/- 0.0254	ID = 0,635+/- 0.0254	ID = 0,635+/- 0.0254
Shaft Specifications			
Shaft Size, mm (Fr)	2,3 (7,5)	2,3(7,5)	2,3(7,5)
Shaft Waviness, mm	0 to 20	0 to 20	0 to 20

Internal Conductors	Lead	38 AWG Copper Wire	38 AWG Copper Wire	38 AWG Copper Wire
Temperature Conductor	Sensor	Thermo-couple: • 38 AWG Constantan Wire • 38 AWG Copper Wire	Thermo-couple: • 38 AWG Constantan Wire • 38 AWG Copper Wire	Thermo-couple: • 38 AWG Constantan Wire • 38 AWG Copper Wire
Tip Conductor/ Lead Wire Joints		3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper	3 twisted pair of 48 AWG copper
Puller Wire, mm		1524 ± 76,2	1524 ± 76,2	1524 ± 76,2
Compression Coil		0,1016 x 0,1524	0,1016 x 0,1524	0,1016 x 0,1524
Handle Specifications				
Construction		Handle with adjustable deflection knob, tension knob and luer fitting	Handle with adjustable deflection knob, tension knob and luer fitting	Handle with adjustable deflection knob, tension knob and luer fitting
Length, cm		21.03+/- 0.05	21.03+/- 0.05	21.03+/- 0.05
Connector		34 pin	34 pin	34 pin
Handle Temperature, °C		0 to 48	0 to 48	0 to 48
NOTE: Tolerance for these parameters were established using Russian GOST R 50444-				

9 Information of EMC tests (EMC Test Report)

The following tables present data regarding compliance of the CARTO® 3 System with IEC 60601-1-2, 3rd edition.

Table 1 — Guidance and Manufacturer's Declaration – Electromagnetic emissions

Table 2 — The CARTO® 3 System is intended for use in the electromagnetic environment specified below. The customer or the user of the CARTO® 3 System should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The CARTO® 3 System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. See CARTO® 3 System Interference with Other Medical Equipment.
RF emissions CISPR 11	Class A	The CARTO® 3 System is suitable for use in all establishments other than domestic.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Compliance	

WARNING This system (including its peripheral components) is intended for use by healthcare professionals only. This system (including its peripheral components) might cause radio interference or might disrupt the operation of nearby equipment. Install the CARTO® 3 System only in an electrophysiological laboratory with fluoroscopic (X-ray) capabilities.

The CARTO[®] 3 System is intended for use in the electromagnetic environment specified below. The customer or the user of the CARTO[®] 3 System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient or burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output	±2 kV for power supply lines ±1 kV for input/output	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	Mains power quality should be that of typical commercial or hospital environment. If the user of the CARTO [®] 3 System requires continue operation during power mains interruptions, it is recommended that the CARTO [®] 3 System be powered from an uninterruptible power supply (UPS).
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

NOTE UT is the AC mains voltage prior to application of the test level.

Table 13 — Recommended separation distances between portable and mobile RF communications equipment and system with Cardiac Catheters

Recommended separation distances between portable and mobile RF communications equipment and the CARTO® 3 System			
The CARTO® 3 System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CARTO® 3 System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CARTO® 3 System as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter, W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,17\sqrt{P}$	80 MHz to 800 MHz $d = 1,17\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,33\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			

NOTES: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

These guidelines might not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

10 Shelf life

The ThermoCool SF NAV Uni-directional and ThermoCool SF NAV Bi-directional Catheters are validated for a one-year product shelf life. Catheter integrity and performance testing was conducted after sterilization by ethylene oxide (EtO) through three commercial cycles followed by accelerated aging (equivalent to one year real time aging), thermal cycling and simulated transportation conditioning. Acceptance criteria were met for all tests.

11 The protection of the environment

11.1 Requirements for the protection of the environment in the application of medical device

Using ThermoCool SF NAV Uni-directional or ThermoCool SF NAV Bi-directional Catheter does not imply compliance with specific requirements for environmental protection. ThermoCool SF NAV Uni-directional and ThermoCool SF NAV Bi-directional Catheters are

EU RoHS (2011/65/EU) and REACH (1907/2006) compliant, are free of latex and latex handling, and are phthalate-free.

11.2 Utilization method/ Disposal

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

Disposal or destruction of the ThermoCool SF NAV Uni-directional and ThermoCool SF NAV Bi-Directional Catheters shall be in accordance with the classification (Class A – for catheters with damaged packaging, class B – for used devices), rules of collection, use, decontamination, placement, storage, transportation, accounting and utilization of medical waste, established by the authorized federal executive body.

12 Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced)

All Biosense Webster catheters are ethylene oxide (EtO) sterilized by Steris Isomedix Inc., a subcontractor for Biosense Webster. The Steris sterilization process utilizes ethylene oxide and incorporates all phases of sterilization (preconditioning, sterilization and aeration) within the sterilization chamber. This process is referred to as “EO Express” since the use of traditional preconditioning and aeration rooms are not necessary. The validation of Biosense Webster products resulted in one common cycle for all Biosense Webster device families. Sterility assurance was verified to exceed a sterility assurance level (SAL) of 1×10^{-6} . The validation met the requirements of EN 550 and ISO 11135. Residual results for EO and ECH are within allowable limits per ISO 10993-7.

13 The list of national and international regulatory documents/standards, which matches medical device

Table 14 — The list standards

Standards	Title
EN ISO 13485	Medical Device Quality Management System
EN ISO 14971	Application of Risk Management to Medical Devices
IEC 60601-1-2	Medical Electrical Equipment-Part 1: General Requirements for safety. Collateral standard: Electromagnetic Compatibility – Requirements and Tests
IEC 61000-4-3	Electromagnetic compatibility (EMC) Part 4-3: Testing and Measurement Techniques- Radiated, Radio-frequency, Electromagnetic Field Immunity Test
IEC 61000-4-6	Electromagnetic Compatibility (EMC) - Part 4-6: Testing and Measurement Techniques - Immunity to Conducted Disturbances, Induced by Radio-frequency Fields
ANSI/AAMI HF18	Standards for Electrosurgical Devices
BS EN ISO 10555-1	Sterile, Single-Use Intravascular Catheters
EN ISO 11135	Sterilization of Health Care Products - Ethylene Oxide
EN 10993-1	Biological Evaluation of Medical Devices – Part 1: Evaluation and Testing
EN 868-1	Packaging Materials for Devices to be Sterilized

ISTA 1G	Packaged-Products Weighing 150 lb (68 kg) or Less (Random Vibration)
ASTM D4332	Standard Practice for Conditioning Containers, Packages or Packaging Components for Testing
ASTM F2096-02	Standard Test Method for Detecting Gross Leaks in Medical Packaging by Internal Pressurization (Bubble Test)
EN ISO 10993-7	Biological Evaluation of Medical Devices – Part 7: Ethylene Oxide Sterilization Residuals
ANSI/AAMI HF18	Standards for Electrosurgical Devices
ISO 15223	Clean rooms and Associated Controlled Environments
EN 556	Sterilization Of Medical Devices
EN 980	Medical Symbols
EN 1041	Information from Manufacturer
ASTM F1980-02	Standard Guide for Accelerated Aging of Sterile Medical Device Packages

14 Customer Service

Contact for all questions, including complaints and service of medical devices the Authorized Representative in Russia:

Johnson & Johnson LLC
 17, Krylatskaya str., bld. 2, Moscow, 121614
 Tel.: (495) 580-77-77 Fax: (495) 580-78-78

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INSTRUCTION FOR USE

THE NAVIGATION CATHETER FOR CARDIAC ABLATION:

THERMOCOOL SMARTTOUCH SF UNI-DIRECTIONAL AND BI-DIRECTIONAL
CATHETERS

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Richard Lauhead
Senior Regulatory Affairs Specialist

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
On May 15, 2017 before me, Richard G. Perez, NOTARY PUBLIC, personally appeared Richard Lauhead, who proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that he executed the same in his authorized capacity and that his signature on the instrument the person, or the entity upon behalf of which the person acted, executed the instrument.

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**I Certify That This Is An
Original Document**

Betty Valverde

District of Columbia: SS

Subscribed and sworn to before me

this 17 day of May, 2017

Kamal S. Ahmad, Notary Public, D.C.

My Commission expires December 14, 2019

THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and The Biosense Webster THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

- **STERILE. Sterilized using ethylene oxide.**
- **For single use only.**
- **Do not resterilize.**
- **Do not use if the package is open or damaged.**

Catheter Description

The THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and The THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter is a multi-electrode luminal catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) current to the catheter tip electrode for ablation purposes.

The catheter shaft measures 7.5 F with 8 F ring electrodes. For ablation, the catheter is used in conjunction with an RF generator and a dispersive pad (indifferent electrode). The catheter has force-sensing technology that provides a real-time measurement of contact force between the catheter tip and the heart wall.

The Uni-Directional Navigation catheter has a high-torque shaft with a uni-directional deflectable tip section containing an array of electrodes which includes a 3.5 mm tip dome. All of the electrodes may be used for recording and stimulation purposes. The tip electrode serves to deliver RF current from the RF generator to the desired ablation site. The tip electrode and ring electrodes are made from noble metals.

The Bi-Directional Navigation catheter has a high-torque shaft with a bi-directional deflectable tip section containing an array of electrodes which includes a 3.5 mm tip dome. All of the electrodes may be used for recording and stimulation purposes. The tip electrode serves to deliver RF current from the RF generator to the desired ablation site. The tip electrode and ring electrodes are made from noble metals.

The Uni-Directional Navigation catheter incorporates a thermocouple temperature sensor that is embedded in the 3.5 mm tip electrode. Tip deflection is controlled at the proximal end by a handpiece in which a piston slides; a thumbknob on the piston controls piston travel. When the thumbknob is pushed forward, the tip is deflected (curved).

When the thumbknob is pulled back, the tip straightens. The high-torque shaft also allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site. Currently, the available curves for the uni-directional version of the THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter include D, F and J.

The Bi-Directional Navigation catheter incorporates a thermocouple temperature sensor that is embedded in the 3.5 mm tip electrode. A Rocker Lever is used to deflect the tip. The high torque

shaft also allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site. Additionally, a variety of curve types are available in symmetric or asymmetric combinations, providing two 180° opposed, single planed curves. Currently, the available curves for the THERMOCOOL SMARTTOUCH® SF Bi-Directional Navigation Catheters include DD, FF, JJ, DF and FJ.

At the proximal end of the Uni-Directional Navigation and Bi-Directional Navigation catheters, a saline input port with a standard Luer fitting terminates from the open lumen. This saline port serves to permit the injection of normal saline to irrigate the tip electrode. During ablation, heparinized normal saline is passed through the internal lumen of the catheter and through the tip electrode, to irrigate and cool the ablation site as well as the electrode tip. An irrigation pump should be used to control the saline irrigation. The Uni-Directional Navigation and Bi-Directional Navigation catheter interfaces with standard recording equipment and a compatible RF generator via accessory extension cables with the appropriate connectors.

This Uni-Directional Navigation and Bi-Directional Navigation catheters features a location sensor embedded in the tip section that transmits location and contact force information to the CARTO 3 Navigation System. An appropriate reference device is required for location reference position purposes. For use in mapping procedures and for information on appropriate reference devices, and for further description of the operation of the CARTO® 3 System, refer to the User Manual for your CARTO® 3 Navigation System.

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Indications

The THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and The THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter is indicated for use in cardiac electrophysiological mapping (stimulation and recording) and, when used in conjunction with a radiofrequency generator, for cardiac ablation.

The THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and The THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter provides a real-time measurement of contact force between the catheter tip and heart wall, as well as location information when used with CARTO 3 Navigation System.

Contraindications

Do not use this device:

1. If the patient has had a ventriculotomy or atriotomy within the preceding eight weeks because the recent surgery may increase the risk of perforation.
2. In the patient with a myxoma or an intracardiac thrombus as the catheter could precipitate an embolus.
3. In patients with prosthetic valves as the catheter may damage the prosthesis.
4. In the coronary vasculature due to risk of damage to the coronary arteries.

5. In patients with an active systemic infection because this may increase the risk of cardiac infection.
6. Via the transseptal approach in a patient with an interatrial baffle or patch because the opening could persist and produce an iatrogenic atrial shunt.
7. Via the retrograde trans-aortic approach in patients who have had aortic valve replacement.
8. With a long sheath or short introducer < 8.5 French in order to avoid damage to the catheter shaft.

Warnings and Precautions

1. Do not use excessive force to advance or withdraw the catheter when resistance is encountered during catheter manipulation through the sheath.
2. Do not manually pre-shape the distal shaft of the catheter by applying external forces intended to bend or affect the intended shape or curve of the catheter.
3. The contact force reading is for information only and is not intended to replace standard handling precautions.
4. The catheter must be warmed up as specified prior to use. If the catheter has not reached a steady state condition, there is potential for a zero-offset drift to occur which could result in an inaccurate contact force reading.
5. Always zero the contact force reading following insertion into the patient or when moving the catheter from one chamber of the heart to another. Ensure the catheter is not in contact with heart tissue prior to zeroing. Refer to the User Manual for your CARTO 3 System for instructions on how to zero the contact force reading.
6. The contact force reading might become inaccurate if the contact force sensor (located between the first and second ring electrode) comes into close proximity with a ferrous material, such as the braided shaft of another catheter. If extreme fluctuations in force are observed, ensure the catheter's contact force sensor is not in close proximity with another catheter's shaft, check zero on the catheter and, if necessary, remove and inspect the catheter.
7. To ensure proper operation of the contact force sensor, all four ring electrodes of the catheter tip must protrude from the distal tip of the guiding sheath.
8. When applying high lateral force during mapping and RF application, the user should monitor the contact force Dashboard and vector display on the CARTO 3 screen to ensure that contact force measurements remain within the accurate range. Refer to the Error Messages and Alerts section of the CARTO 3 System Instructions for Use for System-related alert messages and indications related to inaccurate force readings.
9. Do not use the temperature sensor to monitor tissue temperature. The temperature sensor located within the tip section of the catheter does not reflect either electrode-tissue interface or tissue temperature due to the cooling effects of the saline irrigation of the electrode. The temperature display of the RF generator is the temperature of the cooled electrode, not tissue temperature. The

temperature sensor is used to verify that the irrigation flow rate is adequate. Before initiating the application of RF current, a decrease in electrode temperature confirms the onset of saline irrigation of the ablation electrode. Monitoring the temperature from the electrode during the application of RF current ensures that the irrigation flow rate is being maintained.

10. Do not rely on electrode temperature rise to determine if tissue heating is occurring during RF energy delivery as bench and animal studies showed no significant electrode temperature rise during RF ablation.

11. It is important to carefully follow the power titration procedure as specified in the instructions for use. Too rapid an increase in power during ablation may lead to perforation caused by steam pop.

12. This catheter may damage the prosthetic tricuspid valve of a patient if the catheter is accidentally advanced through the valve.

13. The patient who has had a prior atrial flutter ablation procedure may be at greater risk for perforation and/or pericardial effusion with the use of this catheter system.

14. In accordance with your hospital's protocol, monitor the patient's fluid balance throughout the procedure to avoid fluid volume overload. Some patients may have factors that reduce their ability to handle the volume overload, making them susceptible to developing pulmonary edema or heart failure during or after the procedure. Patients with congestive heart failure or renal insufficiency, and the elderly are particularly susceptible. Prior to the procedure, always identify the patient's risk of volume overload.

15. The safety of discontinuing anticoagulation therapy following catheter ablation of atrial fibrillation has not been established; anticoagulation therapy in such patients should be administered in accordance with the ACC/AHA/ESC Guidelines for the Management of Patients with Atrial Fibrillation.

16. To avoid thromboemboli, intravenous heparin should be used when entering the left heart during ablation. Follow the clinical guideline and common practice for optimal post-procedural anticoagulation.

17. The safety and effectiveness of radiofrequency ablation for the treatment of atrial fibrillation in patients with significant left ventricular dysfunction, advanced heart failure, substantial left atrial enlargement, and structural heart disease have not been established.

18. The catheter has not been shown to be safe at electrode temperatures above 40°C. Verify the CATHETER SELECTION KNOB on the compatible RF generator is on the "TCool SF" or Thermocool SF similar option and ensure that the maximum temperature is set at 40°C.

19. Implantable pacemakers and implantable cardioverter/defibrillator (ICDs) may be adversely affected by RF current. It is important to have temporary external sources of pacing and defibrillation available during ablation and to temporarily reprogram the pacing system to minimum output or OFF mode to minimize the risk of inappropriate pacing. Exercise extreme caution during ablation when in close proximity to atrial or ventricular permanent leads; program

the ICD to the OFF mode during the ablation procedure; and, perform complete implantable device analysis on all patients after ablation.

20. Patients undergoing septal accessory pathway ablation are at risk for complete AV block which requires the implantation of a permanent pacemaker. Patients who experience inadvertent complete AV block as a result of RF ablation may also require permanent pacing.

21. During the trans-aortic approach, adequate fluoroscopic visualization is necessary to avoid placement of the catheter in the coronary vasculature. Intracoronary placement of the ablation catheter, RF energy application, or both have been associated with myocardial infarction.

22. Minimize X-ray exposure during the procedure. Catheter ablation procedures present the potential for significant X-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the X-ray beam intensity and duration of the fluoroscopic imaging. Catheter ablation should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps have been taken to minimize this exposure. Careful consideration must therefore be given for the use of the device in pregnant women.

23. Do not expose catheter to organic solvents such as alcohol.

24. Do not autoclave the catheter.

25. Do not immerse proximal handle or cable connector in fluids; electrical performance could be affected.

26. Do not scrub or twist the distal tip electrode during cleaning.

27. Inspect irrigation saline for air bubbles prior to its use in the procedure. Air bubbles in the irrigation saline may cause emboli.

28. Purge catheter and irrigation tubing with heparinized normal saline.

29. Electrophysiology catheters and systems are intended for use only in X-ray shielded rooms due to electromagnetic compatibility requirements and other hospital safety guidelines.

30. Do not attempt to operate the THERMOCOOL SMARTTOUCH® SF Uni-Directional Navigation Catheter and THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter or the RF generator prior to completely reading and understanding the applicable instructions for use.

31. Cardiac ablation procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory. Appropriate clinical instruction in the use of the THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter should also be completed.

32. The long-term risks of protracted fluoroscopy and creation of RF induced lesions have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children. Furthermore, the risk/benefit in asymptomatic patients has not been studied.

33. When using the THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter with conventional systems (using fluoroscopy to determine catheter tip location), or with the CARTO 3 Navigation System, careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered. The firmness of the braided tip dictates that care must be taken to prevent perforation of the heart. The contact force reading is for information only and is not intended replace standard handling precautions.

34. Always pull the thumbknob back when using the THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter or place the Rocker Lever in the neutral position to when using THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter straighten the catheter tip before insertion or withdrawal of the catheter.

35. Always maintain a constant heparinized normal saline infusion to prevent coagulation within the lumen of the catheter.

36. When RF current is interrupted for either a temperature or an impedance rise (the set limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum, if present. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode or may damage the contact force sensor. A significant change in the baseline reading after cleaning might indicate a damaged contact force sensor. Make sure the irrigation holes are not plugged prior to re-insertion.

37. Apparent low power output, high impedance reading, or failure of the equipment to function correctly at normal settings may indicate faulty application of the indifferent electrode(s) or failure of an electrical lead. Do not increase power before checking for obvious defects or misapplication of the indifferent electrode or other electrical leads.

38. Read and follow the indifferent electrode manufacturer's instructions for use; the use of indifferent electrodes that meet or exceed ANSI/AAMI requirements (AAMI IEC 60601-2-2), is recommended.

39. The THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter are intended for use with compatible RF generators, a compatible irrigation pump, CARTO 3 System and Biosense Webster cables, and other appropriate interface cables and connectors. Use of a compatible irrigation pump is recommended to assure proper irrigation flow rate.

40. The THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter has been shown to create larger lesions than standard non-irrigated RF ablation catheters. Care should be taken when ablating near structures such as the sino-atrial and atrioventricular nodes.

41. The sterile packaging and catheter should be inspected prior to use. Do not use if the packaging or catheter appears damaged.

42. The catheter is sterilized with ethylene oxide gas and should be used by the “Use By” date on the device package. Do not use the catheter if past the “Use By” date.
43. The THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter are intended for single patient use only.
44. Do not use near MRI equipment since movement or heating of the catheter may occur and the image on the display may become distorted.
45. Use both fluoroscopy and electrogram data to monitor catheter advancement and reduce risk of tissue injury.
46. The THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter used in conjunction with a RF generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the catheter and indifferent electrode, particularly when operating the catheter. During energy delivery, the patient should not be allowed to come in contact with grounded metal surfaces.
47. The risk of igniting flammable gases or other materials is inherent in electrosurgery. Precautions must be taken to restrict flammable materials from the electrosurgical suite.
48. Electromagnetic interference (EMI) produced by the THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter, when used in conjunction with an RF generator during normal operation, may adversely affect the performance of other equipment.
49. Electrodes and probes for monitoring and stimulating devices can provide paths for high frequency current. The risk of burns can be reduced but not eliminated by placing the electrodes and probes as far away as possible from the ablation site and the indifferent electrode. Protective impedances may reduce the risk of burns, and permit continuous monitoring of the electrocardiogram during energy delivery.
50. The temperature sensor measures electrode tip temperature, not tissue temperature. The temperature displayed on the RF generator is for the cooled electrode only and does not represent tissue temperature. If the RF generator does not display temperature, verify that the appropriate cable is plugged into the RF generator. If temperature still is not displayed, there may be a malfunction in the temperature sensing system that must be corrected prior to applying RF power.
51. The temperature measurement accuracy of the THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter and THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter, as with any temperature measurement electrophysiology catheter, is largely determined by the temperature accuracy specification of the RF generator used. Please consult the user manual for the RF generator to be used for the temperature accuracy specification.
52. Before use, check irrigation ports are patent by infusing heparinized normal saline through the catheter and tubing.
53. Regularly inspect and test reusable cables and accessories.

- 54. When ablating near adjacent anatomical structures, take precautions to minimize collateral damage to the adjacent structures.
- 55. When ablating near the esophagus (along the posterior wall of the left atrium), take precautions to avoid injuring the esophagus, including appropriately reducing RF power.
- 56. When ablating near the phrenic nerve, take precautions to avoid injuring the phrenic nerve, including appropriately reducing RF power and pacing to identify the proximity of the nerve.

RF Ablation

For RF ablation the catheter is connected to the CARTO® 3 System Patient Interface Unit (PIU), which connects to the RF generator. For setup procedures refer to the User Manual for your CARTO 3 Navigation System. For proper RF generator interface, use only a Biosense Webster or compatible interface cable. To complete the electrical circuit, an indifferent electrode must be connected to the indifferent electrode input on the RF generator. Verify that circuit impedance prior to RF ablation is within expected parameters. Verify that the RF generator displays a temperature not above 37°C after the catheter is inserted into the patient and before applying RF power.

RF Generator Operation

Refer to the applicable RF generator manual for proper connection of the catheter to the generator and for detailed instructions as to generator operation for RF ablation. RF ablation application parameters will vary depending on the ablation site, the specific conditions present in each procedure and the RF generator control circuitry. Based on data obtained from prior animal and clinical studies, recommended RF application parameters are provided below in the “Directions for Use” and in Table 1. Always monitor temperature and impedance rise when using the THERMOCOOL SMARTTOUCH® SF Uni-Directional Navigation Catheter and THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter.

Table 1.

RECOMMENDED RF APPLICATION PARAMETERS		
	ATRIAL ABLATION	VENTRICULAR ABLATION
Power Range	15 W to 30 W*	31W to 50 W
Temperature Monitoring	<40°C**	<40°C**
Irrigation Flow Rate During RF Application	8 ml/min	15 ml/min
Application Time	30 to 120 seconds	60 to 120 seconds

* Power levels exceeding 30 Watts may be used when transmural lesions cannot be achieved at lower energy levels. For power settings > 30 Watts, the recommended irrigation flow rate is 15 ml/min.

** The temperature displayed on the RF generator does not represent tissue temperature or electrode tissue interface temperature.

Additional Recommendation:

For isthmus dependent flutter ablation, power applications exceeding 30 watts and up to 50 watts should only be used if conduction block cannot be achieved at lower power levels.

Sterilization/“Use By” Date

This catheter has been sterilized with ethylene oxide gas. Product and package testing have been conducted to support the “Use By” date printed on the product labels. **DO NOT USE** after the “Use By” date. This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

Storage

Store in a cool, dry, dark place. Storage temperature should be between 5 and 25°C (41 and 77°F).

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How Supplied THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter

- The THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter is supplied STERILE (EtO).
- The catheter has a 7.5 F shaft with 8 F ring electrodes and a usable length of 115 cm.
- Currently, the available curves for the THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheters include DD, FF, JJ, DF, and FJ (Figure 1).
- Additional catheter accessory devices are provided separately.


Figure 1

Packaging of THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter

The THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter is supplied STERILE. The catheter is secured in a two-piece thermoform tray and placed into a Tyvek/Nylon film pouch, sealed, and placed in a box. Both the pouch and thermoform tray are sterile unless the package is damaged or opened.

Directions for Use

THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter

Please refer to the User Manuals for the CARTO 3 System, irrigation pump and RF generator for instructions on connecting and operating these systems in conjunction with the THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter. Use appropriate Biosense Webster accessory cables to connect the THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter to the appropriate accessory equipment.

1. Using aseptic technique, remove the catheter from the package and place in a sterile work area. Inspect the catheter carefully for electrode integrity and overall condition.
2. Create a vascular access in a large central vessel using aseptic techniques.
3. In order to prevent damage to the catheter tip, use the insertion tube supplied with the catheter to advance or retract the catheter through the hemostasis valve of the sheath. After insertion, slide the insertion tube back toward the handle.
4. To verify compatibility between the sheath and catheter, advance the catheter through sheath prior to insertion. Any sheath < 8.5 French is contraindicated.
5. The catheter is connected to the recording equipment and/or the RF generator through the CARTO 3 System using the appropriate interface cables. Connect the catheter to the Patient Interface Unit (PIU) via the appropriate Biosense Webster cable. Connect the PIU to the generator via the appropriate Biosense Webster cable. Connect the PIU to the appropriate recording and mapping systems, including the CARTO 3 Navigation System, with appropriate interface cables. Use only Biosense Webster interface cables. Connect the irrigation pump tubing to the Luer fitting of the catheter. A 3-way stopcock may also be used. Connect the irrigation pump to a room temperature, heparinized (1 IU heparin/ml) normal saline bag using standard safe hospital practices. To complete the electrical circuit, connect an indifferent electrode to the indifferent electrode input on the generator.

6. Flush the catheter and tubing per standard technique to ensure purging of trapped air bubbles and to verify that the irrigation holes are patent.
7. Start continuous irrigation at a flow rate of 2 ml/min.
8. Insert the THERMOCOOL SMARTTOUCH SF Uni-Directional Navigation Catheter via the entrance site using the insertion tube and an appropriately sized introducer sheath. Advance the catheter to the area under investigation.
9. Verify that the "TCool SF" option is selected on the RF generator. When this option is chosen, the generator defaults to the safety parameters established for the THERMOCOOL SMARTTOUCH SF Catheter.
10. Use both fluoroscopy and electrograms (EGM) to aid in proper positioning.
11. In order to achieve optimal force reading accuracy and stability, allow the catheter to warm up for 2 minutes after connection to the CARTO 3 System, prior to use of the force feedback feature.
12. Zero the contact force reading following insertion into the patient. The tip electrode and all four ring electrodes must be outside of the sheath so that the force sensor is inside the body. Ensure the catheter tip is not in contact with tissue by evaluating the location on fluoroscopy and the CARTO 3 System, the EGM amplitude, and catheter movement. Variations in the force reading at the same rate as the cardiac or respiration cycle may indicate contact with cardiac structures. Once these markers indicate the tip is not in contact, the reading can be zeroed. Refer to the User Manual for your CARTO 3 System for instructions on how to zero the contact force reading.
13. Zero the contact force reading when moving the catheter from one chamber of the heart to another or upon re-insertion.
14. The catheter tip can be deflected to facilitate positioning by using the thumbknob to vary tip curvature. Pushing the thumbknob forward causes the catheter tip to deflect; when the thumbknob is pulled back, the tip straightens.
15. For RF application and flow settings refer to Table 1.
16. Recommendation for irrigation **when using Stockert's RF generators**: Increase the irrigation to high flow rate starting up to 5 seconds before the onset of RF energy delivery and maintaining this higher flow rate until 5 seconds after termination of the energy application. For power levels up to 30 W, a high flow rate of 8 ml/min should be used. For power levels between 31-50 W, a high flow rate of 15 ml/min should be used. Do not use catheter without irrigation flow.
17. The application of RF energy must not be initiated until the increase in irrigation flow rate is confirmed by a minimum of 2°C decrease in tip electrode temperature.
18. Monitor the catheter tip temperature throughout the procedure to ensure adequate irrigation. The peak temperature should not exceed 40°C during RF energy delivery. Note: the displayed temperature represents the temperature of the electrode only, not the temperature of the tissue.
19. Start a procedure at 15-20 W. After 15 seconds, power may be increased by 5-10 W increments as needed, until a transmural lesion is achieved, defined by > 80% reduction in unipolar atrial

electrogram amplitude, or emergence of double potentials of equal and low amplitude. It is recommended that power not exceed 50 W when the catheter is parallel to the tissue and 35 W if the catheter is perpendicular to the tissue. The duration of each RF ablation should not exceed 120 seconds. Dragging the catheter to the next location is permissible during the 120-second energy application. RF current may be reapplied to the same or alternate sites using the same catheter.

20. RF current may be reapplied to the same or alternate sites using the same catheter. However, in the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip electrode inspected for coagulum before RF current is reapplied. To remove any coagulum, if present, a sterile gauze pad dampened with sterile saline may be used to gently wipe the tip section clean; do not scrub or twist the tip electrode as damage to the tip electrode bond may occur and loosen the tip electrode or damage may also occur to the contact force sensor and affect measurement accuracy. Prior to reinsertion, ensure that the irrigation holes are not plugged as follows:

- Fill a 1 or 2 ml syringe* with sterile saline and attach to the stopcock or sidearm.

- Carefully inject the saline from the syringe into the catheter. A stream of fluid should be visible from the tip of the catheter.

- Repeat steps a and b, if necessary until the holes are cleared.

- Flush catheter and tubing per standard technique to ensure purging of trapped air bubbles and to verify that the irrigation holes are patent.

- The catheter can now be introduced into the patient.

- Zero catheter following re-insertion into patient.

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Directions for Use

The THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter.

Please refer to the User Manuals for the CARTO 3 System, irrigation pump and RF Generator for instructions on connecting and operating these systems in conjunction with the THERMOCOOL SMARTTOUCH® SF Bi-Directional Navigation Catheter. Use appropriate Biosense Webster accessory cables to connect the THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter to the appropriate accessory equipment.

1. Using aseptic technique, remove the catheter from the package and place in a sterile work area. Inspect the catheter carefully for electrode integrity and overall condition.

2. Create a vascular access in a large central vessel using aseptic techniques.

3. In order to prevent damage to the catheter tip, use the insertion tube supplied with the catheter to advance or retract the catheter through the hemostasis valve of the sheath. After insertion, slide the insertion tube back toward the handle.

4. To verify compatibility between the sheath and catheter, advance the catheter through sheath prior to insertion. Any sheath < 8.5 French is contraindicated.
5. The catheter is connected to the recording equipment and/or the RF generator through the CARTO 3 System using the appropriate interface cables. Connect the catheter to the Patient Interface Unit (PIU) via the appropriate Biosense Webster cable. Connect the PIU to the generator via the appropriate Biosense Webster cable. Connect the PIU to the appropriate recording and mapping systems, including the CARTO 3 Navigation System, with appropriate interface cables. Use only Biosense Webster interface cables. Connect the irrigation pump tubing to the luer fitting of the catheter. A 3-way stopcock may also be used. Connect the irrigation pump to a room temperature, heparinized (1 IU heparin/ml) normal saline bag using standard safe hospital practices. To complete the electrical circuit, connect an indifferent electrode to the indifferent electrode input on the generator.
6. Flush the catheter and tubing per standard technique to ensure purging of trapped air bubbles and to verify that the irrigation holes are patent.
7. Start continuous irrigation at a flow rate of 2 ml/min.
8. Insert the THERMOCOOL SMARTTOUCH SF Bi-Directional Navigation Catheter via the entrance site using the insertion tube and an appropriately sized introducer sheath. Advance the catheter to the area under investigation.
9. Verify that the "TCool SF" option is selected on the RF generator. When this option is chosen, the generator defaults to the safety parameters established for the THERMOCOOL SMARTTOUCH SF Catheter.
10. Use both fluoroscopy and electrograms to aid in proper positioning.
11. In order to achieve optimal force reading accuracy and stability, allow the catheter to warm up for 2 minutes after connection to the CARTO 3 System, prior to use of the force feedback feature.
12. Zero the contact force reading following insertion into the patient. The tip electrode and all four ring electrodes must be outside of the sheath so that the force sensor is inside the body. Ensure the catheter tip is not in contact with tissue by evaluating the location on fluoroscopy and the CARTO 3 System, the EGM amplitude, and catheter movement. Variations in the force reading at the same rate as the cardiac or respiration cycle may indicate contact with cardiac structures. Once these markers indicate the tip is not in contact, the reading can be zeroed. Refer to the User Manual for your CARTO 3 System for instructions on how to zero the contact force reading.
13. Zero the contact force reading when moving the catheter from one chamber of the heart to another or upon re-insertion.
14. Use the Rocker Lever to deflect the catheter tip (Figure 2). When the lever is pulled back from neutral, the tip will deflect relative to the direction of rotation. The amount of deflection is relative to the amount of lever rotation. When the lever is pushed forward, the tip will deflect in the opposite direction. To straighten the tip, return the Rocker Lever to neutral position.



Figure 2

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15. For RF application and flow settings refer to Table 1.

16. Recommendation for irrigation **when using Stockert's RF generators:**

Increase the irrigation to high flow rate starting up to 5 seconds before the onset of RF energy delivery and maintaining this higher flow rate until 5 seconds after termination of the energy application. For power levels up to 30 watts, a high flow rate of 8 ml/min should be used. For power levels between 31-50 watts, a high flow rate of 15 ml/min should be used. Do not use catheter without irrigation flow.

17. The application of RF energy must not be initiated until the increase in irrigation flow rate is confirmed by a minimum of 2°C decrease in tip electrode temperature.

18. Monitor the catheter tip temperature throughout the procedure to ensure adequate irrigation. If temperature increases to 40°C during RF energy delivery, power delivery should be interrupted. The irrigation system must be rechecked prior to restarting RF delivery. Note: the displayed temperature represents the temperature of the electrode only, not the temperature of the tissue.

19. Start a procedure at 15-20 Watts. After 15 seconds, power may be increased by 5-10 W increments as needed, until a transmural lesion is achieved, defined by > 80% reduction in unipolar atrial electrogram amplitude, or emergence of double potentials of equal and low amplitude. It is recommended that power not exceed 50 W when the catheter is parallel to the tissue and 35 W if the catheter is perpendicular to the tissue. The duration of each RF ablation should not exceed 120 seconds. Dragging the catheter to the next location is permissible during the 120-second energy application. RF current may be reapplied to the same or alternate sites using the same catheter.

20. RF current may be reapplied to the same or alternate sites using the same catheter. However, in the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip electrode inspected for coagulum before RF current is reapplied. To remove any coagulum, if present, a sterile gauze pad dampened with sterile saline may be used to gently wipe the tip section clean; do not scrub or twist the tip electrode as damage to the tip electrode bond may occur and loosen the tip electrode or damage may also occur to the contact force sensor and affect measurement accuracy. Prior to reinsertion, ensure that the irrigation holes are not plugged as follows:

If irrigation hole occlusion occurs:

- a) Fill a 1 or 2 ml syringe* with sterile saline and attach to the stopcock or sidearm.
- b) Carefully inject the saline from the syringe into the catheter. A stream of fluid should be visible from the tip of the catheter.
- c) Repeat steps a and b, if necessary until the holes are cleared.
- d) Flush catheter and tubing per standard technique to ensure purging of trapped air bubbles and to verify that the irrigation holes are patent.
- e) The catheter can now be introduced into the patient.
- f) Zero catheter following re-insertion into patient.

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Adverse Reactions

A number of serious adverse reactions have been documented for catheter ablation procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death. The following complications were also noted to have occurred during prior studies or have been reported in the literature:

- Catheterization/catheter procedure related: vascular bleeding/local hematomas, thrombosis, arteriovenous (AV) fistula, pseudoaneurysm, thromboembolism and vasovagal reactions, cardiac perforation, pericardial effusion/tamponade, thrombi, air embolism, arrhythmias and valvular damage, pneumothorax and hemothorax, pulmonary edema, hypoxia, pleural effusion, acute respiratory distress syndrome (ARDS), congestive heart failure, aspiration pneumonia, pneumonia, asthmatic attack, hypotension, implantable cardioverter-defibrillator (ICD) lead malfunction, anemia, thrombocytopenia, disseminated intravascular coagulation, epistaxis, systemic infection, urinary tract infection, sedation induced apnea, sedation induced CO₂ retention with lethargy and cholecystitis.
- RF related: chest pain/discomfort, ventricular tachyarrhythmia, transient ischemic attack (TIA), cerebrovascular accident (CVA), complete heart block, coronary artery spasm, coronary artery thrombosis, coronary artery dissection, cardiac thromboembolism, pericarditis, cardiac perforation/tamponade, valvular damage and increased phosphokinase level.
- Unrelated to device or procedure: urinary retention, transient extremity numbness, Parkinson's disease and gastrointestinal diverticulosis.

DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.

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2. Information of manufacturers of the medical device (name, location address, mailing address, phone number, address of place of production)

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4. Technical specifications of the medical device

4.1. General description

The Biosense Webster THERMOCOOL SMARTTOUCH SF Catheters are a steerable, multi-electrode, luminal catheters with a deflectable tip designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) current to the catheter tip dome electrode for ablation purposes. The catheter shaft measures 7.5 F with 8 F ring electrodes. For ablation, the catheter is used in conjunction with an RF generator and a dispersive pad (indifferent electrode). The catheter has force sensing technology that provides a real-time measurement of contact force between the catheter tip and the heart wall.

The catheter has a high-torque shaft, with either a uni-directional or bi-directional braided deflectable tip section, containing an array of ring electrodes. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site. All ring electrodes are manufactured from platinum-iridium. The tip dome electrode is manufactured from palladium-platinum and serves to deliver RF current from the RF generator to the desired ablation site. The catheter incorporates a temperature sensor that is embedded in the 3.5 mm tip dome electrode. The 3.5 mm tip dome features a 56-hole irrigation pattern as shown in Figure 1.

Figure 1: D-1347-XX-S and D-1348-XX-S tip section

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This catheter features a location sensor embedded in the tip section that transmits location and contact force information to the CARTO 3 System. An appropriate reference device is required for location reference position purposes. The catheter interfaces with standard recording equipment and a compatible RF generator via accessory extension cables with the appropriate connectors.

THERMOCOOL SMARTTOUCH SF Uni-Directional Catheter


Figure 2 Photo of THERMOCOOL SMARTTOUCH SF Uni-Directional Catheter

The Uni-directional catheter's tip deflection is shown in Figure 3 and is controlled at the proximal end by a hand piece in which a piston travels forward and back, controlled by a thumbknob. When the thumbknob is pushed forward, the tip is deflected (curved). When the Thumbknob is pulled back, the tip straightens. The high-torque shaft also allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site. The uni-directional catheter provides 180° of deflection (minimum) and is offered in D, F and J curves.


Figure 3 Uni-Directional Curve Profiles

An overview of the THERMOCOOL SMARTTOUCH SF Uni-Directional Catheter is provided in Figure 4.


Figure 4 Overview of THE THERMOCOOL SMARTTOUCH SF Catheter with uni-directional tip deflection.

THERMOCOOL SMARTTOUCH SF Bi-Directional Catheter


Figure 5 Photo of THERMOCOOL SMARTTOUCH SF Bi-Directional Catheter

The Bi-Directional Catheter's tip deflection is shown in Figure 6 and is controlled by a rocker lever. The high-torque shaft also allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site. Additionally, a variety of curve types are available in symmetrical or asymmetrical combinations, providing two 180° (minimum) opposed, single planed curves. The bi-directional catheter is offered in DD, FF, JJ, DF and FJ curves.


Figure 6 Bi-Directional Curve Profiles

An overview of the THERMOCOOL SMARTTOUCH SF Bi-Directional Catheter is provided in Figure 7


Figure 7 Overview of the THERMOCOOL SMARTTOUCH SF Catheter with bidirectional tip deflection.

4.2. Versions / Options of MD

1. ThermoCool SmartTouch SF Uni-Directional Catheter:
 - Curve type D; electrodes configuration 1-6-2 mm;
 - Curve type F; electrodes configuration 1-6-2 mm;
 - Curve type J; electrodes configuration 1-6-2 mm.
2. ThermoCool SmartTouch SF Bi-Directional Catheter:
 - Curve type D-D; electrodes configuration 1-6-2 mm;
 - Curve type D-F; electrodes configuration 1-6-2 mm;
 - Curve type F-F; electrodes configuration 1-6-2 mm;
 - Curve type F-J; electrodes configuration 1-6-2 mm;
 - Curve type J-J; electrodes configuration 1-6-2 mm.

Table 3 — Options of MD

<u>Catheter</u>	<u>Product code</u>	<u>Configuration</u>	<u>Manufacturing Part Number</u>
Thermocool smarttouch SF Uni-Directional Catheter	D134701	Curve type D; electrodes configuration 1-6-2 mm	D-1347-XX-S
	D134702	Curve type F; electrodes configuration 1-6-2 mm	
	D134703	Curve type J; electrodes configuration 1-6-2 mm	
	D134801	Curve type D-D; electrodes configuration 1-6-2 mm	
	D134802	Curve type F-F; electrodes configuration 1-6-2 mm	

Thermocool smarttouch SF Bi-Directional Catheter	D134803	Curve type J-J; electrodes configuration 1-6-2 mm	D-1348-XX-S
	D134804	Curve type F-J; electrodes configuration 1-6-2 mm	
	D134805	Curve type D-F; electrodes configuration 1-6-2 mm	

a. The possibility and ways of integrating with other MD

The THERMOCOOL SMARTTOUCH SF Catheter is connected to the CARTO 3 System in order to conduct an Electrophysiology procedure. Table below lists devices which are manufactured or distributed by Biosense Webster and also used during a procedure:

Table 4 — Devices used during a procedure

Component Name	Part Number
CARTO 3 System with SMARTTOUCH software module	FG-5400-00
Interface Cable	D-1286-01-S D-1286-03-S
Stockert RF Generator, or smartablate RF Generator	M-5463-01 M-4900-XX
COOLFLOW Irrigation Pump and COOLFLOW tubing set Or Smartablate Irrigation Pump and Smartablate Irrigation Tubing	M-5491-02 D-1233-01-S M-4900-XX D-1320-01-S
CARTO External Reference Patches	D-1283-02
PREFACE Guiding Sheath (or equivalent)	301803A 301803M 301803MS 301803P 301805M

Commercially available devices, not manufactured by Biosense Webster, are also required to conduct an electrophysiology procedure. These include, but are not limited to, the following list of equipment

- Grounding Pad (indifferent electrode)
- Pacing Stimulator
- Electrophysiology (EP) recording equipment
- Electrocardiogram (ECG) leads
- Fluoroscopy/X-ray System
- Cardiac Defibrillator
- Intracardiac Ultrasound (investigator preference)

A connectivity diagram is depicted in Figure 8


Figure 8 Connectivity Diagram

b. The materials from which it is made or its medical device basic parts, biocompatibility

Table 5 — Materials

Description	Material
Materials of following catheters: - Thermocool smarttouch ST SF NAV Uni-Directional Catheter - Thermocool smarttouch ST SF NAV Bi-Directional Catheter	
Tip/ Shaft Adhesive	Polyurethane (PU) P-9231
Peek Housing	PEEK Tube PEEK HOUSING, 8FR, TRITON Sub : PEEK TUBING,.082IDX.092OD
Irrigation Pathway	
Luer Hub	Polycarbonate Plastic Polycarbonate Resin Blue or White - HP3NREU-8H9D273
Irrigation Tubing Material	Polyimide 75 D and Pebax 2553-SA-01
Tip	
Outer Coat Tip	Pallethane 80A
Core	55D/80A Pallethane
Biosensor Housing	Nitinol SE508ELI/SE508/SE506/SE510
Tip Electrode	Palladium/ Platinum 90%(89.6%min) PT/10% (±.25%) IR WIRE BAR

Ring Electrode		Platinum 90%/ Iridium 10% Nobel Met ASTM:B684-97 P/TR 8007 contain P/TR 8007-01	
ACL Ring Electrode		Platinum 80%/ Irridium 20% - 90%(89.6%min) PT/10% (±.25%) IR WIRE BAR	
Temperature Sensor (Thermocouple)		Polyimide Sleeve 75D, Epoxy Encapsulant Constantan/ Cooper (THERMAL CLASSIFICATION 155C)	
Shaft			
Shaft		Pebax2553-SA-01, Vestamid (75D/55D/75D)	
Braiding		Stainless Steel 304v	
Insulation		SPN (SINGLE POLY NYLON) (THERMAL CLASSIFICATION 155C)	
Puller Wire		Stainless Steel 304v	
Handle			
Thermocool smarttouch ST SF NAV Uni- Directional Catheter		Thermocool smarttouch ST SF NAV Bi- Directional Catheter	
Description	Material	Description	Material
Handle	Delrin/ Acetal ASTM-D6778 POM 0211	Rocker Lever	Deflection Knob material: Nylon Zytel®42A NC010 w/ 30% glass and 15% PTFE filler
Barrel Extension	Delrin/ Acetal ASTM-D6778 POM 0211	Handle	Polycarbonate LEXAN™ Resin HP3NREU w/ 30% glass and 15% PTFE filler
Piston	Delrin/ Acetal ASTM-D6778 POM 0211	Connector to handle bond joint	Polyurethane (PU) P- 9231 Polycarbonate LEXAN™ Resin HP3NREU w/30% Glass and 5% PTFE filler
Connector to handle bond joint	Polyurethane (PU) P-9231		

Table 6 — Biocompatibility Testing Results for the Platium / palladium alloy contained within the tip electrode of the catheter

Test	Result
Hemolysis Test (ASTM (F756) Method – Direct Contact	Pass
MEM Elution Using L-929 Mouse Fibroblast Cells (ISO) (Cytotoxicity)	Pass
Intracutaneous Irritation Test (ISO)	Pass

Acute Systemic Injection Test (ISO)	Pass
Materials Mediated Rabbit Pyrogen	Pass
Guinea Pig Maximization Sensitization Test	Pass
Complement Activation C3a Assay and SC5b-9 Assay	Pass
Partial Thromboplastin Time (PTT)	Pass
Platelet and Leukocyte Counts	Pass
Thrombosis (In vivo) – 2 Dog	Pass

3. Main parameters and characteristics of the medical device

Table 7 — Main parameters

1. Thermocool smarttouch SF Uni-Directional Catheter			
Characteristics	D134701	D134702	D134703
General			
Usable Length, cm	115 ± 3	115 ± 3	115 ± 3
Mass, g	76.0 +/-0.4	76.0 +/-0.4	76.0 +/-0.4
Outer diameter, mm	0 to 2.794	0 to 2.794	0 to 2.794
Puller wire material length	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2
Application Time, sec	Not exceed 120	Not exceed 120	Not exceed 120
Power Settings During RF Application, W	15 to 50	15 to 50	15 to 50
Flow Rate, ml/min	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)
	15 ml/min-5% +15% (> 30 W)	15 ml/min-5% +15% (> 30 W)	15 ml/min-5% +15% (> 30 W)
Pre-curve °	≤ 15°	≤ 15°	≤ 15°
Full length, cm*	146 ± 3	146 ± 3	146 ± 3
Mass, g	76.0 +/-0.4	76.0 +/-0.4	76.0 +/-0.4
Luer fitting			
Inner diameter, mm	4.2926 +/- 0.02286	4.2926 +/- 0.02286	4.2926 +/- 0.02286

Length, mm	20.066 +/- 0.25	20.066 +/- 0.25	20.066 +/- 0.25
Major Outer diameter, mm	7.747 +0.000 -0.076	7.747 +0.000 -0.076	7.747 +0.000 -0.076
Side arm Cord length, mm	215,9 +/- 3,175	215,9 +/- 3,175	215,9 +/- 3,175
Tip Specifications			
Tip electrode length, mm	3,6703 +/- 0,5085	3,6703 +/- 0,5083,5	3,6703 +/- 0,5083,5
Number of Irrigation Holes	56	56	56
Irrigation Hole Diameter, mm	0,0889 +/- 0.0127	0,0889 +/- 0.0127	0,0889 +/- 0.0127
Number of Lumens in the Deflectable Tip	4	4	4
Deflectable Tip Lumen Dimensions, mm	2 x 0,4826 +/- 0.0254 2 x 0,889 +/- 0.0254	2 x 0,4826 +/- 0.0254 2 x 0,889 +/- 0.0254	2 x 0,4826 +/- 0.0254 2 x 0,889 +/- 0.0254
Braiding	Braided: 16 strands, 45 +/-1picks per 2,54 cm (1")	Braided: 16 strands, 45 +/-1picks per 2,54 cm (1")	Braided: 16 strands, 45 +/-1picks per 2,54 cm (1")
Tip Electrode Length, mm	3,6703 ± 0.0508	3, 6703 ± 0.0508	3, 6703 ± 0.0508
Number of Ring Electrodes	3	3	3
Ring Electrode Spacing, mm	1- 6- 2	1- 6- 2	1- 6- 2
DC Resistance, ohms	0 to 10	0 to 10	0 to 10
DC Isolation, µA	0 to less than 1	0 to less than 1	0 to less than 1
RF Impedance, ohms	0 to less than 10	0 to less than 10	0 to less than 10
	0 to less than 10	0 to less than 10	0 to less than 10
Isolation and Phase Angle Between Circuits	-86° to -90°	-86° to -90°	-86° to -90°

RF Leakage Current, mA	0 to 252	0 to 252	0 to 252
Puller Wire Anchor Location	Anchored at the distal deflectable tip	Anchored at the distal deflectable tip	Anchored at the distal deflectable tip
Inactive puller wire outside the polymide tube, mm	$1,0 \pm 0,5$	$1,0 \pm 0,5$	$1,0 \pm 0,5$
Curve Types	D	F	J
Deflection Angle Range	0-150°	0-150°	0-150°
Tip Stiffness, N	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647
Buckle Force, N	0 to less than 0,5648	0 to less than 0,5648	0 to less than 0,5648
Side Force, N	Greater than 0,03922 to infinity	Greater than 0,03922 to infinity	Greater than 0,03922 to infinity
Pull Entire Catheter, N	Greater than 17,7929 to infinity	Greater than 17, 7929 to infinity	Greater than 17, 7929 to infinity
Irrigation Pathway Specifications			
OID = 0,6858 +/- 0.0254	OID = 0,6858 +/- 0.0254	OID = 0,6858 +/- 0.0254	OID = 0,6858 +/- 0.0254
Shaft Specifications			
Shaft Waviness, mm	0 to 20	0 to 20	0 to 20
Internal Lead Conductors	6x 40 AWG Copper Wire	6x 40 AWG Copper Wire	6x 40 AWG Copper Wire
Temperature Sensor Conductor	Thermo-couple: 46 AWG Constantan Wire • 46AWG Copper Wire	Thermo-couple: 46 AWG Constantan Wire • 46AWG Copper Wire	Thermo-couple: 46 AWG Constantan Wire • 46AWG Copper Wire
Tip Conductor/Lead Wire Joints	6 twisted pair of 44 AWG copper (1 pair external to cable shielding)	6 twisted pair of 44 AWG copper (1 pair external to cable shielding)	6 twisted pair of 44 AWG copper (1 pair external to cable shielding)
Compression Coil	0.2286 ID x 0.4318 OD	0.2286 ID x 0.4318 OD	0.2286 ID x 0.4318 OD

Handle Specifications			
Construction	Handle with adjustable deflection knob, tension knob and luer fitting	Handle with adjustable deflection knob, tension knob and luer fitting	Handle with adjustable deflection knob, tension knob and luer fitting
Length, cm	21.03+/- 0.05	21.03+/- 0.05	21.03+/- 0.05
Length of the piston, cm	6.985 +/- 0.0762	6.985 +/- 0.0762	6.985 +/- 0.0762
Connector	34 pin	34 pin	34 pin
Handle Temperature, °C	0 to 48	0 to 48	0 to 48

2. ThermoCool SmartTouch SF Bi-Directional Catheter

Characteristics	D134801	D134802	D134803	D134804	D134805
General					
Usable Length, cm	115 ± 3	115 ± 3	115 ± 3	115 ± 3	115 ± 3
Mass, g	89,8± 0.7	89,8± 0.7	89,8± 0.7	89,8± 0.7	89,8± 0.7
Outer diameter, mm	2,7 ± 0,1	2,7 ± 0,1	2,7 ± 0,1	2,7 ± 0,1	2,7 ± 0,1
Puller Wire Material length, mm	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2	1524 ± 76,2
Application Time, sec	Not exceed 120	Not exceed 120	Not exceed 120	Not exceed 120	Not exceed 120
Power Settings During RF Application, W	15 to 50	15 to 50	15 to 50	15 to 50	15 to 50
Flow Rate, ml/min	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)	8 ml/min-5% +15% (≤ 30 W)
	15 ml/min-5% +15%	15 ml/min-5% +15%	15 ml/min-5% +15%	15 ml/min-5% +15%	15 ml/min-5% +15%

	(> 30 W)	(> 30 W)	(> 30 W)	(> 30 W)	(> 30 W)
Pre-curve, °	0 to 15°	0 to 15°	0 to 15°	0 to 15°	0 to 15°
Full length, cm *	146 ± 3	146 ± 3	146 ± 3	146 ± 3	146 ± 3
Mass, g	89.8 ± 0.7	89.8 ± 0.7	89.8 ± 0.7	89.8 ± 0.7	89.8 ± 0.7
Luer fitting					
Inner diameter, mm	4.2926 ± 0.02286	4.2926 ± 0.02286	4.2926 ± 0.02286	4.2926 ± 0.02286	4.2926 ± 0.02286
Major Outer diameter, mm	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762	7.747 +0.000 -0.0762
Length, mm	20.066 ± 0.254	20.066 ± 0.254	20.066 ± 0.254	20.066 ± 0.254	20.066 ± 0.254
Side arm Cord length, mm	76.2 ± 2.54	76.2 ± 2.54	76.2 ± 2.54	76.2 ± 2.54	76.2 ± 2.54
Tip Specifications					
Tip electrode length, mm	3,6703 ± 0,5083,5	3,6703 ± 0,5083,5	3,6703 ± 0,5083,5	3,6703 ± 0,5083,5	3,6703 ± 0,5083,5
Number of Irrigation Holes	56	56	56	56	56
Irrigation Hole Size, mm	0,0889 ± 0.0127	0,0889 ± 0.0127	0,0889 ± 0.0127	0,0889 ± 0.0127	0,0889 ± 0.0127
Number of Lumens in the Deflectable Tip	4	4	4	4	4
Deflectable Tip Lumen Dimensions, mm	2 x 0,4826±/- 0.0254 2 x 0,889±/- 0.0254	2 x 0,4826±/- 0.0254 2 x 0,889±/- 0.0254	2 x 0,4826±/- 0.0254 2 x 0,889±/- 0.0254	2 x 0,4826±/- 0.0254 2 x 0,889±/- 0.0254	2 x 0,4826±/- 0.0254 2 x 0,889±/- 0.0254
Braiding	Braided: 16 strands, 45 +/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45 +/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45 +/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45 +/-1 picks per 2,54 cm (1")	Braided: 16 strands, 45 +/-1 picks per 2,54 cm (1")
Tip Electrode	3,6703 ±	3,6703 ±	3,6703 ±	3,6703 ±	3,6703 ±

Length, mm	0.0508	0.0508	0.0508	0.0508	0.0508
Number of Ring Electrodes	3	3	3	3	3
Ring Electrode Spacing, mm	1- 6- 2	1- 6- 2	1- 6- 2	1- 6- 2	1- 6- 2
DC Resistance, ohms	0 to 10	0 to 10	0 to 10	0 to 10	0 to 10
DC Isolation, μ A	0 to less than 1	0 to less than 1	0 to less than 1	0 to less than 1	0 to less than 1
RF Impedance, ohms	0 to less than 10	0 to less than 10	0 to less than 10	0 to less than 10	0 to less than 10
	0 to less than 25	0 to less than 25	0 to less than 25	0 to less than 25	0 to less than 25
Isolation and Phase Angle Between Circuits	-86° to -90°	-86° to -90°	-86° to -90°	-86° to -90°	-86° to -90°
RF Leakage Current, mA	0 to 252	0 to 252	0 to 252	0 to 252	0 to 252
Puller Wire Anchor Location	Anchored at the distal deflectable tip	Anchored at the distal deflectable tip	Anchored at the distal deflectable tip	Anchored at the distal deflectable tip	Anchored at the distal deflectable tip
Inactive puller wire outside the polyimide tube, mm	1,0 $\pm$ 0,5	1,0 $\pm$ 0,5	1,0 $\pm$ 0,5	1,0 $\pm$ 0,5	1,0 $\pm$ 0,5
Curve Types	DD	FF	JJ	FJ	DF
Deflection Angle Range	0-180°	0-180°	0-180°	0-180°	0-180°
Tip Stiffness, N	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647	0 to less than 0,2647
Buckle Force, N	0 to less than 0,5648	0 to less than 0,5648	0 to less than 0,5648	0 to less than 0,5648	0 to less than 0,5648
Side Force, N	Greater than	Greater than	Greater than	Greater than	Greater than

	0,03922 to infinity	0,03922 to infinity	0,03922 to infinity	0,03922 to infinity	0,03922 to infinity
Pull Entire Catheter, N	Greater than 17,7929 to infinity	Greater than 17, 7929 to infinity	Greater than 17, 7929 to infinity	Greater than 17, 7929 to infinity	Greater than 17, 7929 to infinity
Irrigation Pathway Specifications					
Injection Pressure, kPa	Greater than 552 to infinity	Greater than 552 to infinity	Greater than 552 to infinity	Greater than 552 to infinity	Greater than 552 to infinity
Irrigation Tubing Dimensions, mm	ID = 0,6858+/- 0.0127 ID = 0,4064+/- 0.0127	ID = 0,6858+/- 0.0127 ID = 0,4064+/- 0.0127	ID = 0,6858+/- 0.0127 ID = 0,4064+/- 0.0127	ID = 0,6858+/- 0.0127 ID = 0,4064+/- 0.0127	ID = 0,6858+/- 0.0127 ID = 0,4064+/- 0.0127
Shaft Specifications					
Shaft Waviness, mm	0 to 20	0 to 20	0 to 20	0 to 20	0 to 20
Internal Lead Conductors	6x 40 AWG Copper Wire	6x 40 AWG Copper Wire	6x 40 AWG Copper Wire	6x 40 AWG Copper Wire	6x 40 AWG Copper Wire
Temperature Sensor Conductor	Thermo-couple: 46 AWG Constantan Wire • 46AWG Copper Wire	Thermo-couple: 46 AWG Constantan Wire • 46AWG Copper Wire	Thermo-couple: 46 AWG Constantan Wire • 46AWG Copper Wire	Thermo-couple: 46 AWG Constantan Wire • 46AWG Copper Wire	Thermo-couple: 46 AWG Constantan Wire • 46AWG Copper Wire
Tip Conductor/ Lead Wire Joints	6 twisted pair of 44 AWG copper (1 pair external to cable shielding)	6 twisted pair of 44 AWG copper (1 pair external to cable shielding)	6 twisted pair of 44 AWG copper (1 pair external to cable shielding)	6 twisted pair of 44 AWG copper (1 pair external to cable shielding)	6 twisted pair of 44 AWG copper (1 pair external to cable shielding)
Compression Coil	0.2286 ID x 0.4318 OD	0.2286 ID x 0.4318 OD	0.2286 ID x 0.4318 OD	0.2286 ID x 0.4318 OD	0.2286 ID x 0.4318 OD
Handle Specifications					
Construction	Handle with	Handle with	Handle with adjustable	Handle with	Handle with

	adjustable deflection knob, tension knob and luer fitting	adjustable deflection knob, tension knob and luer fitting	deflection knob, tension knob and luer fitting	adjustable deflection knob, tension knob and luer fitting	adjustable deflection knob, tension knob and luer fitting
Length, cm	21.03+/- 0.05	21.03+/- 0.05	21.03+/- 0.05	21.03+/- 0.05	21.03+/- 0.05
Connector	34 pin	34 pin	34 pin	34 pin	34 pin
Handle Temperature, °C	0 to 48	0 to 48	0 to 48	0 to 48	0 to 48
* Tolerance for these parameters were established using Russian GOST R 50444-92					

Information of EMC tests (EMC Test Report)

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The CARTO® 3 System is intended for use in the electromagnetic environment specified below. The customer or the user of the CARTO® 3 System should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The CARTO® 3 System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. See CARTO® 3 System Interference with Other Medical Equipment.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker	Compliance	The CARTO® 3 System is suitable for use in all establishments other than domestic.

The CARTO® 3 System is intended for use in the electromagnetic environment specified below. The customer or the user of the CARTO® 3 System should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
emissions IEC 61000-3-3		

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Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient or burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output	±2 kV for power supply lines ±1 kV for input/output	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CARTO [®] 3 System requires continued operation during power mains interruptions, it is recommended that the CARTO [®] 3 System be powered from an uninterruptible power supply (UPS).
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

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Table 9 — Recommended separation distances between portable and mobile RF communications equipment and system with Cardiac Catheters

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- These guidelines might not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

5. Shelf life

The Thermocool SmartTouch SF Uni-Directional and Thermocool SmartTouch SF Bi-Directional Catheters are validated for a one-year product shelf life. Catheter integrity and performance testing was conducted after sterilization by ethylene oxide (EtO) through three commercial cycles followed by accelerated aging (equivalent to one-year real time aging), thermal cycling and simulated transportation conditioning. Acceptance criteria were met for all tests.

6. Packaging

The THERMOCOOL SMARTTOUCH SF Uni-Directional Catheters are packaged in a mounting card. A shaft protector is taped together with labels to form a catheter support that is not part of the sterile barrier. The catheter is loaded into the mounting card and a shaft protector. The loaded mounting card is placed in a Tyvek/ nylon pouch, heat-sealed forming a single sterile barrier. The sealed pouch is then placed into a unit carton with an IFU insert.

Figure 9 Photo of primary packaging

Figure 10 Photo of secondary packaging

The THERMOCOOL SMARTTOUCH SF Bi-Directional Catheters packaged in a two-piece thermoform tray with a handle and shaft component. The two-piece thermoform tray is snapped together to form a catheter support tray that is not part of the sterile barrier. The catheter is loaded into the tray. The loaded tray is placed in a Tyvek/nylon pouch, heat-sealed forming a single sterile barrier. The sealed pouch is then placed into a unit carton with an IFU insert.


Figure 11 Photo of primary packaging

Figure 12 Photo of secondary packaging

Table 10 — Packaging materials

The THERMOCOOL SMARTTOUCH SF Uni-Directional Catheter	
Catheter Shaft PETG thermocool Tray	Pebax, Vestamid Pebax 7FR; VESTAMID,CMPND W/20%BS-80% 75D
Nylon/ Tyvek Pouches	Tyvek coated with CR27 adhesive

	Sealed to perfecseal 0,15 mm (6 mil) ULDPE Nylon. (65.0" x 6.0") 1651mm x 152,4 mm
Unit Carton	Carton, (4 5/32" x 1 3/8" x 59 5/8") 105,57 mm x 34,9 mm x 1514,5 mm
Mounting Card Components	Mounting card
The THERMOCOOL SMARTTOUCH SF Bi-Directional Catheter	
Catheter Shaft Thermoform Tray	7.62 mm (0.30") PETG with Glidex
Handle Tray	Bi- Directional/qwikstar Handle
Nylon Pouches	1073B Tyvek coated with CR27 adhesive sealed to perfecseal 0,15 mm (6 mil) ULDPE Nylon. (65.0" x 6.0") 1651mm x 152,4 mm
Foam block	(1.5" x .8" x .8") 38,1mm x 20,32 mm x 20,32 mm
Unit Carton	(4" x 1.5" x 60") 101.6 mm x 38,1 mm x 1524 mm

7. Information on labelling

Figure 13 Example of labelling for The THERMOCOOL SMARTTOUCH SF Uni-Directional Catheter D134701 and The THERMOCOOL SMARTTOUCH SF Bi-Directional Catheter D134801.

Table 11 — List of symbols

Symbol	Definition	Symbol	Definition
	Sterilized Using Ethylene Oxide		Keep Dry
	Do Not Re-Use		Use-By Date
	Caution		Batch Code
	Consult Instructions for Use		Catalog Number
	Do not use if package is damaged		Contents: 1
	Do not use if package is opened.		Manufacturer
	Keep Away from Sunlight		Authorized Representative in the European Community
	Date of Manufacture		Temperature Limit
	Electrodes		Spacing
	Curve Type.		Thermocool smarttouch SF Catheter
	Navigational Catheter		Compatible with Carto 3 EP Navigation System
	Compatible with Carto XP EP Navigation System		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Pin Connector		

8. Conditions of transportation and storage

Table 12 — Storage

Environment Settings	Values
Air temperature, °C	(+ 5- + 25) ± 2 °C
Air humidity (relative), non-condensing, %	(20-60%) $\pm 2\%$
Atmosphere pressure, hpa	101 ± 1 hPa

Table 13 — Transportation

Environment Settings	Values
Air temperature, °C	(+5 - +25) ± 2 °C
Air humidity (relative), non-condensing, %	(20-60%) $\pm 2\%$
Atmosphere pressure, hpa	101 ± 1 hPa

9. The protection of the environment

a. Requirements for the protection of the environment in the application of medical device

Using Thermocool Smarttouch SF Uni-Directional and Thermocool Smarttouch SF Bi-Directional Catheters do not imply compliance with specific requirements for environmental protection. Thermocool Smarttouch SF Uni-Directional and Thermocool Smarttouch SF Bi-Directional Catheters are EU RoHS (2011/65/EU) and REACH (1907/2006) compliant, are free of latex and latex handling, and are phthalate-free.

b. Utilization method/ Disposal

For disposing of Thermocool Smarttouch SF Uni-Directional and Thermocool Smarttouch SF Bi-Directional Catheters required to provide the mandatory implementation of all the necessary procedures in accordance with local applicable laws and regulations.

Disposal or destruction of the Thermocool Smarttouch SF Uni-Directional and Thermocool Smarttouch SF Bi-Directional Catheters shall be in accordance with the classification (Class A – for catheters with damaged packaging, class B – for used devices), rules of collection, use, decontamination, placement, storage, transportation, accounting and utilization of medical waste, established by the authorized federal executive body.

10. Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced)

All Biosense Webster catheters are ethylene oxide (eto) sterilized by Steris Isomedix Inc., a subcontractor for Biosense Webster. The Steris sterilization process utilizes ethylene oxide and incorporates all phases of sterilization (preconditioning, sterilization and aeration) within the sterilization chamber. This process is referred to as EO Express since the use of traditional preconditioning and aeration rooms are not necessary. The validation of Biosense Webster products resulted in one common cycle for all Biosense Webster device families. Sterility assurance was verified to exceed a sterility assurance level (SAL) of 1×10^{-6} . The validation met

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This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

Apostille

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1. District of Columbia, United States of America
2. This public document has been signed by KAMAL S AHMAD
3. acting in the capacity of NOTARY PUBLIC IN AND FOR THE DISTRICT OF COLUMBIA
KAMAL S. AHMAD, NOTARY PUBLIC IN AND FOR THE
DISTRICT OF COLUMBIA
4. bears the seal/stamp of _____

CERTIFIED

5. at Washington, D.C.
6. the 6, NOVEMBER 2017 day of _____
7. by Secretary of the District of Columbia
8. No. 458156
9. Seal/Stamp
10. Signature:

Karen C. Vaughan
Secretary of the District of Columbia

To Federal Service on Surveillance in Healthcare
«Roszdravnadzor»

Place of issue Irvine, CA

date of issue 11-2-2017


CONFIDENTIAL

We Biosense Webster, Inc., located at: 33 Technology Drive Irvine, California 92618 USA, hereby informs about clarification of shelf life for medical device «The cardiac catheter for electrophysiology diagnostics ThermoCool», namely:

- 3-year for Catheter ThermoCool SF NAV Uni-directional, Catheter ThermoCool SF NAV bi-directional, Catheter ThermoCool SF unidirectional with curve visualization NAV, Catheter ThermoCool SF bi-directional with curve visualization NAV;

We certify authenticity of the attached document: Addition to Instruction for Use: The cardiac catheter for electrophysiology diagnostics: ThermoCool SF NAV uni-directional and bi-directional catheters (Registration dossier № RD-18771/44153 from 01.08.2017 for the medical product «The cardiac catheter for electrophysiology diagnostics Thermocool»)

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Name JOANNE WHITEMAN
TITLE Manager, Regulatory Affairs

11-2-2017

Date

I Certify That This Is An
Original Document

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State of California)

County of Orange)

On November 2nd, 2017 before me Steve Oberlander, NOTARY PUBLIC, personally appeared Joanne Whiteman, who proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that she executed the same in her authorized capacity and that by her signature on the instrument the person, or the entity upon behalf of which the person acted, executed the instrument.

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