Siemens Healthcare GmbH,
Henkestr. 127
91052 Erlangen,
Germany

Addendum to the User Manual for Russia

Mobile X-ray Surgical Unit in variants with
accessories: Arcadis Orbic, Arcadis Orbic
3D, Arcadis Varic, Arcadis Avantic

This document is valid for the following product: ARCADIS Orbic

3.

Contr

For 1

con-

for u

medi

Prior

and a

and c

Possibility

For r

catar

nessi

possi

Risks

Risks

neces

needs

when

3.

The d

The d...

and i

Curre

Ambi

Tem

Rel.

Rare

Baro

- Visualization of implants;
- Visualization of surgical instruments.

3.2. Contraindications and Possible side effects:

Contraindications

For this product there are currently no known contra-indications. All residual risks are communicated to the user by warnings and cautions in the operator manual. The final decision for use of the medical device in a certain application is made by the physician based on his/her medical knowledge and the risks of the case involved.

Prior to the examination, it must be confirmed by a physician that an examination is permissible and checked whether increased precautions are necessary.

Possible side effects

For radiation emitting products known side effects of X-ray applications are erythematosis, cataracts, permanent epilation and delayed skin necrosis. The literature considers further possible side effects of radiation applied to humans.

Risks concerning increased risk of cancer are content of the physicians' medical education. The necessity of an X-ray examination in relation to these risks must be considered by the physician when prescribing such a type of examination.



3.3. Conditions for use:

The digital ARCADIS Orbic X-ray system is a mobile system that is designed for use in the ER, OR, and in radiology and endoscopy departments of hospitals, clinics, and outpatient practices. Current country-specific requirements must be complied with.

Ambient conditions for operation:

Temperature range	+15 °C to +35 °C
Rel. humidity	35% to 75% (non condensation)
Barometric pressure	700 hPa to 1060 hPa

3.4. System's overview


The ARCADIS Orbic consists of a C-arm system and a monitor trolley

(1) C-arm system with a 23 cm image intensifier and single-focus tube with generator

(2) Monitor trolley with keyboard, mouse, USB port, two rotatable TFT displays, DVD drive and memory for 60,000 images


- (1) A
- (2) C
- (3) L
- (4) H
- (5) C
- (6) I
- (7) S

C-arm system

- (1) ARCADIS electronic unit
- (2) Control and display panel
- (3) Lifting column
- (4) Horizontal support arm
- (5) C-arm
- (6) Image intensifier with integrated TV camera
- (7) Single tank with X-ray tube unit and integrated collimator

Monitor trolley

- (1) Radiation indicator
- (2) Monitors rotatable by 180°
Optional: folding monitors
- (3) Control panel with storage drawers, keyboard and mouse
- (4) Grips: for neat cable storage
- (5) Drawer for storage of CDs/DVDs, Quick Guide etc.
- (6) DVD R/W drive
- (7) Wheels with cable deflectors
- (8) Monitor column
Optional: Height adjustable
- (9) Storage shelf with USB connection; On/Off switch of ARCADIS Orbic
Optional: Connection for MP3 player including loudspeaker on the rear of the monitor trolley
- (10) Printer bay
- (11) Charge state of uninterruptible power supply (UPS)
- (12) Central locking brake

4. Classification of the Medical Device:

4.1. *Type of the medical device in accordance with the nomenclature for classification:*

Classification according to Annex IX, MDD

Class	IIB
Classification in accordance with Russia Nomenclature classification:	
113880	

4.2. *Class of the Medical Device in accordance with IEC 60601-1:*

Protection against electric shock: Class 1, no applied part according to IEC 60601-1

Operation mode: Continuous operation.

4.3. *Type on contact with the human body:*

The Medical Device is non-sterile, nonexpendable.

Device surfaces that can come into contact with patients are evaluated and proven to be biocompatible according to ISO 10993:2009.

4.4. *IPX classification*

System	IPX0 - No protection
Standard footswitch	IPX8 protected against immersion beyond 1 m
Multifunctional footswitch	IPX8 protected against immersion beyond 1 m
Hand switch	IP54 dust protected, protected against splashing of water

4.5. Laser positioner:

Integrated I.I. laser aimer:

Laser light Class: **1M** (according to IEC 60825-1:2007).

Single-tank laser targeting device:

Laser light Class: Class **2** pursuant to IEC 60825-1:2007.

Horizontal laser light localizer:

Laser light Class: Class **2** pursuant to IEC 60825-1:2007.

4.6. Software classification:

The Medical Device's Software safety Class according to IEC 62304:2006 is **Class B**.

SW version: VC10A

Date of release for the SW: July 2008.

4.7. Waste class

Class A.

General data

Power requirements	100 V, 110 V, 120 V, 127 V, 200 V, 230 V, 240 V ($\pm 10\%$); 50/60 Hz (± 1 Hz). An on-site 16 A/characteristic C overcurrent release (trip breaker) is recommended by DIN VDE 0100-710 for supply network voltages of 230 V/240 V. Country-specific requirements apply for the fuse values to be used for other supply network voltages.
Nominal rating	20 A to 127 V ~, 15 A from 200 V ~ corresponds to the nominal value of the slow-blow fuse in the power input of the product
Internal line impedance	$R_i < 0.3 \text{ Ohm}$ for 100 - 127 V ~ $R_i < 0.8 \text{ Ohm}$ for 200 - 240 V ~
Power consumption	Maximum 2.5 kVA

Weight of the main components

Monitor trolley with 2 monitors, approx. (with PC and UPS)	190 kg
C-arm system approx.	348 kg

Current/voltage values

Voltage (V AC)	Long – term current consumption (A)	Short – term current consumption (A)
100	17	29
110	16	26
120	15	24
127	14	23
200	9	15
230	8	13
240	8	12

Generator

V ($\pm 10\%$); Release (trip) 10-710 for use values	Pulsed Fluoroscopy	Min. pulse width 7 ms Pulse rate 15 f/s/23 mA Max. pulsed output 2.3 kW
ends to the power input	Fluoroscopy	40 kV to 110 kV/0.2 mA to 15.2 mA (max. 1000 W) Max. 8.3 mA/66 kV Max. 15.5 mA/69 kV (in Push Mode at max. 15 sec) Max. 110 kV/5 mA Max. 110 kV/9 mA (in Push Mode at max. 15 sec)
	Digital Radiography	40 kV to 110 kV/0.2 mA to 23 mA (max. 1000 W) Max. 14.3 mA/70 kV Max. 110 kV/9.1 mA
	Tolerances	kV $\pm 10\%$ (measured using spectrometric kV method) mA $\pm 8\% \pm 0.1$ mA (measured in rectified high-voltage circuit) Fluoroscopic time 1 digit (6 s) $\pm 5\%$
	Max. power	1 kW (100 kV/9.9 mA) DR 1 kW (70 kV/14.3 mA) CFC
	Accuracy of radiation output for DR	Variation coefficient < 0.05
	Accuracy of automatic dose rate control	$\pm 20\%$

X-ray tube unit*

SIREPHOS single-tank high-frequency generator	Inverter frequency 15 kHz up to 26 kHz
Inherent filtration	3 mm Al eq. + 0.1 mm Cu
X-ray tube	Stationary anode, focal spot nominal value 0.6

The X-ray tube and image receptor are geometrically aligned such that the emitted radiation line (radiation axis) is always perpendicular to the plane of the image intensifier.

Max. output and focuses in operating modes

	CFC/SUB/ ROADMAP/PFC (≤ 2 f/s)	PFC (> 2 f/s)	DR (k = 1)	DR (k > 1)
tube voltage	40 – 110 kV	40 – 110 kV	40 – 110 kV	40 – 110 kV
tube current	0.2 – 15.2 mA	1 – 23 mA	1 – 23 mA	0.2 – 14.3 mA
exposure time	n.a.	7 – 24 ms	7 – 42 ms	n.a.
pulsed output	n.a.	max. 2.3 kW	max. 1.8 kW	n.a.
clean power	max. 550 W max. 1000 W for CFC and PFC	max. 550 W max. 1000 W for	n.a.	n.a.

	high contrast	PFC high contrast		
--	---------------	-------------------	--	--

C-arm

Orbital movement manual	190° (± 95°)
C-arm angulation	± 190°
C-arm horizontal movement	20 cm
C-arm immersion depth	73 cm
C-arm swivel range	± 10°
C-arm vertical movement	40 cm, motorized
SID	98 cm
Tube-image receptor distance	78 cm

Extension for 3D mode

Motor drive	
Motor-driven orbital movement through	190° for online reconstruction of isotropic 3D data sets
Rotation time	30 seconds with 50 images/scan 60 seconds with 100 images/scan

Image intensifier

SIRECON 23-2int. 1k	Nominal diameter 23 cm (9")
Format switch-over	23 cm/15 cm (9"/6")
Size of I.I. input field	21.5 cm
Detective quantum efficiency (DQE)	≥ 61% according to IEC 61262

Scattered radiation grid at I.I. input

Round grid	17/70 f _{0.100}
System attenuation factor	m = 1.5

Collimator system

Collimator system	Iris diaphragm for concentric collimation and semi-transparent slot diaphragm for collimation with unlimited rotation
-------------------	---

Imaging chain

Camera with CCD Sensor	1024 (H) x 1024 (V)
Aspia imaging system	Features, e.g. digital image rotation +360 °

	monitors	
	FT color displays	Resolution 1280 x 1024 (pixels)
	screen diagonal	48 cm (19")
	luminance	Typically 180 cd/m ² /max. 250 cd/m ²
	fw displays	Resolution 1280 x 1024 (pixels)
	screen diagonal	48 cm (19")
	luminance	Typically 400 cd/m ² /max. 600 cd/m ²
	dose measuring chamber (DAP meter)	
	technology	Ionization chamber
	active area	Diameter 90 mm, area 63.2 cm ²
	DAP resolution	0.01 μGym^2 equivalent to cGycm ²
	Maximum measurable DAP	3000.00 μGym^2 equivalent to cGycm ²
	measuring inaccuracy	< 1% under constant pressure and temperature
tropic 3D	energy range	40 - 125 kV $\pm$ 6%
	attenuation equivalent (inherent filtration)	< 0.5 mm Al equivalent
	dose rate linearity	Better than $\pm$ 2%
	measuring readiness	10 seconds after power on
	accuracy of dose	< 25%
	display/deviation of displayed from measured dose area product	
	Work environment	Temperature range: +10 °C to +70 °C Rel. humidity: 15% to 75% (without condensation) Barometric pressure: 700 hPa to 1060 hPa

Integrated I.I. laser aimer

	laser class	1M (IEC 60825-1:2014)
	laser type	Semiconductor laser (laser diode)
	wave length	635 nm
	color	Red
rent slot	effective power	$\leq$ 0.2 mW

Single-tank laser targeting device

	laser class	2 (IEC 60825-1:2014)
	laser type	Semiconductor laser (laser diode)
	wave length	655 nm
	color	Visible red
	effective power	< 1 mW

Horizontal laser light localizer

Laser class	2 (IEC 60825-1:2014)
Optical setup	Line generating optics
Beam angle	20° to 25°
Wave length	540 - 700 nm
Output power on the module	0.5 mW
Voltage supply	1.5 volt battery operation (4 R6 (AA) batteries)
Current consumption	80 mA
Battery type	round cell R6 (AA) Alternatively primary cells R6/LR6 (1.5 V) or NiCd or NiMH batteries
Service life	At least 20 hours of continuous operation
Aut. time limitation	1 min. (typical)
Case material	Aluminum
Weight	Approx. 1.6 kg
Temperature range	Up to 37.5°C

Standard footswitch

Dimensions	213 mm x 110 mm x 39 mm
Weight	1.5 kg
Length of the cable (if applicable)	4.5 m

Multi-function foot switch

Dimensions	378 mm x 194 mm x 47 mm
Weight	4.5 kg
Length of the cable (if applicable)	4.5 m

DVI video splitter

Outputs	2 single-link DVI-D; female connectors; up to 1920 x 1200 @ 60 Hz or 1080 p
Conditions for use	The displays or monitors used must be compatible with the SXGA standard 1280 x 1024 at 60 Hz. No galvanic separation, therefore the connected monitors must comply with IEC 60601-1-1

Horizontal laser light localizer

Laser class	2 (IEC 60825-1:2014)
Optical setup	Line generating optics
Beam angle	20° to 25°
Wave length	540 - 700 nm
Output power on the module	0.5 mW
Voltage supply	1.5 volt battery operation (4 R6 (AA) batteries)
Current consumption	80 mA
Battery type	round cell R6 (AA) Alternatively primary cells R6/LR6 (1.5 V) or NiCd or NiMH batteries
Service life	At least 20 hours of continuous operation
Aut. time limitation	1 min. (typical)
Case material	Aluminum
Weight	Approx. 1.6 kg
Temperature range	Up to 37.5°C

Standard footswitch

Dimensions	213 mm x 110 mm x 39 mm
Weight	1.5 kg
Length of the cable (if applicable)	4.5 m

Multi-function foot switch

Dimensions	378 mm x 194 mm x 47 mm
Weight	4.5 kg
Length of the cable (if applicable)	4.5 m

DVI video splitter

Outputs	2 single-link DVI-D; female connectors; up to 1920 x 1200 @ 60 Hz or 1080 p
Conditions for use	The displays or monitors used must be compatible with the SXGA standard 1280 x 1024 at 60 Hz. No galvanic separation, therefore the connected monitors must comply with IEC 60601-1-1

Video splitter live

Outputs	Output for connecting an external live/reference monitor. Implemented as a splitter solution for connection with a 15-pin VGA plug
Conditions for use	The displays or monitors used must be compatible with the SXGA standard 1280 x 1024 at 60 Hz. No galvanic separation, therefore the connected monitors must comply with IEC 60601-1-1

Video splitter Ref.

Outputs	Output for connecting an external live/reference monitor. Implemented as a splitter solution for connection with a 15-pin VGA plug
Conditions for use	The displays or monitors used must be compatible with the SXGA standard 1280 x 1024 at 60 Hz. No galvanic separation, therefore the connected monitors must comply with IEC 60601-1-1

Printers

	Sony Printer UP-971	Sony Printer UP-991
Voltage, V	100 bis 240 V AC	100 bis 240 V AC
Frequency, Hz	50/60 Hz	50/60 Hz
Current	2.9 A to 1.2 A	2.9 A to 1.2 A
Type of printing	Thermal paper	Thermal paper
Resolution	325 dpi	325 dpi
Grayscale number	256 bit processing	256 bit processing
Color	black-	black-
Size of paper	210 mm	210 mm

Dimensions, mm	316 x 132.5 x 265 mm	316 x 132.5 x 265 mm
Mass, kg	7 kg	7 kg

265 mm	Motion sound system	
	Motion MIDI interface (Aux In)	3.5 mm plug
	Range of sound	20 – 20.000 Hz

LAN-Client RoW NG

Data transfer	
Ethernet transfer rate	BASE-T max. 100 MBit/s
Wireless transfer rate	max. 150 Mb/s
Radiated power	typically: 100 mW, max. 200 mW
Supported wireless standards	IEEE 802.11a IEEE 802.11b IEEE 802.11e IEEE 802.11g IEEE 802.11h IEEE 802.11i IEEE 802.11n
Supported standards for power supply	IEEE 802.3at Type 1 (802.3af) and 802.3at Type 2 (power supply via Ethernet)

6. Notes concerning electromagnetic compatibility (EMC)

Medical electrical equipment is subject to special precautions regarding EMC. These systems must be installed and put into service according to the EMC information provided in the accompanying documents.

Portable and mobile RF communications equipment can affect medical electrical equipment.

6.1. Guidelines and manufacturer's declaration – Electromagnetic emissions

This system is intended for use in the electromagnetic environment as specified below. The installer or the user of the system should ensure that it is operated in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions according to CISPR 11	Group 1	The system uses RF energy exclusively for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions according to CISPR 11	Class A	The system is not suitable for use in domestic establishments and establishments directly connected to a public power supply network which supplies buildings used for domestic purposes.
Harmonic emissions according to IEC 61000-3-2	Not applicable	The system is a professionally used unit with a total rated power larger than 1 kW. There are no limit values for this rated power.
Voltage fluctuations/ flicker emissions according to IEC 61000-3-3	Complies with the specified regulation	

6.1.1. GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

The system is intended for use in the electromagnetic environment as specified below. The customer or the user of the system should ensure that it is operated in such an environment.

Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic materials, the relative humidity must be at least 30%.
Electrical fast transient disturbances/bursts according to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge according to IEC 61000-4-5	± 1 kV differential-mode voltage ± 2 kV common-mode voltage	± 1 kV differential-mode voltage ± 2 kV common-mode voltage	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines according to IEC 61000-4-11	$< 5\% U_T$ * ($> 95\%$ dip in U_T) for 0.5 cycles $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles	Not tested	As the product has a > 1 kVA power consumption, these tests are not required.
Voltage dips, short interruptions and voltage variations on power supply input lines according to IEC 61000-4-11	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 seconds	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the system requires continued operation during power mains interruptions, it is recommended that the system be powered from an uninterruptible power supply.
Power frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.

Remark: U_T is the AC mains voltage before application of the test level.

6.1.2. GUIDANCE AND MANUFACTURER'S DECLARATION — ELECTROMAGNETIC INTERFERENCE IMMUNITY

The system is intended for use in the electromagnetic environment as specified below. The customer or the user of the system should ensure that it is operated in such an environment.

Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF telecommunications equipment should be used no closer to any part of the system (including cables) than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance
Conducted RF disturbances according to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = 1.2\sqrt{P}$
Radiated RF disturbances according to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ for 80 MHz to 800 MHz
			$d = 2.3\sqrt{P}$ for 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Remark 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Remark 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radios, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to the fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device or system is used exceeds the applicable RF compliance level above, the device or system should be observed to verify normal operation. Should unusual performance features be observed, additional measures (such as change in orientation or change of site of the system) may be necessary.			
^b Over the frequency range of 150 kHz to 80 MHz the field strength should be less than 3 V/m.			

6.1.3. RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE ANGIOGRAPHIC SYSTEM

The system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the system as recommended below, according to the maximum output power of the communications equipment.


Rated maximum output power of transmitter [W]	Safe distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Remark 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Remark 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

7. Scheme of the Production line.

System assembly

- Assembling of the system components (e.g. floorstand, ceilingstand, table, tube collimator, generator etc.) according to the customer order
- Torque wrench and loctite (sealing-wax) are for assembling, where applicable
- Check of shock-, temperature- and humidity sensor, where applicable
- Complete system cabling (screwing and plugging only)

System testing

*X-Ray devices only

System start up (Configuration & adjustment)

- SW installation
- Licenses Installation
- System configuration and startup
- Components Adjustment
- Adjust & check system movements (e.g. stand, table, ...)*
- Tube conditioning*
- Detector calibration*
- System dose adjustment*
- Collimator, detector and tube adjustment*
- Functional check of options

Final testing

- Perform CRST (Clinically relevant simulation testing)*
- Order system documentation
- Perform IQAP (Monitor adjustment, dose control, Image quality test)*
- Perform special legal requirements

Final work

- Perform customer configuration
- Save system data on CD/DVD
- Create system backup
- Apply system label
- Complete customer documents (depending on customer order)

Quality check

- System check by Designated Q-Person SCM (includes a check of functionality, a check of completeness of customer order, test protocols, customer documentation and system labeling)


System disassembly

- Move the system to transport position for applicable components
- Disconnect system cabling
- Assemble transport racks
- Prepare and provide system components for packing


Hand over to logistic provider


8. The list of components and accessories ARCADIS Orbic.

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a continuous 1kV imaging chain. A 23 cm image intensifier with format switchover and a 2.3 kW high-voltage generator ensure an optimal fluoroscopic

acquisition speed of up to 8 f/s is achieved

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



Monitor column with motorized height adjustment and vertically rotatable monitors (180 degrees) for flexible positioning of the TFT displays with integrated cable routing and fold-up function for transport and park position.
Monitor column with monitors that swivel vertically through 210 degrees (-30° to +180°) for flexible positioning of the TFT displays with integrated cable routing.




	(UK)	English (GB / British) layout	
4. Измеритель дозы с измерительной камерой.	Dose measuring chamber # ARC	System integrated dose measuring chamber for displaying the dose area product or air kerma value. The cumulative dose area product is displayed for the current patient on the examination card and saved under the patient data. The cumulated dose is automatically transferred to a radiation summary report and can be retrieved at any time. For each patient a cumulative value is saved in the patient database.	



5. Комплект кабелей и приспособлений для подключения к сети.	Dual IT-port NG	Network hub for simultaneous wired connection of the ARCADIS C-arm system with an IT network and another network-capable system, such as standalone navigation (also required for connecting a W-LAN module).	
6. Педальный переключатель, варианты исполнения: - стандартный; - многофункциональный.	Standard footswitch # ARC Mult. funct. f. ARC Orbic 3D Opt	Standard footswitch for radiation release Ergonomic multifunctional footswitch for the control of radiation functions.	

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









13. Интерфейс подключения навигационной системы.	NaviLink 2D	Integrated navigation interface for connecting navigation systems from other manufacturers with digital DICOM 3.0 image transfer to the navigation system 2D Image transfer		
--	-------------	---	--	--



16. Принтер формат.	Sony Printer UP-971AD Sony Printer UP-991AD	Thermal printer for economical and fast image documentation (gray-scale printing) on thermal paper and film	
17. Инсталляционный комплект для принтера.	Installation kit for printers	Kit is intended for the installation of the printers	


			
18. Ограничитель расстояния фокус-кожа.	DHHS-Spacer # Orb	Single-tank spacer for increasing the minimum source-skin distance to 30 cm.	


19. Кабель заземления.	Grounding cable # ARC	With two clips	
20. Лазерный указатель на усилитель изображения.	Integrated I.I. laser aimer # ARC	I.I. laser aimer integrated in the image intensifier housing for exact positioning without radiation with undertable position of the single tank, with remote control via the control panel at the C-arm.	



21. Моноблочный лазерный указатель.	Single-tank laser target.device#ARC	Laser targeting device integrated in the single tank	
22. Контрольный комплект системы 3-х мерных изображений (3D).	VRT	3D VRT (Volume Rendering Technique) is an additional enhancement of the comprehensive 3D functions with ARCADIS Orbic 3D.	


23.	Интерфейс подключения навигационной системы для 3-х мерного изображения.	NaviLink 3D	Integrated navigation interface for connecting navigation systems from other manufacturers with digital DICOM 3.0 image transfer to the navigation system. 2D and 3D Image transfer
24.	Механический замок (блокиратор) для предупреждения включения рентгеновского излучения		Mechanical lock secured with a key for disabling / enabling the release of X-radiation. The lock only disables radiation release and thus protects against improper use. Use of the system as a viewing workstation and for image processing is possible without the key. Two keys are included in the scope of delivery. Individual keys cannot be reordered. Mandatory in the UK and Ireland, otherwise optional.


25. Калибровочный фантом 3D, варианты исполнения: - стандартный; - расширенной функциональности.	Control kit 3D # Orb NORMI 3D	The control kit contains auxiliary devices for performing the 3D-specific test steps required as part of the acceptance and constancy test. Test object for quality control and acceptance tests of digital volume tomography X-ray units acc. to DIN 6868-150. Includes cylindrical cavities (each 4x with 1.3/1/0.9/0.8/0.7/0.6/0.5 mm) for testing the spatial resolution.	
26. Программный пакет наложения 3D изображений.	SW option 3D image fusion	The 3D Image Fusion package enables spatial alignment and visualization of image data of one patient where image data has been generated at different points in time by different modalities.	

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
30. Акустическая система Emotion.	Sound System EMotion	Integrated Sound System EMotion with interface for operating external (MP3) audio equipment via Aux-in jack (Aux input) including amplifier and stereo loudspeaker for high-fidelity sound.	
--	-----------------------------	--	---

9. Labeling information.

ARCADIS Orbic System label	
Until mid. of July 2016 (also for refurbished systems)	Starting mid. of July 2016.
SIEMENS  (1P) Model No. 01234567  (S) Serial No. 1001 Siemens AG, Wittelsbacherplatz 2, DE- 80333 Muenchen Made in Germany CE 0123 2002 System IVK	SIEMENS ARCADIS Orbic  (01) 0405686900025 GTIN (240) 10143407 MODEL (422) 276 ORIGIN DE (21) XXXXX SERIAL 2017 0123 SYSTEM IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany

The position of the RF beam of the activated WLAN is shown in a warning label.

(1) Warning label: Non-ionizing electromagnetic field

(2) Warning: Read Operator Manual

Artwork of the MD labeling in Russian

Additionally, on the medical devices which are delivered to Russia, the Russia Specific label is applied with the following content (please note that Labeling mentioned below is applied only for the systems which are delivered to Russia not earlier than May 2020):

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10.The list of International Standards and Directives applied to the Medical Device.

Standard organization	Reference number and date	Title of Standard
IEC	60601-1:2005	Medical electrical equipment - Part 1. General requirements for basic safety and essential performance
IEC	60601-1-1: 2000	Medical electrical equipment – Part 1-1: General requirements for safety - Collateral standard: Safety requirements for medical electrical systems
IEC	60601-1-2: 3rd 2007	Medical Electrical Equipment – Part 1-2: General requirements for basic safety and essential performance -Collateral Standard: Electromagnetic Compatibility - Requirements and Tests
IEC	60601-1-3: 3rd 2008	Medical electrical equipment – Part 1-3: General requirements for basic safety and essential performance -Collateral Standard: Radiation protection in diagnostic X-ray equipment
IEC	60601-1-4: 1996 +A1 1999	Medical electrical equipment – Part 1-4: General requirements for safety; 4th collateral standard: Programmable electrical medical systems
IEC	60601-1-6:2010	Medical electrical equipment. Part 1-6. General requirements for basic safety and essential performance - Collateral standard: Usability
IEC	60601-2-7:1998	Medical electrical equipment - Part 2-7: Particular requirements for the safety of high-

		voltage generators of diagnostic X-ray generators
IEC	62366: 2007	Medical devices - Application of usability engineering to medical devices
IEC	14971: 2012	Medical devices - Application of risk management to medical devices
IEC	62304 : 2006 + AC:2008	Medical device software - Software life cycle processes
IEC	60601-2-28: 3rd Edition 2010	Medical electrical equipment – Part 2-28: Particular requirements for the basic safety and essential performance of X-ray tube assemblies for medical diagnosis
IEC	60601-2-54: 2009	Medical electrical equipment – Part 2-54: Particular requirements for the basic safety and essential performance of X-ray equipment for radiography and radioscopy
IEC/PAS	61910-1: 2007	IEC PAS 61910-1 First Edition 2007-07, Medical electrical equipment - Radiation dose documentation – Part 1: Equipment for Radiography and Radioscopy
IEC	62220-1-3	Medical electrical equipment - Characteristics of digital X-ray imaging devices - Part 1-3: Determination of the detective quantum efficiency - Detectors used in dynamic imaging
IEC	60825-1: 2007	Safety of laser products - Part 1: Equipment classification, and requirements
DICOM Standard	2009	Digital Imaging and Communications in Medicine
ISO	14971:2007	Medical devices – Application of risk management to medical devices

11. Disposal.

	The use of this symbol indicates that this product should not be treated as household waste. By ensuring that this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste-handling of this product. For more detailed information concerning the return and recycling of this product, please consult the supplier from whom you purchased the product
---	--

Concerning Batteries (Directive 2006/66/EC) The directive 2006/66/EC requires separate collection and appropriate disposal of spent batteries. This product contains batteries that are not intended to be replaced by the user. Replacement of those batteries will usually be done during regular maintenance or service by service staff who can also arrange proper disposal.

If you want to remove the product from service, take into consideration that public legal directives may contain special regulations regarding disposal of this equipment. In order to ensure that these legal regulations are complied with and to avoid potential environmental hazards which may be caused by the disposal of your system, please consult Siemens Customer Service.

- Batteries and packaging material must be disposed of in an environmentally safe manner according to national regulations.
- The sterile single-use covers must be disposed of in accordance with national regulations or the rules of the hospital.
- For further information about the disposal of the product, please refer to our service documents.

12. Packaging.

The packaging consists of wood, nails, cardboard, paper and plastic. The packing list indicating names, quantities, packaging date, and packer reference number is inserted in the shipping box.

The following items appear on the label affixed to each package:

- Manufacturer trademark;
- Device name;
- Manufacturing date.

The following items appear on the label made according to the manufacturer's working drawings and affixed to each multiple shipping box:

- Medical device manufacturer name;
- Net and gross weight;
- Box dimensions

Transport marking. The shipping box is marked with handling symbols which have the following meanings: "This side up", "Keep away from moisture", "Fragile, handle with care!".

Packaging labeling

	"Fragile" Indicates that medical device can get broken or damaged if not handled carefully.
	"This side up"
	"Keep away from moisture" Indicates that the medical device has to be protected from moisture.

C-arm packaging process


- Packaging and transportation of C-arm is performed after it's positioning on transport platform.
- Than special support unit is mounted which is secured by brackets. Also other parts of C-arm are secured by brackets.
- Bracing cables and mounting brackets secure box with printer and documentation.
- The packing list is inserted in the shipping box. Package is covered by protective layer and secured by yellow ramps.

Monitor trolley packaging process.

- After all wheel lock ups unlock monitor trolley is rolled on transportation platform and secured by special foot blocks.
- Wheels are turned along the line and secure them with brakes.
- Foot blocks positioned at right and left sides provides secure of medical device on transportation platform.
- Secured monitor and keyboard are covered with protective layer which protects them from damage during transportation.

- The packaging process is completed with securing side panel boards and labeled cover up (see beginning p.9).

Real photos of the packaging:

13. Maintenance.

13.1. Regular maintenance

To ensure the safety of patient, operation personnel and others, tests must be performed **regularly** to maintain the safety and proper functioning of the system.

If more frequent inspections and/or maintenance are required by national regulations, compliance with these regulations is essential.

All parts of this system is which wear could lead to a hazard must be inspected by trained technical personnel and replaced if necessary.

Maintenance work can be performed by Siemens Customer Service for a charge.

13.2. Service contacts

Regular inspections should be a part of the annual maintenance performed by Siemens Service under a service contract. We therefore recommend that you take out a service contract.

It is also possible to take out a service contract based on system usage.

Please contact your authorized local representative in this matter.

The name and address of the local representative in Russia:

Siemens Healthcare LLC
Dubininskaja street 96,
115093 Moscow, Russia
Tel.: +7 495 737 12 52
www.siemens-healthineers.ru

14. Service life.

Our medical engineering products are designed for operation under average conditions in accordance with the Operator Manual for a useful life of **10** (ten) years. If the products are operated for longer than this time, additional checks and possibly repairs beyond the usual maintenance procedures may be necessary in order to ensure the functional integrity and operating safety of these products. Please contact Siemens about these measures early enough!

15. Warranty.

Siemens guarantees that the Equipment has no defects caused by installation and, as of the date of risk transfer, is compliant with the approved specifications and has no defects in materials and production within 14 (fourteen) months as of the Equipment delivery date. Warranty terms and conditions for components made by the third parties and supplied along with the Equipment shall be detailed in the respective Order.

Warranty period for operable Equipment components supplied to replace the defected ones shall be 6 (six) months as of their delivery date. Total warranty period for defective components repaired or installed by Siemens shall not exceed 14 months as of the Equipment delivery date.

Defective Equipment (its part), Spare parts, high-vacuum components, other types of Products replaced within the warranty period shall be handed over to Siemens to be further returned to the manufacturer no later than 14 (fourteen) days as of the date of signing the documents (certificate/order form) provided by Siemens, confirming that the respective defect was eliminated. Title to such spare part shall be transferred to Siemens.

16. Complaints

The complaint may be submitted only for those issues which have not been the subject of the device acceptance procedure performed in accordance with the terms of the contract.

In case the system does not operate properly or does not perform commands as described in operating manual contact manufacturer representative service center in RF.

Siemens Healthcare Limited Liability Company (Siemens Healthcare LLC)

115093, Moscow,

Dubininskaya str. 96

Phone: (495) 737-12-52;

www.siemens-healthineers.ru


Siemens Healthcare GmbH

This document has been printed by the manufacturer in Germany.

Siemens Healthcare GmbH

Siegfried Quinger
Head of Quality
Business Area Advanced Therapies


Ali-Ihsan Can
Head of Regulatory Affairs
Business Area Advanced Therapies

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Berlin, 27. Oktober 2021

I, the undersigned notary public Dr. Peter Engel, admitted and resident in Berlin (Germany), hereby **certify** the conformity of the above photocopy with the original presented to me.

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