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ETHICON, LLC

475 C Street Los Frailes Industrial Park
Guaynabo, PUERTO RICO 00969, USA

Manager, Regulatory Affairs Global Strategy Execution

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This 28th day of October, 2021

Keith Petersen
Notary Public

INSTRUCTION OF USE FOR Medical Device

Monofilament Surgical suture with a needle for knotless tissue fixation
STRATAFIX™ Symmetric PDS™ Plus antibacterial coated in variations

*Version for the realization on the territory of
Russian Federation*

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NAME OF PRODUCT

Monofilament Surgical suture with a needle for knotless tissue fixation
STRATAFIX™ Symmetric PDS™ Plus antibacterial coated in variations:

- M2 (3/0) 45 cm Prime Reverse Cutting Needle PS-1 prime
- M2 (3/0) 45 cm Prime Reverse Cutting PS-2 prime
- M3.5 (0) 45 cm Reverse Cutting Needle OS-6
- M4 (1) 45 cm Reverse Cutting Needle OS-6
- M4 (1) 45 cm Blunt Point Ethiguard Needle CTXB
- M4 (1) 45 cm Blunt Point Ethiguard Needle CTB-1
- M3.5 (0) 45 cm Blunt Point Ethiguard Needle CTXB
- M4 (1) 45 cm Taperpoint Needle CTX
- M4 (1) 45 cm Taperpoint Needle CT-1
- M4 (1) 45 Taperpoint Needle CT
- M3.5 (0) 45 cm Taperpoint Needle CT-1
- M3.5 (0) 45 cm Taperpoint Needle CTX
- M3 (2/0) 45 cm Taperpoint Needle CT-1
- M3.5 (0) 45 cm Taperpoint Needle CT
- M3.5 (0) 45 cm Taperpoint Needle CT-2
- M3 (2/0) 45 cm Taperpoint Needle CT-2
- M3 (2/0) 45 cm Taperpoint Needle SH
- M2 (3/0) 45 cm Taperpoint Needle SH
- M3.5 (0) 60 cm Reverse Cutting Needle OS-6
- M4 (1) 60 cm Reverse Cutting Needle OS-6
- M3.5 (0) 60 cm Reverse Cutting Needle OS-8
- M4 (1) 60 cm Reverse Cutting Needle OS-8
- M3.5 (0) 60 cm Blunt Point Ethiguard Needle CTB-1
- M4 (1) 60 cm Blunt Point Ethiguard Needle CTB-1
- M3.5 (0) 60 cm Blunt Point Ethiguard Needle CTXB
- M4 (1) 60 cm Blunt Point Ethiguard Needle CTXB
- M2 (3/0) 15 cm Taperpoint Needle RB-1
- M3.5 (0) 30 cm Taperpoint Needle CT
- M2 (3/0) 15 cm Taperpoint Needle SH

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MEDICAL DEVICE PHYSICAL MANUFACTURER

ETHICON, INC.

Calle Durango No. 2751 Lote Bravo, Ciudad Juarez Chihuahua, MEXICO C.P. 32575, Mexico

INTENDED USE

The device is indicated for general soft tissue approximation.

INDICATIONS

The suture STRATAFIX™ Symmetric PDS™ Plus Devices are indicated for general soft tissue approximation where use of an absorbable suture is appropriate.

CONTRAINDICATIONS

The suture STRATAFIX™ Symmetric PDS™ Plus, being absorbable, is not to be used where prolonged (beyond six weeks) approximation of tissue under stress is required and is not to be used in conjunction with prosthetic devices (i.e., heart valves or synthetic grafts).

The device should not be used in patients with known allergic reactions to IRGACARE MP (triclosan).

DESCRIPTION OF MEDICAL DEVICE

STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device is an antibacterial (polydioxanone) monofilament, synthetic absorbable device prepared from the polyester, poly (p-dioxanone). The empirical molecular formula of the polymer is (C₄H₆O₃)_x. The device contains IRGACARE®* MP (triclosan), a broad spectrum antibacterial agent, at no more than 2360 µg/m. STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device is dyed with D&C Violet No. 2.

STRATAFIX™ Symmetric PDS™ Plus Device consists of an absorbable thread with unidirectional anchors, equipped with a surgical needle at one end and a

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fixation tab at the other. The anchors and fixation tab design allows for tissue approximation without the need to tie surgical knots.

The cross-sectional area of the STRATAFIX™ Symmetric PDS™ Plus Device core is comparable to the cross-sectional area of USP suture with the same size designation.

The USP and EU Pharmacopoeia sizes of the STRATAFIX™ Symmetric PDS™ Plus Device are further defined in Table 1.

Table 1. Diameter Comparison

USP SUTURE SIZE DESIGNATION	EU PHARMACOPOEIA SUTURE SIZE (Metric / Ph. Eur.) DESIGNATION	STRATAFIX™ Symmetric PDS™ Plus DEVICE CORE CROSS-SECTIONAL AREA	
		USP	Metric / Ph. Eur.
1	4	1	4
0	3.5	0	3.5
2-0	3	2-0	3
3-0	2	3-0	2

STRATAFIX™ Symmetric PDS™ Plus Device straight tensile strength is stronger than the knot tensile strength for a USP and EU Pharmacopoeia polydioxanone suture of the equivalent size as shown in Table 2.

Table 2. Tensile Strength Comparison

USP SUTURE SIZE DESIGNATION	EU PHARMACOPOEIA SUTURE SIZE (Metric / Ph. Eur.) DESIGNATION	SUTURE MINIMUM KNOT TENSILE STRENGTH	
		USP (lbf)	Metric / Ph. Eur. (N)

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Polydioxanone has been found to be nonallergenic, nonpyrogenic, and elicits only a slight tissue reaction during absorption.

USAGE CONDITION

This medical device is intended for professional use in the medical and health care institutions.

AREA OF USE

Surgery

POTENTIAL CONSUMER

Medical Device is intended to be used by professional doctor.

WARNINGS

The safety and effectiveness of the suture STRATAFIX™ Symmetric PDS™ Plus have not been established in:

- neural tissue,
- cardiovascular tissue,
- vaginal mucosa tissue,
- microsurgery, or in ophthalmic tissue.

Under certain circumstances, notably orthopedic procedures, immobilization by external support may be employed at the discretion of the surgeon.

Care should be undertaken to avoid leaving barbed suture ends adjacent to the peritoneum in extra-peritoneal tissue closure.

Avoid contacting STRATAFIX™ Device and associated needles with other materials (e.g. surgical gauze, drapes, etc.) in the surgical field to prevent ensnaring on the barbs. If the barbs catch, carefully pull the material in the opposite direction of the needle to disengage it from the barbs.

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Small bowel obstruction (SBO): including volvulus, bowel infarction, and significant morbidity, have been reported due to barbs or barbed suture ends hooking onto adjacent small bowel and/or mesentery, such as in peritoneal closure.

Be sure to read the Instructions for Use section for proper fixation tab seating.

Do not pull on fixation tab after seating it.

Do not resterilize/reuse.

Reuse of this device (or portions of the device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.

Discard open packages and unused devices.

Do not use after expiration date.

STRATAFIX™ Symmetric PDS™ Plus devices are intended for single patient use only.

PRECAUTIONS

The suture STRATAFIX™ Symmetric PDS™ Plus has not been studied for use in interrupted suture patterns.

Subcuticular devices should be placed as deep as possible in order to minimize the erythema and induration normally associated with absorption.

The suture STRATAFIX™ Symmetric PDS™ Plus has not been studied for ligating vessels

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or luminal structures.

As with any device, care should be taken to avoid damage to the strand when handling. Avoid the crushing or crimping action of surgical instruments, such as needle holders and forceps.

Vaginal mucosal devices remaining in place for extended periods may be associated with localized irritation and should be removed when necessary.

Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.

Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse events associated with wound closure (including wounds closed with synthetic absorbable sutures), include:

- wound separation or dehiscence,
- failure to provide adequate wound support in closure of sites where expansion, stretching, or distention occurs, or in elderly, malnourished or debilitated patients or in patients suffering from conditions that may delay wound healing,
- infection,
- minimal acute inflammatory tissue reaction,
- device extrusion and delayed absorption in tissue with poor blood supply,
- calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs,
- transitory local irritation at the wound site, as well as
- allergic reaction to IRGACARE MP (triclosan).

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Broken needles may result in extended or additional surgeries or residual foreign bodies.

Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

MRI SAFETY INFORMATION

STRATAFIX™ Symmetric PDS™ Plus Device is MR Safe.

HOW SUPPLIED

The STRATAFIX™ Symmetric PDS™ Plus Device is available as sterile, monofilament, dyed (violet) strands in sizes from 3-0 to 1 (metric sizes 2-4) in a variety of lengths and needles.

The STRATAFIX™ Symmetric PDS™ Plus Device is available in boxes of one dozen and instruction for use 1 piece.

RISK ANALYSES

Risks associated with the release and uses of the devices have been evaluated as per Franchise Procedure for Medical Device Risk Management Plan.


This procedure guides identification of hazards early in the development process to eliminate or reduce them as far as possible. Where hazards still may exist, controls are established in the design and manufacturing processes to assure that risks are acceptable. Compliance to stated claims and specifications are assured by verification and validation. Design controls demonstrate that the devices perform as intended.

Risks associated with the devices continue to be monitored throughout the product lifecycle, Post-Market Surveillance Procedure.

INSTRUCTIONS FOR USE

1. SEATING THE FIXATION TAB

Seat the fixation tab in intact tissue above the apex and away from the incision

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
- a. Take the first pass directly above or adjacent to the apex in a direction away from the incision (Figure 1a shows directly above the apex).

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- b. Pull the suture STRATAFIX™ Symmetric PDS™ Plus through the tissue to gently seat the fixation tab against the tissue. The fixation tab should be seated above the tissue plane and be visible (Figure 1b). Do not exert additional force on the fixation tab or device.
- c. Moving toward the apex of the incision, take a pass in the intact tissue perpendicular to the initial pass to lock the stitch (Figure 1c). Multiple passes are acceptable. Pull gently on the device to take up any slack.

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


2. CLOSING THE INCISION

- a. Proceed with a continuous suturing pattern to close the incision, taking apposing bites on either side of the wound in standard fashion (Figure 2a).
- b. To achieve the desired approximation and tension, gently pull on the device with each tissue passage.
- c. Continue approximation along the length of the incision taking care not to overtighten the device. Retain approximately 3" (8cm) of the device for completing reverse stitching to secure the terminal end.

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- а. To complete and secure the closure, take two passes in the reverse direction across the incision (Figure 3a). Then gently pull on the free end of the device and cut flush with the surface of the tissue.

If removal of the device is required, cut the core of the device near the fixation tab and carefully pull the device in the direction of the running stitch to disengage the anchors. Cutting the core of the device in multiple places may be required to fully remove the device from the tissue. Retrieve the fixation tab and discard.

DEIMPLANTATION

N/A Absorbable Suture

TECHNICAL CHARACTERISTICS

Table 3. Table of functional and safety characteristics.

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Two important characteristics describe the in vivo performance of absorbable devices: first, tensile strength retention, and second, the absorption rate (loss of mass).

The suture STRATAFIX™ Symmetric PDS™ Plus has been formulated to minimize the variability of these characteristics and to provide wound support through an extended healing period.

Data obtained from an animal implantation study show that the suture STRATAFIX™ Symmetric PDS™ Plus is completely absorbed by 210 days post implantation. The results of implantation studies indicate the percentage of original strength that is retained (Table 1)

Table 4. Original Strength depending on time of implantation Retained

Using in vitro studies, STRATAFIX™ Symmetric PDS™ Plus Device with IRGACARE MP (triclosan) has been shown to inhibit colonization of the device by Staphylococcus aureus, Staphylococcus epidermidis, methicillin-resistant S. aureus, methicillin-resistant S. epidermidis, Escherichia coli, and Klebsiella pneumoniae, which are microorganisms known to be associated with surgical site infections. Animal studies have demonstrated that STRATAFIX™ Symmetric PDS™ Plus Device inhibits bacterial colonization of the device after direct in vivo challenge with bacteria.

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MATERIALS USED FOR THE MEDICAL DEVICE

Suture

Polydioxanone (C₄H₆O₃)_x IRGACARE MP dyed with D&C Violet No. 2 60725, Suture material manufactured by Ethicon, Inc., Cornelia, Georgia, USA
Antibacterial agent - triclosan IRGACARE MP, manufacturer BASF Grenzach GmbH, Germany, GMP certificate №DE_BW_01_GMP_2018_0014.

Needle

Stainless steel AISI grade 420*, Supplied by Carpenter Technology Corp. Reading, PA
Ethalloy Stainless steel*, manufactured by ATI Allvac, USA.
Needle coverage – Silicone (CAS Number 01330-20-7 63148-57-2) or Cerberus, manufactured by Gelest Inc. or Siltech MED-ET1 (Siltech Corporation)

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INFORMATION ABOUT MATERIALS OF ANIMAL/HUMAN ORIGIN, PHARMACEUTICAL SUBSTANCES

This medical product is not manufactured using animal tissue in accordance with Commission Regulation (EC) 722/2012, and is not used with animal tissues or by-products of animal origin within the framework of Commission Regulation (EC) 722/2012.

This medical device does not contain and is not used in any processes with derivatives of human blood (HBD).

This medical product does not contain or is not used in any processes with medicinal products.

This medical product contains a broad spectrum bactericidal agent IRGACARE MP (triclosan), which can not be attributed to medicinal products, since its concentration in this medical product is not more than 2360 µg / m.

The BSI consultation process with the Medical Evaluation Board, CBG-MEB, regarding this medical product was performed in 2007 according to the consultation procedure number: CP # 050. This MEB report is saved in the BSI file.

STERILIZATION

The Suture STRATAFIX™ Symmetric PDS™ Plus is a single-use device provided sterile and sterilized by Ethylene Oxide (EO) to a sterility assurance level (SAL) of 10^{-6} .

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ACCESSORIES

The device in this file does not have accessories included.

No special accessories are required to utilize the suture STRATAFIX™ Symmetric PDS™ Plus.

COMPATIBILITY INFORMATION

The device is utilized with standard surgical instruments, such as graspers/forceps and cutters.

DATA ABOUT PACKAGING AND LABELING PACKAGING

The suture STRATAFIX™ Symmetric PDS™ Plus is wound into a plastic tray packaging which is placed into peel open foil package.

Twelve (12) pouches, a user guide and an IFU are placed into a box. The boxes are wrapped with a film and tear tape is placed during the process to ease the opening of the boxes.

LABELING

Labelling is on the individual primary packaging:

- On the paper lid of the tray
- On the foil package



Labelling is on the group packaging:

- carton box

Labelling is on transport packaging.

For circulation of the product in the Russian Federation a local sticker in the Russian language is added to the original labelling (into group packaging).

Labeling of this medical device includes the following symbols (primary, secondary packaging):

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ENVIRONMENTAL REQUIREMENTS IN THE APPLICATION OF THE MEDICAL DEVICE

This medical device is in use, transportation and storage, has no negative effects on humans and the environment.

PRODUCT EXPIRATION DATE, CRITERIA OF UNSUITABLE CASES OF MEDICAL USE

Shelf life of the product of two (2) years.

Do not use if the packaging is open or damaged.

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Do not use after expiration date.

After expiration date this medical device is to be discarded.

RECYCLING

Discard open packaging and unused devices.

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Dispose the product and its residual elements or waste items in accordance with local laws and regulations.

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After being used the products are potentially infectious waste and are discarded as Class B waste.

Products with expired shelf life and damaged packaging belong to Class A medical waste.

OPERATION, TRANSPORTATION AND STORAGE CONDITIONS

OPERATION CONDITIONS (IN ACCORDANCE WITH THE RUSSIAN LEGISLATION)

The conditions in the rooms where the medical device is used must meet the following requirements:

Temperature from +15°C to 50°C (75% non-condensing humidity equivalent to 30C) when utilized immediately after opening sterile packaging.

STORAGE AND TRANSPORTATION CONDITIONS (IN ACCORDANCE WITH THE RUSSIAN LEGISLATION)

Storage conditions: Store at -60°C up to +30°C

Shipping conditions: Do Not Exceed -60°C or 50°C

Air humidity: packaging in moisture proof foil

Air pressure: N/A

WARRANTY

Manufacturer guarantees to the end user absence of defects and conformity to product specifications for the entire shelf life of the product under the conditions of adherence to storage, transportation and package integrity requirements. There are no warranties that extend beyond the description on the face hereof. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use.

LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS / STANDARDS WHICH THE MEDICAL DEVICE COMPLIES WITH

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