**NOTARIAL CERTIFICATE**

I, **CHOUNG YOUN SIN**, a Notary Public officiating at Klang, Malaysia hereby certify that the **INSTRUCTION FOR USE OF MEDICAL DEVICE SURGICAL NATURAL LATEX GLOVES "BENOVY" STERILE POWDER-FREE** and other documents annexed hereto were presented to me and are to my personal knowledge which information I verily believe to be true belongs to that of **TERANG NUSA (MALAYSIA) SDN. BHD. (COMPANY NO: 224197-U)** of No. 2, Jalan 8, Pengkalan Chepa 2 Industrial Zone, 16100 Kota Bharu, Kelantan, Malaysia.

**IN FAITH AND TESTIMONY WHEREOF I, CHOUNG YOUN SIN** have hereunto set my hand and affixed my seal of office at Klang aforesaid on 22<sup>nd</sup> October 2021.

A handwritten signature in black ink, appearing to read "Choung Youn Sin".

CHOUNG YOUN SIN  
Notary Public  
Klang, Selangor

Attorney General approval reference number PN/189/97 dated 9<sup>th</sup> January 2020 and valid until 13<sup>th</sup> February 2022 under Public Act 1959

TE  
SD  
2, J  
Ind  
Kel

S

VE  
Sem

AP  
Se

|                                                                                                                                           |                                |            |            |
|-------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------|------------|------------|
| <b>TERANG NUSA (MALAYSIA)<br/>SDN. BHD.</b><br>2, Jalan 8, Pengkalan Chepa 2,<br>Industrial Zone, 16100 Kota Bharu,<br>Kelantan, Malaysia | <b>INSTRUCTION FOR<br/>USE</b> | Revision : | 00         |
|                                                                                                                                           |                                | Date :     | 22/06/2021 |

## INSTRUCTION FOR USE

of medical device

### Surgical natural latex gloves “BENOVY” sterile powder-free

VERIFIED BY  
Senior Manager, QA


Ms. Rahayu Bt. Mat Zin

APPROVED BY  
Senior General Manager

Mr. Leong Chew Mun

TERANG NUSA (MALAYSIA) SDN. BHD.  
(Company No. 224197-U)  
2, JALAN 8, PENGKALAN CHEPA 2  
INDUSTRIAL ZONE, 16100 KOTA BHARU,  
KELANTAN, MALAYSIA.  
TEL: +609- 774 7171  
FAX: +609- 774 7757

NOTARISED BY  
  
**CHOUNG YOUN SIN**  
NOTARY PUBLIC  
41050 KLANG,  
SELANGOR.

22 OCT 2021

**Name of the medical device**

**"Surgical natural latex gloves "BENOVY" sterile powder-free"**

in design versions:

1. Surgical natural latex gloves "BENOVY" sterile powder-free with polymer coating microtextured anatomical with curved fingers disposable (beaded cuff), color: natural, size: 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0.
2. Surgical natural latex gloves "BENOVY" sterile powder-free with polymer coating microtextured anatomical with curved fingers disposable with puncture indication in the double glove system (beaded cuff), color: inner glove – green; outer glove – non-glare tan, size: 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0.
3. Surgical natural latex gloves "BENOVY" sterile powder-free with multilayered polymer coating microtextured anatomical with curved fingers disposable with the applicability in the puncture indication system (beaded cuff), color: green, size: 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0.
4. Surgical natural latex gloves "BENOVY" sterile powder-free with multilayered polymer coating microtextured anatomical with curved fingers disposable (beaded cuff), color: natural, size: 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0.
5. Surgical natural latex gloves "BENOVY" sterile powder-free with multilayered polymer coating enhanced tactile sensitivity, including microsurgery, microtextured anatomical with curved fingers disposable (beaded cuff), color: non-glare tan, size: 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0.
6. Surgical natural latex gloves "BENOVY" sterile powder-free with multilayered polymer coating high-strength, including orthopaedics, microtextured anatomical with curved fingers disposable (beaded cuff), color: non-glare brown, size: 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0.

**Name and legal address of the Manufacturer of the medical device**

**Manufacturer: TERANG NUSA (MALAYSIA) SDN. BHD.**

2, Jalan 8, Pengkalan Chepa 2, Industrial Zone, 16100 Kota Bharu, Kelantan, Malaysia

**Data on authorized representative of the Manufacturer**

**"BENOVY" Company Limited**

Russia, 194292, St. Petersburg, 1, 5-y Verkhniy pereulok, litera J, office 205

Tel.: +7 (812) 309-16-32; [www.benovy-company.ru](http://www.benovy-company.ru), e-mail: [benovy@benovy-company.ru](mailto:benovy@benovy-company.ru)

**The scope of delivery of the medical device**

The possible scopes of delivery of the medical device:

- 10 individual packages (pairs) in a consumer package;
- 25 individual packages (pairs) in a consumer package;
- 50 individual packages (pairs) in a consumer package;
- 100 individual packages (pairs) in a consumer package;
- 200 individual packages (pairs) in a consumer package.

**Purpose of the medical device**

The gloves are intended for single use in protecting medical personnel's hands and creating a barrier against bacterial and/or viral contamination between the health worker and the patient during surgical interventions and other sanitary actions that require compliance with sterility conditions.

### **Intended users of the medical device**

Medical personnel, involved in surgical intervention, as well as performing any other sanitary actions that require compliance with aseptic requirements.

### **Scope of the medical device**

Scope of the gloves – all kinds of surgical interventions, general surgery, microsurgery, orthopedics, traumatology, emergency operations and operations with an increased risk of infection, as well as research institutes and veterinary medicine. Gloves are mainly used as a two-sided barrier to protect the patient and the health worker from bacterial and/or viral contamination.

### **Operational principle of the medical device**

The operational principle of the medical device is to create a physical barrier between the skin surfaces of the medical personnel's hands and the surgical field or the other surfaces of the patient, that require compliance with sterility conditions. Thus, the medical device prevents bacterial and/or viral contamination between the medical personnel and the patient and increases the effectiveness of aseptic conditions during surgical interventions, protecting the health worker from possible infection during surgery and other sanitary actions.

The medical device is manufactured without powdering, which eliminates the possible development of a number of postoperative complications associated with the formation of the adhesions in the patient's surgical wound, and also prevents a possible allergic reaction in a health worker.

Medical device design versions allow covering almost all areas of general and specialized surgery, including microsurgery, orthopedics, traumatology, and operations with an increased risk of infection.

The microtextured surface of all design versions of the medical device contributes to the reliable grip of medical instruments and thus increases the probability of a successful surgical outcome.

### **Indications for use of the medical device**

The indications for use of the medical device are the necessity for surgical interventions and other sanitary actions that require compliance with sterility conditions. Gloves represent the personal protective equipment for medical personnel's hands. The use of gloves in healthcare organizations is a mandatory requirement.

### **Contraindications of the medical device**

The contraindication of the medical device is a presence of allergic reaction to natural latex gloves in a patient or a health worker.

### **Side effects**

Health workers may observe the following adverse effects of using the gloves: non-allergic contact dermatitis, allergic contact dermatitis, and immediate allergic reactions accompanied by hives, pruritus and systemic allergic reactions.

Allergic reactions when using the gloves are most often associated with a latex allergy, individual sensitization of a health worker to the chemical components of latex. The manifestation of dermatitis and other skin irritations is possible as a result of improper hand care and harsh conditions of hand skin isolation from natural air ventilation.

If skin irritation occurs, stop using the gloves and consult a doctor.

### Instruction for use of the medical device


Before using the gloves health workers carry out hygienic treatment of their hands having performed surgeon's hand treatment technology.

It is necessary to strictly follow the rules for donning gloves so as not to violate the sterility of their work surface.

Open the individual outer package of sterile gloves (away from the sterile table), remove ("shake out") from it the inner wallet with gloves on a sterile surface. Then open the inner wallet with sterile forceps or aseptically treated hands. Initially, the glove is donned on the dominant (right) hand. Grab the edge of the right glove cuff turned inside out with the thumb and forefinger of the left hand and put it on the right hand without touching the outer sterile surface of the glove. The cuff remains inverted. Next, place fingers of the right hand, on which the glove is already worn, under the cuff lapel of the left glove and don it on the left hand, without touching the skin and the wrong side of the cuff. Spread the cuff of the left glove on the gown in circular motions. The cuff of the glove should overlap the cuff of the gown sleeve by more than 5-10 cm. Only then spread the inverted cuff edge of the right glove in the same circular movements. Pull on the cuff with your four fingers and finally adjust the gloves, ensuring their comfortable use.

Visually the process of surgical gloves donning is shown in Fig.1.

Figure 1. Surgical gloves donning algorithm

There are features of work in surgical gloves.

Hands with sterile gloves should not be lowered below the sterile surfaces of the instrument and operating table. It is not recommended to treat gloves with alcohol-containing and other antiseptic agents before the operation and during surgery, as they increase the porosity/permeability of the glove material and create conditions for overcoming the protective barrier by microorganisms. To maintain an optimal level of protection, the gloves must be changed every 60-120 min, depending on the nature of the surgical intervention. Before donning new gloves, it is necessary to treat the hands with an alcohol-containing antiseptic.

Gloves should be replaced in the following situations:

- after damage (puncture);
- after defect detection;
- after an accidental electric shock from an electrosurgical instrument;
- if any liquid gets (absorbs) under the glove;
- upon appearance of stickiness;
- upon sensation of "glove juice" appearance;
- in the transition from a "dirty" surgical stage to a "clean".


If one glove is damaged, both gloves must be changed. When changing gloves after removing them from the hands, the hands should be treated with a skin antiseptic. Don a new pair of gloves on dry hands.

Remove used gloves carefully, avoiding splashes from the surfaces of the glove, which could cause microbial contamination of the hands and the environment.

To remove the glove, it is necessary to take the left glove cuff from the outside by fingers of the right hand in the glove. Remove the glove from the left hand turning it inside out and leave it holding the lapel in the right hand. Then place the fingers of the free left hand in the glove cuff on the right hand and pull it off similarly to the first one. Having turned the glove inside out, use it as a bag for another glove. Immerse both gloves in a container with a disinfectant solution for disinfection or in a disposable container for temporary waste storage for subsequent centralized disinfection/neutralization.

Visually the process of surgical gloves removing is shown in Fig.2.

Figure 2. Surgical gloves removal algorithm

After finishing the work and gloves removal, it is necessary to carry out hygienic hand treatment. Disposable gloves after use must be decontaminated/neutralized as a medical waste of the appropriate class. To prevent dryness and adverse effects on the skin, it is recommended for health workers to always use professional nourishing hand cream.

### **Application conditions of the medical device**

- Gloves must be used as per the instruction on the medical device.
- The use of the device does not require special qualification or special education.
- Medical gloves are used in medical or healthcare facilities indoors, in a dry place, away from direct sunlight, at a temperature of 15°C to 25°C and relative humidity from 65% up to 85%.
- The device is not intended for reuse.
- Do not use in case of size mismatch.
- For individual use only.

### **Shelf life of the medical device**

The established shelf life of surgical gloves is 5 years.

### **Potential risk for intended use of the medical device:**

class 2a

### **Warnings and precautions**

- The device is intended for single use only.
- For external use only (non-invasive device).
- Every time it is necessary to wash hands before donning or changing gloves, as well as after removing them.
- It is necessary to choose gloves of the right size and depending on the nature of the performed manipulations – they should be easily donned and should not hamper movements or squeeze the hands.
- Touching is permissible only on the limited area of the glove corresponding to the wrist (upper edge of the cuff).
- Sterile. Do not use if individual package is damaged.
- Sterilized by radiation method.
- If the integrity of the individual package is broken, the gloves must be disposed of as General Wastes (Class A Wastes according to the regulatory guidelines of the Russian Federation SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of urban and rural settlements, for water facilities, drinking water and drinking water supply, atmospheric air, soils, living quarters, the operation of industrial, public premises, the organization and implementation of sanitary and anti-epidemic (preventive) measures").
- Do not wet and protect from moisture.
- Avoid exposure to direct sunlight or bright artificial light that has a high UV content.
- Store in a place protected from sunlight. Storage rooms should not contain any equipment capable of producing ozone, such as mercury lamps or high voltage equipment, causing electric sparks or inaudible electric discharges. Combustible gases and organic vapors should also be excluded from the premises.
- The gloves contain natural rubber latex that can cause an allergic reaction.
- The gloves do not guarantee 100% hand protection against punctures, cuts, and damage caused by exposure to chemicals.

### Technical specification of the medical device

**Form.** Surgical gloves have curved fingers, with division into the right and left hands (anatomical), beaded cuff.

**Surface.** The gloves have a microtextured surface over the entire area of palm.

**Dimensions.** Depending on the size and design version of the medical device, tables 1 and 2 indicate the width, length, and thickness of the surgical disposable sterile gloves.

Table 1. Dimensions of the medical device (Design versions Nos. 1 - 5)

| Size | Length $l$ , mm, min | Width $W$ , mm | Thickness, mm, min |      |      |
|------|----------------------|----------------|--------------------|------|------|
|      |                      |                | finger             | palm | cuff |
| 5.5  | 280                  | $72 \pm 4$     | 0.13               | 0.13 | 0.10 |
| 6.0  |                      | $77 \pm 5$     |                    |      |      |
| 6.5  |                      | $83 \pm 5$     |                    |      |      |
| 7.0  |                      | $89 \pm 5$     |                    |      |      |
| 7.5  |                      | $95 \pm 5$     |                    |      |      |
| 8.0  |                      | $102 \pm 6$    |                    |      |      |
| 8.5  |                      | $108 \pm 6$    |                    |      |      |
| 9.0  |                      | $114 \pm 6$    |                    |      |      |

Table 2. Dimensions of the medical device (Design version No. 6)

| Size | Length $l$ , mm, min | Width $W$ , mm | Thickness, mm, min |      |      |
|------|----------------------|----------------|--------------------|------|------|
|      |                      |                | finger             | palm | cuff |
| 5.5  | 280                  | $72 \pm 4$     | 0.20               | 0.20 | 0.15 |
| 6.0  |                      | $77 \pm 5$     |                    |      |      |
| 6.5  |                      | $83 \pm 5$     |                    |      |      |
| 7.0  |                      | $89 \pm 5$     |                    |      |      |
| 7.5  |                      | $95 \pm 5$     |                    |      |      |
| 8.0  |                      | $102 \pm 6$    |                    |      |      |
| 8.5  |                      | $108 \pm 6$    |                    |      |      |
| 9.0  |                      | $114 \pm 6$    |                    |      |      |

The thickness of the bead for all design versions of the medical device is not more than 2.5 mm.

**Physico-mechanical characteristics.** The medical device is manufactured according to the uniform requirements for physico-mechanical characteristics for all design versions. The values of physico-mechanical characteristics before and after accelerated aging of the medical device are presented in table 3.

Table 3. The values of physico-mechanical parameters of all design versions of the medical device

| Characteristics                     | Before accelerated aging | After accelerated aging |
|-------------------------------------|--------------------------|-------------------------|
| Tensile strength, MPa, min          | 24                       | 18                      |
| Elongation at break, %, min         | 750                      | 560                     |
| Stress at 500% elongation, MPa, max | 5.5                      | —                       |
| Force at break, N, min              | 12.5                     | 9.5                     |
| Force at 300% elongation, N, max    | 2.0                      | —                       |

**Inspection levels and acceptable quality level.** For the purposes of quality control, gloves are selected and inspected according to the requirements of ISO 2859-1:1999 standard "Sampling procedures for inspection by attributes – Part 1: Sampling schemes indexed by acceptance quality limit (AQL) for lot-by-lot inspection".

Inspection levels and acceptable quality levels AQL should correspond to the values given in table 4.

Table 4. Inspection levels and acceptable quality levels AQL of the medical device

| Parameter                         | Defect related to parameter                               | Inspection level | AQL |
|-----------------------------------|-----------------------------------------------------------|------------------|-----|
| Watertightness                    | Presence of holes                                         | G1               | 1.5 |
| Dimensions                        | Mismatch of length, width and thickness                   | S-2              | 4.0 |
| Physico-mechanical parameters     | Mismatch of parameters before and after accelerated aging | S-2              | 4.0 |
| Residual powder concentration     | Exceeds maximum value                                     | N=5              | —   |
| Water-extractable protein content | Exceeds recommended maximum value                         | N=3              | —   |
| Antigenic protein content         | Exceeds recommended maximum value                         | N=1              | —   |

**Watertightness.** Gloves must be watertight. The acceptable quality level (AQL) is  $\leq 1.5$ .

**Powder content.** In powder-free gloves, the residual powder content does not exceed 2 mg per one glove.

**Water-extractable protein content in latex gloves.** The content of water-extractable proteins is not more than  $200 \mu\text{g}/\text{dm}^2$ .

**Antigenic protein content.** The content of antigenic proteins is not more than  $10 \mu\text{g}/\text{dm}^2$ .

**Resistance to viral penetration.** The medical device is resistant to viral penetration.

### Transportation, storage and disposal conditions of the medical device

#### Transportation conditions

Transportation of the device can be carried out by all types of covered vehicles as per the transportation rules applied to this type of transport. During transportation the storage conditions must not be violated.

#### Storage conditions

The devices must be stored in a dry place, protected from direct sunlight or bright artificial light having a high ultraviolet content, at a temperature of  $0^\circ\text{C}$  to  $25^\circ\text{C}$  and relative air humidity – from 65 % to 85%.

The devices should be handled with caution when they are stored at temperature below  $15^\circ\text{C}$ , as during careless handling, they can be stiff and sensitive to deformation. After storage at a low temperature, gloves should be warmed up to  $30^\circ\text{C}$  before use throughout all device volume.

Storage rooms should not contain ozone-generating equipment, such as mercury lamps or high voltage electrical equipment, generating excessive sparks or inaudible electric discharges.

The devices must be stored away from sources of ignition. Combustible gases and organic vapors should be excluded from the storage facilities.

### **Disposal conditions**

The used medical device belongs to the category of Infectious Wastes (Class B Wastes according to the regulatory guidelines of the Russian Federation SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of urban and rural settlements, for water facilities, drinking water and drinking water supply, atmospheric air, soils, living quarters, the operation of industrial, public premises, the organization and implementation of sanitary and anti-epidemic (preventive) measures").

The collection, temporary storage, and removal of waste should be performed as per the medical waste management scheme adopted by the organization that carries out healthcare and/or pharmaceutical practice.

### **Manufacturer's warranties**

The manufacturer guarantees the preservation of all parameters and characteristics of the gloves throughout the shelf life (5 years) upon the condition of package integrity, observance of storage, transportation and operation conditions in accordance with the instruction and acceptable quality level  $AQL \leq 1.5$ .

### **List of national standards, applied by manufacturer (producer) of the medical device**

|                  |                                                                                                                                                                                                         |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| EN 455-1:2000    | Medical gloves for single use. Part 1. Requirements and testing for tightness                                                                                                                           |
| EN 455-2:2015    | Medical gloves for single use. Part 2. Requirements and testing for physical properties                                                                                                                 |
| EN 455-3:2015    | Medical gloves for single use. Part 3. Requirements and testing for biological evaluation                                                                                                               |
| EN 455-4:2009    | Medical gloves for single use. Part 4. Requirements and testing for shelf life determination                                                                                                            |
| EN 556-1:2002    | Sterilization of medical devices. Requirements for medical devices to be designated "sterile". Part 1. Requirements for terminally sterilized medical devices                                           |
| EN 868-5:2019    | Packaging materials and systems for medical devices which are to be sterilized. Part 5: Heat and self-sealable pouches and reels of paper and plastic film construction – Requirements and test methods |
| ISO 2859-1:1999  | Sampling procedures for inspection by attributes – Part 1: Sampling schemes indexed by acceptance quality limit (AQL) for lot-by-lot inspection                                                         |
| ISO 10012:2003   | Measurement management systems – requirements for measurement processes and measuring equipment                                                                                                         |
| ISO 11137-1:2015 | Sterilization of health-care products – Radiation – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices                             |

|                   |                                                                                                                                                                                |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ISO 11137-2:2015  | Sterilization of health-care products – Radiation – Part 2: Establishing the sterilization dose                                                                                |
| ISO 11137-3:2017  | Sterilization of health-care products – Radiation – Part 3: Guidance on dosimetric aspects of development, validation and routine control                                      |
| ISO 11139:2018    | Sterilization of health care products – Vocabulary of terms used in sterilization and related equipment and process standards                                                  |
| ISO 10993-1:2018  | Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process                                                                     |
| ISO 10993-10:2010 | Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization                                                                                |
| ISO 10993-12:2012 | Medical devices. Biological evaluation of medical devices – Part 12: Sample preparation and reference materials                                                                |
| ISO 11607-1:2019  | Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems                                        |
| ISO 11607-2:2019  | Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes                                              |
| ISO 11737-1:2018  | Sterilization of medical devices – Microbiological methods – Part 1: Determination of a population of microorganisms on products                                               |
| ISO 11737-2:2019  | Sterilization of medical devices – Microbiological methods – Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process     |
| ISO 13485:2016    | Medical devices – Quality management systems – Requirements for regulatory purposes                                                                                            |
| ISO 14971:2007    | Medical devices – Application of risk management to medical devices                                                                                                            |
| ISO 15223-1:2016  | Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements                                        |
| ASTM D412         | Standard Test Methods for Vulcanized Rubber and Thermoplastic Elastomers – Tension                                                                                             |
| ASTM D573         | Standard Test Method for Rubber – Deterioration in an Air Oven                                                                                                                 |
| ASTM D3577        | Standard Specification for Rubber Surgical Gloves                                                                                                                              |
| ASTM D3767        | Standard Practice for Rubber – Measurement of Dimensions                                                                                                                       |
| ASTM D5151        | Standard Test Method for Detection of Holes in Medical Gloves                                                                                                                  |
| ASTM D5712        | Standard Test Method for Analysis of Aqueous Extractable Protein in Natural Rubber and its Products Using the Modified Lowry Method                                            |
| ASTM D6124        | Standard Test Method for Residual Powder on Medical Gloves                                                                                                                     |
| ASTM F88          | Standard test method for seal strength of flexible barrier materials                                                                                                           |
| ASTM F1671        | Standard test method for resistance of materials used in protective clothing to penetration by blood-borne pathogens using Phi-X174 bacteriophage penetration as a test system |

**Instruction for use of the medical device “Surgical natural latex gloves “BENOVY” sterile powder-free” is released for the first time (Rev: 00 (June 2021)).**

In case of an undesirable event showing the signs of an adverse event (incident), it is necessary to notify thereof the manufacturer or an authorized representative of the manufacturer in the Russian Federation.

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