
CODMAN NEURO

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January 9, 2020

Letter of Declaration


As a responsible representative of **Codman & Shurtleff, Inc.**, I, **Vivian Perez**, hereby certify and attest that the attached document, IFU Codman Microsensor Basic Kit - Variant I - REF 626631, is a true, exact, unaltered copy of the original document issued by **Codman & Shurtleff, Inc.**, located at 325 Paramount Drive, Raynham, MA 02767-0350, USA.

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
INSTRUCTIONS FOR USE
Medical Device
CODMAN MICROSENSOR Kit
Variant I. CODMAN MICROSENSOR Basic Kit
REF 626631

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Instructions for use
CODMAN MICROSENSOR Basic Kit

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STERILE EO

Description

The CODMAN MICROSENSOR Basic Kit, catalog no 626631, consists of the CODMAN MICROSENSOR ICP Transducer and a 14-gauge Tuohy needle with stylet. See Figure 1.

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The CODMAN MICROSENSOR Basic Kit is designed for use with the CODMAN Cranial Hand Drill, catalog no. 82-6607. The drill facilitates access to the intraparenchymal area. The drill is also available as a component of the CODMAN Cranial Access Kit, catalog no. 82-6612, 82-6614, and 82-6616.

Intended Use

CODMAN MICROSENSOR Basic kit is designed for use for direct measurement of intracranial pressure.

Indications

Use of the CODMAN MICROSENSOR Basic Kit is indicated when direct ICP monitoring is required. The kit is indicated for use in both subdural and intraparenchymal pressure monitoring applications only.

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Exposure to electrostatic discharge (ESD) energy could damage this device. High levels of ESD could damage the electronic components and cause the transducer to be rendered inaccurate or inoperable. Take all precautions to reduce the buildup of electrostatic charge during the use of this product and avoid touching the transducer connector pins, which are identified with the ESD symbol. (Refer to *Electrostatic Discharge (ESD) Information* section). Take care when tying sutures onto the sensor. Tying sutures too tightly can collapse the wall of the sensor body, causing damage to internal wires.

The transducer must be zeroed at atmospheric pressure prior to implantation. The transducer tip must remain wet during the zeroing process. Do not submerge the tip of the transducer or catheter vertically in a deep pool or cup of sterile water/sterile saline. Doing so will impose a hydrostatic pressure on the transducer diaphragm that is higher than atmospheric zero, resulting in an inaccurate zero reference.

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Do not forcibly pull or jerk the transducer catheter.

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patient with an implanted CODMAN

MICROSENSOR. Failure to adhere to the Conditions for Safe Use may result in serious injury to the patient.

The CODMAN MICROSENSOR is MR Conditional

MRI safety information:

Non-clinical testing demonstrated that the CODMAN MICROSENSOR is MR Conditional. A patient implanted with this device can be safely scanned in an MR system which meets or is operated under the following conditions:

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Preparation for the MRI procedure:

1. Immediately prior to entering the MRI environment, verify that the CODMAN MICROSENSOR is functioning properly. DO NOT perform an MRI procedure if the CODMAN MICROSENSOR is damaged or otherwise not functioning properly.
2. Disconnect all cables and patient monitoring devices attached to the CODMAN MICROSENSOR prior to transporting the patient into the MRI environment. DO NOT bring the ICP EXPRESS Monitor, cables or other accessories such as the trocar or stylet into the MRI environment.
3. Special positioning of the CODMAN MICROSENSOR is required to ensure patient safety during the MRI procedure. The CODMAN MICROSENSOR must be placed in a specific geometry to minimize the potential for excessive heating of the sensor tip. Coil the tubing of the CODMAN MICROSENSOR near the base of the electrical connector into 5 or 6 loops approximately 6cm in diameter and center on top of the patient's head (see Figure 7). Do not perform MRI with the CODMAN MICROSENSOR in a "straight line" configuration (i.e., uncoiled). Failure to follow this guideline can result in serious injury to the patient.
4. Insert a dry gauze pad at least 1cm thick between the CODMAN MICROSENSOR electrical connector with coiled tubing and the patient's scalp. Secure in place using tape (see Figure 7). Use care when removing the tape to prevent damage to the CODMAN MICROSENSOR.

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All components have been tested and were determined to be nonpyrogenic, **except** for the electrical connector of the sensor and the silicone tubing used for packaging, which are not tested. The device is packaged using a combination of recyclable and non-recyclable materials. Recycle or dispose of all packaging waste in accordance with hospital procedures and regulations.

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1. Connect the transducer to a compatible CODMAN ICP Monitor using an appropriate CODMAN interface cable. Refer to instructions for use provided with the cable for sterilization information.

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3. Zero and calibrate the external patient monitor according to the instructions provided with the CODMAN ICP monitor, as well as the external patient monitor manufacturer's instructions.

4. Prepare to zero the catheter by laying the tip of the catheter flat in a shallow pool of sterile water or sterile saline. The accompanying sterile blister package has a marked well that is suitable for this procedure. Pour sufficient sterile water or sterile saline into the well, then lay at least a 5cm section of the catheter horizontally just under the surface of the sterile water or sterile saline.

CAUTION: Do not submerge the tip of the catheter vertically in a deep pool or cup of sterile water or sterile saline. Doing so will impose a

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hydrostatic pressure on the sensor diaphragm that is higher than atmospheric pressure, resulting in an inaccurate zero reference.

5. While keeping the tip of the catheter flat and still in the sterile water or sterile saline, zero the catheter according to the instructions provided with the ICP Monitor.

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CAUTION: The sensor tip must remain still during the zeroing process.

Motion of the sensor may be interpreted by the CODMAN ICP Monitor as a fluctuating ICP signal which will prevent the sensor zeroing process from successfully completing.

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Measuring Subdural Pressure

1. Following craniotomy and bone flap removal connect and zero the transducer. Refer to *Connecting and Zeroing the Transducer*.
2. Choose the burr hole through which the transducer will be placed and bevel the edge on the side the transducer will exit to facilitate removal.

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3. Use the Tuohy needle to tunnel under the scalp from the craniotomy site to the desired transducer exit site.

4. Remove the Tuohy needle stylet and thread the transducer from the tip of the needle until the approximate length for desired placement exits from the hub. See Figure 2.

CAUTION: The inner edges of the Tuohy needle are sharp; exercise caution when threading the catheter.

5. Gently remove the needle and estimate the length of the transducer from the tip to the first bend.

6. Fold back the transducer once completed at the desired bend site to leave a kink in the catheter. Verify that the kink is on the opposite side of the transducer from the transducer sensing element, as shown in Figure 3.

7. Place the tip of the transducer on the brain tissue under the dura opposite the beveled burr hole. The kink must be placed at the bottom of the burr hole so that the transducer sensing element faces the cortex, as shown in Figure 4. **CAUTION: To ensure accurate ICP measurements, position the sensing element of the transducer tip towards the cortex during subdural pressure monitoring.**

CAUTION: To minimize movement artifacts of the bone flap, place the transducer under the intact skull.

8. Close and suture down the dura following standard neurosurgical procedure.

9. Replace the bone flap and close the scalp incision.

10. Secure the catheter to the scalp. To provide additional strain relief, make a small loop with the catheter and tie the loop down.

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Measuring Intraparenchymal Pressure

NOTE: The CODMAN Cranial Hand Drill, catalog no. 82-6607, is recommended for this procedure.

1. Connect and zero the transducer. Refer to Connecting and Zeroing the MICROSENSOR Transducer.
2. Perform craniotomy and retraction procedures required to expose the skull. Put the 2.7 mm drill bit into the drill chuck and make a drill hole through the outer table. **CAUTION: Proceed gently through the inner table with care to avoid injury to the dura or parenchyma.**
3. Carefully bevel the bone at the incision site towards the side from which the catheter will exit.
4. Use the Tuohy needle to make a cruciate puncture in the dura.
5. Use the Tuohy needle with stylet in place to tunnel under the scalp from the incision site to the exit site.
6. Remove the stylet and thread the transducer from the tip of the needle until approximately twice the placement length exits from the hub. See Figure 5. **CAUTION: The inner edges of the Tuohy needle are sharp; exercise caution when threading the transducer.**
7. Gently remove the needle and estimate the length of the transducer from the tip to where it will bend upon exiting the skull.
8. Fold the transducer once completely at the desired bend site to leave a kink in the transducer.
9. Place the tip of the transducer into the parenchyma through the puncture in the dura until the kink is at the edge of the hole, as shown in Figure 6.
10. Carefully pull back the excess transducer.
11. Remove the retractor, verify hemostasis in the insertion area, and then suture the incision site closed.
12. Secure the catheter to the scalp. For additional strain relief, make a small loop with the catheter and tie the loop down.
13. Close and dress the incision site.

Package contents

One consumer package included:

1. CODMAN MICROSENSOR Basic Kit for intracranial pressure measurement consisting of:
 - ICP Transducer – 1 pcs
 - Tuohy needle – 1 pcs
 - Stylet – 1 pcs
2. Instructions for use – 1 pcs

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- ICP Transducer

Technical characteristics

Sensing element	Strain gauge silicon microchip
Probe usable length	100 cm nominal
Probe tip diameter	1,3 mm maximum
Probe tubing diameter	0,8 mm maximum
Probe tubing length	105 $\pm$ 2,5 cm
Sensor housing dimensions	51 x 17 x 9 mm $\pm$ 0,1 mm
Functional pressure range	-50 mmHg to +250 mmHg (- 6666 to +33330 Pa)
Functional overpressure range without damage	-700 mmHg to +1250 mmHg (- 93324 to +166650 Pa)
Input/output impedance	1000 ohms nominal
Zero drift	No greater than 5 mmHg over 30 days
Output signal (sensitivity)	5 μ V/V/mmHg nominal
Frequency response	Greater than 200 Hz
Markings of probe tubing	Thin markings (16 places) dimensions 0,6 mm $\pm$ 0.1 mm thick Thick markings centered at 5 cm intervals from the tip (10 markings) dimensions 1,1 mm $\pm$ 0.1 mm thick. One mark @5 cm, two marks @10 cm, three marks @15 cm and four marks @20 cm. Spacing between thick markings to be 0,5 mm $\pm$ 0.1 mm.
Length from the tip to the first mark	9.7 $\pm$ 1 mm.
Pull force	Joint/bonds – 2,2 N Epoxy cap/ball – 1 N

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- Tuohy needle

Technical characteristics

Ø _{ext}	2,11± 0,03 mm (14-gauge)
Ø _{int}	1,83± 0,03 mm
Length of hub	34.29± 0,38 mm
Length of needle	86,36± 1,52 mm
Wing dimensions	26,92± 0,38 x 12,19± 0,38 x 3,048± 0,13 mm
Markings thick	10 mm, the marks are located every 10 mm from the end of the needle
Type of connector	Female luer
Hardness	52-54 RC
Roughness	0,4 µm

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Environmental Specifications

For non-implantable portion of device

Operating temperature range	5 °C to 45 °C
Operating humidity range	30% to 90% relative humidity (non-condensing)
Operating atmospheric pressure range	700 millibar to 1060 millibar

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For implantable portion of device

This portion intended to be used in the body's internal environment.

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Transportation and Storage Conditions

Temperature range	5 °C to 50 °C
Humidity range	25% to 80% relative humidity (noncondensing)
Atmospheric pressure range	700 millibar (hPa) to 1060 millibar (hPa)

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Sterility

This device is intended for SINGLE USE ONLY; DO NOT RESTERILIZE.
Sterilization method used: Ethylene Oxide.

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The CODMAN MICROSENSOR ICP Transducer is designed for use with the ICP

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EXPRESS Monitor, catalog no. 82-6634 (117 VAC) and 82-6635 (230 VAC) (Russian Registration Certificate N RZN 2013/952 dated August 07, 2013) or any other suitable CODMAN pressure monitoring interface or device.

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Quality Systems	
ISO 13485: 2003	Medical devices -- Quality management systems -- Requirements for regulatory purposes
Risk	
ISO 14971: 2012	Risk Management for Medical Devices
Product Safety and Performance Standards	
ASTM F2052-14	Standard Test Method for Measurement of Magnetically Induced Displacement Force on Medical Devices in the Magnetic Resonance Imaging
ASTM F2182-11a	Standard Test Method for Measurement of Radio Frequency Induced Heating On or Near Passive Implants During Magnetic Resonance Imaging
ASTM F2119-07	Standard Test Method for Evaluation of MR Image Artifacts from Passive Implants
ASTM F2213-06	Standard Test Method for Measurement of Magnetically Induced Torque on Medical Devices in the Magnetic Resonance Environment
ASTM F2503-13	Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment
IEC 60601-1	Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2	Collateral Standard: Electromagnetic disturbances - Requirements and tests
IEC 60601-1-6	Collateral standard: Usability
IEC 60601-1-9	Collateral Standard: Requirements for environmentally conscious design

EN ISO 14630: 2012	Non-active surgical implants. General requirements (components other than the MICROSENSOR probe)
Clinical	
MedDev 2.7.1: 2009	Guidelines on medical devices-Evaluation of clinical data. A guide for manufacturers and Notified Bodies
Biocompatibility	
ISO 10993- 1: 2009/AC2010	Biological Evaluation of Medical Devices: Part 1: Evaluation and testing with a risk management process
Packaging	
ISO 11607: 2009 Part 1	Packaging for terminally sterilized medical devices Requirements for materials, sterile barrier systems and packaging systems
ISO 11607:2009 Part 2	Packaging for terminally sterilized medical devices Validation requirements for forming, sealing and assembly processes
Labeling	
EN 980: 2008	Graphic Symbols for use in labeling of medical devices
EN 1041:2008	Information Supplied by Manufacturer of Medical Devices
EN ISO 15223-1:2012	– Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied
IEC TR 60878:2003	– Graphical symbols for electrical equipment in medical practice
Sterilization	
ISO 11135: 2007	Validation and Routine Control of Ethylene Oxide Sterilization
ISO 10993-7:2008	Biological Evaluation of Medical Devices – Part 7: Ethylene Oxide Sterilization Residuals
ISO 11737-2:2009	Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
ANSI ST72:2011	Bacterial endotoxins — Test methods, routine monitoring, and alternatives to batch testing

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Warranty

Codman & Shurtleff, Inc., warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.** Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.


Figure 1

A. Tuohy Needle with Stylet

B. MICROSENSOR ICP Transducer with cm depth markings

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
Figure 3

CAUTION: Avoid direct contact with the sensing element of the transducer.

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
Figure 4
Subdural Placement under Intact Skull
A.Dura
B.Subdural Space

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
Figure 7 Transducer configuration to ensure patient safety during MRI

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
January 9, 2020

Letter of Declaration


As a responsible representative of **Codman & Shurtleff, Inc.**, I, **Vivian Perez**, hereby certify and attest that the attached document, Instructions For Use (IFU) Codman Microsensor Plastic Skull Bolt Kit – 121119 - Variant II - Ref 626632, is a true, exact, unaltered copy of the original document issued by **Codman & Shurtleff, Inc.**, located at 325 Paramount Drive, Raynham, MA 02767-0350, USA.

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INSTRUCTIONS FOR USE
Medical Device
CODMAN MICROSENSOR Kit
Variant II. CODMAN MICROSENSOR Plastic Skull Bolt Kit
REF 626632

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**Instructions for Use
CODMAN MICROSENSOR Plastic Skull Bolt Kit**

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IMPORTANT INFORMATION

Please Read Before Use

**CODMAN MICROSENSOR
Plastic Skull Bolt Kit**

REF 626632

STERILE EO

Rx Only

Description

The CODMAN MICROSENSOR Plastic Skull Bolt Kit, catalog no. 626632, consists of the CODMAN MICROSENSOR ICP Transducer with a Tuohy Borst adapter mounted on the catheter, one Adult and one Pediatric Skull Bolt with spacing washers, and a Disposable ICP Screw Handle for facilitating bolt removal (see Figure 1). The CODMAN MICROSENSOR ICP Transducer is a catheter with a microminiature strain gauge sensor mounted at one end and an electrical connector at the other end. It is designed for use with the ICP EXPRESS Monitor, catalog no. 82-6634 (117 VAC) and 82-6635 (230 VAC) or any other suitable CODMAN pressure monitoring interface or device.

The single lumen plastic skull bolt is designed to be used with the CODMAN Cranial Hand Drill, catalog no. 82-6607. The Cranial Hand Drill is also available as a component of the CODMAN Cranial Access Kit, catalog no. 82-6612, 82-6614, and 82-6616.

CAUTION: Bolt fit can be compromised if used with other drills.

The skull bolts are preassembled with one thrust washer and one sealing washer. Skull bolts are supplied with spacing washers, male LUER-LOK caps, and obturators. See Figure 2.

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Page 4 of 31

Intended Use

CODMAN MICROSENSOR Plastic Skull Bolt Kit is designed for use for direct measurement of intracranial pressure.

Indications

Use of the CODMAN MICROSENSOR Plastic Skull Bolt Kit is indicated when direct intracranial pressure (ICP) monitoring is required. The kit is indicated for use in both subdural and intraparenchymal pressure monitoring applications.

Contraindications

Insertion of the skull bolt is contraindicated in children less than one year of age.

Note: This kit is not designed, sold, or intended for any use except as indicated. This kit is not designed, sold, or intended for use as a therapeutic device.

Conditions of Use

The medical device is intended for use at the standard conditions of hospital operating room environment. All procedures should be performed by a qualified neurosurgeon using standard surgical procedure. Neurosurgeon is responsible for the implementation of appropriate steps and procedures to avoid infection and complications.

Warnings

Installation of the skull bolt must be performed with the bolt held within 10

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degrees of the perpendicular to the incision site. Installing at an angle can result in fracture of the device or in an inadequate seal between the washer and the skull.

Do not remove the thrust washer or the sealing washer that are preassembled on the skull bolt, as this can result in an inadequate seal between the washer and the skull.

Excessive torque applied to the skull bolt during insertion can cause breakage. This device is designed to break off at the area between the hexagonal, if excessive torque is applied. The lower portion of the screw can then be removed, using the screw handle.

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Before conducting an MRI procedure on a patient with an implanted CODMAN MICROSENSOR Transducer, read the MRI Safety section. Failure to read and strictly adhere to these guidelines can result in serious injury to the patient.

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It is essential to maintain strict sterile technique during bolt application and transducer placement

Exposure to electrostatic discharge (ESD) energy could damage this device. High levels of ESD could damage the electronic components and cause the transducer to be rendered inaccurate or inoperable. Take all

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precautions to reduce the buildup of electrostatic charge during the use of this product and avoid touching the transducer connector pins, which are identified with the ESD symbol. (Refer to *Electrostatic Discharge (ESD) Information* section). Bolt fit can be compromised if used with other drills.

The transducer must be zeroed at atmospheric pressure prior to implantation. The transducer tip must remain wet during the zeroing process.

Do not submerge the tip of the transducer or catheter vertically in a deep pool or cup of sterile water/sterile saline. Doing so will impose a hydrostatic pressure on the transducer diaphragm that is higher than atmospheric zero, resulting in an inaccurate zero reference.

The transducer can be damaged if exposed to pressures over 1250 mmHg (166 650 Pa).

Do not forcibly pull or jerk the transducer catheter. Excessive force can result in fracture of the transducer or unintended withdrawal from the skull bolt

Do not expose the transducer to solvents or cleaning agents, including alcohol; these may cause damage leading to inaccurate ICP measurements.

Read all instructions included with the monitoring display device prior to use.

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MRI Information

MR Conditional

Read and understand this document in its entirety prior to performing a Magnetic Resonance

Imaging Procedure on a patient with an implanted CODMAN MICROSENSOR. Failure to adhere to the Conditions for Safe Use may result in serious injury to the patient.

The **CODMAN MICROSENSOR** is MR Conditional

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- Static magnetic field of 1.5 and 3 tesla only.
- Maximum spatial gradient magnetic field of 1000 G/cm (10 T/m).
- Maximum gradient field slew rate of 170 T/m/s.
- Horizontal cylindrical bore MRI scanner.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode), for 15 minutes of continuous scanning.
- Special positioning of the CODMAN MICROSENSOR is required to ensure patient safety during the MRI procedure (see "PREPARATION FOR THE MRI PROCEDURE" below for specific instructions).
- WARNING: Do not bring the ICP EXPRESS monitor, cables or other accessories such as Tuohy needles, trocar or stylet into the MRI suite.
- WARNING: Do not use Transmit / Receive or Transmit-only RF Head coils. Only use Transmit / Receive RF Body coil or Transmit RF Body coil / Receive-only RF Head coil.
- WARNING: Do not scan a patient with an elevated body temperature.

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Preparation for the MRI procedure:

- 1.Immediately prior to entering the MRI environment, verify that the CODMAN MICROSENSOR is functioning properly. DO NOT perform an MRI procedure if the CODMAN MICROSENSOR is damaged or otherwise not functioning properly.
- 2.Disconnect all cables and patient monitoring devices attached to the CODMAN MICROSENSOR prior to transporting the patient into the MRI environment. DO NOT bring the ICP EXPRESS Monitor, cables or other accessories such as the trocar or stylet into the MRI environment.
- 3.Special positioning of the CODMAN MICROSENSOR is required to ensure patient safety during the MRI procedure. The CODMAN MICROSENSOR must be placed in a specific geometry to minimize the potential for excessive heating of the sensor tip. Coil the tubing of the CODMAN MICROSENSOR near the base of the electrical connector into 5 or 6 loops approximately 6 cm in diameter and center on top of the patient's head (see Figure 8). Do not perform MRI with the CODMAN

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MICROSENSOR in a “straight line” configuration (i.e., uncoiled).

Failure to follow this guideline can result in serious injury to the patient.

4. Insert a dry gauze pad at least 1 cm thick between the CODMAN MICROSENSOR electrical connector with coiled tubing and the patient’s scalp. Secure in place using tape (see Figure 8). Use care when removing the tape to prevent damage to the CODMAN MICROSENSOR.

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- a. Maximum spatial gradient magnetic field of 1,000 G/cm (10 T/m). The highest SG magnetic field is commonly located off-axis, at a side wall, and near the opening of the bore of the scanner. Please refer to MRI manufacturer’s published value and location of the peak SG that is accessible to the patient.
- b. Maximum gradient field slew rate of 170 T/m/s.
- c. Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode), for 15 minutes of continuous scanning.

Electrostatic Discharge (ESD) Information

CAUTION: Exposure to electrostatic discharge (ESD) energy could damage this device. High levels of ESD could damage the

electronic components and cause the transducer to be rendered inaccurate or inoperable. Take all precautions to reduce the buildup of electrostatic charge during the use of this product.

- Provide patient grounding; e.g., grounding straps on gurneys
- Avoid the use of materials that could generate ESD during patient movement and transport; e.g., nylon transfer boards with bedding

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- Before touching the patient, caretakers should discharge ESD buildup by touching a grounded metal surface, such as a bedrail

Provide patient grounding; e.g., grounding straps on gurneys

Avoid the use of materials that could generate ESD during patient movement and transport; e.g., nylon transfer boards with bedding

Before touching the patient, caretakers should discharge ESD buildup by touching a grounded metal surface, such as a bed rail. It is recommended that all hospital personnel in contact with these devices receive an explanation of the ESD symbol and training in ESD precautionary procedures. Training should include, at a minimum, an introduction to electrostatic discharge, when and why it occurs, precautionary measure, and the damage that can be done to electronic components if touched by a user who is electrostatically charged.

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Codman Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or re-sterilization, after a single patient use. These devices are intended to come into contact with the central nervous system and the ability does not currently exist to destroy possible contaminants such as Creutzfeldt-Jakob Disease. Reuse can also compromise device performance and any usage beyond the design intent of this single-use device can result in unpredictable use hazards or loss of functionality.

Codman & Shurtleff will not be responsible for any product that is re-sterilized, nor

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accept for credit or exchange any product that has been opened but not used.

As long as the individual package is not damaged or opened, the product is sterile. All components have been tested and were determined to be nonpyrogenic, except for the electrical connector of the sensor, the silicone tubing used for packaging, and the disposable ICP handle, which are not tested. The device is packaged using a combination of recyclable and non-recyclable materials. Recycle or dispose of all packaging waste in accordance with hospital procedures and regulations.

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Connecting and Zeroing the transducer

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1.Connect the transducer to the CODMAN ICP Monitor using an appropriate CODMAN interface cable. Refer to instructions for use provided with the cable for sterilization information.

2.If applicable, connect the CODMAN ICP Monitor to an available pressure channel on an external patient monitor using a CODMAN Patient Monitor Interface Cable.

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3.Zero and calibrate the external patient monitor according to the instructions provided with the CODMAN ICP Monitor, as well as the external patient monitor manufacturer's instructions.

4.Prepare to zero the catheter by laying the tip of the catheter flat in a shallow pool of sterile water or sterile saline. The accompanying sterile blister package has a marked well that is suitable for this procedure. Pour sufficient sterile water or sterile saline into the well, then lay at least a 5 cm section of the catheter horizontally just under the surface of the sterile water or sterile saline.

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5.While keeping the tip of the catheter flat and still in the sterile water or sterile saline, zero the catheter according to the instructions provided with the ICP EXPRESS Monitor.

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CAUTION: The sensor tip must remain still during the zeroing process. Motion of the sensor may be interpreted by the CODMAN ICP Monitor as a fluctuating ICP signal which will prevent the sensor zeroing process from successfully completing.

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6. If available, record the three-digit zero reference number provided by the CODMAN ICP Monitor. Mark this number on the catheter connector housing and the patient's chart for future reference.

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Installing the Skull Bolt

1. Connect and zero the transducer. Refer to Connecting and Zeroing the MICROSENSOR Transducer.

2. Perform craniotomy and retraction procedures required to expose the skull. Use the 5.8 mm drill bit when drilling hole.

CAUTION: Bolt fit can be compromised if used with other drills. Make a cruciate puncture in the dura to establish an assured path between the bolt and the intraparenchymal area.

3. Select adult or pediatric bolt. Note: if necessary, adjust the seating depth of the bolt using the spacing washers provided. Two washers can be used with the adult skull bolt and one can be used with the pediatric skull bolt to adjust for variations in skull thickness in patients. Position spacing washers prior to insertion of the bolt. The range of adjustment available to the surgeons:

Pediatric: 3.5 mm to 5 mm

Adult: 5.0 mm to 10 mm

WARNING: Do not remove the thrust washer and sealing washer that are preassembled on the skull bolt, as this

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can result in a fracture of the device or in an inadequate seal between the washer and the skull.

Note: Do not remove the obturator from the skull bolt.

4. Put the skull bolt into position and place the screw handle on top of the bolt as shown in Figure 4. Turn the screw handle clockwise until the spacing washer beneath the screw flange rests close against the outer table of the skull. Remove the screw handle.

WARNING: Excessive torque applied to the skull bolt during insertion can cause breakage. This device is designed to break off at the area between the hexagonal, if excessive torque is applied. The lower portion of the screw can then be removed using the screw handle. To remove the bolt, turn the screw handle counterclockwise.

5. Remove the obturator from the skull bolt.

6. Guide the transducer through the bolt until the tip is in the subdural or parenchymal space, as shown in Figure 5.

7. Irrigate the channel with non-bacteriostatic sterile saline.

8. Loosen the Tuohy-Borst adapter from the catheter body by twisting the blue adapter counterclockwise. Slide the Tuohy-Borst adapter into place as shown in Figure 6. Turn the blue adapter clockwise to secure.

9. Remove the retractor and close the skin incision in one layer using interrupted sutures, as shown in Figure 7. Apply appropriate dressing to the wound site.

Note: If desired, the male LUER-LOK cap may be used to occlude the bolt outlet when monitoring is discontinued.

Removing the Skull Bolt

1. Remove all connections from the skull bolt.
2. Remove the dressing around the bolt and one or two sutures closest to the stem.

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3. Using the sterile screw handle provided, turn the screw counterclockwise until free.

CAUTION: Make certain that all the washers have been removed from the incision area before closing.

4. Apply pressure with sterile gauze until hemostasis is achieved. Close the incision again in one layer with uninterrupted sutures. Apply a dry sterile dressing.

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Package contents

One consumer package included:

1. CODMAN MICROSENSOR plastic Skull Bolt Kit for intracranial pressure measurement consisting of:
 - ICP Transducer – 1 pcs
 - Tuohy-Borst adapter – 1 pcs
 - Skull Bolt Adult – 1 pcs
 - Skull Bolt Pediatric – 1 pcs
 - Luer-lok Cap – 2 pcs
 - Disposable ICP Handle – 1 pcs
 - Spacing Washer – 3 pcs
2. Instructions for use – 1 pcs
3. Patient labels – 6 pcs

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Functional overpressure range without damage	-700 mmHg to +1250 mmHg (- 93324 to +166650 Pa)
Input/output impedance	1000 ohms nominal
Zero drift	No greater than 5 mmHg over 30 days
Output signal (sensitivity)	5 μ V/V/mmHg nominal
Frequency response	Greater than 200 Hz
Markings of probe tubing	Thin markings (16 places) dimensions 0,6 mm $\pm$ 0.1 mm thick
	Thick markings centered at 5 cm intervals from the tip (10 markings) dimensions 1,1 mm $\pm$ 0.1 mm thick. One mark @5 cm, two marks @10 cm, three marks @15 cm and four marks @20 cm. Spacing between thick markings to be 0,5 mm $\pm$ 0.1 mm.
Length from the tip to the first mark	9.7 $\pm$ 1 mm.
Pull force	Joint/bonds – 2,2 N Epoxy cap/ball – 1 N

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- Tuohy-Borst adapter

Technical characteristics

Housing length	23,62 ± 0,51 mm
Housing type of connectors	At inlet – male luer At outlet - tapered adapter
Head Length	11,94 ± 0,13 mm
Øint	3,05 ± 0,18 mm
Type of connector	male luer

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- Skull Bolt, Adult

Technical characteristics

Bolt	
Length	25,44 ± 0,38 mm
Length of the threaded portion	10,41 ± 0,13 mm
Diameter	10,92 ± 0,13 mm
Type of connector	Female luer
Obturator	
Total length	35 mm
Puncturer length	21 mm
Male luer length	8 mm
Type of connector	Male luer

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- Skull Bolt, Pediatric

Technical characteristics

Bolt	
Length	24,63±0.38 mm
Length of the threaded portion	6,60±0,13 mm
Diameter	10,92±0,13 mm
Type of connector	Female luer
Obturator	
Total length	32,64 mm
Puncturer length	21 mm
Male luer length	8 mm
Type of connector	Male luer

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Methods and Detergents for Disinfection and Pre-cleaning

Codman Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or resterilization, after a single patient use

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Sterility

This device is intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Sterilization method used: Ethylene Oxide.

Installation and Maintenance Requirements

Codman Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or re-sterilization, after a single patient use

Possible Combinations with Other Equipment

The CODMAN MICROSENSOR ICP Transducer is designed for use with the ICP EXPRESS Monitor, catalog no. 82-6634 (117 VAC) and 82-6635 (230 VAC) (Russian Registration Certificate N RZN 2013/952 dated August 07, 2013) or any other suitable CODMAN pressure monitoring interface or device.

The CODMAN MICROSENSOR Kits are designed for use with the CODMAN Cranial Hand Drill, catalog no. 82-6607. The drill facilitates access to the intraparenchymal area. The drill is also available as a component of the CODMAN Cranial Access Kit, catalog no. 82-6612, 82-6614, and 82-6616.

If desired, attach the drainage port of the ventricular catheter (for catalog no. 626633 and 626653) to an appropriate external drainage system such as the Codman EDS 3 (catalog no. 82-1731).

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**The list of documents/standards
applicable to the medical device**

Quality Systems	
ISO 13485: 2003	Medical devices -- Quality management systems -- Requirements for regulatory purposes
Risk	
ISO 14971: 2012	Risk Management for Medical Devices
Product Safety and Performance Standards	
ASTM F2052-14	Standard Test Method for Measurement of Magnetically Induced Displacement Force on Medical Devices in the Magnetic Resonance Imaging
ASTM F2182-11a	Standard Test Method for Measurement of Radio Frequency Induced Heating On or Near Passive Implants During Magnetic Resonance Imaging
ASTM F2119-07	Standard Test Method for Evaluation of MR Image Artifacts from Passive Implants
ASTM F2213-06	Standard Test Method for Measurement of Magnetically Induced Torque on Medical Devices in the Magnetic Resonance Environment
ASTM F2503-13	Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment
IEC 60601-1	Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2	Collateral Standard: Electromagnetic disturbances - Requirements and tests
IEC 60601-1-6	Collateral standard: Usability
IEC 60601-1-9	Collateral Standard: Requirements for environmentally conscious design
EN ISO 14630: 2012	Non-active surgical implants. General requirements (components other than the MICROSENSOR probe)
Clinical	
MedDev 2.7.1: 2009	Guidelines on medical devices-Evaluation of clinical data. A guide for manufacturers and Notified Bodies
Biocompatibility	
ISO 10993- 1: 2009/AC2010	Biological Evaluation of Medical Devices: Part 1: Evaluation and testing with a risk management process
Packaging	
ISO 11607: 2009 Part 1	Packaging for terminally sterilized medical devices Requirements for materials, sterile barrier systems and packaging systems
ISO 11607:2009 Part 2	Packaging for terminally sterilized medical devices Validation requirements for forming, sealing and assembly processes
Labeling	
EN 980: 2008	Graphic Symbols for use in labeling of medical devices
EN 1041:2008	Information Supplied by Manufacturer of Medical Devices
EN ISO 15223-1:2012	– Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied
IEC TR 60878:2003	– Graphical symbols for electrical equipment in medical practice
Sterilization	
ISO 11135: 2007	Validation and Routine Control of Ethylene Oxide Sterilization

ISO 10993-7:2008	Biological Evaluation of Medical Devices – Part 7: Ethylene Oxide Sterilization Residuals
ISO 11737-2:2009	Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
ANSI ST72:2011	Bacterial endotoxins — Test methods, routine monitoring, and alternatives to batch testing

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Warranty

Codman & Shurtleff, Inc. warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.** Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the label hereof.

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Figure 6
A. To tighten
B. Blue Adapter


Figure 8. Transducer configuration to ensure patient safety during MRI

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
CODMAN NEURO

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

January 9, 2020

Letter of Declaration


As a responsible representative of **Codman & Shurtleff, Inc.**, I, **Vivian Perez**, hereby certify and attest that the attached document, Instructions For Use (IFU) Codman Microsensor Ventricular Catheter with Tuohy-Borst Adapter Kit Variant III - Ref 626633, is a true, exact, unaltered copy of the original document issued by **Codman & Shurtleff, Inc.**, located at 325 Paramount Drive, Raynham, MA 02767-0350, USA.


Vivian Perez, MBA, RAC
Associate Director Regulatory Affairs

State of Florida
County of Miami-Dade

On this 9th day of January, 2020, before me, the undersigned notary public, personally appeared **Vivian Perez**, proved to me through satisfactory evidence of identification, which were Personal Knowledge, to be the person whose name is signed on the preceding or attached document, and who swore or affirmed to me that the contents of the document are truthful and accurate to the best of her knowledge and belief.


Codman & Shurtleff, Inc.
325 Paramount Drive, Raynham, MA 02767-0350 • (508) 880-8100 • www.codman.com

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INSTRUCTIONS FOR USE
Medical Device
CODMAN MICROSENSOR Kit
*Variant III. CODMAN MICROSENSOR Ventricular Catheter with
Tuohy-Borst Adapter Kit
REF 626633*

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Instructions for Use
CODMAN MICROSENSOR Ventricular Catheter with
Tuohy-Borst Adapter Kit

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IMPORTANT INFORMATION

Please Read Before Use

CODMAN MICROSENSOR
Ventricular Catheter with Tuohy-
Borst Adapter Kit

REF 626633

STERILE EO

Rx Only

Description

The CODMAN MICROSENSOR Ventricular Catheter with Tuohy-Borst Adapter Kit, catalog no. 626633, consists of the CODMAN MICROSENSOR ICP Transducer, a 35 cm ventricular catheter with female LUER-LOK connector, a Tuohy-Borst adapter with drainage port, a split trocar, a 10-gauge ventricular needle and a 15 cm stylet. See Figure 1.

The CODMAN MICROSENSOR ICP Transducer is a catheter with a microminiature strain gauge sensor mounted at one end and an electrical connector at the other end. It is designed for use with the ICP EXPRESS Monitor, catalog no. 82-6634 (117VAC) and 82-6635 (230 VAC) or any other suitable CODMAN pressure monitoring interface or device.

Intended Use

CODMAN MICROSENSOR Ventricular Catheter with Tuohy-Borst Adapter Kit is designed for use for direct measurement of intracranial pressure.

Indications

Use of the CODMAN MICROSENSOR Ventricular Catheter with Tuohy-Borst Adapter Kit is indicated when direct intracranial pressure monitoring is required. The kit is indicated for use in intraventricular pressure

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monitoring and cerebrospinal fluid (CSF) drainage applications.

Contraindications

Ventriculostomy is contraindicated in patients with coagulopathy or active infection in the area of the catheter. Use of the CODMAN MICROSENSOR Ventricular Catheter with Tuohy-Borst Adapter Kit is contraindicated in children less than one year of age.

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Before conducting an MRI procedure on a patient with an implanted CODMAN MICROSENSOR Transducer, read the *MRI Information* section. Failure to read and strictly

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adhere to these guidelines can result in serious injury to the patient.

Precautions

Inspect the sterile package carefully. Do not use if:

- the package or seal appears damaged,
- contents appear damaged, or
- the expiry date has passed.

Avoid direct contact with the sensing element at the transducer tip. Care must be taken at all times during handling of the transducer to protect the tip from impact. Damage could result.

Do not hit the transducer tip with the stylet. Damage could result.

It is essential to maintain strict sterile technique during ventriculostomy, and subsequent transducer placement.

Exposure to electrostatic discharge (ESD) energy could damage this device. High levels of ESD could damage the electronic components and cause the transducer to be rendered inaccurate or inoperable. Take all precautions to reduce the buildup of electrostatic charge during the use of this product and avoid touching the transducer connector pins, which are identified with the ESD symbol. (Refer to *Electrostatic Discharge (ESD) Information* section).

Silicone is highly electrostatic and, as a result, attracts airborne particles and surface contaminants that could produce tissue reaction. Use extreme to prevent the ventricular catheter from coming in contact with towels, drapes, talc, or any other lint or granular surface, as well as airborne particles. Silicone has a low cut and tear resistance; therefore, care should be exercised when placing ligatures so as not to tie them too tightly. The use of stainless-steel ligatures on silicone is not advised. In addition, tying sutures too tightly can collapse the wall of the sensor body, causing damage to internal wires. The transducer must be zeroed at atmospheric pressure prior to implantation. The transducer tip must remain wet during the zeroing process.

Do not submerge the tip of the transducer or catheter vertically in a deep pool or cup of sterile water/sterile saline. Doing so will

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impose a hydrostatic pressure on the transducer diaphragm that is higher than atmospheric zero, resulting in an inaccurate zero reference. The transducer can be damaged if exposed to pressures over 1250 mmHg (166 650 Pa). Do not forcibly pull or jerk the transducer catheter. Do not expose the transducer to solvents or cleaning agents, including alcohol; these may cause damage leading to inaccurate ICP measurements. Read all instructions included with the monitoring display device prior to use.

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Adverse Events

The following Adverse Events may occur with the use of the CODMAN MICROSENSOR:

- Hemorrhage*
- Infection
- Subcutaneous CSF leakage
- Neurological sequelae

*Subarachnoid, intracerebral, or extracerebral hemorrhage may occur at the site of transducer placement (either skull, cortical, or dural areas). Testing of the blood clotting factor should be conducted on patients before insertion.

MRI Information

MR Conditional

Read and understand this document in its entirety prior to performing a Magnetic Resonance Imaging Procedure on a

patient with an implanted CODMAN MICROSENSOR. Failure to adhere to the Conditions for Safe Use may result in serious injury to the patient. The **CODMAN MICROSENSOR** is MR Conditional

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MRI SAFETY INFORMATION:

Non-clinical testing demonstrated that the **CODMAN MICROSENSOR** is MR Conditional. A patient implanted with this device can be safely scanned in an MR system which meets or is operated under the following conditions:

- Static magnetic field of 1.5 and 3 tesla only.
- Maximum spatial gradient magnetic field of 1000 G/cm (10 T/m).
- Maximum gradient field slew rate of 170 T/m/s.
- Horizontal cylindrical bore MRI scanner.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of
- 2.0 W/kg (Normal Operating Mode), for 15 minutes of continuous scanning.
- Special positioning of the CODMAN MICROSENSOR is required to ensure patient safety during the MRI procedure (see "PREPARATION FOR THE MRI PROCEDURE" below for specific instructions).
- WARNING: Do not bring the ICP EXPRESS monitor, cables or other accessories such as Tuohy needles, trocar or stylet into the MRI suite.
- WARNING: Do not use Transmit / Receive or Transmit-only RF Head coils. Only use Transmit / Receive RF Body coil or Transmit RF Body coil / Receive-only RF Head coil.
- WARNING: Do not scan a patient with an elevated body temperature.

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MRI-Related Heating

Under the scanning conditions defined above, the CODMAN MICROSENSOR is expected to produce a maximum temperature rise of less than 2°C after 15 minutes of continuous scanning. The effects of scanning beyond 15 minutes are undetermined.

Artifact Information

In non-clinical testing, the maximum artifact size was seen on the gradient echo pulse sequence at 3T and extends to a zone approximately 2mm relative to the size and shape of the CODMAN MICROSENSOR.

Preparation for the MRI procedure:

1. Immediately prior to entering the MRI environment, verify that the CODMAN MICROSENSOR is functioning properly. DO NOT perform an MRI procedure if the CODMAN MICROSENSOR is damaged or otherwise not functioning properly.
2. Disconnect all cables and patient monitoring devices attached to the CODMAN MICROSENSOR prior to transporting the patient into the MRI environment. DO NOT bring the ICP EXPRESS Monitor, cables or other accessories such as the trocar or stylet into the MRI environment.
3. Special positioning of the CODMAN MICROSENSOR is required to ensure patient safety during the MRI procedure. The CODMAN MICROSENSOR must be placed in a specific geometry to minimize the potential for excessive heating of the sensor tip. Coil the tubing of the CODMAN MICROSENSOR near the base of the electrical connector into 5 or 6 loops approximately 6 cm in diameter and center on top of the patient's head (see Figure 7). Do not perform MRI with the CODMAN MICROSENSOR in a "straight line" configuration (i.e., uncoiled). Failure to follow this guideline can result in serious injury to the patient.
4. Insert a dry gauze pad at least 1cm thick between the CODMAN MICROSENSOR

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electrical connector with coiled tubing and the patient's scalp. Secure in place using tape (see Figure 7). Use care when removing the tape to prevent damage to the CODMAN MICROSENSOR.

5. Do not exceed the following MRI parameters during imaging:

- a. Maximum spatial gradient magnetic field of 1,000 G/cm (10 T/m). The highest SG magnetic field is commonly located off-axis, at a side wall, and near the opening of the bore of the scanner. Please refer to MRI manufacturer's published value and location of the peak SG that is accessible to the patient.
- b. Maximum gradient field slew rate of 170 T/m/s.
- c. Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode), for 15 minutes of continuous scanning.

Electrostatic Discharge (ESD) Information

CAUTION: Exposure to electrostatic discharge (ESD) energy could damage this device. High levels of ESD could damage the electronic

components and cause the transducer to be rendered inaccurate or inoperable. Take all precautions to reduce the buildup of electrostatic charge during the use of this product.

- Provide patient grounding; e.g., grounding straps on gurneys
- Avoid the use of materials that could generate ESD during patient movement and transport; e.g., nylon transfer boards with bedding
- Before touching the patient, caretakers should discharge ESD buildup by touching a grounded metal surface, such as a bedrail

It is recommended that all hospital personnel in contact with these devices receive an

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explanation of the ESD symbol and training in ESD precautionary procedures. Training should include, at a minimum, an introduction to electrostatic discharge, when and why it occurs, precautionary measure, and the damage that can be done to electronic components if touched by a user who is electrostatically charged.

Avoid touching the connector pins, which are identified with the ESD symbol, before following ESD precautionary procedures. Avoid touching the transducer tip (sensing element) at all times.

How Supplied

This device is intended for SINGLE USE ONLY; DO NOT RESTERILIZE.

Codman Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or re-sterilization, after a single patient use. These devices are intended to come into contact with the central nervous system and the ability does not currently exist to destroy possible contaminants such as Creutzfeldt- Jakob Disease. Reuse can also compromise device performance and any usage beyond the design intent of this single-use device can result in unpredictable use hazards or loss of functionality.

Codman & Shurtleff will not be responsible for any product that is re-sterilized, nor accept for credit or exchange any product that has been opened but not used.

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As long as the individual package is not damaged or opened, the product is sterile. All components have been tested and were determined to be nonpyrogenic, **except** for the electrical connector of the sensor. The device is packaged using a combination of recyclable and non-recyclable materials. Recycle or dispose of all packaging waste in accordance with hospital procedures and regulations.

Environment protection requirements

Medical device during use, transportation and storage has no negative effects on humans and the environment.

Medical Device Recycling or Disposal

The device is packaged using a combination of recyclable and non-recyclable materials. Recycle or dispose of all packaging waste in accordance with hospital procedures and regulations.

Unused product should be recycled in accordance with the rules for the treatment of medical waste of class A.

Used product should be recycled in accordance with the rules for the treatment of medical waste of class B

Connecting and Zeroing the Transducer

CAUTION: The transducer must be zeroed at atmospheric pressure before implantation.

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1. Connect the transducer to the CODMAN ICP Monitor using an appropriate CODMAN interface cable. Refer to instructions for use provided with the cable for sterilization information.

2. If applicable, connect the CODMAN ICP Monitor to an available pressure channel on an external patient monitor using a CODMAN Patient Monitor Interface Cable.

CAUTION: Use CODMAN Patient Monitor Interface Cables only with the patient monitors for which they are specifically designed and designated.

3. Zero and calibrate the external patient monitor according to the instructions provided with the CODMAN ICP Monitor, as well as the external patient monitor manufacturer's instructions.

4. Prepare to zero the catheter by laying the tip of the catheter flat in a shallow pool of sterile water or sterile saline. The accompanying sterile blister package has a marked well that is suitable for this procedure. Pour sufficient sterile water or sterile saline into the well, then lay at least a 5-cm section of the catheter horizontally just under the surface of the sterile water or sterile saline.

CAUTION: Do not submerge the tip of the catheter vertically in a deep pool or cup of sterile water or sterile saline. Doing so will impose a hydrostatic pressure on the sensor diaphragm that is higher than atmospheric pressure, resulting in an inaccurate zero reference.

5. While keeping the tip of the catheter flat and still in the sterile water or sterile saline, zero the catheter according to the instructions provided with the ICP EXPRESS Monitor.

CAUTION: The sensor tip must remain wet during the zeroing process.

CAUTION: The sensor tip must remain still during the zeroing process. Motion of the sensor may be interpreted by the ICP Monitor as a fluctuating ICP signal which will prevent the sensor zeroing process from completing successfully.

6. If available, record the three-digit zero reference number provided by the ICP Monitor. Mark this number on the catheter connector housing and the patient's chart for future reference.

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General Surgical Procedure

The following is a general guide for informational purposes only. The surgeon may wish to alter details in accordance with his or her own clinical experience and medical judgment. The CODMAN Cranial Access Kit is recommended for this procedure.

и Connect and zero the transducer. Refer to *Connecting and Zeroing the MICROSENSOR Transducer*.

2. Incise the scalp to expose the skull and then make a twist drill hole using a 5.8 mm drill bit.
3. Gently bevel the burr hole on one side to allow the catheter to exit without a sharp angulation.
4. Puncture the dura using a needle or cautery.
5. Use the split trocar to tunnel the ventricular catheter under the scalp from the craniotomy site to any other convenient site located at least 3 cm from the burr hole.
6. Remove the stylet from the split trocar.
7. Insert the ventricular catheter through the split trocar from the exit site to the burr hole wound, as shown in Figure 2.
8. Remove the split trocar from under the scalp.
9. Verify that the transducer tip is located at the first large hole at the distal end of the ventricular catheter, as shown in Figure 8. If necessary, loosen the Tuohy-Borst adapter cap (turn counterclockwise) and adjust transducer position as required.

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10. Firmly tighten the Tuohy-Borst adapter cap (turn clockwise) and LUER-LOK connections. Verify that the transducer tip remains in position as shown in Figure 8. Readjust if necessary.

11. Insert the tip of the 15 cm stylet through the middle large hole at the distal end of the ventricular catheter. See Figure 3.

CAUTION: Do not hit the transducer tip with the stylet. Damage could result.

Accidental contact with the transducer can be detected by a simultaneous transient waveform response on the monitor trace and/or an increase in the numeric display. In such cases, do not use the transducer.

12. Advance the stylet and catheter at a right angle to the skull into the lateral ventricle at a depth of approximately 7 cm. See Figure 4.

Note: Some surgeons may prefer to first use the 10-gauge ventricular needle to locate the ventricle, as shown in Figure 5, and then pass the stylet and catheter through the tract thus created.

13. To verify ventricular placement, remove the cap from the side drain port of the Tuohy-Borst adapter to allow CSF to flow through the catheter.

14. Hold the catheter in place and gently withdraw the stylet.

15. Press down on the ventricular catheter and pull any excess catheter away from the incision site.

16. Secure the catheter to the scalp with sutures at the exit site. Additional sutures may be placed through the holes in the female LUER-LOK connector flange.

17. Close and dress the incision site.

18. If desired, attach the drainage port of the ventricular catheter to a ventricular drainage system.

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Alternate Tunneling Technique

The following is an alternate technique for tunneling the ventricular catheter under the scalp without the transducer preloaded. The surgeon may wish to alter details in accordance with his or her clinical experience and medical judgment.

1. Carefully detach the ventricular catheter from the white female LUER-LOK connector and Tuohy-Borst adapter.
2. Gently pull the Tuohy-Borst adapter and transducer assembly completely out of the catheter.
3. Connect and zero the transducer as described in *Connecting and Zeroing the MICROSENSOR Transducer*.
4. Use either a split trocar (included) or a standard tunneling trocar (not included), see Figure 6 to tunnel the ventricular catheter under the scalp from the craniotomy site to any other convenient site located at least 3 cm from the burr hole.
5. At this point, the catheter can be cut to the desired length. Leave enough excess catheter at the incision site to permit cannulation.
6. Thread the tip of the transducer back through the proximal end of the catheter as shown in Figure 7. Sterile saline may be used to irrigate the inner channel of the catheter to facilitate transducer advancement. Continue to advance the transducer until the tip can be seen through the first hole at the distal end of the catheter, as shown in Figure 8.
7. Reattach the white female LUER-LOK connector and Tuohy-Borst adapter to the

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ventricular catheter as shown in Figure 9.

Finger-tighten all LUER-LOK connections.

8. Firmly tighten the Tuohy-Borst adapter cap (turn clockwise), verifying that the transducer tip remains in position as shown in Figure 8. Readjust if necessary.

9. Introduce the ventricular catheter into the ventricle as described in *General Surgical Procedure*.

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Package contents

One consumer package included:

1. CODMAN MICROSENSOR ventricular catheter Kit with Tuohy-Borst adapter for intracranial pressure measurement consisting of:
 - ICP Transducer - 1 pcs
 - ventricular catheter - 1 pcs
 - Tuohy-Borst adapter with drainage port - 1 pcs
 - ventricular needle - 1 pcs
 - Stylet - 1 pcs
 - Split trocar - 1 pcs
2. Instructions for use - 1 pcs
3. Patient labels - 6 pcs

Specifications

The tolerances of the characteristics are $\pm 5\%$, unless otherwise specified.

- ICP Transducer

Technical characteristics

Sensing element	Strain gauge silicon microchip
Probe usable length	100 cm nominal

Probe tip diameter	1,3 mm maximum
Probe tubing diameter	0,8 mm maximum
Probe tubing length	105 ± 2,5 cm
Sensor housing dimensions	51 x 17 x 9 mm ± 0,1 mm
Functional pressure range	-50 mmHg to +250 mmHg (- 6666 to +33330 Pa)
Functional overpressure range without damage	-700 mmHg to +1250 mmHg (- 93324 to +166650 Pa)
Input/output impedance	1000 ohms nominal
Zero drift	No greater than 5 mmHg over 30 days
Output signal (sensitivity)	5 µV/V/mmHg nominal
Frequency response	Greater than 200 Hz
Markings of probe tubing	Thin markings (16 places) dimensions 0,6 mm ± 0.1 mm thick
	Thick markings centered at 5 cm intervals from the tip (10 markings) dimensions 1,1 mm ± 0.1 mm thick. One mark @5 cm, two marks @10 cm, three marks @15 cm and four marks @20 cm. Spacing between thick markings to be 0,5 mm ± 0.1 mm.
Length from the tip to the first mark	9.7 ± 1 mm.
Pull force	Joint/bonds – 2,2 N Epoxy cap/ball – 1 N

Note: All performance specifications based on 5 VDC excitation voltage

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- Tuohy-Borst adapter with drainage port

Technical characteristics

Length of adapter cap (head)	11,94±0,13 mm
Length of adapter body	23,11±0,13 mm
Length of the drainage port tube	38,10±8,89 mm
Outer diameter of drainage port tube	3,18±0,13 mm
Internal diameter of the drainage port tube	1,91±0,13 mm
Length of drainage port connector	19,05±0,38 mm

Type of connectors	Male luer (adapter) Female luer (drainage port)
--------------------	--

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Warranty

Codman & Shurtleff, Inc. warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.** Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.

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

January 9, 2020

Letter of Declaration


As a responsible representative of **Codman & Shurtleff, Inc.**, I, **Vivian Perez**, hereby certify and attest that the attached document, Instructions For Use (IFU) Codman Microsensor Metal Skull Bolt - Kit Variant IV - Ref 626638, is a true, exact, unaltered copy of the original document issued by **Codman & Shurtleff, Inc.**, located at 325 Paramount Drive, Raynham, MA 02767-0350, USA.

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
INSTRUCTIONS FOR USE Medical Device CODMAN MICROSENSOR Kit *Variant IV. CODMAN MICROSENSOR Metal Skull Bolt Kit REF 626638*

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**Instructions for Use
CODMAN MICROSENSOR Metal Skull Bolt Kit**

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STERILE EO

Description

The CODMAN MICROSENSOR Metal Skull Bolt Kit, catalog no. 626638, consists of the CODMAN MICROSENSOR ICP Transducer, Skull Bolt with integrated compression cap and optional washer, Obturator/Dura Pierce, and Drill Bit with depth guide and hex wrench. See Figure 1.

The CODMAN MICROSENSOR ICP Transducer is a catheter with a microminiature strain gauge pressure sensor mounted at one end and an electrical connector at the other end. It is designed for use with the ICP EXPRESS Monitor, catalog no. 82-6634 (117 VAC) and 82-6635 (230 VAC) or any other suitable CODMAN pressure monitoring interface or device.

The CODMAN MICROSENSOR Metal Skull Bolt Kit is designed for use with the CODMAN Cranial Hand Drill, catalog no. 82-6607. The drill facilitates access to the intraparenchymal area. The drill is also available as a component of the CODMAN Cranial Access Kit, catalog no. 82-6612, 82-6614, and 82-6616.

The skull bolt is preassembled with a washer for optional use in adjusting the depth of the bolt and a compression cap for securing the transducer in place. An obturator/dura pierce is also provided for clearing the passage for the transducer and piercing the dura. See Figure 2.

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Intended Use

CODMAN MICROSENSOR Metal Skull Bolt Kit is designed for use for direct measurement of intracranial pressure.

Indications

Use of the CODMAN MICROSENSOR Metal Skull Bolt Kit is indicated when direct intracranial pressure (ICP) monitoring is required. The kit is indicated for use in both subdural and intraparenchymal pressure monitoring applications.

Contraindications

Use of the skull bolt is contraindicated in children less than one year of age.

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This kit is not designed, sold, or intended for use as a therapeutic device.

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Warnings

Installation of the skull bolt must be performed with the bolt held within 10° degrees of the perpendicular to the incision site. Installing at an angle can result in fracture of the device or in an inadequate seal between the washer and the skull.

Excessive torque applied to the skull bolt during insertion can cause breakage.

The use of electrosurgical equipment; e.g., monopolar, bipolar, diathermy, can cause damage to the CODMAN MICROSENSOR Transducer and/or the ICP EXPRESS

MONITOR. This could lead to permanent or temporary disabling of either device.

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Exposure to electrostatic discharge (ESD) energy could damage this device. High levels of ESD could damage the electronic components and cause the transducer to be rendered inaccurate or inoperable. Take all precautions to reduce the buildup of electrostatic charge during the use of this product and avoid touching the transducer connector pins, which are identified

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Preparation for the MRI procedure:

- 1.Immediately prior to entering the MRI environment, verify that the CODMAN MICROSENSOR is functioning properly. DO NOT perform an MRI procedure if the CODMAN MICROSENSOR is damaged or otherwise not functioning properly.
- 2.Disconnect all cables and patient monitoring devices attached to the CODMAN MICROSENSOR prior to transporting the patient into the MRI environment. DO NOT bring the ICP EXPRESS Monitor, cables or other accessories such as the trocar or stylet into the MRI environment.
- 3.Special positioning of the CODMAN MICROSENSOR is required to ensure patient safety during the MRI procedure. The CODMAN MICROSENSOR must be placed in a specific geometry to minimize the potential for excessive heating of the sensor tip. Coil the tubing of the CODMAN MICROSENSOR near the base of the electrical connector into 5 or 6 loops approximately 6 cm in diameter and center on top of the patient's head (see Figure 9). Do not perform MRI with the CODMAN MICROSENSOR in a "straight line" configuration (i.e., uncoiled). Failure to follow this guideline can result in serious injury to the patient.
4. Insert a dry gauze pad at least 1 cm thick between the CODMAN MICROSENSOR

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electrical connector with coiled tubing and the patient's scalp. Secure in place using tape (see Figure 9). Use care when removing the tape to prevent damage to the CODMAN MICROSENSOR.

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It is recommended that all hospital personnel in contact with these devices receive an explanation of the ESD symbol and training in ESD precautionary procedures. Training should include, at a minimum, an introduction to electrostatic discharge, when and why it occurs, precautionary measure, and the damage that can be done to electronic components if touched by a user who is electrostatically charged.

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How Supplied

2 This device is intended for SINGLE USE ONLY; DO NOT RESTERILIZE.

Codman Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or re-sterilization, after a single patient use. These devices are intended to come into contact with the central nervous system and the ability does not currently exist to destroy possible contaminants such as Creutzfeldt- Jakob Disease. Reuse can also compromise device performance and any usage beyond the design intent of this single-use device can result in unpredictable use hazards or loss off functionality.

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As long as the individual package is not damaged or opened, the product is sterile. All components have been tested and were determined to be nonpyrogenic, **except** for the electrical connector of the transducer, the hex wrench, the dura pierce, the bolt assembly, and the silicone tubing used for packaging, which are not tested.

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4. Prepare to zero the catheter by laying the tip of the catheter flat in a shallow pool of sterile water or sterile saline. The accompanying sterile blister package has a marked well that is suitable for this procedure. Pour sufficient sterile water or sterile saline into the well, then lay at least a 5 cm section of the catheter horizontally just under the surface of the sterile water or sterile saline.

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6. If available, record the three-digit zero reference number provided by the ICP Monitor. Mark this

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3. If desired, adjust the seating depth of the bolt using the washer provided.

4. Condition the compression cap and sensor locking mechanism before installation of the bolt. First, hand-tighten the compression cap by turning it clockwise. Then loosen the cap by turning it counterclockwise the same or a lesser amount.

5. Put the skull bolt in position and screw manually in a clockwise direction until properly seated. See Figure 3.

CAUTION: Do not loosen the cap during installation of the bolt.

WARNING: Excessive torque applied to the skull bolt during insertion can cause breakage. Installation of the skull bolt must be performed with the bolt held within 10° of perpendicular to the incision site. Installing at an angle may result in fracture of the device.

6. Use the obturator/dura pierce to clear the pathway through the bolt and to pierce the dura. See Figure 4.

CAUTION: Do not plunge the obturator. Carefully push the obturator until the dura has been pierced.

7. Irrigate the channel with non-bacteriostatic sterile saline. See Figure 5.

8. Guide the transducer through the bolt until the transducer tip is in the subdural or parenchymal space, as shown in Figure 6.

9. Turn the compression cap firmly in a clockwise direction to secure the transducer. See Figure 7.

WARNING: If the skull bolt compression cap is not initially tightened properly, transducer movement and/or leakage can occur. If during the implantation period the patient is moved or becomes active, check the security of the transducer in the skull bolt compression cap and tighten, if necessary.

10. Remove retractor and close the skin incision in one layer using interrupted sutures, as shown in Figure 8. Apply appropriate dressing to wound site.

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Removing the Skull Bolt

1. Loosen the compression cap by turning in a counterclockwise direction.
2. Carefully withdraw the transducer from the bolt assembly.
3. Manually unscrew the bolt in a counterclockwise direction until the bolt is free.

CAUTION: Make certain that the depth washer, if used, has been removed from the incision area before closing.

4. Apply pressure with sterile gauze until hemostasis is achieved. Close the incision again in one layer, with uninterrupted sutures. Apply a dry sterile dressing.

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Package contents

One consumer package included:

1. CODMAN MICROSENSOR metal Skull Bolt Kit for intracranial pressure measurement consisting of:
 - ICP Transducer – 1 pcs
 - Skull Bolt – 1 pcs
 - compression cap – 1 pcs
 - Spacing Washer – 1 pcs
 - Obturator/Dura Pierce – 1 pcs
 - Drill Bit with depth guide – 1 pcs
 - hex wrench – 1 pcs
2. Instructions for use – 1 pcs
3. Patient labels – 6 pcs

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	One mark @5 cm, two marks @10 cm, three marks @15 cm and four marks @20 cm. Spacing between thick markings to be 0,5 mm $\pm$ 0.1 mm.
Length from the tip to the first mark	9.7 $\pm$ 1 mm.
Pull force	Joint/bonds – 2,2 N Epoxy cap/ball – 1 N

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Quality Systems	
ISO 13485: 2003	Medical devices -- Quality management systems -- Requirements for regulatory purposes

Risk	
ISO 14971: 2012	Risk Management for Medical Devices
Product Safety and Performance Standards	
ASTM F2052-14	Standard Test Method for Measurement of Magnetically Induced Displacement Force on Medical Devices in the Magnetic Resonance Imaging
ASTM F2182-11a	Standard Test Method for Measurement of Radio Frequency Induced Heating On or Near Passive Implants During Magnetic Resonance Imaging
ASTM F2119-07	Standard Test Method for Evaluation of MR Image Artifacts from Passive Implants
ASTM F2213-06	Standard Test Method for Measurement of Magnetically Induced Torque on Medical Devices in the Magnetic Resonance Environment
ASTM F2503-13	Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment
IEC 60601-1	Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2	Collateral Standard: Electromagnetic disturbances - Requirements and tests
IEC 60601-1-6	Collateral standard: Usability
IEC 60601-1-9	Collateral Standard: Requirements for environmentally conscious design
EN ISO 14630: 2012	Non-active surgical implants. General requirements (components other than the MICROSENSOR probe)
Clinical	
MedDev 2.7.1: 2009	Guidelines on medical devices-Evaluation of clinical data. A guide for manufacturers and Notified Bodies
Biocompatibility	
ISO 10993- 1: 2009/AC2010	Biological Evaluation of Medical Devices: Part 1: Evaluation and testing with a risk management process
Packaging	
ISO 11607: 2009 Part 1	Packaging for terminally sterilized medical devices Requirements for materials, sterile barrier systems and packaging systems
ISO 11607:2009 Part 2	Packaging for terminally sterilized medical devices Validation requirements for forming, sealing and assembly processes
Labeling	
EN 980: 2008	Graphic Symbols for use in labeling of medical devices
EN 1041:2008	Information Supplied by Manufacturer of Medical Devices
EN ISO 15223-1:2012	– Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied
IEC TR 60878:2003	– Graphical symbols for electrical equipment in medical practice
Sterilization	
ISO 11135: 2007	Validation and Routine Control of Ethylene Oxide Sterilization
ISO 10993-7:2008	Biological Evaluation of Medical Devices – Part 7: Ethylene Oxide Sterilization Residuals
ISO 11737-2:2009	Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
ANSI ST72:2011	Bacterial endotoxins — Test methods, routine monitoring, and alternatives to batch testing

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
January 9, 2020


Letter of Declaration

As a responsible representative of **Codman & Shurtleff, Inc.**, I, **Vivian Perez**, hereby certify and attest that the attached document, Codman Microsensor Ventricular Catheter Kit – 121119- Variant V - REF 626653, is a true, exact, unaltered copy of the original document issued by **Codman & Shurtleff, Inc.**, located at 325 Paramount Drive, Raynham, MA 02767-0350, USA.

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
INSTRUCTIONS FOR USE Medical Device

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**Instructions for Use
CODMAN MICROSENSOR Ventricular Catheter Kit**

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STERILE EO

Description

The CODMAN MICROSENSOR Ventricular Catheter Kit, catalog no 626653, consists of a 38 cm ventricular catheter with an integrated stylet and intracranial pressure (ICP) sensor, a catheter anchoring clip, and a 7-gauge tunneling trocar. See Figure 1.

The ventricular catheter is a ventricular drainage catheter with a microminiature strain gauge pressure sensor integrated into the catheter wall. It is designed for use with the CODMAN ICP EXPRESS Monitor, catalog no. 82-6634 (117 VAC) and 82-6635 (230 VAC) or any other suitable pressure CODMAN monitoring interface or device.

Intended Use

CODMAN MICROSENSOR Ventricular Catheter Kit is designed for use for direct measurement of intracranial pressure.

Indications

Use of the CODMAN MICROSENSOR Ventricular Catheter Kit is indicated when direct intraventricular pressure monitoring is required. The kit is indicated for use in ICP monitoring and cerebrospinal fluid (CSF) drainage applications.

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Ventriculostomy is contraindicated in patients with coagulopathy or active infection in the area of the catheter. Use of the Ventricular Catheter Kit is contraindicated in children less than one year of age.

Note: This kit is not designed, sold, or intended for any use except as indicated.

Conditions of Use

The medical device is intended for use at the standard conditions of hospital operating room environment.

All procedures should be performed by a qualified neurosurgeon using standard surgical procedure. Neurosurgeon is responsible for the implementation of appropriate steps and procedures to avoid infection and complications.

Warnings

Take extreme care to avoid damage to the dura and underlying cerebrum.

The use of electrosurgical equipment; e.g., monopolar, bipolar, diathermy, can cause damage to the CODMAN MICROSENSOR Transducer and/or the ICP EXPRESS MONITOR. This could lead to permanent or temporary disabling of either device.

Before conducting an MRI procedure on a patient with an implanted CODMAN MICROSENSOR Transducer, read the *MRI Information* section. Failure to read and strictly adhere to these guidelines can result in serious injury to the patient.

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Avoid direct contact with the sensing element at the transducer tip. Care must be taken at all

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times during handling of the transducer to protect the tip from impact. Damage could result.

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Connections could loosen during shipping and handling.

Exposure to electrostatic discharge (ESD) energy could damage this device. High levels of ESD could damage the electronic components and cause the transducer to be rendered inaccurate or inoperable. Take all precautions to reduce the buildup of electrostatic charge during the use of this product and avoid touching the transducer connector pins, which are identified with the ESD symbol. (Refer to Electrostatic Discharge (ESD) Information section).

Silicone is highly electrostatic and, as a result, attracts airborne particles and surface contaminants that could produce tissue reaction. Use care to prevent the ventricular catheter from coming in contact with towels, drapes, talc, or any other linty or granular surface, as well as airborne particles.

Silicone has a low cut and tear resistance; therefore, care should be exercised when placing ligatures so as not to tie them too tightly. The use of stainless-steel ligatures on silicone is not advised. In addition, tying sutures too tightly can collapse the wall of the sensor body, causing damage to internal wires.

The transducer must be zeroed at atmospheric pressure prior to implantation.

The transducer tip must remain wet during the zeroing process.

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Preparation for the MRI procedure:

1. Immediately prior to entering the MRI environment, verify that the CODMAN MICROSENSOR is functioning properly. DO NOT perform an MRI procedure if the CODMAN MICROSENSOR is damaged or otherwise not functioning properly.
2. Disconnect all cables and patient monitoring devices attached to the CODMAN MICROSENSOR prior to transporting the patient into the MRI environment. DO NOT bring the ICP EXPRESS Monitor, cables or other accessories such as the trocar or stylet into the MRI environment.
3. Special positioning of the CODMAN MICROSENSOR is required to ensure patient safety during the MRI procedure. The CODMAN MICROSENSOR must be placed in a specific geometry to minimize the potential for excessive heating of the sensor tip. Coil the tubing of the CODMAN MICROSENSOR near the base of the electrical connector into 5 or 6 loops approximately 6 cm in diameter and center on top of the patient's head (see Figure 7). Do not perform MRI with the CODMAN MICROSENSOR in a "straight line" configuration (i.e., uncoiled). Failure to follow this guideline can result in serious injury to the patient.
4. Insert a dry gauze pad at least 1 cm thick between the CODMAN MICROSENSOR

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It is recommended that all hospital personnel in contact with these devices receive an explanation

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of the ESD symbol and training in ESD precautionary procedures. Training should include, at a minimum, an introduction to electrostatic discharge, when and why it occurs, precautionary measure, and the damage that can be done to electronic components if touched by a user who is electrostatically charged. Avoid touching the connector pins, which are identified with the ESD symbol, before following ESD precautionary procedures. Avoid touching the transducer tip (sensing element) at all times.

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As long as the individual package is not damaged or opened, the product is sterile.

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All components have been tested and were determined to be nonpyrogenic, except for the electrical connector of the sensor, and the trocar guard and silicone tubing used for packaging, which are not tested.

The device is packaged using a combination of recyclable and non-recyclable materials.

Recycle or dispose of all packaging waste in accordance with hospital procedures and regulations.

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The device is packaged using a combination of recyclable and non-recyclable materials.

Recycle or dispose of all packaging waste in accordance with hospital procedures and regulations.

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1. Connect the transducer to the CODMAN ICP Monitor using an appropriate CODMAN interface

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cable. Refer to instructions for use provided with the cable for sterilization information.

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General Surgical Procedure

The following is a general guide for informational purposes only. The surgeon may wish to alter details in accordance with his or her own clinical experience and medical judgment. The CODMAN Cranial Access Kit is recommended for this procedure. This device is not designed, sold, or intended for use as a therapeutic device.

1. Connect and zero the catheter. Refer to Connecting and Zeroing the MICROSENSOR Transducer.
2. Incise the scalp to expose the skull, and then make a twist drill hole using a 5.8 mm drill bit.
3. Gently bevel the burr hole on one side to allow the catheter to exit without a sharp angulation.
4. Puncture the dura using a needle or cautery.
5. Inset the tip of the catheter into the end of the clear tubing of the tunneling trocar until the tubing is filled. See Figure 2.
6. Use the tunneling trocar to tunnel the ventricular catheter under the scalp to the craniotomy site. The point of entry can be any convenient location at least 3 cm from the burr hole. See Figure 3. Due to the rigidity of the integrated stylet, it is advisable to limit tunneling to 7 cm in order to have sufficient catheter length for insertion into the ventricle.
7. Continue to advance the trocar until the trocar and the clear tubing exit the scalp at the burr hole. Remove the tunneling trocar and tubing from the catheter.
8. Pull the tunneled catheter toward the craniotomy site until the Y-connector almost touches the initial tunneling entrance. This will provide adequate catheter length for insertion.

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9. Holding the catheter at a right angle to the skull, advance the catheter into the lateral ventricle to the desired depth. See Figure 4.

10. To verify ventricular placement, observe CSF flowing through the catheter.

11. Holding the catheter as shown in Figure 5, bend the catheter to the side opposite the two black arrows to expose approximately 3 cm of the catheter stylet (See Figure 6). The black arrows indicate the portion of the catheter that has been slit for stylet removal.

12. Holding the catheter firmly in place, gently withdraw the stylet. See Figure 7.

13. Remove any excess catheter from the burr hole by gently pulling back on the catheter when it exits the scalp incision.

14. Place the catheter clip over the catheter at the exit site; suture the clip to the scalp. Additional sutures may be placed through the holes in the catheter Y-connector.

15. Close and dress the incision.

16. If desired, attach the drainage port of the ventricular catheter to an appropriate external drainage system such as the Codman EDS 3 (catalog no.82-1731)

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Package contents

One consumer package included:

1. CODMAN MICROSENSOR Ventricular Catheter Kit for intracranial pressure measurement consisting of:
 - Ventricular Catheter with an integrated Stylet and ICP Sensor – 1 pcs
 - catheter anchoring clip – 1 pcs
 - tunneling trocar – 1 pcs
2. Instructions for use – 1 pcs
3. Patient labels – 6 pcs

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- Ventricular Catheter with an integrated Stylet and ICP Sensor

Technical characteristics

ICP Sensor	
Sensing element	Strain gauge silicon microchip
Probe usable length	100 cm nominal
Probe tip diameter	1,3 mm maximum
Probe tubing diameter	0,8 mm maximum
Probe tubing length	105 $\pm$ 2,5 cm
Sensor housing dimensions	51 x 17 x 9 mm $\pm$ 0,1 mm
Functional pressure range	-50 mmHg to +250 mmHg (- 6666 to +33330 Pa)
Functional overpressure range without damage	-700 mmHg to +1250 mmHg (- 93324 to +166650 Pa)
Input/output impedance	1000 ohms nominal
Zero drift	No greater than 5 mmHg over 30 days
Output signal (sensitivity)	5 μ V/V/mmHg nominal
Frequency response	Greater than 200 Hz
Note: All performance specifications based on 5 VDC excitation voltage	
Markings of probe tubing	Thin markings (16 places) dimensions 0,6 mm $\pm$ 0.1 mm thick

	Thick markings centered at 5 cm intervals from the tip (10 markings) dimensions 1,1 mm $\pm$ 0.1 mm thick. One mark @5 cm, two marks @10 cm, three marks @15 cm and four marks @20 cm. Spacing between thick markings to be 0,5 mm $\pm$ 0.1 mm.
Length from the tip to the first mark	9.7 $\pm$ 1 mm.
Pull force	Joint/bonds – 2,2 N Epoxy cap/ball – 1 N
Ventricular Catheter	
Catheter outside diameter	3.5 mm maximum
Catheter usable length	38 cm nominal
Marking	every 10 mm, addition of the marker is after every 50 mm (1-2-3-4 token ring)
Stylet	
Length	200 mm
Hardness	52-54 RC
Roughnes	0,4 μ m

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alteration, such as disassembly, cleaning or re-sterilization, after a single patient use

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Sterilization method used: Ethylene Oxide.

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The CODMAN MICROSENSOR ICP Transducer is designed for use with the ICP EXPRESS Monitor, catalog no. 82-6634 (117 VAC) and 82-6635 (230 VAC) (Russian Registration Certificate NRZN 2013/952 dated August 07, 2013) or any other suitable CODMAN pressure monitoring interface or device.

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The list of documents/standards applicable to the medical device

Quality Systems	
ISO 13485: 2003	Medical devices -- Quality management systems -- Requirements for regulatory purposes
Risk	
ISO 14971: 2012	Risk Management for Medical Devices
Product Safety and Performance Standards	
ASTM F2052-14	Standard Test Method for Measurement of Magnetically Induced Displacement Force on Medical Devices in the Magnetic Resonance Imaging
ASTM F2182-11a	Standard Test Method for Measurement of Radio Frequency Induced Heating On or Near Passive Implants During Magnetic Resonance Imaging
ASTM F2119-07	Standard Test Method for Evaluation of MR Image Artifacts from Passive Implants
ASTM F2213-06	Standard Test Method for Measurement of Magnetically Induced Torque on Medical Devices in the Magnetic Resonance Environment
ASTM F2503-13	Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment
IEC 60601-1	Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2	Collateral Standard: Electromagnetic disturbances - Requirements and tests
IEC 60601-1-6	Collateral standard: Usability
IEC 60601-1-9	Collateral Standard: Requirements for environmentally conscious design
EN ISO 14630: 2012	Non-active surgical implants. General requirements (components other than the MICROSENSOR probe)
Clinical	
MedDev 2.7.1: 2009	Guidelines on medical devices-Evaluation of clinical data. A guide for manufacturers and Notified Bodies
Biocompatibility	
ISO 10993- 1: 2009/AC2010	Biological Evaluation of Medical Devices: Part 1: Evaluation and testing with a risk management process
Packaging	
ISO 11607: 2009 Part 1	Packaging for terminally sterilized medical devices Requirements for materials, sterile barrier systems and packaging systems
ISO 11607:2009 Part 2	Packaging for terminally sterilized medical devices Validation requirements for forming, sealing and assembly processes
Labeling	
EN 980: 2008	Graphic Symbols for use in labeling of medical devices
EN 1041:2008	Information Supplied by Manufacturer of Medical Devices
EN ISO 15223-1:2012	– Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied
IEC TR 60878:2003	– Graphical symbols for electrical equipment in medical practice
Sterilization	
ISO 11135: 2007	Validation and Routine Control of Ethylene Oxide Sterilization
ISO 10993-7:2008	Biological Evaluation of Medical Devices – Part 7: Ethylene Oxide Sterilization Residuals
ISO 11737-2:2009	Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
ANSI ST72:2011	Bacterial endotoxins — Test methods, routine monitoring, and alternatives to batch testing
Система управления качеством	
ISO 13485: 2003	Системы менеджмента качества – медицинские изделия – системные требования для целей регулирования

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Warranty

Codman & Shurtleff, Inc. warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

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