2009/APR/21 at 4:34 p.m. Doc number: M934622A001 Rev. A [NO17]

**Medtronic****DLP®**

Arteriotomy Cannula
 Canule d'artériotomie
 Arteriotomiekánüle
 Cânula de arteriotomia
 Arteriotomiecanule
 Cannula per arteriotomia
 Arteriotomikanyyle
 Arteriotomiakanyyl
 Arteriotomikanyl
 Κάнуλα αρτηριτομής
 Arteriotomikanyle
 Kaniula do kardioplegii
 Cânula de arteriotomia
 Arteriyotomi Kanülü
 Kanyla pro arteriotomii
 Arteriotomiás kanül
 Артериотомическая канюля

Josh
 Johnson
 on

Digitally
 signed by
 Josh Johnson

Date:
 2021.08.03
 12:07:04
 -05'00'

Varina Phil 24 Aug 2021

VARINA V. PHILAVONG
 NOTARY PUBLIC - MINNESOTA
 MY COMMISSION EXPIRES 01/31/26

As a registered Notary Public for the State of Minnesota,
 I hereby certify that this is a true and accurate copy of
 the original document. The original is filed at
 Medtronic, Inc.

Instructions for Use • Mode d'emploi • Gebrauchsanweisung • Instrucciones de uso
 Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning • Käyttöohjeet
 Bruksanvisning • Οδηγίες χρήσης • Brugsanvisning • Instrukcja użytkowania
 Instruções de utilização • Kullanım Talimatları • Návod k použití • Használati útmutató
 Инструкция по эксплуатации

[USA] Caution: Federal law (USA) restricts this device to sale by or on the
 order of a physician.

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Arteriotomy Cannula

Product Description

The IMA (Internal Mammary Artery) Cannula consists of a plastic female slip fit connector and 26-gauge stainless steel tubing that ends in a 1-mm bulb. The Arteriotomy Cannulae consist of a plastic winged female luer lock connector with a length of polyurethane tubing that ends in a radiopaque bulb tip.

Sterile, nonpyrogenic, single use.

Indications for Use

These cannulae are intended for the perfusion of cardioplegia solution directly into a cardiac bypass graft or artery, or retrograde flushing of the free internal mammary artery with a vasodilator. These devices are intended for use up to six hours or less.

Contraindications

This device is not intended for use except as indicated above.

Warnings

These cannulae may not be detectable on x-ray in some cases due to the variances in anatomical densities and x-ray techniques.

Precautions

Note: Proper surgical procedures and techniques are the responsibility of the medical professional. The described procedure is furnished for information purposes only. Surgeons must evaluate the appropriateness of the procedure based on their own medical training and experience, and the type of surgical procedure.

These cannulae are designed and intended for single use only. DO NOT REUSE.

Adverse Effects

This device, as do all extracorporeal blood system devices, has possible side effects which include, but are not limited to, infections, blood loss, thrombus formation, embolic events, and dislocation. Vessel damage and complications at the puncture site may occur if the Instructions for Use are not followed.

Instructions for Use

1. Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the cannula onto the sterile field utilizing an aseptic technique.
2. Connect the cannula to the cardioplegia solution delivery tubing or a pressure sensing syringe.
3. Prime the cannula to remove any air before perfusing the bypass graft or artery.
4. Carefully insert the cannula into the target vessel using caution to avoid any trauma to the internal intima.
5. Slowly increase the flow to the desired rate or the maximum line pressure prescribed by the surgeon. When using a pressure sensing syringe, caution must be exercised so as not to over pressurize or over inflate the artery or graft.
6. When the surgical procedure is completed, dispose of the cannula according to hospital policy. DO NOT REUSE.

The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE ARTERIOTOMY CANNULA, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining

portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE ARTERIOTOMY CANNULA, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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