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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

00552193

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of BLUE SAIL MEDICAL CO., LTD. on the annexed DOCUMENT is genuine.

China Council
of International Trade

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Authorized
Signature: Wei Shanshan

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证书查询网址 Website for verifying the certificate: <http://www.raccpit.com/validata.html>

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bluesail+

No.:IFU_MDD_Nitrile Gloves_1.1
released: 11.10.2021

INSTRUCTION FOR USE

Diagnostic examination disposable gloves nitrile powder-free
textured non-sterile

Approved
Blue Sail Medical Co., Ltd.
QA Department Manager: Robin Liu

/Signature/
Stamp * Blue Sail Medical Co., Ltd. *

Blue Sail Medical Co., Ltd.
No.21, Qingtian Road, Qilu Chemical Industrial Park,
255414 Zibo, Shandong Province, PEOPLE'S REPUBLIC OF CHINA;
Tel.: +86 0533 7871098; e-mail: info@bluesail.cn

revision 1.1
released: 11.10.2021

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1. NAME OF MEDICAL DEVICE:

«Diagnostic examination disposable gloves nitrile powder-free textured non-sterile»

"Diagnostic examination disposable gloves nitrile powder-free textured non-sterile" in the composition:

1. Diagnostic examination disposable gloves nitrile powder-free textured non-sterile, sizes XS/ S/ M/ L/ XL

2. Instruction for use

Note: The following product designation may be used further in this document:

The gloves

The product

2. INFORMATION ABOUT THE PRODUCER, MANUFACTURER, AND PLACE OF MANUFACTURE OF THE MEDICAL DEVICE

2.1. Producer

Blue Sail Medical Co., Ltd., No.21, Qingtian Road, Qilu Chemical Industrial Park, 255414 Zibo, Shandong Province, PEOPLE'S REPUBLIC OF CHINA.

2.2. Designer

Blue Sail Medical Co., Ltd., No.21, Qingtian Road, Qilu Chemical Industrial Park, 255414 Zibo, Shandong Province, PEOPLE'S REPUBLIC OF CHINA.

2.3. Place of manufacture

Blue Sail Medical Co., Ltd., No.21, Qingtian Road, Qilu Chemical Industrial Park, 255414 Zibo, Shandong Province, PEOPLE'S REPUBLIC OF CHINA.

3. PURPOSE, INDICATIONS, SCOPE OF THE MEDICAL DEVICE

3.1. Purpose of the medical device

The product is intended to be used as a two-way barrier to protect the hands of medical personnel from various contamination during medical research, diagnostic and therapeutic procedures, as well as to prevent cross-infection of patients and staff.

The gloves are intended for medical personnel, patients and visitors in inpatient, outpatient medical and preventive institutions, polyclinic conditions, public organizations and other social institutions, as well as at home.

3.2. Scope of a medical device

Inpatient, outpatient medical and preventive institutions, polyclinics, public organizations and other social institutions, as well as at home.

3.3. Indications for use

The product is for single use individual protection of the skin of the hands during medical research, diagnostic and therapeutic procedures in order to protect the patient and medical worker from mutual infection.

3.4. Contraindications:

Individual intolerance (hypersensitivity to the components of the product).

The residual risks associated with the use of a medical device after the realization of measures to minimize the risk are small in comparison with the benefits of using the product.

The residual risks are acceptable for the following hazards:

- non-compliance with the rules and established conditions of storage, transportation or operation;
- re-use of an already used product;
- use of the product after the end of its service life;
- violation of the established procedure for the disposal of used products;
- inattentive reading of the instructions for use;
- deliberate violation of the established order of application;
- insufficiently informative and accessible instructions for the use of the product;
- rupture of the product during sharp manipulations during putting on and operation as a result of excessively strong impact on the product of the consumer;

- manufacturing defects.

3.5. Possible side effects:

Allergic reactions to certain components of gloves are possible on the skin of people who are prone to allergic reactions, especially against the background of increased sweating.

4. DIRECTION FOR USE

1. Open the group package of gloves, remove one pair of gloves/

ATTENTION! Do not use the product if mechanical damage, ruptures, contamination is detected.

2. Consistently put the gloves on first on one hand, then on the other hand.

3. At the end of the application, grasp the glove at the base of the palm with the fingers of the other hand. Pull off the glove and hold it with the other hand. Place the fingers of your free hand in the cuff of the glove on the other hand and pull it off in the same way as the first. After turning the glove inside out, use it as a cover for another glove.

4. After use, dispose of in a Class B waste container.

5. SAFETY PRECAUTIONS

1. Nails should be cut short and neatly. The use of artificial nails is not allowed.

2. Jewellery and watches should be removed when gloves are used.

3. Gloves are intended solely for one-time use in medical research, diagnostic and therapeutic procedures in order to protect the patient and researcher from mutual infection.

4. Re-use of used gloves is not allowed.

5. Processing with alcohol-containing and aggressive chemical solutions is prohibited, because these solutions destroy the protective layer of gloves.

6. It is not allowed to use the same pair of gloves when moving from one patient to another or when moving from one manipulation to another.

7. Shouldn't use gloves in the following cases:

- after detecting a visible defect;

- after a possible or really occurred puncture;

- if the glove begins to swell, stretch, stick, or loosely sits on the hand.

8. The time of operation (application) of the gloves should not exceed 2-3 hours.

9. Immediately discontinue use if reactions are suspected. Discontinue use if reddening or irritation occurs and have a consultation with medical specialist.

6. DESCRIPTION OF MEDICAL DEVICE

The constant characteristics of gloves are the following:

- material of manufacture: synthetic nitrile rubber;
- repeatability: single-use;
- external surface: textured;
- internal surface processing: powder-free;
- cuff design: with a roller on the edge;
- form: non-anatomical, universal for both hands

The varieties characteristics of gloves are the following::

- color: black, white, blue, blue purple, cobalt blue;
- size: XS/ S/ M/ L/ XL

7. DESCRIPTION OF THE MANUFACTURE MATERIALS

The gloves are made of a composite material consisting of a base material and auxiliary additives. The composition of the manufacture material is shown in Table 1.

Table 1. Composition of the manufacture material


Name of the material	Concentration (by weight%)	Function
Synthetic Nitrile Latex KNL830	96,98%	The basis of the composite material
Zinc oxide	0,51%	Activator
Lump sulfur	0,51%	Substance forming cross-links
Titanium dioxide	0,81%	Bleaching agent
Zinc Dibutyl Dithiocarbamate Powder	0,34%	An antioxidant Suppression of oxidation
The bottled reaction product of p-cresol and dicyclopentadiene	0,26%	Accelerating additive
Potassium hydroxide	0,56%	Accelerating additive
The pigments used for gloves are listed below:	0,04%	Dyes
For blue: pigment NBR-008		
For white: pigment NBR-001		
For black: pigment NBR Black 7		
For blue purple: pigment NBR BLUE 7212		
For cobalt blue: pigment NBR BLUE 7214		
Calcium Stearate Dispersion	-	Coagulant

The main material of the product: synthetic nitrile latex KNL830 - hypoallergenic polymer rubber material. It does not leak moisture and provides increased protection from splashes of acids, alkalis, oils and other aggressive substances. It does not contain heavy metals that can harm human health. The remaining components: zinc oxide, sulfur, titanium dioxide, accelerating additives, antioxidant, potassium hydroxide and pigments, calcium stearate are non-toxic, do not contain heavy metals that can harm human health.

The product does not contain medicines, materials of animal or human origin.

Type of contact with the human body: The product has short-term contact with intact skin and mucous membranes of the human body..

8. TECHNICAL SPECIFICATION

1. Thickness measurement location

Distance:	mm
AND	13 +/- 13
B	33 +/- 5
C	48 +/- 9, (note 1)

Note 1: Distance C defines the approximate center of the palm relative to the size of the gloves.

2. Length measurement location

Distance:	mm
D	From the tip of the second finger to the outer edge of the cuff

3. Width measurement location

Distance:	mm
E	Measured at the level between the base of the index finger and thumb.

The main sizes of gloves are shown in Table 2.

Table 2. Glove sizes

The size	Width mm	*Length, mm not less	** Thickness according to the drawing, not less than, mm		
			Cuff	Palm	Finger
XS	≤80	220	0,08	0,11	0,11
S	80±10	220			
M	95±10	230			
L	110±10	230			
XL	≥110, но не 130	230			

* maximum permissible length (regardless of size) not more than 270 mm

** maximum allowable thickness (cuff, palm and finger) not more than 2.0 mm

Gloves with a round roller on the edge of the cuff, roller diameter no more than 2.5 mm.

The strength characteristics of the gloves are presented in Table 3.

Table 3. Strength characteristics

Characteristic	Value
Breaking force to accelerated aging, N, not less	7,0
Elongation at break to accelerated aging, %, not less	500
Breaking force after 3 months from the date of manufacture, N, not less	6,0
Elongation at break after 3 months from the date of manufacture, %, not less	400
Breaking force after accelerated aging, N, not less	6,0
Elongation at break after accelerated aging, %, not less	400

Powder-free gloves have a residual amount of dusting agent < 2.0 mg per glove in accordance with EN 455-3.

Gloves are airtight, level of quality of tightness (AQL) no more than 1,5

The weight of one glove is not more than 10 g (for all sizes)

9. COMPLETENESS OF DELIVERY

The scope of delivery of gloves should include:

- gloves of the same color and size in a group/transport package;
- instructions for use.

10. LABELING AND PACKAGING OF THE PRODUCT.

10.1. Description of the marking

Marking on consumer packaging of gloves contains:

- name of the product;
- symbol "Producer", name and location of the producer;
- location of the place of manufacture;
- the name, location of the authorized representative of the producer;
- symbol "Don't contain natural rubber latex" and the inscription: "Does not contain natural rubber latex";
- symbol and the inscription: "The product is non-sterile";
- symbol "Don't re-use" and the inscription: "Don't re-use! Disposable gloves!";
- symbol "Keep away of sunlight";
- symbol "Consult instructions for use";
- symbols "Keep dry"; "Temperature limit";
- symbol "Don't use if package is damaged";
- the number and date of registration certificate;
- reference to standard GOST R 52239-2004;
- symbol "Catalogue number" and catalogue number;
- symbol "Batch code" and batch code;
- symbol «Date of manufacture» and the date (month, year) of manufacture;
- symbol «Use by date» and the expiration date (month, year);
- size and quantity of gloves;
- information on the color of the gloves;
- barcode (if available);
- information about the conditions of storage.

10.2. Description of symbols

Table 4. Symbols used in labeling

Symbol	Title of symbol
	Manufacturer
	Date of manufacture
	Temperature limit
	Consult instructions for use
	Non-sterile
	Don't re-use

	Catalogue number
	Batch code
	Use by date
	Don't use if package is damaged
	Keep away of sunlight
	Keep dry
	Don't contain natural rubber latex

10.3. Packaging:

One consumer package contains 50 pieces (25 pairs) or 100 pieces (50 pairs) of gloves.

Material of manufacture of consumer packaging: coated double-layer cardboard.

Dimensions of a consumer cardboard box (Length x Width x Height, mm, ± 10 mm): 250x130x80.

Package weight with gloves (50 pieces): no more than 500 g (regardless of the version) Package weight with gloves (100 pieces): no more than 1000 g (regardless of the version).

The box allows you to maintain the quality of gloves during transportation and storage.

One transport package contains 10 consumer packages.

Material of manufacture of transport packaging: goffered cardboard.

Dimensions of the transport cardboard box (Length x Width x Height, ± 10 mm.): 500x300x400.

Package weight, gross: no more than 10.0 kg.

11. DISPOSAL REQUIREMENTS

The used product may pose an infectious, microbial, or environmental hazard due to possible infection with infectious agents: viruses, bacteria, and microbial-infested biological fluids. The products used for their intended purpose are disposed of in accordance with the requirements of SanPiN 2.1.3684-21 for medical waste of class B.

Unused products with an expired shelf life are disposed of as medical waste of class A according to SanPiN 2.1.3684-21. The products are recyclable for secondary use. The rules of delivery and processing should be clarified with the regional processors of recyclable materials.

The used product at home is disposed of as ordinary solid household waste after packaging into plastic bag according to SanPiN 2.1.3684-21 in accordance with the territorial waste management scheme.

12. REQUIREMENTS FOR TRANSPORTATION, STORAGE AND OPERATING CONDITIONS

12.1. Transportation conditions

Products must be transported only in packaging by all types of closed vehicles and according to the rules of cargo transportation applicable to each type of transport. Temperature during transportation - from -50 ° C to +50 ° C. Relative humidity: up to 80%.

After transportation at sub-zero temperatures, unpack, having stood for 2 hours at room temperature.

12.2. Storage conditions

The products must be stored in dry, ventilated storage areas, protected from direct sunlight and atmospheric influences, and located at a distance of at least 1 m from heating devices.

The storage temperature of the products should be from +5 ° C to +40 ° C. Relative humidity: up to 80%

12.3. Operating conditions

The product must be operated at a temperature of +10 ° C to +35 ° C and a relative humidity of up to 80%.

13. REQUIREMENTS FOR CLEANING, DISINFECTION AND STERILIZATION

The products are delivered in a non-sterile form and are not subject to cleaning, sterilization and disinfection.

14. MAINTENANCE

Gloves are for single use products and are not subject to maintenance.

15. LIST OF ACCESSORIES

Product is supplied without accessories.

16. PRODUCER'S WARRANTY

The product is designed, manufactured, tested and packaged in compliance with all standards.

The company provides a warranty for this product until the expiration date, provided that its packaging is not broken, opened or damaged.

The producer guarantees that the product meets all the requirements of the technical documentation, subject to the conditions of use, transportation and storage.

Warranty period of storage of gloves is 5 years from the date of manufacture.

17. INTERNATIONAL STANDARDS USED IN PRODUCTION

- COUNCIL DIRECTIVE 93/42/EEC of 14 June 1993 on medical devices, Medical Device Directive (MDD) Directive 2007/47.
- ISO 15223-1, Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
- ISO 13485, Medical devices - Quality management systems - Requirements for regulatory purposes
- ISO 14971, Medical devices - Application of risk management to medical devices
- ASTM D 6319-05 Standard specification for nitrile examination gloves for medical application
- ISO 10993-1-2011 Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and testing
- ISO 10993-5, Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
- ISO 10993-10, Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
- EN 455-1, Medical gloves for single use - Part 1: Requirements and testing for freedom from holes
- EN 455-2, Medical gloves for single use - Part 2: Requirements and testing for physical properties
- EN 455-3, Medical gloves for single use - Part 3: Requirements and testing for biological evaluation
- EN 455-4, Medical gloves for single use - Part 4: Requirements and testing for shelf life determination

18. COMPLAINT

Complaints are subject to:

- products with defects detected during operation, storage and transportation during the established warranty period;
- products with non-compliance of packaging, packaging, labeling, and completeness with the requirements of the accompanying documentation.
- information about side effects and adverse events during the use of gloves

No complaints are made:

- after the expiration of the warranty period for the product;
- if the detected defects were the result of non-compliance by the consumer with the conditions and rules of use, storage and transportation.

Complaints should be sent to the address of the authorized representative of the producer in the territory of the Russian Federation, LLC «TDE».

19. AUTHORIZED REPRESENTATIVE OF THE PRODUCER IN THE RUSSIAN FEDERATION

LIMITED LIABILITY COMPANY «Torgoviy Dom EVRASIA» (LLC «TDE»), 125080, Moscow, Volokolamsk road, 2, floor 23, apt. I, room 6
tel.: +7 (495) 744 32 82; e-mail: info@mediok.ru

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