**WORLD WORK SRL**Via Del Progresso, 47 - 36054 Montebello Vicentino
(VI), Italy

Tel.: +39 0444 574297

Email: info@worldwork.it

APPROVED BY

General Manager

Mrs. Daniela Baldissera


36054 MONTEBELLO VICENTINO (VI)

Partita IVA 01679160244

Tel. +39 0444/574297 - Fax +39 0444/37054

INSTRUCTION FOR USE
of medical device**Dental disposable saliva ejectors DENTAL FORMULA®**

Rev: 00 (June 2021)

1. NAME

Dental disposable saliva ejectors DENTAL FORMULA®

design versions: Dental disposable saliva ejectors DENTAL FORMULA® with removable tip, length 15 cm, color: clear, light blue, blue, lilac, fuchsia, orange, green, mint green.

Manufacturer and producer of the medical device

Manufacturer WORLD WORK SRL
Via Del Progresso, 47 - 36054 Montebello Vicentino (VI), Italy
Address of the manufacturing facility Via Del Progresso 47/49, Municipality of Montebello Vicentino (VI), Italy

Authorized representative of the Manufacturer on the territory of the Russian Federation

"BE Plus" Company Limited
"BE Plus" Co. Ltd.
Russia, 194292, St. Petersburg, ext. ter. c. municipal district Sergievskoe, 12, 6-y Verkhniy per, litera A, part of office 5.3.1. Tel.: +7 (812) 309-36-62, e-mail: beplus@beplus.spb.ru

The scopes of delivery of the medical device: - 250 pcs in a consumer package;
- 100 pcs in a consumer package;
- 40 pcs in a consumer package;
- 15 pcs in a consumer package;
- 5 pcs in a consumer package.

2. PURPOSE

The saliva ejectors are designed to remove liquids and suspensions from the patient's oral cavity during the dental procedures as part of an aspiration system.

The saliva ejectors can be used by different specialists involved during the dental procedures. The use of the saliva ejectors does not require special qualification or special education.

3. FUNCTIONAL SPECIFICATION

The operational principle of the medical device is to supply negative pressure generated by an aspiration system to various parts of the patient's oral cavity without damaging the tissue of the mucous membrane. Thus, the medical device ensures the dryness of the surgical field, the safety, and comfort of the patient, preventing the swallowing of potentially hazardous dental solutions and suspensions, as well as the leakage of saliva from the patient's oral cavity, whose swallowing function is impaired during the procedure.

Scope of the ejectors – using as part of the aspiration system for all types of dental procedures in medical or healthcare facilities.

4. RISK OF USING, CONTRAINDICATIONS AND SIDE EFFECTS

Potential risk of using of the saliva ejectors – Class 2a.
Contraindications and side effects of the medical device use are not identified.

The shelf life of the saliva ejectors is 5 years.

5. TECHNICAL SPECIFICATION

The cannula of the ejector is a phthalate-free PVC tube and is reinforced with carbon steel wire throughout the length of the tube. The filter element is a removable tip made of phthalate-free PVC or low-density polyethylene.

Tube is manufactured in the following colors: clear / light blue / blue / lilac / fuchsia / orange / green / mint green.

Tip: clear.

Total ejector length (tube + tip): 150 ± 10 mm.

Ejector (tube) outer diameter: 6.5 ± 2 mm.

6. APPLICATION PROCEDURE

Before use, check the tightness of the tip on the cannula using a slight force. The tip should fit snugly at the end of the cannula.

Install the ejector into the saliva ejector hose adapter, if any, or directly into the suction circuit of the pneumatic or hydraulic system. Shape the ejector according to the patient's oral cavity. Bend the ejector smoothly, avoiding sharp blocking kinks that reduce the efficiency of the saliva ejector.

If it is necessary to make a strong bend of the ejector, make it in two stages. Using the ball of your thumb as a stop, make a smooth first fold, then make the second one, rotating the ejector around its axis.

7. APPLICATION CONDITIONS

The ejectors must be used as per the instruction on the medical device. The use of the device does not require special qualification or special education.

The saliva ejectors are used in medical or healthcare facilities indoors, in a dry place, away from direct sunlight, at a temperature of 15°C to 35°C and relative humidity from 45% up to 85%.

The device is intended for single use only. For individual use only. When using the device, it is necessary to protect it from mechanical damage. Avoid sharp blocking kinks in the ejector.

8. TRANSPORTATION, STORAGE AND DISPOSAL CONDITIONS OF THE MEDICAL DEVICE

The transportation of the device may be carried out by all types of covered vehicles in the original package as per the transport regulations, which are applied to this type of transport. During transportation the storage conditions must not be violated.

The devices must be stored indoors at a temperature of 0°C to 35°C away from heating appliances under conditions that exclude exposure to aggressive environments (acidic, alkaline, etc.), as well as flammable and hot liquids. The devices must be protected from direct exposure to sunlight and high humidity (relative air humidity up to 85%).

An unused medical device that has become useless, including due to the expiration date, belongs to the category of General waste (Class A waste according to the regulatory guidance of the Russian Federation SanPiN 2.1.3684-21).

The used medical device belongs to the category of Infectious waste (Class B waste

according to the regulatory guidance of the Russian Federation SanPiN 2.1.3684-21). Collection, temporary storage, and waste disposal should be carried out following the medical waste management scheme adopted in the given organization carrying out medical and/or pharmaceutical activities.

9. MANUFACTURER'S WARRANTIES

The manufacturer guarantees the safety of all parameters and characteristics of the saliva ejectors for 5 years from the manufacturing date, subject to the conditions of transportation, storage and application set forth in the instruction.

10. LIST OF NATIONAL STANDARDS, APPLIED BY MANUFACTURER (PRODUCER) OF THE MEDICAL DEVICE

UNI CEI EN 1041: Information supplied by the manufacturer of medical devices
2014

EN ISO 10993-1:2018 Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process

ISO 10993-5:2009 Medical devices. Biological evaluation of medical devices — Part 5: Tests for in vitro cytotoxicity

ISO 10993-10:2010 Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization

ISO 13485:2016 Medical devices – Quality management systems – Requirements for regulatory purposes

UNI CEI EN ISO 14971:2015 Medical devices – Application of risk management to medical devices

UNI CEI EN ISO 15223-1:2017 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements

CE European Conformity. Requirements

93/42/EEC Medical Device Directive

Instruction for use of the medical device: Dental disposable saliva ejectors DENTAL FORMULA® is released for the first time (Rev: 00 (June 2021)).

In case of an undesirable event showing the signs of an adverse event (incident), it is necessary to notify thereof the manufacturer or an authorized representative of the manufacturer in the Russian Federation.

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