“APPROVE”

CEO

(position)

Farideh Bischoff


(name)

[Signature]

(signature)

Stamp

User Manual

[Signature] 05/17/2022

ZipThaw®202 rapid plasma thawing system

FDA
CLEARED

Fremon Scientific, Inc.
5726 La Jolla Blvd., Suite 304
La Jolla, CA 92037
USA

Authorized representative of the manufacturer in the Russian Federation for complaints and communications
(name, address, telephone, e-mail)

Individual entrepreneur Ratnesh Prashant
Apt. 42, 23/16 Sovetskaya Street, Yaroslavl, Yaroslavl region, Russia, 150003
+7-910-664-28-62
Prashant Ratnesh shreyayar@yahoo.com

TABLE OF CONTENTS

TABLE OF CONTENTS	3
TABLE OF FIGURES	5
INTRODUCTION	6
3.1 Indications for use	6
3.2 Important!	6
3.3 Terms and abbreviations	7
DESCRIPTION OF THE DEVICE	8
4.1 ZipThaw®202 system modules	8
Thawing chamber module (TCM)	8
4.1.2 System control processor (SCP)	8
4.1.3 Chamber module controller (CMC).....	8
4.1.4 ZIPSleeve™	8
4.1.5 Power supply	8
4.2 ZipThaw®202, front view	10
4.3 ZipThaw®202, back view.....	10
4.4 ZipSleeve®	10
4.5 Operating principle.....	10
4.6 Temperature default settings	10
4.7 Marking and symbols	11
4.7.1 Back panel	11
4.7.2 Symbols used in this User Manual.....	11
4.8 Compliance with the Rules and Regulations of the US Federal Communications Commission, Title 47 C.F.R. Part 18	12
CONTRA INDICATIONS	12
DEVICE SAFETY	13
OPERATION	13
7.1 Placing of the device and power supply connection	13
7.2 Door position	14
7.3 Turning on	14
7.3.1 Turning on failure	14
7.3.2 Failure of internal diagnostic device	15
7.3.3 System recovery after reset	15
7.3.4 Successful power up	16
7.4 Description of home screen	17
7.5 Settings	17
7.5.1 Description of keys	18
7.5.2 Language preference	18

7.5.3 Selection of temperature unit	18
7.5.4 Site name input	19
7.5.5 Date and time setting	19
7.5.6 Sound on	20
7.5.7 Adding new user.....	20
7.5.8 Edit user list	21
7.6 Setting frequency for replacement of heating cushion	22
7.7 Start operation	23
7.8 Input of data for thawing	23
7.8.1 Input of bar code data	23
7.8.2 Manual data input	24
7.8.3 Override of all data input for processing speed-up	24
7.9 START-UP OF THAWING CYCLE	25
7.9.1 Loading a bag with frozen plasma into ZipSleeve®	25
7.9.2 Loading ZipSleeve® unit into the chamber	26
7.9.3 Choosing MELT or THAW options.....	27
7.9.4 Start	27
7.9.5 Intended stop of thawing	28
7.9.6 Thawing using the second chamber B.....	28
7.10 Completion of thawing cycle	29
7.11 Waiting for thawing cycle	30
DATA EXPORT 31	
8.1 Export of data structure	32
POWER OFF 33	
9.1 Standby mode	33
9.2 Off mode	33
MESSAGE FOR USER 35	
CLEANING 43	
11.1 Cleaning frequency	43
11.2 Cleaning and disinfection of external surface	43
11.3 Screen cleaning	43
11.3 Cleaning of spill tray	43
WARRANTY 43	
CUSTOMER SUPPORT SERVICE 44	
SPECIFICATIONS 45	
PROTECTION AGAINST HARMFUL INTERFERENCE 46	
15.1 EMC safety	46

LIST OF FIGURES

FIGURE 1.	ZIPThaw® 202, back view	9
FIGURE 2.	ZIPThaw® 202, back view	9
FIGURE 3.	Appearance of ZipSleeve®	10
FIGURE 4.	ZIPThaw® 202 label	11
FIGURE 5.	COMPLIANCE WITH TITLE 47 C.F.R., PART 18	12
FIGURE 6.	Open door pop-up error screen	14
FIGURE 7.	ZIPThaw® 202 activation screen	14
FIGURE 8.	System error pop-up screen	15
FIGURE 9.	Chamber error pop-up screen	15
FIGURE 10.	Chamber error home screen	15
FIGURE 11.	Home screen	16
FIGURE 12.	Home screen	17
FIGURE 13.	Entering the settings screen	17
FIGURE 14.	Settings screen	18
FIGURE 15.	Language preferences	18
FIGURE 16.	Temperature units.....	19
FIGURE 17.	Site name input.....	19
FIGURE 18.	Date and time setting	20
FIGURE 19.	Sound on	20
FIGURE 20.	Adding new user screen	20
FIGURE 21.	Adding new user	21
FIGURE 22.	Heating cushion replacement screen	22
FIGURE 23.	Cushion replacement reminder screen.....	23
FIGURE 24.	Chamber selection on home screen.	23
FIGURE 25.	Bar code input screen	24
FIGURE 26.	Manual data input screen	24
FIGURE 27.	Error of data input into mandatory fields and Override of data input	25
FIGURE 28.	Home screen after data input (chamber A).....	26
FIGURE 29.	ZipSleeve® orientation	26
FIGURE 30.	Device door open.....	27
FIGURE 31.	Replacement of ZipSleeve® in accordance with the screen of number of uses	27
FIGURE 32.	Home screen with loaded chamber A – start parameters; Melt (until ice free state) or Thaw	28
FIGURE 33.	Thawing chamber A screen	29
FIGURE 34.	Prompt of intended stop of thawing cycle	29
FIGURE 35.	Screen of chamber A thawing cycle, chamber B readiness screen ...	30
FIGURE 36.	Chamber A thawing cycle completed, chamber B loaded and active	30
FIGURE 37.	Waiting screen.....	31
FIGURE 38.	Data export screen	32
FIGURE 39.	Data load screen	32
FIGURE 40.	First screen of device turning off prompt	34
FIGURE 41.	DEVICE TURNING OFF PROMPT SCREEN	34

INTRODUCTION

This document contains the operating instructions and technical details of ZipThaw®202 system based on system requirements that were obtained from several sources including customer preferences, marketing requirements and previous generations of systems.

These instructions are intended for use by healthcare professionals at hospitals.

Name of medical device

ZipThaw®202 rapid plasma thawing system, consisting of:

1. ZipThaw®202 blood plasma thawer - 1 pc;
 2. Power cord - 1 pc;
 3. Consumables for ZipThaw®202 blood plasma thawer: ZipSleeve® polymer bags - from 1 to 40,000 pcs
 4. User manual - 1 pc;
 5. Warranty certificate - 1 pc.
- (hereinafter referred to as device or medical device)

3.1 INDICATIONS FOR USE

The FMS ZipThaw®202 is a thawing device intended for the following applications:

- thawing of Fresh Frozen Plasma (FFP)
- thawing of Plasma Frozen within 24 hours after Phlebotomy (PF24)

CAUTION!

The device cannot be used for heating of previously thawed cryoprecipitated antihemophilic factor, platelets and granulocytes.

This device is not intended to thaw, among others, red blood cells, hematopoietic progenitor cells or cord blood, except for research use.

3.2 IMPORTANT!

The following icons, symbols, and words are used in these instructions to highlight warnings and special recommendations:

 ATTENTION!	Attention! warning message warns the user of potential danger to persons associated with improper use of the device.
 CAUTION!	Caution! warning message warns the user of a potential problem with the device that could pose risk to property, the device itself, or its essential functions

These instructions are intended for use by healthcare professionals.

When using the device, be sure to follow the instructions and take into account safety information contained in this user manual. This will eliminate risk to patients and users and prevent damage to equipment.

FreMon Scientific assumes no responsibility for damage resulting from failure to follow these instructions.

ATTENTION!

These instructions for use should be kept throughout the service life of the device and handed over to any subsequent user.

ATTENTION!

The electromagnetic compatibility (EMC) safety measures are assumed during operation of electrical medical devices. Installation and operation of the device must be carried out in accordance with EMC guidelines contained in these instructions.

3.3 TERMS AND ABBREVIATIONS

Temp.	Temperature	DC	Direct current
VTR	Variable thawing rate	W	Watt
GUI	Graphical User Interface	V	Volt
EIS	Electronic Information Systems	A	Ampere
Min	Minute(s)	HMI	Human Machine Interface
FEC	Frontend controller	Chamber	Place for processing of bag with product
AC	Alternating current	MCS	Modular computer system
CMC	Chamber module controller	TCM	Thawing chamber module
SCP	System control processor	BDD	Built-in diagnostic device: diagnostic program that checks the system health.

DESCRIPTION OF THE DEVICE

ZipThaw®202 system platform has a state-of-the-art modular architecture that provides configuration flexibility.

4.1 ZIPTHAW®202 SYSTEM MODULES

The ZipThaw®202 system consists of universal basic models:

4.1.1 THAWING CHAMBER MODULE (TCM)

This is a universal module that uses a closed heating chamber, mixing mechanism and temperature monitoring system. It can be connected to ZipSleeve® mechanically and wirelessly.

TCM includes two sealed heating cushions, each filled with 300 ml of water-based liquid. Each cushion is attached to a heating element and controlled by the temperature sensor and chamber module controller (CMC). Moreover, each heating element has its own heating sensor for overheating protection.

General overheating protection is achieved through the use of independent temperature sensor that cuts off power in case of cushion overheating.

4.1.2 SYSTEM CONTROL PROCESSOR (SCP)

System control processor (SCP), which controls 4 CMC. **SCP** is a high-level central processing unit that controls and monitors system functions, controls a graphical interface module with touch screen and an external memory system.

4.1.3 CHAMBER MODULE CONTROLLER (CMC)

Chamber module controller: universal chamber control module. It carries out control of one chamber unit (i.e. sensors, mixing, limit switches and heating) and data exchange with SCP for data/status synchronization.

4.1.4 ZIPSLEEVE®

A protective cover with a built-in wireless temperature sensor that measures the temperature of the bag with thawed product and sends a signal to CMC to stop heating as soon as the sensor temperature reaches the "target temperature" of 31 °C (87.8 ° F) for THAW option or 15 °C (59.0 ° F) for MELT option to *ice-free* condition.

4.1.5 POWER SUPPLY

A universal medical AC/DC power supply, which is suitable for international power systems and provides the power required for operation of the modules with 2 pockets.

4.2 ZIPTHAW®202 FRONT VIEW

FIGURE 1. ZipThaw®202, front view

4.3 ZIPTHAW®202, BACK VIEW

FIGURE 2. ZipThaw®202, back view

4.4 ZIPSLEEVE® :


FIGURE 3. Appearance of ZipSleeve®

4.5 OPERATING PRINCIPLE

A bag with frozen plasma is placed inside ZipSleeve®, in a “smart” shell, bag and inserted into the thawing chamber module (TCM).

During the thawing process, the plasma bag is thawed using the heating cushions of the chamber. The system monitors the temperature of the heating cushions to maintain the heating temperature so that the temperature of the bag with product would not exceed 37 °C (98.6 °F).

The temperature of the frozen plasma bag is monitored by ZipSleeve™ temperature sensor (RFID) which signals the system to stop heating as soon as the sensor reaches a “target temperature” of 31 °C (87.8 °F) for THAW option or 15 °C (59.0 °F) for MELT option. Moreover, ZipThaw™ system monitors the mixing mechanism to ensure uniform mixing and heat distribution.

4.6 TEMPERATURE DEFAULT SETTINGS

The final “Target temperatures” for THAW and MELT options are fixed as factory settings. The ZipSleeve® sensor monitors the temperature of the frozen plasma bag and signals the system to stop heating as soon as the target temperature is reached.

Note: The user cannot change the “Target temperature” but can choose between two start options: “Thaw” and “Melt until ice-free”.

4.7 MARKING AND SYMBOLS

4.7.1 Back panel


FIGURE 4. Label of ZipThaw®202

4.7.2 Symbols used in this User Manual:

No.	Symbol	Description
1		CE MARK
2		ATTENTION, PLEASE READ THE ACCOMPANYING DOCUMENTATION
3		WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT CONFORMITY SYMBOL (WEEE)
4		MANUFACTURER
5		AUTHORIZED REPRESENTATIVE IN THE EUROPEAN UNION
6		FUSE 15A/250V
7		SERIAL NUMBER
8		CATALOGUE NUMBER

No.	Symbol	Description
9		THE FEDERAL LAW RESTRICTS THIS DEVICE FOR SALE BY OR ON BEHALF OF A LICENSED PROFESSIONAL.
10		WARNING ABOUT HIGH TEMPERATURE
11		UL MARK

4.8 COMPLIANCE WITH THE RULES AND REGULATIONS OF THE US FEDERAL COMMUNICATIONS COMMISSION, TITLE 47 C.F.R. PART 18

This device complies with Part 18 of FCC Rules and Regulations. The following two conditions must be met during operation:

- 1 This device must not cause harmful interference.
- 2 This device must accept any interference including interference that may cause malfunctioning.

FIGURE 5. COMPLIANCE WITH TITLE 47 C.F.R., PART 18

CONTRA INDICATIONS

Except for the substances described in the section “Purpose” (see Section 3.1), ZipThaw®202 is not intended to thaw, heat, or keep warm any blood products or other substances including but not limited to red blood cells, hematopoietic progenitor cells, cord blood, previously thawed cryoprecipitated antihemophilic factor, platelets or granulocytes. There are no contraindications to thawing or heating blood products.

DEVICE SAFETY

ATTENTION!

- Please read the instructions contained in this user manual and any other user documentation before using the device.
- This device is referred to class 1 electrical equipment.
- Use this device only in accordance with the instructions given in the user manual.
- It is not allowed to exceed the heating temperature and the time limit for heating blood and blood products.
- If the device is faulty, do not use it and contact FreMon Scientific.
- The electrical installations must comply with applicable codes and standards.

- To connect to the power supply, use the power cord supplied by FreMon Scientific only.
- To ensure complete isolation of the device from the power supply, the power cord plug must be removed from the outlet.
- Except for replacing the heating cushions, do not attempt to perform maintenance of the device yourself. Maintenance operations and modifications may only be performed by Fremon Scientific or authorized personnel of the customer.
- Strong magnetic fields may cause malfunction of ZipThaw®202. If electromagnetic interference occurs, increase the distance between ZipThaw®202 and the device causing interference. ZipThaw®202 operates within the range of limit values specified in the standard GOST R IEC 61326-1 and can be affected outside the range of these limit values.
- The use of the product is carried out solely with gloves registered in the Russian Federation in accordance with the established procedure.
- Portable and mobile carrier equipment may interfere with operation of the device.
- Do not use the device near:
 - Flammable materials
 - Flammable mixtures of anesthetic agents with air, oxygen or nitrous oxide with a flashpoint below 50 °C (122 ° F).

OPERATION

7.1 PLACING OF THE DEVICE AND POWER SUPPLY CONNECTION

- There must be enough space to place ZipThaw®202 on a bench or cart.
- Leave at least 10 cm (4 inches) between the back of ZipThaw®202 and the wall for easy access to the main On/Off switch located as part of the outlet.
- Leave at least 40 cm (15.7 inches) of free space on each side for loading and unloading of the bags in and out of the chambers.
- Connect the power cord to the rear connector of ZipThaw®202 and plug the cord into the power outlet.
- Always use the power cord supplied by FreMon Scientific to connect to the main power supply source.

7.2 DOOR POSITION

If any chamber door is open or not closed properly during system operation, the following error message will be displayed:

FIGURE 6. OPEN DOOR POP-UP ERROR SCREEN

Close the chamber door and press "OK", the message will disappear.

7.3 TURNING ON

Switch the power button on the back side of the device to ON position to supply power to the system. (See [Figure 2](#)).

The system runs self-diagnostic and start-up (warm-up) processes, which may take up to 3 minutes.

FIGURE 7. ZIPThaw™ 202 activation screen

7.3.1 TURNING ON FAILURE

If the self-diagnostic and start-up processes fail, an error alarm will appear for 1 second and a system error message will be displayed on the touch screen.

Turn off the power at the back of the device and unplug the power cord for 5 minutes. Turn on the device again. If the problem persists, please contact your FMS distributor or us at info@fremonscientific.com, and seek assistance through our website <https://fremonscientific.com>!

***System Error! FMS Technical Support Required**

FIGURE 8. SYSTEM ERROR POP-UP SCREEN

7.3.2 Failure of internal diagnostic device

During self-diagnostic tests, the system will determine if both chambers are working or one or both chambers are faulty and must be disabled. The system will disable the faulty chamber(s) but will allow full operation of the device with one chamber running.


The following error message will appear:

* Chamber B Error, Use Chamber A

FIGURE 9. CHAMBER ERROR POP-UP SCREEN

If you press OK, the error message will disappear and the next screen will be displayed, on which the faulty chamber will be disabled and will be displayed in grey.

* Select Chamber

Error – Inactive

Chamber A

FIGURE 10. CHAMBER ERROR HOME SCREEN

▪ 7.3.3 System recovery after reset

If the system enters unexpectedly the Reset mode by default, the home screen will appear after a successful internal diagnostics, and the system will return to pre-thawing state, see clause 7.3.2.

Note: The user decides to start a new thawing cycle using these bags or not.

To start a new thawing cycle, go to clause 7.3.2.

▪ 7.3.4 Successful power up

After a successful preheating cycle the home screen will appear. A query will appear on the screen: “**Select Chamber**” as shown below.


FIGURE 11. HOME SCREEN

7.4 DESCRIPTION OF HOME SCREEN


FIGURE 12. Home screen

- | | | | |
|---|---|---|--------------------------------------|
| 1 | FreMon logo + site name | 6 | Basic action keys. |
| 2 | Data export. | 7 | Instructions for the active chamber. |
| 3 | Settings. | 8 | ZipSleeve® usage indicator. |
| 4 | Date and time. | | |
| 5 | Work area A + B is the main area of user interaction. | | |

7.5 SETTINGS

Pressing the “**Settings**” icon will open the settings access screen. **The following screens are intended for the device administrator and technical support only.**

The administrator password is required to enter the settings screen. The administrator will have a default password (45253, factory default value) that can be changed in the settings screen.

To enter the **Technical Screens**, press the **Settings** icon and enter number “2” as access code. These screens will be useful for troubleshooting and setting up replacement of a heating cushion (see section 7.6 below).


FIGURE 13. Entering the settings screen

Note: If the administrator password is forgotten, only a certified professional can reset it.

After entering the correct **administrator password** and pressing **OK**, the following settings screen will appear:


FIGURE 14. Settings screen

7.5.1 Description of keys

HOME

When pressing this key a return to the home screen without saving will occur.

DEFAULT SETTINGS

Pressing this key will return to default settings.

SAVE


Pressing this key will save the settings and return to home screen.

ENTER key.

7.5.2 Language preferences

Tapping the displayed flag key opens a drop-down menu with additional flag images allowing the user to select the language represented by the flag.

English
Russian
French

Figure 15. Language preferences

When you tap the selected flag, the drop-down menu will close, the language will change, and the selected flag will be displayed.

Note: The default language is English represented by the US flag icon.

7.5.3 Selection of temperature unit

The administrator can select the temperature unit (for example, Celsius or Fahrenheit) by pressing the corresponding toggle button.

Figure 16. Temperature units

7.5.4 Site name input

To enter a site name, touch the screen in the New Site field to select it and use the keyboard to enter the site name. The typed name will appear in the header at the top left next to FreMon Scientific logo.


Figure 17. Site name input

Note: This field is limited to 25 characters.

7.5.5 Date and time setting

Touch the Date field to set the date. A drop-down menu will open. Select the date format you want (for example, DD/MM/YY or MM/DD/YY), then use the on-screen keyboard to enter the current date.

Touch the Time field to set the time. The time format is 24:00. Use the keyboard to enter the current time.


Figure 18. Date and time setting

Note: If the information entered is incorrect the field will be marked in red

7.5.6 Sound on

To enable sound alerts, turn Sound on or off by sliding the corresponding switch.

Figure 19. Sound on

7.5.7 Adding new user

Touching the Users field allows entering a new username by typing it in the New name field and ID number in the ID/bar code field using the on-screen keyboard, or by scanning the new user ID badge.


FIGURE 20. Adding new user screen

Pressing the ID badge/bar code turns on the built-in bar code reader (i.e. you will see the red light on the bar code reader light up) and the user can either scan the ID or enter it using the on-screen keyboard.

Press for new user  confirmation - the new name and ID will be added to the list.

7.5.8 Edit user list

- Check Users (i.e. tap) and a drop-down list will appear.
- To remove a username from the list, press X to the right of the user ID.
- Press ENTER to add a new user.
- When finished editing the entire screen, press SAVE to save the input data.

- Press **HOME** key button to return to home screen.


FIGURE 21. Adding new user screen

1.6 SETTING FREQUENCY FOR REPLACEMENT OF HEATING CUSHION


Figure 22. Heating cushion replacement screen

A reminder schedule for replacement of the heating cushions is set up on this screen. It is configured by the user depending on the expected level of ZipThaw202 usage. The manufacturer may replace the cushions at its discretion in advance at least once a year. They can be replaced earlier, that is, at intervals of 3 or 6 months or with maximum cycles depending on the higher intensity of use. The relevant warning screen is just a reminder. The use of the cushion is not limited by the software. The manufacturer recommends replacement as soon as possible. You can continue using the cushions as may be necessary as long as they are not damaged or leaking.

- Date of the last cushion change: it sets the current date as the date of the last cushion change.
- Next Replacement Date: it sets the current date plus the designated replacement interval (see below).
- Replacement interval: Maximum number of cycles = the expected maximum allowable number of cycles before replacing the heating cushions.
- Replacement interval: 3 months, 6 months or 12 months depending on the intended method of use
- Cycle Usage A: the actual number of cycles for chamber A since the last change of the heating cushion.
- Cycle Usage B: the actual number of cycles for chamber B since the last change of the heating cushion.
- RESET key: it resets the dates for replacement interval and number of usage cycles for chambers A and B

Figure 23. Reminder for cushion replacement “Error alert”

7.7 START OPERATION

Select chamber A or B in the home screen by touching the screen area as shown below.


FIGURE 24. Select the chamber on the home screen

Selecting a chamber will open the Insert Data Screen.

7.8 INPUT OF DATA FOR THAWING

Data can be entered:

- Using built-in bar code reader.

► Manual data input.

7.8.1 Input of bar code data

Place the plasma bag before the bar code reader and aim the red beam at the top bar code field (1) of the plasma bag label. As soon as you hear a beep, move to the next bar code field (for example, 2, 3, 4), one at a time. Make sure that the correct data have been entered, especially the donor identification number (DIN) (i.e. Field marked with 1). The user ID should be entered by scanning the user ID or manually as described below (7.7.2).


FIGURE 25. Bar code, data input screen

7.8.2 MANUAL DATA INPUT

To enter data manually, touch the relevant field: The following screen will appear:

**Mandatory fields for data input
are marked with "*"**
(e.g. User Name/ID, Donation
ID)

FIGURE 26. Manual data input screen

- Touching the field will open the keyboard. Pressing **OK** will insert the data, close the keyboard and move to the next field.
- User ID: Up to 25 characters when scanning. 4 characters are allowed (user ID).
- DIN: Up to 16 characters.
- ABO/RhD: 4 characters.

- ▶ Product code: 8 characters.
- ▶ Expiration date - 10 characters.

7.8.3 OVERRIDE OF DATA INPUT

After data input

- Press **OK** key to accept data input and return to home screen.
- Press **CANCEL** to cancel and return to home screen.

If **OK** key was pressed but the required field was not filled in, the system will mark the empty field so that the user could enter the data and a pop-up screen with the following warning will be displayed. "Data is Incorrect or Missing, Continue?" If you want to continue without entering any data, press **YES**. If you want to return and fill in the missing data press **NO**.


FIGURE 27. Error of data input into mandatory fields and Override of data input

If there are no data input errors, or if you decide to continue without entering data, and if the **OK** key was pressed, the system will return to home screen by default, as shown below:


FIGURE 28. Home screen after data input (chamber A)

Load the frozen plasma unit into ZipSleeve. Section 7.8.1 below.


7.9 START-UP OF THAWING CYCLE

7.9.1 Loading a bag with frozen plasma into ZipSleeve®

- 1 Place the frozen plasma bag inside the ZipSleeve™.
- 2 Insert the frozen bag with the label facing the transparent side of the ZipSleeve™.

Note:

Once the ZipThaw™ system verifies that a valid ZipThaw™ unit is placed inside the device (i.e. the number of uses is less than or equal to 8), this will allow the device to start a thawing cycle.

Note:
The maximum number of uses
of the ZipSleeve™ is 8

FIGURE 29. ZipSleeve™ orientation.

7.9.2 Loading ZipSleeve® unit into the chamber

1. Open the chamber door and place the loaded ZipSleeve® inside the chamber.
2. Align 3 slots on the top of the ZipSleeve® with the 3 holders inside the chamber.

Note: Note that the ZipSleeve® can only be loaded in the orientation indicated by the slots on the rigid plastic plank at the top of the ZipSleeve®.


FIGURE 30. Device door open

3. Close the chamber door.

After closing the chamber door, the ZipThaw® system will check that a suitable ZipSleeve is placed inside the device (i.e. less than or equal to 8 uses) before activating the START key.

If the ZipSleeve™ unit has been used 8 times, the system will display the following screen:


FIGURE 31. Replace ZipSleeve* in accordance with the screen of number of uses.

Replace the ZipSleeve* as described above (7.8.1, 7.8.2).

After closing the chamber door with a suitable ZipSleeve* the following screen will be displayed in the graphical user interface:


FIGURE 32. Home screen with loaded chamber A - start parameters; Melt (until ice free state) or Thaw

7.9.3 CHOOSING MELT OR THAW OPTIONS

Choose the shortened Melt cycle, which stops the process for a liquid product **without ice** at the product temperature of 15 °C, or select Thaw with a stop point of 31 °C for a completely thawed product at the transfusion temperature of 30 to 37 °C. Melt option is used to prepare “thawed liquid plasma”, which can be then moved to the refrigerator for further storage at 1-6 °C or carried to the patient's bed for infusion via a bedside blood warmer. The Melt option is often used in STAT situations where the desired operation is intended to produce the maximum number of liquid plasma units for much higher performance. It is useful for creation and maintenance of the working stock of cooled liquid plasma.

7.9.4 Start-up

Press one of two start options MELT or THAW to start the thawing process. The following graphical interface will be displayed:


FIGURE 33. Thawing chamber A screen.

When the thawing cycle starts:

- The graphical interface of the bag will turn **ORANGE**
- Plasma temperature parameters, elapsed thawing time chart are displayed.
- The LED indicator will turn **ORANGE**

Note: Mixing will start automatically as soon as the thawing cycle starts.

7.9.5 INTENDED STOP OF THAWING

In some cases, the thawing cycle may be stopped at a confirmed liquid stage of the substance (i.e. without ice crystals). When the temperature in the chamber reaches + 15 °C/ 59 °F, the thawed substance (i.e. plasma) will be in a liquid state. To interrupt the thawing cycle at any stage, press **STOP** key. The following message screen will appear:


FIGURE 34. Prompt of intended stop of thawing cycle

Press **OK** to confirm the stop of the thawing cycle or press **Cancel** to continue the thawing cycle.

7.9.6 CURRENT THAWING USING THE SECOND CHAMBER B

The ZipThaw™ system allows thawing in both chambers simultaneously but independently of each other. Just press the chamber icon B, the screen will refresh:



FIGURE 35. Screen of chamber A thawing cycle and chamber B readiness screen

Repeat clauses 7.8.1-7.8.4 for chamber B.

7.10 COMPLETION OF THAWING CYCLE

Upon completion of the chamber A thawing cycle the following screen will appear:

31°

- The bag icon will turn **GREEN**.
- A single beep will sound for 1 second.
- The chamber LED turns green

FIGURE 36. Chamber A thawing cycle completed and chamber B loaded and active

Open the chamber door and remove the ZipSleeve® from the chamber, then remove the thawed plasma bag.

Note: If the product remains inside the chamber, the system will start counting the time after the end of the thawing cycle, warn the user (i.e. a sound alarm will appear and the LED indicator will blink) once and again for 1 second after 5 minutes. The operator **must remove** the thawed product from the ZipThaw202 device as soon as the cycle is completed. **Close the chamber door immediately after removal!**

7.11 ENDING WAITING FOR THAWING CYCLE

The ZipThaw™ system has a cycle wait mechanism that stops the thawing cycle in the cases when the thawing cycle does not reach the specified temperature (for example, 31 °C).

Note: The waiting time is currently set to 60 minutes for each thawing cycle. The system will display the following message:


FIGURE 37. Waiting screen

Press **OK**, open the chamber door and remove the ZipSleeve® from the chamber, then remove the thawed plasma bag. *Be sure to close the chamber door immediately!*

DATA EXPORT

At the end of the thawing process or before starting the thawing process, the recorded data can be exported to a USB drive. Press "Export Data" icon, the following graphical interface will appear.


FIGURE 38. Data export screen

After entering the screen, the user scans/enters own user ID (the bar code will light up upon entering the screen and will close after scanning or after the user enters own ID and presses enter). The export data key will be inactive until the identification number is entered.

The user connects then the USB drive and presses **Export Files**. During the downloading process, a pop-up window with the progress bar will appear. *(When choosing a suitable USB drive, make sure that its body is not too wide for the port, which may prevent the USB drive from being fully inserted into the port opening.)*


FIGURE 39: Data load screen

When finished, exit the screen by pressing **EXIT**.

The device will export the copies of all records stored in the device (2000 records can be stored). The user cannot sort the records stored in the device.

The ZIPThaw™ device keeps the records of the last 2000 thawing cycles. When the storage space on the device is full, the oldest record will be deleted when you attempt to create a new record.

8.1 EXPORT OF DATA STRUCTURE

The exported file will have the following fields. The file will be in the format of comma separated values.

No.	Log No. (from 1 to 1999)
Date	Thawing date (MM / DD / YYYY)
start time	Thawing start time
completion time	time (when the cycle is stopped or finished) Format: 24:00
User ID	trawing completion time (when the cycle is stopped or finished) Format: 24:00
Bag ID	Bag ID as specified in the data insertion screen
duration	From the thawing cycle start to end (33 °C/green icon)
extra time after completion	bag kept inside the chamber after the thawing cycle
initial bag temperature	Temperature of the first thawed bag
final bag temperature	Temperature of the last thawed bag
error 01	Errors occurred during the thawing cycle (corresponding chamber)
error 02	
error 03	
error 04	
Bag temp. at 00:30	
Bag temp. at 01:00	
Bag temp. at 01:30	
.....	The temperature is measured twice a minute until the bag is removed.
.... Bag temp. at 60:00	

SHUTDOWN

9.1 STANDBY MODE

To turn off the system, press the (standby)  button located on the front panel (see [Figure 1](#)).

The following prompt will appear on the screen:


FIGURE 40: First screen of device turning off prompt

- ▶ Press CANCEL to cancel and return to the previous screen.
- ▶ To continue shutdown, press OK, the following pop-up message will appear:


FIGURE 41: Device turning off prompt screen

- ▶ Press **CANCEL** to cancel the action and return to the previous screen.
- ▶ Press **YES** to confirm standby mode.

The system will switch to standby mode.

Pressing the standby button will wake up the system and initiate the power up process as described above (see [Section 7.3](#)).

9.2 OFF MODE

To turn off the ZipThaw®202 completely (i.e. to power off), press the power switch on the back panel of the device.

MESSAGE FOR USER

The ZipThaw®202 system will prompt the user with the following operational and error messages:

PHASE	MESSAGE	DESCRIPTION	SOUND ALARM	SYSTEM RESPONSE	ACTION REQUIRED
POWER UP SCREEN	Chamber A Door is Open. Please Close Chamber Door	Warning: Chamber Door A is open.	NO	Pop-up warning message that the chamber door is open.	1. Close the door to resume the power up sequence. 2. The pop-up message will disappear OR press OK.
	Chamber B Door is Open. Please Close Chamber Door	Warning: Chamber Door B is open.	NO	Pop-up warning message that the chamber door is open.	1. Close the door to resume the power up sequence. 2. The pop-up message will disappear OR press OK.
	System Error, FMS Technical Support Required	Error message: The error was found during internal testing of the system.	YES	System error pop-up message. Operational stop.	1. Turn off the power on the back panel of the device. 2. Unplug the power cord for 30 seconds. 3. Turn on the device again.
	Chamber A Error, Use Chamber B	Error message: An error was found during internal testing of the system: Chamber A malfunction.	YES	System error pop-up message. Chamber A will be disabled.	1. Pressing OK will close the pop-up window. 2. Follow the instructions for technical support!
	Chamber B Error, Use Chamber A	Error message: An error was found during internal testing of the system: Chamber B malfunction.	YES	System error pop-up message. Chamber B will be disabled and the graphical interface of chamber will be changed.	1. Pressing OK will close the pop-up window. 2. Follow the instructions for technical support!
SELECT THE CHAMBER SCREEN	Select Chamber	System prompt: Indication of the user's choice to use the available chamber.	NO	No	Follow the operating instructions
PHASE	MESSAGE	DESCRIPTION	SOUND ALARM	SYSTEM RESPONSE	ACTION REQUIRED

DATA LOAD SCREEN	Enter Required Fields	System prompt: The user must enter the required data for the thawing cycle.	NO	System pop-up message.	Pressing OK will close the pop-up window. Follow the operating instructions.
	Incorrect ID, Please Re-Enter ID	System warning: When the wrong ID was entered.	NO	System pop-up message.	Pressing OK will close the pop-up window. Follow the operating instructions
	Data is Incorrect or Missing, Please Recheck Data	System prompt: If incorrect bar codes are entered OR required fields are not completed.	NO	System pop-up message.	Pressing OK will close the pop-up window. Follow the operating instructions
PLASMA BAG INSERT SCREEN	Place Plasma Bag in Chamber A and Press Start	System prompt: Operating instructions	NO	System prompt.	Follow the operating instructions
	Place Plasma Bag in Chamber B and Press Start	System prompt: Operating instructions	NO	System prompt.	Follow the operating instructions
	ZipSleeve Reached Max. Uses, please replace ZipSleeve	System prompt: It indicates that ZipSleeve has reached the maximum number of uses.	NO	1. System pop-up message. 2. Operational stop. 3. An indication to change the user interface (not specified).	1. Pressing OK will close the pop-up window and allow the user to use the new ZipSleeve. 2. Follow the instructions for use and replace ZipSleeve.
	Door A is Open, Please Close Chamber Door	Warning: Chamber A door is closed incorrectly.	NO	1. System pop-up message. 2. The system will disable the thawing operation of chamber A until the chamber door is closed properly.	1. Pressing OK will close the pop-up message. 2. Close chamber A door correctly.
	Door B is Open, Please Close Chamber Door	Warning: Chamber B door is closed incorrectly.	NO	1. System pop-up message. 2. The system will disable the thawing operation in chamber B until the chamber	1. Pressing OK will close the pop-up message. 2. Close chamber B door correctly.

THAWING DURING THE PROCESS				door is closed properly.	
	No ZipSleeve Detected in Chamber A, Please Insert ZipSleeve	Warning: It indicates that the ZipSleeve has not been inserted into chamber A.	NO	1. System pop-up message. 2. Operational pause of chamber A.	1. Pressing OK will close the pop-up message. 2. Follow the operating instructions and insert the ZipSleeve into the intended chamber.
	No ZipSleeve Detected in Chamber B, Please Insert ZipSleeve	Warning: It indicates that the ZipSleeve has not been inserted into the intended chamber B.	NO	1. System pop-up message. 2. Operational pause of chamber B.	1. Pressing OK will close the pop-up message. 2. Follow the operating instructions and insert the ZipSleeve into the intended chamber.
	Chamber A Failure, Technical Support Required	Error message: It indicates that a failure occurred in chamber A after starting the thawing cycle.	YES	1. System pop-up message. 2. The system turns off chamber A after 3 tests. 3. The system allows chamber B to continue operation.	1. Pressing OK will close the pop-up message. 2. Follow the instructions for technical support!
	Chamber B Failure, Technical Support Required	Error message: It indicates that a failure occurred in chamber B after starting the thawing cycle.	YES	1. System pop-up message. 2. The system turns off chamber B after 3 tests. 3. The system allows chamber A to continue operation.	1. Pressing OK will close the pop-up message. 2. Follow the instructions for technical support!
	Chamber Door A Opened Resume Thawing?	Warning: It indicates that the door of chamber A was opened during the thawing cycle.	YES	1. System pop-up message. 2. The system suspends thawing in chamber A.	1. Close the chamber door. 2. Pressing YES will resume thawing. 3. Pressing NO will interrupt the thawing cycle.

	Chamber Door B Opened Resume Thawing?	Warning: It indicates that the door of chamber B was opened during the thawing cycle.	YES	1. System pop-up message. 2. The system suspends thawing in chamber B.	1. Close the chamber door. 2. Pressing YES will resume thawing. 3. Pressing NO will interrupt the thawing cycle.
	Thawing Exceeded Max. Time! Please Remove the Bag From Chamber A/B	Error message: It indicates that the system has interrupted the thawing cycle due to inability to thaw the bag to the specified temperature within the pre-set time.	YES	1. System pop-up message. 3. The system stops the thawing process in the specified chamber.	4. Press OK, open the chamber door and remove the ZipSleeve™ from the chamber, then remove the thawed plasma bag.

POST THAWING CYCLE	Chamber A Door Opened, Thawing Cycle Stopped!	Error message: This means that the door of chamber A has been left open for more than 5 minutes.	YES	1. System pop-up message. The system stops thawing in chamber A.	1. Pressing OK will close the pop-up message. 2. The graphical user interface of chamber A will return to the Select Chamber screen.
	Chamber B Door Opened, Thawing Cycle Stopped!	Error message: This means that the door of Chamber B has been left open for more than 5 minutes.	YES	1. System pop-up message. 2. The system stops thawing in chamber B.	1. Pressing OK will close the pop-up message. 2. The graphical interface of chamber B will return to the Select Chamber screen.
ADMINISTRATOR SCREEN	Password Incorrect	Error message: It indicates that an incorrect password was entered.	YES	1. System pop-up message.	1. Pressing OK will close the pop-up message. 2. Re-enter the correct password.
	Incorrect Data, Please Re-Enter Data	Warning: It indicates that the two entered administrator passwords do not match	NO	1. System pop-up message. 2. Incorrect fields will be marked in red.	1. Pressing OK will close the pop-up message. 2. Re-enter the correct answer.
	Incorrect Data, Please Re-enter Data -OK	Warning message: It indicates that the date and time were filled in incorrectly.	NO	1. System pop-up message. 2. Incorrect fields will be marked in red.	3. Pressing OK will close the pop-up message. 4. Re-enter the correct date and time.

	The Device Will Go to Default Settings All Users Will be Erased. Are You Sure?	Warning message: When pressing the default settings key.	NO	System pop-up message.	- Pressing YES will delete the changed settings and return to the Home screen. - Pressing NO will close the pop-up window and return to the settings/administrat or screen.
	Exit the Screen Without Saving. Are You Sure?	Warning message: When pressing Home after changing the settings, the data will not be saved.	NO	System pop-up message.	- Pressing YES will delete the changed settings and return to the default settings. - Pressing NO will close the pop-up window and return to the screen.

EXPORT DATA	Scan/Enter Your ID, then Insert USB and Export	System prompt: when entering the Export Data screen.	NO	None	
	Incorrect ID, Re-Enter ID -OK	Error message: It indicates that an incorrect ID was entered.	YES	System pop-up message.	- Pressing OK will close the pop-up window. - Re-enter password.
	Downloading...	System prompt: During downloading	NO	The message disappears automatically upon completion of download	None
	No USB Media Identified, Please Insert USB Media	Warning message: It indicates that the user has not connected a USB media to the port.	NO	System pop-up message.	Pressing OK will close the pop-up window.
	Downloading Failure	Warning message: It indicates that the data download has failed.	NO	System pop-up message.	Pressing OK will close the pop-up window.
	Media Memory Full	Warning message: It indicates that the USB key is full	NO	System pop-up message.	Pressing OK will close the pop-up window.
SYSTEM TURNING ON/OFF	(*) Turning OFF ZIPThaw Are You Sure?	Warning message: It indicates that ON/OFF button has been pressed.	NO	System pop-up message.	Pressing NO will remove the pop-up window and return to the previous system status.

					2. Pressing YES will turn off the system.
	(**) System Shutdown! Thawing Cycle Will be Stopped	Warning message: It indicates that ON/OFF button was pressed while the thawing cycle was in progress.	YES	1. System pop-up message (**). 2. The thawing cycle is paused. 3. IF the user selects YES: a message will appear (*).	1. Pressing NO will remove the pop-up window and resume thawing. 2. When pressing YES, a pop-up message will appear (*). 3. Pressing NO will remove the pop-up window and return to the previous system status. 4. Pressing YES will turn off the system.

CLEANING

11.1 CLEANING FREQUENCY

Cleaning should be performed on daily basis and as needed if there is visible contamination on the surfaces or if you can see any liquid in the overflow tray.

11.2 CLEANING AND DISINFECTION OF EXTERNAL SURFACE OUTSIDE CLEANING AND DISINFECTION

To clean the ZipThaw®202 system, use disinfectant alcohol wipes or spray that both cleans and disinfects as well as has effect on a wide range of pathogens including bacteria, yeast/mold, and viruses.

11.3 SCREEN CLEANING

To clean the touch screen of the product, use lens wipes containing water, isopropyl alcohol, detergent, and antistatic agent.

11.4 CLEANING OF SPILL TRAY

The spill trays should be emptied, cleaned and disinfected as needed, for example, if a ZipSleeve® plasma bag is leaking inside or if blue liquid is leaking from the heating cushion. If blue liquid is visible outside the heating cushion, do not use the chamber with a leaking cushion and call to service division.

The spill trays can be accessed through the back panel of the ZipThaw®202 by opening the drain cover and removing the disposable trays.

Note:

Liquid inside the spill tray bags indicates a leak in the heating cushion that requires technical support.

ZipSleeve® bags can be used up to 8 times. If moisture appears on them between uses, it can be wiped off with a dry napkin, if blood has leaked on the bag it must be replaced with a new one.

WARRANTY

Three-year full replacement warranty. See the Warranty document enclosed with each device.

CUSTOMER SUPPORT SERVICE

For maintenance and repair:

Contact your authorized representative in the Russian Federation, email us at shreyayar@yahoo.com or visit our website <https://fremonscientific.com> for assistance.

Follow the instructions for return specified on the Return Merchandise Authorization (RMA) form issued by the Product Management Department. The device must be cleaned and disinfected before being sent for maintenance or disposal.

Medical device operating conditions

from 10 °C to 32°C (from 50°F to 86°F)

from 20% to 70% without condensation

SPECIFICATIONS

Thawing technology:	Electric power, disposable bag
Typical capacity of thawed bag:	ZipSleeve® bags
Dimensions, cm (ZipSleeve® bags)	15x29 ± 1 (with packaging 19x32.5 ± 1)
Temperature of heating cushions:	
Standby mode:	36.5°C±0.5°C
Thawing:	36.5°C±0.5°C
ZipSleeve® sensor shutdown temperature:	33°C
Average thawing time (minutes), single or double bag	100 ml - ≤ 10 250 ml - ≤ 16 450 ml - ≤ 40
Mixing method:	Electronically controlled mechanical movement
Mixing parameters:	Cycle: 2.2 c/sec. Amplitude: ± 32mm
Visual display/User interface:	7" Touch screen with graphical user interface
System:	Main CPU: VAR-SOM-SOLO/Dual: Freescale i.MX6 • High performance 1.0 GHz single and dual core Cortex-A9™ processor • 1GB DDR3, 1GB SLC NAND and eMMC 64GB • Full HD 1080p with video encoding/decoding option • Vivante CPU providing 2D/3D acceleration Chamber control: STMicroelectronics, STM32L073, MCU 32BIT 192KB
Interface and network:	USB • USB 2.0 OTG • USB 2.0 Host
Memory bank:	Internal for 2000 records of thawing cycle
Time required to reach full operation	5-10 sec. • Consumption: 100-240 V AC, maximum rated current 9A Fuse 15A/250V
Heating power:	Max. 400W, each chamber
Electrical characteristics (V/Hz/A):	100-240 VAC, 50/60Hz 9Amps
Dimensions: [WxDxH]	cm: 47.1 x 27.6 x 41.6 ± 1 (including top handle) inches: 18.5 x 10.9 x 16.4 ± 1 (including top handle)
Weight: [Lb/kg]	13.6 ± 1kg/16.85 ± 1 (with packaging)
Portable:	Yes
Cushions (inside ZipThaw®202 blood plasma thawer)	The cushions are made of plastic, the cushions are filled with deionized water with blue dye.

The level of radio interference and acoustic noise during operation must not exceed 50 dB.

The external surfaces of the device must not have damaged protective decorative coatings. The surface of the devices must not have cracks, bubbles, chips, cavities and other defects that worsen their marketable condition and technical characteristics. The metal and non-metal parts of the device must be made of the material, which is resistant to corrosion or have reliable protection against it.

The device must provide continuous operation during the entire service life of the device.

The device must be mechanically strong and remain functional during and (or) after exposure to mechanical loads:

Description of	Frequency	Accelerati	Impact duration, ms
----------------	-----------	------------	---------------------

mechanical load	band, Hz	on, m/s ²	
Vibration	1-60	29.4	-
Multiple impacts	-	78.4	5-10
Impact stress during transportation	-	147	10-15

The device must remain functional during and (or) after exposure to climatic factors:

- air temperature from +10 °C to + 32 °C;
- relative humidity from 20 to 70%.

Average service life is 10 years. Power cord length 2.5 m (± 10%).

Ingress protection: housing IP52

The device must comply with GOST R IEC 62366. For electromagnetic compatibility, the device must comply with GOST R IEC 61326-1.

Audio-visual alarm: in case of system errors, an error audio alarm will appear for 1 second and a system error message will be displayed on the touch screen. (see USER MESSAGES).

Software version V01.31 dated 11.09.2020 and class A according to GOST R IEC 62304.

Service life data

Service life - 10 years.

Transportation and storage requirements

The devices are transported by all types of covered vehicles except for unheated aircraft compartments and as per the rules approved in accordance with the established procedure for each transport mode.

Temperature during storage and transportation:

(in original packaging): from -20 °C to 70 °C (from 4 °F to 158 °F)

Humidity during storage (non-condensing): 30% - 90% without condensation

Requirements for safe destruction and disposal

Disposal of the instrument must be carried out in accordance with applicable legislation. The used instruments are referred to class A waste and shall be disposed in accordance with SanPiN 2.1.3684-21. The electronic and electrical components of the device must be disposed separately from household waste. The used bags are referred to class B waste and shall be disposed in accordance with SanPiN 2.1.3684-21. Disposal of the instrument after the end of its service life shall include: disassembly, cleaning from residues of the working medium and contamination, drying and sorting of materials.

Maintenance

Following the maintenance procedures is recommended for safe operation. Regular maintenance is required for safe operation of the instrument. Regular maintenance is carried out by qualified engineers according to the schedule indicated in the service manual. Description and frequency of regular maintenance is given in the user manual supplied with equipment. The frequency of inspections is indicated on the basis of the assumption that the equipment is in use for 12 hours a day.

Maintenance plan:

Check/calibration	Period
Appearance of the instrument	Condition based

Check for physical damage to the instrument	
Checking for physical damage to LCD screen	Condition based
Checking of the operation interface	Condition based
Checking of the power cord	Once a month
Checking of the fuse	Once a year
Checking of the cushion	Before every use replacement once a year

PROTECTION AGAINST HARMFUL INTERFERENCE

15.1 EMC SAFETY

- The device was tested and found to comply with the limits for medical devices as per standard GOST R IEC 61326-1. These limits are designed to provide the required protection against harmful interference at default installation in a clinical setting. This device generates, uses and can radiate radio frequency energy. Installation and use with breach of the instructions may cause harmful interference to other devices nearby. However, there is no guarantee that interference with other devices, which can be determined by turning the device on and off, is caused by this tool.

The user is recommended to attempt to remove the interference using one or more of the following methods:

- Reorient or relocate the receiving device.
- Increase the distance between devices.
- Connect the device to a power that has not been used before.
- Contact FreMon Scientific service personnel for assistance.

Interference with the device can be caused by portable and mobile RF communications equipment. In case of interruption, avoid close proximity of such a device.

- Using the system with any accessories, sensors, or cables other than those specified may result in increase in RADIATION or decrease in INTERFERENCE IMMUNITY compared to specified ones.
- The system must not be used next to other equipment or placed on top of other equipment. If it is required to be located next to or on top of other equipment, the system should be monitored to ensure that it is operating normally in the configuration it will be used in.

Recommendations and declaration of the manufacturer - electromagnetic emission

The ZipThaw®202 system is intended for use in the electromagnetic environment specified below. The buyer or user of ZIPThaw® 202 must make sure that it is used in such an environment.

Emission test	Conformity	Electromagnetic environment - recommendations
RF emission CISPR 11	Group 1	The ZIPThaw® 202 system uses RF energy only for its internal functions. Therefore, its RF emission is very low and unlikely to cause any interference with nearby electronic equipment. The ZIPThaw® 202 system is suitable for use at all institutions other than home use and can be used in residential premises and those that are directly connected to the public low-voltage grid supplying buildings used for domestic purposes provided for compliance with the following warning: Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or disrupt operation of nearby equipment. Mitigation measures such as reorientation or relocation of ZIPThaw® 202 or shielding of the location may be required.
RF emission CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	In conformity	

Recommendations and declaration of the manufacturer - electromagnetic interference immunity

The ZIPThaw® 202 system is intended for use in the electromagnetic environment specified below. The buyer or user of ZIPThaw® 202 must make sure that it is used in such an environment.

Immunity tests	Control level as per IEC 61326-1	Conformity level	Electromagnetic environment - recommendations
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	The floor must be wooden, concrete or ceramic. If the floor is covered with synthetic material, the relative humidity must be at least 30%
The electrical fast transient processes/impulse interference as per IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 2 kV for input/output lines	The quality of the mains must comply with the standard for commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV in differential mode ± 2 kV in common mode	± 1 kV in differential mode ± 2 kV in common mode	The quality of the mains must comply with the standard for commercial or hospital environment.

Voltage dips, short interruptions and voltage fluctuations on input power supply lines according to IEC 61000-4-11	$<5\% U_T$ (dip $> 95\% U_T$) for 0.5 cycle $40\% U_T$ (dip $60\% U_T$) for 5 cycles $70\% U_T$ (dip $30\% U_T$) for 25 cycles $<5\% U_T$ (dip $> 95\% U_T$) for 5 sec.	Dip $> 95\%$ for 10 ms Dip 60% for 100 ms Dip 30% for 500 ms Dip 95% for 5000 ms	The quality of the mains must comply with the standard for commercial or hospital environment. If the user of ZIPThaw® 202 requires continuous operation during a power outage, it is recommended to connect ZIPThaw® 202 to uninterruptible power supply.
Power frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	The power frequency magnetic fields must be at the levels peculiar to a typical location in a typical commercial or hospital environment.

NOTE. - U_T is the AC mains voltage before the test level is applied.

Recommendations and declaration of the manufacturer - electromagnetic interference immunity

The ZIPThaw® 202 system is intended for use in the electromagnetic environment specified below. The buyer or user of ZIPThaw® 202 must make sure that it is used in such an environment.

IMMUNITY TESTS	Control level According to IEC 61326-1	Conformity level	Electromagnetic environment - recommendations
RF test performed as per IEC 61000-4-6 Level of RF emission as per IEC 61000-4-3	3 V rms From 150 kHz to 80 MHz 3 V/m From 80 MHz to 2.7 GHz	[3] V [3] V/m	The portable and mobile RF telecommunications equipment should not be used closer to any part of ZIPThaw®202, including cables, than the recommended distance calculated from the equation applicable to the frequency of the transmitter. Recommended safety distance From 800 MHz to 2.5 GHz, where P is the maximum power output of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended distance in meters (m). The field density from fixed radio transmitters, as

			determined by the electromagnetic site survey,^a must be less than the conformity level in each frequency band.^b Interference may occur near equipment marked with the following symbol: 
--	--	--	--

NOTE 1. At 80 MHz and 800 MHz, the higher frequency band applies.

NOTE 2: These recommendations may apply not in all situations. Propagation of electromagnetic waves is affected by absorption and reflection from structures, objects and people.

^a The field density from fixed transmitters such as base stations for radio telephones (cellular/wireless) and land mobile stations, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted with high accuracy in theoretical terms. To evaluate the electromagnetic environment generated by fixed RF transmitters, the electromagnetic site survey should be considered. If the measured field density at the location where ZIPThaw™ 202 is used exceeds the above applicable RF compliance level, the [MEDICAL EQUIPMENT or MEDICAL SYSTEM] must be checked to verify normal operation. If malfunctions are observed in operation, the additional measures may be required, such as reorientation or relocation of ZIPThaw™ 202.

^b The field density strengths in the frequency band from 150 kHz to 80 MHz must be less than [3] V/m.

Recommended safety distances between portable and mobile RF communications equipment and ZipThaw 202			
ZIPThaw® 202 is intended for use in electromagnetic environment in which emitted RF interference is controlled. The buyer or user of ZIPThaw® 202 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and ZIPThaw® 202 as recommended below in accordance with the maximum output power of the communications equipment.			
Rated maximum output power of the transmitter [W]	Safety distance depending on the frequency of the transmitter [m]		
	From 150 kHz to 80 MHz	From 80 MHz to 800 MHz	From 800 MHz to 2.5 GHz
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33
For the transmitters rated for maximum power output not indicated above, the recommended safe distance d in meters (m) can be calculated using the equation applicable to the frequency of the transmitter, where P is the maximum power output of the transmitter in watts (W) according to the specifications of the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency band applies.			
NOTE 2: These recommendations may apply not in all situations. Propagation of electromagnetic waves is affected by absorption and reflection from structures, objects and people.			

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Copyright © FreMon Scientific, Inc.