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SECRETARY OF STATE
STATE OF INDIANA

Holli Sullivan
Secretary of State


APOSTILLE

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1. Country: United States of America
2. This public document has been signed by *Elissa Robinson*
3. Acting in the capacity of notary public in & for *Monroe County*
4. And bears the seal/stamp of notary public in & for the State of Indiana

CERTIFIED

5. at Indianapolis, Indiana
6. This *Second* day of *May*, 2022
7. by the Secretary of State of Indiana
8. No. A2022-0502022705
9. Seal/Stamp:
10. Signature:

Holli Sullivan
Indiana Secretary of State

This document was revised April 12, 2021.

This apostille is NOT for use anywhere within the United States of America.


INSTRUCTION FOR USE ENDOVASCULAR WIRE GUIDES

VASCULAR MICRO WIRE GUIDE

INTENDED USE

Endovascular Wire Guides are used for percutaneous endovascular access.

1

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The TPMG micro wire guide has a stainless steel shaft and a soft distal spring coil portion with a radiopaque platinum tip.

NOTE: Not all TPMG micro wire guides are hydrophilic coated.

COATING

All TPMG micro wire guides with product codes ending on -HC are hydrophilic coated with AQ®.

AQ is a bio-compatible hydrophilic coating which becomes very lubricious when wet, greatly reducing the coefficient of friction of the surface of the micro wire guide. AQ hydrophilic coating is applied to the distal end of the micro wire guide. This facilitates control of the proximal portion of the micro wire guide during procedural use. Refer to label for further information.

INTENDED USE

The micro wire guide is intended for diagnostic and interventional procedures.

CONTRAINDICATIONS

None known.

WARNINGS

- Possible allergic reactions should always be considered.
- Insert micro wire guide from the flexible end. (Inserting from the rigid end may cause damage to tissue or device).
- Manipulation of products requires fluoroscopic control.

PRECAUTIONS

- Manipulation of micro wire guides always requires fluoroscopic visualization.
- End hole size and length of the micro catheter must be taken into consideration to ensure proper fit between micro wire guide and micro catheter.
- Altering tip configuration of the micro wire guide may result in damage.
- Withdrawal or manipulation of micro wire guide through a needle tip may result in damage or breakage.
- The products are intended for use by physicians trained and experienced in diagnostic and interventional techniques.
- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Never use an AQ hydrophilic coated micro wire guide through a metal cannula. Destruction of the coating may result. A plastic outer sheath needle is recommended if the micro wire guide is used for initial placement.
- Always keep the surface of the AQ hydrophilic coated micro wire guide wetted for optimum results.
- When an AQ hydrophilic coated micro wire guide has been removed from a micro catheter, wipe it with a sterile gauze pad moistened with saline solution or heparinized saline solution and place it in the holder again. Fill the holder with saline solution. This will maintain the lubricious state of its surface for subsequent procedural use.
- The AQ hydrophilic coated micro wire guides are very slippery when wet. Always maintain control of the micro wire guide when manipulating it through the catheter.

POTENTIAL ADVERSE EVENTS

Potential adverse events are including:

- Allergic reaction
- Arrhythmia or cardiac arrest
- Cardiac perforation
- Death

- Embolism
- Foreign body/wire fracture
- Hemorrhage or hematoma
- Infection or sepsis
- Organ injury
- Pseudoaneurysms
- Stroke or other neurological event
- Vascular access site complications
- Vascular thrombus/thromboembolism
- Vessel occlusion
- Vessel spasm
- Vessel trauma, dissection, perforation, rupture or injury
- Wire entrapment/entanglement
- Wound complications

Instructions for activating hydrophilic coating

1. Attach the syringe with heparinized saline solution or sterile water to the Luer lock fitting of the micro wire guide holder.

2. Inject enough solution to wet the micro wire guide surface entirely (Fig. 1). This will activate the AQ hydrophilic coating, making the micro wire guide surface very lubricious.

3. Standard micro wire guide use may now be initiated. An initial flush of heparinized saline through the catheter lumen will enhance initial movement of the micro wire guide through the catheter.

NOTE: If the surface of the micro wire guide becomes dry after removal from the holder, wetting with additional heparinized saline or sterile water will renew the hydrophilic effect.

INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the micro wire guide from its outer packaging and place in the sterile field.

2. Rinse or wipe the micro wire guide with heparinized saline or sterile water prior to use.

3. Carefully insert the distal tip of the micro wire guide into the micro catheter.

4. Under fluoroscopy and while maintaining position of the micro catheter advance the micro wire guide to the targeted site.

NOTE: Under fluoroscopy, observe all micro wire guide movement in the vessel. The micro wire guide should be advanced only when visualizing the tip of the micro wire guide fluoroscopically.

NOTE: If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take action necessary to relieve the resistance. Advancement and withdrawal of the micro wire guide should be performed slowly and carefully.

5. With the micro wire guide secured in place, advance the micro catheter or other therapeutic device to the target site.

NOTE: Always maintain a continuous flush while exchanging micro wire guides to prevent air entry into the system.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred. Do not use the device if damage has occurred or if the sterilization barrier has been damaged or broken; instead, return the product to Cook.

REFERENCES


These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

Microcoils™-18 Coil Pusher

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Microcoils™-18 Coil Pusher is a .018 inch (0,46 cm) diameter PTFE coated stainless steel mandril shaft with a 49 cm very soft distal spring coil tip, a 1 cm radiopaque marker is placed 1 cm from the distal tip.

The product also contains a wire guide inserter, and a pin vise handle for manipulation (Fig 1). Refer to label for further information.

INTENDED USE

The Microcoils™-18 Coil Pusher is intended for placement of .018 inch (0,46 cm) microcoils through a microcatheter with an ID of .018 to .025 inches (0,46-0,64 cm)

CONTRAINDICATIONS

None known.

WARNINGS

- Possible allergic reactions should be considered.
- Manipulation of products requires fluoroscopic control.
- To determine correct catheter position, perform an angiogram prior to embolization.
- Never advance the coil pusher outside the tip of the microcatheter after the embolization coil has been placed. This might displace the coil and damage the vessel.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.
- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Prior to introduction of the coil, flush the microcatheter with saline.

INSTRUCTIONS FOR USE

1. For coil loading instructions, please see IFU supplied with the coil.
2. Introduce the coil pusher through the wire guide inserter into the catheter.
3. Slowly advance the coil pusher under fluoroscopic control until the coil has reached the target lesion at the distal tip of the catheter. The coil must be advanced completely

outside the microcatheter.

NOTE: Never advance the coil pusher outside the microcatheter tip.

4. To facilitate coil pusher control the radiopaque marker (1 cm from the tip) should remain inside the microcatheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. The coil pusher may be used to deliver multiple coils in one procedure. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

AMPLATZ ULTRA STIFF WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Amplatz Ultra Stiff Wire Guide is a PTFE-coated stainless steel wire guide available in a straight or J-tip configuration. The suffix –AUS1 indicates a 3 cm flexible tip and –AUS2 a 7 cm flexible tip.

INTENDED USE

The wire guide is intended for diagnostic and interventional procedures.

CONTRAINDICATIONS

None known

WARNINGS

- Possible allergic reactions should always be considered.
- Wire guide stiffness should be considered.
- Insert wire guide from the flexible end. (Inserting from the rigid end may cause damage to tissue or device).
- If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take appropriate steps to relieve the resistance. Advancement and withdrawal of the wire guide should be performed slowly and carefully.

PRECAUTIONS

- Manipulation of wire guides requires fluoroscopic visualization.
- End hole size and length of the catheter must be taken into consideration to ensure proper fit between wire guide and catheter.
- Altering tip configuration of the wire guide may result in damage.
- Withdrawal or manipulation of wire guide through a needle tip may result in breakage.
- The products are intended for use by physicians trained and experienced in diagnostic and interventional techniques.
- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

POTENTIAL ADVERSE EVENTS

Potential adverse events are including:

- Allergic reaction
- Arrhythmia or cardiac arrest
- Cardiac perforation
- Death
- Embolism
- Foreign body/wire fracture
- Hemorrhage or hematoma
- Infection or sepsis
- Organ injury
- Pseudoaneurysms
- Stroke or other neurological event
- Vascular access site complications
- Vascular thrombus/thromboembolism
- Vessel occlusion
- Vessel spasm
- Vessel trauma, dissection, perforation, rupture or injury
- Wire entrapment/entanglement
- Wound complications

INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the wire guide from its outer packaging and place in the sterile field.
2. Rinse or wipe the wire guide with heparinized saline or sterile water prior to use.
3. Carefully insert the distal tip of the wire guide.
4. Under fluoroscopy and while maintaining position of the peripheral vascular access catheter advance the wire guide to the targeted site.

NOTE: Under fluoroscopy, observe all wire guide movement in the vessel. The wire guide should be advanced only when visualizing the tip of the wire guide fluoroscopically.

NOTE: If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take action necessary to relieve the resistance. Advancement and withdrawal of the wire guide should be performed slowly and carefully.

5. With the wire guide secured in place, advance the percutaneous catheter or other therapeutic device to the target site.

NOTE: Always maintain a continuous flush while exchanging wire guides to prevent air entry into the system.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

Do not use the device if damage has occurred or if the sterilization barrier has been damaged or broken; instead, return the product to Cook.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

LUNDERQUIST Extra Stiff WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Lunderquist Extra Stiff Wire Guide (TS(C)MG-/-LES) is a PTFE -coated stainless steel wire guide. Lengths of ≥ 260 cm have 4 cm of tip flexibility, including an inner gold coil for enhanced visibility. Lengths of < 260 cm the distal tip have 4 or 7 cm of tip flexibility, and the TS(C)MG-/LES1 has a 3, 5, or 7 cm flexible tip.

The Lunderquist Extra Stiff DC Wire Guide (TSCMG-/(E)-LESDC) is a PTFE -coated stainless steel wire guide with a double curved tip design. The double curved distal tip has 4 cm of tip flexibility, including an inner gold coil for enhanced visibility. LES-DC has a primary/secondary curve radius of 75/15 mm, and E-LESDC has a primary/secondary curve radius of 55/15 mm.

The wire guide is used both to assist in anatomical access for other devices (not included), and to support the delivery of medical devices. The wire guide is introduced into the target vessel; other devices, such as a sheath, catheter, stent, or endovascular graft can then be passed over the wire guide to be positioned or manipulated within the vascular system.

INTENDED USE

The device is intended for complex diagnostic and interventional procedures where increased body, flexibility, and low surface friction of the wire guide is needed.

CONTRAINDICATIONS

None known.

WARNINGS

- Possible allergic reactions should always be considered.
- Wire guide stiffness should be considered.
- Insert wire guide from the flexible end. (Inserting from the rigid end may cause damage to tissue or device).
- If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take appropriate steps to relieve the resistance to avoid the risk of vessel perforation. Advancement and withdrawal of the wire guide should be performed slowly and carefully.

PRECAUTIONS

- Manipulation of wire guides always requires fluoroscopic visualization.
- End hole size and length of the catheter must be taken into consideration to ensure proper fit between the wire guide and catheter.
- Altering the tip configuration of the wire guide may result in damage.
- Withdrawal or manipulation of the wire guide through a needle tip may result in damage or breakage.
- The products are intended for use by physicians trained and experienced in diagnostic and interventional techniques.
- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

POTENTIAL ADVERSE EVENTS

Potential adverse events are including:

- Allergic reaction
- Arrhythmia or cardiac arrest
- Cardiac perforation

- Death
- Embolism
- Foreign body/wire fracture
- Hemorrhage or hematoma
- Infection or sepsis
- Organ injury
- Pseudoaneurysms
- Stroke or other neurological event
- Vascular access site complications
- Vascular thrombus/thromboembolism
- Vessel occlusion
- Vessel spasm
- Vessel trauma, dissection, perforation, rupture or injury
- Wire entrapment/entanglement
- Wound complications

INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the wire from its outer packaging and place it in the sterile field.
2. Attach a syringe with heparinized saline solution or sterile water to the Luer Lock fitting of the wire guide holder.
3. Inject enough solution to fill the holder and bathe the wire in solution.
4. Carefully insert the distal tip of the wire guide in the peripheral vascular access catheter.
5. Under fluoroscopy, maintain the position of the peripheral vascular access catheter and advance the wire guide to the targeted site.

NOTE: Under fluoroscopy, observe all wire guide movement in the vessel. Do not advance the wire guide if you cannot visualize its tip fluoroscopically.

NOTE: If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take action necessary to relieve the resistance to avoid the risk of perforation. Be sure to advance and withdraw the wire guide slowly and carefully.

6. With the wire guide secured in place, advance the therapeutic device to the target site.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred. Do not use the device if damage has occurred or if the sterilization barrier has been damaged or broken; instead, return the product to Cook.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

APPROACH® HYDRO ST MICRO WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Approach Hydro ST Micro Wire Guides are 0.014 inches (0.36 mm) in diameter and are available in a variety of lengths and configurations

Approach Hydro ST incorporates a stainless steel shaft with distal stainless steel and platinum coils and a hydrophilic floppy tip, used for general peripheral vascular interventions.

INTENDED USE

Approach Hydro ST Wire Guides are intended for use in facilitating delivery of percutaneous catheters into the peripheral vasculature.

CONTRAINDICATIONS

None known

WARNINGS

- This wire guide is a delicate instrument and should be handled carefully; forceful angulation should be avoided.
- The wire guide should be advanced only when visualizing the tip of the wire guide fluoroscopically. The wire guide should not be torqued without evidence of corresponding movement of the distal tip. **NOTE:** Torquing against resistance may cause vessel trauma or wire guide damage which may lead to device fracture.
- Refer to the labeling of the percutaneous catheter or other therapeutic devices for contraindications and potential complications associated with use.
- Operate these wire guides carefully to avoid injury to the blood vessel, observing information in these instructions.
- Approach Hydro ST wire guides should never be used in conjunction with or through a metal needle cannula.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Avoid manipulating or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape the coating or shear the wire guide. A catheter, introducer sheath or vessel dilator should replace the needle as soon as the wire guide has been inserted in the vessel.
- Excessive tightening of the torque device (if provided) onto the wire guide may result in abrasion of the coating on the wire.
- Always flush the Hydro ST before use.
- Always maintain a wet surface on the hydrophilic portion of the Hydro ST for optimal results.

INSTRUCTIONS FOR USE

1. Flush the wire guide holder by attaching or pressing a syringe with heparinized saline or sterile water to the fitting of the wire guide holder and injecting enough solution to wet the wire guide surface entirely.
2. Carefully remove the wire guide from the holder.
3. If desired, wire guide tip may be carefully shaped using standard tip shaping techniques. Do not use a shaping instrument with a sharp edge.

Page 9 of 53

4. Insert a wire guide insertion tool through the valve assembly of the guiding sheath or guiding catheter.
5. Carefully insert the distal tip of the wire guide through the insertion tool.
6. Attach the torque device (if provided) to the wire guide.
7. Under fluoroscopy and while maintaining position of the peripheral vascular access catheter, advance the wire guide to the targeted site.

NOTE: Under fluoroscopy, observe all wire guide movement in the vessels. The wire guide should not be torqued without evidence of corresponding movement of the distal tip. Further torquing against resistance may cause vessel trauma or wire guide damage which may lead to device fracture. The wire guide should be advanced only when visualizing the tip of the wire guide fluoroscopically.

NOTE: If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take action necessary to relieve the resistance. Advancement and withdrawal of the wire guide should be performed slowly and carefully.

8. With the wire guide secured in place, advance the percutaneous catheter or other therapeutic device to the target site.

NOTE: Always maintain a continuous flush while exchanging wire guides to prevent air entry into the system.

NOTE: Prior to re-insertion of the exchanged wire guide, verify that the distal tip of the cardiovascular device is free within the vessel. If the tip is against the wall, vessel trauma may occur upon re-introduction of the wire guide.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

Disposable Reuter Tip-Deflecting Wire Guides

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Disposable Reuter Tip-Deflecting Wire Guide is a handle/wireguide assembly that allows for controlled deflection of the wire guide tip within a catheter lumen. The handle is constructed of plastic, with a plastic side-arm adapter. The wire guides are constructed of stainless steel, with an optional heparin coating, and have a tip-deflected curve radius of 5 or 10 mm.

INTENDED USE

Disposable Reuter Tip-Deflecting Wire Guides are intended for use in curving or deflecting catheter tips for pulmonary angiography, selective angiography, translumbar aortography, bronchography, repositioning of central venous catheter tips and other vascular and non-vascular applications. The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

None known

WARNINGS

- Exercise caution when using a Disposable Reuter Tip-Deflecting Wire Guide. This wire guide does not have a soft tip; therefore, to avoid damaging the vessel, **do not advance the distal tip of the wire guide beyond the distal tip of the catheter.**
- If the side-arm fitting is not used, exercise extra caution to ensure that wire guide tip does not extend beyond catheter tip.
- Do not use excessive force when deflecting a Disposable Reuter Tip-Deflecting Wire Guide. Keep fingers and thumb in rings while deflecting wire guide. Throw length is designed to help limit excessive use of force in order to prevent damage to handle and wire guide.
- Do not attempt deflection of wire guide until it is inserted into a catheter. Such deflection may result in damage to tip of wire.
- When inserted in catheter, wire guide tip should be near the endhole, but should not be lodged in the catheter tip orifice. If the wire guide tip is lodged in the tip orifice of the catheter, flushing of the catheter tip is ineffective.
- Table-top test the handle, wire guide and catheter assembly prior to use. **Confirm that the wire guide does not protrude through or obstruct the catheter endhole. Ensure that the wire guide, when deflected, cannot exit through a catheter sideport.**
- Always monitor manipulation using fluoroscopic guidance.

PRECAUTIONS

- Always use a Disposable Reuter Tip-Deflecting Wire Guide with a cross-sectional diameter at least one size larger than the endhole diameter of the catheter with which it is paired (e.g., use a .038 inch (0.97 cm) wire guide with a .035 inch (0.89 cm) endhole catheter). This will help prevent wire guide tip from inadvertently exiting catheter tip during manipulation.
- A catheter with sideports may be used with the Disposable Reuter Tip-Deflecting Wire Guide. However, care should be taken to ensure that the wire guide cannot accidentally exit through a catheter sideport.

INSTRUCTIONS FOR USE

NOTE: Throughout duration of procedure, periodically use side-arm fitting to flush system with heparinized saline solution to prevent clot formation on catheter or deflecting wire guide surface.

1. After catheter is properly positioned in vessel, perform a wire guide exchange for the Disposable Reuter Tip-Deflecting Wire Guide/Handle assembly. Attach the male Luer lock on the sidearm fitting to the catheter hub.
2. To deflect the wire guide and catheter, insert thumb and fingers in proper rings and slowly but firmly pull fingers back to deflect tip.

WARNING: Exercise caution when using a Reuter Tip-Deflecting Wire Guide. The distal tip of the wire guide should not be advanced beyond the distal tip of the catheter. Damage to the vessel wall can occur.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

ROADRUNNER® UNIGLIDE™ HYDROPHILIC WIRE GUIDE

CAUTION: Federal law (USA) restricts this device to the sale and use by or on the order of a physician.

DEVICE DESCRIPTION

Utilizing proprietary processes, the Roadrunner UniGlide Hydrophilic Wire Guide is constructed from a steerable, metallic core with a polymer coating. A hydrophilic coating is applied over the radiopaque polymer jacket. Wire guide Length, Diameter, Coating, Tip, and Core configurations are indicated on the product label.

Standard Shaft Configuration

Stiff Shaft Configuration

Angled Tip Configuration

Straight Tip Configuration

INTENDED USE

The Roadrunner UniGlide Hydrophilic Wire Guide is intended for use in facilitating the delivery of percutaneous catheters into the peripheral vasculature.

CONTRAINDICATIONS

None known

WARNINGS

- The hydrophilic wire guide may slide entirely into the catheter, sheath introducer, vessel dilator or other device because of its low sliding friction. To prevent this, keep at least 5 cm of the wire protruding from the device fitting at all times.
- To prevent possible tissue damage, care should be taken when manipulating a device over a wire guide during the device's placement and withdrawal. If resistance is felt during device placement, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of resistance cannot be determined, remove the wire guide and device as a unit to prevent possible damage and/or complications.
- When using a wire guide, potential exists for thrombus formation or emboli, arterial or venous wall damage and/or plaque dislodgment. The physician should be familiar with the literature concerning the complications of angiography.
- Cook does not recommend a particular technique for the use of this wire guide. The steps contained in the following directions are for information purposes only. Each physician should evaluate their appropriateness according to individual patient condition and his or her medical training and experience.


PRECAUTIONS

- Do not resterilize and/or reuse this device.
- Carefully read all instructions prior to use. Observe all warnings and cautions. Failure to do so may result in complications.
- Avoid manipulating or withdrawing the hydrophilic wire guide back through a metal needle or cannula. A sharp edge may scrape the coating or shear the wire guide. A catheter, introducer sheath or vessel dilator should replace the needle as soon as the wire guide has been inserted in the vessel.
- It is recommended that a plastic torque device be used to handle the hydrophilic wire guide. Use of a metal torque device may damage the wire.

INSTRUCTIONS FOR USE

Preparation

1. Before removing the hydrophilic wire guide from its dispenser, inject sterile heparinized saline

solution into the Luer lock hub end of the dispenser.

2. Inject enough solution to fill the dispenser coil. This will completely cover the wire guide surface and activate the hydrophilic coating.
3. Remove the hydrophilic wire guide from its dispenser by gently withdrawing the wire's tip.
4. If the hydrophilic wire guide cannot easily be removed from its dispenser, inject more heparinized saline solution into the dispenser and then try again.

Wire Guide Use

1. Fill catheter or other device with heparinized saline solution before and during use to ensure smooth movement of the hydrophilic wire guide within the device.
2. When wet with saline solution or blood, the hydrophilic wire guide will become lubricious. Use of sterilized gauze moistened with heparinized saline solution facilitates handling the wire.
3. Insert the wire guide into the device and advance to desired position.

NOTE: If movement of the wire within the device becomes diminished, remove the wire guide and reactivate the hydrophilic coating by wetting its entire surface with heparinized saline solution.

After Use

1. Excess blood can be removed from the wire guide by wiping the surface once or twice with gauze moistened with heparinized saline solution. **DO NOT** use dry gauze as this may damage the wire surface resulting in increased resistance when the wire is reinserted into the device.
2. Use of alcohol, antiseptic solutions or other solvents must be avoided, because they may adversely affect the surface of the hydrophilic wire guide.

Standard Shaft Configuration

Stiff Shaft Configuration

Angled Tip Configuration

Straight Tip Configuration

3. After cleaning the wire, replace it by its proximal end into its dispenser filled with heparinized saline solution. This wire guide can only be used during the same procedure on the same patient.
4. Upon reinsertion into the device, if the wire guide does not move with initial ease, exchange it for a new UniGlide™ Hydrophilic Wire Guide.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

HIWIRE™ HYDROPHILIC WIRE GUIDE

Federal law (USA) restricts this device to the sale and use by or on the order of a physician.

PACKAGE CONTENTS:

1 - HiWire™ Hydrophilic Wire Guide

NOTE: Wire guide Length, Diameter, Coating, Tip, and Core configurations are indicated on the product label.

DESCRIPTION:

Utilizing proprietary processes, the HiWire™ Hydrophilic Wire Guide is constructed from a steerable, metallic core with a polymer coating. A hydrophilic coating is applied over the radiopaque polymer jacket.

INTENDED USE:

To facilitate the placement of devices during diagnostic and interventional procedures.

STORAGE:

Store in a cool, dry place.

Standard Shaft Configuration


Stiff Shaft Configuration

Angled Tip Configuration

Straight Tip Configuration

PREPARATION FOR USE:

1. Before removing the hydrophilic wire guide from its dispenser, inject sterile heparinized saline solution into the Luer lock hub end of the dispenser.
2. Inject enough solution to fill the dispenser coil. This will completely cover the wire guide surface and activate the hydrophilic coating.
3. Remove the hydrophilic wire guide from its dispenser by gently withdrawing the wire's tip.
4. If the hydrophilic wire guide cannot easily be removed from its dispenser, inject more heparinized saline solution into the dispenser and then try again.

DIRECTIONS FOR USE:

1. Fill catheter or other device with heparinized saline solution before and during use to ensure smooth movement of the hydrophilic wire guide within the device.
2. When wet with saline solution or blood, the hydrophilic wire guide will become lubricious. Use of sterilized gauze moistened with heparinized saline solution facilitates handling the wire.
3. Insert the wire guide into the device and advance to desired position.

NOTE: If movement of the wire within the device becomes diminished, remove the wire guide and reactivate the hydrophilic coating by wetting its entire surface with heparinized saline solution.

AFTER USE:

1. Excess blood can be removed from the wire guide by wiping the surface once or twice with gauze moistened with heparinized saline solution. **DO NOT use dry gauze** as this may damage the wire surface resulting in increased resistance when the wire is reinserted into the device.
2. Use of alcohol, antiseptic solutions or other solvents must be avoided, because they may adversely affect the surface of the hydrophilic wire guide.
3. After cleaning the wire, replace it by its proximal end into its dispenser filled with heparinized saline solution. This wire guide can only be used during the same procedure on the same patient.
4. Upon reinsertion into the device, if the wire guide does not move with initial ease, exchange it for a new HiWire™ Hydrophilic Wire Guide.

WARNINGS:

1. The hydrophilic wire guide may slide entirely into the catheter, sheath introducer, vessel dilator or other device because of its low sliding friction. To prevent this, keep at least 5 cm of the wire protruding from the device fitting at all times.
2. To prevent possible tissue damage, care should be taken when manipulating a device over a wire guide during the device's placement and withdrawal. If resistance is felt during device placement, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of

resistance cannot be determined, remove the wire guide and device as a unit to prevent possible damage and/or complications.

3. When using a wire guide, potential exists for thrombus formation or emboli, arterial or venous wall damage and/or plaque dislodgment. The physician should be familiar with the literature concerning the complications of angiography.

4. Cook does not recommend a particular technique for the use of this wire guide. The steps contained in the preceding directions are for information purposes only. Each physician should evaluate their appropriateness according to individual patient condition and his or her medical training and experience.

CAUTIONS:

- Prior to use, inspect for damage. If damaged, DO NOT USE.
- This device is sterilized by ethylene oxide gas and intended for one-time use only. Do not use any unit if its package is opened or damaged. Do not resterilize and/or reuse.
- Carefully read all instructions prior to use. Observe all warnings and cautions. Failure to do so may result in complications.
- Avoid manipulating or withdrawing the hydrophilic wire guide back through a metal needle or cannula. A sharp edge may scrape the coating or shear the wire guide. A catheter, introducer sheath or vessel dilator should replace the needle as soon as the wire guide has been inserted in the vessel.
- It is recommended that a plastic torque device be used to handle the hydrophilic wire guide. Use of a metal torque device may damage the wire.

ROADRUNNER® PC WIRE GUIDES WITH AQ® HYDROPHILIC COATING

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Roadrunner PC Wire Guide is coated with AQ (a biocompatible hydrophilic coating), which when activated becomes lubricous. The wire guide is made of a nitinol core with distal platinum inner spring coil tip, which is completely covered by a bismuth-impregnated polymer jacket

	Taper Length (cm)
Firm	10.5
Firm LT	18.5
Nimble	11.5
Nimble Floppy	9

INTENDED USE

The Roadrunner PC Wire Guide is used for catheter positioning and exchange in diagnostic and interventional procedure, exclusive of the coronary arteries.

CONTRAINDICATIONS

None known

WARNINGS

- This product is not intended for use in the coronary arteries.
- This product is a delicate instrument. Avoid forceful angulation.
- Avoid manipulation or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape or shear material from the wire guide.
- Altering the tip's configuration or curve manually may damage the wire guide.
- Excessive tightening of a torque device may abrade the coating on the wire guide.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of devices in angiographic procedures should be employed.
- Use medical imaging when manipulating the wire guide. Do not advance or torque the wire guide without visual evidence of the corresponding movement of the distal tip.
- When using the wire guide with another device, consider the end-hole size and the length of the device in order to ensure a proper fit between the wire guide and the device.
- Always keep the surface of the wire guide wet during use for optimal results.

INSTRUCTIONS FOR USE

Activating the AQR Hydrophilic Coating


Fig. 1

1. Flush the wire guide holder by attaching or pressing a syringe with saline, heparinized saline solution, or sterile water to the fitting or opening of the wire guide holder. Inject enough solution to wet the wire surface entirely. This will activate the AQ® coating.
2. Carefully remove the wire guide from the holder.
3. Insert a wire guide insertion tool through the valve assembly of the guiding sheath or guiding catheter. Insert the distal tip of the wire guide through the insertion tool.
4. Attach a torque device to the wire guide (if provided).
5. Standard wire guide endovascular techniques may now be employed.
NOTE: An initial flush of heparinized saline solution through the catheter lumen will enhance initial movement of the wire guide through the catheter.
6. When the wire guide is removed from a catheter, wipe it with a sterile gauze pad moistened with saline, heparinized saline solution, or sterile water and place it back into the original holder. Fill the holder with heparinized saline solution to maintain the lubricous state of the surface of the product for subsequent use.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

Roadrunner ® extra support Wire Guide

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Roadrunner Extra Support Wire Guide features a nitinol shaft and a short, angled platinum coil tip.

INTENDED USE

The Roadrunner Extra Support Wire Guide is intended for use in facilitating delivery of percutaneous catheters into the cardiovascular system.

CONTRAINDICATIONS

None known

WARNINGS

- The wire guide should be advanced only when visualizing its tip fluoroscopically. Do not torque wire guide without evidence of corresponding movement of distal tip.

- Refer to labeling of percutaneous cardiovascular access catheter for contraindications and potential complications associated with use.
- Carefully read all instructions prior to use. Failure to observe all warnings and precautions may result in complications.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in cardiovascular system access procedures. Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.
- Withdrawal or manipulation of distal spring coil portion of wire guide through needle tip may result in breakage.
- This wire guide is a delicate instrument and should be handled carefully; forceful angulation should be avoided.

INSTRUCTIONS FOR USE

1. Attach pin vise to wire guide.
2. Under fluoroscopy and while maintaining position of cardiovascular access catheter, advance wire guide to targeted site.
WARNING: Observe all wire guide movement in vessels under fluoroscopy. Do not torque wire guide without evidence of corresponding movement of distal tip. Further torquing against resistance may cause vessel trauma or wire guide damage, which may lead to device fracture.
WARNING: Wire guide should be advanced only when visualizing wire guide tip fluoroscopically. If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take action necessary to relieve resistance. Advancement and withdrawal of the wire guide should be performed slowly and carefully.
3. Upon reaching desired location, advance distal tip of wire guide distal to lesion, placing a more supportive section of wire guide at targeted site.
4. Confirm location of wire guide in two planes to assure placement of distal tip in desired vessel and not within a side branching vessel.
5. With wire guide secured in place, advance cardiovascular access catheter to target site.

Wire Guide Exchange

1. Remove pin vise from wire guide.
2. Open hemostatic valve and flush line.
WARNING: Always maintain a continuous flush while exchanging wire guides to prevent air entry into the system.
3. Withdraw wire guide under fluoroscopy slowly and gently to prevent trauma. **NOTE:** After removing wire guide from catheter, wipe with sterile gauze pad moistened with heparinized saline solution or sterile water and place back into holder.
4. Close hemostatic valve, flush line, and prepare a new wire guide.
WARNING: Prior to re-inserting exchanged wire guide, verify that distal tip of cardiovascular device is free within vessel. If tip is against vessel wall, vessel trauma may result upon re-introduction of wire guide.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the

product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.


Reuter Tip Deflecting Wire Guide Handle and Wire Guides

Caution: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Reuter Tip Deflecting Wire Guides are constructed of stainless steel and are available with an optional heparin coating. The wire guides have a tip-deflected curve radius of 5 or 10 mm.

The Reuter Tip Deflecting Wire Guide Handle is a metal handle with a replaceable plastic side-arm rotating adapter for catheter flushing.

INTENDED USE

Reuter Tip Deflecting Wire Guides are intended for use with the Reuter Tip Deflecting Handle assembly for curving or deflecting catheter tips for pulmonary angiography, selective angiography, translumbar aortography, bronchography, repositioning of central venous catheter tips and other vascular and non-vascular applications.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

None known

WARNINGS

- The Reuter Tip Deflecting Wire Guide Handle is shipped non-sterile. Sterilize prior to use.
- Exercise caution when using a Reuter Tip Deflecting Wire Guide. The distal tip of the wire guide should not be advanced beyond the distal tip of the catheter. Damage to the vessel wall can occur.

- Do not use excessive force when deflecting a Reuter Tip Deflecting Wire Guide. Excessive use of force can result in damage to the wire guide.
- Do not attempt deflection of wire guide until it is inserted into a catheter. Such deflection may result in damage to tip of wire.
- Table-top test the handle, wire guide and catheter assembly prior to use. Confirm that the wire guide does not protrude through or obstruct the catheter endhole. Ensure that the wire guide, when deflected, cannot exit through a catheter sideport.

PRECAUTIONS

- Use new plastic side-arm fitting at the beginning of each procedure. The order number for the fitting is **TBY-RA (Global Part Number G02666)**.
- Special care should be taken when selecting a catheter for use with a Reuter Tip Deflecting Wire Guide. The inner lumen of the catheter must be larger than the outer dimension of the wire guide.
- Always monitor manipulation using fluoroscopic guidance.
- A catheter with sideports may be chosen for use with the Reuter Tip Deflecting Wire Guide. However, care should be taken to ensure that the deflecting wire guide tip cannot accidentally exit through a catheter sideport.
- Read and follow the Product Recommendations presented in this IFU.

TIP DEFLECTING WIRE GUIDE

TIP DEFLECTING WIRE GUIDE HANDLE

5 mm TIP DEFLECTED

10 mm TIP DEFLECTED

PRODUCT RECOMMENDATIONS


Fig. 1

Assembly

WARNING: The Reuter Tip Deflecting Wire Guide Handle is shipped non-sterile. Sterilize prior to use.

NOTE: Referring to Fig. 1, ensure all parts are present and in working condition.

1. After opening and inspecting the handle and deflecting wire guide, insert proximal end of deflecting wire guide into large hole in front end of handle.
2. Move ring slide grip forward toward deflecting wire guide until end of deflecting wire guide can pass through small slot on side of ring slide grip. Tighten wire guide lock screw on ring. Slide to seat proximal end of deflecting wire guide.
3. Insert deflecting wire guide into sealing cap end of side-arm fitting. Slide side-arm fitting on deflecting wire guide until sealing cap is well seated in hole in front end of handle.
4. Anchor side-arm fitting to handle by tightening lock screw.

CAUTION: Use new plastic side-arm fitting at the beginning of each procedure.

CAUTION: The sealing feature on the side-arm fitting does restrict movement of the deflecting wire guide. It is suggested that the sealing cap be tightened only slightly to allow free movement of the deflecting wire guide and still prevent leakage through seal.


Fig. 2

Care of Reuter Tip Deflecting Wire Guide Handle

1. Clean handle in the same manner as surgical instruments.
2. Use gas or autoclave to sterilize.
3. The plastic side-arm fitting must be replaced after each use. Do not attempt to reuse the fitting. The order number for the fitting (with rubber "O" ring) is TBY-RA (Global Part Number G02666).

INSTRUCTIONS FOR USE

CAUTION: Throughout duration of procedure, periodically use side-arm fitting to flush system with heparinized saline solution to prevent clot formation on catheter or deflecting wire guide surface.

1. After catheter is properly positioned in vessel, perform a wire guide exchange for the Reuter Tip Deflecting Wire Guide/Handle assembly. Attach the male Luer lock on the side-arm fitting to the catheter hub.


Fig. 3

2. To deflect the wire guide and catheter, insert thumb and fingers in proper rings and slowly but firmly pull. (Figure 3)

WARNING: Exercise caution when using a Reuter Tip Deflecting Wire Guide. The distal tip of the wire guide should not be advanced beyond the distal tip of the catheter. Damage to the vessel wall can occur.

HOW SUPPLIED

Wire Guides: Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

Wire Guide Handle: Supplied non-sterile in peel-open packages. Sterilize prior to use. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

TORQ-FLEXR WIRE GUIDE TFE Coated

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The TFE coated Torq-Flex Wire Guide is etch marked in five centimeter increments from the distal tip.

INTENDED USE

The TFE coated Torq-Flex Wire Guide is intended for placement of Turbo-Flo® Over-The-Wire Peripherally Inserted Central Venous Catheter.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in venous access techniques. Standard techniques for placement of peripherally inserted central venous catheters should be employed.
- When exchanging a Twin Turbo-Flo® double lumen catheter, ensure that the wire guide is inserted into the hub marked "LARGE" (round lumen).

PRODUCT RECOMMENDATIONS

This wire guide is recommended for use only with Turbo-Flo® peripherally inserted central venous catheters.

INSTRUCTIONS FOR USE

Turbo-Flo Over-The-Wire Peripherally Inserted Central Venous Catheter placement using the Micropuncture® Peel-Away® Introducer technique with the graduated TFE coated Torq-Flex Wire Guide.


Fig. 1

1. After prepping the access site, introduce a thinwall percutaneous entry needle into the vessel. **NOTE:** The use of ultrasound is helpful to determine suitability for vessel access and patency. The EchoTip® marking on the needle is used to locate the tip of the needle during access into the vessel. (Fig. 1)


Fig. 2

2. Using fluoroscopic guidance, introduce the wire guide through the needle; advance the wire guide 15 to 20 cm into the vessel. (Fig. 2)


Fig. 3

3. Leaving the wire guide in place, withdraw the needle. If necessary, enlarge the puncture site with a #11 scalpel blade. (Fig. 3)


Fig. 4

4. Introduce a sheath-introducer assembly which has an internal diameter of .018 inch (0.46 cm) over the wire guide. With a twisting motion, advance the assembly into the vessel. (Fig. 4)


Fig. 5

5. Using fluoroscopic control, determine the correct catheter length by advancing the wire guide to the desired catheter tip location. Once the wire guide tip is in proper position, and using the proximal portion of the wire that is external to the patient, measure the distance on the wire guide from the triple etch marks to the puncture site. (The distal 60 cm of the wire guide is marked in 5 cm increments, with the triple etch marks positioned proximally at 60 cm.) (Fig. 5) Trim the Turbo-Flo catheter to the appropriate (60 cm minus X cm) length. A Turbo-Flo catheter's original length is 60 cm.

6. Continue catheter placement as designated by institutional protocol.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

Fixed Core Wire Guides

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Fixed Core Wire Guides are made with a stainless steel inner mandril and stainless steel outer coil. They are available in a variety of diameters, lengths, and tip configurations. Some devices are provided with a green PTFE coating to ease coaxial use with other devices. Refer to the product label for product specifications.

INTENDED USE

Fixed Core Wire Guides are intended to facilitate the placement of devices during diagnostic and interventional procedures.

CONTRAINDICATIONS

None known

WARNINGS

- This product is not intended for use in the coronary arteries.
- This product is a delicate instrument. Avoid forceful angulation.
- Avoid manipulating or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape or shear material from the wire guide.
- Altering the tip's configuration or curve manually may damage the wire guide.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of devices in angiographic procedures should be employed.

- Do not attempt to torque the wire guide.
- Use medical imaging when you manipulate the wire guide. Do not advance or manipulate the wire guide without visual evidence of the corresponding movement of the distal tip.
- When you use the wire guide with another device, consider the end-hole size and the length of the device in order to ensure a proper fit between the wire guide and the device.
- Do not advance or withdraw a wire guide when resistance is encountered, as a perforation could occur.
- Use only the flexible end as the initial inserted end. Inserting from the rigid end may cause damage to tissue or the device.
- Inspect the wire guide for kinks or damage prior to use.

INSTRUCTIONS FOR USE

1. Flush the wire guide holder by attaching a syringe with heparinized saline or sterile water to the fitting of the wire guide holder. Inject enough solution to wet the wire guide surface entirely.

NOTE: If flushing through the wire guide holder is not possible, remove the wire guide from the holder and place it in a bowl of heparinized saline or sterile water, or wet the wire guide surface over the entire length using gauze moistened with heparinized saline solution.

2. Carefully remove the wire guide from the holder.

3. If needed, insert a wire guide insertion tool (provided) through the valve assembly or hub of the guiding sheath or other interventional device. Insert the tip of the wire guide through the insertion tool.

4. Standard wire guide techniques may now be employed.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for more information.

Mandril Wire Guides

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Mandril Guides are made with a core mandril of stainless steel or nitinol material that runs the entire length of the product. Attached at the distal end is a coil of stainless steel, platinum, or palladium of varying length. Mandril Wire Guides come in a variety of lengths, diameters, and tip configurations. Some devices are coated with TFE coating to ease coaxial use with other devices. Some devices are provided with two usable ends. Refer to the product label for product specifications.

INTENDED USE

Mandril Wire Guides are used to facilitate the placement of devices during diagnostic and interventional procedure.

CONTRAINDICATIONS

None known

WARNINGS

- This product is not intended for use in the coronary arteries.
- This product is a delicate instrument. Avoid forceful angulation.
- Avoid manipulating or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape or shear material from the wire guide.
- Altering the tip's configuration or curve manually may damage the wire guide.
- Excessive tightening of a torque device may abrade the coating on the wire guide.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of devices in angiographic procedures should be employed.
- Use medical imaging when you manipulate the wire guide. Do not advance or torque the wire guide without visual evidence of the corresponding movement of the distal tip.
- When you use the wire guide with another device, consider the end-hole size and the length of the device in order to ensure a proper fit between the wire guide and the device.

INSTRUCTIONS FOR USE

1. Flush the wire guide holder by attaching a syringe with heparinized saline or sterile water to the fitting of the wire guide holder. Inject enough solution to wet the wire guide surface entirely.

NOTE: If flushing through the wire guide holder is not possible, remove the wire guide from the holder and place it in a bowl of heparinized saline or sterile water, or wet the wire guide surface over the entire length using gauze moistened with heparinized saline solution.

2. Carefully remove the wire guide from the holder.

3. If needed, insert a wire guide insertion tool (provided) through the valve assembly or hub of the guiding sheath or other interventional device. Insert the tip of the wire guide through the insertion tool and advance the wire guide to the desired location.

4. Attach a torque device to the wire guide (if provided).

5. Standard wire guide techniques may now be employed.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for more information.

Medical devices, which are not required to have an Instruction for Use according to the EU Regulations.

Products: Endovascular Wire Guides

The Vascular Wire Guides are intended to assist in inserting catheters, introducer sheaths and other medical devices in vascular applications.

Vascular wire guides designed for specific intended uses include the Roadrunner® Wire Guides, Roadrunner® Extra Support Wire Guides and the Cook Incorporated Micro-wires identified as the Approach® Hydro ST Micro Wire Guides and Roadrunner® UniGlide™ Hydrophilic Wire Guides.

The following statements regarding device specific intended uses are provided in the following sentences:

- Roadrunner® Wire Guides are specifically intended for use in complex diagnostic and interventional procedures.
- The Roadrunner® UniGlide™ Hydrophilic Wire Guide is intended for use to facilitate the delivery of percutaneous catheters into the peripheral vasculature.
- The Roadrunner® Extra Support Wire Guides are intended for use in complex diagnostic and interventional procedures where extra support is needed for catheter exchange or manipulation of devices and in facilitating delivery of percutaneous catheters into the cardiovascular system.
- The Cook Approach® Micro Wire Guides (Approach® Hydro ST) Wire Guides are intended for use in facilitating delivery of percutaneous catheters into the peripheral vasculature.

Hereinafter the Vascular Wire Guides will be referred to according to their respective device categorization when describing specific attributes of each particular device. These general categorizations are based upon differences in materials, device attributes and specifications. There are devices categorizations are as follows: Fixed Core, Moveable Core, Mandril, Roadrunner®, Approach® Wire Guides and the Alight™ Superselective Microwire Guide.

Fixed Core Wire Guides

Fixed Core Wire Guides are designed with a tapered fixed core and a tapered mandril wire guide for a gradual transition to a flexible tip. These wire guides are used for vascular applications but may also be used at times for non-vascular applications. Fixed Core Wire Guides are manufactured with stainless steel and stainless steel with PTFE coating, and can be manufactured with a distal straight and/or curved distal tip. These wire guides are manufactured with safety wire, a tapered fixed core construction, and with a gradual or blunt transition to the distal tip. The taper of the flexible tip is designed relative to the curve radius. Fixed Core Wire Guides may also be manufactured with large diameters, and with a variety of curve radii. The Fixed Core Wire Guides are manufactured in

diameter measurements within the range of 0.015-0.052 inches (0.038 – 0.13 cm), and varying lengths ranging from 15-400 centimeters. Distal tip curve radius of these fixed core wire guides may fall within a range of 2-15 millimeters or may follow a template to form various curves. All device diameters, lengths and curve radius are specific to the prefix in the product catalogue number.

Fixed Core Wire Guides may also be manufactured with different length tapered mandril wires, allowing for a tapered long flexible tip (Newton Wire Guide), manufactured with different degrees of coil relaxation in the distal tip (Bentson Wire Guides), manufactured with an increased diameter of the wire guide and stiffer mandril wire (Amplatz Wire Guides), manufactured with a straight and a curved flexible tip (Double Flexible Tipped Wire Guide), and manufactured with a heavy duty stainless steel core. In addition, specific distal curve formations (Safe-T-J® Multiple Curved Wire Guide) and distal curve radius (Rosen Curved Wire Guide) are also manufactured as a Fixed Core Wire Guides.

Safe-T-J® Fixed Core Wire Guides

The Safe-T-J® Fixed Core Wire Guides are manufactured from stainless steel as a straight wire guide (SF) or a PTFE coated stainless steel wire guide (TSF) with a flexible tip. These wire guides are manufactured in diameters ranging from 0.015 to 0.038 inches (0.038 – 0.13 cm) and in lengths ranging from 15 to 80 centimeters.

The Safe-T-J® Fixed Core Wire Guides are also manufactured from stainless steel with a distal curve configuration (SCF) or as a curved wire guide with PTFE coating (TSCF). These wire guides are manufactured in diameters ranging from 0.015 to 0.052 inches (0.038 – 0.13 cm) and in lengths ranging from 15 to 260 centimeters. The distal curves for the uncoated wire guides are designed with a 2 or 3 millimeter radius, in contrast to the PTFE coated wires which are manufactured with a 1.5, 3, 7.5 or 15 millimeter curve radius. The available length measurements and curve radius are specific to the wire guide diameter.

PTFE coated stainless steel curved Safe-T-J® Fixed Core Wire Guides are also manufactured with a larger mandril (THSCF) wire for a stiffer body. These wire guides are manufactured in diameters ranging between 0.025 to 0.038 inches (0.064 – 0.097 cm) and lengths from 40 to 260 centimeters and distal curve radius of 1.5, 3, 7.5 and 15 millimeters.

Product catalogue number prefixes of devices are used to identify the basic construction features of the device. The letters SF represents a Stainless Steel Fixed Core Wire Guide. When a wire guide is manufactured with a distal curve radius a C is used to represent the presence of a curve (SCF). Many wire guides are manufactured with green PTFE coating and the letter T is used to indicate this characteristic, for example TSF, a Stainless Steel Fixed Core Wire Guide with PTFE coating. Wire guides manufactured with a Heavy duty core use the letter H in the product catalogue number prefix to indicate this construction feature, such as THSF. The following methodology is used to identify Safe-T-J® Fixed Core Wire Guides.

Table. Methodology for Safe-T-J® Fixed Core Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Distal Curve Radius (mm)	Additional Identifier
SF	18	50		

SCF	18	40	2	
TSF	35	145		
TSCF	35	60	3	PNS
TSCF	52	55	3	PTS

The additional identifiers used in the product catalogue numbers for some of the Safe-T-J[®] Fixed Core Wire Guides in the previous table are used to represent specific characteristics of the finished device. Some additional identifiers are used to indicate the set of which the wire guide is a component of, such as PTS, Percutaneous Tracheostomy Set and PNS, Percutaneous Nephrostomy Drainage Set. Some additional identifiers are used to indicate a change from the standard packaging. The letters SGH indicates that the wire guide is packaged inside a Spiral Wire Guide Holder instead of the standard wire guide holder.

Newton Wire Guide

Newton Wire Guides can be designed with a ten centimeter tapered core (TSFNA) or a fifteen centimeter tapered core (TSFNB). All Newton Wire Guides are designed with a long tapered flexible tip for negotiating tortuous anatomical strictures. Newton Wire Guides are manufactured using PTFE coated stainless steel coils, a tapered mandril wire and safety wire. Newton Wire Guides manufactured as a straight wire are manufactured in a diameter of 0.035 inches (0.089 cm) and lengths ranging from 145 to 260. The Newton Wire Guides are used for vascular applications but may also be used for non-vascular applications.

Newton Wire Guides are also manufactured with a distal curve radius and designed with a ten centimeter tapered core (TSCFNA) or a fifteen centimeter tapered core (TSCFNB). These Newton wire guides are manufactured in diameters ranging from 0.018 to 0.038 inches (0.046 – 0.097 cm) and in lengths ranging from 145 to 260 centimeters. The curved Newton wire guide is also designed with a distal curve radius of 3 15 millimeters.

These wire guides are also manufactured with a larger diameter core with a fifteen centimeter tapered core (THSFNB) and have all the characteristics of the standard Newton Wire Guides, with the exception of a stiffer body. These Newton Wire Guides are manufactured in a 0.035 inch (0.089 cm) diameter and a length of 260 centimeters.

The product catalogue number prefixes for the Newton Wire Guides also identify the basic construction features of the device. The letters TSFNA represents a PTFE coated Stainless Steel Fixed Core Wire Guide with a Newton ten centimeter tapered core (A). The letters TSFNB represents a PTFE coated Stainless Steel Fixed Core Wire Guide with a Newton fifteen centimeter tapered core (B). These wire guides manufactured with a curve (C) are represented by the letters TSCFNA and TSCFNB. In a similar fashion Newton Wire Guides are manufactured with a Heavy duty core and are represented by the letters THSFNB.

Table. Methodology for Newton Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Distal Curve Radius (mm)	Additional Identifier
TSFNA	35	145		
TSFNB	35	145		

Bentson Wire Guide

Bentson Wire Guides are designed with a floppy flexible tip (TSFB) and are manufactured with PTFE coated stainless steel coils, a large diameter core, a tapered mandril wire and a safety wire. These wires are manufactured without a set distal curve, but the degree of coil relaxation provides a 20 centimeter flexible tip with increased floppiness over the distal 6 centimeters. Bentson Wire Guides are manufactured in 0.018, 0.025, and 0.035 inch (0.046, 0.064, 0.097 cm) diameters and lengths of 80 to 260 centimeters. The Bentson Wire Guides are used for vascular applications and may be used for non-vascular applications. Bentson Wire Guides manufactured in a length of 260 centimeters are typically used as an exchange wire guide during the medical procedure.

Bentson Wire Guides may also be manufactured with a distal curve radius (TSCFB) of 1.5 or 3 millimeters. These Bentson Wire Guides are manufactured in 0.035 inch (0.089 cm) diameters and lengths of 145 centimeters.

Bentson Wire Guides constructed as a Bentson Plus[®] (TSFBP) wire are designed with a flat wire outer coil, enabling for a larger inside diameter and a larger mandril wire. The construction of this wire provides for a stiffer body, with the same distal tip characteristics as other Bentson Wire Guides. Wire guides with the Bentson Plus[®] design are manufactured in a 0.035 inch (0.089 cm) diameter and a length of 145, 180 and 260 centimeters.

Bentson Wire Guides are also manufactured with a moveable inner core (TSMB). This wire guide is designed with flat wire and a stainless steel cannulated tapered mandril. The moveable core is designed with a coons taper allowing for a gradual flexible tip transition. Although the Bentson Wire Guide with a moveable core is manufactured without a designated distal curve, it is designed with a very floppy distal six centimeters of coil allowing for an atraumatic tip.

The moveable core Bentson wire guides are manufactured in a 0.035 inch (0.089 cm) diameter and a length of 145 centimeters.

The product catalogue number prefixes for the Bentson Wire Guides also identify the basic construction features of the device. The letters TSFB represents a PTFE coated Stainless Steel Fixed Core Wire Guide with a Bentson flexible floppy distal tip. The letters TSCFB represents a PTFE coated Stainless Steel Curved Fixed Core Wire Guide with a Bentson flexible floppy tip. The letters TSMB identifies a PTFE coated Stainless Steel Moveable Core Bentson Wire Guide and the letters TSFBP refers to the PTFE coated Stainless Steel Fixed Core Bentson Plus[®] Wire Guide.

Table. Methodology for Benston Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Distal Curve Radius (mm)	Additional Identifier	Additional Identifier
TSFB	18	180			
TSCFB	35	145	3		
TSFBP	35	260			
TSMB	35	145			

Amplatz Wire Guides

Amplatz Wire Guides are manufactured with varying degrees of stiffness in the shaft and are labelled as still, ultra stiff and extra stiff. These wire guides may be used for interventional cardiovascular and catheter exchanges or may be used for drainage applications and urology procedures. Amplatz Wire Guides are designed with an increased inner diameter of wire guide coil allowing for the insertion of a stiffer mandril wire. The Amplatz Wire Guides are manufactured from PTFE coated stainless steel.

Amplatz Stiff Wire Guides

Amplatz Stiff Wire Guides (THSF-/ASG) are designed with an increased inner diameter of wire guide coil allowing for the insertion of a stiffer mandril wire. The Amplatz Stiff Wire Guides are manufactured with an 8 centimeter flexible tip length and a stiff shaft. These wire guides are manufactured with PTFE coating applied to the surface of the wire guide and are designed with a stiff shaft and a gradual transition to a very flexible distal tip. The Amplatz Stiff Wire Guides are manufactured in 0.035 inch (0.089 cm) diameters and lengths of 80, 145, 180 and 260 centimeters. Amplatz Stiff Wire Guides are only manufactured as a straight wire guide. These wire guides may be used for interventional cardiovascular and catheter exchanges or may be used for drainage applications and urology procedures. Amplatz Stiff Wire Guides manufactured in a length of 260 centimeters are typically used as an exchange wire guide during the medical procedure.

The following methodology is used to identify the Amplatz Stiff Wire Guides. The product catalogue number prefix THSF identifies the device design as a PTFE coated Heavy core Stainless steel Fixed core wire guide.

Table. for Amplatz Stiff Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Additional Identifier
THSF	35	260	ASG
THSF	35	145	ASG

The additional identifier ASG in the previous table identifies the wire guide as an Amplatz Stiff Wire Guide.

Amplatz Ultra Stiff Wire Guides

Amplatz Ultra Stiff Wire Guides (THSF-/AUS) are also designed with an increased inner diameter of wire guide coil allowing for the insertion of a stiffer mandril wire. The Amplatz Ultra Stiff Wire Guides are manufactured designed with an 8 centimeter flexible tip length. The Amplatz Ultra Stiff Wire Guides are manufactured with PTFE coating applied to the surface of the wire guide in 0.035 inch (0.089 cm) diameters and lengths of 80, 145 and 180 centimeters. These devices are manufactured with a straight flexible tip.

Amplatz Ultra Stiff Wire Guides are also manufactured in these two diameters and identical lengths with a distal curve radius of 3 millimeters (THSCF-/AUS). These wire guides may be used for interventional cardiovascular and catheter exchanges or may be used for drainage applications and urology procedures.

The following methodology is used to identify the Amplatz Ultra Stiff Wire Guides. The product catalogue number prefix THSF identifies the device design as a PTFE coated Heavy core Stainless steel Fixed core wire guide, and the catalogue number prefix THSCF identifies the device design as a PTFE coated Heavy core Stainless steel Curved Fixed core wire guide.

Table. Methodology for Amplatz Ultra Stiff Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Distal Curve Radius (mm)	Additional Identifier
THSF	35	180		AUS
THSF	35	145		AUS
THSCF	35	180	3	AUS

The additional identifier AUS in the previous table identifies the wire guide as an Amplatz Ultra Stiff Wire Guide.

Amplatz Extra Stiff Wire Guides

Amplatz Extra Stiff Wire Guides (HSF-/AES) are designed with an increased inner diameter using an extra stiff mandril wire tapered for a flexible distal tip. These wire guides are manufactured in 0.038 inch (0.097 cm) diameters and lengths of 40 centimeters as straight wire guides.

Amplatz Extra Stiff Wire Guides are also manufactured with PTFE coating applied to the surface of the wire guide (THSF-/AES) and as a straight wire. These wire guides are manufactured in 0.025, 0.035 and 0.038 inch (0.064, 0.089, 0.097 cm) diameters and lengths of 80, 100, 145, 180, 260, 300, and 400 centimeters. Amplatz Extra Stiff Wire Guides manufactured in a length of 260, 300 or 400 centimeters are typically used as an exchange wire guide during the medical procedure. A modification of this wire guide (THSF-/AES-ST) is designed with a shorter distal taper which provides for a stiffer tip.

Amplatz Extra Stiff Wire Guides are also manufactured with PTFE coating applied to the surface of the wire guide (THSCF-/AES) with a distal curve. These wire guides are manufactured in 0.025, 0.032, 0.035 and 0.038 inch (0.064, 0.089, 0.097 cm) diameters and lengths of 40, 80, 145, 180, and 260 centimeters with a 3, 7.5 or 15 millimeter curve radius. Similar wire guides are manufactured with a 0.028 inch (0.071 cm) diameter in a length of 40 centimeters with a 2 millimeter curve radius.

The following methodology is used to identify the Amplatz Extra Stiff Wire Guides. The product catalogue number prefix THSF identifies the device design as a PTFE coated Heavy core Stainless steel Fixed core wire guide, and the catalogue number prefix THSCF identifies the device design as a PTFE coated Heavy core Stainless steel Curved Fixed core wire guide.

Table. Methodology for Amplatz Extra Stiff Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Distal Curve Radius (mm)	Additional Identifier	Additional Identifier
THSCF	35	145	3	AES	
THSF	35	180	7.5	AES	
THSF	35	300		AES	ST

The additional identifier AES identifies the device as an Amplatz Extra Stiff Wire Guide.

Amplatz Extra Stiff Whisker Wire Guide

Amplatz Extra Stiff Whisker Wire Guides are manufactured with PTFE coating applied to the surface of the wire guide (THSF/-AESW) and as a straight wire. The Amplatz Extra Stiff Whisker Wire Guide is manufactured in 0.035 inch (0.089 cm) diameters and lengths of 260 centimeters. These wire guides are designed with an increased diameter of the wire guide coil allowing for an extra stiff mandril wire. Amplatz Extra Stiff Whisker Wire Guides are manufactured with a very short taper length of the mandril wire which provides added support at the distal tip.

The following methodology is used to identify the Amplatz Extra Stiff Whisker Wire Guides. The product catalogue number prefix THSF identifies the device design as a PTFE coated Heavy core Stainless steel Fixed core wire guide.

Table. Methodology for Amplatz Extra Stiff Whisker Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Additional Identifier
THSF	35	260	AESW

The additional identifier AESW is used to identify the device as an Amplatz Extra Stiff Whisker Wire Guide.

Curved Amplatz Apex Catheter Exchange Wire Guide

The Curved Amplatz Apex Catheter Exchange Wire Guide is manufactured from PTFE coated stainless steel coils, a stainless steel heavy duty cannulated mandril wire and a safety wire. The Curved Amplatz Apex Catheter Exchange Wire Guide is designed with an 18 centimeter flex length and a double curve formed to a specific template. This wire guide is manufactured with a larger inner core similar to the Amplatz wire guides and is manufactured in 0.035 inch (0.089 cm) diameters and a length of 260 centimeters.

The following methodology is used to identify the Curved Amplatz Apex Catheter Exchange Wire Guide. The product catalogue number prefix THSCF identifies the device design as a PTFE coated Heavy core Stainless steel Curved Fixed core wire guide.

Table. Methodology for Curved Amplatz Apex Catheter Exchange Wire Guide

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Distal Curve Formation
THSCF	35	260	APEX-1

Rosen Curved Wire Guide

Rosen Curved Wire Guides are manufactured from PTFE coated stainless steel and are designed with a heavy duty centerless ground mandril wire, a two centimeter flexible tip, safety wire, and a Rosen 1.5 millimeter curve radius. The Rosen 1.5 millimeter curve radius is characterized by an inward curve of the distal tip and formed to a specific template.

Rosen Curved Wire Guides are manufactured in 0.035 inch (0.089 cm) diameters and lengths of 80, 100, 145, 180, 220 and 260 centimeters. Rosen Curved Wire Guides are used for renal angioplasty catheter exchange.

The following methodology is used to identify the Rosen Curved Wire Guides. The product catalogue number prefix THSCF identifies the device design as a PTFE coated Heavy core Stainless steel Curved Fixed core wire guide.

Table. Methodology for Rosen Curved Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Distal Curve Radius (mm)
THSCF	35	145	1.5-Rosen
THSCF	35	80	1.5-Rosen
THSCF	35	220	1.5-Rosen

Double Flexible Tipped Wire Guide

Double Flexible Tipped Wire Guides are manufacture with PTFE coated stainless steel coils (C-TDOC). These devices are manufactured diameters of 0.035 (0.089 cm) inches and a length of 50 centimeters. The curve radius of the distal tip for these devices is 3 millimeters.

The following methodology is used to identify the Double Flexible Tipped Wire Guides. The letters TDOC identifies the device as a PTFE coated Double Flexible Tipped Wire Guide with a Straight (O) and a Curved (C) end.

Table. Methodology for Double Flexible Tipped Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Proximal Curve Radius (mm)	Distal Curve Radius (mm)	Additional Identifier
C-TDOC	35	50	0	3	

Heavy Double Flexible Tipped Wire Guides

Double Flexible Tipped Wire Guides may also be manufactured using a heavy duty mandril using stainless steel or PTFE coated stainless steel coils. Heavy Double Flexible Tipped Wire Guides (C-HDOC) is manufactured with flexibility of both ends for the dual purpose of interchangeable use. These wire guides are also constructed of a tapered mandril wire and a safety wire and are designed with a straight flexible tip and a curved flexible tip. These wire guides are manufactured from stainless steel and have an inner mandril which is tapered on both ends and a safety-wire. Double Flexible Tipped Wire Guides are used for access and catheterization procedures.

A modification of the Heavy Double Flexible Tipped Wire Guides is manufactured as an Amplatz extra stiff wire with a flat safety wire, and a tapered mandril with a specified relaxation of the coils (C-HDOC/-HTS). This modification of the Heavy Double Flexible Tipped Wire Guide is manufactured with a 0.032 inch (0.082 cm) diameters and a length of 60 centimeters. This device is also designed with a straight flexible tip and a curved tip with a 2.5 millimeter curve radius. These

wire guides are manufactured with stainless steel, offering a high tensile strength due to the specific application of the device, often used with a somewhat larger diameter catheter.

Heavy Double Flexible Tipped Wire Guides are also manufactured with PTFE coated stainless steel coils (C-TDOC). These devices are manufactured in identical fashion to the stainless steel double flexible tipped wire guides with the exception of a larger mandril wire which offers a stiffer shaft. PTFE coated Heavy Double Flexible Tipped Wire Guides are manufactured in 0.035 inch (0.089 cm) diameters and lengths of 50 centimeters. The curve radius of the distal tip for these devices is 3 millimeters.

The following methodology is used to identify the Heavy Double Flexible Tipped Wire Guides. The product catalogue number prefix HDOC identifies the device design as a Heavy Double Flexible Tipped Wire Guide with a Straight (O) and a Curved (C) end.

Table. Methodology for Heavy Double Flexible Tipped Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Proximal Curve Radius (mm)	Distal Curve Radius (mm)	Additional Identifier
C-HDOC	32	60	0	2.5	HTS
C-TDOC	35	50	0	3	

In the previous table the additional identifier HTS indicates that this particular wire guide is manufactured from stainless steel with a High Tensile Strength.

Moveable Core Wire Guides

Moveable Core Wire Guides are used to negotiate tortuous vessels. Moveable Core Wire Guides are designed with an increased inner diameter of wire guide coils and the use of a Tefcor polytetrafluoroethylene (PTFE) coated moveable core. Moveable Core Wire Guides are manufactured with PTFE coated stainless steel, and can be manufactured with a straight or curved distal tip. Moveable core wire guides are manufactured with a handle coil, safety wire, and a tapered moveable core with a gradual transition to the distal tip. Moveable core wire guides are designed with a five centimeter handle on the moveable core. The handle allows for manipulation of the distal tip or distal curve radius. The tapered distal end of the wire allows for a gradual flexible tip transition and is designed relative to the curve radius.

Moveable Core Wire Guides are manufactured with a Tefcor rounded end mandril wire and a safety wire. Moveable Core Wire Guides are manufactured with PTFE coating (TMT) applied to the wire guide coils. These Moveable Core Wire Guides are manufactured in 0.035 inch (0.089 cm) diameters and lengths of 145 centimeters.

PTFE coated Moveable Core Wire Guides manufactured with a Tefcor rounded end mandril wire is also manufactured with a distal curve (TCMT). These Moveable Core Wire Guides are manufactured in 0.032, 0.035 and 0.038 inch (0.081, 0.089, 0.097 cm) diameters and lengths of 80, 125 and 145 centimeters. These wire guides are manufactured with a 1.5, 3 or 7.5 millimeter curve radius.

Moveable wire guides are manufactured with a tapered tefcor mandril and are designed with the features of an Amplatz wire guide are also manufactured with a 3 millimeter distal curve (TCMTNA-/-AMC) configuration. This wire guide is manufactured with a 0.035 inch (0.089 cm) diameter and a length of 145 centimeters.

The following methodology is used to identify the Moveable Core Wire Guides. The letters TMT identifies the device as a wire guide manufactured with a PTFE coating and a Moveable Tefcor core. Coated wire guides manufactured with a distal curve are represented by the letters TCMT; PTFE coated Curved Moveable Tefcor. In addition, these devices designed with a longer flexible tip are identified with the letters TMTNA (PTFE Moveable Tefcor Newton A) and is identified with the letters TCMTNA when manufactured with a distal Curve.

Table. Methodology for Moveable Core Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Distal Curve Radius (mm)	Additional Identifier
TMT	35	145		
TCMT	32	80	3	SGH
TMTNA	35	145		
TCMTNA	35	145	3	AMC

The additional identifier in the previous table, SGH identifies a change from the standard packaging for this wire guide, and designates that the wire is placed inside a Spiral Guide Holder. The additional identifier AMC indicates that the device is manufactured with an Amplatz Moveable Core.

Mandril Wire Guides

Mandril Wire Guides are typically manufactured with stainless steel and can be manufactured with a platinum, stainless steel or nitinol distal tip. The shaft of the wire guide is usually manufactured from stainless steel, with the distal end of the wire guide being constructed of a spring coil tip over the tapered mandril wire. Some specialty mandril wire guides are manufactured utilizing a PTFE coated nitinol shaft, with a distal platinum tip. The length of the distal spring coil portion of any mandril wire guide is specific to the product catalogue number. Mandril Wire Guides are generally available in diameters within the range of 0.015 to 0.018 inches (0.039-0.046 cm) and are manufactured in varying lengths ranging from 30 to 480 centimeters. A template is used to form the specific distal end curve. All device diameters, length and curve type are specific to the prefix in the product catalogue number.

Mandril Wire Guides are constructed with a mandril shaft with the distal portion of the wire guide designed with a spring coil construction over the tapered mandril. The distal spring coils allow for flexibility at the distal tip. The mandril shaft provides rigidity comparable to a 0.89 millimeter diameter wire guide. Although similar, different types of mandril wire guides are designed with slightly different features. Mandril Wire Guides include Cope Mandril Wire Guides Torq-Flex[®] Wire Guides and Cor-Flex[™] Wire Guides.

Cope Mandril Wire Guides

Cope Mandril Wire Guides are manufactured with either a stainless steel mandril (PMG-/COPE) or a nitinol mandril wire (PMG-/COPE-NT). The distal spring coil portion of both types of Cope Mandril Wire Guides is constructed of 7 centimeter segment of platinum coils. Cope Mandril Wire Guides are used for complex diagnostic interventional procedures.

Cope Mandril Wire Guides manufactured with a stainless steel mandril (PMG-/COPE), and platinum coils are manufactured with 0.015, 0.018 and 0.018sp inch (0.039-0.046 cm) diameters. The 0.018SP diameter wire guides are slightly smaller in diameter than the 0.018 inch (0.046 cm) diameter guides. The distal platinum coils are formed to a template to create the COPE curve and flexible distal tip. These wire guides are manufactured in a length of 40, 60, 80, 100, 125 and 145 centimeters.

Cope Mandril Wire Guides manufactured with a nitinol mandril wire (PMG-/COPE-NT), and platinum coils are manufactured with 0.018 and 0.018SP inch (0.046 cm) diameters and in lengths of 25, 30, 40, 50, 60, 80, 100, 125, and 145 centimeters. The inner mandril wire has a long taper and allows for a stiffer wire. The coil length for this wire guide measures approximately 6.5 centimeters. These wire guides are also curved to match the COPE template.

A modification of the Cope Mandril Wire Guides manufactured with a nitinol mandril wire and platinum coils device is manufactured with a shorter taper of the mandril and allows for more flexibility (PMG-/COPE-NT-ST) of the coils. The Cope Mandril Wire Guides designed with a shorter tip are characterized by a shorter coil length ranging from 3 to 4 centimeters. These wire guides are manufactured with a 0.018SP diameter and lengths of 60, 80, 100, and 180 centimeters.

The following methodology is used to identify the Cope Mandril Wire Guides. The product catalogue number prefix PMG identifies the device as a Platinum Mandril Guide.

Table. Methodology for Cope Mandril Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Additional Identifier (Distal Tip)	Mandril Wire (NT=Nitinol)	Additional Identifier
PMG	15	100	COPE		
PMG	18	125	COPE		
PMG	18	145	COPE	NT	
PMG	18SP	50	COPE	NT	
PMG	18SP	100	COPE	NT	ST
PMG	18SP	80	COPE	NT	ST

The additional identifier ST indicates that the wire guide is manufactured with a Short Tip with approximately half the distal coil length as similar wire guides.

Torq-Flex® Wire Guides

Torq-Flex® Wire Guides are designed for torque control and are manufactured with a stainless steel mandril and may have stainless steel coils (STF). Torq-Flex® Wire Guides may also be designed

with PTFE coated stainless steel mandril (THTF) and platinum coils. These wire guides are used for small vessel procedures, peripheral angioplasty and studies of distal branches in the vasculature.

Torq-Flex[®] Wire Guides manufactured with stainless steel mandril and stainless steel coils (STF-AUST) are designed with a distal spring coil tip measurement of 8 centimeters. These wire guides are manufactured in 0.018 inch (0.046 cm) diameters and lengths of 40, 60 and 80 centimeters.

The Torq-Flex[®] Wire Guide identified with the product catalogue number THTF-18SP-130-AUST-ST is manufactured with a 0.018sp (0.046 cm) diameter and a length of 130 centimeters. However, this wire guide is designed with a shorter distal taper to create a stiffer distal tip.

The following methodology is used to identify the Torq-Flex[®] Wire Guides. The product catalogue number prefix STF identifies the device as a Stainless steel Torq-Flex[®] Wire Guide. The letters THTF identifies the devices as a PTFE coated Torq-Flex[®] Wire Guide with platinum coils.

Table. Methodology for Torq-Flex[™] Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Additional Identifier
STF	18	40	AUST
STF	18	60	AUST
THTF	18SP	130	AUST-ST

The additional identifier AUST identifies the device as the Australian modification of the initial device and is designed with a shorter platinum tip. The additional identifier AUST-ST identifies the Australian modification as being manufactured with a Short Taper at the distal tip allowing for a stiffer tip than similar devices and is also designed with the shorter platinum tip.

Cor-Flex[™] Wire Guides

Cor-Flex[™] Wire Guides are designed with a stainless steel mandril wire, safety wire and stainless steel distal spring coils. These devices have a safety wire soldered onto the tapered end of the mandril wire. The Cor-Flex[™] Wire Guides may also be manufactured as a double-ended (CFWD) Cor-Flex mandril wires. These wire guides are used for small vessel procedures, peripheral angioplasty and studies of distal branches or tight lesions in the vasculature. All the Cor-Flex[™] Wire Guides are manufactured with a 0.018 inch diameter.

Double-ended (CFWD) Cor-Flex[™] Wire Guides manufactured in lengths of 30 centimeters and have a straight flexible and a curved end with a 2 millimeter curve radius. The manufacturing of a double ended wire guide requires a safety wire soldered to both ends of the tapered mandril wire.

Table. Methodology for Cor-Flex[™] Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Secondary Curve Radius (mm)	Distal Curve Radius (mm)	Distal Coils (P=Platinum)	Additional Identifier
CFWD	18	30	0	2		

Micropuncture Nitinol Mandril Wire Guides

Micropuncture Nitinol Mandril Wire Guides are manufactured with a nitinol mandril with a distal spring coil portion constructed of a 5 centimeter segment of palladium coils. Micropuncture Nitinol Mandril Wire Guides are used for complex diagnostic interventional procedures.

Micropuncture Nitinol Mandril Wire Guides manufactured with a nitinol mandril and palladium coils are manufactured with 0.018 inch (0.046 cm) diameters. The distal palladium coils are formed to a template to create the flexible distal tip. (PMG-/MPIS) These wire guides are manufactured in a length of 40, 60 and 80 centimeters.

Micropuncture Nitinol Mandril Wire Guides are also manufactured with a distinct distal curve formation (PMG-/MPIS-C). These wire guides are manufactured in a 0.018 inch (0.046 cm) diameter and lengths of 40, 60 and 80 centimeters. The distal tip curve is also formed to a template; however for this device the curve is set. The formation of the distal tip is similar to the Cope Mandril Wire Guides.

The following methodology is used to identify the Micropuncture Nitinol Mandril Wire Guides. The product catalogue number prefix PMG identifies the device as a Palladium Mandril Guide.

Table. Methodology for Micropuncture Nitinol Mandril Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Additional Identifier	Distal Tip Formation
PMG	18	40	MPIS	
PMG	18	60	MPIS	
PMG	18	80	MPIS	
PMG	18	40	MPIS	C
PMG	18	60	MPIS	C
PMG	18	80	MPIS	C

The additional identifier MPIS indicates that the wire guide is manufactured as a compatible wire guide with a Micropuncture Introducer Set. The letter C associated with the distal tip formation indicates that the distal tip of this wire guide has a set Curve, in contrast to similar wire guides that are only designed with distal tip flexibility.

Roadrunner® Wire Guides

Roadrunner® Wire Guides are manufactured with a nitinol core and a radiopaque platinum distal spring coil tip. These devices have an outer polyurethane sleeve with Slip-Coat™ hydrophilic coating. The length of the distal taper of the spring coil portion of the wire guide is specific to the product catalogue number. Roadrunner® Wire Guides are manufactured with a short distal taper (RFSPC) with a flexible shaft or a standard distal taper (RFPC) with a flexible shaft. In addition, Roadrunner® Wire Guides are manufactured with a standard distal taper and a stiff shaft (RPC), or a long distal taper and a stiff shaft (RLPC). The distal tip curve of these wire guides is generally of a straight or angled configuration. All device diameters, length and curve types are specific to the

prefix in the product catalogue number. All Roadrunner® wire guides are supplied with slip-coating on the outer surface of the polymer jacket.

Roadrunner® Wire Guides are manufactured with a short distal taper (RFSPC) with a flexible shaft are manufactured in 0.035 and 0.038 inch (0.089- 0.097 cm) diameters and lengths of 145, 180 and 260 centimeters. These wire guides are manufactured from a polymer coated nitinol core secured to a platinum coils at the distal tip. These Roadrunner® Wire Guides are characterized by a nimble and angled distal tip. All Roadrunner® Wire Guides manufactured in a length of 260 centimeters are typically used as an exchange wire guide during the medical procedure.

Roadrunner® Wire Guides are manufactured with a standard distal taper (RFPC) with a flexible shaft are manufactured in 0.035 inch (0.089 cm) diameters and lengths of 145, 180 and 260 centimeters. These wire guides are manufactured from a polymer coated nitinol core secured to a platinum coils at the distal tip. These Roadrunner® Wire Guides are characterized by a nimble floppy and angled distal tip, however, can be manufactured with a nimble floppy straight distal tip (RFPC-35-145-0).

Roadrunner® Wire Guides are manufactured with a standard distal taper and a stiff shaft (RPC) are manufactured in 0.035 and 0.038 inch (0.089- 0.097 cm) diameters and lengths of 80, 145, 180 and centimeters. These wire guides are manufactured from a polymer coated nitinol core secured to a platinum coils at the distal tip. These Roadrunner® Wire Guides are characterized by a firm and angled distal tip, however, can be manufactured with a firm straight distal tip (RPC-35-145-0).

A modification of the Roadrunner® Wire Guides manufactured with a standard distal taper and a stiff shaft (RPC-/CONNORS) is designed with the distal 65 centimeters of the wire guide manufactured with a Bentson jacket and the remaining portion of the wire guide PTFE coated. These modifications are manufactured with a 0.035 inch (0.089 cm) diameter and a length of 145, 180 and 260 centimeters with an angled tip or in a length of 260 centimeters manufactured with a 2 millimeter distal curve radius (RPC-35-260-2-CONNORS).

Roadrunner® Wire Guides are manufactured with a long distal taper and a stiff shaft (RLPC) are manufactured in 0.035 inch (0.089 cm) diameters and lengths of 145, 180 and 260 centimeters. These wire guides are manufactured from a polymer coated nitinol core secured to a platinum coils at the distal tip. These Roadrunner® Wire Guides are characterized by an extra firm and angled distal tip

An Olcott Torque Device (OTD-100) is supplied with the Roadrunner® Wire Guides. These devices are used for complex diagnostic and interventional procedures and are designed for torque control.

The following methodology is used to identify the Roadrunner® Wire Guides. The product catalogue number prefixes for the various Roadrunner® Wire Guides identifies the device as a Roadrunner® (Polymer Coated) Wire Guide (RPC, RLPC, RFPC and RFSPC) and also identifies the physical attribute of the wire guide. For example the initial Roadrunner® Wire Guides, RPC were manufactured with a standard taper and a stiff shaft. Changes in the device design created the RLPC Roadrunner® Wire Guides manufactured with a Long taper and a stiff shaft. In like manner, changes to the device created the RFPC Roadrunner® Wire Guides with a standard taper and a Flexible shaft and the RFSPC Roadrunner® Wire Guides manufactured with a Flexible shaft and a Short taper.

Table. Methodology for Roadrunner® Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Distal Curve Configuration (mm)	Additional Identifier
RPC	35	260	2	CONNORS
RLPC	35	260		
RFPC	35	145	0	
RFSPC	35	145		

Since most Roadrunner® Wire Guides are manufactured with an angled distal tip, distal curve configurations indicated by "0" represents straight coils and distal curve configurations indicated by "2" represents a 2 millimeter curve radius in the coil tip. The additional identifier CONNORS indicated that the device is a Connors modification and is manufactured with a distal polymer jacket and a proximal section manufactured with PTFE coating.

Roadrunner® Extra Support Wire Guides

Roadrunner® Extra Support Wire Guides (RSTF) is constructed with a heavy duty PTFE (Dupont Compound Number 850-321) coated nitinol mandril wire. The mandril wire portion of the device is designed with a centerless ground mandril wire. The distal spring coils attached to the distal end of the mandril wire are manufactured from platinum and are designed with an angled tip. The flexible tip portion of the wire guide measures between 2.5 to 4.0 centimeters, depending upon device specification. These wire guides are manufactured with a 0.014 or a 0.018 inch (0.036-0.046 cm) outer diameter and are manufactured in length measurements from 180 to 300 centimeters. All available length measurements are relevant to the outer diameter of the wire guide. A lubricious silicone coating is applied to the outer surface of the nitinol mandril wire to facilitate smooth manipulations.

Each device is supplied with a pin vise assembly (PPV-200) for imparting torque response to the distal tip.

The following methodology is used to identify the Roadrunner® Extra Support Wire Guides. The product catalogue number prefix RSTF identifies the device as the Roadrunner® Extra Support Wire Guide manufactured with a Torque Flex characteristic.

Table. Methodology for Roadrunner® Extra Support Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)
RSTF	14	180
RSTF	14	300
RSTF	18	180
RSTF	18	270
RSTF	18	300

Roadrunner® UniGlide™ Hydrophilic Wire Guide

The Roadrunner® UniGlide™ Hydrophilic Wire Guide is constructed of a nitinol core wire with a black radiopaque polyurethane jacket. The outer surface of the polyurethane jacket is hydrophilic

coated for a lubricious outer surface. These wire guides are manufactured with 0.018, 0.025, 0.035 and 0.038 inch (0.046, 0.064, 0.089, 0.097 cm) diameters and in lengths of 80, 150, 180 and 320 centimeters. The construction of the distal tip of the core wire guide is tapered for flexibility and can be manufactured with a straight or an angled tip configuration.

All finished diameters except the 0.018 inch (0.046 cm) wires are manufactured in both a standard shaft version and a stiff shaft version, where the stiff shaft version has a larger diameter core wire than the respective standard shaft version. The smaller 0.018 inch (0.046 cm) diameter wires are manufactured only with the standard shaft version. All Roadrunner® UniGlide™ Hydrophilic Wire Guides with a standard shaft version are designed with an 11 centimeter tapered distal section, with an additional 3 centimeters constant diameter floppy distal tip, extending from the edge of the tapered section. The Roadrunner® UniGlide™ Hydrophilic Wire Guides manufactured with a stiff core are designed in the same fashion however; these wire guides are manufactured with a 14 centimeter tapered distal section followed by the 3 centimeter diameter floppy distal tip.

A torque device is supplied with the Roadrunner® UniGlide™ Hydrophilic Wire Guides. The torque device is intended for use in complex diagnostic and interventional procedures and is designed for torque control.

The following methodology is used to identify the Roadrunner® UniGlide™ Hydrophilic Wire Guides. The product catalogue number prefixes for the Roadrunner® UniGlide™ Hydrophilic Wire Guides identifies the device by the product label description. The letters HPW represents the device as a Hydrophilic Coated Wire Guide. In addition the Roadrunner® UniGlide™ Hydrophilic Wire Guides may be manufactured with an Angled distal tip (HPWA), or with as a straight wire with a Stiff core (HPWS), or with an Angled tip and a Stiff core (HPWAS) construction.

Table. Methodology for Roadrunner® UniGlide™ Hydrophilic Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Shaft Version	Distal Tip
HPW	18	80	Standard	Straight
HPW	25	150	Standard	Straight
HPWA	35	180	Standard	Angled
HPWA	38	150	Standard	Angled
HPWS	35	150	Stiff	Straight
HPWAS	38	80	Stiff	Angled

Approach® Hydro ST Micro Wire Guides

Approach® Hydro ST Micro Wire Guides (HMW) is constructed using a stainless steel core wire guides coated with green (PTFE) polytetrafluoroethylene. The wire guide is constructed with the core wire tapering toward the distal tip where a stainless steel coil and a platinum coil are affixed using silver solder. A blue polyurethane jacket with hydrophilic coating covers the coiled distal tip of the wire guide.

The distal tip of the Approach® Hydro ST Micro Wire Guide consists of a stainless steel coil (distal) measuring nine centimeters and a platinum coil (proximal) measuring three centimeters. The two coils are joined together by a single-point joint and the coil assembly is secured to the core wire in three locations with silver solder. A thirteen centimeter segment of polyurethane jacket surrounds all

coils and joints. The Approach® Hydro ST Micro Wire Guides are manufactured with a 0.014 inch (0.036 cm) diameter and lengths of 135, 190 and 300 centimeters. The Approach® Hydro ST Micro Wire Guides are also manufactured as a straight wire, but are designed with a stainless steel core that allows shapeability of the distal tip if desired. The distal tip is designed with a flex length of 12 centimeters.

The following methodology is used to identify the Approach® Hydro ST Micro Wire Guides. The product catalogue number prefix HMW for the Approach® Hydro ST Micro Wire Guides identifies the device by product type as the Hydro Micro Wire.

Table. Methodology for Approach® Hydro ST Micro Wire Guides

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Additional Identifier
HMW	14	135	ST
HMW	14	190	ST
HMW	14	300	ST

Additional identifiers used in the product catalogue number for the Approach® Hydro ST Micro Wire Guides; ST identifies both the device as the Approach® Hydro ST Micro Wire Guide and indicates that the micro wire is a straight guide with a Shapeable Tip.

HiWire™ Hydrophilic Wire Guide (HW,HWA, HWAS, HWS) is a device that is constructed from a steerable, metallic core covered with a polymer coating. The wire guide's radiopaque jacket is coated with a hydrophilic coating, increasing the lubricity of the device.

The HiWire Hydrophilic Wire Guide is a wire guide that may be used in general intravascular interventional procedures to ease the placement of other devices. The HiWire™ Hydrophilic Wire Guide is available diameters of 0.018, 0.025, 0.035 and 0.038 inches (0.046, 0.064, 0.089, 0.097 cm) and in lengths of 80, 150, 180, and 260 centimeters. The HiWire™ Hydrophilic Wire Guides are manufactured with an angled or a straight distal tip. In addition, these wire guides are manufactured with a standard shaft or one constructed with a stiffer shaft. The HiWire™ Hydrophilic Wire Guides are designed with a 1 to 1 torque response providing maximum control.

Methodology for HiWire Hydrophilic Wire Guide

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Device Construction
HW	18	80	Standard Shaft with Straight Tip Configuration
HWA	25	150	Standard Shaft with Angled Tip Configuration
HWAS	38	260	Stiff Shaft with Angled Tip Configuration
HWS	35	180	Stiff Shaft with Straight Tip Configuration

Disposable Reuter Tip Deflecting Wire Guides (DTDW) are manufactured with stainless steel wire guides that have a distal dip configuration when deflected. These devices are manufactured with .035 inch (0.089 cm) diameters and lengths of 110 and 145 centimeters. The distal curve when formed is generally manufactured with a ten millimeter curve radius.

The Disposable Reuter Tip Deflecting Wire Guides are manufactured to be 5.5 centimeters longer than the appropriate catheter length so that a flush fitting may be used with the handle assembly and allow for appropriate wire guide length to deflect the distal tip of the catheter. Disposable Reuter Tip Deflecting Wire Guides are manufactured with a plastic handle assembly. With the Disposable Reuter Tip Deflecting Wire Guide, the handle assembly and flushing adapter are provided, as the wire and handle assembly are sold as a unit.

Methodology for Disposable Reuter Tip Deflecting Wire Guides (DTDW)

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Distal Curve Radius (mm)
DTDW	35	110	10
DTDW	35	145	10
DTDW	38	110	10

Reuter Tip Deflecting Wire Guides (TDW) are manufactured the same as the disposable model, however, the handle assembly is not provided with the device, but is sold separately. These devices are constructed stainless steel. The Reuter tip Deflecting Wire Guides are manufactured to have a distal curve configuration when deflected. The distal curve when formed is generally manufactured as either a five millimeter or a ten millimeter curve radius. The tip deflecting wire guides are manufactured to be 5.5 centimeters longer than the appropriate catheter length so that a flush fitting may be used with the handle assembly and allow for appropriate wire guide length to deflect the distal tip of the catheter.

Methodology for Reuter Tip Deflecting Wire Guides (TDW)

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Distal Curve Radius (mm)
TDW	25	100	5
TDW	25	125	5
TDW	25	145	10
TDW	35	100	10
TDW	35	100	5
TDW	35	145	10
TDW	35	145	5

Non-Vascular Wire Guides are manufactured with a fixed core. They may also be manufactured with stainless steel or polytetrafluoroethylene (TFE) coated stainless steel, and with a flexible distal tip. These wire guides are designed with a tapered fixed core construction, and with a gradual transition to the distal tip.

The designated categorization of each wire guide type is referenced on the product label, as Lunderquist Exchange Wire Guide, Lunderquist –Ring Torque Wire Guide or Coons Interventional Wire Guide. Complex Fixed Core Guides, such as the Lunderquist Exchange Wire Guide and the Lunderquist-Ring Torque Wire Guide are not manufactured with a safety wire, but contain a mandril

wire which extends from the proximal to the distal tip. In contrast, based upon device design, the Coons Interventional Wire Guide is constructed with both a mandril and safety wire.

Lunderquist-Ring Torque Wire Guide (THG), is constructed of stainless steel or TFE coated stainless steel. This wire guides is manufactured in diameters of 0.035 (0.89 cm) inches, and length measurements of 125 and 145 centimeters. These wire guides have a long tapered mandril wire that extends the length of the device, with anchor solders at specified points at the distal portion of the wire guide. Both ends of the wire guide are welded.

Methodologyfor Lunderquist Ring Torque Wire Guide (THSF)

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Additional Information
THSF	35	125	THG
THSF	35	145	THG

Coons Interventional Wire Guide (THSF) is characterized with having a heavy duty shaft, with a gradual transition to a very flexible distal tip. This wire guide type is manufactured in diameters of .035 and .038 inches (0.089 and 0.97 cm), and length measurements of 80, 100, 125, 145, 180, and 260 centimeters. This wire guide is constructed of a TFE coated stainless steel, tapered mandril wire and a braided safety-wire. The safety-wire is soldered to the mandril wire at a specified distance to allow for a floppy distal tip. Both ends of the wire guide are welded.

Methodologyfor Coons Interventional Wire Guide (THSF)

Product Catalogue Number Prefix	Outer Diameter (X .001")	Length (cm)	Additional Information
THSF	35	100	COONS-SGH
THSF	35	125	COONS-SGH
THSF	35	145	COONS
THSF	35	145	COONS-SGH
THSF	35	180	COONS
THSF	35	260	COONS
THSF	35	80	COONS
THSF	38	145	COONS
THSF	38	145	COONS-SGH
THSF	38	80	COONS

Lunderquist Extra Stiff Wire Guide (TSMG, TSCMG) are comprised of the following products:

1. Straight (TSMG/-LES1) and Curved (TSCMG/-LES1): these are PTFE-coated stainless steel wire guides with a 3 cm, 5 cm or 7 cm flexible tip. The length of the wire guides is 90 cm. The J-curve radius for the curved wire guides is 3 mm.
2. Straight (TSMG/-LES) and Curved (TSCMG/-LES): these are PTFE-coated stainless steel wire guides with a 4 cm or 7 cm flexible tip. For the wire guide lengths 260 cm and 300 cm, the flexible tip includes an inner gold coil for enhanced visibility. The length of the wire guide ranges from 90 cm to 300 cm. The J-curve radius for the curved wire guides is either 3 mm or 7.5 mm.
3. Double Curved (TSCMG/-LECDC and TSCMG/-E-LESDC): these are PTFE-coated stainless steel wire guides with a 4 cm flexible tip that includes an inner gold coil for

enhanced visibility. TSCMG-/-LESDC has a primary/secondary curve radius on 75/15 mm, and TSCMG-/-E-LESDC has a primary/secondary curve radius on 55/15 mm.

Amplatz Ultra Stiff Wire Guides (THSF-/-AUS, THSCF-/-AUS) have less shaft stiffness compared to the Lunderquist Extra Stiff Wire Guides. The increased inner diameter of the wire guide coil allows utilization of an ultra-stiff mandril while maintaining tip flexibility. The Amplatz Ultra Stiff Wire Guide is a PTFE-coated stainless steel wire guide available in a straight or J-tip configuration. The suffix -AUS1 indicates a 3 cm flexible tip and -AUS2 a 7 cm flexible tip.

The Amplatz Ultra Stiff Wire Guides are comprised of the following products:

1. Straight (THSF-/-AUS1, THSF-/-AUS2) and Curved (THSCF-/-AUS1, THSCF-/-AUS2): These are PTFE-coated stainless steel wire guides with a 3 cm (-AUS1) or 7 cm flexible tip (-AUS2). The length of the wire guide ranges from 60 cm to 260 cm. The J-curve radius for the curved wire guides is 3 mm.

Vascular Micro Wire Guides (TPMG) feature a uniquely fabricated stainless steel shaft facilitating torque control and steerability (mandril wire guides). The atraumatic coil tip, made of platinum to enhance visualization, is shapeable and steerable. The spring coil portion is AQ® hydrophilic coated for the TPMG-/-CEP-HC, TPMG-/-CIR-HC and TPMG-/-SAG-HC wire guides, which when wet, reduces the wire guide surface coefficient of friction.

The Vascular Micro Wire Guide has a stainless steel shaft and a soft distal spring coil portion with a radiopaque platinum tip. Not all Vascular Micro Wire Guides are hydrophilic coated.

All Vascular Micro Wire Guides with product codes ending on -HC are hydrophilic coated with AQ. AQ is a bio-compatible hydrophilic coating which becomes very lubricious when wet, greatly reducing the coefficient of friction of the surface of the wire guide. AQ® hydrophilic coating is applied to the distal end of the wire guide. This facilitates control of the proximal portion of the wire guide during procedural use.

The Vascular Micro Wire Guides are comprised of the following products:

1. The Vascular Micro Wire Guides (TPMG-/-FLX12, TPMG-/-CIRRUS, TPMG-/-UNIC,) and the hydrophilic coated Vascular Micro Wire Guides (TPMG-/-CEP-HC, TPMG-/-CIR-HC, TPMG-/-SAG-HC): these are made of a PTFE-coated stainless steel shaft and has a coil tip made of platinum. The distal spring coil portion is made of platinum (TPMG-/-FLX12, TPMG-/-CIRRUS) or of both platinum and a PTFE-coated stainless steel coil (TPMG-/-UNIC, TPMG-/-CIR-HC, TPMG-/-CEP-HC, TPMG-/-SAG-HC). The length of the distal spring coil portion ranges from 5 to 39.5 cm and is AQ hydrophilic coated for TPMG-/-CEP-HC, TPMG-/-CIR-HC, TPMG-/-SAG-HC. The Vascular Micro Wire Guides are available in diameters within the range of 0.14 inch, 0.16 inch or 0.21 inch (0.36, 0.41, 0.53 cm). The overall length of the wire guides is 135 cm, 150 cm or 180 cm.

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Cook Incorporated and William Cook Europe ApS manufacture all medical products in accordance with Medical Devices Directive 93/42/EEC.

Cook Incorporated manufacture all medical products in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820, Quality System Regulations.

All of Cook facilities are also certified to ISO 13485 (2016) Quality Management System.

SHELF LIFE

The "use by" date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

All Vascular Wire Guides have a five year expiration date with the exception of the Approach[®] Hydro ST (HMW) Wire Guide, Roadrunner[®] Wire Guides (RPC, RFSPC, RFPC and RLPC), Roadrunner[®] UniGlide[™] Hydrophilic Wires (HPW)

The applicable expiration date for the various Vascular Wire Guides is provided in the following table.

Expiration Date - Shelf Life

Device Description	Device Prefix	Expiration Date – Shelf Life
Roadrunner [®] Wire Guides	RPC RFSPC RFPC RLPC	3 Years
Approach [®] Hydro ST Micro Wire Guide	HMW	3 Years or less (expiration date varies as assigned by the vender)
Roadrunner [®] UniGlide [™] Hydrophilic Wire Guide	HPW HPWA HPWS HPWAS	3 Years
HiWire [™] Hydrophilic Wire Guide	HW HWA HWAS HWS	3 Years
Tip Deflecting Wire Guides	DTDW TDW	5 Years
Non-Vascular Wire Guides	THSF	5 Years
Lunderquist Extra Stiff Wire Guide	TSMG TSCMG	5 Years
Amplatz Ultra Stiff Wire Guide	THS(C)F/- AUS	5 Years
Vascular Micro Wire Guides	TPMG/-CEP- HC TPMG/-CIR- HC TPMG/- SAG-HC	3 Years
All Other Wire Guides		5 Years

WARRANTY

An authorized representative of the manufacturer in the territory of the Russian Federation is the ZAO "SHAG".

Page 51 of 53

Russia, 119002, Moscow, Karmanitskiy per., 9.office 501.
phone: (495) 956-13-09;
Fax: (495) 956-1310,
e-mail: info@schag.ru
Manufacturer's warranty binding on quality products.

The manufacturer guarantees the operability and preservation of the sterility of the devices during the entire shelf life from the date of manufacture indicated on the package), under the storage conditions.

STERILIZATION DETAILS

Products falling under the scope of this Technical data Summary are supplied in sterile packages. Cook utilizes an EtO (Ethylene Oxide) sterilizes devices using a validated ethylene oxide (EtO) method to a sterility assurance level (SAL) of 10⁻⁶. EtO was the sterilization method selected for these devices as it is effective and does not alter the physical properties of the device.

STORAGE RULES

These devices should be stored in a dark, dry, cool place; avoid extended exposure to light, avoid sudden temperature changes.

TRANSPORTATION RULES

Transported by all modes of transport in covered vehicles in accordance with the current rules of transportation.

TERMS OF USE

Endovascular Wire Guides used in hospitals.

DELIVERY SET

Endovascular Wire Guides supplied in 1 piece in individual packaging. Individual packages are packed into corrugated cardboard boxes. The boxes are sealed with a control tape with a call to the side walls or around the perimeter.

Labeling

Endovascular Wire Guides supplied in 1 piece, packaging which should be free of signs of damage or opening. If bends, fractures or other defects are found, the product cannot be used.

The product label contains the following information:

- manufacturer name and address of production;
- product name (trade mark);
- serial number;
- a mark on the presence of the CE certificate;
- symbol of sterilization;

- address of the European representative;
- trademark of the manufacturer;
- date of manufacture;
- dimensions;
- picture;
- batch number;
- symbol "No re-use";
- symbol for caution "Attention, see Instructions for Use";
- symbol "Keep dry";
- symbol "Keep away from sunlight";
- shelf life.

manufacturer name and address	
shelf life	
serial number (LOT)	
batch number	
symbol of sterilization	
symbol for keep away from sunlight	
symbol for keep dry;	
No re-use	
symbol for "Consult Instruction for Use"	
symbol for caution "Attention, see Instructions for Use"	


DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste. There is no environmentally harmful effect as the result of the device disposal in accordance with the relevant national guidelines.

 Hernán J Serrano Vargas
 Regulatory Affairs Specialist
 Cook Incorporated

United States of America)
 State of Indiana) ss:
 County of Monroe)

02 May 2022
 Date

Subscribed and sworn to before me this 02ND day of MAY of 2022.

My commission expires on March 5, 2025

 Elissa Robinson, Notary Public
 Commission # 697487

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