Susan Uí Uiginn

5 MARY STREET
GALWAY CITY H91 NXW0
IRELAND

TEL: +353 91 567550
FAX: +353 91 396340
EMAIL: contact@susanuiuginn.ie
WEB: www.susanuiuginn.ie

NOTARY PUBLIC

FOR THE COUNTIES OF GALWAY AND
ADJOINING COUNTIES
COMMISSIONED FOR LIFE

NOTARIAL CERTIFICATE

BE IT KNOWN, that I, Susan Uí Uiginn, Notary Public, duly authorised, admitted and sworn by the Chief Justice of the Supreme Court, Republic of Ireland and commissioned and practising at 5 Mary Street, Galway, H91 NXW0, Ireland as a Notary Public in the County of Galway, Republic of Ireland and its surrounding counties


DO HEREBY CERTIFY that the attached Instruction for Use for the Medical Device EmboTrap Revascularization Device in 2 variations (diameter 5mm working lengths 21mm and 33mm respectively) was signed before me and in my presence by Niall Fox, Director of Regulatory, of Neuravi Limited (Cerevovus) whose Registered Office is Block 3 Ballybrit Business Park, Galway H91 K5YD, Ireland

The Country to which the document is addressed is Russia.

IN FAITH AND TESTIMONY whereof I have set my hand and affixed my seal of office at 5 Mary Street, Galway H91 NXW0, Ireland, this date the 12th day of
of October 2021.

Susan Uí Uiginn, Notary Public

The notarial act applied to the within document extends only to verification of the identity & legal capacity of the appearer(s) and to the due execution of the document where applicable unless otherwise expressly stated in the English Language

Susan Ui Uiginn, Notary Public
Galway, Ireland. Commissioned for Life

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Block 3, Ballybrit Business Park, Galway, Ireland, Ирландия

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
Instruction for Use for the medical device

EmboTrap Revascularization Device in variations:

- EmboTrap II Revascularization Device (diameter 5 mm, working length 21 mm);
- EmboTrap II Revascularization Device (diameter 5 mm, working length 33 mm);

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Adverse reactions and events

Possible complications include: air embolism, hematoma or hemorrhage at puncture site, infection, distal embolization, vascular spasm, neurological deterioration including stroke and death, thrombosis, dissection or perforation, device malfunction or fracture, vascular occlusion, ischemia, intracranial hemorrhage, pseudoaneurysm formation, or other minor or major complications related to cerebral angiography.

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Method of use (Instructions)

PROCEDURAL INFORMATION:

The Device is delivered endovascularly under fluoroscopic guidance in a similar manner to that of other neurovascular clot-retrieval systems. The use of antiplatelet and anticoagulation agents is recommended at the discretion of the treating physician.

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Preparation & Delivery:

- 1) Using standard interventional techniques, access the arterial system and using angiography, determine the location of the occluded vessel.
- 2) Advance an appropriate Guide Catheter, Sheath or Balloon Guide Catheter as close to the occlusion site as possible. Connect a rotating hemostasis valve (RHV) to the proximal end of this catheter and connect to a continuous flush system.
- 3) Select an appropriate Microcatheter (see section «The possibility and ways of integrating with other medical devices»). Connect an RHV to the proximal end of the Microcatheter and connect to a continuous flush system.
- 4) With the aid of a suitable Guidewire (see section «The possibility and ways of integrating with other medical devices»), and using standard catheterization techniques and fluoroscopic guidance, advance the Microcatheter up to and across the occlusion so that the distal end of the Microcatheter is positioned distal to the occlusion.
- 5) Remove the Guidewire from the Microcatheter and, if desired, gently infuse contrast media through the Microcatheter to visualize the distal end of the occlusion.
- 6) Remove the Insertion Tool with the preloaded Device from the packaging hoop, taking care not to accidentally retract or advance the Device from the Insertion Tool while doing so.
- 7) Insert the distal end of the Insertion Tool through the RHV of the Microcatheter and wait until fluid is seen exiting the proximal end of the Insertion Tool, confirming that the Device is flushed. Advance the Insertion Tool until it is fully seated in the hub of the Microcatheter. Fully tighten the RHV to hold the Insertion Tool securely in position.
- 8) Confirm that the Insertion Tool is fully seated in the hub of the RHV before proceeding to advance the Device. Advance the Device until at least half of the shaft length has been inserted into the Microcatheter, at which point the Insertion Tool may be removed.
- 9) Do not discard the Insertion Tool as it will be needed if additional device passes are required during the procedure.

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- 1) Continue to advance the Device until the distal radiopaque tip of the Device approaches the distal region of the Microcatheter.

Note: There are two silver bands positioned on the proximal portion of the Device shaft to assist in minimizing the amount of fluoroscopic exposure required during Device insertion. For example, if using a standard Microcatheter (total length of 155cm and a 7cm RHV), then when the first silver band on the shaft approaches the RHV, the Device tip will be approximately 8cm from the distal end of the Microcatheter. When the second silver band on the Device shaft approaches the RHV the Device tip will be nearing the distal end of the Microcatheter.

- 2) Advance the Device in the Microcatheter and position the EmboTrap within the clot by aligning the proximal markers of the Device with the proximal margin of the clot (Figure 1).
- 3) If the proximal markers are not visible, position within the clot by advancing the Device until the proximal radiopaque coil is approximately 7mm from the proximal margin of the clot. (Figure 2).


Figure 1


Figure 2

- 4) While stabilizing the Device shaft, retract the microcatheter to deploy the Device within the clot. The Device will be fully deployed once the distal tip of the Microcatheter has been retracted to a position over the proximal radiopaque coil of the Device.
- 5) The Device may be left to embed within the clot for 3-5 minutes prior to withdrawal.

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Retrieval:

- 1) Optimum control of the Device and Microcatheter relative positions during retrieval may be achieved by advancing a torque device over the Device shaft and tightening it against the Microcatheter RHV.
- 2) If using a Balloon Guide Catheter inflate the balloon to occlude the vessel as per the Balloon Guide Catheter IFU.
- 3) Withdraw the Device and Microcatheter slowly and carefully as a unit to the Guide Catheter tip while aspirating through the Guide Catheter with a syringe.
- 4) As the device reaches the Guide Catheter apply vigorous aspiration, withdraw the Device and Microcatheter into the Guide Catheter and continue to aspirate until the Device reaches the proximal RHV.

Note: If withdrawal into the Guide Catheter is difficult (as may be the case with a large clot burden) then deflate the balloon (if applicable) and withdraw the Guide, Microcatheter and Device together through the Introducer Sheath.

- 5) Disconnect the RHV from the Guide Catheter and remove the Device, Microcatheter and RHV together from the Guide Catheter.
- 6) Use a syringe to further aspirate the Guide Catheter to ensure it is clear of any thrombus material.

Thrombus Removal and Preparation For Additional Passes:

- The Device may be used for up to three retrieval attempts.
- If an additional pass is to be made with the Device then carefully remove any captured thrombus from the Device and clean the Device in heparinized saline, rubbing gently from proximal to distal to remove any residual thrombus material. Inspect the Device carefully and if any damage or deformity is observed do not use the Device – instead use a new EmboTrap II Revascularization Device for any subsequent passes. If using the same Device, replace the Insertion Tool onto the proximal shaft and retract the Device into the Insertion Tool until it is fully resheathed.
- Do not attempt more than three retrieval attempts in one vessel.

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

EmboTrap II Revascularization Device (diameter 5 mm, working length 33 mm)

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
Note: The indicator band positions are measured from the distal margin of the indicator band to the distal tip of the device while the device is in the collapsed configuration.

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
LIST AND DESCRIPTION OF MD MATERIALS


Medical device contact type: Short-term contact with circulating blood (<24 hours).




The shaft is made of nitinol and coated with PTFE with a Pebax outer layer. Platinum coils are made of platinum alloy. The outer chamber is made of nitinol and contains gold markers, channel is made of nitinol. Dymax MD adhesive for medical devices is used as an adhesive. The introducer is made of polytetrafluoroethylene.

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	Warnings – see instructions for use
	Use only once – Do not re-use
	See instructions for use
	Sterilized using ethylene oxide
	Catalogue Number
	Minimum Microcatheter inner diameter (ID)
	Keep dry
	Lot Number

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INFORMATION ABOUT MEDICINES AND (OR) PHARMACEUTICAL SUBSTANCES INCLUDED IN MD

Devices do not contain medicines and (or) pharmaceutical substances and are not used in combination with medicines and (or) pharmaceutical substances.

INFORMATION ABOUT THE STERILITY OF A MEDICAL DEVICE

All EMBOTRAP models are supplied sterile and are sterilized by 100% Ethylene Oxide using the same validated cycle and method. Minimum Sterility Assurance Level (SAL) of 10^{-6} .

REQUIREMENTS FOR TECHNICAL SUPPORT AND REPAIR OF MD

Not applicable.

CONDITIONS OF TRANSPORTATION AND STORAGE

Table 3 – Transportation conditions

Designation	Range
Temperature Limitation, °C	-35 to +40
Relative Humidity Limitation, without condensation, %	30 - 95
Atmospheric Pressure Limitation, hPa	700-1060

Table 4 – Storage conditions

Designation	Range
Temperature Limitation, °C	+10 to +35
Relative Humidity Limitation, without condensation, %	<75
Atmospheric Pressure Limitation, hPa	860-1060

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MD DISPOSAL RULES AND CONDITIONS

Disposal should be in accordance with generally the requirements SanPiN 2.1.3684-21 and applicable local, state and federal laws and regulations.

If package is damaged or expiration date is exceeded, device is disposed in conformity with the rules of handling medical wastes of Class A.

The used device must be disposed in conformity with the rules of handling medical wastes of Class B.

WARRANTY

The manufacturer guarantees the conformity of EmboTrap Revascularization Device to the requirements of the standards in compliance with operating, storage and transportation conditions of the expiration date.

Warranty storage period: 3 years

NAME AND LEGAL ADDRESS OF MANUFACTURER AND AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION**MANUFACTURER:**

Neuravi Limited, Block 3, Ballybrit Business Park, Galway, Ireland.

MANUFACTURING SITES:

Neuravi Ltd, Block 3, Ballybrit Business Park, Galway, Ireland.
Advant Medical Ltd, Parkmore Business Park West, Galway, Ireland.

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RUSSIA:

Johnson & Johnson LLC,
121614, Moscow, Krylatskaya str., 17/2.
Tel. +7 (495) 580-7777
RA-JNJRU-RZN@ITS.JNJ.com

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Susan Uí Uiginn

MARY STREET
GALWAY CITY H91 NXW0
IRELAND

TEL: +353 91 567550
FAX: +353 91 396340
EMAIL: contact@susanuiuginn.ie
WEB: www.susanuiuginn.ie

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DO HEREBY CERTIFY that the attached Instruction for Use for the Medical Device EmboTrap Revascularization Device in 3 variations (diameter 5 mm, working lengths 22mm and 37mm respectively and diameter 6.5mm, working length 45mm) was signed before me and in my presence by Niall Fox, Director of Regulatory, of Neuravi Limited (Cerevovus) whose Registered Office is Block 3 Ballybrit Business Park, Galway H91 K5YD, Ireland

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- Patients with angiographic evidence of carotid dissection.

Adverse reactions and events

Possible complications include: air embolism, hematoma or hemorrhage at puncture site, infection, distal embolization, new territory embolization, vascular spasm, neurological deterioration including stroke and death, thrombosis, dissection or perforation, device malfunction or fracture, vascular occlusion, ischemia, intracranial hemorrhage, pseudoaneurysm formation, or other minor or major complications related to cerebral angiography.

INFORMATION ON POTENTIAL USERS OF MD AND DIRECTIONS FOR USE

Directions for use

The device is used in the internal environment of the body.

Table 1 – Environmental conditions

Environmental conditions	
Temperature	+32°C to +42°C
Humidity	100%
Pressure	700-1060 hPa

Potential users of MD

EmboTrap Revascularization Device must be used in compliance of the rules and precautions notified in the Instructions for Use.

The Device should only be used by physicians trained in neurointerventional catheterization and the treatment of ischemic stroke.

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Indication of possibilities and features of medical device use in people with implanted medical devices, pregnant women, breastfeeding women, children, adults with chronic diseases

EmboTrap has not been evaluated for use in people with implanted medical devices, pregnant women, breastfeeding women, children, or adults with chronic diseases.

EmboTrap does not have any special features intended for these specific patient populations. There is limited or no data on the safety or performance of EmboTrap in pregnant women, breastfeeding women, or children.

Risks associated with interactions between EmboTrap and previously stented vessels (existing implanted Med Device) has been considered in the device risk assessments and appropriate warnings / precautions are provided in the labeling.

DESCRIPTION OF MAIN FUNCTIONAL ELEMENTS OF MD AND PRINCIPLE OF OPERATION

General description of functional elements of MD and principle of operation

The Device comprises of a retrievable, self-expanding Nitinol shaped section mounted on the distal end of a PTFE-coated, tapered Nitinol shaft. The Device is supplied pre-loaded within an Insertion Tool.

It is designed to restore blood flow in the neurovasculature by removing thrombus in patients experiencing ischemic stroke. The EmboTrap Revascularization Device is supplied sterile and is intended for single-use only by physicians trained in neuro-interventional catheterization and the treatment of ischemic stroke.

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PROCEDURAL INFORMATION:

The Device is delivered endovascularly under fluoroscopic guidance in a similar manner to that of other neurovascular clot-retrieval systems. Physicians are advised to select devices which are appropriately sized for the estimated occlusion length and location in the vasculature in accordance with the information presented in Table 2.

The use of antiplatelet and anticoagulation agents is recommended at the discretion of the treating physician.

Preparation & Delivery:

1. Using standard interventional techniques, access the arterial system and using angiography, determine the location of the occluded vessel.
2. Advance an appropriate Guide Catheter, Sheath or Balloon Guide Catheter as close to the occlusion site as possible. Connect a rotating hemostasis valve (RHV – see section «The possibility and ways of integrating with other medical devices») to the proximal end of this catheter and connect to a continuous flush system.

Note: Depending on tortuosity and/or occlusion location the use of an Intermediate Catheter may be considered to provide additional support for ease of retrieval.

3. Select an appropriate microcatheter (see section «The possibility and ways of integrating with other medical devices»). Connect an RHV to the proximal end of the microcatheter and connect to a continuous flush system.
4. With the aid of a suitable guidewire (see section «The possibility and ways of integrating with other medical devices»), and using standard catheterization techniques and fluoroscopic guidance, advance the microcatheter up to and across the occlusion so that the distal end of the microcatheter is positioned distal to the occlusion.
5. Remove the guidewire from the microcatheter and, if desired, gently infuse contrast media through the microcatheter to visualize the distal end of the occlusion.

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6. Remove the Insertion Tool with the preloaded Device from the packaging hoop, taking care not to accidentally retract or advance the Device from the Insertion Tool while doing so.
7. Insert the distal end of the Insertion Tool through the RHV of the microcatheter and wait until fluid is seen exiting the proximal end of the Insertion Tool, confirming that the Device is flushed. Advance the Insertion Tool until it is fully seated in the hub of the microcatheter. Fully tighten the RHV to hold the Insertion Tool securely in position.
8. Confirm that the Insertion Tool is fully seated in the hub of the RHV before proceeding to advance the Device. Advance the Device until at least half of the shaft length has been inserted into the microcatheter, at which point the Insertion Tool may be removed.
9. Do not discard the Insertion Tool as it will be needed if additional device passes are required during the procedure.

Positioning and Deployment:

1. Continue to advance the Device until the distal radiopaque tip of the Device approaches the distal region of the microcatheter. Care should be taken not to advance the distal tip of the Device out of the microcatheter.

Note: There are two silver bands positioned on the proximal portion of the Device shaft to assist in minimizing the amount of fluoroscopic exposure required during Device insertion. For example, if using a standard microcatheter (working length of 155cm and a 7cm RHV), then when the first silver band on the shaft approaches the RHV, the Device tip will be approximately 5cm to 8cm from the distal end of the microcatheter (when using the longest or shortest models, respectively). When the second silver band on the Device shaft approaches the RHV the Device tip will be nearing the distal end of the microcatheter.

2. Advance the Device in the microcatheter and position the EM-BOTRAP III such that the clot lies within the working length of the device, i.e. between the most proximal and most distal radiopaque markers.

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3. Where possible, and if appropriate, position the Device proximally such that the proximal radiopaque markers align with the proximal margin of the clot (Figure 3).


Figure 3

4. While stabilizing the Device shaft, retract the microcatheter to deploy the Device within the clot. The Device will be fully deployed once the distal tip of the microcatheter has been retracted to a position over the proximal radiopaque coil of the Device.
5. The Device may be left to embed within the clot for 3-5 minutes prior to withdrawal.

Retrieval:

1. Optimum control of the Device and microcatheter relative positions during retrieval may be achieved by advancing a torque device over the Device shaft and tightening it against the microcatheter RHV.
2. If using a balloon guide catheter inflate the balloon to occlude the vessel as per the balloon guide catheter IFU.
3. Withdraw the Device and microcatheter slowly and carefully as a unit to the guide catheter tip while aspirating through the guide catheter with a syringe.
4. As the device reaches the guide catheter apply vigorous aspiration, withdraw the Device and microcatheter into the guide catheter and continue to aspirate until the Device reaches the proximal RHV.

Note: If withdrawal into the guide catheter is difficult (as may be the case with a large clot burden) then deflate the balloon (if applicable) and withdraw the guide, microcatheter and Device together through the introducer sheath.

5. Disconnect the RHV from the guide catheter and remove the Device, microcatheter and RHV together from the guide catheter
6. Use a syringe to further aspirate the guide catheter to ensure it is clear of any thrombus material.

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• The Device may be used for up to three retrieval attempts.

• If an additional pass is to be made with the Device then carefully remove any captured thrombus from the Device using heparinized saline, rubbing gently from proximal to distal until all residual thrombus material is removed. Inspect the Device carefully and if any damage or deformity is observed do not use the Device – instead use a new EMBOTRAP III Revascularization Device for any subsequent passes. If using the same Device, replace the Insertion Tool onto the proximal shaft and retract the Device into the Insertion Tool until it is fully resheathed. Do not retract the device into the Insertion Tool against significant resistance –ascertain the cause of the resistance and if required use a new device for any subsequent passes.

• Do not attempt more than three retrieval attempts in one vessel.

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The possibility and ways of integrating with other medical devices

Note: other medical devices of any manufacturer that meet the conditions described in this section can be used.

Microcatheter: The Device should be introduced through a microcatheter with an internal diameter (ID) of 0.021"- 0.027".


Balloon Guide Catheter / Guide Catheter: The Device should be used in conjunction with a Neurovascular Guide or Sheath. The guide/sheath chosen is to be compatible with the selected microcatheter.

Additional Accessories: The following devices may be selected for use during the procedure. Only accessory devices which are compatible with the selected microcatheter and guide catheter combination should be used.

- Intermediate Catheter / Aspiration Catheter
- Rotating Hemostasis Valve (RHV)


Note: RHV chosen should be appropriate for use with the EMBOTRAP III Insertion Tool, i.e. the insertion tool is required to seat fully into the microcatheter hub and be secured firmly by the RHV. See *Preparation and Delivery* for further information.

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EmboTrap III Revascularization Device (diameter 5 mm, working length 37 mm)

PTFE-Coated Nitinol Shaft w/PEBAJ Jacket

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
EmboTrap III Revascularization Device (diameter 6,5 mm, working length 45 mm)

PTFE-Coated Nitinol Shaft w/PEBAX Jacket

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


7. DEFINITIONS OF SYMBOLS USED ON THE LABEL

	Warnings – see instructions for use
	Use only once – Do not re-use
	See instructions for use
	Sterilized using ethylene oxide
	Catalogue Number
	Minimum Microcatheter inner diameter (ID)
	Keep dry
	Lot Number
	Use by

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Service life – N/A

HOW SUPPLIED

EmboTrap Revascularization Device:

№	Name	Qty*
	<i>Consumer package contains:</i>	
1	Instruction for Use Booklet / IFU	1
2	Individual package (Tyvek Pouch):	1
2.1	Packaging Hoop:	1
2.1.1	Device (EmboTrap Revascularization Device)	1

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- Do not withdraw the Device against significant resistance. Assess the cause of resistance using fluoroscopy and if necessary, advance the microcatheter over the Device to resheath or partially resheath to aid withdrawal. Resheathing a Device with captured clot may result in loss of clot and distal embolization.
- This Device is for use in one patient only. Do not resterilize or reuse in more than one patient. Reusing the Device in another patient could result in compromised Device performance, cross infection and other safety related hazards.
- Do not use any EMBOTRAP III Revascularization Device more than 3 times.
- Do not perform more than 3 retrieval attempts in any one vessel.
- Do not withdraw the deployed Device through a previously stented vessel. If it becomes necessary to cross a deployed stent to access a clot, first ensure that the stent can be crossed with the Guide or Intermediate Catheter. The Device should then be fully pulled back into Guide or Intermediate Catheter prior to retrieval through the stent.
- Users shall take all necessary precautions to limit X-radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible. The risks of angiographic and fluoroscopic X-ray radiation include alopecia, burns ranging in severity from skin reddening to ulcers, cataracts and delayed neoplasia that increase in probability as procedure time and number of procedures increase.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Serious Incident Reporting:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- a) the death of a patient, user or other person,
- b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- c) a serious public health threat.

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