NOTARY CERTIFICATION BY DOCUMENT CUSTODIAN

I, Ashni De Alwis, of legal age, being first duly sworn, hereby
Name of document custodian

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A h h.

Signature of document custodian (affiant)

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
State of California

County of Orange

Subscribed and sworn to (or affirmed) before me on this 31st day of Oct.
20 17 by Ashni De Alwis, proved
to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Steve Oberlander
Signature (Seal)

INSTRUCTIONS FOR USE

Marksman™ Catheter

h h h

70968-001 Rev. 12/16

This page is intentionally left blank.

TABLE OF CONTENTS

Marksman™ Catheter

English	4
Français	5
Deutsch	6
Italiano	7
Español	8
Svenska	9
Nederlands	10
Português	11
Suomi	12
Dansk	13
Ελληνικά	14
Česky	15
Magyar	16
Русский	17
Polski	18
Türkçe	19
Norsk	20
Slovenčina	21
Română	22
Български	23
한국어	24
عربي	25
Symbol Glossary	26

English
Instructions for Use

Marksman™ Catheter

CAUTION

Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician. This device should be used only by physicians with a thorough understanding of angiography and/or percutaneous neurointerventional procedures. Inspect the catheter prior to use for any kinks or bends. Any catheter damage may decrease the desired performance characteristics.

DESCRIPTION

The Marksman™ Catheter is a variable stiffness, single lumen catheter designed to access small, tortuous vascular areas. The outer surface of the catheter's distal segment is coated with a hydrophilic material to provide lubricity during use. The catheter also incorporates a PTFE liner to facilitate movement of introduction devices passed through its lumen. The Marksman™ Catheter has a radiopaque marker at the distal tip to facilitate fluoroscopic visualization. The distal tip of the catheter is shapeable.

Table with 4 columns: DIMENSIONS, FRENCH, INCH, mm. Rows include Proximal O.D., Distal O.D., Proximal I.D., and Distal I.D.

CONTENTS

- 1 ea Marksman™ Catheter
- 1 ea Shaping Mandrel

INDICATIONS FOR USE

This device is intended for the introduction of interventional devices and infusion of diagnostic or therapeutic agents into the neuro, peripheral, and coronary vasculature.

COMPATIBILITY

Table with 2 columns: DEVICE DESCRIPTION, SIZE. Rows include Guiding Catheter and Guidewire.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to, the following:

- Distal emboli (air, foreign material, tissue, or thrombus)
- Hemorrhage
- Ischemia
- Neurological deficit that could lead to stroke and death
- Hematoma
- Infection
- Pain at insertion site
- Vessel spasm, dissection, perforation or injury

WARNINGS

THIS DEVICE IS INTENDED FOR SINGLE USE ONLY. Discard after one procedure. Structural integrity and/or function may be impaired after reuse or cleaning. The catheter should be manipulated under fluoroscopy only. Do not attempt to move the catheter without observing the resultant tip response. Never advance or withdraw an intraluminal device against resistance until the cause of resistance is determined by fluoroscopy. If the cause cannot be determined, withdraw the catheter. Movement of the micro catheter against resistance may result in damage to the micro catheter, or the vessel. Static pressure for this catheter should not exceed 400 psi. This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

EN

DIRECTIONS FOR USE

- Before removing the catheter from the packaging hoop, flush the hoop with heparinized saline through the luer fitting connected to the end the packaging hoop. Repeat if friction is felt when attempting to remove the catheter.
- Flush the lumen of the catheter with heparinized saline.
- Gently remove the catheter from packaging hoop and inspect to verify that it is undamaged.
- Shaping mandrel:
 - Remove shaping mandrel from card and insert into the distal tip of the catheter.
 - Bend the shaping mandrel to the desired shape.
 - Hold shaping mandrel and catheter tip directly over a steam source for approximately 30 seconds to set shape. Remove from heat source and allow cooling in either air or liquid prior to removing the mandrel.
- Carefully insert an appropriate guidewire into the hub and advance into the catheter lumen.
- Place an appropriate guiding catheter using standard technique. Attach rotating hemostasis valve to the guiding catheter luer connector and maintain a continuous flush.
- Introduce the catheter and wire assembly into the guiding catheter through the hemostasis sidearm adaptor. Tighten the valve around the catheter to prevent backflow but allowing movement of the catheter through the valve.
- Advance the guidewire and catheter to the selected target site by alternately advancing the guidewire and then tracking the catheter over the guidewire.
- Withdraw the guidewire from the catheter. Connect the desired syringe with infusate or interventional device and continue the procedure according to the manufacturer's instructions.

Flow Rate Table and Approximate Values table. Flow Rate Table has columns for 100% Contrast, 50/50 Contrast, and 100% Saline. Approximate Values table has columns for Catheter Working Length, Dead Space, and flow rates at 100 and 300 psi.

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Table with 4 columns: Description, French, Inches, mm. Rows include Proximal outer diameter, Distal outer diameter, Proximal inner diameter, and Distal inner diameter.

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Table with 2 columns: Description of device, Size. Rows include Catheter-guide and Guide wire.

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Table with 4 main columns: Flow rate, 100% Contrast, 50/50 Contrast, 100% Physiological solution. Sub-table with approximate values for working length, dead space, and pressure at different flow rates.

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COPY CERTIFICATION BY DOCUMENT CUSTODIAN

I, PRIYA NAVKAR, of legal age, being first duly sworn, hereby
Name of document custodian

Swear (or affirm) that the attached reproduction of

Markoman Appendix to IFU
Description of documents

Is a true and exact copy of the correct and complete original document.

Priya Navkar
Signature of document custodian (affiant)

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.


State of California

County of Orange

Subscribed and sworn to (or affirmed) before me on this 28th day of Nov.
20 17 by PRIYA NAVKAR — proved

to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Steve Oberlander
Signature (Seal)

Appendix to the Instruction for use

Name of Medical Device

1. Name of the medical device	1. Marksman catheter in various versions:
	1. Marksman catheter, versions: 1.1. The Marksman catheter, 105 cm long, consisting of: - Marksman catheter, length 105 cm - Forming mandrels 1.2. The Marksman catheter, 135 cm long, consisting of: - Marksman catheter, length 135 cm - Forming mandrels 1.3. The Marksman catheter, 150 cm long, consisting of: - Marksman catheter, length 150 cm - Forming mandrels 1.4. The Marksman catheter, 160 cm long, consisting of: - Marksman catheter, length 160 cm - Forming mandrels
2.1 With regards to the manufacturer of the medical device - full name and, if available, abbreviated name, including company name, legal form of the legal entity,	Developer: Micro Therapeutics, Inc. DbA Ev3 Neurovascular 9775 Toledo Way, Irvine, CA 92618, USA +1-949-837-3700 Manufacturer: Micro Therapeutics, Inc. DbA Ev3 Neurovascular

Available upon request

Priya Markan

Appendix to the Instruction for use
Name of Medical Device

location address or last name, first name and (if available) patronymic, details of the identity document, address of residence of the individual entrepreneur, as well as phone numbers and (if available) e-mail address of the legal entity or individual entrepreneur;	9775 Toledo Way, Irvine, CA 92618, USA +1-949-837-3700 Holder of Registration Certificate in Russia (name, address, tel, fax, e-mail, website) Micro Therapeutics, Inc. DbA Ev3 Neurovascular 9775 Toledo Way, Irvine, CA 92618, USA +1-949-837-3700 Addresses of Manufacturing Sites for medical device: 1) Micro Therapeutics, Inc., DBA ev3 Neurovascular, 9775 Toledo Way, Irvine, CA 92618, USA. 2) Medtronic Mexico S.de R.L.de CV, Av. Paseo Cucapah, 10510 El Lago, C.P. 22210, Tijuana, Baja California, MEXICO, Baja California, MEXICO
2.2 With respect to authorized representative of the manufacturer of the medical device - full and, if available, abbreviated name, including company name, legal form of the legal entity, address (location), telephone numbers and (if available) e-mail address of the	Authorized Representative of Manufacturer: "Medtronic" LLC, 123317, Russia, Moscow, Presnenskaya Naberezhnaya., 10 Tel: (495) 580-73-77, Fax: (495) 580-73-78 www.medtronic.ru

Available upon request

Appendix to the Instruction for use
Name of Medical Device

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Available upon request

Appendix to the Instruction for use
Name of Medical Device

7. Description of accessories, medical devices or non-medical devices that are intended for use in combination with the medical device (if applicable)	Формирующий мандрен Forming mandrel		Длина, мм: 76 Диаметр ушка, мм: 6 Толщина стержня, мм: 0,7 Масса, г: 0,09
8. Information on the presence in the medical device of a medicinal product for human use, human and/or animal materials	Not applicable		
9. Information on the procedure of installation, adjustment, calibration and other actions necessary to put the medical device into operation	Not applicable		
10. Requirements for the premises in which the installation of the medical device shall be carried out, as well as	Not applicable		

Available upon request

Appendix to the Instruction for use
Name of Medical Device

requirements for training or qualification of persons who install the medical device (if applicable)	
11. Information necessary to verify the correct installation of the medical device and its readiness for safe operation, including:	Not applicable
a. Maintenance and its frequency, including cleaning and disinfection of the medical device;	Not applicable
b. List of information provided by the manufacturer of the medical device, keys, access passwords, programs required for installation, adjustment, operation and maintenance of the medical device;	Not applicable

Available upon request

Appendix to the Instruction for use
Name of Medical Device

c. List of consumables (components, reagents), as well as the procedure for their use and replacement;	Not applicable
d. Need for calibration to ensure the proper and safe operation of the medical device during its lifetime;	Not applicable
e. Risk minimization methods associated with the installation, calibration or maintenance of the medical device;	Not applicable
f. Information on installation, adjustment, calibration and other actions necessary for putting the medical device into operation and its correct operation	Not applicable

Available upon request

Appendix to the Instruction for use
Name of Medical Device

(application);																	
g. Information on the list of the main characteristics for the operation (application) of the medical device, transportation and storage conditions (for example, temperature and humidity, lighting and other characteristics)	The medical device is intended for use under standard conditions. A medical device is allowed to be used only by specialists who have passed possession and have relevant work experience. Operating conditions in the health facility: <table border="1" data-bbox="539 448 1195 635"> <tr> <td>t°</td><td>10 to 40 C°</td></tr> <tr> <td>Pressure</td><td>700-1060 gPa</td></tr> <tr> <td>Humidity</td><td>30 -75 %</td></tr> </table> Operating conditions in the internal environment of the body: <table border="1" data-bbox="539 759 1195 959"> <tr> <td>t°</td><td>32 to 42 C°</td></tr> <tr> <td>Pressure</td><td>700-1060 gPa</td></tr> <tr> <td>Humidity</td><td>100 %</td></tr> </table> Storage Conditions Store with other sterile devices in normal operating room temperatures in a dry place and in a manner that protects the integrity of the package and the sterile barrier. Rotate inventory so that device is used before the expiry date on the package. <table border="1" data-bbox="539 1206 1195 1335"> <tr> <td>Air temperature, °C</td><td>+10 - +32</td></tr> <tr> <td>Air humidity (relative), non-</td><td>< 75</td></tr> </table>	t°	10 to 40 C°	Pressure	700-1060 gPa	Humidity	30 -75 %	t°	32 to 42 C°	Pressure	700-1060 gPa	Humidity	100 %	Air temperature, °C	+10 - +32	Air humidity (relative), non-	< 75
t°	10 to 40 C°																
Pressure	700-1060 gPa																
Humidity	30 -75 %																
t°	32 to 42 C°																
Pressure	700-1060 gPa																
Humidity	100 %																
Air temperature, °C	+10 - +32																
Air humidity (relative), non-	< 75																

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	condensing, %	
	Atmosphere pressure, hPa	860 - 1060
	Transport Conditions	
	Air temperature, °C	-35 +40
	Air humidity (relative), non-condensing, %	30 - 95
	Atmosphere pressure, hPa	700 - 1060
h. List of national standards applied by the manufacturer	International Standards:	
	Number	Standard Title
	ISO 594-1:1986E	Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment -- Part 1: General requirements
	ASTM F 2052-06	Standard Test Method for Measurement of Magnetically Induced Displacement Force on Passive Implants in the Magnetic Resonance Environment
	ASTM F 2096-04	Standard Test Method for Detecting Gross Leaks in Porous Medical Packaging by Internal Pressurization (Bubble Test)
	ASTM F 2182-02a	Standard Test Method for Measurement of Radio Frequency Induced Heating Near Passive Implants During MR Imaging
	ASTM F 2119-07	Standard Test Method for Evaluation of MR Imaging Artifacts from Passive Implants
	ASTM F 2503-05	Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	ISO 15223-1:2012	Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	EN 1041:2008	Information supplied by the manufacturer of medical devices
	EN ISO 14971:2012	Medical devices – Application of risk management to medical devices
	EN ISO 10993-1 :2009	Biological evaluation of medical devices – Part 1: Evaluation and testing
	EN ISO 10993-3:2003	Biological Evaluation of Medical Devices Part 3: Tests for Genotoxicity, Carcinogenicity and Reproductive Toxicity
	EN ISO 10993-4:2002; A1:2006	Biological Evaluation of Medical Devices-Part 4: Hemocompatibility
	EN ISO 10993-5:2009	Biological Evaluation of Medical Devices -Part 5: Tests for In Vitro Cytotoxicity
	EN ISO 10993-6:2007	Biological Evaluation of Medical Devices -Part 6: Implantation
	EN ISO 10993-7:2008/Cor. 1:2009	Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals
	EN ISO 10993-10:2002; A1:2006	Biological Evaluation of Medical Devices -Part 10: Tests for Irritation and Delayed- type Hypersensitivity

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	EN ISO 10993-11:2006	Biological Evaluation of Medical Devices -Part 11: Test for Systemic Toxicity
	EN ISO 10993-12:2007	Biological Evaluation of Medical Devices -Part 12: Sample Preparation and Reference Materials
	ISO 10993-13:2010	Biological evaluation of medical devices – Part 13: Identification and quantification of degradation products from polymeric medical devices
	ISO 10993-17:2002	Biological evaluation of medical devices – Part 17: Establishment of allowable limits for leachable substances
	ISO 10993-18:2005	Biological evaluation of medical devices – Part 18: Chemical characterization of materials
	ISO 15223-1:2012	Medical devices - Symbols to be used with Medical Device Labels, Labeling and Information to be Supplied –Part 1: General Requirements
12. Information on the sterile condition of the medical device, method of its sterilization and the procedure for dealing with a sterile package violation (if the medical device is delivered sterile), or if the medical device is not sterile, indication of the	Microcatheters are sterilized using medical equipment, using the standard ethylene oxide sterilization cycle (EtO), which results in a guaranteed level of sterility (SAL) of 10-6. ANSI / AAMI / ISO 11135-1: 2007 and ANSI / AAMI / ISO 11135-2: 2008. The microcatheters are sterilized according to the parameters listed in the table below:	
	Sterility Validation Method.	ANSI/AAMI/ISO 11135-1 2007, Sterilization of health care products – Ethylene oxide – Requirements for development, validation and routine control of a sterilization process for medical devices. Overkill approach.
	EO Residuals	ISO 10993-7 2008 Biological evaluation of medical devices

Available upon request

Appendix to the Instruction for use
Name of Medical Device

need for sterilization before use		– Part 7: Ethylene oxide sterilization residuals.	
	Sterilization Method	100% ethylene oxide gas	
	Sterility Assurance	10 ⁻⁶	
	Pyrogen Tests	The device is non pyrogenic. Pyrogen testing is conducted on a lot-to-lot basis using the “Guideline for Validation of LAL Test as an End-Product Endotoxin Test for Medical Devices, (FDA 1978).”	
	Contract Sterilization	Sterigenics, Inc., Los Angeles, CA	
13. Information on the procedure for the preparation of the medical device for its reuse, including cleaning, disinfection, packaging and, if necessary, method	Not applicable		

Available upon request

Appendix to the Instruction for use
Name of Medical Device

of re-sterilization (if the medical device is reusable), as well as criteria for the unfitness of the medical device for use	
14. Information necessary for the identification of medical devices in order to obtain a safe combination, and information on known restrictions on concomitant use of medical devices (for medical devices intended for use with other medical devices and/or supplies)	Not applicable
15. Information on the nature, type, and (if necessary) the intensity and distribution of radiation (electromagnetic, ionizing, other) of the medical device and methods for	Not applicable.

Available upon request

Appendix to the Instruction for use
Name of Medical Device

the protection of consumers and third parties from unintentional radiation during the operation of the medical device (if the medical device creates a dangerous or potentially hazardous radiation level when used as intended)	
16. Information on precautions taken in case of:	WARNINGS • Infusion pressure with this device should not exceed 600 PSI (4137 kPa). Pressures in excess of 600 PSI (4137 kPa) may result in catheter rupture, which may result in patient injury. • If flow through the micro catheter becomes restricted, do not attempt to clear the device by high pressure infusion. Either remove the micro catheter to determine the cause of the obstruction or replace it with a new micro catheter. Excessive pressure may cause catheter rupture, which may result in patient injury. • Never advance or withdraw an intraluminal device against resistance until the cause of the resistance is determined by fluoroscopy. Excessive force against resistance may result in damage to the device or vessel perforation or other patient injuries. • This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance. Warnings Shaping Mandrel is not intended for use in the human body. Use only a steam source to shape the catheter tip. Do not use other heat sources. Prior to use, inspect the catheter tip for any damage that may have resulted from shaping. Do not use a catheter that has been damaged in any way. Damaged catheters may rupture causing vessel trauma or tip detachment during steering maneuvers.

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	PRECAUTIONS • Ensure embolic material compatibility with catheter prior to use. • Prior to use, carefully examine the micro catheter and packaging to verify that they have not been damaged during shipment. • Prior to use, all accessory devices and agents should be fully prepared according to the manufacturer's instructions. • Always monitor infusion rates when using the micro catheter. • When injecting contrast for angiography, ensure that the micro catheter is not kinked or occluded. • The micro catheter has a lubricious hydrophilic coating on the outside of the catheter. It must be kept hydrated in order to be lubricious. This can be accomplished by attaching the Y-connector to a continuous saline drip. • When the infusion catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the micro catheter without observing the resultant tip response. • Because the micro catheter may be easily advanced into narrow, selective vasculature, repeatedly verify (for example by slightly withdrawing the catheter) that the micro catheter has not been advanced so far as to interfere with its removal.
a. Failures of the medical device, malfunction or deviations in the functioning that may affect the safety of the medical device, including those determined by external signs;	Transient Use – medical device temporary contacts human body. Limited duration of contact – less than 24 hours (during complicated surgery with several stages) Device is for single use only Device is intended for use on central circulatory system.
b. Effects on the operation of the medical device of external factors associated with the use of the medical device in combination with	Not applicable

Available upon request

Appendix to the Instruction for use
Name of Medical Device

other medical devices and/or equipment, or such predictable factors as external electromagnetic fields, electrostatic discharges, radiation (electromagnetic, ionizing, other), atmospheric pressure and its differences, humidity and air temperature;	
c. The risk of electromagnetic interference created by the medical device for other medical devices, equipment and communication equipment during the conduct and evaluation of the results of diagnosis, treatment or its intended use (for example, electromagnetic	Not applicable

Available upon request

Appendix to the Instruction for use
Name of Medical Device

radiation from a medical device that affects other equipment);	
17. Warnings and/or precautions taken by the consumer in the application of the medical device containing a medicinal product for human use, human and/or animal material, materials that are carcinogenic, mutagenic or toxic, whose possible release or washing-out leads to sensitization, allergic reaction or negatively affects the reproductive function	Not applicable. The medical device does not contain materials of animal or human.
18. Warnings and/or precautions taken by the consumer when disposing of the medical device,	For disposing of the microcatheters is required to provide the mandatory implementation of all the necessary procedures in accordance with local applicable laws and regulations.

Available upon request

Appendix to the Instruction for use
Name of Medical Device

accessories and consumables used together with it (if any), including information about the infectious, microbial, environmental or physical hazard of the medical device	
19. Information on the circumstances which should prompt the consumer to consult a health professional	A medical device used only by qualified medical personnel who have undergone specialized training and have an appropriate level of surgical experience.
20. Information on the initial release or the latest revision of the operational documentation;	Not applicable
21. Procedure and conditions for the disposal or destruction of a medical device.	<u>After use, this product may be a potential biohazard.</u> Dispose the medical device in accordance with the requirements of the national legislation on environment and health and/or hospital Unused, expired devices are classified as medical waste class «A» Catheters extracted from the patient's body are medical waste Class B, that is epidemiologically hazardous waste.

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	Waste Class B subjects to mandatory disinfection (disinfection)/ disposal. The choice of the method of disinfection/ neutralization is determined by the capabilities an organization engaged in medical and/ or pharmaceutical activity and performed in the development of the treatment scheme of medical waste.
--	---

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