

This is to certify that the attached Instruction For Use For Medical Device For In Vitro Diagnostics: BD FACSPresto cartridge with reagents for direct calculation of the absolute and relative number of CD4 lymphocytes by fluorescence photomicroscopy in samples of human venous and capillary whole blood is a true and accurate copy in the possession of Becton, Dickinson and Company, BD Biosciences, 2350 Qume Drive, San Jose, California, 95131, USA

Cristina Renteria
Sr. International Affairs Specialist

February 21, 2022
Date

I, Cristina Renteria, hereby affirm that the attached reproduction of Instruction For Use For Medical Device For In Vitro Diagnostics: BD FACSPresto cartridge with reagents for direct calculation of the absolute and relative number of CD4 lymphocytes by fluorescence photomicroscopy in samples of human venous and capillary whole blood is a true, correct and complete copy of the document in the possession of Becton, Dickinson and Company, BD Biosciences.

Cristina Renteria
Becton, Dickinson and Company, BD Biosciences
2350 Qume Dr., San Jose, CA 95131


A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California
County of Santa Clara

Subscribed and sworn to (or affirmed) before me on this 21 day of February 2022, by Cristina Renteria, proved to me on the basis of satisfactory evidence to be the person who appeared before me.

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INSTRUCTION FOR USE
FOR MEDICAL DEVICE FOR IN VITRO DIAGNOSTICS:
**BD FACSPresto cartridge with reagents for direct calculation of the
absolute and relative number of CD4 lymphocytes by fluorescence
photomicroscopy in samples of human venous and capillary whole
blood**

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
Becton Dickinson and Company,

BD Biosciences

Senior International Regulatory Affairs Specialist

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Developer and legal manufacturer:

Becton, Dickinson and Company, BD Biosciences
2350 Qume Drive, San Jose, California, 95131, USA

Manufacturing site:

Becton Dickinson and Company, BD Biosciences
2350 Qume Drive, San Jose, California, 95131, USA

Authorized representative of the manufacturer in Russia:

Becton Dickinson Vostok Limited Liability Company
24/27 Sadovaya-Samotyochhnaya Street, 2nd floor,
premises I, office 18, 127051 Moscow, the Russian
Federation.
Tel./Fax: +7 (495) 7758582, +7 (495) 77585833

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CD4 counts and CD4 percentages (%CD4) have been used to evaluate the immune status of patients diagnosed with, or suspected of developing, immune deficiencies. The CD4 T cell percentage of total lymphocytes shows less variation in all ages group.^{31,33}

The human immunodeficiency virus specifically targets and binds the receptors on CD4 T cells.^{29,31,34} A central feature in the pathogenesis of HIV is the relentless destruction of CD4 T cells, either directly or indirectly, resulting in the depletion of CD4-bearing T lymphocytes in peripheral blood.³⁵ This loss of CD4 T cells results in failure of an HIV-specific immune response, recall antibody response, and finally, non-specific immune response in the AIDS stage.²⁹

The amount of CD4 T cells are counted in peripheral blood of HIV infected individuals to measure the degree of immune deficiency in order to assess disease progression, to guide therapeutic decisions by assessing eligibility for antiretroviral treatment (ART) and antimicrobial prophylaxis, and immune restoration in patients receiving ART.^{29,31,33-36}

In 2017, the WHO redefined terms now referred as advanced HIV infection based on the use of viral load.³⁸

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CD4 T cell percentages of lymphocytes and absolute CD4 T cell counts provide the single effect prognostic factors for estimating the risk of AIDS and death in HIV infected persons³³

and provides the best direct measurement of a patient's immune status and their risk for opportunistic infections^{30,31,36,37} and AIDS-defining neoplasms³¹

The absolute number and percentage of CD4 T lymphocytes are the cellular parameters most closely associated with HIV disease progression and patient prognosis.⁴

CD4 T cell counting is also useful for applications unrelated to HIV, such as for follow-up and diagnostic evaluation of primary cellular immune deficiencies, including severe combined immunodeficiency and immune reconstitution after hematopoietic transplantation.³⁵

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6. Principles of the procedure

The BD FACSPresto cartridge¹ contains dried fluorochrome-conjugated antibody reagents. When blood reacts with the reagents, the antibodies in the reagents bind to the surface antigens on the lymphocytes and monocytes. After the incubation period, the cells are analyzed on BD FACSPresto instrument. The software, which is embedded into BD FACSPresto instrument, identifies the cell populations of interest and calculates CD4 absolute counts, CD4 percentages of lymphocytes.

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7. Reagent composition

The cartridge contains CD4 PE -Cy5/CD3 APC /CD45RA APC/CD14 PE dried antibody reagents. The dried antibody reagents include inert ingredients such as buffer, bovine serum albumin (BSA), and ProClin as a preservative. CD3 APC and CD45RA APC enumerate total lymphocytes. Lymphocytes labeled

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with CD4 PE-Cy5 are designated CD4⁺lymphocytes. CD14 PE identifies monocytes, which are excluded from the analysis.

The CD4 antigen^{11,12} (55 kDa)¹³ is present on the helper/inducer Tlymphocyte subset^{14,15} CD3⁺CD4⁺, which consists of 28% to 58%¹⁶ of lymphocytes in normal peripheral blood,¹² and is present in low density on the surface and in the cytoplasm of monocytes. The CD4 antibody, clone SK3,¹¹ is derived from the hybridization of mouse myeloma cells with spleen cells from BALB/c mice immunized with human peripheral blood T lymphocytes.

The CD4 antibody is composed of mouse IgG1 heavy chains and kappa light chains. The CD3 antibody, clone SK7, is derived from the hybridization of NS-1 mouse myeloma cells with spleen cells from BALB/c mice immunized with human thymocytes. CD3 is composed of mouse IgG1 heavy chains and kappa light chains.

The CD45RA antigen is present on approximately 50% of CD4⁺ T lymphocytes, approximately 75% of CD8⁺ T lymphocytes, and on essentially all B lymphocytes and natural killer (NK) lymphocytes.¹⁷ The helper/inducer Tlymphocyte subset expresses the phenotype CD4⁺CD45RA⁺.¹⁷ The CD45RA antigen is expressed on naïve T lymphocytes. Antigen density decreases upon *in vitro* activation.¹⁸ A selective loss of the CD4⁺CD45RA⁺ subset during active multiple sclerosis has been demonstrated.^{17,19}

The CD45RA antibody, clone HI100,²⁰ is derived from the hybridization of mouse myeloma cells with spleen cells isolated from mice immunized with human whole blood cells (WBCs).

The CD14 antigen is present on the majority of normal peripheral blood monocytes.²¹ The CD14 antibody recognizes a human monocyte/macrophage antigen of 55-kDa.²² The CD14 antibody, clone MφP9, is derived from the hybridization of Sp2/0 mouse myeloma cells with spleen cells from BALB/c mice immunized with peripheral blood monocytes from a patient with rheumatoid arthritis. The CD14 antibody is composed of mouse IgG2b heavy chains and kappa light chains.

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8. Package contents

Each box contains 100 cartridges and 100 pipets as well as one IFU.

9. Storage and Handling

Store the cartridge:

- In its original foil pouch. Do not use the cartridge if the pouch has been opened for more than 30 minutes.
- At 4°C–31°C (39°F–88°F).
- In 10%–95% non-condensing humidity.
- Until the expiration date. Do not use the cartridge after the expiration date on the package.
- Incubate the cartridge at 10°C–40°C.

10. Cartridge QC

Cartridge Quality Control (QC) uses immobilized antibodies. The instrument verifies that the reagent is present and that there is sufficient sample in the cartridge.

Cartridge QC runs automatically at every cartridge run.

11. Patient specimens

The assay is designed to be used only with peripheral whole blood collected by venipuncture into EDTA tubes or by finger stick.

Capillary²³ or venous blood samples are transferred directly into the cartridge and incubated. Samples are run on the instrument after incubation.

Follow these guidelines for handling your samples:

- Do not dilute whole blood before adding it to the cartridge.

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- Do not refrigerate the whole blood specimen before sample preparation.
- Store whole blood collected in EDTA tubes at 20°C–25°C (68°F–77°F) up to 24 hours before applying the blood to the cartridge.
- Minimize exposure of the cartridges to light.
- Do not remove the channel protector on the cartridge until just before you insert the cartridge into the instrument.

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12. Reference intervals

The reference intervals for the BD FACSPresto cartridge shown in Table 1 are representative for hematologically normal adults.

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Table 1. Representative reference intervals.

Subset	Gender	N	Mean	Reference interval
CD4 ²⁷	Male	77	811	462–1,306 cells/ μ L
	Female	83	866	440–1,602 cells/ μ L
%CD4 ²⁷	Male	77	41	29%–54%
	Female	83	44	32%–55%

We recommend that laboratories and other users establish their own reference intervals for their patient populations using the BD FACSPresto system to reflect potential sources of variability, such as patient gender, race, age, and preparation techniques.

13. Expected results. Performance characteristics of BD FACSPresto system

Performance of the BD FACSPresto cartridge (the cartridge) was established by the testing at the BD Biosciences laboratories in San Jose, CA, USA and at one clinical laboratory in Kisumu, Kenya, Africa.

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Accuracy Evaluation

Absolute counts of CD4-positive cells, the percentage of CD4 positive cells in the lymphocyte population in whole blood from HIV-infected patients were determined using the BD FACSPresto system. Results were compared with results from the BD Tritest CD3 FITC/CD4 PE/CD45 PerCP reagent, in BD Trucount tubes, on the BD FACSCalibur flow cytometer, with the BD FACS Loader, using BD Multiset software.

Whole blood samples were collected at the clinical laboratories. Regression statistics reported in Table 2 and Table 3 indicate the results are substantially equivalent.

Details are shown in Figure 1 through Figure 4.

Table 2. Method comparison for venous blood

Parameter	N	R ²	Slope	Intercept	Range
CD4	189	0.98	0.97	7.37	55-2478 cells/ μ L
%CD4	188	0.96	1.03	0.13	5.06-53.77%

Table 3. Method comparison for capillary blood

Parameter	N	R ²	Slope	Intercept	Range
CD4	162	0.97	1.03	13.47	69-2474 cells/ μ L
%CD4	161	0.96	1.02	-0.26	5.9-56.2%

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/ Figure 1. Scatter plot with weighted Deming fit ($y = 0.97x + 7.37$) for venous blood (CD4)

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/ Figure 2. Scatter plot with Deming fit ($y = 1.03x + 0.13$) for venous blood (%CD4)

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Precision Evaluation

Twenty replicates were created by one operator in one day, using two concentration levels of process controls, one lot of cartridges, and one instrument. Results are shown in Table 4.

Estimates of precision were determined at one site, BD Biosciences, using CD4 cellular controls. Two replicates of each CD4 control (normal and low) were analyzed in each run, and two runs were performed per day for a total of 21 days. Three different instruments and three cartridge lots with three different operators were used, each for seven of the 21 days. Coefficients of variation (CVs) and standard deviations (SDs) are provided for CD4 absolute counts and CD4 percentages for within-run precision in Table 5 and Table 6, respectively.

Table 4. Repeatability study

	Level	Mean	CV	SD ^a	N
CD4	Low	155.10	5.78		20
	Normal	926.60	2.59		20
%CD4	Low	12.77		0.73	20
	Normal	43.99	1.53		20

a. Standard deviations (SDs) are reported instead of coefficients of variation (CVs) for %CD4 low level.

Table 5. Within-run and within-device precision CD4 absolute counts

	Low control (CV)	Normal control (CV)
Within-run	6.79	2.18
Within-device	6.79	3.30

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Table 6. Within-run and within-device precision CD4 percentages

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Linearity

Linearity was assessed in triplicate measurements of multiple concentrations of CD4+ cells across the reportable range of the assay for CD4 absolute counts, total lymphocytes, on the instrument.

Results are linear in the CD4 range (50–4,000 cells/μL) and absolute lymphocyte range (200–10,000 cells/μL).

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Analytical Sensitivity

The limit of blank (LoB) for CD4 absolute count was determined to be 10 CD4 cells/μL and the limit of detection (LoD) was determined to be 22 CD4 cells/uL.

Therefore, analytical sensitivity of the BD FACSPresto system can be defined by a LoD of 22 CD4+cells/uL.

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CD4 Absolute Count Limit of Quantitation Determination

The limit of quantitation (LoQ) of the CD4 assay for the BD FACSPresto system is determined to be 35 CD4 cells per uL.

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14. Interfering Conditions

Table 8 lists the substances that were tested for interference with the reagents in the cartridge. Testing for interference was performed in accordance with EP7²⁸.

There was no detectable interference at the following concentrations.

Table 8. Interfering substances

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During working with patient blood samples must be observed the following limitations:

WARNING Do not use previously fixed and stored samples. Whole blood specimens refrigerated before staining can give incorrect results. Specimens from patients taking immunosuppressive drugs can yield poor resolution.²⁴ Do not test hemolyzed specimens.

WARNING The reagents contain antibodies of mouse and rat origin.

WARNING All biological specimens and materials coming in contact with them are considered biohazards. Handle as if capable of transmitting infection^{25,26} and dispose of with proper precautions in accordance with federal, state, and local regulations. Never pipette by mouth. Wear suitable protective clothing, eyewear, and gloves.

WARNING Each cartridge is for single use only. Use one cartridge per specimen.

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18. Procedure

Before you begin

- Open the cartridge package. Write the patient ID on the cartridge. Face the inlet port up. You must use the cartridge within 30 minutes of opening the package.

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

To prepare a capillary blood sample:

1. Prepare for puncture by increasing blood flow to the fingers.

Ask the patient to extend his or her arm towards the floor and shake and squeeze his or her hand. If the patient's hand is cold, use a warming device or run the hand under warm water.


2. Firmly squeeze the base of the patient's fingers.
3. Select either the middle or ring finger and clean the fingertip with an alcohol. Let the alcohol dry.
4. Place the tip of the lancet on the circled area of the finger, as shown in the following figure. Push the lancet to puncture the finger.

5. Wipe away the first drop of blood with the sponge.
The initial drop might contain tissue fluids that may dilute the sample.
6. Wait until a large drop of blood forms.
7. Apply the blood to the inlet port of the cartridge.

Hold the cartridge by its ridges only. Squeeze and hold the finger to maintain blood flow until blood reaches the top of the inlet port.

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
8. Make sure the cartridge is level, with the barcode side up, at all times. Make sure that blood appears in the part of the channel not covered by the channel protector, next to the containment zone.

9. Apply the sponge and a bandage to the patient's finger to stop excess bleeding.
10. Close the cartridge cap securely.
11. Make sure both latches on the cap snap closed.

Caution: Biological!

If necessary, use a cloth dampened with bleach diluted to 0.5% sodium hypochlorite concentration to clean excess blood outside the containment zone. Be careful to not contaminate the inlet port or introduce any debris into the underside of the cartridge. Do not smear, contaminate, or damage the barcode. If you drop the cartridge into a contaminated area, discard the cartridge into a biohazardous waste container and start over with a new cartridge.


12. Set the on-board timer.
13. Make sure the fill indicator is full.
14. Place the cartridge, barcode side up, on the workstation.

To analyze samples, follow the BD FACSPresto Instrument Instructions for Use

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To prepare a venipuncture sample:

1. Invert the tube 10 times to mix the contents well.
2. Use the pipet to obtain the specimen.
Use one pipet per specimen.
3. Gently squeeze the bulb on the pipet to form a drop of blood on the tip of the pipet.
4. Carefully dispense the specimen into the inlet port.

Hold the cartridge by its ridges only. Make sure the blood reaches the top of the inlet port. If necessary, gently squeeze the bulb on the pipet to dispense more blood. Make sure the cartridge is level, with the barcode side up, at all times. Make sure that blood appears in the part of the channel not covered by the channel protector, next to the containment zone.

5. Discard the pipet into a biohazardous waste container.
6. Close the cartridge cap securely. Make sure both latches on the cap snap closed.

Caution: Biological! If necessary, use a cloth dampened with bleach diluted to 0.5% sodium hypochlorite concentration to clean excess blood outside the containment zone. Be careful to not contaminate the inlet port or introduce any debris into the underside of the cartridge. Do not smear, contaminate, or damage the barcode. If you drop the cartridge into a contaminated area, discard the cartridge into a biohazardous waste container and start over with a new cartridge.

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7. Set the on-board timer.

See Setting the on-board timer (page 26) in the IFU for "Portable cell analyzer BD FACSPresto for counting CD4 lymphocytes by fluorescence photomicroscopy in samples of human whole blood (BD FACSPresto™ Near-Patient CD4 Counter), with accessories".

8. Make sure the full indicator is full.

9. Place the cartridge, barcode side up, on the work.

To analyze samples, follow the BD FACSPresto Instrument Instructions for Use

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19. Technical description

BD Presto cartridge

The BD FACSPresto cartridge is a unit-dose disposable device that contains all the reagents that are needed to test one sample or specimen. Use one cartridge per specimen.

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Disposable pipette

Disposable Pipette are used for collection samples from tube by medical personal.

There is not contact with human body. Non-Sterile. For each samples must be used a new pipette.

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Population and demographic aspects of the use

Device applicable to HIV-infected male and female subjects over 2 years and normal subjects free of hematological abnormalities.

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Potential users

Potential users of the in vitro medical device can be health counselor, medical laboratory technician, medical laboratory assistant, nurse are capable of operating BD FACSPresto instrument and possess all the required knowledge on collecting blood samples and handling biohazards.

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Indications for use

An indication for the use of a medical device is the need to determinate of absolute count and percentage of CD4 lymphocytes for a patient.

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Utilization rules




Separate collection of waste at the end of life as required by local rules and laws. Dispose of in accordance with the applicable country regulation.

User of medical device is liable for recycling of all wastes (sample materials), accumulated in course of measurements, in compliance with legislative and local norms, as well as with protocols of the medical institution.

It is recommended to recycle the system in compliance with Public health regulations SanPiN 2.1.3684-21.

At the end of the shelf life, not opened pouch of the cartridge and pipettes can be disposed of as Class B waste.

Cartridges, patient samples and wastes from the various stages of sample handling, as well as all components of a medical device that have come into contact with biological samples, must be disposed of using disinfectants in accordance with national requirements for the safe handling of pathogens and parasitic pathogens and class B medical waste management.

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1. Giorgi JV, Hultin LE. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. *Clinical Immunology Newsletter*. 1990;10:55-61.
2. Schmidt RE. Monoclonal antibodies for diagnosis of immunodeficiencies. *Blut*. 1989;59:200-206.
3. Dalgleish AG, Beverley PC, Clapham PR, Crawford DH, Greaves MF, Weiss RA. The CD4 (T4) antigen is an essential component of the receptor for the AIDS retrovirus. *Nature*. 1984;312:763-767.
4. Fahey JL, Taylor JM, Detels R, et al. The prognostic value of cellular and serologic markers in infection with human immunodeficiency virus type 1. *N Engl J Med*. 1990;322:166-172.
5. Lewis DE, Puck JM, Babcock GF, Rich RR. Disproportionate expansion of a minor T cell subset in patients with lymphadenopathy syndrome and acquired immunodeficiency syndrome. *J Infect Dis*. 1985;151:555-559.
6. Ohno T, Kanoh T, Suzuki T, et al. Comparative analysis of lymphocyte phenotypes between carriers of human immunodeficiency virus (HIV) and adult patients with primary immunodeficiency using twocolor immunofluorescence flow cytometry. *Tohoku J Exp Med*. 1988;154:157-172.
7. Stites DP, Casavant CH, McHugh TM, et al. Flow cytometric analysis of lymphocyte phenotypes in AIDS using monoclonal antibodies and simultaneous dual immunofluorescence. *Clin Immunol Immunopathol*. 1986;38:161-177.
8. Henry DH, Beall GN, Benson CA, et al. Recombinant human erythropoietin in the treatment of anemia associated with human immunodeficiency virus (HIV) infection and zidovudine therapy. Overview of four clinical trials. *Ann Intern Med*. 1992;117:739-748.
9. Firnhaber C, Smeaton L, Saukila N, et al. Comparisons of anemia, thrombocytopenia, and neutropenia at initiation of HIV antiretroviral therapy in Africa, Asia, and the Americas. *Int J Infect Dis*. 2010; 14:e1088-1092.
10. Levine AM, Berhane K, Masri-Lavine L, et al. Prevalence and correlates of anemia in a large cohort of HIV-infected women: Women's Interagency HIV Study. *J Acquir Immune Defic Syndr*. 2001;26:28-35.
11. Bernard A, Boumsell L, Hill C. Joint report of the first international workshop on human leucocyte differentiation antigens by the investigators of the participating laboratories. In: Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman SF, eds. *Leucocyte Typing*. New York, NY: Springer-Verlag; 1984:9-108.
12. Evans RL, Wall DW, Platsoucas CD, et al. Thymusdependent membrane antigens in man: inhibition of cell-mediated lympholysis by monoclonal antibodies to TH2 antigen. *Proc Natl Acad Sci USA*. 1981;78:544-548.
13. Ledbetter JA, Evans RL, Lipinski M, Cunningham-Rundles C, Good RA, Herzenberg LA. Evolutionary conservation of surface molecules that distinguish T-lymphocyte helper/inducer and cytotoxic/suppressor subpopulations in mouse and man. *J Exp Med*. 1981;153:310-323.
14. Engleman EG, Benike CJ, Glickman E, Evans RL. Antibodies to membrane structures that distinguish suppressor/cytotoxic and helper T lymphocyte subpopulations block the mixed leukocyte reaction in man. *J Exp Med*. 1981;154:193-198.
15. Kotzin BL, Benike CJ, Engleman EG. Induction of immunoglobulin-secreting cells in the allogeneic mixed leukocyte reaction: regulation by helper and suppressor lymphocyte subsets in man. *J Immunol*. 1981;127:931-935.
16. Reichert T, DeBruyère M, Deneys V, et al. Lymphocyte subset reference ranges in adult Caucasians. *Clin Immunol Immunopath*. 1991;60:190-208.
17. Sobel RA, Hafler DA, Castro EA, Morimoto C, Weiner HL. The 2H4 (CD45R) antigen is selectively decreased in multiple sclerosis lesions. *J Immunol*. 1988;140:2210-2214.
18. Serra HM, Krowka JF, Ledbetter JA, Pilarski LM. Loss of CD45R (Lp220) represents a post-thymic T cell differentiation event. *J Immunol*. 1988;140:1435-1441.
19. Rose LM, Ginsberg AH, Rothstein TL, Ledbetter JA, Clark EA. Selective loss of a subset of T helper cells in active multiple sclerosis. *Proc Natl Acad Sci USA*. 1985;82:7389-7393.
20. Schwinzer R. Cluster report: CD45/CD45R. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:628-634.
21. Bernstein ID, Self S. Joint report of the myeloid section of the Second International Workshop on Human Leukocyte Differentiation Antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human Myeloid and Hematopoietic Cells*. Vol 3. New York, NY: Springer-Verlag; 1986:1-25.

22. Goyert SM, Ferrero E. Biochemical analysis of myeloid antigens and cDNA expression of gp55 (CD14). In: McMichael AJ, Beverley PC, Ceballos S, et al, eds. *Leucocyte Typing III: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1987:613-619.
23. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2008. CLSI document GP42-A6.
24. Giorgi JV. Lymphocyte subset measurements: significance in clinical medicine. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:236-246.
25. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI document M29-A3.
26. Siegel JD, Rhinehart E, Jackson M, Chiarello L, Healthcare Infection Control Practices Advisory Committee. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. <http://www.cdc.gov/ncidod/dhqp/pdf/isolation2007.pdf>.
27. Zeh C, Amornkul PN, Inzaule S, et al. Populationbased biochemistry, immunologic and hematological reference values for adolescents and young adults in a rural population in western Kenya. *PLOS One*. 2011;6:1-10.
28. McEnroe RJ, Burritt MF, Powers DM, et al. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI document EP7-A2.
29. World Health Organization WHO. Laboratory guidelines for enumerating CD4 T lymphocytes in the context of HIV/AIDS. In: Regional Office for South-East Asia; 2009.
30. Stein DS, Korvick JA, Vermund SH. CD4+ lymphocyte cell enumeration for prediction of clinical course of human immunodeficiency virus disease: a review. *J Infect Dis* 1992;165:352-363.
31. Turner BJ, Hecht FM, Ismail RB. CD4+ T-lymphocyte measures in the treatment of individuals infected with human immunodeficiency virus type 1. A review for clinical practitioners. *Arch Intern Med*. 1994;154:1561-1573.
32. Wilfert CM, et al. Quality standard for the enumeration of CD4+ lymphocytes in infants and children exposed to or infected with human immunodeficiency virus. *Clin Infect Dis*. 1995;21:S134-135.
33. O'Gorman MR, Zijenah LS. CD4 T cell measurements in the management of antiretroviral therapy -- A review with an emphasis on pediatric HIV-infected patients. *Cytometry B Clin Cytom*. 2008;74:S19-26.
34. Barnett D, Walker B, Landay A, Denny TN. CD4 immunophenotyping in HIV infection. *Nat Rev Microbiol*. 2008;6:S7-15.
35. Kestens L, Mandy F. Thirty-five years of CD4 T-cell counting in HIV infection: From flow cytometry in the lab to point-of-care testing in the field. *Cytometry B Clin Cytom*. 2017;92:437-444.
36. Ford N, Meintjes G, Vitoria M, Greene G, Chiller T. The evolving role of CD4 cell counts in HIV care. *Curr Opin HIV AIDS*. 2017;12:123-128.
37. Achhra AC, Petoumenos K, Law MG. Relationship between CD4 cell count and serious long-term complications among HIV-positive individuals. *Curr Opin HIV AIDS*. 2014;9:63-71.
38. World Health Organization WHO. Guidelines for managing advanced HIV and rapid initiation of antiretroviral therapy. 2017.

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BD, логотип BD и все остальные товарные знаки являются собственностью компании Becton, Dickinson and Company. © BD, 2016