### Certification for Information in Dossier

The submitted dossier for "Tubes for Dual flow pump for arthroscopy" for the Russian application contains the following:

- Instructions for use of **Day set for Dual Flow Pump**
- Instructions for use of **Patient set with check valve for dual flow pump**
- Instructions for use of **Intermediary tubing with check valve for Dual Flow Pump**

I certify that the documents which are submitted are true and accurate to the best of my knowledge. In addition, HEMODIA SAS takes responsibility for any legal obligations related to this submission.

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
## ENGLISH

### INSTRUCTIONS FOR USE

**INDICATIONS FOR USE:** The device for daily use delivers the irrigation fluid from the bag to the intermediary irrigation tube. This device is used in conjunction with the Intermediary irrigation tube or Patient Set (Ref. 11010078 or Ref. 11010077) and the Dual Flow Pump (Ref. 11010072).

**DESCRIPTION:** Daily irrigation tube DAY SET for Dual Flow Pump.

**SUPPLY DETAILS:** Sterile, if the package is unopened and undamaged. Re-usable for the surgical day when used with Ref. 11010078 or Ref. 11010077. See warnings below and appropriate Dual Flow Pump user manual (Document 104072).

#### INSTRUCTIONS FOR USE:

1. Open using aseptic technique.
2. Open the outer package.
3. Close all three clamps and connect the spikes to the irrigation fluid bags.
4. Slip the green connector from the tubing into the green bracket.
5. Place and center the tubing around the roller head.
6. Pull the tube tight.
7. Place the pressure chamber (vertically) in holder, close the pressure arm against the tube and place the transducer of the tubing onto the CPC™ connector of the pump.
8. Open one of the clamps connected to the irrigation fluid bags.
9. Close the left cover ("INFLOW") on the pump.
10. Press on the "Run/Stop" button of the Dual Flow Pump. Blinking green light will turn to solid green light.
11. Press the "Fill Chamber" button, several times, until chamber reaches 1/3 full.

**WARNING:** The pressure chamber must remain vertical during the procedure.

**WARNING:** The DAY SET tube can be re-used for multiple procedures during the same surgical day only when used in conjunction with Ref. 11010078 or Ref. 11010077.

**WARNING:** Sterility is preserved only if proper setup and tube changing procedures have been followed according to the Dual Flow Pump instruction manual.

#### WARNING:

1. Do not use if the package is damaged or opened.
2. Refer to and read the Dual Flow Pump instruction manual (Document 104072) for proper use.
3. Re-sterilization:

This product has not been designed to be re-sterilized. Reprocessing may lead to changes in material characteristics such as deformation and material degradation which may impact the strength of the device and compromise device performance. These risks may potentially affect patient safety.

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# INSTRUCTIONS FOR USE

**PATIENT SET**  
**WITH CHECK VALVE**  
**FOR DUAL FLOW PUMP**  
**REF. 11010077**

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







**INDICATIONS FOR USE:** The device is intended to be used in conjunction with the Dual Flow Pump (Ref. 11010072) during arthroscopic procedure and is composed of an intermediary tube (Ref. 11010078) and a suction tubing. The intermediary tube delivers the irrigation fluid from the DAY SET tube (Ref. 11010076) to the arthroscope. The suction tubing is used to remove the fluid in the joint via the cannula or a shaver.

**DESCRIPTION:** PATIENT SET made up of intermediary & suction tube sets for Dual Flow Pump.

**SUPPLY DETAILS:** Sterile, single use. Sterile if the package is unopened and undamaged.

**INSTRUCTIONS FOR USE:**

1. Open using aseptic technique.
2. Open the outer package.
3. The scrub nurse removes the both tube sets from the package.
4. Remove the smaller tube set (INTERMEDIARY) from the small package.
5. Ensure that all Luer locks are tight.
6. The scrub nurse hands the tubing end with the blue cap to the circulating nurse.
7. Remove the blue cap from the end of the INTERMEDIARY tubing (Fig I. D) and the DAY SET tubing (Fig I. C).
8. Immediately connect the INTERMEDIARY tubing (Fig II. D) and the DAY SET tubing (Fig II. C). Tighten the nut to lock the both Luer connectors.
9. Press on the "Run/Stop" button of the Dual Flow Pump. Blinking green light will turn to solid green light.  
Open the clamp under the fill chamber to prime tubing.
10. Once the tubing has been primed, close off clamp.
11. The scrub nurse takes the large tube set (suction tubing) and hands the tubing end marked "WASTE" to the circulating nurse.
12. Insert the "WASTE" end into the waste collection device.
13. The circulating nurse slips the orange connector of the tubing into the orange notch.
14. Place and center the tubing around the pump head.
15. Pull the tube tight, place it into the housing provided for this purpose and close the pressure arm against the tube.
16. Insert the right tubing "OUTFLOW CANNULA" into the right side of the pinch valve. Press the "Pinch Valve" button of the Dual Flow Pump, then insert the tubing marked "SHAVER" into the left side of the pinch valve.
17. Pass the tubes, after the pinch valve, into the hook provided for this purpose.
18. Close the right cover ("OUTFLOW") on the pump.
19. The scrub nurse plugs the tubing end marked "OUTFLOW CANNULA" onto the cannula. Plug the tubing end marked "SHAVER" onto the shaver device.

**END of the operation:**

Close the clamp under the pressure chamber.

Disconnect the white Luer lock (Fig III. E) under the check valve and discard tubing. The check valve will remain connected to the DAY SET tubing until the start of the next operation.

**START of the next operation:**

1. The circulating nurse removes the old check valve by disconnecting the used check valve (FIG IV. D) from the DAY SET tubing (FIG IV. C).
2. Immediately connect the blue Luer lock (Fig. V. D) of the new INTERMEDIARY tubing to the blue connector (FIG V. C) of the DAY SET tubing as per step 7 in the Instructions for use.

**WARNING:**

1. Do not use if the package is damaged or opened.
2. Refer to and read the Dual Flow Pump instruction manual (Document 104072) for proper use.
3. This product is for single-use only. It has not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as deformation and material degradation which may impact the strength of the device and compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
4. Never attach one check valve to another.

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


**INDICATIONS FOR USE:** The device is intended to be used in conjunction with the Dual Flow Pump (Ref. 11010072) and delivers the irrigation fluid from the DAY SET tube (Ref. 11010076) to the arthroscope during arthroscopic procedures.

**DESCRIPTION:** INTERMEDIARY tube set for Dual Flow Pump.

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**INSTRUCTIONS FOR USE:**

1. Open using aseptic technique.
2. Open the outer package.
3. The scrub nurse removes the tube set from the package.
4. Ensure that all Luer locks are tight.
5. The scrub nurse hands the tubing end with the blue cap to the circulating nurse.
6. Remove the blue cap from the end of the INTERMEDIARY tubing (Fig I. D) and the DAY SET tubing (Fig I. C).
7. Immediately connect the INTERMEDIARY tubing (Fig II. D) to the DAY SET tubing (Fig II. C). Tighten the nut to lock the both Luer connectors.
8. Press on the "Run/Stop" button of the Dual Flow Pump. Blinking green light will turn to solid green light. Open the clamp under the fill chamber to prime tubing.
9. Once the tubing has been primed, close off the clamp.

**END of the operation:**

Close the clamp under the pressure chamber

Disconnect the white Luer-lock (Fig III. E) under the check valve and discard tubing. The check valve will remain connected to the DAY SET tubing until the start of the next operation.

**START of the next operation:**

1. The circulating nurse removes the old check valve by disconnecting the used check valve (FIG IV. D) from the DAY SET tubing (FIG IV. C).
2. Immediately connect the blue Luer-lock (Fig. V. D) of the new INTERMEDIARY tubing to the blue connector (FIG V. C) of the DAY SET tubing as per step 7 in the Instructions for use.

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