APPROVAL

PhysIOL
Organization

Allée des Noisetiers, 4, 4031 Liège, BELGIUM
Address

Research and Development Director
Occupation

Christophe PAGNOULLE
Surname, Name

12th April 2023
Date DD.MM.YYYY


/ /
Signature

INSTRUCTION FOR USE

for a medical device

Surgical ophthalmological instrument for implantation of intraocular lenses
edition in accordance with the request dated 31.03.2023 No. 93880/2

The undersigned notary, member of
VAN HALTEREN, Associated Notaries, at
1000 Brussels, rue de Ligne 13, do hereby
certify that the signature / set on this
document is/are stated to be the
signature / of Christophe Pagnouille
Brussels, May 11th 2023

PhysIOL
Science Po
Allée des Noisetiers, 4, 4031 Liège, Belgium
Tél. +32 (0)4 341 40 00

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1. Name of a medical device

Surgical ophthalmological instrument for implantation of intraocular lenses, the set:

- a single use injector 1.2.3. Premium, model 226NG – 1 pc.;
- label - 11 pcs.;
- instruction for use - 1 pc.

Hereinafter — injector, single use injector 1.2.3. Premium, injector 1.2.3 Premium, device, 1.2.3. Premium, product, medical device, device.

2. Manufacturer and place of production

Public limited company PhysiOL, Belgium

Short name of a legal entity: PhysiOL, Belgium

Address: Allée des Noisetiers, 4, 4031 Liège, Belgium

Phone: +32 (0)4 361 05 49

Email: info@physiol.be

Place of production: PhysiOL, Allée des Noisetiers, 4, 4031 Liège, Belgium

3. Authorized manufacturer representative in the Russian Federation

RegMed prof. LLC

Address: 119607, Moscow, Ramenki str., 9, building 3, flat 505

Phone: +7(495)120-11-59, +7(495)120-11-64


Email: info@regmed.biz

4. Intended use

Designed to be used in combination with the preloaded PhysiOL posterior chamber intraocular lenses (IOLs) for implantation in the capsular bag after phacoemulsification during cataract surgery.

5. Potential consumer

The injector must be used by a skilled health care professional qualified in this type of surgery. As part of the risk analysis related to the product, the end user is an ophthalmic surgeon, trained and operating according to the good clinical practices.

6. Scope of application

The scope of application is ophthalmology, in the conditions of inpatient and outpatient medical institutions.

7. Medical indications

Surgical ophthalmological instrument for implantation of intraocular lenses is intended to be used in combination with a PhysIOL preloaded intraocular lens (IOL) contained in the preloaded cartridge for capsular bag implantation after phacoemulsification during cataract surgery.

8. Device performance

Surgical ophthalmological instrument for implantation of intraocular lenses has been designed with the aim to:

- Ensure consistent and reproducible IOL injection
- Limit the corneal incision size
- Eliminate IOL manual handling during lens loading and injection.


9. Contraindications

Contraindications have not been identified when used in accordance with the intended use indicated by the manufacturer.

10. Undesirable Side Effects

As with any surgical procedure, potential side effects accompanying cataract or ophthalmic implant surgery may include but are not limited to the following:

- Ocular infection (endophthalmitis, microbial keratitis);
- Inflammatory reaction (e.g. uveitis, Toxic anterior segment syndrome (TASS), hypopyon);
- Corneal edema;
- Corneal endothelial damage;
- Cystoid Macular edema (CME);
- Secondary surgical intervention (include, but are not limited to: lens repositioning, lens replacement, vitreous aspiration, iridectomy for pupillary block, wound leak repair, and retinal detachment repair);

- Loss of visual acuity;
- Hyphema;
- Vitreous Loss;
- Wound Leak (positive Seidel);
- Pain;
- Eye redness;
- Irritation;
- Elevated IOP
- Small tear in Descemet's membrane;
- Posterior capsular opacification..


11. Precautions, warnings, and recommendations

11.1. Precautions for storage and handling conditions

- For optimal functionality of the injector, do not store it at a temperature above 45°C or below 2°C.
- The SINGLE USE INJECTOR 1.2.3. PREMIUM is for single use only. Do not re-use this product to avoid the risk of cross-contamination.
- Do not resterilize.
- Do not use if the packaging is damaged.
- The injectors have an expiry date and may not be used past this date. PhysIOL makes no warranty for a device used past its expiry date.

11.2. Recommendations

- You must always follow the instruction for use. Any deviation may result in hazardous situation for the patient and damage the implant.
- All stages for preparing and using the 1.2.3. PREMIUM should follow each other in a continuous and quick workflow to minimize the risks of underperformance of the injector. The timeframe of injection procedure should be reduced to the maximum.

- For optimal injector performance, it is recommended that the OVD (Ophthalmic viscoelastic device), the IOL and the injector are stored at room temperature for a few hours.
- Reuse of the device is contraindicated as it could lead to (including but not limited to) device failure, cross-contamination and patient injury or illness.
- The disconnection of both the cartridge and the SINGLE USE INJECTOR 1.2.3. PREMIUM prior implantation is contraindicated as performance of the injector is no longer guaranteed.
- Corneal incision not less than 2.2 mm.

NOTE: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user is established.

12. Description of a medical device

12.1. Description

The SINGLE USE INJECTOR 1.2.3. PREMIUM is a non-implantable medical device aimed at transient, single use only. The materials used during the SINGLE USE INJECTOR 1.2.3. PREMIUM manufacture can present specific interactions with surrounding electrical systems, such as Magnetic Resonance Imaging (MRI) technology. The SINGLE USE INJECTOR 1.2.3. PREMIUM should be used within the scope of IOL implantation only.

12.2. Compatibility

As the single use injector 1.2.3. Premium is only designed to be used in combination with the preloaded PhysIOL IOLs (MICRO+ A(Y) 123, the hydrophobic Family 2 (MIC-GFY) MICROPURE 123 and ISOPURE 123 preloaded IOLs), the use with other preloaded IOLs is contraindicated.

Note: This IFU is intended to inform the user to the proper use of the SINGLE USE INJECTOR 1.2.3. PREMIUM only. Please refer to the preloaded IOL IFU for any information related to the use of the preloaded IOL.

13. Main technical features of the device

The overall characteristics of the injector are shown in Figure 1. Dimensional tolerance $\pm 10\%$.



Figure 1

Weight ($\pm 10\%$)	6.2 g
Injection Force (F)	$1,5\text{ N} < F < 25\text{ N}$

14. List and description of medical device materials

The table below shows device materials.


Name of the device and its components		Material, brand	Manufacturer	Type of contact
Injector	Housing	Styrene-butadiene copolymer K-RESIN Colorant: blue transparent 5821	ADMajoris, France	Indirect contact with a doctor through gloves
	Plunger	Polyoxymethylene POM ULTRAFORM 2320 Naturel	BASF, Germany	Indirect contact with a doctor through gloves.
	Silicone tip	Silicone elastomer d MED-4020 with color masterbatch MED-4102-4 (orange)	STERNE SAS, France	Short-term contact with the soft tissues of the eye (cornea, eyeball).
	Spring mechanism	Stainless steel AISI 302 HS	RDO SPRING, FRANCE	No contact
	Stopper	Polypropylene PPH 7060	Total petrochemicals, Belgium	Indirect contact with a doctor through gloves


15. Handling and injection

Before each use, examine the unopened box of both the injector and the preloaded IOL for integrity and expiration date. Check if the information on the carton box label is consistent with the information present on the secondary blister and self-adhesive labels (e.g. model and serial number).


1. The external surface of the injector blister is not sterile. Open the blister containing the injector, then drop the injector onto a sterile field.
2. For the use of the preloaded IOL refer to its IFU.
3. Whilst keeping the preloaded IOL container vertical in its cartridge, peel off the container lid.
4. Grasp the body of the injector (blue) by its central section, WITHOUT EXERTING PRESSURE ON THE PLUNGER, to prevent expulsion of the deformable pillow from the injector. The injector features a transparent stop catch at the level of the plunger.

5. Assemble the injector onto the preloaded cartridge: hold the injector firmly by its body (blue) in a vertical position with one hand, and the preloaded IOL container in the cartridge in the other hand, resting on a table. Then, “clip” the injector onto the visible part of the cartridge. Exert increasing pressure until you hear the “clip,” which indicates both elements have been assembled. If you do not hear the “clip”, there is a possibility that the connection of both elements could not be secured. In the event you do not hear the “clip”, first remove the assembled device from the container. Secondly, replace the assembled device vertically into the container and proceed once again to the “clipping”. The cartridge/injector non-return safety clip is then engaged and locked, ensuring proper solidarization of both the injector and the cartridge.
6. Whilst keeping the system in the container, push the plunger down until it abuts the safety catch. This operation aims to bring the lens into a pre-folded shape in the folding zone of the cartridge.
7. Keep the plunger in this position for about 3 seconds.
8. Allow the plunger to return in a controlled manner to its original position under the spring action. During this operation, keep one finger on the plunger to prevent the inappropriate, uncontrolled release of the plunger.
9. Completely release the safety catch of the injector by a twist motion, rotating it with one hand whilst keeping the injector in the container with the other hand. Prior to injection, always keep the cartridge with the IOL immersed in the container to avoid IOL or cartridge deterioration.
10. Remove the assembled system from the container and check that the cartridge is properly locked onto the injector. The non-return lock clip of the cartridge should be located just behind both lateral marks of the injector body. For correct positioning of the IOL at the moment of the injection, the inscription “PhysIOL” on the blue injector body must be FACING UPWARDS.
11. Rinse abundantly the IOL before implantation with sterile solution (NaCl 0.9%) or Balanced Saline Solution (BSS): introduce the non-piercing injection cannula into the small hole on the body of the injector and leave the fluids in the cartridge to flow onto the operating field for 2 to 3 seconds, cartridge pointing downwards. Afterward, inject a generous amount of viscoelastic into the same hole.
12. Start to inject the IOL immediately by exerting a constant, controlled, and regular pressure on the plunger. When the first leading haptics of the IOL are out of the cartridge, release the plunger a few millimeters to ensure that both trailing haptics are not captured by the silicone pillow.
13. As with conventional lens injection procedures, it is possible that the trailing haptics will have to be repositioned inside the capsular bag with the help of a suitable positioning hook.
14. Carefully remove all traces of the viscoelastic solution by irrigation/aspiration, particularly in the space between the lens and the posterior face of the bag.


NOTE: (1) No maintenance is required for the use of the device as it is single use only.

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



16. Usage, transportation and storage conditions

16.1. Instruction for use

Keep for several hours before use at temperature greater than or equal to 17 °C.
Use at temperatures between 17 °C and 45 °C.

16.2. Shipping and Storage

Store and transport at a temperature from 2 °C to 45 °C with a relative humidity of no more than 70 %. Keep for several hours before use at a temperature greater than or equal to 17 °C.

Storage conditions must be observed during shipment and at any intermediate stages of storage.

PhysIOL is not responsible for use of the SINGLE-USE 1.2.3 PREMIUM INJECTOR, whose storage, use, or transport conditions were violated.

17. Expiration date

Valid for 5 years.

18. Symbols on labels

	Manufacturer, address, name.
	Do not re-use
	Refer to the instruction for use
	Refer to the instruction for use
	Avoid exposure to sunlight
	Protect from moisture
	Do not re-sterilize
	Best before date

LOT	Lot code
STERILEEO	ethylene oxide sterilization
CE	the compliance with the European Directive
	The sterilization indicator

19. Packaging

The injector in a blister is packed into a carton.

Should any part of the packaging be damaged, do not use the device and send it back to the supplier. The expiration date is indicated on the outside of the injector package. Do not use the injector after its expiration date.

20. Information on the substances contained in the medical device

The device does not contain medicines, drugs or pharmaceutical substances.

21. Information on materials of animal and (or) human origin

The SINGLE USE INJECTOR 1.2.3. PREMIUM has been verified for its biocompatibility and is deemed safe for use as a transient surgical invasive device. None of the materials composing the injector include the following:


- Medicinal substances, including a human blood or plasma derivative;
- Tissues or cells or their derivative, of human origin;
- Tissues or cells or their derivatives, of animal origin.

22. Manufacturer warranty

PhysIOL guarantees compliance with the declared properties and features of the medical device, subject to the conditions of handling, transportation and storage.

PhysIOL shall not be liable for any injury or damage suffered by a patient as a result of:

- Improper prescription or model selection,
- Improper device handling,
- Surgical technique applied,

- Any iatrogenic error caused by the implanting surgeon,
- Off-label use of the product,
- Use of expired devices
- Improper compliance to the information contained in this IFU.

PhysIOL makes no expressed or implied warranties in connection with the resale of its medical devices. The user should properly follow those instructions for use as a general recommendation.

23. Description of the sterilization method


The medical device is sterilized in accordance with ISO 11135:2014 (low-temperature ethylene oxide sterilization).

24. Maintenance and repair

This medical device is single used. The device is beyond repair and maintenance.

25. The list of applicable national and international standards

Standard No.	Title
Directives and standards	
Directive 93/42/EEC	
Quality management systems	ISO 9001: 2015 ISO 13485: 2016
Sterilization	EN ISO 11135:2014
Package, labeling	EN 1041+A1:2013 EN ISO 11607-1:2018 EN ISO 11607-2:2018
Biosafety	EN ISO 10993-7:2008 ISO 10993-5:2009 ISO 10993-10:2010 ISO 10993-11:2017
Risk analysis	EN ISO 14971:2012

26. Utilization

The devices that have come into contact with blood and/or other biological liquids must be utilized and/or disposed of in accordance with SRN 2.1.3684-21 (Class B).

Unused devices (which have not come into contact with blood and/or biological liquids), including those that have expired, are subject to utilization and/or disposal in accordance with SRN 2.1.3684-21 (Class A).

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(Research and Development Director)

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