Ophthalmic stand VS-1000, VS-2000

User manual

Viscope Inc.

Add: No.86, Qingyang Road, Changzhou city, Jiangsu Province.

Tel: +86 519 85966517

Contents:

1. Introduction
2. The requirement to the personnel operating this product
3. Description, process of assembly and use order to destination
 - 3.1 Main technical characteristics
 - 3.2 Product structure
 - 3.3 Assembly process
 - 3.4 Operating procedure to destination
 - 3.5 Wiring Diagram
4. Maintenance and repair
5. Storage conditions and operation
6. Transportation conditions
7. Utilization order
8. Cleaning and disinfection
9. Guarantee certificates of the producer

1. Introduction

The product is intended for carrying out a wide range of ophthalmologic researches and is basic system for placement of the ophthalmologic devices intended for these purposes (creates space for their placement).

Is the basic module of equipment of the clinics which are profiling on diagnostics and restoration of broken sight.

Basic elements of installation — a chair of the patient, a rack with a little table for devices and a lamp of the general lighting.

The combine rack with a mobile little table creates space for placement:

Slit lamp of any firm manufacturer

Autorefkeratometr or keratometer of any firm manufacturer

Chart projector of any firm manufacturer

Foropter of any firm manufacturer

Manual ophtalmoscope of any firm manufacturer

2. The requirement to the personnel operating this product

Specially trained and certified personnel is allowed to independent operation of products of medical equipment only 18 years, suitable on a condition of qualification to performance of the specified works aren't younger.


Before the admission to work the personnel has to pass introduction and primary personnel safety notification with display of safe and rational working methods with registration in instructing magazines.

In rooms where the medical equipment is constantly operated, have to be hung out in a place available to the personnel, taking into account norms of production sanitation, the safety rule in which personnel actions in case of accidents, fires, electric traumas have to be accurately formulated.

Modification of the equipment can be executed only by specialists in service of the manufacturing company or other authorized officers.

3. Description, process of assembly and use order to destination

The ophthalmologic combine is used in work by doctors ophthalmologists for sight research, for selection of points, etc. For installation use in other purposes responsibility is born by the user.

3.1 Main technical characteristics:

VS 1000

Chair type	Motorize
Chair color	Black
Angle of chair rotation	$\pm 60^\circ$
Maximum chair lifting	500-700 mm
Arm angle of foropter rotation	$\pm 30^\circ$
Table-top size	840 mm x 405 mm
Maximum load of a chair	300 kg
Angle of table-top rotation	90°
Back inclination	165°
Height of lifting of a head restraint	150 mm
Electric potential	110/220B 50Hz-60Hz
Fluorescent lamp	12V
Dimensions	1700 mm x 1740 mm x 2005 mm
Weight	175 kg
Working temperature	$-30^\circ\text{C} \sim +50^\circ\text{C}$
Humidity	Less 90 %
Safety lock	5A

VS 2000

Chair type	Motorize
Chair color	Wight
Angle of chair rotation	$\pm 60^\circ$
Maximum chair lifting	500-700 mm
Arm angle of foropter rotation	$\pm 30^\circ$
Table-top size	840 mm x 405 mm
Maximum load of a chair	300 kg
Angle of table-top rotation	90°
Back inclination	165°
Height of lifting of a head restraint	150 mm
Electric potential	110/220B 50Hz-60Hz
Fluorescent lamp	12V
Dimensions	1400 mm x 1000 mm x 1900 mm
Weight	120 kg
Working temperature	$-30^\circ\text{C} \sim +50^\circ\text{C}$
Humidity	Less 90 %
Safety lock	5A

3.2 Product structure

VS 1000

1. Lamp
2. Lamp pipe holder.
3. Core for fastening of a foropter
4. Arm for fastening of a foropter
5. Chair
6. Table-top
7. System of horizontal shift of a table-top
8. Arm of fastening of a table-top
9. Control unit
10. Control panel
11. Central block
12. Handle of blocking of rotation of a chair
13. Handle of an inclination of a chair
14. The socket under a food cable
15. Electric lifting mechanism
16. Additional control unit of a workplace
17. Sockets of power supply of devices
18. System of an emergency stop of the lifting mechanism of a chair
19. System of stabilization of tension and uninterrupted food
20. Sockets for charging of the additional equipment (2 sockets)

VS 2000

1. Lamp
2. Lamp pipe holder.
3. Core for fastening of a foropter
4. Arm for fastening of a foropter
5. Short pipe holder
6. Chair
7. Table-top
8. System of horizontal shift of a table-top
9. Arm of fastening of a table-top
10. Control unit
11. Control panel
12. Central block
13. Handle of blocking of rotation of a chair
14. Handle of an inclination of a chair
15. The socket under a food cable
16. Electric lifting mechanism
17. Additional control unit of a workplace
18. Sockets of power supply of devices
19. System of an emergency stop of the lifting mechanism of a chair

20. System of stabilization of tension and uninterrupted food
21. Sockets for charging of the additional equipment (3 sockets)

Accessories (VS 1000 and VS 2000):

1. Support under feet
2. The holder for a foroptor
3. The holder for a projector of signs
4. Bolts and keys six-sided for assembly
5. Safety locks electric
6. Pedals for remote control of the lifting mechanism
7. Maintenance instruction


3.3 Assembly process

First of all, open the packing box and read the instruction manual carefully, then check the parts and tools.

Put the main part on a suitable place. Make sure that the mainframe is stable.

Install the chair

Take out the backrest, unscrew the bolts and nuts on the chair frame, align the three holes on the back and their corresponding holes on the chair frame, fix the back to chair frame reliably with the bolts and nuts(see Fig.1-4).

(Fig 1-4)

Then connect the chair up/down switch line to the motor cover(see Fig.2-4), take out the headrest and fix it to the back reliably by fastening the knob (see Fig.2-4)


(Fig 2-4)

1. Install the short post


Open the back cover through loosening the two screws on it, take out the short post and make the power line of the projector go through the hole on the mainframe that close to chair, then fit the short post on the mainframe and fix it to half fastening by using the screw (M10*100) and gasket (see Fig .3-4).

(Fig 3-4)

2. Install the bend lamp tube and post.

Take out the bend lamp tube and make the power line of the lamp go through the long post, align the bend lamp tube and the holes on the top of the long post, fit the bend lamp tube on the post through fastening the 4 screws, in this step, you should keep the lamp tube balance (see Fig.4-4), then install the long post with the same way of installing the short post (see Fig.4-4)


(Fig 4-4)

3. Put the cannula in the screws between the posts, then fix it to half fastening with the hexagon screw (M8*120), adjust the post' positions to make them vertical to the ground. after this ,fasten the corresponding screws, connect he power lines of the projector and lamp to corresponding positions of the circuit board

4. Install the phoropter arm

Take out the phoropter arm and hang it onto the short post, stick the sliding cover into the hole on the top of the rocker, then install the horizontal sliding rod and tighten it with the screw. at last, fix the phoropter arm locker to the phoropter arm rocker, make sure that phoropter arm can run softly and comfortably (see Fig.5-4)

(Fig 5-4)

5. Install the chart projector tray
6. Fix the chart projector onto the tray with prepared screws and plug in.
7. Hang the view tester onto the horizontal sliding rod, tighten the knob then adjust the view tester to a suitable position.
8. Place the auto refractometer or other related instruments on the tabletop, connect it to the outlet below the with the prepared power line.
9. Power on to test the wire connection, check the keys' function of the control panel, if the connection is correct.

Note:

1. By default food on the device moves always when the workplace is included. To regulate food on the control panel, remove a short closing ring (see an electric circuit).
2. When the table-top is turned on 20 degrees, food is switched off when the provision of a table-top comes back back, food proceeds. If you separate contacts on a payment, food will always move, irrespective of table-top turn (see an electric circuit).

3.4 Operating procedure to destination

1. Include a workplace, having pressed the button [ON/OFF] on the control panel.
2. Adjust at the necessary height a forofter, using the handle of blocking of an arm for a forofter.
3. Adjust chair height, using buttons [Up] and [Down] on the control panel. You also can use buttons on the control panel under a table-top (fig. 23) or on the back party of a chair (fig. 24). The workplace is equipped with the torzmozny device from below table-tops (fig. 25). When it works, the chair stops lifting.

Fig. 23

Fig. 24


Fig. 25

4. If on a table-top there are 2 devices, you can move a table-top to the left or to the right, a little her having pressured.

Note: At device change, food switches over for 12 V. For tension control, switch yellow wires for tension change (see an electric circuit).

5. If additional lighting is necessary, you can include the top lighting lamp, having pressed the button [Light] and to adjust lighting having turned the spherical handle.
6. For use of a projector of signs, press the button [Projector] to include it.
7. After carrying out measurements, return an arm for a forofter to a starting position and switch off a workplace.

3.5 Wiring Diagram

4. Maintenance and repair

Service

At the correct operation, the ophthalmologic combine doesn't demand special service. In regular service there is no need. In case of malfunction, contact qualified personnel or our distributors. We recommend to carry out equipment service each three years, for observance of security measures and maintenance of normal functioning of the ophthalmologic combine.

Cleaning and disinfection

The correct cleaning has to be made by a rag with soap solution (not the damp!) Plastic surfaces can be disinfected.

Attention:

During cleaning of the ophthalmologic combine, it is impossible to allow liquid hits inside or on keys. Before cleaning, always disconnect a workplace from a network.

The key of lighting can be not included.

Before again including a workplace of the doctor, be convinced that it completely dried.

Wipe all painted parts and an upholstery of sitting by damp fabric. You can add a little detergent in water. Liquid or detergent shouldn't get in a workplace of the doctor. workplace.

Here are some frequent failures or something happened possibly when you operate this stand. You can examine the facility or repair it according to following list.

Failures	Possible reasons	Repair suggestion
1. The indicator light has no light when power on.	The fuse is blown The plug is not inserted well The outlet has no power	Check and change Insert the plug well Test and confirm
2. The indicator light has light but it doesn't work yet.	The plug of control console is not inserted well	Check and insert it well
3. The machine has power but doesn't work	The control console is broken The voltage doesn't meet the required	Check and change Test and confirm
4. The projector has no power	The circuit is connected wrongly The power line of projector is short-circuit	Check the circuit and fix the breakdown Change the power line

5. The lamp has no light	The lamp is damaged The circuit is connected wrongly	Change the lamp Check the circuit and fix the breakdown
6. The two DC supply power on the board is not prepared	The circuit is connected wrongly The plug is loose	Check the circuit and fix the breakdown Insert the plug well
7. The chair can go down but can't rise	The up inhibiting switch under the table is damaged	Change the switch
8. The chair goes down automatically	The screw on the electromagnetic valve is loose	Fix the screw
9. When going up or down, there's a big noise with the motor	Someplace is leaked The temperature is low and the density of the oil is high	Check and change It will be normal within ten minutes or so after it is started

5. Storage conditions and operation

It is necessary to avoid a condition of the environment provided below at operation and storage of the equipment:

- Where the equipment is affected by water vapor.
- Where the equipment is affected by sunshine.
- At frequent changes of temperature (normal temperature for equipment operation in the range - 30oC ~ + 50oC, and humidity of air less than 90%).
- Where any heating devices are used close to the equipment.
- At high humidity and where there are heat removal and/or ventilation problems.
- Where the equipment is exposed to excessive shake-ups and vibrations.
- Where the equipment can be affected by chemical or combustible substances.


6. Transportation conditions

It is necessary to avoid places at transportation where ambient temperature falls below -25 °C or above 40 °C. Humidity has to be supported in limits between 30% and 75%. It is necessary to avoid conditions in which the equipment is affected by strong blows or vibrations.

The ophthalmologic combine transport covered transport in covered vehicles according to requirements of GOST P 50444-92 and the rules of transportation of goods existing on each type of transport.

Transportation conditions – under the terms of storage of 5 GOST 15150-69.

Ophthalmologic combines have to be stored in warehouses of the supplier and the consumer in storage conditions of 1 GOST 15150-69.

7. Utilization order

Any equipment at utilization has to be separated from other waste. Utilization of electric devices and the electronic equipment is carried out according to existing national legislation and other regulations of the country in which territory this equipment is operated.

8. Cleaning and surface disinfection

Cleaning has to be made by a rag with soap solution (not the damp!) Plastic surfaces can be disinfected. For disinfection of surfaces, moisten fabric with not abrasive cleaner and accurately wipe surfaces.

ATTENTION: Don't use the caustic or abrasive cleaners.

9. Guarantee certificates of the producer

Warranty period – 12 months.

The guarantee doesn't extend:

- On the damages caused by non-compliance with service regulations
- On the damages caused by use of a product for for which it isn't intended.
- On the damages caused by negligent storage, transportation by the owner or the trading organization.
- On the damages caused by the wrong leaving.
- On the mechanical damages caused by the consumer.

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