20

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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of SAC (XIAMEN) MEDICAL TECHNOLOGY CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature: Zhu Chaoyang

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(Date: Jan. 05, 2023)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Date: 8th of November, 2021

«Approved»

Wu Xiaohong (signature)
Wu Xiaohong,
International sales director
SAC (Xiamen) Medical Technology Co.Ltd.

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE

**Single Use Sterile Pulse Flushing System for Wound Care,
model SAC WPWS-1000**

Edition 2.

1 Name of the medical device

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000, consisting of:

1. Main body handle – 1 pc.
2. Short debridement nozzle – 1 pc.
3. Long debridement nozzle – 1 pc.
4. Irrigation line – 1 pc.
5. Drainage line – 1 pc.
6. Battery box – 1 pc.
7. Flow controller – 1 pc.
8. Instruction for use – 1 pc.

Hereinafter referred to as the “Pulse Flushing System”, “system”, “device”.

2 Information about the manufacturer, developer and authorized representative of the manufacturer

Manufacturer

SAC (Xiamen) Medical Technology Co., Ltd.
Room 401, No. 218, Houxiang Road, Haicang District, Xiamen, Fujian Province, 361022, P.R. China.

Telephone: +86 592 6087101

Fax: +86 592 6087101

E-mail: info@double-medical.com

Developer

SAC (Xiamen) Medical Technology Co., Ltd.
Room 401, No. 218, Houxiang Road, Haicang District, Xiamen, Fujian Province, 361022, P.R. China.

Telephone: +86 592 6087101

Fax: +86 592 6087101

E-mail: info@double-medical.com

Authorized representative of manufacturer in the Russian Federation

Zdorovie LLC

Legal address: 115054, Moscow, Dubninskaya street, 57, building 2, unit IV office 9,10

Business address: 115054, Moscow, Dubninskaya street, 57, building 2, unit IV office 9,10

Tel.: 8 (495) 134-34-42

E-mail: sales@trauma-msk.ru

3 The device purpose with the indication of a potential customer

3.1.Device purpose

Intended for treating wounds, bone marrow cavity and cleansing soft tissues in orthopedic surgery for locomotor apparatus injuries.

3.2.Potential customers

This device must be used only by qualified specialists, surgical specialization doctors (surgeons, orthopedists, traumatologists).

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000

4 Description of the device

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000, is composed of the main body handle, debridement nozzle, irrigation line, drainage line, flow controller and battery box. By the vibration of the pulse flow, bacteria and dirt in the wound detach from tissues, which makes it easy to clean thoroughly. It can clean extensive wounds, no easy-clean crack, and deep wounds effectively.

5 Indications, contraindications and possible side effects

Indications:

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000 is used for locomotor apparatus injuries, non-union of fractures and other diseases associated with orthopedic surgery; for hip and knee joint replacement, joint restoration, osteomyelitis and other diseases associated with joint surgery; for spinal cord injuries, spinal tumor, scoliosis, degenerative cord diseases.

Contraindications:

It must not be used for the following patients:

- Severe diabetes;
- Patients with hypoproteinemia;
- Patients with severe hemorrhagic shock.

Possible side effects: not detected.

6 Precautions

- 6.1. This device is for single-use, reuse is not allowed. Reuse may cause infection to wounds.
- 6.2. This device is sterilized. If the sealed bag is damaged before use, do not use it.
- 6.3. Do not use in the place with flammable and explosive anesthetic or gas.
- 6.4. Do not immerse the mainframe in the fluid, to avoid the system corruption.
- 6.5. When using, the flushing fluid tube's hanging position should be more than 20 cm higher than the mainframe.
- 6.6. Turn off the flow controller when suspending the surgery.
- 6.7. Do not use the product near the high frequency equipment, for instance, high frequency electrical knife, mobile phones. Otherwise, it would cause the risk of malfunction.
- 6.8. Do not directly administer to spinal nerves during operation.
- 6.9. This product is disposable and cannot be repeatedly sterilized and used. The use time is no more than 30 min.
- 6.10. The rinsing solution is physiological saline, provided by the hospital.
- 6.11. Before using this product, please read the instructions thoroughly to avoid improper operation due to which the product cannot achieve the expected results.

7 Usage procedure

1. Peel out the pack and remove all the accessories.
2. Select an appropriate nozzle, connect it to the main body assembly and lock it. The short debridement nozzle is used for cleaning large-area wounds; long debridement nozzle is used for cleaning deep wounds.
3. Attach the drainage line to the suction pump.
4. Attach the irrigation line to the irrigation source.
5. Select the low speed or high speed fluid delivery.
6. Use a flow controller to turn rinse solution on or off.

There are three modes:

“Run Stop”: the device stops operating in the current mode.

“Slow Speed”: select this mode, when the wound requires a large flushing force. The flow at reduced gear must be not less than 300 ml / min, and the water flow frequency must be not less than 3 times per minute.

“Fast Speed”: select this mode, when the wound requires a small flushing force. The flow rate at high gear must be not less than 500 ml / min, and the water flow frequency must be not less than 4 times per minute.

8 Technical characteristics.

Table 1

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000 consisting of	
Main body handle	
The device overall dimensions	
Width (± 2) mm	140
Height (± 3) mm	190
Weight (± 10) g	230 g
Static axial load of connections N	≤ 15 N
Short debridement nozzle	
Connecting diameter, D1($\emptyset$) (± 2), mm	12mm
The diameter of a protective screen against splashing, D2($\emptyset$) (± 2), mm	34mm
Length, L1 (± 5), mm	40mm
Length, L2 (± 5), mm	114mm
Weight (± 2) g	17,4g
Long Debridement nozzle	
Connecting diameter, D1($\emptyset$) (± 5), mm	30mm
Tube diameter, D2($\emptyset$) ($\pm 0,5$), mm	9mm
Length, L1 (± 5), mm	40mm

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000

Length, L2 (± 5), mm	205mm
Weight (± 2) g	15,2g
Irrigation line	
Diameter, D($\varnothing$) (± 1), mm	9 mm
Length, L (± 5), cm	290 cm
Maximum pressure MPa (no leakages at room temperature)	0,05 MPa
Drainage line	
Diameter, D($\varnothing$) (± 1), mm	9mm
Length, L (± 5), cm	300 cm
Maximum pressure MPa (no leakages at room temperature)	0,05 MPa
Flow controller	
Overall dimensions length x width x height (± 2)	25x38x17,5mm
Liquid aspiration appliance (provided by the hospital)	The negative pressure value of the liquid aspiration appliance: 100-300 mm•h
Wound flushing solution (provided by the hospital)	0.9% NACL saline
Battery	
The device overall dimensions length x width x height (± 2) mm	195x75 x21mm
Weight (± 10) g	299 g
Voltage	12 Vdc
Power supply elements, type LR6	1.5 V * 8AA without recharging
Input power	12 VA
Operating mode	short-term
Operating time, min	≤ 2
Interval, min	≥ 5
System connection	fixed wire
Length, L ($\pm 0,1$), m	2.70 m
Working part	Type B
Ingress protection rating	IPX0
The device total weight (± 50), g	1010 g

8.1. Performance Specifications.

8.1.1. Irrigation line and drainage line can tolerate pressure of 0.05 MPa without leaks at room temperature.

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000

- 8.1.2. Each junction can sustain static axial tension not less than 15 N.
- 8.1.3. There is no backflow or leakage of liquid in the tube when the flow controller is under a negative pressure of 200 mmHg.
- 8.1.4. The high gear flow should be no less than 500 ml/min, and water flow frequency should be no less than 4 times per minute; the low gear flow should be no less than 300 ml / min, and water flow frequency should be no less than 3 times per minute.
- 8.1.5. Flushing pressure: 0.310-0.703 kg / cm.
- 8.1.6. Total extractable metal content: Sample test solution and blank control solution of the same batch, the total content of lead, zinc, tin and iron is ≤ 1 ug/mL, cadmium ≤ 0.1 ug/mL.
- 8.1.7. pH: The difference between the pH of the sample test solution and blank solution of the same batch is ≤ 1.5 .
- 8.1.8. Reducing substance: The difference in consumption of potassium permanganate titration solution (0.002 mol / L) between the sample test solution and the blank control solution shall not exceed 2.0 mL.
- 8.1.9. Depending on the perceived mechanical influences, the product belongs to portable, but not intended for use when carrying and moving of the device within a stationary room. A medical device is used only within a medical institution.
- 8.1.10. Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000 is a device with an internal power supply.

8.2 Characteristics and parameters

8.2.1 Tube system sealing performance

Each end of the irrigation line and drainage line should be connected to the negative pressure source which is 0.05 MPa. Immerse the other end of the tube into the water. There should be no leakage when the clamp is closed.

8.2.2 Connecting firmness

The connection part of the tube should bear the separate force which is no less than 15 N.

8.2.3 Flow controller

It should open and close the tube effectively. When closed and under the condition of 200 mmHg negative pressure, fluid in the tube should no flow back or leak.

8.2.4 Work noise

When functioning properly, the product shall not produce more than 95 dB of noise when measured at a distance of 2 m with a standard psophometer.

8.2.5 Performance requirements

The product has high gear flow and low gear flow mode, and its performance requirements are listed in Table 2, as below:

Table 2. Performance Requirements of High gear and Low gear

Gear	Flow requirements	Hydrodynamic frequency
High gear	≥ 500 ml/min	≥ 4 times/s
Low gear	≥ 300 ml/min	≥ 3 times/s



9 Interpretation of symbols used in the marking of the consumer and transportation

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000

Packaging

Table 3. Symbols and Designations

Symbol	Interpretation	Symbol	Interpretation
	Manufacturer		Manufacture date
	Ethylene oxide sterilization		Prohibition of reuse
	Catalogue number (part number)		Lot number
	Compliance with European standards		Please, refer to the instructions for use
	Do not use if the packaging is damaged		Use before
	Temperature limit		Humidity range
	Atmospheric pressure limitation		Caution!
	Type B applied part		The accumulator battery must not be disposed of in undesignated areas
	Fragile. Handle with care (for transportation packaging)		The right position of the packaging is this side up (for transportation packaging)
	Keep away from moisture (for transportation packaging)		Direct current (for the battery)

10 Requirements for maintenance and repair of the medical device

The devices do not require maintenance and repair..

11 Sterilization

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000 shall be sterilized with ethylene oxide.

The information on the sterilization process validation is represented in a separate document "EO Sterilization validation report: (see Annex 3 of this document).

EO Sterilization is an outsourcing process, which is controlled according to the control procedure of outsourcing processes.

Sterilization company of the device:

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000

Surgaid Medical (Xiamen) Co., Ltd.

3F Comprehensive building, No. 218 Building, Houxiang road, Haicang District, Xiamen, 361026, P.R. China.

Tel: +86-592-6880677

e-mail: mail@surgaid.com

The devices must not be re-sterilized. The re-sterilization may disrupt the structural integrity of devices, create a risk of contamination of devices, lead to infection of a patient and / or to cross-infection, including the transfer of infectious diseases from one patient to another, as well as become a cause of injury or death of a patient.

The manufacturer carries out the sterilization process validation. The manufacturer has evaluated the medical device safety and performance. The manufacturer has established and applies a quality management system for the design engineering and development, manufacture and distribution of the Pulse Flushing System.

The verification and validation procedures create a basis for the preparation of the documentation confirming that the results of engineering the device comply with the project specifications, and that the device is suitable for its intended use.

In "EO Sterilization validation report" the device sterilization method, preparation and procedure are indicated. All procedures are repeated (maximum - twice for one batch) under the control of the quality department of the device manufacturer. The conditions for carrying out the sterilization stage have the parameters, which values are clearly formulated with an indication of the established ranges and errors.

12 Transportation, storage and operation conditions

Transportation

Do not throw, do not expose to pressure, do not crumple during transportation. The device must be well protected from the influence of rain and snow. The number of tiers, when stacking, must be less than five.

Transportation conditions:

Temperature: -40 to +70°C

Relative humidity: 10-80%

Atmospheric pressure: 1000-1060 hPa

Storage

The device should be stored in a clean, dry, ventilated location. Store away from corrosive materials.

Storage conditions:

Temperature: +18 to +28°C

Relative humidity: 10-80%

Atmospheric pressure: 1000-1060 hPa

Usage

The devices are used in medical institutions.

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000

Environmental operation conditions of the medical device:

Temperature: +18 to +28°C

Relative humidity: 10-80%

Atmospheric pressure: 1000-1060 hPa

13 Disposal

The medical devices "Single Use Sterile Pulse Flushing System for Wound Care" after using according to their intended purpose belong to epidemiologically hazardous waste (devices contaminated with blood and / or other biological fluids of patients).

In this case, the medical organization carries out only the collection and temporary storage of the above waste. The collection is carried out into the single-use bags or containers with the colour markings established for hazardous medical waste indoors or in places for temporary storage of the medical waste until subsequent transportation to the disposal site.

The transportation of the medical waste from the medical institutions and their disposal are carried out by specialized organizations holding a license to handle epidemiologically hazardous medical waste.

Packaging, defective devices and unused devices with expired shelf life belong to epidemiologically safe waste (non-hazardous waste that does not contact with biological fluids of patients) that is close in composition to municipal solid waste (MSW), with the exception of the battery, which must be disposed of through special licensed organizations.

The collection of the non-hazardous waste is carried out in reusable vessels or disposable bags (with the exception of bags of the colour established for the hazardous waste). The disposable bags are located inside the reusable containers. The waste collection vessels must be marked and installed on a special site (indoors). The reusable containers after emptying are subject to washing and disinfecting. The procedure for washing and disinfecting the reusable containers is determined in accordance with the waste management scheme of each particular organization.

The transportation and disposal is carried out taking into account the requirements established for this territory, in accordance with the sanitary legislation of a particular country applied for the maintenance of the territories of populated areas and the treatment of the waste of manufacture and consumption.

14 Shelf life

The device shelf life is 2 years from the date of sterilization.

Do not use upon the shelf life expiration.

15 List of international standards that the medical device complies with

Directive 93/42/EEC, EN 980, EN ISO 1041, EN IEC 62366, IEC 60601-1, EN IEC 60601-1, EN 60601-1-2, EN ISO 10993-1, EN ISO 10993-2, EN ISO 10993-5, EN ISO 10993-10, EN ISO 10993-12, EN ISO 10993-18, ISO 11737-2, EN 550, EN ISO 13485, EN ISO 14971, EN ISO 15223-1, MEDDEV 2.7.1, rev. 4, MEDDEV 2.12.1, rev. 8, MEDDEV 2.12.2, rev. 2, ISO 11607

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE

Single Use Sterile Pulse Flushing System for Wound Care, model SAC WPWS-1000

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