Ref: 20220876.01 IG


LEGALISATION

Seen for legalisation by me, Anneke Heijnen, Civil-Law Notary, practising in Katwijk, the signature of: Mr Mohammed Aldirawi, residing at 2272 XD Voorburg, Leidschendam-Voorburg, Rozentuin 74, born in Gaza on the fourth day of July nineteen hundred and seventy-six, identifying himself with his passport with number 17FV17947, issued in Brussels on the ninth day of August two thousand seventeen,

in this respect acting in his capacity as Holder of a Power of Attorney with representative powers of **Mentor Medical Systems B.V.**, a private company with limited liability (besloten vennootschap met beperkte aansprakelijkheid), incorporated and existing under the Laws of the Netherlands, having its corporate seat in Leiden, with offices at Zernikedreef 2, 2333 CL Leiden, registered with the Traderegister of the Chamber of Commerce under number 28062039, and as such authorised to validly represent Mentor Medical Systems B.V., according to an extract of the commercial register of the Chamber of Commerce, which is attached to this document.

The aforementioned statement only relates to the signature of Mr Mohammed Aldirawi, aforementioned, and his capacity as Holder of a Power of Attorney of the aforementioned company Mentor Medical Systems B.V., and does not relate to the contents of the document on which his signature was endorsed.

Signed in Katwijk, the Netherlands, on the twenty-second day of April two thousand twenty-two.

Business Register extract Netherlands Chamber of Commerce

CCI number 28062039

Page 1 (of 3)

Legal entity

RSIN	802689954
Legal form	Besloten Vennootschap (comparable with Private Limited Liability Company)
Name given in the articles	Mentor Medical Systems B.V.
Corporate seat	Leiden
First entry in Business Register	31-05-1994
Date of deed of incorporation	16-05-1994
Date of deed of last amendment to the Articles of Association	29-12-2009
Issued capital	EUR 45.380,00
Paid-up capital	EUR 45.380,00
Filing of the annual accounts	The annual accounts for the financial year 2020 were filed on 22-12-2021. The annual accounts have not been adopted.

Company

Trade name	Mentor Medical Systems B.V.
Company start date	22-11-1993
Activities	SBI-code: 32502 - Manufacture of medical instruments and supplies (no dental laboratories)
Employees	140

Establishment

Establishment number	000015832236
Trade name	Mentor Medical Systems B.V.
Visiting address	Zernikedreef 2, 2333CL Leiden
Telephone number	+31717513600
Fax number	+31715217422
Internet address	www.mentormedical.nl
Email address	info@mentormedical.nl
Date of incorporation	22-11-1993
This legal person has conducted the establishment since	16-05-1994
Activities	SBI-code: 32502 - Manufacture of medical instruments and supplies (no dental laboratories)

For further information on activities, see Dutch extract.

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This extract has been certified with a digital signature and is an official proof of registration in the Business Register. You can check the integrity of this document and validate the signature in Adobe at the top of your screen. The Chamber of Commerce recommends that this document be viewed in digital form so that its integrity is safeguarded and the signature remains verifiable.

2022-04-22 12:07:29

Business Register extract

Netherlands Chamber of Commerce

CCI number 28062039

Page 2 (of 3)

Employees 140

Sole shareholder

Name Mentor B.V.
Visiting address Zernikedreef 2, 2333CL Leiden
Registered under CCI number 64647323
Sole shareholder since 21-12-2015 (registration date: 21-12-2015)

Board member

Name Draaijer, Matthijs Geert
Date of birth 05-05-1972
Date of entry into office 28-04-2020 (registration date: 07-05-2020)
Title Director
Powers Solely/independently authorised

Authorised representatives

Name Henkes, Marcel
Date of birth 08-08-1961
Date of entry into office 23-11-2001
Title Inkoper
Contents of power of attorney There are other restrictions. See Dutch extract.
Commencement (present) power of attorney 10-11-2014

Name Tuinier, Erwin
Date of birth 08-05-1965
Date of entry into office 21-06-2012 (registration date: 25-06-2012)
Title Controller
Contents of power of attorney Full power of attorney

Name Mos, Jacoba Alberta
Date of birth 07-07-1966
Date of entry into office 29-09-2017 (registration date: 20-10-2017)
Contents of power of attorney There are other restrictions. See Dutch extract.

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This extract has been certified with a digital signature and is an official proof of registration in the Business Register. You can check the integrity of this document and validate the signature in Adobe at the top of your screen. The Chamber of Commerce recommends that this document be viewed in digital form so that its integrity is safeguarded and the signature remains verifiable.

2022-04-22 12:07:29

KVK

Business Register extract Netherlands Chamber of Commerce

CCI number 28062039

Page 3 (of 3)

Name	Claes, Hilde Bernadette J.
Date of birth	08-07-1967
Date of entry into office	29-09-2017 (registration date: 25-10-2017)
Contents of power of attorney	There are other restrictions. See Dutch extract.

Name	Aldirawi, Mohammed
Date of birth	04-07-1976
Date of entry into office	02-06-2020 (registration date: 29-06-2020)
Title	Gevolmachtigde
Contents of power of attorney	Full power of attorney

Name	Koopman, Cees Willem
Date of birth	18-08-1986
Date of entry into office	03-06-2021 (registration date: 04-06-2021)
Contents of power of attorney	There are other restrictions. See Dutch extract.

Extract was made on 22-04-2022 at 12.07 hours.

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This extract has been certified with a digital signature and is an official proof of registration in the Business Register. You can check the integrity of this document and validate the signature in Adobe at the top of your screen. The Chamber of Commerce recommends that this document be viewed in digital form so that its integrity is safeguarded and the signature remains verifiable.

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Quality Operations Director
Должность/Occupation

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Instruction for use:

«Becker tissue Expanders/Breast Implants with Siltex shell in the following versions»/

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revision surgery may be necessary including implant removal or replacement. The surgeons in the European Community (EC) member states and Australia are required to provide each prospective patient with the Mentor patient brochure titled: “Gel-Filled Breast Implant Surgery: Making an Informed Decision.” The purpose of this brochure is to assist patients in making informed decisions about breast augmentation and breast reconstructive surgery. A Patient Signature Form is provided in the back of the patient brochure for this product. For patients in the EC member states and Australia, the patient brochure must be read, understood and the Patient Signature Form must be signed by the patient prior to surgery.

It is the responsibility of the individual surgeon to decide the best method by which to counsel a patient prior to surgery. Mentor relies upon the surgeon to advise the patient of all potential complications and risks associated with the use of breast implants.

After surgery, surgeons in the EC member states and Australia are to provide the patient with the Patient ID Card with information on the implant(s) used. Labels are supplied with the implant to be applied to the Patient ID Card. Patients should be advised to carry the patient card to facilitate medical care in case of emergency.

Following surgery, the surgeon should inform the patient of necessary postoperative office visits and the need to continue to consult a physician for routine examinations to detect for breast cancer. Patients should be advised to inform a physician or surgeon of the presence of an implant if any surgery of the breast area is scheduled, and to consult a physician or pharmacist before the use of topical medicines such as steroids in the breast area. In addition, the patient should be instructed to contact her surgeon should she experience any problem related to her breast implants.

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Risks specific to breast implant surgery

Capsular Contracture

Capsular contracture is the most common side effect of breast implants. To accept the implant, a surgical pocket behind the breast is made somewhat larger than the implant itself. Normally a healing scar forms an envelope around the implant, which, on occasion, will shrink sufficiently to squeeze the implant, producing varying degrees of firmness. At its worst, the implant can feel hard, be painful and/or distorted. This can occur soon after surgery or years later and may be unilateral, bilateral or asymmetric. Surgical release or excision of the scar is often successful, but recurrence is not uncommon. The cause of the contracture phenomenon is poorly understood. In the past, closed disruption of the scar by squeezing the breast was common, but this is rarely practiced today.

Capsular contracture is graded in severity on a scale of I to IV by Baker classification. Calcification of the capsule can also occur. Calcification is a phenomenon that is occasionally seen with long-term scarring especially if there is irritation such as tight burn scars that cross-joints.

Calcified capsules may require removal if the patient wishes relief from her contracture but otherwise seem to be harmless. Small foci of calcification are commonly seen anywhere in the breast parenchyma. They can usually be identified as benign by the radiologist but on occasion may require biopsy to rule out a malignancy.

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Rupture/Deflation of the Implant

Breast implants may not last a lifetime. Although rupture can occur at any time following implantation, a number of studies examining rupture of current generation, single-lumen, CE marked, round and shaped silicone gel-filled breast implants of various manufacturers over time (using MRI evaluations) have consistently reported similar results indicating the average expected lifetime exceeds 10 years. The cumulative incidence of rupture at 5 years (identified by physician exam rather than MRI) for MENTOR Becker tissue Expanders/Breast Implants with Siltex shell (double-lumen breast implants) was 6% among 1,428 primary reconstruction patients, 9% among 490 revision reconstruction patients, and 11% among 60 revision augmentation patients (all with follow-up at 5 years).

While the silicone material itself has not been shown to biodegrade, the shell may rupture due to wear and tear, or direct injury. If the implant shell is ruptured, the escaping gel is usually contained by the scar envelope in the surgical pocket (intracapsular) and may be undetectable except by MRI (silent rupture), which is about 85% effective in detecting rupture. If the scar envelope is torn, the gel can be driven into the local tissue planes and breast tissue (extracapsular). Most of the escaped gel remains in the immediate environment of the breast but on rare occasions it has been reported to migrate down the arm, into nerve sheaths or into the abdominal wall. Ultrasound, mammography and physical examination may also diagnose these ruptures, which have escaped the scar envelope. Most of the reported cases occurred in the more fragile, thinner shell devices implanted in the late 1970s.

Rupture should be suspected if there is a change in character of the device such as a new, persistent unilateral burning sensation or a change in

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softness, texture or shape of the implant. Because of the silent nature of most ruptures and the difficulty of diagnosis without surgical exploration, the true incidence is unknown. Current products have thicker and stronger shells and more cohesive gel contents. Caution should be used when comparing expected or actual rupture rates of current devices to historical incidences, especially when, as is often the case, the brand, vintage and type of device is unknown. Explantation and/or replacement may be indicated if the implant fails, especially if it is seen in the breast parenchyma as it could be confused with or mask a tumor.

Causes of implant rupture include, but are not limited to: damage from surgical instruments, intraoperative or postoperative trauma, excessive stresses or manipulations as may occur during daily routines such as vigorous exercise, contact athletics, routine manual massage, intimate physical contact and from compression required during mammography.

Complications of Tissue Expansion

- Tissue thinning.
- Sloughing of poorly vascularized tissue.
- Closed, postoperative hematoma, manifested by enlargement, tenderness and discoloration leading, if untreated, to extrusion of the device.
- Undue pressure on the tissue located over the device or trauma to surrounding tissues which may lead to venous thrombosis, the breakdown of skin over the device and subsequent extrusion. Deflation or removal of the device may be necessary for tissue repair.

Changes in Nipple and Breast Sensation/Breast Pain

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Any surgery that undermines the skin of the breast can result in hypersensitivity or loss of sensation in the nipple-areolar region. These changes can vary in degree and may be temporary or permanent. Changes in nipple/breast sensation may, on occasion, affect sexual response or comfort while nursing. These changes are believed to be a result of nerve damage or stretching of the nerves from the surgery. There is no specific treatment for this condition.

Most women undergoing augmentation or reconstruction with a breast implant will experience some breast and/or chest pain postoperatively. While this pain normally subsides in most women as they heal after surgery, it can become a chronic problem in other women. Chronic pain can be associated with hematoma, migration, infection, and implants that are too large or capsular contracture. Sudden severe pain may be associated with implant rupture.

Interference with Mammography in Detection of Cancer/Calcium Deposits

As silicone is opaque to x-rays, an implant may theoretically interfere with the early detection of cancer by mammography as it may obscure part of the breast. Newer techniques of breast compression improve the amount of breast that can be visualized. Alternatively, most surgeons feel that the device may improve the detection of tumors by palpation. While of considerable theoretical concern, delayed detection due strictly to the presence of an implant has not been reported. Women at high risk of developing breast cancer should consider getting implants with caution. Since the breast is compressed during mammography, it is possible for the implant to rupture, but this is rare and should not deter a woman from regular, routine mammographic screening. Before the mammography exam, women should inform the technologist that they have implants.

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Calcium deposits are seen occasionally in old scars anywhere in the body and this is true of the implant capsule. This usually does not occur until years after implant surgery. Benign calcifications are also commonly seen on mammography in otherwise normal breast parenchyma even in breasts that have never been operated on. These benign calcium deposits usually have a different x-ray appearance than the calcifications that signal a malignancy. An expert radiologist can usually determine if a calcium spot is benign or malignant but occasionally a biopsy may be necessary to rule out malignancy. There is no evidence that these deposits occur more or less frequently in women with implants than those who have not received implants. After many years, some patients may develop a thin layer of calcium in the scar capsule that surrounds the implant. This is almost always associated with capsular contracture but otherwise causes no known problem.

Extrusion of the Implant/Interruption of Wound Healing

Skin necrosis and/or sloughing resulting in implant exposure may occur from undue tension of the tissues overlying the implant, the overdosage of steroids placed in the implant pocket or surgical or external trauma.

This is quite rare in the augmentation patient but occurs occasionally in the more challenging tissue environment of reconstruction after mastectomy, local tissue injury and/or radiation. The areas of the scar, especially after

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radiation to that area, seem to be the most vulnerable to disruption. Extruded implants may need to be surgically removed.

Wrinkling of the Implant/Dissatisfaction with Cosmetic Results/Asymmetry

Visible and/or palpable wrinkles in an implant are related to the thinness of the overlying tissue cover, the degree of capsular contracture and the texturing of the implant shell surface. Traditional smooth-walled gel devices rarely demonstrate wrinkles. Some have described camouflaging of wrinkling with autologous fat transfer superficial to the implant. Surgical error, preexisting asymmetry or deformity, keloid formation of the incisional scar, vagaries of time, weight gain or loss, pregnancy and nursing can all contribute to an immediate or late poor aesthetic result. With time, most breasts, with or without implants, become ptotic to some degree. Asymmetry is usually related to the inability to totally correct for preexisting disparity between the two breasts. It can also be attributed to surgical error, asymmetric contracture, or rupture of the implant.

Excessive buildup of collagen at the incision site during the healing process causes some patients to develop scars of cosmetic concern. Keloid scars, which do not respond well to treatment, often extend beyond the edges of the original scars and can continue to enlarge over time. Hypertrophic scars are generally confined to the original site and respond well to scar revision treatment, which may include steroid injections to break down the collagen or surgery to revise the position, direction or line of the scar.

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Possible Reactions to Silicone

This text contains a brief summary of information from the medical literature. Mentor recognizes that the information contained within this text is highly technical. However, medical ethics and practice dictate that the physician must be an intervening party between the manufacturer of prescriptive medical devices and the patient. The issue of the possible relationship between silicone (and other implantable materials) and various diseases has been the subject of significant scientific debate. Concerns include immunological and neurological disorders, carcinogenicity and connective tissue disorders.

Despite numerous anecdotal reports of established and newly described connective tissue diseases attributed to the device, multiple epidemiological studies have very consistently shown no significant association between silicone breast implants and any established or new connective tissue or other immunological diseases. Four prestigious multispecialty scientific expert panels have reviewed the published literature on this topic, specifically as it relates to silicone breast implants, and have issued lengthy reports on their findings. These expert panels include the Independent Review Group (commissioned by the Chief Medical Officer of the U.K.), the National Science Panel (appointed by Judge Pointer for MDL 926), the IOM and the Scientific and Technical Options Assessment (STOA) Programme (commissioned by the European Parliament). These four panels have uniformly concluded that there was no discernible evidence of causal association or positive risk ratio between exposure to silicone breast implants and recognized or new autoimmune or connective tissue diseases. Essentially the same finding was reported in a 2011 updated review of the epidemiology literature published in the peer-reviewed medical literature in 2011.

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Coincident neurological problems such as multiple sclerosis and amyotrophic lateral sclerosis (ALS) have occurred in a small number of breast implant patients. Subsequent studies have not demonstrated a relationship between any neurological disease and breast implants.

Several studies have been conducted to determine the carcinogenic risk of breast implants and no evidence of increased risk of cancer has been demonstrated. Continuing assessment of both known and possible risks associated with breast implant surgery is ongoing.

No credible reports on birth defects or other reproductive effects in humans associated with implantation of silicone breast implants of any type have been identified in the literature. Recent studies sponsored by Mentor provide further evidence that silicone materials used in breast implants do not cause adverse reproductive effects in experimental animals.

Although any breast surgery, including breast implantation, could theoretically interfere with the adequacy of a woman's milk supply, many women with breast implants have nursed their babies successfully. It is known that any breast surgery such as breast biopsy can affect the quantity of milk produced. In recent years, the question has been raised regarding the potential transfer of silicone into the breast milk of women with silicone breast prostheses and possible effects on the health of the breast-fed children. However, more recent studies have provided strong evidence of

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the lack of association between silicone breast implants and adverse effects in breast-fed children. The American Academy of Pediatrics has stated that there is no reason why a woman with implants should refrain from nursing. The European Committee on Quality Assurance and Medical Devices in Plastic Surgery (EQUAM) concluded in 2000 that silicone gel-filled breast implants do not adversely affect pregnancy, fetal development, breast-feeding, nor the health of breast-fed offspring.

The report sponsored by the IOM, “Safety of Silicone Breast Implants,” released in July 1999 states that women with silicone breast implants are no more likely than the rest of the population to develop cancer, immunologic diseases, or neurological problems. The committee also concluded that there is no evidence that mothers with implants pass silicone on to infants when breast-feeding.

Anaplastic Large Cell Lymphoma

Based on information reported to global regulatory agencies and found in medical literature, an association has been identified between breast implants and the development of anaplastic large cell lymphoma (ALCL), a type of non-Hodgkin’s lymphoma. Women with breast implants have a very small but increased risk of developing Breast Implant Associated ALCL (BIA-ALCL) in the fluid or scar capsule adjacent to the implant, with documented potential for local, regional, and distant spread of the cancer with mortality reported in rare cases.

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BIA-ALCL has been reported globally in patients with an implant history that includes Mentor's and other manufacturers' breast implants with various surface properties, styles, and shapes. Most of the cases in the literature reports describe a history of the use of textured implants. Several journal articles explore the risk factors for BIA-ALCL, including the varied methods used to create surface texture of the implant and the role of biofilm in causing disease, among others.

You should consider the possibility of BIA-ALCL when a patient presents with a late onset seroma (at least 1 year and on average 7 to 10 years following implantation), breast mass adjacent to the implant capsule or lymphadenopathy. Work-up of delayed enlargement of a previously implanted breast includes an ultrasound or MRI to look for the presence of seroma or lymphoma, Diagnosis starts at presentation, with needle aspiration of as much of the fluid as possible (minimum 50 mL), using ultrasound for implant displacement and protection. A suspicious mass requires needle or open tissue biopsy. Specimens should be sent for cell morphology by cytology, CD30 immunohistochemistry and ALK (anaplastic lymphoma kinase), and flow cytometry for evaluation, quantification, and characterization of T cells within the specimen. Hematopathology consultation at a tertiary cancer center is strongly encouraged to establish or exclude a diagnosis of BIA-ALCL. Once the diagnosis of BIA-ALCL has been established, physicians are strongly encouraged to consult with a multi-disciplinary team including oncologists, pathologists, surgical oncologists, and plastic surgeons. The United States National Comprehensive Cancer Network (NCCN) guidelines for

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treatment of established BIA ALCL includes the removal of the implant and the entire surrounding implant capsule in addition to associated masses and lymph node involvement.

Please report all confirmed cases of BIA-ALCL per the details provided in the ADVERSE EVENT REPORTING.

Please note for latest diagnostic and treatment options refer The National Comprehensive Cancer Network (NCCN) guidelines

Gel Bleed

The gel in an implant consists of large three-dimensional mesh-like molecules that constitute about 20% of the total weight of the gel. The spaces within these molecules are filled with medical grade viscous silicone fluid. This viscous fluid is similar to the materials available in many products including antifatulence medication, available without prescription for infants and adults. A small amount of this material can diffuse or bleed through the shell of the implant. The major portion of this

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material stays on the implant wall. A small amount moves into the scar capsule where it is gradually picked up by certain “scavenger” cells of the body’s immune system, called macrophages. Normally, these cells try to destroy foreign material such as bacteria. However, if the material (such as silicone) cannot be destroyed, it is carried to the lymph glands by the macrophages. As noted in the section on Possible Reactions to Silicone, the STOA report published in 2000 concluded: “Studies do not point to an association between silicone implants and serious health risks, such as cancer and connective tissue diseases.”

Granulomas

It is possible for a granuloma to form around a tiny amount of silicone. Although these lumps are noncancerous, they can be difficult to distinguish from cancerous lumps without being removed (biopsied) and examined. Granulomas, if large or suspicious for malignancy, may need to be biopsied or surgically removed and examined.

Other Possible Reactions

- Thrombosed veins, resembling large cords, have temporarily developed in the area of the prosthesis and have resolved without surgical or medical therapy.
- Pain from an improperly sized and/or placed implant, such as from compression of nerves or interference with muscle movement, may occur.
- Hypertrophic scarring has been reported.
- The prosthesis may become difficult to explant if the degree of tissue adhesion is significant.

Instructions and Precautions for Implant Removal

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due to normal body motion, unintentional sharp instrument damage during surgery, and physical stress such as that caused by car accidents or strenuous activity. Rupture rates can vary depending upon the type of breast implant surgery performed. In MENTOR clinical studies, patients undergo MRI screening on a regular basis to assess whether rupture of the implant has occurred. It is important to note that most patients had no symptoms related to their reported implant rupture.

Becker Expander/Breast Implant devices are designed with both a saline and a gel filled lumen, so the functional lifetime of the device can be impacted by deflation or rupture.

Deflation can occur if there is a tear in the device shell (inner shell) that contains the saline

Rupture can occur if there is a tear in the device shell (outer shell) that contains the gel

Suspected ruptures are ruptures that are suspected based on the MRI imaging findings, but which have not been confirmed by surgical explantation and physical examination of the implant, usually because the patient chose not to have the implant removed.

Confirmed ruptures are those verified by examining the implant once removed from the patient.

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Of the overall ruptures (those confirmed as well as those suspected but not confirmed) reported for each device, the proportion of confirmed ruptures was 56% (10 of 18 suspected or confirmed ruptures).

5-year suspected or confirmed implant rupture rates and deflation rates for Becker Expander/Breast Implants for Reconstruction Patients:

Suspected or confirmed	7.2%
Deflation	1.3%

Smaller sample sizes in the Becker Implants clinical study limit data accuracy.

Becker Expander/Breast Implant devices have data out to 5 years at this time with less than a 10% rupture/deflation rate for reconstruction patients.

- Many of the changes to a patient's preoperative breast following implantation are irreversible (cannot be undone). If a patient later chooses to have her implant(s) removed, native breast volume loss, unacceptable dimpling, puckering, wrinkling, or other cosmetic changes of the breast may be an unavoidable consequence without reimplantation.

- Breast feeding is safe following implantation however implants may affect a patient's ability to produce a sufficient supply of milk for breast-feeding. Also, breast implants will not prevent a patient's breast from sagging after pregnancy.

- With breast implants, routine screening mammography may require additional images and radiation exposure to be performed optimally

Patient ID Card

Enclosed with each gel-filled breast implant is a Patient ID Card. To complete the Patient ID Card, stick one Patient Record Label for each implant on the back of the Patient ID Card. Patient Record Labels are located on the internal product packaging attached to the label. If a Patient Record Label is unavailable, the lot number, catalog number and

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Becker 25 tissue Expander/Breast Implant with Siltex shell

- o Gel volume 25% nominal implant size.
- o Cohesive I (standard) gel filling material.
- o Indicated for temporary overexpansion.

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Becker 35 tissue Expander/Breast Implant with Siltex shell

- o Gel volume 35% nominal implant size.
- o Cohesive II (moderate) gel filling material.
- o Indicated for temporary overexpansion.

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The Mentor Becker Expander Implants contain the following main components: Outer Shell, Inner Shell, a Valve, Gel between the two shells, Laser Button, sealing Washers and Dipcoat. Textured outer shells have an additional layer that is textured. The shell shape, dimensions and volume vary and define the various styles and sizes of Becker Expander Implants. The Laser Button identifies the style and size. The Dipcoat seals the gel fill hole. The other components are identical among the Becker Expander Implants, with the gel volume varying to match the shell volume. Figure below, following, depicts the main components of the Mentor Becker Expander Implants, using the round implant as an example.


Рисунок 2– Описание компонентов МИ «Расширитель-имплантат грудной тканевый Becker 25 с оболочкой Siltex - Siltex Round Becker 25 Expander/Breast Implant Cohesive I» / Figure 2 – Description of Becker 25 tissue Expander/Breast Implant with Siltex shell - Siltex Round Becker 25 Expander/Breast Implant Cohesive components

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Table 2 – Technical characteristics of Becker 25 tissue Expander/Breast Implant with Siltex shell - Siltex Round Becker 25 Expander/Breast Implant Cohesive I

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Table 12 - TRUE-LOCK Clamp Specifications

Inner Diameter (A)	4.5 mm ± 0.1 mm
Length (B)	4.6 mm ± 0.07 mm
Outer Diameter (C)	5.5 mm ± 0.1 mm

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Table 14 - Stainless Steel Connector Specifications

Inner diameter (A)	12.7 mm ± 0.1 mm
Length (B)	11.6 mm ± 1.0 mm
Outer diameter (C)	23.1 mm ± 0.07 mm
Taper	7° ± 1°

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Technical characteristics of accessories

Accessories of the medical device (non sterile demo samples) are identical to the Becker tissue Expanders/Breast Implants with corresponding volume in terms of design features and the materials of manufacturing. The table below consists of Becker Tissue Expander/Breast Implant on the left side and the non sterile demo sample with corresponding volume on the right side.

Table 15. Compliance of medical device and its accessories

Medical device	Accessory with identical design, technical characteristics and materials
Becker 25 Tissue Expander/Breast Implant with Siltex shell - Siltex Round Becker 25 Expander/Breast Implant Cohesive I: 300 ml	Becker 25 Tissue Expander/Breast Implant with Siltex shell - Siltex Round Becker 25 Expander/Breast Implant Cohesive I Demo - non sterile demo sample, volume 300 ml.
Becker 25 Tissue Expander/Breast Implant with Siltex shell - Siltex Round Becker 25 Expander/Breast Implant Cohesive I: 350 ml	Becker 25 Tissue Expander/Breast Implant with Siltex shell - Siltex Round Becker 25 Expander/Breast Implant Cohesive I Demo - non sterile demo sample, volume 350 ml.
Becker 35 Tissue Expander/Breast Implant with Siltex shell - Siltex Contour Profile Becker 35 Expander/Breast Implant Cohesive II: 325 ml	Siltex Contour Profile Becker 35 Expander/Breast Implant - Siltex Contour Profile Becker 35 Expander/Breast Implant Cohesive II Demo - non sterile demo sample, volume 325 ml.
Becker 35 Tissue Expander/Breast Implant with Siltex shell - Siltex Contour Profile Becker 35 Expander/Breast Implant Cohesive II: 365 ml	Siltex Contour Profile Becker 35 Expander/Breast Implant - Siltex Contour Profile Becker 35 Expander/Breast Implant Cohesive II Demo - non sterile demo sample, volume 365 ml.

с оболочкой Siltex - Siltex Contour Profile Becker 35 Expander/Breast Implant Cohesive II: 365 мл	Siltex - Siltex Contour Profile Becker 35 Expander/Breast Implant Cohesive II Demo - нестерильный демонстрационный образец, объем 365 мл	Becker 35 Tissue Expander/Breast Implant with Siltex shell - Siltex Contour Profile Becker 35 Expander/Breast Implant Cohesive II: 400 ml	Siltex Contour Profile Becker 35 Expander/Breast Implant - Siltex Contour Profile Becker 35 Expander/Breast Implant Cohesive II Demo - non sterile demo sample, volume 400 ml.
Расширитель-имплантат грудной тканевый Becker 35 с оболочкой Siltex - Siltex Contour Profile Becker 35 Expander/Breast Implant Cohesive II: 400 мл	Расширитель-имплантат грудной тканевый Becker 35 с оболочкой Siltex - Siltex Contour Profile Becker 35 Expander/Breast Implant Cohesive II Demo - нестерильный демонстрационный образец, объем 400 мл		

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








DATA ABOUT LABELING OF MD AND ITS PACKAGING

Table 17-Definitions of symbols used on the label

Left Breast	Left Breast is the location of the implanted breast implant
Right Breast	Right Breast is the location of the implanted breast implant
Date of Implant	Date of the implant surgery
Patient Name	Patient Name
Surgeon Name	Surgeon Name
Address	Surgeon's Address
Phone	Surgeon's Telephone Number
GTIN	Global Trade Item Number
REF	Catalogue number
LOT	Batch code
QTY 1	Quantity One
STERILE !	Sterilized using steam or dry heat

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	Sterilized using irradiation
	Sterilized using ethylene oxide
	Serial number
	Use-by date
	Caution
	Do not reuse
	Diameter, Projection
	Width, Height, Projection
	Manufacturer
	CE-mark
	Recycling code 2 High-density polyethylene

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- Only one prosthesis should be implanted per breast. Mentor recommends against the stacking of implants, one upon the other. The devices have not yet been tested for this use and the integrity of the implants cannot be guaranteed as the materials may abrade and wear. Such abnormal stress may result in weakening or deflation/rupture of the prostheses.
- Do not insert or attempt to repair a damaged or altered prosthesis.
- The action of drugs (examples: antibiotics and steroids) in contact with the prosthesis has not been tested by the manufacturer and their use cannot be recommended. Each physician who chooses to use chemotherapeutic drugs with this prosthesis must assure compatibility of the drug with the silicone elastomer.
- In vitro testing has demonstrated that even low concentrations of povidone-iodine solution placed within the breast implant will compromise implant integrity in the long term. Therefore, we recommend that no povidone-iodine solution or other antibacterial, antiseptic or cleaning agent be added to the injection media. If a cleaning solution is to be used within the implant space, the site should be carefully rinsed to remove the residual solution.
- Do not introduce or make injections of drugs or other substances into the implant. Injections through the implant shell will compromise the product's integrity, causing it to leak while in use and eventually deflate and/or rupture.
- Preoperative evaluation of the implant design, size, and implant site should include allowances for adequate tissue coverage. Pressure, force, tension and other stresses to which the implant site will be susceptible must be considered.
- Excessive inflation of the device may result in tissue necrosis/thrombosis.
- Final Expansion Volume should not be less than the Minimum Recommended Volume or more than the Maximum Recommended Volume. Underfilled prostheses may buckle, fold or wrinkle causing crease/fold failure of the device and subsequent deflation and/or rupture can occur. Inflation beyond the Maximum Recommended Volume may also cause crease/fold failure or shell rupture.
- Sepsis, hemorrhage or thrombosis may result from the placement of any foreign object in the body.
- The use of microwave diathermy in patients with breast implants has been reported to cause tissue necrosis, skin erosion and extrusion of the implant. Its use in patients with breast implants is not recommended.

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- The patient should be made aware that any abnormal stress or trauma to the breast could result in rupture of the prosthesis. Women should be advised to see their physician immediately if they suspect that their implant has ruptured.

- Mentor strongly recommends against the use of forceful external stress (such as closed capsulotomy) to treat capsule firmness. Mentor is not responsible for the structural integrity of the implant should the surgeon elect to perform such a procedure. If the physician uses this technique, several complications may occur: hematoma, displacement of the implant and/or shell rupture. The physician should inform the patient of these potential complications and of alternatives to the procedure. Capsule firmness must not be treated by overexpansion of the device. Such abnormal stress or trauma to the breast and the prosthesis could result in rupture of the prosthesis.

- The American College of Radiology has stated that mammography may be less effective on implanted breasts and may interfere with early detection of breast cancer. The mammographer should be trained and experienced with the most current radiologic techniques and equipment. This may increase cost and radiation exposure to the patient. Patients should inform the mammographer that they have breast implants and should also be instructed how to distinguish the prosthesis from normal or abnormal breast tissue during self-examinations for breast cancer.

- Mentor has not tested the in vivo effects of radiation therapy on tissues of patients who have breast implants; however, the literature suggests that radiation therapy may increase the likelihood of experiencing capsular contracture.³ In addition, the literature notes the following about the effects of radiation therapy on implanted breasts: "(a) when the implanted breast was free of fibrotic changes, radiotherapy produced acceptable results, (b) whenever feasible, 45 Gy/5 weeks seemed preferable over higher doses, (c) irradiation immediately after the reconstructive surgery appeared to produce poorer cosmetic results."⁴ The decision regarding the use of radiation therapy following breast implantation should be made by the surgeon and the radiation oncologist.

- The use of the periumbilical approach in placement of the implant has not been established and is not recommended.

- Careful hemostasis is important to prevent postoperative hematoma formation. Should excessive bleeding persist, it is recommended that the device not be implanted until the bleeding

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is controlled.

- If a physician treats a hematoma or serous fluid accumulation by aspiration, or if a biopsy or lumpectomy is performed, care must be taken to avoid damaging the implant. These procedures present possible risk of implant puncture.

- The incidence of extrusion of the prosthesis has been shown to increase when the prosthesis has been placed in injured areas: scarred, heavily irradiated or burned tissue, crushed bone areas or where severe surgical reduction of the area has been previously performed.

- Excessive fibrous capsular formation or contracture may occur around any implant placed in contact with soft tissues. The incidence and severity of this occurrence may increase if postoperative local hematoma or injection occurs.

- The physician should use personal discretion when deciding to use these prostheses in patients who exhibit psychological instability.

- Surgical implantation of a breast implant may interfere with the ability to breast-feed. The Institute of Medicine concluded that there is limited evidence that implantation, especially through a peri-areolar incision, may interfere with lactation and breast-feeding. However, it should be noted that previous breast reconstruction surgery, such as mastectomy, may be the initial case of this interference.

- Patients with connective tissue disease (CTD) may experience an increased risk of wound dehiscence, infection and bleeding (likely due to their ongoing medical management) that may require further treatment.⁸

- Some patients with breast implants have reported experiences of neurological and/or rheumatological diseases. Mentor is currently collecting data to further understand these potential risks and their possible association with breast implants.

- A patient should seek evaluation by a plastic surgeon to evaluate for BIA-ALCL if they experience delayed onset seroma in the breast (delayed means developing more than 12 months after implantation), a new breast mass or lymphadenopathy.

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Precautions

It is the responsibility of the surgeon to advise prospective patients or their representatives, prior to surgery, of the possible complications associated with the use of this product.

- Preexisting infection should be treated and resolved before implantation of the prosthesis.

- It is possible that bubbles may form in the silicone gel as a result of the manufacturing or sterilization process. These bubbles will not detract from the safety or efficacy of the prosthesis, and will diffuse and dissipate of their own accord.

- Any surgeon performing augmentation or reconstructive mammaplasty with implants should be familiar with the currently available techniques for measuring the patient, determining implant size and performing surgery.

- The dual self-sealing valve of the Becker tissue Expanders/Breast Implants with Siltex shell family of prostheses is unique and may be unfamiliar to the surgeon. The fill tube is inserted into the prosthesis at the time of manufacture and should be handled carefully to avoid accidental dislodgment from its prepositioned location. Do not hold the device by its fill tube.

- The silicone elastomer shell, fill tube and injection dome may be easily cut by a scalpel or ruptured by excessive stress, manipulation with blunt instruments or penetration by a needle. Subsequent deflation and/or rupture will result. All prostheses should be carefully inspected for structural integrity prior to and during implantation.

- Meticulous care must be exercised in handling, connecting and implanting the device.

- Any subsequent surgical procedures in the area of the implant should be undertaken with extreme caution as damage to the implant could occur. In the event that an implant is damaged, it must be removed.

- Each device should be checked for patency prior to surgery and continuously monitored throughout the surgical procedure to ensure the structural integrity of the device is not compromised in any way. This prosthesis should not be implanted following any modifications to its original design. A prosthesis which has been damaged, or on which repairs or modifications have been attempted, should not be implanted. A standby prosthesis should be available at the time of surgery.

- When removing the fill tube and injection dome, the fill tube should be removed first. Grasp the fill tube beyond the connector to prevent separation of

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the injection dome from the fill tube. Do not exert sudden or undue tension on the fill tube during removal. Avoid instrument damage to the fill tube which may result in tube breakage, retraction of tube into the pocket and subsequent deflation and/or rupture of the device.

- Tissue ingrowth can occur when using the TRUE-LOCK Connector. Surgeons should anticipate the need to dissect the capsule prior to removing the fill tube and injection dome. Grasp beyond the connector and remove the fill tube before taking out the injection dome.

- Do not contact the device with disposable, capacitor-type cautery instruments as damage to the outer shell of the prosthesis may result.

- The fill tube which connects the implant to the injection dome should be carefully sized to avoid kinks. Careful attachment of the fill tube to the connector is important to prevent separation. Failure of the device to inflate may be due to kinking of the fill tube, leakage, separation of the components or injections which do not penetrate the injection dome.

- Extreme care should be taken when connecting the fill tube to the connector. The tube is easily damaged with surgical instrumentation (e.g., forceps contact), and their use should be avoided.

- Surgeons should ensure themselves of the position of the injection dome prior to adding or withdrawing fluid.

- Potential for contamination exists when fluid is added to or removed from the device. Use aseptic technique in the introduction of saline into the implant; a single-use, sterile saline container is recommended.

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Precautions for winged catheter use:

Use care not to puncture the expansion envelope.

- Replace set per CDC guidelines

1. Remove end protector cap from female luer of Winged Infusion Set.

2. Connect female luer of Winged infusion Set to male luer of delivery system

3. Remove protective sheath from needle of the Winged Infusion Set and insert needle of the Winged Infusion Set and insert needle into injection reservoir according to the instructions for the tissue expander.

Additional Precautions for Becker tissue Expanders/Breast Implants with Siltex shell

- Avoid too small an incision. A larger incision than is normally used for smooth-shelled implants may be required to facilitate insertion and to avoid damage to the device. A device which is damaged during insertion may result in postoperative deflation and/or rupture.

- Mentor recommends the surgeon consider the size of the implant and the firmer nature and higher profile of the SILTEX shell when choosing optimum incision size and surgical approach. Certain surgical approaches may cause higher stresses on the device during implantation.

Peri-operative Infection and Inflammation Source Reduction Practices

The following practices could potentially reduce peri-operative bacterial and particulate based infection and inflammation:

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- Perform a layered closure to separate skin/wound bacteria from the implant cavity.

- Avoid, where possible, pocket dissection that communicates with breast ducts and glands

- Avoid unnecessary adjacent tissue trauma and maintain meticulous hemostasis in the creation of the implant pocket.

- Minimize intra-operative device exposure, skin contact, and pocket placement trials.

- Administer preoperative intravenous antibiotic prophylaxis.

- Lint, dust, talc, surgical glove powder, drape and sponge lint, fingerprints, skin oils and other surface contaminants deposited on the device by improper handling may cause foreign body reactions. Strict adherence to clean, aseptic techniques should be maintained to prevent contamination and possible complications. Surgical instruments and gloves should be rinsed clean of any impurities before handling the device.

MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION

Non-clinical testing has demonstrated the Mentor with Injection Domes and Stainless Steel Connection Systems (provided with the Becker Expander/Implant) are MR Conditional. A patient with this device can be safely scanned in an MRI system meeting the following conditions:

- Static magnetic field of 3.0 T or 1.5 T
- Maximum spatial field gradient of 1900 gauss/cm (19 T/m)
- Maximum MRI system reported, whole body averaged specific absorption rate (SAR) of 4.0 W/kg (Normal or first level mode)

Under the scan conditions defined above, non-clinical testing results indicate the Mentor and Becker implant with Injection Domes and Stainless Steel

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guaranteed. Furthermore, the integrity of the device cannot be guaranteed due to the risk of damage to the implant. The established shelf life of the products is compromised and thus null and void if compliance with the single-use only indication is not followed. Sterility, safety and efficacy cannot be assured for damaged devices. In the event the product becomes contaminated, contact your local Mentor representative.

For more detailed information see Validation Sterilization Report № 100479294.

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REQUIREMENTS FOR TECHNICAL SUPPORT AND REPAIR OF MD

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CONDITIONS OF TRANSPORTATION AND STORAGE

Table 18 – Transportation conditions

Designation	Range
Temperature	-20 °C to +50°C
Relative humidity	not more 75%
Barometric pressure	700-1060 hPa

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Table 19 – Storage conditions

Designation	Range
Temperature	+10°C to +25°C
Relative humidity	not more 75%
Barometric pressure	700-1060 hPa

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MD DISPOSAL RULES AND CONDITIONS

Unused products which had no contact with the internal environment,

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RA-JNJRU-RZN@ITS.JNJ.com

MANUFACTURER:

Mentor Medical Systems B.V., Zernikedreef 2, 2333 CL Leiden, The Netherlands.

MANUFACTURERS SITES:

Mentor Medical Systems B.V., Zernikedreef 2, 2333 CL Leiden, The Netherlands.

B.Braun Medical Inc., 824 Twelfth Avenue, Bethlehem, PA 18018, USA

Mentor, 3041 Skyway Circle North Irving, TX 75038, USA

**AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN
RUSSIA:**

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designed to join the preinserted or pre-attached fill tube of the device to the injection dome.

The Injection Domes and Connector Systems are supplied in a sterile and nonpyrogenic double-wrap packaging system. The package is sterilized by gamma radiation. The double-wrap system facilitates the preferred method of sterile product transfer from the circulating area to the sterile field. Sterility cannot be guaranteed if the double-wrap packaging system has been damaged. Sterility, safety and efficacy cannot be assured if the packaging has been opened or the devices have been damaged.

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1. The TRUE-LOCK system consists of a straight connector with clamps. The connector holder is utilized to facilitate placement of the tubing into the cage of the connector. The TRUE-LOCK connector does not require a suture tie.

2. The Stainless Steel Connector (see figure 14) does require suture material tied around the tube and connector to secure the connection.

INSTRUCTIONS FOR USE

Stainless Steel Connector Systems

The physician may trim the tubing to size and then join it to the injection dome using the connector supplied. Care should be taken to tailor the tubing length so that the tubing will not kink or shorten. Nonabsorbable sutures should be tied around the tubing and connector (Figure 5) to secure the connection.


Figure 14 - Securely tie both tubings with nonabsorbable sutures.

It is important to securely suture the fill tube both distally and proximally to the connector so the entire fill tube assembly will be removed when the injection dome is removed from the patient. Care must be taken to secure the tubing to the connector with ligatures in such a manner as to avoid cutting or occluding the tubing or connector. The use of forceps or hemostats to aid in the connection and suture tying process is contraindicated as tubing or connector damage may lead to deflation of the device.

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TRUE-LOCK connector system

The TRUE-LOCK Connection System consists of a straight connector with clamps. The connector holder is utilized to facilitate placement of the tubing into the cage of the connector.

Injection Domes

Micro Dome

1. The Micro Injection Dome is intended for use when diminished dermal palpability and projection of a fill dome are required. It is suggested that the Micro Injection Dome be placed in a relatively superficial location to allow ease of identification and access during subsequent filling procedures. Inflation is accomplished by using pyrogen-free, sterile sodium chloride solution for injection. Use a 21-gauge (or finer) standard or butterfly 12° bevel needle. Extreme care should be taken to puncture only the center of the top surface of the Micro Injection Dome at an angle perpendicular $\pm 30^\circ$ to the top surface (Figure 15). This dome is designed to withstand up to 10 total injections.

Standard Dome

2. The Standard Injection Dome is larger in diameter and height than the Micro Injection Dome. It is suggested that the Standard Injection Dome be placed in a relatively superficial location to allow ease of identification and access during subsequent filling procedures. Inflation is accomplished by using pyrogen-free, sterile sodium chloride solution for injection. Use a 21-gauge (or finer) standard or butterfly 12° bevel needle. Extreme care should be taken to puncture only the center of the top surface of the Standard Injection Dome at an angle perpendicular $\pm 30^\circ$ to the top surface (Figure 3). This dome is designed to withstand up to 20 total injections.

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
Figure 15.


Figure 16.

1. Trim tubes to apply clamp Trim tubes squarely and slide on clamp with tapered end first:

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
Figure 17.

2. Slide connector onto the connector holder.


Figure 18.

3. Insert tube into the connector cage and advance tube to the center wall. (This should be verified visually.)

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
Figure 19.

4. Slide clamp over connector until it also meets the center wall. Remove connector holder.


Figure 20.

5. Place second clamp on tube and insert tube into connector. Advance tube to the center wall. Slide second clamp over connector until it also meets the center wall.

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
Figure 21.

For more information, see the instructions for use of the medical device «Becker tissue Expanders/Breast Implants with Siltex shell in the following versions».

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How supplied

Stainless steel and TRUE-LOCK connectors are double packaged sterile. The first sterile packaging is a DuPont Tyvek Uncoated bag Tyvek1073B (DuPont), a Mylar film (DuPont). The second sterile layer (general bag containing sterile packaging of connectors and injection ports) is a Tyvek Uncoated 1073B bag (manufactured by DuPont), sealed with Mylar branded (manufactured by DuPont)

The injection ports are packaged in sterile packaging - Tyvek Uncoated 1073B (DuPont) sealed with Mylar (DuPont).

Packaging material of Winged infusion set and Safsite injection Site and Cap – Surlyn/Eva Film (manufacturer Amcor) and medical grade paper Uncoated 1059B Breathable Tyvek (manufacturer DuPont).

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