Sworn to and subscribed
before me this

26th day of March, 2021

MATEUSZ F. LESZCZAK
NOTARY PUBLIC OF NEW JERSEY
Commission # 50085913
My Commission Expires 08/10/2022

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INSTRUCTION OF USE FOR Medical Device

Surgical suture unidirectional with needle for knotless
tissue fixation STRATAFIX™ Spiral PDS™ Plus
in configurations

*Version for the realization on the territory of
Russian Federation*

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NAME OF PRODUCT

Surgical suture unidirectional with needle for knotless tissue fixation STRATAFIX™ Spiral PDS™ Plus in configurations:

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DESIGN OWNER

ETHICON, INC., PO Box 151, Route 22 West, Somerville, New Jersey 08876-0151, USA

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INTENDED USE

The STRATAFIX™ Spiral PDS™ Plus Device is indicated for use in soft tissue approximation where the use of absorbable sutures is appropriate.

USAGE CONDITION

This medical device is intended for professional use in the medical and health care institutions.

AREA of USE

Surgery

POTENTIAL CONSUMER

Surgical Suture STRATAFIX™ Spiral PDS™ Plus is intended to be used by professional doctor, who is aware about surgical procedures, techniques and materials, as well as how to work with absorbable surgical materials.

INDICATIONS

The STRATAFIX™ Spiral PDS™ Plus Device is indicated for use in soft tissue approximation where the use of absorbable sutures is appropriate.

CONTRAINDICATIONS

The STRATAFIX™ Spiral PDS™ Plus Device is not to be used where prolonged (beyond six weeks) approximation of tissues under stress is required and is not to

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be used in conjunction with or for fixation of prosthetic devices (e.g., heart valves or synthetic grafts) that are nonabsorbable in nature. This device, being absorbable, should not be used where extended approximation of tissue under stress is required, such as in fascia. The STRATAFIX™ Spiral PDS™ Plus Device should not be used in patients with known allergic reactions to IRGACARE MP (triclosan).

WARNINGS

- The fixation loop of the STRATAFIX™ Spiral PDS™ Plus Device should be deployed by users familiar with surgical procedures, techniques and tissues and the use of absorbable devices before employing the STRATAFIX™ Spiral PDS™ Plus Device for wound closure, as risk of wound separation may vary with the site of application and the suture material used.
- Physicians should consider the quantity and quality of tissue in which the loop will be anchored.
- The use of this device may be inappropriate in highly vascularized tissue or fragile tissue that cannot withstand potential further tightening/constriction within the loop.
- Care should be undertaken to avoid leaving barbed suture ends adjacent to the peritoneum in extra-peritoneal tissue closure.
- Avoid contacting STRATAFIX™ Device and associated needles with other materials (e.g. surgical gauze, drapes, etc.) in the surgical field to prevent ensnaring on the barbs. If the barbs catch, carefully pull the material in the opposite direction of the needle to disengage it from the barbs.
- Physicians should consider the in vivo performance (under ACTIONS section) when selecting a device for use in patients. The use of this device may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing.
- As with any foreign body, prolonged contact of any device with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation.
- As an absorbable device, STRATAFIX™ Spiral PDS™ Plus may act transiently as a foreign body.
- Acceptable surgical practice should be followed for the management of contaminated or infected wounds.
- Small bowel obstruction (SBO): including volvulus, bowel infarction, and significant morbidity, have been reported due to barbs or barbed suture ends

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hooking onto adjacent small bowel and/or mesentery, such as in peritoneal closure.

- As this is an absorbable device material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distention or which may require additional support.
- The safety and effectiveness of the STRATAFIX™ Spiral PDS™ Plus Device have not been established for use in fascial closures (including abdominal wall, thoracic and extremity fascial closures), gastrointestinal anastomoses, cardiovascular tissue, neural tissue, osseous tissue, tendinous tissue, ophthalmic surgery or for use in microsurgery.
- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.
- Discard open packages and unused devices. Discard packaging that is damaged or opened.
- Do not use after expiration date.

PRECAUTIONS

- The STRATAFIX™ Spiral PDS™ Plus Device contains a loop end and barbs to anchor tissues and does not require knots to secure the device. Tying of knots on the barbed section of the device will damage the barbs and potentially reduce their effectiveness. For the device to function properly, the STRATAFIX™ Spiral PDS™ Plus Device must first be anchored in robust tissue using the fixation loop. Then subsequently engage the barbs into the tissue in a standard closure pattern. When completing placement of the barbed segment in subcuticular tissue, take at least one additional pass in the reverse direction to lock the STRATAFIX™ Spiral PDS™ Plus Device in place. In all other tissue layers, take at least two passes in the reverse direction to lock the STRATAFIX™ Spiral PDS™ Plus Device in place.

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- When using the STRATAFIX™ Spiral PDS™ Plus Device subcutaneously, the device should be placed as deeply as possible in order to minimize erythema and induration normally associated with absorption.
- Care should be taken to avoid damage when handling. Avoid crushing or crimping the suture material with surgical instruments, such as needle holders and forceps. Do not attempt to remove memory in the polymer by running fingers down the suture material as this can damage the barbs.
- Infections, erythema, foreign body reactions, transient inflammatory reactions and in rare instances wound separation are typical or foreseeable risks associated with any suture and hence are also potential complications associated with the STRATAFIX™ Spiral PDS™ Plus Device.
- Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds. To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.
- Discard used needles in designated biohazard "sharps" containers.

ADVERSE REACTIONS

Adverse events associated with wound closure (including wounds closed with synthetic absorbable sutures), include wound separation or dehiscence, failure to provide adequate wound support in closure of sites where expansion, stretching, or distention occurs, or in elderly, malnourished or debilitated patients or in patients suffering from conditions that may delay wound healing, infection, minimal acute inflammatory tissue reaction, device extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site, as well as allergic reaction to IRGACARE MP (triclosan). Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

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MRI SAFETY INFORMATION

STRATAFIX™ Spiral PDS™ Plus Device is MR Safe.

DESCRIPTION

The STRATAFIX™ Spiral PDS™ Plus Device is an antibacterial monofilament, synthetic absorbable device consisting of dyed (violet) polyester, poly(p-dioxanone), the empirical molecular formula of which is $(C_4H_6O_3)_x$.

The device contains IRGACARE MP (triclosan), a broad spectrum antibacterial agent, at no more than 2360 µg/m. The pigment for the violet dye is D&C Violet No. 2. Polydioxanone polymer has been found to be nonallergenic, nonpyrogenic and elicits only a slight tissue reaction during absorption.

The STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device, Variable Loop Design consists of barbed suture material, armed with a surgical needle on one end and a fixation loop at the opposite end. The STRATAFIX™ Spiral PDS™ Plus Device is designed to anchor with a closed loop at one end and a unidirectional barbed section on the other end. The STRATAFIX™ Spiral PDS™ Plus Device barbs are oriented in one direction to allow tissue approximation without the need to tie surgical knots. While the formation of barbs in the STRATAFIX™ Spiral PDS™ Plus Device reduces the tensile strength relative to non-barbed suture material of the same size, tying of knots in non-barbed suture materials also reduces their effective strength. For this reason, the strength of the STRATAFIX™ Spiral PDS™ Plus Device can be compared to USP knot strength of non-barbed sutures. Additionally, USP designations for diameter are used to describe the STRATAFIX™ Spiral PDS™ Plus Device suture material after barbing, except for minor variation in suture diameter with a maximum overage of 0.1 mm.

The USP and EU Pharmacopoeia sizes of the STRATAFIX™ Spiral PDS™ (poly[p-dioxanone]) Plus Device are further defined in Table 1.

Table 1. The USP and EU Pharmacopoeia sizes of the STRATAFIX™ Spiral PDS™ Plus Device

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TENSILE STRENGTH

STRATAFIX™ Spiral PDS™ Plus Device straight tensile strength is stronger than the knot tensile strength for a USP and EU Pharmacopoeia polydioxanone device of the equivalent size as shown in Table 2.

Table 2. Tensile Strength Comparison

USP DEVICE SIZE DESIGNATION	EU PHARMACOPOEIA DEVICE SIZE (Metric/Ph.Eur.) DESIGNATION	Suture Minimum Knot Tensile Strength	
		USP (kgf)	EU Ph. (N)
1	4	5,08	50,8
0	3,5	3,90	39,0
2-0	3	2,68	26,8
3-0	2	1,77	17,5
4-0	1,5	0,95	9,5

ACTIONS

The results of implantation studies in animals using the STRATAFIX™ Spiral PDS™ Plus Device indicate that for sizes 3-0 and larger, approximately 80% of its original strength remains after two weeks and four weeks of implantation. At six weeks post implantation, approximately 40% to 70% of its original strength is retained. For size 4-0, approximately 67% of its original strength remains after two weeks, and 50% at four weeks of implantation. At six weeks post implantation, approximately 37% of its original strength is retained.

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Data obtained from implantation studies show that absorption of the STRATAFIX™ Spiral PDS™ Plus Device is essentially complete within 210 days (Table 3).

Table 3. Percentage preservation of the suture strength from the original depending on the number of days after implantation

DAYS IMPLANTATION	APPROXIMATE ORIGINAL STRENGTH REMAINING (%) Sizes 3-0 and larger	APPROXIMATE ORIGINAL STRENGTH REMAINING (%) Size 4-0
14 days	80%	67%
28 days	80%	50%
42 days	40-70%	37%

Using in vitro studies, STRATAFIX™ Spiral PDS™ Plus Device with IRGACARE MP (triclosan) has been shown to inhibit colonization of the device by *Staphylococcus aureus*, *Staphylococcus epidermidis*, methicillin-resistant *S. aureus*, methicillin-resistant *S. epidermidis*, *Escherichia coli*, and *Klebsiella pneumoniae*, which are microorganisms known to be associated with surgical site infections. Animal studies have demonstrated that STRATAFIX™ Spiral PDS™ Plus Device inhibits bacterial colonization of the device after direct in vivo challenge with bacteria.

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RISK ANALYSES

The evaluation of risk management was conducted in accordance with the International Standards Organization EN ISO 14971:2012, entitled Medical Devices – Application of risk management to medical devices.

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LIST OF MATERIALS ANIMALISTIC or HUMAN ORIGATION

STRATAFIX™ Spiral PDS™ Plus Suture doesn't contain any materials animalistic and /or human origination

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INSTRUCTION FOR USE

STRATAFIX™ Spiral PDS™ Plus Suture is designed to be used in continuous suture patterns. STRATAFIX™ Spiral PDS™ Plus Device is intended to be used without anchoring knots to begin or terminate the device line.

Use as required per surgical procedure.

1. Secure the fixation loop portion to robust tissue by taking a bite in the designated tissue, passing the needle through the loop and pulling tautly around the anchoring tissue.
2. After STRATAFIX™ Spiral PDS™ Plus Device is anchored into the tissue, grasp the strand and approximate tissue to the desired tension. Bites, or passes through the tissue, can be taken in a continuous

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suturing technique to approximate the wound. Care should be taken to utilize the STRATAFIX™ Spiral PDS™ Plus Device on the barbed segments only. Do not attempt to approximate wounds using the non-barbed segment near the needle, as the barbs are required for successful wound approximation with the STRATAFIX™ Spiral PDS™ Plus Device.

3. а) To complete the closure in subcuticular closure, take at least one pass in the reverse direction. Then pass the needle across the incision and take a split thickness bite perpendicular to the incision and exit the skin.

б). To complete the closure in all other tissue layers, take at least two passes in the reverse direction to complete the closure.

4. Then gently pull on the free end of the device and cut flush with the surface of the tissue.

To remove the entire STRATAFIX™ Spiral PDS™ Plus Device: Cut the STRATAFIX™ Spiral PDS™ Plus Device at the end anchored with the loop, and then pull the distal end to remove the device.

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MATERIALS**Suture**

Polydioxanone (C₄H₆O₃)_x IRGACARE MP dyed with D&C Violet No. 2 60725, Suture material manufactured by Ethicon, Inc., Cornelia, Georgia, USA
 Antibacterial agent - triclosan IRGACARE MP, manufacturer BASF Grenzach GmbH, Germany, GMP certificate №DE_BW_01_GMP_2018_0014.

Needle

Stainless steel AISI grade 420*, Supplied by Carpenter Technology Corp. Reading, PA, USA or
 Stainless steel AISI grade 4310*, Supplied by Zapp Precision Metals GmbH Schwerte, Germany or
 Ethalloy Stainless steel*, manufactured by ATI Specialty Materials, USA.
 Needle coverage – Silicone (CAS Number 01330-20-7 63148-57-2), Supplied by Nusil Technology, USA or
 Multipass (polydimethylsiloxane silicone with the addition of wax powders), Supplied by MicroPowders Inc., USA and Nusil Technology, USA or
 Cerberus, manufactured by Gelest Inc., USA or Siltech MED-ETI (Siltech Corporation, Canada).

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ACCESSORIES

The device in this file does not have accessories included.

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STERILIZATION

The STRATAFIX™ Spiral PDS™ Plus Device is sterilized by ethylene oxide gas. Do not resterilize. Do not use if package is opened or damaged. Discard opened, unused sutures. Do not use after expiration date.

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HOW SUPPLIED

The STRATAFIX™ Spiral PDS™ Plus Device is available sterile. Each STRATAFIX™ Spiral PDS™ Plus Device is single armed with various needle sizes. Device is placed into a plastic tray packaging which is placed into peel open foil package. Twelve (12) pouches and an IFU are placed into a carton box. The boxes are wrapped with a film and tear tape is placed during the process to ease the opening of the boxes.

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CRITERIA OF UNSUITABLE CASES OF MEDICAL DEVICES USE

Do not use of the packaging is open or damaged. Discard open packages and unused devices.

EXPIRATION DATE

Shelf life of the product is 24 months.

DEIMPLANTATION

N/A, Absorbable Suture.

In case suture removal follow the 4 paragraph of Instruction of Use section of the current Instruction.

STORAGE, TRANSPORTATION AND OPERATION CONDITIONS

The operation conditions in the rooms where the medical device is used must meet the following requirements:

Temperature from -10°C to 50°C (75% non-condensing humidity equivalent to 30°C) when utilized immediately after opening sterile packaging.

Storage: -60°C до +30°C

transportation conditions: -60°C до +50°C

air humidity is N/A (packaged in moisture proof foil).

WARRANTY

Manufacturer guarantees to the end user absence of defects and conformity to product specifications for the entire shelf life of the product under the conditions of adherence to storage, transportation and package integrity requirements. There are no warranties that extend beyond the description on the face hereof. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use.

RECYCLING

Discard open packaging and unused devices.

Discard used needles in "sharps" containers.

Dispose the product and its residual elements or waste items in accordance with local laws and regulations.

After being used the products are potentially infectious waste.

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The device has no negative impact on humans and the environment during use, transportation and storage.

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