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JD

True Copy

I hereby certify that the following is a true copy of the original.


Hamburg, 29. November 2021

W

Heiko Zier, notary public

Approved by

Sarah Stöttinger
Senior Regulatory Affairs Specialist
Global Brand, Innovation and Sustainability
BSN medical GmbH

 **essity**
BSN medical GmbH
Schützenstraße 1-3
20253 Hamburg, Germany

 **essity**
BSN medical GmbH
Quickbornstraße 24
20253 Hamburg, Germany

BSNmedical
an Essity company
BSN medical GmbH
Quickbornstr. 24
20253 Hamburg

Instruction for use medical devices

Medical adhesive tapes

Manufactured by
BSN medical GmbH, Germany

1. Name of MD.

Medical adhesive tapes

Models:

Leukopor[®], Leukoflex[®], Leukofix[®]

Available sizes: 1,25 cm x 5 m, 2,5 cm x 5 m, 5 cm x 5 m

Leukosilk[®], Leukoplast[®]

Available sizes: 1,25 cm x 4,6 m, 2,5 cm x 4,6 m, 5 cm x 4,6 m

Leukoplast[®] skin sensitive

Available sizes: 1,25 cm x 2,6 m, 2,5 cm x 2,6 m

2. Manufacturer and developer: BSN medical GmbH

Postal address:

Quickbornstrasse 24

20253 Hamburg

Germany

Juristic address:

Schützenstrasse 1-3

22761 Hamburg

Germany

Production site:

BSN medical GmbH

Heykenaukamp 10

21147 Hamburg

Germany

For Leukoplast[®]

BSN medical SAS

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CS 90022

72320 Vibraye

France

3. Intended use of MD

Medical fixing adhesive tapes are intended for short-term use on intact skin for fixing dressings of all types and sizes, as well as for fixing tubes, catheters, probes, cannulae, electrodes and other measuring devices.

Unsterile, 1st risk class MD.

4. Indications for use of MD

Fixing dressings of all types and sizes, tubes, catheters, probes, cannulas, electrodes and other measuring devices.

5. Con

6. Adv

Beside
risks a

7. Pre

Risk m
toward
intende
Leuko

8. Inf

Medic

9. Fun

Leuko

Test

Adhes

Air pe

MVTI

Antim

Sorpti

Functi

Water

Thick

Surfac

// MV

Leuko

Test

Adhes

Air pe

MVTI

Antim

Sorpti

Functi

Water

Thick

Surfac

5. Contraindications

Leukoplast®: Known intolerance to latex or colophony

6. Adverse reactions

Besides the possibility of intolerance to latex or colophony concerning the Leukoplast® product, no risks and adverse events are known for the self-adhesive medical tapes.

7. Precautions

Risk management revealed two risks requiring evaluation: Skin reactions (e.g. by allergic reactions towards latex or colophony or macerations caused by high humidity) and wrong usage (e.g. not intended use).

Leukoplast®: Tape contains natural rubber latex and colophony.

8. Information about potential consumers of MD

Medical adhesive tapes are used by healthcare professionals and patients

9. Functional characteristics of MD and indications for use

Leukosilk®

Test	Limit	Test method
Adhesion	1.5 – 2.9 N/cm or 150 - 290 N/m	JBS 40000083
Air permeability	N/A	N/A
MVTR	>300 g/(m ² *24h	JBS.40000101
Antimicrobial activity	N/A	N/A
Sorption capacity of wound pad	N/A	N/A
Functional pad wetting time, s	N/A	N/A
Water vapor permeability	N/A	N/A
Thickness	0,21 ± 0,04 mm	DIN EN ISO 9073-2
Surface density	166 ± 25 g/m ²	DIN EN 14079

// MVTR = Water vapor permeability//

Leukopor®

Test	Limit	Test method
Adhesion to steel	1.2 – 2.6 N/cm or 120 - 260 N/m	JBS 40000083
Air permeability	> 1 [cm ³ /cm ² /s]	JBS.40000050
MVTR	>500 g/(m ² *24h	JBS.40000101
Antimicrobial activity	N/A	N/A
Sorption capacity of wound pad	N/A	N/A
Functional pad wetting time, s	N/A	N/A
Water vapor permeability	N/A	N/A
Thickness	0,16 ± 0,04 mm	DIN EN ISO 9073-2
Surface density	80 ± 10 g/m ²	DIN EN 14079

Leukofix®

Test	Limit	Test method
Adhesion to steel	1.6 – 3 N/cm or 160 - 300 N/m	JBS 40000083
Air permeability	> 5 [cm ³ /cm ² /s]	JBS.40000050
MVTR	>500 g/(m ² *24h	JBS.40000101
Antimicrobial activity	N/A	N/A
Sorption capacity of wound pad	N/A	N/A
Functional pad wetting time, s	N/A	N/A
Water vapor permeability	> 500 g/m ² 24h	JBS 40000101
Thickness	0,17 ± 0,04mm	DIN 53370
Surface density	130 ± 25 g/m ²	DIN EN 14079

Leukoflex®

Test	Limit	Test method
Adhesion to steel	1.5 – 3.1 N/cm or 150 - 310 N/m	JBS 40000083
Air permeability	N/A	N/A
Antimicrobial activity	N/A	N/A
Sorption capacity of wound pad	N/A	N/A
Functional pad wetting time, s	N/A	N/A
Water vapor permeability	N/A	N/A
Thickness	0.1 mm +/- 0.02 mm	DIN 53370
Surface density	95 ± 25 g/m ²	DIN EN 14079

Leukoplast® skin sensitive

Test	Limit	Test method
Adhesion to steel	0.15 – 0.6 N/cm or 15 - 60 N/m	JBS 40000083
Air permeability	N/A	N/A
MVTR	>300 g/(m ² *24h	JBS.40000101
Antimicrobial activity	N/A	N/A
Sorption capacity of wound pad	N/A	N/A
Functional pad wetting time, s	N/A	N/A
Water vapor permeability	N/A	N/A
Thickness	0,27 ± 0,07mm	DIN EN ISO 9073-2
Surface density	250 ± 50 g/m ²	DIN EN 14079

Leukoplast®

Test	Limit	Test method
Adhesion to steel	2.5 – 4 N/cm 250 - 400 N/m	JBS 40000083
Air permeability	>1 [cm ³ /cm ² /s]	JBS.40000050
MVTR	>500 g/(m ² *24h)	JBS.40000101
Antimicrobial activity	N/A	N/A
Sorption capacity of wound pad	N/A	N/A
Functional pad wetting time, s	N/A	N/A

Water vapor permeability	N/A	N/A
Thickness	0.32 mm +/- 0.064 mm	DIN EN ISO 9073-2
Surface density	220 ± 50 g/m ²	DIN EN 14079

MVTR: Moisture Vapour Transmission Rate

11. Technical characteristics of MD

Measurement point, mm	Model of MD		
	Leukopor [®] , Leukofix [®] , Leukoflex [®] 1,25 cm x 5m	Leukopor [®] , Leukofix [®] , Leukoflex [®] 2,5 cm x 5m	Leukopor [®] , Leukofix [®] , Leukoflex [®] 5 cm x 5m
Width of MD, mm	12,5 ± 1,5	25 ± 1,5	50 ± 1,5
Length of MD, m	5000 ± 150	5000 ± 150	5000 ± 150

Measurement point, mm	Model of MD		
	Leukosilk [®] , Leukoplast [®] 1,25 cm x 4,6 m	Leukosilk [®] , Leukoplast [®] 2,5 cm x 4,6 m	Leukosilk [®] , Leukoplast [®] 5 cm x 4,6 m
Width of MD, mm	12,5 ± 1,5	25 ± 1,5	50 ± 1,5
Length of MD, m	4600 ± 150	4600 ± 150	4600 ± 150

Measurement point, mm	Model of MD	
	Leukoplast [®] skin sensitive 1,25 cm x 2,6 m	Leukoplast [®] skin sensitive 2,5 cm x 2,6 m
Width of MD, mm	12,5 ± 1,5	25 ± 1,5
Length of MD, m	2600 ± 150	2600 ± 150

12. How to apply MD

Unroll a strip and cut it off.

Then fix the dressing or medical device with it.

13. Handling conditions

Environmental conditions	
Temperature (°C)	+15°C - + 30°C
Humidity (%)	Max. 65%

14. Labeling of the MD and consumer pack.

Medical device contains:

- name of MD
- name and address of manufacturer
- Lot number
- expiration date

Graphical symbols on the MD

Symbol	Explanation
	Product conformity symbol for quality requirements in the EU
	Products meet the requirements of environmental standards.

Snapring contains:

- name of MD
- name and address of manufacturer
- Article number (REF)
- dimensions of MD
- barcode
- «Latex or Latex free» (for Leucoplast/Leukoplast skin sensitive)
- logo «essity»

Graphical symbols on snapring of MD

Symbol	Explanation
	Products meet the requirements of environmental standards.
	Product conformity symbol for quality requirements in the EU
	Avoid exposure to sunlight (for Leukoplast skin sensitive and Leukosilk)
	Keep away from moisture (for Leukoplast skin sensitive and Leukosilk)
	Manufacturer
REF	Article number
	Medical device
	EU registered trade mark

Graphical symbols on folding box (consumer group pack)

Symbol	Explanation
	Product conformity symbol for quality requirements in the EU
	Manufacturer
	Expiration date

15. I

Med
direc
20°C
The
envi
is re

Here

Stor
Temp
Hum

16. S
For 1
expir
in ac

For 1
the p
cond

17. U
Med
biolo
epid
acco

18. M

Med
patie

	Batch code
	Products meet the requirements of environmental standards

15. Information about transportation of MD and storage conditions

Medical adhesive tapes are recommended to store in a dry and clean place without exposure to direct sunlight. The temperature during transport shall not be higher than 30°C and fall not under -20°C for one week.

The product must be transported in transport packaging that protects the product from damage, in an environment that does not allow excessive moisture and direct sunlight. No special pressure regime is required during storage and transportation of the product.

Here are storage ranges:

Storage conditions	
Temperature (°C)	+15°C - +30°C
Humidity (%)	Max. 65%

16. Shelf life

For medical adhesive tapes (Leukosilk®, Leukoflex®, Leukofix®, Leukopor® and Leukoplast®) the expiration date is set within 5 years from the production date indicated on the package, when stored in accordance with the specified conditions.

For medical adhesive tape Leukoplast® skin sensitive the expiration date is set within 3 years from the production date indicated on the package, when stored in accordance with the specified conditions.

17. Utilization

Medical adhesive tapes are intended for single use only. After use, the products may pose a potential biological hazard and must be disposed of in accordance with legal requirements as epidemiologically hazardous waste (hazard class B). Handling and disposal must be performed in accordance with good medical practice and applicable local, state, and federal laws and regulations.

18. Manufacturer's warranty

Medical adhesive tapes warranty covers material, function and workmanship when used on one patient only. Reuse is prohibited.

19. Authorized representative in Russia

Authorized representative in Russia:

LLC Essity
Krzhizhanovskogo str.14 bld.3.
117218 Moscow
Russia

Manufacturer and developer: BSN medical GmbH

Postal address:
Quickbornstrasse 24
20253 Hamburg
Germany

Juristic address:
Schützenstrasse 1-3
22761 Hamburg
Germany

Production sites:

BSN medical GmbH
Heykenaukamp 10
21147 Hamburg
Germany

BSN medical SAS
Rue du Millenaire
CS 90022
72320 Vibraye
France

20. Lis

1. ISO

2. DIN

3. EN I

4. DIN

5. DIN

6. ISO

7. MEI

20. List of applicable national standards

1. ISO 13485:2016	Medical devices - Quality management systems – Systems requirements for regulatory purposes This international standard specifies requirements for a quality management system that can be used by an organization for the design and development, production, installation and servicing of medical devices, and the design, development, and provision of related services.
2. DIN EN ISO 14971	Risk management This international standard is a framework for effective management by the manufacturer of the risks associated with the use of medical devices. It includes a procedure to investigate, using available information, the safety of a medical device, including accessories, by identifying hazards and estimating the risks associated with the device.
3. EN ISO 10993-1	Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management.
4. DIN EN 62366-1	Medical devices – Part 1: Application of usability engineering to medical devices
5. DIN EN 1041	Information supplied by the manufacturer of medical devices This standard specifies the information to be supplied by a manufacturer for different categories of medical devices, as required by the relevant EEC Directives. It does not specify the language to be used for such information.
6. ISO 15223-1	Medical devices – Symbols to be used with medical device labels, labeling and information supplied. This international standard identifies symbols conventionally used to convey information essential for proper use to the user and others for safe and effective use of medical devices.
7. MEDDEV 2.7.1 rev. 4	Guidelines on Medical Devices, Clinical Evaluation

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