APPROVED BY:

General Manager of TOP GLOVE SDN. BHD.


(Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru,
41050 Klang, Selangor D. E., Malaysia)Tel.: + 603-3392 1992, sales@topglove.com.my

/ Noor Akilah bt Saidin /

MF.

INSTRUCTION FOR USE FOR
MEDICAL DEVICE:Medical diagnostic (examination) disposable non-sterile
gloves BENOBY® specialized

NOTARISED BY

CHOUNG YOUN SIN
NOTARY PUBLIC
41050 KLANG,
SELANGOR.

15 APR 2022

*Initial issue of operational documentation (instruction for use for) medical device: 01/11/2020*

2022

TABLE OF CONTENTS

1 Data on the medical device	3
1.1 Name of the medical device	3
1.2 Manufacturer (producer) of the medical device	3
1.3 Purpose of the medical device	3
1.4 Scope of the medical device	3
1.5 Application conditions of the medical device	3
1.6 Intended users of the medical device	4
1.7 Indications for use of the medical device	4
1.8 Contraindications for use of the medical device	4
1.9 Operating principle of the medical device	4
1.10 Shelf life of the medical device	6
2 Instruction for use of the medical device	7
3 Warning and precautions	8
4 Transportation, storage and disposal conditions of the medical device	9
4.1 Transportation conditions of the medical device	9
4.2 Storage conditions of the medical device	9
4.3 Disposal conditions of the medical device	9
5 Cleaning and disinfection methods of the medical device	9
6 Specifications of the medical device	10
6.1 Description of the medical device	10
6.2 Physical specifications of the medical device	11
6.2.1 Physical and mechanical performance of the medical device	11
6.2.2 Inspection levels and acceptable quality level	12
7 Information on completeness and packaging of the medical device	13
8 Scope of delivery of the medical device	13
9 List of national standards, applied by manufacturer (producer) of the medical device	14
10 Possible adverse reactions	15
11 Circumstances, in which the consumer must consult healthcare professional	15
12 Manufacturer's warranty	15
13 Data on Authorized representative of the manufacturer	15

1 Data on the medical device

1.1 Name of the medical device

Medical diagnostic (examination) disposable non-sterile gloves BENOBY® specialized,
in hardware versions:

1. Medical diagnostic (examination) disposable non-sterile gloves BENOBY® Dual Color nitrile powder free chlorinated finger textured two-colored with two-layer structure for additional protection against punctures and cuts (beaded cuff), color: black&orange, orange&black, size: XS, S, M, L, XL.
2. Medical diagnostic (examination) disposable non-sterile gloves BENOBY® Grip Plus nitrile powder free chlorinated finger textured for better grip (beaded cuff), color: blue, black, size: XS, S, M, L, XL.
3. Medical diagnostic (examination) disposable non-sterile gloves BENOBY® Wave Grip nitrile powder free chlorinated with wave texture over the entire surface for better grip (beaded cuff/elongated beaded cuff), color: blue, black, orange, size: XS, S, M, L, XL.
4. Medical diagnostic (examination) disposable non-sterile gloves BENOBY® Chemotherapy nitrile powder free chlorinated finger textured resistant to chemical and cytostatic agents/for chemotherapy (beaded cuff), color: sky-blue, cornflower, size: XS, S, M, L, XL.
5. Medical diagnostic (examination) disposable non-sterile gloves BENOBY® Chemotherapy nitrile powder free chlorinated high density textured resistant to chemical and cytostatic agents/for chemotherapy (elongated beaded cuff), color: blue, size: XS, S, M, L, XL.

1.2 Manufacturer (producer) of the medical device

Manufacturer (producer) of the medical device: TOP GLOVE SMD. BHD., Malaysia, Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor D.E., Malaysia.

Address of the manufacturing facility of the medical device: Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor D.E., Malaysia.

1.3 Purpose of the medical device

The gloves are intended for single use in protecting medical personnel's hands and creating a barrier against bacterial and/or viral contamination between the health worker and the patient during examination/treatment, as well as in case of possible allergies caused by latex and powder.

1.4 Scope of the medical device

Gloves are used during medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces, providing hand hygiene, as also protection of patient and healthcare professional from cross-contamination.

1.5 Application conditions of the medical device

Medical disposable gloves are used in medical or healthcare facilities indoors in a dry place, away from direct sunlight, at temperature from 10°C to 30°C and a relative humidity up to 80%, in accordance with the instruction for use on medical device.

Gloves for external use only (the device is non-invasive).

1.6 Intended users of the medical device

Medical personnel, who carry out medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces and other activities, which require compliance with the requirements for protection of patient/healthcare professional from cross-infection.

1.7 Indications for use of the medical device

The indications for use of the medical device are the necessity of conduction of medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces.

1.8 Contraindications for use of the medical device

The contraindication to the medical device is a presence of allergic reaction to gloves made from synthetic latex (nitrile) in a patient or a medical worker.

1.9 Operating principle of the medical device

The basis of nitrile diagnostic disposable gloves limits the movement of infectious substances, increases the effectiveness of providing antibacterial protection in a medical facility by protecting a medical worker from possible infections when examining / treating a patient, as well as contributes to the reliable capture of medical instruments and increases the safety of work.

Due to the differences in the method of surface treatment, its texture, glove thickness, cuff length, and additional additives in the latex, the medical device allows you to capture almost all areas of diagnostics and therapy, including dentistry, medical, sanitary, laboratory and other medical procedures and researches.

All hardware versions of the medical device are powder-free - chlorinated. The gloves are treated with chlorine, which makes the surfaces of the medical device smooth and eases putting on and removing, as well as prevents the gloves from sticking together.

Dual Color Gloves. The gloves have a two-layered structure, they are made of two-layered nitrile, the outer layer of the glove has a contrasting color in relation to the inner layer of a glove, in order to provide the visual control of the integrity of a glove before putting it on, when working with the known infected patients and the patients at risk.

Grip Plus Gloves. The Grip Plus technology, due to the sticky texture on the fingers, allows the use of gloves when working with thin and small medical instruments, as well as in cases when it is extremely important to hold the instrument in the hand, in order not to let it move.

The improved features of the textured surface on the fingers provides the better holding of tools, both dry and wet ones.

Description of Grip Plus Technology. For the production of gloves, a related ceramic mould of hands is used (see Image № 1), its fingers' areas are coated with a special compound which makes the fingertips of the gloves become shiny and sticky (see image № 2).

The coating composition consists of alumina oxide (CAS № 1344-28-1), silicon dioxide (CAS № 14808-60-7), manufactured by Concord Chemicals Corporation Sdn. Bhd., Malaysia.

Image № 1.

Image № 2.

Wave Grip Gloves. Their significant wavy texture throughout the glove surface provides the perfect grip in wet and oily environments. The extended cuff option provides the improved protection of a person's wrist and forearm against possible contamination from infectious medical materials.

Chemotherapy Gloves (fourth (4) hardware version). The gloves for chemotherapeutic manipulations are resistant to the permeability of chemotherapy drugs (cytotoxic agents). The gloves have been tested for cytotoxic permeability in accordance with standard ASTM D 6978 (the results are presented in Table 1).

Table № 1.

№	Name of Chemotherapeutic Drug	Concentration, mg/ml	Detection Time for Diffusion, min.
1	Carmustin (BCNU)	3.3	Average time of 15.47 min.
2	Cisplatin	1.0	No break up to 240 min.
3	Cyclophosphamide (Cytosan)	20.0	No break up to 240 min.
4	Dacarbazine (ABVD)	10.0	No break up to 240 min.
5	Doxorubicin hydrochloride	2.0	No break up to 240 min.
6	Etoposide (Toposar)	20.0	No break up to 240 min.
7	Fluorouracil	50.0	No break up to 240 min.
8	Paclitaxel (Taxol)	6.0	No break up to 240 min.
9	Thiotepa	10.0	Average time of 3.39 min.
10	Methotrexate	25.0	No break up to 240 min.
11	Mitomycin C	0.5	No break up to 240 min.
12	Vincristine sulfate	1.0	No break up to 240 min.

Gloves Chemotherapy (fifth (5) hardware version). The gloves for chemotherapeutic manipulations are resistant to the permeability of chemotherapy drugs (cytotoxic agents). The gloves have been tested for cytotoxic permeability in accordance with standard ASTM D 6978 (the results are presented in Table 1).

Besides, the gloves possess the feature of high density (High Density). They are characterized by the increased wall thickness and increased resistance to damage. The high density gloves are used in conditions of high infection risk (risk of contamination), while performing the procedures where there is a risk of mechanical damage to the gloves, thus, the reliable hand protection is required.

The extended cuff option provides the improved protection of a person's wrist and forearm against possible contamination from infectious medical materials.

Table № 2.

№	Name of Chemotherapeutic Drug	Concentration, mg/ml	Detection Time for Diffusion, min.
1	Carmustin (BCNU)	3.3	Average time of 31.9 min.
2	Cisplatin	1.0	No break up to 240 min.
3	Cyclophosphamide (Cytosan)	20.0	No break up to 240 min.
4	Dacarbazine (ABVD)	10.0	No break up to 240 min.
5	Doxorubicin hydrochloride	2.0	No break up to 240 min.
6	Etoposide (Toposar)	20.0	No break up to 240 min.
7	Fluorouracil	50.0	No break up to 240 min.
8	Paclitaxel (Taxol)	6.0	No break up to 240 min.
9	Thiotepa	10.0	Average time of 136.2 min.

1.10 Shelf life of the medical device

The established shelf life of the medical nitrile diagnostic gloves is 5 years.

2 Instruction for use of the medical device

Preliminary hygienic treatment of your hands is required before use of the gloves (treat hands with skin antiseptic or wash with soap and water). Gloves are put on dry hands only.

How to put the gloves on:

- 1) Take a glove from the box. It is recommended to touch the bounded surface of the glove corresponding to the wrist (upper edge of the cuff).
- 2) Pull on the first glove.
- 3) Take the second glove with your bare hand, touching the bounded surface of the glove corresponding to the wrist.
- 4) To avoid touching the skin of the forearm with a gloved hand, rotate the outer surface of the glove to be put on with folded fingers in the glove so that you could pull the glove over the other hand.
- 5) Having put the gloves on, it is required to touch only those objects provided for by the indications for use and conditions of the use of the gloves.

How to take the gloves off:

- 1) To take a glove off, pick it up at the level of the wrist without touching the skin of the forearm and pull it off your hand, turning it inside out.
- 2) Holding the removed glove with the gloved hand, it is necessary to place the fingers of your bare hand between the glove and wrist. Remove the second glove by rolling it down the hand and put into the first glove.
- 3) After use, the gloves should be disposed of, in accordance with the recommendations related to the disposal of medical waste. Wash the hands thoroughly.

3 Warning and precautions

- The device is intended for single use;
- The medical device is non-sterile and does not require any sterilization;
- For individual use only;
- Every time it is necessary to perform hygienic treatment of your hands before putting on or changing gloves, as well as after removing the gloves;
- It is necessary to choose the gloves of the right size – they should be easily put on and should not impede movement or squeeze the hands. Do not use in case of size mismatch;
- It is recommended to touch only those areas of the glove that are designed for wrists (top edge of the cuff);
- It is necessary to change the gloves in case of any defect, puncture, damage, contamination, fluid contact with skin under the glove, as well as stickiness;
- If one glove is damaged, change both ones;
- Depending on the type of the performed manipulation, to keep an optimal level of protection, it is recommended to change gloves every 60-120 minutes;
- Do not wet and keep away from moisture;
- It is necessary to put the gloves on dry hands only;
- Avoid exposure to direct sunlight;
- Gloves do not guarantee 100% hand protection against punctures, cuts and damage caused by exposure to chemicals;
- It is recommended to consult healthcare professional in the event of dermal hypersensitivity reactions when using the device.

4 Transportation, storage and disposal conditions of the medical device

4.1 Transportation conditions of the medical device

The transportation of the device may be made by all means of transport (sea, land, railway or air transport) according to the rules of transportation, which are applied to such means of transport. During transportation, the storage conditions should not be violated.

4.2 Storage conditions of the medical device

The devices should be stored in a dry place protected from direct sunlight, within a temperature range from 10°C to 30°C and a relative humidity up to 80%. The devices should be stored away from the sources of ignition.

4.3 Disposal conditions of the medical device

Disposal of the medical device in the Russian Federation is carried out by the medical organization as class B waste in accordance with SanPiN 2.1.3684-21 'Sanitary and epidemiological requirements for the maintenance of urban and rural settlements, for water facilities, drinking water and drinking water supply, atmospheric air, soils, living quarters, the operation of industrial, public premises, the organization and implementation of sanitary and anti-epidemic (preventive) measures' and approved guideline of the organization.

The unused medical gloves or gloves with the expired shelf life should be disposed as class A waste in accordance with SanPiN 2.1.3684-21 'Sanitary and epidemiological requirements for the maintenance of urban and rural settlements, for water facilities, drinking water and drinking water supply, atmospheric air, soils, living quarters, the operation of industrial, public premises, the organization and implementation of sanitary and anti-epidemic (preventive) measures' and approved guideline of the organization.

5 Cleaning and disinfection methods of the medical device

Disposable non-sterile examination gloves are a disposable medical device. Reuse of the medical device is a gross violation of application methods.

6 Specifications of the medical device

6.1 Description of the medical device

Form of the medical device. Examination gloves have straight fingers and are ambidextrous (non-anatomic).

The glove's cuff may be elongated (regardless of size, the length of the gloves is $300 \text{ mm} \pm 10 \text{ mm}$) or may have a standard length (length indicators are presented in Table № 3). Rolled cuff.

Treatment of the surface of the medical device. Medical gloves are powder-free. Powder-free gloves undergo only one treatment, i.e. chlorination with chlorine.

Texture of the external surface of the medical device. Depending on the version, the gloves have rough textured surface on the fingers or along the entire surface, as well as the significant wave texture throughout the entire glove surface.

Sizes of the medical device. Depending on the size, Table 3 gives the values of width, length and thickness of gloves.

Table 3 – Sizes of the medical device.

Value name	Size				
	XS	S	M	L	XL
Width W (Figure 1), mm	≤ 80	80±10	95±10	110±10	≥110
Length l (Figure 1), mm, min	220		230		
Thickness (at points shown in Figure 2), mm, min					
- smooth surface					
- textured surface					
Thickness of the cuff, mm, max	2,50				


Figure 1 – Points for measuring the width and length of the glove


Figure 2 – Points for measuring the wall thickness of the glove

6.2 Physical specifications of the medical device

6.2.1 Physical and mechanical performance of the medical device

Physical and mechanical parameters of the medical device are represented in Tables 4, 5.

Table 4 - Values of physical and mechanical characteristics of the medical device before accelerated ageing.

Value of physical and mechanical characteristics before accelerated ageing:		
Elongation at break, %, min	Force at break, N, min	Conditional strength, MPa, min
500	7,0	14

Table 5 - Values of physical and mechanical characteristics of the medical device after accelerated ageing.

Value of physical and mechanical characteristics after accelerated ageing		
Elongation at break, %, min	Force at break, N, min	Conditional strength, MPa, min
400	6,0	14

6.2.2 Inspection levels and acceptable quality level

For the purposes of quality control medical gloves are selected and examined in accordance with standard ISO 2859-1.

Table 6 – Control Levels and Acceptable Quality Level (AQL).

SPECIFICATION	CONTROL LEVEL	AQL
Dimensions (width, length, thickness)	S-2	4,0
Impermeability	G-1	2,5
Elongation at break and force at break (before and after accelerated ageing)	S-2	4,0

7 Information on completeness and packaging of the medical device

The package provides protection from the exposure to mechanical and climatic factors during transportation and storage.

Consumer packaging. Gloves shall be packed in boxes made of thick cardboard. Consumer packaging of the gloves provides protection from external mechanical damage.

Depending on the scope of delivery, the consumer packaging shall contain a certain number of gloves of the same size, one or several colors within the same version. When packing gloves of the same version in two or more colors, a package with an internal partition shall be used in order to separate the gloves along with several perforations for separate carrying of these gloves.

Transportation package. Gloves in consumer package are packed in boxes made of corrugated board. Carton boxes are strong enough and provide safety of devices during transportation and storage.

8 Scope of delivery of the medical device

Gloves can be delivered in the following scopes of supply:

- 2000 gloves, i.e. 1000 pairs in a consumer package;
- 400 gloves, i.e. 200 pairs in a consumer package;
- 200 gloves, i.e. 100 pairs in a consumer package;
- 100 gloves, i.e. 50 pairs in a consumer package;
- 50 gloves, i.e. 25 pairs in a consumer package;
- 40 gloves, i.e. 20 pairs in a consumer package;
- 30 gloves, i.e. 15 pairs in a consumer package;
- 20 gloves, i.e. 10 pairs in a consumer package;
- 10 gloves, i.e. 5 pairs in a consumer package;
- 4 gloves, i.e. 2 pairs in a consumer package;
- 2 gloves, i.e. 1 pairs in a consumer package.

9 List of national standards, applied by manufacturer (producer) of the medical device

Table 7 – List of standards applied by the Manufacturer.

Index	Name
ISO 9001:2008	Quality Management Systems. Requirements.
ISO 13485:2016	Medical Devices. Quality Management Systems. Requirements for Regulatory Purposes.
ISO 14971:2007	Medical Devices - Application of Risk Management to Medical Devices.
ISO 15223-1:2016	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. Part 1. General Requirements.
ISO 11193-1	Single-Use Medical Examination Gloves. Part 1. Specification for Gloves Made of Rubber Latex or Rubber Solution.
ISO 2859-1	Statistical Methods. Sampling procedures for Inspection by Attributes.
ASTM D 6319	Standard Specification for Nitrile Examination Gloves for Medical Application.
ASTM D 6978	Uniform Procedure for Assessing the Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

10 Possible adverse reactions

Healthcare providers may observe developments of type IV allergies associated with an allergic and sensitizing reaction to chemical components of gloves made from synthetic latex (nitrile). Moreover, the manifestation of dermatitis and other skin irritations is possible as a result of improper hands skincare and harsh conditions of hand skin isolation from natural air ventilation. If skin irritation occurs, stop using gloves and consult a doctor.

11 Circumstances, in which the consumer must consult healthcare professional

The patient is recommended to consult healthcare professional in the event of dermal hypersensitivity reactions, when using the device.

12 Manufacturer's warranty

The Manufacturer (producer) guarantees preservation of all the parameters and characteristics of the medical device taking into account Acceptable Quality Level (AQL), as well as the quality of the medical gloves within all shelf life (5 years) from the date of manufacture upon the condition of package integrity, observance of transportation, storage and operation conditions in accordance with the instruction.

13 Data on Authorized representative of the manufacturer

Please contact the Authorized representative of the manufacturer for any and all purposes, related to circulation of the medical device on the territory of the Russian Federation:

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