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GETINGE *

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Manufacturer:
Maquet Cardiovascular LLC
45 Barbour Pond Drive Wayne, NJ 07470, USA

Susan Eichler-Huston

SUSAN EICHLER-HUSTON
NOTARY PUBLIC
STATE OF NEW JERSEY
MY COMMISSION EXPIRES APRIL 18, 2027

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2022
New Jersey, USA


Figure 1.
AXIUS Blower/Mister

1. Atraumatic Tip
2. Malleable Shaft
3. Flow Volume Controller
4. Hand Piece
5. IV Luer Lock Connection
6. Pinch Clamp
7. Roller Clamp
8. IV Spike and Chamber
9. Filter
10. Braided Tubing

INSTRUCTIONS FOR USE

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions.

DESCRIPTION:

The AXIUS Blower/Mister is comprised of a handpiece and a malleable shaft which enables positioning of the device near the targeted site. Connectors on the proximal end of the device allow for connection to standard gas and saline sources. The device permits both an adjustable flow of gas and a mist of saline to be deposited at the targeted site.

INDICATIONS:

The AXIUS Blower/Mister is intended to clear an anastomotic site for improved visibility.

WARNINGS AND PRECAUTIONS:

- **Single Use Only.** Do not re-sterilize. Single use medical devices can have extremely small and/or structurally complex components. It can be difficult to remove blood, tissue, body fluids and pathogens from the device. In addition, the cleaning and re-sterilization process may compromise the materials and functionality of the device.
- Contains Phthalates, Bis(2-ethylhexyl) Phthalate (DEHP). Device is not intended for use in children or pregnant or lactating women.
- Do not use if the product sterilization barrier or its packaging is compromised.
- DO NOT set flow at a rate greater than 8 L/min (0.13 L/s).
- Use caution when moving the tip of the AXIUS Blower/Mister closer than 3 cm to the surgical site. Do not allow tip to contact tissue.
- DO NOT USE OXYGEN WITH THIS DEVICE.
- This product is for single use only. DO NOT RESTERILIZE. DO NOT REUSE.
- When removing the device from the package and preparing the device for use take care not to stretch or damage between the Hand Piece (Figure 1;#4) and the IV Luer Lock Connection (Figure 1;#5). Doing so may damage the saline line resulting in poor or no flow of saline to the anastomotic site.

INSTRUCTIONS FOR USE:

1. Connect the braided (0.64 cm) tubing to a regulated source of CO₂. It is recommended that the regulator be set to no more than 50 psi (344 kPa). Initially, adjust the flow control to 5 L/min (0.08 L/s). DO NOT EXCEED 8 L/min (0.13 L/s).
2. Close the pinch clamp on the IV tubing. Aseptically connect the IV spike to a new bag of sterile saline and place in a pressure cuff. Inflate the cuff to approximately 150 mm Hg (20 kPa). Open the pinch clamp and adjust sterile saline flow to 1–5 (mL/min) using the roller clamp on the IV tubing.
3. Adjust the irrigation mist as needed by making adjustments to either the roller clamp on the IV tubing or to the flow volume controller on the AXIUS Blower/Mister hand piece.
4. Bend the malleable shaft to the desired shape to easily access the surgical site.
5. Hold the tip 5 cm–15 cm from the surgical site to be visualized. Adjust the distance as needed to clear the site. Use caution when moving the tip closer than 3 cm to the surgical site. DO NOT ALLOW TIP TO CONTACT TISSUE.


HOW SUPPLIED:

The AXIUS Blower/Mister is supplied **STERILE** and **NON-PYROGENIC** in an unopened and undamaged package. The AXIUS Blower/Mister has been sterilized using gamma irradiation and is for single use only. **DO NOT RESTERILIZE. DO NOT REUSE.**

***CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.*

WARRANTY:

MAQUET Cardiovascular LLC (MAQUET) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond MAQUET's control directly affect the instrument and the results obtained from its use. MAQUET's obligation under this warranty is limited to the repair or replacement of this instrument and MAQUET shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. MAQUET neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **MAQUET assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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MAQUET is a Trademark of MAQUET GmbH & Co. KG
All rights reserved.