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


INSTRUCTION FOR USE

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices

Revision 2

2019

1. Name of the medical devices

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices, versions:

1. Rolls EUROSTERIL without folds (50 mm x 200 m, 75 mm x 200 m, 100 mm x 200 m, 150 mm x 200 m, 200 mm x 200 m, 250 mm x 200 m, 300 mm x 200 m, 400mmx200m).
2. Rolls EUROSTERIL gusseted (200 mm x 50 mm (fold) x 100 m, 250mmx50mm (fold) x100m, 250 mm x65 mm (fold) x100m).
3. Pouches EUROSTERIL (50 mmx250mm, 75mm x 250 mm, 100mm x250 mm)
4. Pouches EUROSTERIL (Self-adhesive) (90 x 250 mm, 140 x 260 mm, 190 x 330 mm).

Hereinafter – devices, rolls and/or pouches EUROSTERIL.

2. Designer, manufacturer and authorized representative information

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Manufacturer's authorized representative on the territory of the Russian Federation

«PROTECO», JSC

Legal address 192102, Saint-Petersburg, Salova street, 57, building1, letter "I" and 196128, Saint-Petersburg, Varshavskaya street, 5/2, office 401, business center "Varshavsky"

Phone +7 (812) 635-88-90; Email: info@protecodent.ru

INSTRUCTION FOR USE

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices

3. Intended purpose, indications and potential consumers of the medical devices

Rolls and pouches EUROSTERIL are intended for packaging of medical devices that are subject to sterilization in order to preserve their sterility until opening or during the established shelf life.

Rolls and pouches EUROSTERIL are recommended for the following consumers:

- Personnel in health care institutions, medical and preventive treatment institutions of various profiles and accreditation levels, including obstetric in-patient departments, physiotherapy departments, maternity hospitals, children and day in-patient departments, polyclinics, dental clinics and offices, hospitals, ambulatories, dispensaries, medical and obstetric centers, organs transplantation centers, medical units and medical stations, first-aid stations, anatomic pathology departments, health resorts, recreation resorts, rehabilitation centers, social welfare institutions, specialized medical centers, donor centers and blood transfusion centers, clinical, microbiological, biochemical, bacteriological, virological, serological and other laboratories, etc.;
- Personnel of companies manufacturing sterile medical devices; pharmaceutical, chemical, microbiological industries;
- Personnel of other epidemiologically significant facilities the activity of which requires medical devices sterilization and further storage in sterile condition in accordance with the existing hygienic and sanitary regulations and anti-epidemic rules, normative legal documents.

4. Contraindications, precautions, side effects associated with the intended use of the medical devices

No contraindications.

Read instructions for use prior to application. Do not use expired devices. Use of expired devices is prohibited.

Rolls and pouches EUROSTERIL are classified as single use devices!

Prior to packaging, rolls or pouches shall be inspected for integrity. Use of damaged devices is prohibited!

Availability of sterilization indicator and its conformity to the color specified in the instruction shall be checked. Pouches without an indicator or with an indicator of non-conforming color are not permitted for use.

Sterilization of packed devices is performed in the sterilizers approved for use according to the procedure established in the Russian Federation, in accordance with the modes

INSTRUCTION FOR USE

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices

regulated by the existing documents, as well as according to the operation manual for a certain-type sterilizer.

Possible side effects associated with medical device intended use are not registered.

5. General information, design and operation principle.

Assortment of products manufactured by EURONDA S.p.A. is represented in the form of combined rolls (with and without folds) and pouches, including self-adhesive pouches, of different sizes (see table 1), resistant to steam and gas sterilization (with the use of ethylene oxide).

Rolls without folds and pouches for sterilization are intended for packaging of small-thickness devices. Rolls with folds are used for larger devices.

Rolls and pouches are made of transparent blue film and paper. Polyester and polypropylene film (transparent side) and water-repellant medical paper (non-transparent side) are connected by 3-strip thermal seam for increased strength (2-strip seam for rolls 50 mm x 200 m and pouch EUROSTERIL 50mm×250 mm). The sealed edge of a pouch has a cutout for fingers, which facilitates opening of packages when taking sterilized devices out from them. Sealed zone is free of any texts in order to prevent printing ink penetration into the package. Heat-resistant glue with 0.07 mm thickness closed with protective paper strip is applied on the open edge of self-adhesive pouches, on the internal side of the protruding paper part. Rolls and pouches are sealed by sealing or thermo sealing machines of different types, permitted to use in a prescribed manner according to instructions for use for every type of machine in particular.

Pouch corners are additionally sealed in order to prevent accumulated dust in the open package. Closing of devices is not difficult, which allows reducing medical staff expenses for packaging. 1st class indicators are applied on the paper side of the bag, they change color upon steam or ethylene oxide sterilization (see the section "Indicators"), which allows to distinguish between processed and unprocessed devices.

Rolls and pouches EUROSTERIL:

- Permeable for relevant sterilizing agents, allowing to sterilize packaged devices;
- Non-permeable for microorganisms when the rules of package closing, storage conditions and dates are observed;
- Preserve integrity (including tightness of seams) and appearance (except for indicator color change) after sterilization by the appropriate method;
- Manufactured with the use of transparent film materials, ensuring easy identification of sterilized devices;

INSTRUCTION FOR USE

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices

- Moisture-proof.
- Sterilized pouches and rolls must be handled as less as possible, kept possibly in closed cabinets in a clean environment, at a temperature between 8° and 35°C and a relative humidity of between 30 and 70%. If maintained as described, the devices will keep the contents sterile for a maximum 6 months and no later than the indicated expired date.

6. Sizes

Table 1.

Version	Width, mm	Length, m	Fold, mm
Rolls EUROSTERIL without folds	50 ± 5	200 ± 5 %	-
	75 ± 5	200 ± 5 %	-
	100 ± 5	200 ± 5 %	-
	150 ± 5	200 ± 5 %	-
	200 ± 5	200 ± 5 %	-
	250 ± 5	200 ± 5 %	-
	300 ± 5	200 ± 5 %	-
	400 ± 5	200 ± 5 %	-
Rolls EUROSTERIL gusseted	200 ± 5	100 ± 5 %	50 ± 6
	250 ± 5	100 ± 5 %	50 ± 6
	250 ± 5	100 ± 5 %	65 ± 6
Pouches EUROSTERIL	50 ± 5	250 ± 5 mm	-
	75 ± 5	250 ± 5 mm	-
	100 ± 5	250 ± 5 mm	-
Pouches EUROSTERIL (self-adhesive)	90 ± 5	250 ± 5 mm	-
	140 ± 5	260 ± 5 mm	-
	190 ± 5	330 ± 5 mm	-

INSTRUCTION FOR USE

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices

7. Drawings of the devices

Fig. 1 Pouches EUROSTERIL (self-adhesive)




Fig. 2 Rolls EUROSTERIL without folds

INSTRUCTION FOR USE

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices


Fig. 3. Rolls EUROSTERIL gusseted.

INSTRUCTION FOR USE

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices

Fig. 4. Pouches EUROSTERIL

applied on the protruding paper part, fold the protruding part along the fold line indicated on the pouch, press with the adhesive layer to the package surface, pressing in the direction from the center to the edges; check adhesion zone for air bubbles. Rolls and pouches (except self-adhesive pouches) are sealed by a thermo sealing device.

- Packaged devices are grouped into baskets, cassettes, trays and plates. Packages are placed according to the principle film-to-film, paper-to-paper. The package shall be arranged with its transparent part facing the tray or plate surface.
- Sterilization of packed devices is performed in the sterilizers approved for use according to the procedure established in the Russian Federation, in accordance with the modes regulated by the existing documents, as well as according to the operation manual for a certain-type sterilizer.
- When loading a sterilizer, it is necessary to observe correctness of packaged devices placement in the sterilization chamber: do not allow packages contact with the walls, door (cover), follow the loading norm; do not put them one on another, and do not put baskets, cassettes, trays and plates with packaged devices on the upper shelf.
- Sterilizing chamber shall not be loaded by more than 2/3 of the volume, so that steam and gas could freely circulate between and inside the packages.
- In the absence of any violations, sterilization packages are opened aseptically, from the non-applied part of the device; then devices are placed on the "sterile table" or immediately used according to their intended purpose.

Sterility of devices is ensured by:

1. Sterilization process:

- Steam sterilization: at 121 °C during 20 minutes or at 134 °C during 5 minutes.
- Ethylene oxide sterilization: at 40 °C during 24 h.

2. Rolls and pouches sealing process:

- Tension or displacement of the device shall be avoided during sealing.
- Sealing area shall be clean and dry.
- In order to avoid folds and corrugations, it is required to strain and hold the edge to be sealed until it is completely placed into the sealing part of the machine.
- Avoid shutdowns of a sealing machine during the sealing process.

INSTRUCTION FOR USE

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices


	Version	Quantity of devices in a consumer package	Quantity of devices/consumer packages in transport package
Rolls EUROSTERIL without folds	50 mm × 200 m	-	8 rolls
	75 mm × 200 m	-	5 rolls
	100 mm × 200 m	-	4 rolls
	150 mm × 200 m	-	2 rolls
	200 mm × 200 m	-	2 rolls
	250 mm × 200 m	-	1 roll
	300 mm × 200 m	-	1 roll
	400mm × 200 m		1 roll
Rolls EUROSTERIL gusseted	200 mm × 50 mm (fold) × 100 m	-	2 rolls
	250 mm × 50 mm (fold) × 100 m		1 roll
	250 mm × 65 mm (fold) × 100 m		1 roll
Pouches EUROSTERIL	50mm × 250 mm	500 pcs	1 box
	75mm × 250 mm	500 pcs	1 box
	100mm × 250 mm	500 pcs	1 box
Pouches EUROSTERIL (self-adhesive)	90 mm × 250 mm	200 pcs.	9 boxes
	140mm × 260 mm	200 pcs.	6 boxes
	190 mm × 330 mm	200 pcs.	5 boxes

12. Description of symbols used in labeling of the package

Symbol	Description	Symbol	Description
	Date of production		Manufacturer
	Protect from moisture		Temperature range
	Lot number		Conformity to European standards

INSTRUCTION FOR USE

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices

Symbol	Description	Symbol	Description
	Read instructions for use		Conformity mark upon declaring conformity in the system GOST R
	Keep away from sunlight		For single use only
	Use before		Humidity range
	"Fragile" handle carefully		"Do not use if package is damaged"

13. Requirements for maintenance and repair of the medical devices.

Devices are easy to apply, intended for single use and do not require any special training or qualification of the user. They do not require any maintenance or repair.

14. Storage and transportation conditions.

Storage. Rolls and pouches EUROSTERIL shall be stored in consumer and transport packages away from direct sunlight, in a dry and clean place, at the temperature from +5 °C to +30 °C and relative humidity from 5% to 50 %.

In case of partial use of packages from the transport box, further storage of remaining packages (within the established shelf life period) shall be effected in the same transportation box, which shall be closed tightly.

Transportation. Transportation of devices packed into transport package shall be effected by all types of roofed vehicles with observance of all precautions specified on the transport package and in accordance with the cargo transportation rules effective for a certain type of vehicle. Transportation conditions are identical to the storage conditions: at the temperature from +5 °C to +30 °C and relative humidity from 5% to 50 %.

15. Sterilization, disinfection

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices are non-sterile and intended for single use only.

INSTRUCTION FOR USE

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices

16. The procedure of disposal and destruction of medical devices.

Medical devices shall be disposed in accordance with the medical procedures adopted in health care institutions and with local, state or federal regulations. Choice of methods of safe disposal and destruction of devices depends on the capacity and profile of the healthcare institution, existence of wastes decontamination/detoxification installations, method of wastes decontamination/detoxification adopted on the administrative territory (incineration, removal to landfills, disposal).

Medical devices Rolls and pouches EUROSTERIL for sterilization and storage of medical devices, consumer package, defective devices and expired devices are classified as epidemiologically safe wastes similar in composition to household wastes.

17. Shelf life

Medical devices shelf life is calculated from the date of production is 5 years provided that the storage and transportation conditions are observed.

Storage period of sterilized material is not more than 6 months provided that all the requirements are observed.

Storage period of sterilized material shall not exceed the shelf life of the medical device Rolls and pouches EUROSTERIL for sterilization and storage of medical devices.

INSTRUCTION FOR USE

Rolls and pouches EUROSTERIL for sterilization and storage of medical devices

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