JD 2407

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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

01134441

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of LAKONG MEDICAL
DEVICES CO., LTD on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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
Micro applicator dental disposable EuroType

Outlet options:

- | | |
|---|----------------------|
| 1. Micro applicator dental disposable EuroType LK4000R | - 100 pieces / pack; |
| 2. Micro applicator dental disposable EuroType LK4000F | - 100 pieces / pack; |
| 3. Micro applicator dental disposable EuroType LK4000S | - 100 pieces / pack; |
| 4. Micro applicator dental disposable EuroType LK4000XS | - 100 pieces / pack; |
| Instruction for use | - 1 piece. |

2. Information about the manufacturer of a medical device.

Manufacturer: Lakong Medical Devices, Co. Ltd

No 13, Jingci Road, Shilong, Dongguan, Guangdong, 523326 China

Phone +86-769-81807881, E-mail: sales@lakongdental.com web: www.lakonggroup.com

The authorized representative in the Russian Federation:

«Eurotype Rus» LLC

119334, Moscow, inner territorial urban municipal district Donskoy, 5th Donskoy driveway, 21B, building 10, 5th floor, room 1, apartment 35

Ph.: + 7 (495) 478-82-13, E-mail import@nika-dent.ru

3. Purpose of a medical device.

Dental disposable micro applicator EuroType (hereinafter referred to as micro applicator or device) is intended for applying dental materials. For single use only.

The device is a portable hand-held dental instrument. Consists of a handle with a tip / brush, on the surface / bristles of which the applied material is held; the device does not look like a syringe or plunger.

The device is for single use.

The device is non-sterile, not sterilizable.

3.1. Indications to use: application of a dental filling material (for example, amalgam, composite) to the tooth / teeth of a patient in dentistry, including for applying:

- hemostatic solutions;
- solutions for the indication of dental plaque;
- sealants;
- Dycal accessories;
- etching agents;
- binding agents;
- pads for the oral cavity;
- fluoride varnishing

3.2. Side effect: no.

3.3. Contraindications: individual intolerance of materials.

3.4. Potential users: the qualified medical personnel in health care organizations.

3.5. Route of administration: use according to the instruction for use.

3.6. Application conditions: in the organizations of health care by qualified specialists.

3.7. Precautionary measures:

ATTENTION! The device is intended for single use only!

ATTENTION! Reuse may lead to cross-infection.

ATTENTION! Medical personnel must wear gloves.

4. General parameters and properties of a medical device.

4.1. Description of the device


An applicator is a single-use device.

The micro applicator consists of a plastic rod and a nylon brush. The surface of the plastic rod is clean, even, colorfast, free of burrs and stains, non-deformable. The nylon brush is even, clean and not plucked, arranged in a sphere shape.

The neck of the micro applicator bends at different angles and does not break.

"Micro applicator dental disposable EuroType" in outlet options:

Drawing 1 – Micro applicator dental disposable EuroType LK4000R, general view

Drawing 2 – Micro applicator dental disposable EuroType LK4000F, general view

Drawing 3 – Micro applicator dental disposable EuroType LK4000SF, general view


Drawing 4 – Micro applicator dental disposable EuroType LK4000XS, general view

4.2. Technical properties

5.	Properties	Outlet option			
		LK4000R	LK4000F	LK4000S	LK4000XS
1	Color	Green	Yellow	White	Black
2	Length, mm	103			
3	Diameter, mm	3			
4	Diameter of a head (working part), mm	1		0.5	
5	Head length (working part), mm	3.5	2.5	2.0	4.0
6	Weight,	0.60			
7	Distance to a notch, a, mm	17			
8	Tensile strength of a plastic core, not less, than, MPa	7.0			
9	Overall dimensions of a tube, diameter x height, mm	4.5 x 110			
10	Tube gross weight, no more then	100			
11	Number of micro applicators in a tube, piece	100			
12	Overall dimensions of a box, length x width x height, mm	100 x 100 x 115			
13	Box weight gross, no more than	400			
14	Quantity of tubes in a box, piece	4			

Permissible tolerance of parameters $\pm 10\%$.

4.3. Requirements to materials

The materials contacting to a human body:

1. Polypropylene of the "Kayshun" brand, manufactured by "Gaishun Plastic" (China), dyed with a green colorant of the "Kayshun" brand, manufactured by "Gaishun Plastic" (China) - a plastic core of the disposable EuroType LK4000R micro applicator;
2. Polypropylene of the "Kayshun" brand, production of "Gaishun Plastic" (China) painted by yellow colorant of the "Kayshun" brand, production of "Gaishun Plastic" (China) - a plastic core of the micro applicator of disposable EuroType LK4000F;
3. Polypropylene of the "Kayshun" brand, production of "Gaishun Plastic" (China) - a plastic core of the micro applicator of disposable EuroType LK4000S;
4. Polypropylene of the "Kayshun" brand, production of "Gaishun Plastic" (China) painted by black colorant of the "Kayshun" brand, production of "Gaishun Plastic" (China) - a plastic core of the micro applicator of disposable EuroType LK4000XS;
5. Nylon (polyamide) of the "Kayshun" brand, production of "Gaishun Plastic" (China) - a nylon brush of the micro applicator of disposable EuroType LK4000R; nylon brush of the micro applicator of disposable EuroType LK4000F; nylon brush of the micro applicator of disposable EuroType LK4000S; nylon brush of the micro applicator of disposable EuroType LK4000XS.

Medical personnel must wear gloves.

Type of contact with the patient's body: Short-term contact with damaged tissues of the tooth (dentin) of the patient.

In manufacturing process of a device are not used: materials of animal and (or) human origin, medical devices;

5. Route of administration.

- 5.1. Check the integrity of the packaging (tube) and the expiration date of the device.
 - 5.2. Specify the type of the device.
 - 5.3. Put on gloves after pretreating hands with antiseptic solutions. Gently, open the package, remove the applicator.
 - 5.4. Use directly for the preparation of solutions (eg etching agents, fluoride varnishing agents, etc.).
 - 5.5. Use directly for precise application in hard-to-reach areas.
 - 5.6. Execute necessary manipulations.
- ATTENTION!** To exclude reuse, the doctor must break the micro applicator along a notch after use!
- 5.7. Dispose of the applicator in the prescribed manner..

6. Safety requirements.

The device is intended for single use only!

Reuse may lead to cross-infection.

ATTENTION! To exclude reuse, the doctor must break the micro applicator along a notch after use! Do not use the devices if the packaging has been opened or damaged.

The symbols used on a device and its packing

Symbol	Description
	Manufacturer Specifies the manufacturer of a medical device
	Date of production Date when was the device is made
	Use until Specifies date after which expiration the device should not be applied or be used
	Code of a consignment of goods Code of a series of a device (or consignments of goods)
	Ban on repeated use The medical device is intended for single use or for use on one patient during one procedure
	Temperature limit Specifies borders of temperature range within which the medical device can be applied without damage of its safety

Symbol	Description
	Not to use at packing damage Specifies, that in case of damage or opening of packing the medical device cannot be used
	"To throw out in a garbage can" The sign means that this packing should be thrown out in a casket

7. Cleaning and disinfection.

This medical device is delivered unsterile, for single use. Additional cleaning and disinfection before use is not required.

8. Sterilization.

This medical device is delivered unsterile.

The device is not subject to sterilization.

9. Conditions of use, transportation and storage.

This medical device can be transported by all types of transport in top-covered vehicles in accordance with the applicable rules for the carriage of goods for each type of transport.

During operation, the products are resistant to the influence of biological fluid.

During use, transportation and storage, this medical device is resistant to temperatures ranging from -5°C to $+55^{\circ}\text{C}$.

10. Life span (shelf life) and an indication of the prohibition of the use of a medical device after the end of life span (shelf life).

The shelf life of a medical device is 5 years from the date of manufacture.

Not to use a device after the expiration date.

11. Procedure of utilization of a medical device.

11.1. Packing is utilized to places of collecting household garbage.

Note – For the Russian Federation, the packaging of a new device is disposed of in accordance with SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, residential premises, operation of industrial, public premises, organization and conduct of sanitary anti-epidemic (preventive) measures.

11.2. Used devices are decontaminated and disposed of, in accordance with the requirements of local organizations and disposal regulations in the area. For more information on the dates and locations, where disposal takes place, please, contact your local authority or waste disposal company.

Note – For the Russian Federation, the disposal of devices used in the infectious ward of a medical facility is carried out in the manner prescribed by SanPiN SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, residential premises, operation of industrial, public premises, organization and conduct of sanitary anti-epidemic (preventive) measures. for class B waste (epidemiologically hazardous waste).

11.3. Recycling and disposal of unused devices after the expiration date is carried out in accordance with the requirements for the disposal of epidemiologically safe waste.

Note – For the Russian Federation, in order to prevent environmental pollution, all waste generated during the disposal of devices is subject to mandatory collection with subsequent disposal, in accordance with the established procedure and in accordance with the current requirements and norms of local regulatory documents, including in accordance with SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, for water bodies, drinking water and drinking

water supply, atmospheric air, soils, residential premises, operation of industrial, public premises, organization and conduct of sanitary anti-epidemic (preventive) measures, for Class A waste.

11.4. Since the regulations for disposal may differ in different countries, if necessary, contact the manufacturer's authorized representative in the Russian Federation Eurotype Rus Limited Liability Company (Eurotype Rus LLC), 119334, Moscow, inner territorial urban municipal district Donskoy, 5th Donskoy driveway, 21B, building 10, 5th floor, room 1, apartment 35 E – mail: import@nika-dent.ru (tel. +7 (495) 478-82-13).

12. Maintenance and repair.

It is not provided, a device of single use.

13. Warranty obligations.

Lakong Medical Devices, Co. Ltd (Лакун Медикал Девайсес Ко, Лтд) assumes a warranty for the medical device to be free from defects and conform to the manufacturer's declared characteristics during the shelf life, subject to the conditions specified in the instructions for use.

The manufacturer is not responsible for damage or accidents resulting from misuse, user error or external emergencies. For device quality questions, please contact an authorized representative of the manufacturer.

13. Claim.

In the event of a defect, the developer and manufacturer Lakong Medical Devices, Co. Ltd located at No 13, Jingci Road, Shilong, Dongguan, Guangdong, 523326, China guarantees the quality of the device. In case of found defect, please, contact your authorized representative in the Russian Federation.

Authorized representative of the developer and manufacturer in the Russian Federation - Limited Liability Company "Eurotype Rus" (LLC "Eurotype Rus"), 119334, Moscow, inner territorial urban municipal district Donskoy, 5th Donskoy driveway, 21B, building 10, 5th floor, room 1, apartment 35 E – mail: import@nika-dent.ru (tel. +7 (495) 478-82-13) - accepts complaints.

Violations of the device, received during operation and careless handling are not considered as defects.

14. List of the international regulatory documents.

- ISO 9001:2015 "Quality management system. Requirements";
- ISO 13485:2016 "Medical devices. Quality management systems. Regulation requirements";
- ISO 15223-1:2016 "Medical devices. Symbols used on labels of medical devices when labeling and in the information provided. Part 1. General requirements";
- ISO 14155:2011 "Clinical research on medical devices for humans. Appropriate Clinical Practice";
- ISO 14971:2019 "Medical devices. Application of risk management to medical devices";
- ISO 10993-1:2018 "Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing within the Risk Management Process";
- ISO 10993-10:2010 "Medical devices. Assessment of the biological effect of medical devices. Part 10. Investigations of irritant and sensitizing effects".

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