1221

MAQUET

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Manufacturer:
Maquet Cardiovascular LLC
45 Barbour Pond Drive Wayne, NJ 07470, USA

Patel

June 15, 2022

/ Signature

Rachana Patel
Regulatory Affairs Specialist II
Maquet Cardiovascular LLC

stamp

Susan M Eichler-Huston

SUSAN M. EICHLER-HUSTON
NOTARY PUBLIC
STATE OF NEW JERSEY
MY COMMISSION EXPIRES APRIL 18, 2027

2022
New Jersey, USA

Ultima Activator II Drive Mechanism

ENGLISH


Figure 1.

Ultima Activator II Drive Mechanism and AccessRail Platforms (Sold Separately)

1. Ultima Activator II Drive Mechanism
2. Platform Release Latches
3. AccessRail Platforms
4. Suture Holders
5. Drive Handle

INSTRUCTIONS FOR USE

Carefully read all Instructions prior to use. Observe all warnings and precautions noted throughout these instructions.

CAUTION: Federal (U.S.A.) law restricts this device to sale, distribution and use by or on the order of a physician.

DESCRIPTION

The Ultima Activator II Drive Mechanism (P/N UA-5001) is a reusable device designed for use with the AccessRail Platform (Sold Separately). The AccessRail Platform along with the Ultima Activator II Drive Mechanism, is designed to create surgical access to, and direct visualization of, the thoracic cavity through a sternotomy incision.

INDICATIONS

The Ultima Activator II Drive Mechanism in combination with the Ultima AccessRail Platform is used to spread the sternum, providing access and direct visualization to the thoracic cavity.

WARNINGS AND PRECAUTIONS

- Physicians should be properly trained to perform cardiac surgical procedures with this device prior to use.

- As with all surgical retractors, care should be taken to use only as much retraction as necessary to provide adequate access and visualization.
- Take care not to adjust or move the AccessRail Platform while pericardial sutures are engaged in the platform.
- The Ultima Activator II Drive Mechanisms are supplied non-sterile and must be cleaned and sterilized prior to each use. Never use saline solution for cleaning surgical instruments.
- Steam sterilization is recommended to sterilize the Ultima Activator II Drive Mechanism. Other methods of sterilization are not recommended.

PREPARATION FOR USE

Thoroughly clean and sterilize the Drive Mechanism prior to each use. Follow recommended cleaning and sterilization instructions as described in this section.

The Ultima Activator II Drive Mechanism is manufactured from stainless steel. Proper care and handling of the instrument is required to ensure that the product maintains its appearance, corrosion resistance, and proper function over time.

Prior to use, sterilize the Ultima Activator II Drive Mechanism by cleaning and autoclaving the mechanism. MAQUET Cardiovascular provides the following recommended cleaning, sterilization and handling parameters for the Ultima Activator II Drive Mechanisms.

REPROCESSING INSTRUCTIONS:

General Information:

1. Maquet Getinge Group reusable medical devices are supplied ~~non-~~ **sterile**.
2. These devices must be thoroughly cleaned and sterilized prior to each use.
3. Reusable medical devices should be handled carefully to avoid damage.
4. Reusable medical devices that are used for sterile procedures or in sterile parts of the body should be sterilized and stored in a protective container/basket that is perforated.
5. Do not cool a reusable medical device after sterilization by rapidly exposing it to cool air or liquid. Sudden temperature changes may cause damage to medical devices.
6. Do not use immediate use steam sterilization (IUSS) for sterilization of the Maquet Getinge Group reusable medical devices. An IUSS cycle can introduce sudden temperature changes which can cause damage to some materials used to manufacture some medical devices.
7. For disassembly instructions see the IFU for the specific device.

Special Accessories:

1. There are no special accessories required for cleaning of any of the Maquet Getinge Group reusable medical devices.
2. A soft-bristled, non-metallic brush can be used for removing gross debris from the reusable medical devices.
3. Any cloth used for cleaning should be lint-free.

Cleaning Agents (ENZYMATIC DETERGENTS):

Enzymatic detergents should meet the following characteristics. The customer should evaluate the detergent(s) used at the facility to ensure that these criteria are met.

1. The enzymatic formulation contains ingredients for the removal of protein and lipid based soils. Additional ingredients that will remove other types of soil are acceptable as long as the detergent will remove proteins and lipids.
2. The detergent contains surfactants.
3. The solution has a neutral pH, between 6.0 and 8.5, and is compatible with device substrates including: stainless steels, thermoplastic polymers, elastomer, glass, soft metals and medical grade silicone.
4. Enzymatic detergents with a higher pH are acceptable for cleaning in an automated washer.
 - a. In an automated washer the diluted pH can be 6.0 to 12.
5. Examples of detergents that meet the acceptable characteristics include:

- a. GETINGE POWERCON TRIPLE ENZYME CONCENTRATE DETERGENT (used in the cleaning validations that were performed by Maquet Getinge Group)
- b. METREX METRIZYME
- c. ECOLAB OPTIPRO MULTI-ENZYMATIC LOW FOAMING DETERGENT
- d. STERIS POLYSTICA ULTRA CONCENTRATE ENZYMATIC
- e. BORER MEDICAL DECONEX TWIN ZYME

Cleaning Agents (ALKALINE DETERGENTS):

Alkaline detergents should meet the following requirements. The customer should evaluate the detergent(s) used at the facility to ensure that these criteria are met.

1. The solution has a pH of 8.5 – 12.5 diluted, 13.0 – 14.0 concentrate.
2. Dissolves and disperses clinical soils.
3. Effective for removal of protein, blood, and fat.
4. Compatible with device substrates including: stainless steels, thermoplastic polymers, elastomers, glass, soft metals (e.g. aluminum) and medical grade silicone.
5. Examples of detergents that meet the acceptable characteristics include:
 - a. GETINGE POWERCON ALUMINUM SAFE CONCENTRATE DETERGENT (used in the cleaning validations performed by Maquet Getinge Group)
 - b. STERIS POLYSTICA ULTRA CONCENTRATE ALKALINE DETERGENT
 - c. NEODISHER MEDICLEAN FORTE
 - d. BORER MEDICAL DECONEX TWIN PH10


Water Quality:

1. Water quality can affect the efficacy of cleaning agents. Depending on water quality more or less detergent will be needed.
2. Water hardness can be a concern because mineral deposits from the water can be left on device materials and this can impact decontamination of the device.
3. Purified water is recommended for the final rinse for most cleaning processes.
4. Inadequate water quality can contribute to residue left on a medical device which could cause problems with functioning of the device, e.g., film on the lens optics of an endoscope.
5. AAMI TIR34 *Water for the reprocessing of medical devices* and the detergent manufacturer should be consulted for more information on determining water quality and detergent dosage.


POINT-OF-USE PROCESSING:

Disassembly: UA-5001 Ultima Activator II Drive Mechanism

1. The drive mechanism disassembles into three (3) pieces: blade, retracting mechanism and handle.
2. Grasp the drive handle and turn counter-clockwise until the moveable retracting mechanism slides off the end of the serrated blade.
3. Remove drive handle from moveable retracting mechanism.

UA-5001 Ultima Activator II Drive Mechanism – Assembled

UA-5001 Ultima Activator II Drive Mechanism – Disassembled

Remove Gross Soil/Pre-Cleaning:

1. Rinse device by immersion in a basin of water or rinsing under running water.
2. Use a syringe to flush hard to reach areas including grooves, crevices, recesses and mated or near-mated surfaces.
3. If necessary, remove gross debris by brushing the outside of the device with a soft-bristled non-metallic brush or wiping with a lint-free cloth.
4. Spray the device with Getinge Clean Renuzyme Foam Spray (or equivalent enzymatic foam spray), ensuring all surfaces are thoroughly covered with foam for a minimum of three (3) minutes.
5. Prior to cleaning, rinse the device with tap water for a minimum of two (2) minutes, using a syringe and/or soft bristled non-metallic brush to assist in the removal of gross soil and debris, including the Renuzyme foam (or equivalent). The rinse may be performed before, or after transport to processing area.

Maintain Moisture:

1. Keep the device(s) moist until the time that cleaning begins.
2. Place a damp towel on top of the device(s); or
3. Place device(s) in a humidification bag or container.

Transport to Processing Department/Area:

1. Medical devices should be confined and contained during transport.
2. Transport device(s) to the processing area as soon as possible after performing the point-of-use processing.
3. Soiled devices should be transported from the point-of-use to the processing area in a puncture-proof container that completely covers the contaminated device(s).
4. Containers used for transportation of contaminated medical devices should be appropriately identified as biohazard.
See OSHA 29 CFR 1910.1030.

Choose one of the following cleaning methods for reprocessing of the device:

1. Manual Cleaning
2. Ultrasonic Cleaning or
3. Automated Washer Cleaning.

MANUAL CLEANING:**Soak:**

1. Prepare a detergent bath using a neutral pH enzymatic detergent. Follow the detergent manufacturer's instructions for detergent concentration and solution temperature.
2. Soak the device in the detergent solution for 10 minutes.
3. Use a syringe to flush all hard to reach areas of the device, including grooves, crevices, and recesses and mated or near-mated surfaces.

4. While the device(s) is soaking brush all surfaces of the device(s) with a soft-bristled non-metallic brush for 4 minutes.
5. Ensure that all outer surfaces of the device are brushed.
6. Clean device under water to prevent aerosolization of contaminants.

Rinse:

1. Rinse the device(s) under running warm tap water for 2 minutes.
2. Use a syringe to flush all hard to reach areas including grooves, crevices, recesses and mated or near-mated surfaces.
3. Rinse the device(s) with purified water for 30 seconds.
4. Dry devices with a clean, dry, lint-free cloth, or filtered compressed air.

ULTRASONIC CLEANING:

If the device was not disassembled at the point-of-use, the device will need to be disassembled and the pre-cleaning instructions will need to be followed here. Go to "POINT-OF-USE-PROCESSING", "Disassembly" and "Remove Gross Soil/Pre-Cleaning" sections and follow the instructions there.

Sonication:

1. Prepare a detergent bath using a neutral pH enzymatic detergent. Follow the detergent manufacturer's instructions for detergent concentration and solution temperature.
2. Flush all cannulas, lumens and crevices with detergent solution using a syringe, pipette or water jet.
3. Sonicate the device in the ultrasonic bath for 10 minutes.
4. Rinse device(s) under running tap water for 2 minutes.
5. Use a syringe, pipette or water jet to flush all lumens, grooves and crevices.
6. Rinse device(s) with purified water for 30 seconds.
7. Dry devices with a clean, dry, lint-free cloth, or filtered compressed air.

AUTOMATED WASHER CLEANING:**Soak:**

1. Prepare a detergent bath using a neutral pH enzymatic detergent. Follow the detergent manufacturer's instructions for detergent concentration and solution temperature.
2. Soak device(s) for 10 minutes.
3. While device(s) are soaking flush all surfaces of the device with detergent solution using a syringe, water jet or water pistol.
4. Rinse under running tap water for 2 minutes.

Automated Washer:

1. Place the device(s) into a washer basket to contain the device(s) during the automated washer process.
2. Process device(s) through the automated washer with the following parameters.

Cycle/Phase	Minimum Time	Temperature Set Point	Type of Detergent and Dose
Pre-wash (tap water)	5 minutes	Cold tap water	N/A
Enzyme Wash	5 minutes	Hot tap water	• Neutral pH enzymatic detergent • Consult detergent manufacturer for dosage to use.
Enzyme Rinse	30 seconds	Hot tap water	N/A
Wash 1	5 minutes	55°C (minimum)	• Either a mild alkaline detergent or a neutral pH enzymatic detergent can be used. • Mild alkaline detergent • Neutral pH enzymatic • Consult detergent manufacturer for exact dosage to use for each detergent.
Rinse 1 (tap water)	2 minutes	43°C (minimum)	N/A
PURW Rinse	30 seconds	43°C (minimum)	N/A
Dry	7 minutes	85°C (maximum)	N/A

3. Exact dosage will depend on water quality. The detergent manufacturer should be consulted to assist in analyzing water quality and setting of detergent dosage.

Thermal Disinfection/Automated Washer (for OUS users only):

1. If thermal disinfection is desired a thermal disinfection step can be programmed after Rinse 1 in the automated washer program. The purified water rinse should follow the thermal disinfection rinse.

2. Medical devices that will be sterilized should be thermally disinfected with an $A_0 = 600$ to protect employees who handle the device(s) between cleaning and sterilization. The device must be sterilized prior to use.
3. Depending on the level of disinfection desired use one the following parameters for the thermal rinse cycle/phase.

Cycle/Phase	Motor Speed	Recirculation Time	Water Temperature	Desired A ₂ Value	Type of Detergent and Dose
1) Thermal Rinse	High or Low	5 minutes	90°C	3000	N/A
2) Thermal Rinse	High or Low	10 minutes	80°C	600	N/A
3) Thermal Rinse	High or Low	1 minute	90°C	600	N/A
4) Thermal Rinse	High or Low	30 seconds	93°C	600	N/A
5) Thermal Rinse	High or Low	2 minutes 30 seconds (3 minutes)	93°C	3000	N/A

INSPECTION:

- After cleaning, devices should be visually inspected for flaws, damage, debris, detergent residue, and completeness. Automated Washer cleaning may cause cosmetic blemishes, which do not affect product performance or biological safety, when the retractor is processed according to the instructions for use. See ANSI/AAMI ST79 for more information on visual inspection of medical devices after cleaning.
- If debris or detergent residue is found on the device, the device should be returned to the decontamination area and re-cleaned until a visually clean endpoint has been reached.
- Remove the device from service if there is any sign of damage. A damaged device should be repaired or replaced.

STERILIZATION

STERILIZATION PROCESSES THAT ARE NOT RECOMMENDED

The following sterilization processes are not recommended.

- Gravity Displacement Steam Sterilization
- Immediate Use Steam Sterilization (IUSS)
- STERIS SYSTEM 1E
- Ethylene Oxide Sterilization
- STERRAD Sterilization
- STERIS VPRO Sterilization
- STERIZONE VP4 Sterilization

General Requirements:

- The UA-5001 Ultima Activator II Reusable Drive Mechanism has been validated for steam sterilization.
- When a device can be sterilized with steam it is preferable that steam sterilization is used.
- If the device can be disassembled then the device should be disassembled for sterilization.
- Adequate cleaning must be performed prior to any sterilization process.

STEAM STERILIZATION:

- Place device(s) in a container or basket appropriate for steam sterilization.
- If the device(s) is part of an instrument set, the set should not weigh more than 25 pounds (11 kilograms).
- Wrap the container/basket with the device(s) in standard sterilization wrap.
- Load the sterilizer chamber in accordance with the sterilizer manufacturer's instructions on loading of the sterilizer.
- Choose appropriate cycle in accordance with the information below and then close the sterilizer door.
- Initiate the cycle.
- Cycles
 - Prevacuum
 - Conditioning: 3-4 pulses
 - Exposure time: 4 minutes
 - Exposure temperature: 132°C
 - Dry time: 20 minutes minimum
 - Cool down time: 1 hour
 - Prevacuum (For OUS only)
 - Conditioning: 3-4 pulses
 - Exposure time: 3 minutes
 - Exposure temperature: 134°C
 - Dry time: 20 minutes minimum
 - Cool down time: 1 hour
 - Prevacuum (For OUS only)
 - Conditioning: 3-4 pulses
 - Exposure time: 3 minutes
 - Exposure temperature: 135°C

iv. Dry time: 20 minutes minimum

v. Cool down time: 1 hour

- Prevacuum - The following cycle utilizes the parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding Creutzfeldt-Jakob Disease or Transmissible Spongiform Encephalopathy¹. This Cycle is not for use in the United States and is not to be used for inactivation of prions.

- Conditioning: 3-4 pulses
- Exposure time: 18 minutes
- Exposure temperature: 134°C
- Dry time: 20 minutes minimum
- Cool down time: 1 hour

INSTRUCTIONS FOR USE

- Prior to use, close the Ultima Activator II Drive Mechanism by turning the drive handle clockwise.
- Assemble the AccessRail Platform onto the Ultima Activator II Drive Mechanism, matching the correct left and right platforms to the correct drive attachment side. Ensure that the drive and platforms are fully engaged.
- Ensure that the AccessRail Platform blades are properly seated on the sternum.
- Once the AccessRail Platform is in place, begin spreading the sternum to the desired opening by slowly turning the drive handle counter-clockwise.
- If using sutures, slide the sutures into the suture holder slots. Engage only one suture strand per slot to ensure proper hold.
- Release engaged suture from the platform by concurrently pulling back and up on the suture while pulling the suture through the holder slot.
- Perform Surgical procedure.
- To remove the Ultima Activator II System, fully close the system by turning the drive handle clockwise. Gently remove the system from incision.
- To remove AccessRail Platform from the Ultima Activator II Drive Mechanism, turn release latches on the platform in the direction of the arrows and pull the platform away from the Ultima Activator II Drive Mechanism.

HOW SUPPLIED

The Ultima Activator II Drive Mechanisms are supplied NON-STERILE and must be steam sterilized prior to use.


WARRANTY

MAQUET Cardiovascular LLC (MAQUET) warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond MAQUET's control directly affect the instrument and the results obtained from its use. MAQUET's obligation under this warranty is limited to the repair or replacement of this instrument and MAQUET shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. MAQUET neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. MAQUET assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

¹ WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies

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MAQUET is a Trademark of MAQUET GmbH & Co. KG
All rights reserved.