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Abbott Laboratories
100 Abbott Park Road
Abbott Park, Illinois, 60064, USA


I, Mary Caren Murawski, Associate Director, Abbott Diagnostics Division Regulatory Affairs, Abbott Laboratories, hereby certify that the attached documents that keep the number of pages indicated below are the Insert. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: 
Mary Caren Murawski
Associate Director, Regulatory Affairs
Abbott Diagnostics Division

Date: 16 November 2017

Abbott Laboratories
100 Abbott Park Road
Abbott Park, Illinois, 60064, USA


State of Illinois, County of Lake
United States of America

Signed before me on November 16, 2017 (Date)

By, Mary Caren Murawski, Associate Director, Regulatory Affairs, Abbott Diagnostics Division.

(Seal)

(Signature of Notary Public)

Multiconstituent Calibrator

FOR USE WITH
ARCHITECT

en
Multiconstituent Calibrator
1E65
307292/R09
S1ET50

Read Highlighted Changes: Revised August 2017.

This package insert contains information for the use of Multiconstituent Calibrator on the ARCHITECT cSystems.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

INTENDED USE

For use in the calibration of the Albumin, Calcium, Cholesterol, Creatinine, Glucose, Iron, Lactic Acid, Magnesium, Phosphorus, Total Protein, Triglyceride, Urea Nitrogen, and Uric Acid assays.

CONTENTS / MATERIALS PROVIDED

REF 1E65-05 Multiconstituent Calibrator

CAL 1+2 3 x 5 mL

CAL VALUE MEDIA 1

Multiconstituent Calibrator is prepared from a human-based matrix containing the following analytes: albumin, calcium, cholesterol, creatinine, glucose, iron, lactic acid, magnesium, phosphorus, total protein (made from human sourced albumin material), triglyceride, urea nitrogen (urea), and uric acid. Sodium azide is present as a preservative. Lot-specific calibrator values are listed in the Multiconstituent Calibrator Value Sheet, packaged with the calibrator. Users can import the values electronically (refer to the INSTRUCTIONS FOR USE section of this package insert). The last two digits of the lot numbers can vary.

When converting the units for an assay, the following items must be edited, if applicable, using the appropriate conversion factor.

- Calibrator Concentration
- Qualitative/interpretation Range
- Linear Range
- Retest rules ranges
- Constituent assay ranges for calculated assays
- Reference Range
- Panic/Extrema Range
- Expected QC Values

STANDARDIZATION

Multiconstituent Calibrator is prepared and standardized as described in the table below.

Analyte	REF	Reference Material	Reference Method
Albumin	7D53, 7D54	ERM-DA470	Visible spectrometry
Calcium	3L79	NIST SRM 956	ID-ICP-MS
Cholesterol	7D62	SPRL Total Cholesterol Verification Set	Abell-Kendall procedure (CDC reference method)
Creatinine (Serum/Plasma)	3L81	NIST SRM 987	ID-LC/MS*
Creatinine (Urine)	3L81	NIST SRM 914	Gravimetric
Glucose	3L82	NIST SRM 985	ID-GC/MS
Iron	6K95	NIST SRM 3126	Gravimetric
Lactic Acid	9D89, 9P18	Reagent Grade Lactic Acid	Titration
Magnesium**	3P88	NIST SRM 956	ID-ICP-MS
Phosphorus	7D71	NIST SRM 2186-I	Gravimetric
Total Protein	7D73	NIST SRM 927	Biuret
Triglyceride	7D74	ACS Grade Glycerol	Gas Chromatography
Urea Nitrogen***	7D75	NIST SRM 912	Differential scanning calorimetry
Uric Acid	3P99	NIST SRM 913	Gravimetric

* Use IDMS-traceable MDRD equation with factor 175 to estimate GFR.

** Available Ex-US only.

*** 7D75 Urea Nitrogen assay reports concentrations of urea nitrogen (mg/dL) in conventional units and concentrations of urea (mmol/L) in SI units.

ACS—American Chemical Society

ERM—European Reference Material

GFR—Glomerular Filtration Rate

ID-GC/MS—Isotope Dilution-Gas Chromatography Mass Spectrometry

ID-ICP-MS—Isotope Dilution Inductively Coupled Plasma Mass Spectrometry

ID-LC/MS—Isotope Dilution-Liquid Chromatography Mass Spectrometry

MDRD—Modification of Diet in Renal Disease

NIST—National Institute of Standards and Technology

SRM—Standard Reference Material

PRECAUTIONS

Precautions for Users

- **IVD**
- For *In Vitro* Diagnostic Use.
- Do not use components beyond the expiration date.
- Do not mix components from different kit lot numbers.
- **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the CONTENTS/MATERIALS PROVIDED section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that this calibrator and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹ Biosafety Level 2² or other appropriate biosafety practices^{3,4} should be used for materials that contain or are suspected of containing infectious agents.
- The human serum albumin used in the calibrators has been tested and found to be nonreactive for HBsAg, anti-HIV-1/HIV-2, anti-HCV, and syphilis.

The following warning and precaution apply to **CAL 1+2**:

Contains sodium azide.

EUH032 Contact with acids liberates very toxic gas.

P501 Dispose of contents/container in accordance with local regulations.

NOTE: Refer to Section 8 of the ARCHITECT System Operations Manual for proper handling and disposal of products containing sodium azide.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

STORAGE AND STABILITY

Store unopened Multiconstituent Calibrator upright at 2 to 8°C.

Unopened calibrator is stable until the expiration date when stored at 2 to 8°C.

Opened calibrator is stable for 7 days at 2 to 8°C or for 24 hours at 15 to 30°C if kept tightly capped.

PREPARATION OF CALIBRATOR

Multiconstituent Calibrator requires no preparation prior to use.

INSTRUCTIONS FOR USE

Configure Calibrator Values

Required User Access Level: System Administrator

Calibrator values may be configured manually, or configured electronically using the **CAL VALUE MEDIA** packaged with the calibrator or using e-files accessed and imported from www.abbottdiagnostics.com, or Abbott mail.

Ensure the assay numbers on your system match the assay numbers listed in the MCC value sheet. If the system assay number(s) differ from those listed in the value sheet, a status of "No Assay" is generated when importing the calibrator file. You must manually configure that assay.

MANUAL CONFIGURATION

1. Verify the lot number listed on each calibrator carton agrees with the lot number printed on the value sheet.
2. Enter the calibrator values provided in the value sheet into the Configure calibrator set screen.

INSTRUCTIONS FOR USE (Continued)

ELECTRONIC CONFIGURATION

The **CAL VALUE MEDIA** contains the calibrator e-files. Refer to the Import Calibrator Values section of this package insert for instructions.

To Download Calibrator Values from www.abbottddiagnostics.com

1. Access www.abbottddiagnostics.com.
2. Select **Support > Technical Library > Value Assignments > ARCHITECT cSystems > cSystems Calibrators**.
3. The calibrator .pdf (Printable) file and .xml (Downloadable) file options are available; select the .xml (Downloadable) file.
4. Right click and select **Save Target As** to save the calibrator file to a USB flash drive or a CD-ROM. Refer to the Import Calibrator Values section of this package insert.

To Access Calibrator Values in Abbott Mail

Availability of calibrator values from Abbott mail may vary by area and requires ARCHITECT System Software v8.00 or higher AND AbbottLink software v4.00 or higher.

NOTE: Calibrator value sheet(s) can be viewed only from the Abbott mail screen. The calibrator files sent via Abbott mail are stored and available on the hard drive when performing the **Import calibrator data (cSystem)** procedure.

1. Select the Abbott mail icon on the Snapshot screen. The Abbott mail screen displays.
2. Select **F5 - View** or **F7 - Review** to view the calibrator values .pdf (Printable) file. Refer to **Section 2** of the ARCHITECT System Operations Manual for additional information. Refer to the Import Calibrator Values section of this package insert.

Import Calibrator Values

1. From the menu bar Select **System > Configuration > QC-Cal** settings.
2. Select **Calibrator set** from the **QC-Cal** categories list, select **MCC** from the Calibrator sets list, then select **F8 - Configure**.
3. Select the **Lot number list** button and then select **Import**.
4. If using electronic media, insert the CD-ROM or USB flash drive and then select **OK**. Abbott Mail users select **OK** or **Cancel**. Available calibrator files from Abbott mail or electronic media display on the **Import lot file selection** window.
5. Select the calibrator file to be imported; select **Done**.
6. Review the **Status** column in the **Assign assays for calibrator set** window. For any assays with a status other than **OK** or **Assigned - System**, refer to **Import calibrator data (cSystem)** in **Section 2** of the ARCHITECT System Operations Manual.
7. Select **Done** in the **Assign assays for calibrator set** window.
8. Select **Done** to save the calibrator values and exit the screen.
9. Select **OK** if the data were imported from a USB flash drive then remove the drive from the USB port.

Calibration

Calibration is performed by running a water blank and the Multiconstituent Calibrator set. Water for the blank is provided by the instrument.

1. Verify that the correct calibrator values have been entered into the calibration file.
2. Allow calibrator to come to room temperature.
3. Mix bottle five times by gentle inversion.
4. Open bottle, place an appropriate amount of each calibrator in a separate sample cup, and place in the assigned positions.
5. Cap bottle tightly and return to refrigerated storage immediately after use.
6. Perform calibration as indicated in the ARCHITECT System Operations Manual.

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are visible signs of leakage, extreme turbidity, microbial growth, or if calibration does not meet the criteria in the reagent package insert or the ARCHITECT System Operations Manual.

LIMITATIONS OF THE PROCEDURE

Accurate and reproducible results are dependent upon properly functioning instruments and reagents, storage of product as directed, and good laboratory technique.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections: Approved Guideline - Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Key to Symbols

	Assay number
	Calibrators 1 and 2
	Calibrator 1
	Calibrator 2
	Calibrator value media
	Contains sodium azide. Contact with acids liberates very toxic gas.
	Conventional Units
	Authorized Representative in the European Community
	Identifies products to be used together
	Information needed for United States of America only
	In Vitro Diagnostic Medical Device
	Batch code/Lot number
	Manufactured for
	Product of USA
	Reagent
	Catalog number/List number
	SI Units
	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	Waste Electrical and Electronic Equipment

The ARCHITECT cSystem family of instruments consists of c4000, c8000, and c16000 instruments.

ARCHITECT, c4000, c8000, c16000, and cSystem are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott Laboratories
Diagnostics Division
Abbott Park, IL 60064
USA

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

MANUFACTURED FOR

Abbott Laboratories

Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottddiagnostics.com

August 2017

©2011, 2017 Abbott Laboratories

ADDENDUM

to the Instruction for Use "Multiconstituent calibrator" (for Russia)

SHIPPING

Recommended shipping conditions for device correspond to storage conditions provided in the Instructions for Use (2-8°C)

Device is stable until the expiration date indicated on the package.

Device is no longer suitable for use if shipping temperature limits have been violated (GOST R 51088-2013). Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. Manufacturer set the warranty shelf life of device which corresponds to period between the date of manufacture and expiration date that are indicated on the label and/or in accompanying documentation. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required

REPAIRS

Not required

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

GOST R 51088-2013, GOST R 51352-2013 GOST R ISO 23640-2015, GOST R 56894-2018 / GHTF / SG1 / N063:2011, GOST ISO 17511-2011, GOST ISO 18153-2011, GOST R ISO 18113-1-2015, GOST R ISO 18113-2-2015, ISO 9001:2000, ISO 13485:2003

POTENTIAL USER

Professional user

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

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